
Computed Tomography

Detectors & Crystals

Cutting-edge technologies for
Siemens Healthineers CT detectors

siemens-healthineers.com/ufc

Contents

Welcome to Detectors & Crystals	4
Origins of Computed Tomography	6
State-of-the-art detectors since 1996	8
Locations	12
The detector: Core component of every CT system	15
Operating principle: Energy-integrating detectors (EID)	16
CT scintillator material	18
Operating principle: Photon-counting Computed Tomography (PCCT)	21
From sensor to detector	24
Highly innovative front-end electronics	26
Detector resolution and z-coverage	28
Quality characteristics	30
Overview: Producing UFC and photon-counting detectors	34
OEM business for UFC	58
Sustainability and environmental protection	60

Welcome to Detectors & Crystals

Pioneering CT technology since 1996

At Detectors & Crystals, we develop and manufacture ultra-modern X-ray detectors for computed tomography scanners.

Since we began in Forchheim in 1996, our detectors have been at the heart of all Computed Tomography scanners from Siemens Healthineers. Our CT scanners are reliable and very precise high-tech systems in use in medical facilities all over the world.

At Detectors & Crystals, an interdisciplinary team of experts in research and development, production process development, and process management is constantly working on cutting-edge technology. Innovation is part of our DNA as well as our daily mission.

As the first company to offer photon-counting detectors that deliver high image quality with surprisingly small radiation doses, we've revolutionized the field of CT scanning. Our ultra-fast scintillation ceramics are also used in our successful fleet of energy-integrating CT scanners.

We engage our creativity and passion to develop these technologies on a daily basis. With the construction of our new factory in Forchheim, we're making an investment in the future of healthcare technology.

But we're not just active in Forchheim: Learn more about our other locations in this brochure!

The manufacturing depth at our Detectors & Crystals Center is exceptionally high. We're intensively involved in developing the new production processes required for our highly precise and innovative manufacturing techniques.

From the raw material of scintillator ceramics (a high-purity powder made from rare earth elements) to the front-end sensors and the finished CT detector, the entire manufacturing chain is represented here at our Forchheim location.

In our new factory, we'll produce photon-counting sensors made of cadmium telluride, structure and connect them with highly integrated circuits, and install them in our ultra-modern CT detectors.

After intensive testing, the detectors are delivered to the CT manufacturing line, which you can visit here at the Forchheim site.

We cordially invite you to learn more about our unique competence center for detectors on the following pages. It's the only one of its kind in Europe!

Uwe Rückl

VP SHS DI Computed Tomography
Head of Detectors & Crystals

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Origins of Computed Tomography

1895

On November 8, 1895, German physicist Wilhelm Conrad Röntgen made an important discovery in his lab at Würzburg: X-rays. For this breakthrough, he was awarded the first Nobel Prize in Physics in 1901.

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
1895

With X-rays, it became possible for the first time in history to see inside a person's body without having to cut it open. However, anatomical structures often overlapped, and soft tissue couldn't be differentiated.

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1972

In 1972, English electrical engineer Godfrey Hounsfield and radiologist Jamie Ambrose introduced the first CT image from an experimental prototype. Hounsfield was awarded the Nobel Prize in Medicine in 1979.

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1975


The first Computed Tomography scanner from Siemens was the SIRETOM, a head scanner that went into serial production in 1975. This marked the beginning of a wonderful success story. In 2025, we can proudly look back on a 50-year history packed with innovation and technology leadership.

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State-of-the-art detectors since 1996

1996 – Opening of the CT Detector Center

1997 – First UFC solid-state dedector in SOMATOM Plus4

2003 – Development of multislice technology

1997

First UFC solid-state detector with special material for converting X-rays premieres in the SOMATOM Plus4.

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1998

First multislice detector "Volume Zoom" with four slices reduces scan time and increases z-coverage.

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2003


Multislice technology becomes the foundation for detectors with 6, 16, and 64 slices and large z-coverage.

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2005

World's first Dual Source scanner with two X-ray tubes and two detectors supports cardiac scanning without beta blockers.

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
2008 – Matrix detector module

2011 – Sillian tile technology for Stellar detector

2013 – Sillian tile technology

2008

First detector platform to provide tileable detector elements: Matrix technology enables 2 and 4 cm detectors (64 and 128 slices) for Single and Dual Source scanners. It establishes Dual Source technology for everyday use in the SOMATOM Definition Flash.

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2011

Sillian tile technology combines the photo diode and AD converter, resulting in reduced electronic noise and X-ray dose. The Stellar detector is the first fully integrated detector at Siemens Healthineers with this technology.

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2011


Japanese company Acrorad joins the Siemens Healthineers family. Located in Okinawa, Acrorad specializes in the development and production of CdTe-based X-ray detectors. This lays the foundation for the later introduction of photon-counting CT detectors.

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2013

New developments in the evolution of Sillian tile technology increase the manufacturing depth and lead to more reliable production processes and products without fragile Flex PCBs.

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2014 – Stellar module

2014 – In-house production of PDFEE board

2017 – Introduction of the new UFC detector platform

2014

Thanks to its higher resolution, our SOMATOM Force detector facilitates the development of Stellar Infinity Technology (192 slices). It's the first Siemens Healthineers detector with a 3D anti-scatter grid.

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2014

Even better manufacturing depth: The new production process for flip-chip handling enables in-house production of PDFEE boards for Stellar detectors.

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2017

Introduction of the UFC detector platform, the foundation of the new SOMATOM go. platform, using standardized components for 1, 2, and 4 cm Single Source scanners.

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
2019

SOMATOM X.cite is the first product for the new SOMATOM X. platform. The combination of our modular detectors with the Vectron tube results in high image quality in Single Source systems.

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2021 – Photon-counting sensors

2021 – Photon-counting detector module

2023 – Start of construction of the new photon-counting Sensor Fab

2021

NAEOTOM Alpha: The world's first medical CT system with a photon-counting detector. The Quantum Technology used in this system brings clinical advantages at low radiation doses.

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2023

SOMATOM Pro.Pulse is designed to make Dual Source technology more accessible to a broader range of clinical institutions.

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2023

Construction starts on the new crystal growing factory at Forchheim, which will expand the photon-counting capacities in addition to those at Acrorad.

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
2024

Introduction of the NAEOTOM Alpha class based on photon-counting Quantum technology with 4- and 6 cm-wide Dual Source and 6 cm-wide Single Source detectors.

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Locations

Detector development and production

Production modules and detectors

Photon-counting crystal growth and semiconductor manufacturing starting in 2026

Detectors & Crystals Forchheim

Development and manufacturing

The Detectors & Crystals Center has been located in Forchheim since 1996. The detector is the heart of a computer tomograph, and they're developed right here in Forchheim.

The Forchheim location also manages the production of UFC, arrays, and modules and the assembly and testing of both UFC-based detectors and photon-counting detectors.

Beginning in 2026, the sensors required for photon-counting detectors will also be produced in Forchheim. The production process will comprise crystal growth and dicing as well as metallization and lithography processes.

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SSME Shanghai

Module and DMS production

Module manufacturing and the assembly and testing of detectors take place at the Shanghai location.

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Acrorad Okinawa


Developing and producing materials for photon-counting detectors

Acrorad has a complete semiconductor production line. Starting with raw cadmium and tellurium, they grow high-purity crystals that are made into wafers. Surface-treated and pixel structures are created using photolithography. Finally, the sensors are diced, inspected, and sent to Germany.

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The detector: Core component of every CT system

In a Computed Tomography (CT) scanner, an X-ray tube and an X-ray detector rotate around a patient lying on a table in the classic configuration (Single Source). At Siemens Healthineers, we also use a configuration with two X-ray tubes and two detectors (Dual Source) to increase the temporal resolution of the CT system. This makes it better suited for cardiac and other precision imaging.

The main task of the detector is to register the X-rays attenuated by passing through the patient in one projection. Several thousand projections with a resolution in the range of microseconds are typically recorded in each rotation around the patient. Because the detector has anywhere from several hundred thousand to over one million pixels, the CT scan generates large amounts of data that are pre-processed on the detector and transferred to the stationary part of the CT. They're then further processed into reconstructed sectional images ("tomograms") – 3D reconstructions of the patient – and analyzed using specialized post-processing tools, depending on the clinical question, and finally read by the radiologist.

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Operating principle: Energy-integrating detectors (EID)


Before 2021, all CT detectors were designed based on the energy-integrating principle. Today they still represent the majority of detectors produced at our Detectors & Crystals Center. In the future, photon-counting technology – which was first introduced to the market by Siemens Healthineers – will take over a larger and larger share of the worldwide CT market.

Energy-integrating detectors use scintillator materials that convert the X-ray signal attenuated by the patient in two steps. First, scintillation in the detector material generates (visible) light proportional to the energy of the X-rays. At Siemens Healthineers, we call our scintillator material UFC (ultra-fast ceramics). In the second

step, this light is converted into an electrical signal by photo diodes that register all light within a defined time window. This means that energy information is lost in the process as the signals are combined in each projection, and the resulting signal is an energy-weighted cumulative signal that's further processed by front-end ASICs and other electronics.

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Two-step conversion: X-rays → light → electric current

CT scintillator material

Thanks to its fast decay behavior and low afterglow, our UFC scintillator material is ideal for CT systems with high rotational speeds and short integration periods.

The decay behavior of a scintillator is determined by the decay constant, which is 3 μs for our UFC. This means that UFC is ideally suited for measurement intervals (integration periods) down to 100 μs .

In contrast, the afterglow describes the signal behavior of a scintillator material over numerous integration periods and is therefore a decisive parameter for image quality. Our UFC shows an afterglow of 0.0005 percent after 40 ms, making it the ideal detector material for most CT applications on the market

One of the definitive parameters for optimal dose utilization is the efficiency with which a scintillator

material absorbs X-ray radiation. Compared to other materials of the same thickness, our UFC scintillator material has greater X-ray absorption efficiency over a wide energy range. For example, a LuTAG detector needs to be 50 percent thicker in order to be capable of absorbing the same number of X-ray quanta as a UFC detector. Due to their greater thickness and the associated higher aspect ratio, the LuTAG detector's luminous efficiency is decreased.

What impact do these differences have on image quality? Even with short integration periods, our CT scanners don't require an artificial modification of the image quality via correction algorithms. Our fast and efficient scintillator material also enables a short scan time and consistently high image quality. Both of these factors mean a lower radiation dose for the patient.

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

X-ray absorption efficiency
for 1.4 mm detector material

Afterglow behavior
log I/I₀

Single-step direct conversion: X-rays \rightarrow electric current

Operating principle: Photon-counting Computed Tomography (PCCT)

Unlike energy-integrating detectors, the new photon-counting CT (PCCT) detectors directly convert the incoming X-ray quanta into an electric signal. This is realized by using a direct X-ray converter material. At Siemens Healthineers, we've found that CdTe, a semiconductor crystal made of cadmium and tellurium, exhibits optimal properties for application in ultra-modern CT systems.

The incoming X-ray photons produce a charge that's proportional to the X-ray energy. By applying a defined voltage over the semiconductor, these charges can be measured in a read-out ASIC.

This is how individual photons are measured, and this makes it possible to measure the energy information relative to a defined threshold in the ASIC. With the introduction of multiple threshold values, dual or more energy measurements at the detector have become a reality.

In a PCCT detector, the pixels no longer need to be sawed individually: The structure is defined by the electrical field and the read-out ASIC. This allows for smaller pixels and therefore a higher detector resolution.

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
Smaller detector pixel


Smaller detector pixel

Because there's no need for reflecting septa or comb filters, PCCT can deliver ultra high-resolution imaging with no dose penalty

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Eliminate electronic noise


Elimination of electronic noise

Advanced image quality in obese patients and low-dose scans. Stable, reproducible CT numbers (i.e., HU) for quantitative CT at a low dose.

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Intrinsic spectral sensitivity


Intrinsic spectral sensitivity

Established Dual Energy applications available for any scan. Multiple thresholds support research applications (for example, separation of different contrast penalty).

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Equal energy quanta weighting

Equal energy contribution

Contrast-rich images, increased iodine CNR. Increased iodine CNR translates into a potential dose reduction.

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Photon-counting CT direct conversion material

Thanks to its material properties, our cadmium telluride (CdTe) detector material offers a number of technical advantages in an X-ray detector, especially for photon-counting CT (PCCT).

Due to the high atomic number and high density of this material, CdTe can effectively absorb X-rays and gamma rays with energies of more than 10 keV.

CdTe is a direct semiconductor with a band gap of approximately 1.5 eV, which enables efficient absorption of X-rays.

In addition, CdTe can be operated at room temperature thanks to its wide band gap (low leakage currents, less cooling).

Electrical parameters like charge carrier mobility and electrical conductivity can be deliberately altered by intentionally adding foreign atoms (doping).

CdTe is highly resistant to radiation, which makes it appropriate for applications in medical imaging.

What effect do these differences have on image quality?

Photon-counting is the latest generation of CT detectors. Depending on the clinical question, PCCT enables a higher resolution and/or a lower dose, reduce image noise, and simplify spectral imaging.

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



From sensor to detector

- 1 Sensor material
- 2 Read-out ASIC
- 3 Ceramics
- 4 Read-out interface

The detector (also referred to as “DMS” or “data measurement system”) consists of many identical detector modules, each of which contains multiple sensors. Structured scintillator ceramics or directly converted semiconductor materials are connected with a photo diode (EID) and highly precise read-out electronics (ASIC) and positioned on a substrate (for example, ceramics).


These sensors are assembled to form modules using module carrier mechanics, a module backplane PCB, and an anti-scatter grid that eliminates scattered radiation to improve the detected signals. The assembled detector integrates these modules with more read-out and control electronics, power supply sources, and a carrier structure.

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External power supplies
1 Backplanes (power, data)
2 Modules
3 Carrier

Highly innovative front-end electronics



There are extremely stringent requirements for high resolution and precision in the development of detector electronics. Because our CT units are optimized for low radiation exposure for patients and therefore use the lowest possible X-ray dose, the detector signals are in the picoampere range. The detector electronics need to measure and process these signals with very high accuracy in order to generate a high-quality, artifact-free CT image for clinical diagnosis.


The extremely large number of measurement channels – which must work quickly and in parallel with high precision and dynamic range – requires mixed-signal ICs that are specifically developed for this purpose. Only a few chip makers worldwide are able to produce such chips.

Siemens Healthineers has worked with these partners to optimize the technology so that the distance between the photo diode and the AD converter is less than 0.5 mm in energy-integrating detectors. This is a crucial requirement for low noise and disturbances. Siemens Healthineers and their partners have set the standard for photon-counting detectors; they introduced the world's first photon-counting ASIC in a medical CT product that features high spectral accuracy and optimal high-flux capabilities.

The digital measurement data from this large number of sensor pixels is further processed in special logic modules (FPGAs) and transmitted to the image processor on the stationary part of the CT system. A modern CT detector generates a data stream of up to 40 Gbit/s.

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Detector resolution and z-coverage


When working with modern high-performance CT scanners, radiologists expect maximum resolution and fast acquisition of the largest possible regions of the human body, with the physiologically predetermined periods taken into account (for example, resting phase of the heart, breath-hold phase, and contrast medium flow).


The spatial resolution of the detector system is defined by the number of channels and pixels. To ensure high detector efficiency, it's essential that the dead zones separating the pixels be minimized.

Our many years of experience with high-precision manufacturing technologies for structuring detector ceramics allow us to ensure this for our energy-integrating detectors.

The resolution is higher in our photon-counting detectors: The number of pixels in a high-end PCCT detector is about nine times higher than in an EID detector. This is even bigger than the step from HD to UHD/4K TV resolution.

In addition, the seamless placement of the individual detector elements (arrays) in two directions enables the assembly of very large detector areas with high z-coverage, so that large volumes can be covered in a short scan time. This requires modern techniques for tiling the sensors and the corresponding read-out electronics on all sides using innovative packaging solutions and highly precise assembly methods.

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Quality characteristics

High detector efficiency and low dose

The materials in the detector elements absorb the X-ray quanta and convert them into an electrical signal. The better this works, the higher the efficiency of the detector and the lower the required radiation dose.

Short afterglow (in EID)

When the radiation is converted into light, the UFC scintillator ceramic itself starts to glow. The sooner it becomes "dark" again, the sooner it will be ready to detect new X-ray information and convert it into light.


PCCT detectors eliminate the intermediate step to visible light and are therefore free of afterglow effects.

High image resolution

The higher the resolution of the CT images, the more details of the human body can be detected and the more reliable the diagnosis. In a CT detector, just like in a digital camera, the number of pixels per unit of area plays a decisive role, and this is determined by the number of detector channels and rows. High-precision production processes are required to achieve structuring of the detector on a sub-millimeter scale.


Spectral separation (in PCCT)

Today more and more CT applications require not only a gray-scale (density-based) reconstruction of patient attenuation; they also need information about the specific material penetrated by the X-rays. This can be realized using different techniques (source-based or detector-based) – including the new photon-counting technology – to enable functional imaging (for example, separation of different materials inside the body).


Mid- and inner-ear ultra high-resolution acquisition of 0.4 mm slices with a UFC-based system (SOMATOM Force)

Courtesy of Erasmus MC, University Medical Center Rotterdam, The Netherlands

Lung fibrosis assessment with NAEOTOM Alpha, reconstructed 1 mm and 0.2 mm slices

Courtesy of Hospital ICP, Lille, France

A whole-body CT angiography acquired with TwinBeam Dual Energy shows an infra-renal abdominal aortic aneurysm (AAA) 3D rendering and slice images

Courtesy of University of Erlangen-Nuremberg, Germany

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Quality characteristics

Image quality

In addition to the quality of the X-ray tube, the detector system is a major factor that determines the image quality and efficiency of a CT system. It takes maximum-efficiency detector material and precise data processing to produce high-quality CT scans with a low radiation dose for the patient.



Rotation time

Up to 5,000 measurement values per detector element are captured during a single rotation of the CT system. Therefore, extremely fast-reacting detector materials need to be used to enable the separation of the individual measurement values from each other.

Multislice CT and other technologies


High-precision manufacturing technologies and processes as well as integrated high-speed detector electronics make it possible to drastically increase the number of measurement channels so that a large number of layers can be scanned simultaneously with every rotation. Modern multislice CT systems can cover large volumes in a very short scan time and produce very precise images. In Dual Source CTs, the scan speed and the temporal resolution can be doubled, which allows hearts with high and irregular heart rates to be scanned. The introduction of photon-counting detectors that emit even lower radiation doses has created entirely new clinical fields, thanks to their high resolution and functional imaging.

Our modern EID and PCCT detector technologies have the potential to help physicians make their diagnosis at an earlier stage of disease. It's our goal to deliver consistently high image quality at a high resolution while also minimizing the radiation dose, so the patient's radiation exposure can be kept to an absolute minimum.

Coronary CT angiography for the evaluation of coronary artery disease demonstrating NAEOTOM Alpha.Prime's high-resolution imaging at a low radiation dose.

Courtesy of Diagnostikum Graz, Austria

Chest-abdomen-pelvis imaging acquired with NAEOTOM Alpha.Prime makes it possible to examine even large patients (190 kg in this case) with high image quality and an extended field of view.

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Overview: Producing UFC and photon-counting detectors

UFC

Production

Producing the powder
Sintering and sawing the ceramic
Testing the wafer

Array production

Structuring, filling, and testing the arrays

Photon-counting

Ingot production

Raw material and ampoule preparation
Syntheses of pre-products (including polycrystalline feed and solvent)
THM single crystal growth and quality inspection

Wafer production

Wafer slicing with multiwire saw
Multi-step polishing
Final chemical/mechanical polishing (CMP)

Sensor production

Chemical and physical metal deposition
Pixel structuring using the lift-off process
Sensor separation using the dicing process
Optical inspection

UFC module production

Hybridization: Assembling the PDFEE board

Hybridization: Assembling the sensor board

Assembling the UFC module

Photon-counting module production

Hybridization I and II

Hybrid assembly

Sensor board assembly

Assembling the photon-counting module

Module testing

Testing the module

DMS production

Assembling and testing the DMS

Production

Producing the powder

The initial process step is performed in a synthesis plant. In the first step, the ceramic raw materials from the family of rare earth elements are dissolved in water. In the next step, heat is applied, causing microscopically small crystallites to grow from the oversaturated solution. After the water has been filtered out and the residue dried, a powder-like intermediate product is obtained. In the subsequent oven process, the powder is reduced at high temperatures in a gaseous atmosphere to the actual fluorescing material used in the detector. The powder now consists of our special UFC chemical formula.

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Figure 1 Synthesis
Figure 2 Reduction furnace
Figure 3 Rotating evaporator

Sintering and sawing the ceramic

The UFC properties required for CT detectors emerge only after the powder is sintered to create a completely homogeneous ceramic material. This densification is achieved in a hot-pressing process at very high temperatures of more than 1,000° C, with the simultaneous application of a compressive force of more than 30 tons.


After cooling down, the UFC pressed block is cut into wafers on a multiwire saw by an extremely thin diamond wire. The wafers are then ground to their exact final thickness. In the next step, their physical characteristics are tested.

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Production

Testing the wafer

The UFC wafers are tested on an X-ray test bench developed specifically for this purpose. A high-resolution camera is used to measure all of the parameters that are critical for later use in a CT system, including luminous efficiency (brightness), afterglow, and decay behavior. Any impurities or mechanical defects are detected in this same process.

Only the wafers that have met the specifications in all quality tests are approved for the next production step.

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

Figure 1 and 2
X-ray measurement

Array production

Structuring, filling, and testing the arrays

After the UFC wafer has been tested, a reflector layer is laminated onto it. The wafer is then structured and separated into UFC arrays. Each CT product line has its own requirements in terms of the size, number, and layout of the pixels. The required structure is first transferred to the wafer – fully automatically and in two dimensions – using an ultra high-precision saw.

These saws operate so precisely that they can split a human hair lengthwise into several parts. We only tolerate deviations of less than 15 μm . The interstitial zones between the pixels in the array are coated with a special reflecting polymer that prevents the loss of any of the light produced in the scintillator and optically separates the pixels. Conformance with the tight mechanical tolerances is verified in several test steps.

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

Figure 1 Fully automatic saw
Figure 2 Filling equipment
Figure 3 Type of array
Figure 4 Measuring equipment

Ingot production

Raw material and ampoule preparation

Before the first process step, all raw materials and supplies need to be tested for purity and chemically cleaned. The technical requirements for cleanliness are very high in the field of semiconductor production, because even impurities in the ppm range (parts per million) can have a negative impact on the resulting electrical properties.

At least 99.9999 percent (6N) purity of the raw materials must be maintained throughout the entire process chain (coating, filling in the clean room, sealing, synthesis).

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Syntheses of pre-products

Before the actual crystal growth process, several preliminary products must be produced. They include synthesizing the polycrystalline CdTe feed material – the “solvent” that’s essential for THM growth – and a perfect seed crystal from a previous crystal growth process. The general rule for every crystal growth process is “no seed – no crystal,” because it’s the seed that determines the direction of growth and the quality of the new crystal.

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Ingot production

THM single crystal growth and quality inspection

The traveling heater method (THM) is a powerful and versatile technique for growing a wide variety of materials, especially those with a high melting point or that melt non-congruently.


Major advantages of this technique include a lower growth temperature than the maximum solidus temperature of the material, resulting in fewer material impurities, lower decomposition pressures, and more perfect crystals.

After a process time of 12 weeks, all crystals need to be carefully examined for internal and external defects.

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Wafer production

Wafer slicing with multiwire saw

Wafer cutting using a multiwire saw is a precise and efficient method used in the semiconductor and photovoltaic industries to slice monocrystalline ingots.


The multiwire saw consists of multiple ultra-fine steel wires, often coated with diamond abrasives, that are passed through grooves on the guide rollers. The wires are tensioned and run at high speeds. The wires slice through the ingot, simultaneously cutting multiple thin wafers.

The sliced wafers are cleaned to remove any debris or residue from the cutting process. Then they're subjected to a quality inspection that checks for defects and imperfections.

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Wafer production

Multi-step polishing

Wafer polishing is a crucial step in semiconductor manufacturing; it involves smoothing and planarizing CdTe wafers to create a flat, defect-free surface. The wafers coming directly from the slicing process may have surface irregularities, and the initial treatment is intended to smooth out any larger surface defects.

In this multi-step procedure, the wafers are affixed to a carrier plate that's placed upside-down on a polishing plate. The plates move in opposite directions as an abrasive slurry is applied to the wafers to remove surface irregularities and create a relatively smooth surface.

Finally, the wafers are polished using a combination of chemical and mechanical forces (CMP) to achieve a higher degree of evenness and smoothness.

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Sensor production

Chemical and physical metal deposition

Plating is a process in which a thin layer of metal is chemically deposited onto a substrate surface. The substrate is immersed in a chemical bath containing metal ions, which are deposited onto the surface with or without the support of an electrical field.

Sputtering is a widely used technique in the physical vapor deposition (PVD) of metal layers. The process takes place in a vacuum chamber to minimize contamination and allow for precise control of the environment. Inert gas ions generated in a plasma are used to knock out atoms from a high-purity metal target. The metal atoms condense on the substrate and form a thin film: The sputtered layers may be as thin as $1/1,000$ the diameter of a human hair.

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Sensor production


Pixel structuring in the lift-off process

Pixel structuring in the lift-off process is a common technique used in micro-fabrication to create precise patterns on a substrate like a CdTe wafer.

In the first step, the wafer is cleaned and prepared, after which a sacrificial layer (photoresist) is deposited on the surface. The sacrificial layer is patterned using techniques like photolithography or electron beam lithography to create an inverse pattern of the desired structure. To create an electrode, a metal is deposited over the entire surface, filling the patterned areas and covering the remaining sacrificial layer. The sacrificial layer is then removed using a solvent that lifts off the excess metal and leaves behind the desired pattern on the wafer.

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Sensor separation by dicing process & optical inspection


Wafer dicing is a crucial step in semiconductor manufacturing in which individual chips or dies are separated from a wafer. This is typically done in a clean-room environment to prevent contamination. The wafer is usually mounted on a dicing tape that's affixed to a metal frame. The wafer is cut into individual dies or chips using a diamond-coated blade rotating at high speeds.

The dicing process has to be very precise in order to avoid damaging the delicate structures on the wafer.

The individual dies are inspected for quality in an automated procedure and sorted for further processing (final testing).

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Module production

Hybridization: Assembling the PDFEE board

The PDFEE board is a sensor element that combines the essential functionalities and components of the front-end electronics (FEE) with a photodiode (PD) in one chip. This chip allows for PDFEE board assembly without using conventional circuit board technology. The ASICs are placed on a ceramic substrate, soldered in a conventional reflow oven, and stabilized with underfill material. Electrical functionality is tested in the subsequent test phase, and then the PDFEE board is ready to be incorporated into a sensor board.

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
Figure 1 Chip placer
Figure 2 Inspection system

1


2

2

Figure 1 Automatic assembly system
Figure 2 Capillary underfill tool

1

Hybridization: Assembling the sensor board

In sensor board assembly, the PDFEE boards are attached to the UFC arrays using optically transparent adhesives. The components are first aligned in an automated six-axis robotic assembly system and then joined together using UV adhesive. The accuracy of these sensor-controlled assembly systems can reach $\pm 10 \mu\text{m}$.

The gaps between the components are then filled with optically conductive adhesives with the help of capillary processes and the adhesive is cured.

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Photon-counting module production

Hybridization I and II: Hybrid assembly sensor board assembly

The photon-counting sensor board is created by affixing the semiconductor material to an ASIC chip on a ceramic carrier. This is performed on a flip-chip hybridization line. In the first step, the photon-counting semiconductor is soldered with the ASICs. This requires more than 3,500 solder connections per 100 μm of surface area. In the second step, the photon-counting sensor is soldered to the ceramic carrier. The resulting sensor boards have more than 10,000 solder connections.

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


Figure 1 Sensorboard
Figure 2 Chip placer


Assembling the UFC and photon-counting modules

UFC module

Photon-counting module

The structuring of the UFC array defines the spatial resolution and, similar to a digital photo camera, the maximum achievable resolution in the CT image.

The structuring of the photo diode corresponds exactly to the pixel distribution of the UFC array. The array and the photo diode board (PDFEE) are precisely aligned relative to each other in a fully automated production step and then attached.

The anti-scatter grid is positioned above the UFC array. This ensures that disruptive scatter radiation from the side is effectively suppressed, which enhances the contrast in the CT image.

In this case, the resolution is the determining factor for structuring the counting sensors and the corresponding ASICs.

The sensor material and the ASICs are joined together precisely and then affixed to the carrier ceramic. As part of this process, an anti-scatter grid is also attached to the photon-counting sensor board, where it performs the same function as it does for the UFC modules.

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Assembling the UFC module

The sensor boards are now precisely positioned onto the module carriers and glued in place. After that, the collimator (anti-scatter) is mounted onto the sensor board.

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
Figure 1 Module assembly (manual)

Figure 1 Module assembly system (automated)
Figure 2 Sensor bank

Assembling the photon-counting module

The photon-counting module is assembled in a fully automated process using a module assembly system. First, the sensor boards are placed on the module racks with extreme precision and then glued. Then the ASGs and other components are mounted before the final module test.

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Module testing

Testing the module

Each individual module is carefully inspected and evaluated based on various criteria. The test results are recorded in a database, which is used to calculate the optimal module arrangement for each detector. This allows for the targeted selection of modules with similar properties for detector assembly, enabling uniform signal processing and high image quality.

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DMS production

Assembling and testing the DMS

In the final production step, the detector housing is equipped with the presorted modules and completed with a variety of electronic components.

The final test is now conducted under real-life conditions. The assembled detector is installed in a CT test system in order to test each specified characteristic one more time under maximum loads.

The detector isn't delivered to our internal customer, CT Final Assembly, until we're absolutely confident that it meets our highest quality standards.

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

Figure 1 DMS assembly
Figure 2 DMS final test

OEM business for UFC

The Siemens Healthineers UFC scintillator has the ability to provide high-resolution images in fractions of a second – which means optimally reliable detectors. It's not just useful for our medical applications, it's also valuable for industrial applications like security inspections at airports and container ports and for recycling and sorting many different kinds of materials.

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Security inspection

The importance of security in air travel and international freight transportation has grown significantly. Pieces of luggage, parcels, complete truckloads, and even shipping containers are now checked for illegal goods and substances using X-ray radiation. Speed and high resolution are of the utmost importance in these applications.

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Recycling and sorting

When it's equipped with our UFC scintillator material, X-ray technology can be used to identify the ingredients of materials of all kinds, sort them regardless of color and contamination, and recycle them. With this technology, substances can be separated according to their atomic weight and density. This makes it possible to separate the ingredients of a material into different material types.

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Wood industry

Thanks to UFC ceramics, round logs can be classified in such a way that the saw cuts for boards are optimized. The sorting is performed based on customer-specific product quality requirements. A 360-degree optimization for the sawing process using CT scanning allows for continuous, qualitative, and complete reconstruction of round logs in 3D. The appearance, quality, and strength can be assessed on the basis of recognizable material structures. This permits conclusions to be drawn about the impact of sawing on the finished product before the log is actually cut.

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Sustainability and environmental protection

Consumption of recycling

Innovation is the foundation of our success. We're also well aware of our responsibility for environmental protection in every step of production. From the outset, we've consistently pursued the goal of achieving continuous and sustainable improvements to our company's environmental performance through voluntary initiatives.

Over the past several years, we've implemented numerous process improvements that have enabled us to reduce the amount of raw materials required to produce a defined UFC surface by more than two-thirds.

Production of UFC ceramics also takes place in a closed recycling loop, which has a particularly positive impact thanks to the reuse of rare earth elements.

We're already developing comparable resource-saving process approaches for our photon-counting detectors.

We're proud that our products and processes have been certified under the strict EU environmental management regulations since 1996. Our environmental management system also meets the requirements of ISO 14001 and OHSAS 18001.

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
Consumption of resources

Reduction of virgin Gadolinium oxide in the production of CT detector ceramics


*Changes during several combined years shown in grey

Sustainability and environmental protection

Consumption of resources


Rare earth GdOx ^{>5N}

Consumption of raw material

Wastewater


Consumption of deionized water

Recycling of dangerous waste

Rare Earth GdOx ^{>5N}


Re-use of hazardous substance

GdOx in returns, service parts, and hazardous waste

S – Sulfur ^{>6N}

Reuse of hazardous substance

At Siemens Healthineers, we pioneer breakthroughs in healthcare. For everyone. Everywhere. Sustainably. As a market leader, we want to advance a world in which breakthroughs in healthcare create new possibilities with a minimal impact on our planet. We've been pushing the boundaries in medical technology for more than 125 years. By consistently bringing innovations to the market, we enable healthcare professionals to innovate personalized care, achieve operational excellence, and transform the system of care.

With the unique combination of our strengths in patient twinning¹, precision therapy, as well as digital, data, and artificial intelligence (AI), we are well positioned to take on the greatest challenges in healthcare. We will continue to build on these strengths to help overcome the world's most threatening diseases, enable efficient operations, and expand access to care.

Our portfolio, spanning in vitro and in vivo diagnostics to image-guided therapy and cancer care, is crucial for clinical decision-making and treatment pathways. We are committed to improving healthcare access for all, limiting our environmental impact as we pioneer breakthroughs, and engaging our diverse Healthineers to achieve this impact on a global scale.

Motivated by our purpose and guided by our values, we are building an inclusive culture, where we embrace diversity in all its forms. We are a team of 73,000 Healthineers in over 70 countries passionately pushing the boundaries of what is possible in healthcare to help improve the lives of people around the world.

The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

¹ *Personalization of diagnosis, therapy selection and monitoring, aftercare, and managing health.*

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com