

Urologic clinical cases Cios Alpha

Michael Straub, MD
Senior Urologist and Head of Endourology & Stone Center
University Hospital rechts der Isar, Munich, Germany



Urology clinical cases with Cios Alpha

Case 1:
Cystography follow-up control,
15 days after radical prostatectomy



Case 2:
Bilateral nephrostomies for urine
drainage



Case 3:
Bilateral Double J (DJ) stent
insertion to prepare radical
gastrointestinal surgery





Case 1

Cystography follow-up control, 15 days after radical prostatectomy

Cystography

- **Visualize the urinary bladder under X-ray and evaluate volume, morphology, leak tightness, and functionality**
- **During the procedure, the bladder is filled through a catheter with contrast medium, and then voided**
- **Consists of up to four exposures of the bladder:**
 - **before contrast medium is inserted**
 - **filled with contrast medium, a.p.**
 - **filled with contrast medium, lateral**
 - **after contrast medium is voided**



Condition	Radical prostatectomy – the complete removal of the prostate – is the surgical state-of-the-art treatment for prostate cancer. The procedure implies a surgical anastomosis of urethra and bladder. Cystography as a follow-up control checks for tightness or potential leakage of the anastomosis. A leak requires further treatment in order to prevent scarring and inflammation.
Patient (Gender, Age, Size, Weight)	Male, 64 years, 170 cm, 88 kg
Key procedure steps	Take exposures of the bladder in a.p. and lateral projection before and after filling with contrast medium. Assess the leak tightness of the anastomosis. Magnification for precise depiction of the leakage.
Outcome	A small leak, located right dorsal, was detected. Due to the small size of the leak, the urine catheter was removed and no further treatment was required.

Key procedure steps during the case

Prepare the contrast medium and position the FD



Take fluoroscopy scene of the filling phase




View the resulting fluoro sequences and X-ray images



Bladder filled with contrast medium a.p. projection



Bladder filled with contrast medium, lateral projection, depicting a small leak, right dorsal



Magnification with focus on the anastomosis





Case-related system highlights

Motorized **single-touch positioning*** between a.p. and lateral projection – precise and fast

Magnification of the X-ray image at the push of a button, for a better view of the anatomical detail

Brilliant image quality provided by the **flat detector technology** and the **Retina Imaging Chain** of Cios Alpha



*option



*„I was **very happy with the image quality** in this case, the system provided good contrast between the bladder and the surrounding tissue. I was able to detect a minimal leakage of the bladder outlet.”**

Michael Straub, MD

Managing Senior Physician

University Hospital rechts der Isar, Munich, Germany

*The surgeon is employed by an institution that receives/received financial support from Siemens Healthineers for collaborations. The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.


Case 2

Bilateral nephrostomies for urine drainage

General procedure description

Percutaneous nephrostomy

- Placement of a catheter percutaneously into the renal pelvis for urinary drainage
- A nephrostomy is necessary if an obstruction blocks the flow of urine from the kidney to the bladder
- During the procedure an antegrade pyelography is used to control the correct position of the tube



Condition	Prostate cancer can lead to ureter obstructions and subsequently blocked urinary flow. Nephrostomy tubes ensure urine drainage and preserve kidney function. Usually nephrostomy tubes must be exchanged every 4-6 weeks.
Patient (Gender, Age, Size, Weight)	Male, 74 years, 175 cm, 75 kg
Key procedure steps	Replacement of nephrostomy tubes under local anesthesia and in Seldinger technique. Check for correct positioning of the tubes in the renal pelvis with contrast medium and fluoroscopy. Connection of the nephrostomy tubes to an external urinary bag.
Outcome	Bilateral nephrostomy tubes were successfully exchanged. Side effect: small extravasation caused by the stiff guide wire.

Key procedure steps during the case


Place the nephrostomy tube over a guide wire into the kidney




Insert contrast medium through the nephrostomy tube for position control



Control the placement and filling of the nephrostomy tube via antegrade pyelography




First step of Seldinger technique: positioning the guide wire



Antegrade pyelography of the left kidney. Finding: nephrostomy tube correctly placed, small extravasation caused by the stiff guide wire




Antegrade pyelography of the right kidney




Case-related system highlights

The **FD mounted brake controls** allow for simple and convenient operation of the C-arm in the sterile field; here to change from one kidney to the other



Flat detector technology and **Retina Imaging Chain** for brilliant image quality; here to view very fine differentiated imaging of the renal collecting system and the extravasation behind the renal pelvis



Side-by-side viewing for comparison of the left and the right kidney



*„The patient needed bilateral nephrostomy tube exchanges. Normally this requires repositioning of the patient during the case including additional efforts for the urologist and his team. With Cios Alpha I moved the C-arm within seconds very conveniently from the left to the right kidney without moving the patient nor leaving my sterile position.” **

Michael Straub, MD

Managing Senior Physician

University Hospital rechts der Isar, Munich, Germany



*The surgeon is employed by an institution that receives/received financial support from Siemens Healthineers for collaborations. The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.



Case 3

Bilateral Double J (DJ) stent insertion to prepare radical gastrointestinal surgery

General procedure description

DJ stent insertion


- DJ stents are placed between the renal pelvis and the bladder
- Before major surgery, DJ stents are inserted to prevent inadvertent injury and to maintain urine flow in the ureter
- The special shape (double-J) of the catheter ends prevents the catheter from dislocating




Condition	A large colorectal tumor requires open surgery. DJ stents are inserted on both sides as a safety measure in preparation for gastrointestinal surgery.
Patient (Gender, Age, Size, Weight)	Female, 70 years, 165 cm, 68 kg
Key procedure steps	Insertion of DJ stents under endoscopic and fluoroscopic control. Verification and documentation of proper placement of the DJ stent in the kidney and bladder region.
Outcome	Bilateral DJ stents are successfully inserted and the patient is prepared for the upcoming surgical intervention.

Key procedure steps during the case


Prepare the DJ stent on the sterile OR table




Insert the DJ stent in Seldinger technique through the endoscope into the urinary tract




Place the DJ stent under fluoroscopic and endoscopic view




Endoscopy image of the DJ stent passing the right orifice of the bladder



Fluorolooop documentation of the placement of 2nd DJ stent



Placement in the bladder




Case-related system highlights

Imaging of the complete urinary tract thanks to the large **30 cm x 30 cm flat detector**



SmartView provides comfortable side-by-side viewing of X-ray and endoscopic images



With **HD EndoStore**, images from all modalities can be stored in one patient folder for complete case documentation





*„For an endourologist it is very important to have the fluoroscopic and endoscopic images in the same view. I have to work with both imaging modalities simultaneously for stent placement. This is a big step towards **efficient work in the OR.**” **

Michael Straub, MD

Managing Senior Physician

University Hospital rechts der Isar, Munich, Germany

*The surgeon is employed by an institution that receives/received financial support from Siemens Healthineers for collaborations. The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.



*„ After 15 years of experience in endourology I worked with many different C-arms. In my opinion Cios Alpha is very comfortable with respect to image quality and flexibility in movements in the OR. It offers all possibilities that I need in modern endourology.” **

Michael Straub, MD

Managing Senior Physician

University Hospital rechts der Isar, Munich, Germany

*The surgeon is employed by an institution that receives/received financial support from Siemens Healthineers for collaborations. The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.

Disclaimer

Cios Alpha is not commercially available in all countries. Due to regulatory reasons its future availability cannot be guaranteed. Please contact your local Siemens organization for further details.

The statements by Siemens' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

The speaker is employed by an institution that provides Siemens product reference services for compensation pursuant to a written agreement.



Thank you for your enthusiasm!

Siemens Healthineers Headquarters

Siemens Healthcare GmbH

Henkestr. 127

91052 Erlangen

Germany

Phone: +49 9131 84-0

[siemens-healthineers.com](https://www.siemens-healthineers.com)

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this presentation are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

The customers cited are employed by an institution that might provide Siemens product reference services, R&D collaboration or other relationship for compensation pursuant to a written agreement.