

SIEMENS



DICOM Conformance Statement

**syngo®** MR XA30A

# 1 Conformance Statement Overview

syngo® MR conforms to the DICOM Standard and supports the network services as described in Table 1: Network Services and the media services as described in Table 2 - Media Services.

**Table 1: Network Services**

| SOP Classes  | SOP Class UID                 | User of Service (SCU) |      | Provider of Service (SCP) |         |
|--|-------------------------------|-----------------------|------|---------------------------|---------|
| Verification   |                               |                       |      |                           |         |
| Verification   | 1.2.840.10008.1.1             | Yes                   |      | Yes                       |         |
| SOP Classes created by syngo® MR                           |                               |                       |      |                           |         |
|  |                               | Create                | Send | Store                     | Display |
| MR Image Storage   | 1.2.840.10008.5.1.4.1.1.4     | Yes                   | Yes  | Yes                       | Yes     |
| Enhanced MR Image Storage                                  | 1.2.840.10008.5.1.4.1.1.4.1   | Yes                   | Yes  | Yes                       | Yes     |
| MR Spectroscopy Storage                                    | 1.2.840.10008.5.1.4.1.1.4.2   | Yes                   | Yes  | Yes                       | Yes     |
| Enhanced MR Color Image Storage                            | 1.2.840.10008.5.1.4.1.1.4.3   | Yes                   | Yes  | Yes                       | Yes     |
| Secondary Capture Image Storage                            | 1.2.840.10008.5.1.4.1.1.7     | Yes                   | Yes  | Yes                       | Yes     |
| Multi-frame Grayscale Word Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.3   | Yes                   | Yes  | Yes                       | Yes     |
| Multi-frame True Color Secondary Capture Image Storage     | 1.2.840.10008.5.1.4.1.1.7.4   | Yes                   | Yes  | Yes                       | Yes     |
| Grayscale Softcopy Presentation State Storage              | 1.2.840.10008.5.1.4.1.1.11.1  | Yes                   | Yes  | Yes                       | Yes     |
| Raw Data Storage   | 1.2.840.10008.5.1.4.1.1.66    | Yes                   | Yes  | Yes                       | No      |
| Spatial Registration Storage                               | 1.2.840.10008.5.1.4.1.1.66.1  | Yes                   | Yes  | Yes                       | No      |
| Segmentation Storage                                       | 1.2.840.10008.5.1.4.1.1.66.4  | Yes                   | Yes  | Yes                       | Yes     |
| Real World Value Mapping Storage                           | 1.2.840.10008.5.1.4.1.1.67    | Yes                   | Yes  | Yes                       | No      |
| Enhanced SR Storage  | 1.2.840.10008.5.1.4.1.1.88.22 | Yes                   | Yes  | Yes                       | No      |
| Comprehensive SR Storage                                   | 1.2.840.10008.5.1.4.1.1.88.33 | Yes                   | Yes  | Yes                       | No      |
| Encapsulated PDF Storage SOP Class                         | 1.2.840.10008.5.1.4.1.1.104.1 | Yes                   | Yes  | Yes                       | Yes     |
| SOP Classes managed by syngo® MR                           |                               |                       |      |                           |         |
| Computed Radiography Image Storage                         | 1.2.840.10008.5.1.4.1.1.1     | No                    | Yes  | Yes                       | Yes     |
| Digital X-Ray Image Storage - For Presentation             | 1.2.840.10008.5.1.4.1.1.1.1   | No                    | Yes  | Yes                       | Yes     |
| Digital X-Ray Image Storage – For Processing               | 1.2.840.10008.5.1.4.1.1.1.1.1 | No                    | Yes  | Yes                       | Yes     |
| Digital Mammography X-Ray Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.2   | No                    | Yes  | Yes                       | Yes     |
| Digital Mammography X-Ray Image Storage – For Processing   | 1.2.840.10008.5.1.4.1.1.1.2.1 | No                    | Yes  | Yes                       | Yes     |
| CT Image Storage   | 1.2.840.10008.5.1.4.1.1.2     | No                    | Yes  | Yes                       | Yes     |
| Enhanced CT Image Storage                                  | 1.2.840.10008.5.1.4.1.1.2.1   | No                    | Yes  | Yes                       | Yes     |

| SOP Classes  | SOP Class UID                  | User of Service (SCU) |     | Provider of Service (SCP) |     |
|--|--------------------------------|-----------------------|-----|---------------------------|-----|
| Ultrasound Multi-frame Image Storage                       | 1.2.840.10008.5.1.4.1.1.3.1    | No                    | Yes | Yes                       | Yes |
| Ultrasound Image Storage                                   | 1.2.840.10008.5.1.4.1.1.6.1    | No                    | Yes | Yes                       | Yes |
| Multi-frame Single Bit Secondary Capture Image Storage     | 1.2.840.10008.5.1.4.1.1.7.1    | No                    | Yes | Yes                       | Yes |
| Multi-frame Grayscale Byte Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.2    | No                    | Yes | Yes                       | Yes |
| 12-lead ECG Waveform Storage                               | 1.2.840.10008.5.1.4.1.1.9.1.1  | No                    | Yes | Yes                       | No  |
| General ECG Waveform Storage                               | 1.2.840.10008.5.1.4.1.1.9.1.2  | No                    | Yes | Yes                       | No  |
| Ambulatory ECG Waveform Storage                            | 1.2.840.10008.5.1.4.1.1.9.1.3  | No                    | Yes | Yes                       | No  |
| Hemodynamic Waveform Storage                               | 1.2.840.10008.5.1.4.1.1.9.2.1  | No                    | Yes | Yes                       | No  |
| Cardiac Electrophysiology Waveform Storage                 | 1.2.840.10008.5.1.4.1.1.9.3.1  | No                    | Yes | Yes                       | No  |
|  |                                |                       |     |                           |     |
| Color Softcopy Presentation State Storage SOP Class        | 1.2.840.10008.5.1.4.1.1.11.2   | No                    | Yes | Yes                       | No  |
| Pseudo-Color Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.3   | No                    | Yes | Yes                       | No  |
| Blending Softcopy Presentation State Storage               | 1.2.840.10008.5.1.4.1.1.11.4   | No                    | Yes | Yes                       | No  |
| X-Ray Angiographic Image Storage                           | 1.2.840.10008.5.1.4.1.1.12.1   | No                    | Yes | Yes                       | Yes |
| Enhanced XA Image Storage                                  | 1.2.840.10008.5.1.4.1.1.12.1.1 | No                    | Yes | Yes                       | Yes |
| X-Ray Radiofluoroscopic Image Storage                      | 1.2.840.10008.5.1.4.1.1.12.2   | No                    | Yes | Yes                       | Yes |
| Enhanced XRF Image Storage                                 | 1.2.840.10008.5.1.4.1.1.12.2.1 | No                    | Yes | Yes                       | Yes |
| X-Ray 3D Angiographic Image Storage                        | 1.2.840.10008.5.1.4.1.1.13.1.1 | No                    | Yes | Yes                       | Yes |
| Breast Tomosynthesis Image Storage                         | 1.2.840.10008.5.1.4.1.1.13.1.3 | No                    | Yes | Yes                       | Yes |
| Spatial Fiducials Storage                                  | 1.2.840.10008.5.1.4.1.1.66.2   | No                    | Yes | Yes                       | No  |
| Deformable Spatial Registration SOP Class                  | 1.2.840.10008.5.1.4.1.1.66.3   | No                    | Yes | Yes                       | No  |
| Surface Segmentation Storage                               | 1.2.840.10008.5.1.4.1.1.66.5   | No                    | Yes | Yes                       | No  |
| Basic Text SR Storage                                      | 1.2.840.10008.5.1.4.1.1.88.11  | No                    | Yes | Yes                       | No  |
| Procedure Log Storage Storage                              | 1.2.840.10008.5.1.4.1.1.88.40  | No                    | Yes | Yes                       | No  |
| Mammography CAD SR Storage                                 | 1.2.840.10008.5.1.4.1.1.88.50  | No                    | Yes | Yes                       | No  |
| Key Object Selection Document Storage                      | 1.2.840.10008.5.1.4.1.1.88.59  | No                    | Yes | Yes                       | No  |
| X-Ray Radiation Dose SR Storage                            | 1.2.840.10008.5.1.4.1.1.88.67  | No                    | Yes | Yes                       | No  |

| SOP Classes                                     | SOP Class UID                 | User of Service (SCU) |     | Provider of Service (SCP) |     |
|---|-------------------------------|-----------------------|-----|---------------------------|-----|
| Positron Emission Tomography Image Storage      | 1.2.840.10008.5.1.4.1.1.128   | No                    | Yes | Yes                       | Yes |
| RT Image Storage                                | 1.2.840.10008.5.1.4.1.1.481.1 | No                    | Yes | Yes                       | Yes |
| RT Dose Storage                                 | 1.2.840.10008.5.1.4.1.1.481.2 | No                    | Yes | Yes                       | No  |
| RT Structure Set Storage                        | 1.2.840.10008.5.1.4.1.1.481.3 | No                    | Yes | Yes                       | No  |
| RT Beams Treatment Record Storage               | 1.2.840.10008.5.1.4.1.1.481.4 | No                    | Yes | Yes                       | No  |
| RT Plan Storage                                 | 1.2.840.10008.5.1.4.1.1.481.5 | No                    | Yes | Yes                       | No  |
| RT Treatment Summary Record Storage             | 1.2.840.10008.5.1.4.1.1.481.7 | No                    | Yes | Yes                       | No  |
| Transfer (Private SOP Class)                    |                               |                       |     |                           |     |
| Syngo Non-Image Storage                         |                               | Yes                   |     | Yes                       |     |
| Storage Commitment                              |                               |                       |     |                           |     |
| Storage Commitment Push Model SOP Class         | 1.2.840.10008.1.20.1          | Yes                   |     | Yes                       |     |
| Worklist Management                             |                               |                       |     |                           |     |
| Modality Worklist Information Model - FIND      | 1.2.840.10008.5.1.4.31        | Yes                   |     | No                        |     |
| Modality Performed Procedure Step SOP Class     | 1.2.840.10008.3.1.2.3.3       | Yes                   |     | No                        |     |
| Query/Retrieve                                  |                               |                       |     |                           |     |
| Patient Root Q/R Information Model - FIND       | 1.2.840.10008.5.1.4.1.2.1.1   | Yes                   |     | Yes                       |     |
| Patient Root Q/R - Information Model - MOVE     | 1.2.840.10008.5.1.4.1.2.1.2   | Yes                   |     | Yes                       |     |
| Study Root Q/R - Information Model - FIND       | 1.2.840.10008.5.1.4.1.2.2.1   | Yes                   |     | Yes                       |     |
| Study Root Q/R - Information Model - MOVE       | 1.2.840.10008.5.1.4.1.2.2.2   | Yes                   |     | Yes                       |     |
| Patient/Study Only Q/R - Information Model FIND | 1.2.840.10008.5.1.4.1.2.2.1   | Yes                   |     | Yes                       |     |
| Patient/Study Only Q/R - Information Model MOVE | 1.2.840.10008.5.1.4.1.2.3.2   | Yes                   |     | Yes                       |     |
| Print Management                                |                               |                       |     |                           |     |
| Basic Grayscale Print Management Meta           | 1.2.840.10008.5.1.1.9         | Yes                   |     | No                        |     |
| Basic Color Print Management Meta               | 1.2.840.10008.5.1.1.18        | Yes                   |     | No                        |     |
| Basic Film Sesssion                             | 1.2.840.10008.5.1.1.1         | Yes                   |     | No                        |     |
| Basic Film Box                                  | 1.2.840.10008.5.1.1.2         | Yes                   |     | No                        |     |
| Basic Grayscale Image Box                       | 1.2.840.10008.5.1.1.4         | Yes                   |     | No                        |     |
| Basic Color Image Box                           | 1.2.840.10008.5.1.1.4.1       | Yes                   |     | No                        |     |
| Printer   | 1.2.840.10008.5.1.1.16        | Yes                   |     | No                        |     |
| Print Job                                       | 1.2.840.10008.5.1.1.14        | Yes                   |     | No                        |     |
| Presentation LUT                                | 1.2.840.10008.5.1.1.23        | Yes                   |     | No                        |     |

**Table 2 - Media Services**

| Media Storage Application Profile | Write Files (FSC or FSU) | Read Files (FSR) |
|-----------------------------------|--------------------------|------------------|
| Compact Disk - Recordable         |                          |                  |

|                  |     |     |
|------------------|-----|-----|
| STD-GEN-CD       | Yes | Yes |
| AUG-GEN-CD       | Yes | Yes |
| <b>DVD</b>       |     |     |
| AUG-GEN-DVD      | Yes | Yes |
| AUG- GEN-DVD-J2K | Yes | Yes |
| STD-GEN-DVD      | Yes | Yes |
| STD-GEN-DVD-J2K  | Yes | Yes |
| <b>USB</b>       |     |     |
| AUG- GEN-USB-J2K | Yes | Yes |
| STD-GEN-USB-J2K  | Yes | Yes |

**Table 3 - Implementation Identifying Information**

| <b>Name</b>                 | <b>Value</b>           |
|-----------------------------|------------------------|
| Application Context Name    | 1.2.840.100008.3.1.1.1 |
| Implementation Class UID    | 1.3.12.2.1107.5.2      |
| Implementation Version Name | SYNGO_MR_XA30A         |

## 2 Table of Contents

|            |   |           |
|------------|---|-----------|
| <b>1</b>   | <b>Conformance Statement Overview .....</b>                     | <b>2</b>  |
| <b>2</b>   | <b>Table of Contents .....</b>                                  | <b>6</b>  |
| <b>3</b>   | <b>Introduction.....</b>  | <b>8</b>  |
| 3.1        | Scope and Field .....   | 8         |
| 3.2        | Audience .....  | 8         |
| 3.3        | Remarks .....   | 8         |
| 3.4        | Definitions, Terms and Abbreviations .....                      | 8         |
| 3.5        | References .....  | 9         |
| <b>4</b>   | <b>Networking .....</b>   | <b>10</b> |
| <b>4.1</b> | <b>Implementation Model.....</b>                                | <b>10</b> |
| 4.1.1      | Application Data Flow .....                                     | 10        |
| 4.1.2      | Functional Definitions of Application Entities .....            | 13        |
| 4.1.3      | Sequencing of Activities.....                                   | 14        |
| <b>4.2</b> | <b>Application Entity Specification .....</b>                   | <b>18</b> |
| 4.2.1      | Verification AE Specification .....                             | 18        |
| 4.2.2      | Storage AE Specification .....                                  | 20        |
| 4.2.3      | Storage Commitment AE Specification .....                       | 26        |
| 4.2.4      | Query/Retrieve AE Specification .....                           | 29        |
| 4.2.5      | Modality Worklist AE Specification .....                        | 34        |
| 4.2.6      | Modality Performed Procedure Step AE Specification .....        | 40        |
| 4.2.7      | Print AE Specification .....                                    | 43        |
| 4.2.8      | Print SCU Specification .....                                   | 53        |
| <b>4.3</b> | <b>Network Interfaces .....</b>                                 | <b>62</b> |
| 4.3.1      | Physical Network Interface .....                                | 62        |
| 4.3.2      | Additional Protocols .....                                      | 62        |
| 4.3.3      | IPv4 and IPv6 Support.....                                      | 62        |
| <b>4.4</b> | <b>Configuration.....</b>                                       | <b>62</b> |
| 4.4.1      | AE Title/Presentation Address Mapping .....                     | 62        |
| 4.4.2      | Parameters .....  | 63        |
| <b>5</b>   | <b>Media Interchange .....</b>                                  | <b>65</b> |
| <b>5.1</b> | <b>Implementation Model.....</b>                                | <b>65</b> |
| 5.1.1      | Application Data Flow Diagram .....                             | 65        |
| 5.1.2      | Functional definitions of AEs .....                             | 66        |
| 5.1.3      | Sequencing of Real-World Activities .....                       | 66        |
| 5.1.4      | File Meta Information for Implementation Class and Version..... | 67        |
| <b>5.2</b> | <b>AE SPECIFICATIONS.....</b>                                   | <b>67</b> |
| 5.2.1      | Media Storage AE – Specification .....                          | 67        |
| <b>5.3</b> | <b>AUGMENTED AND PRIVATE APPLICATION PROFILES.....</b>          | <b>69</b> |
| 5.3.1      | Augmented Application Profiles.....                             | 69        |
| <b>5.4</b> | <b>MEDIA CONFIGURATION.....</b>                                 | <b>69</b> |

|          |  |                   |
|----------|--|-------------------|
| <b>6</b> | <b><i>Support of Extended Character Sets .....</i></b> | <b><i>70</i></b>  |
| <b>7</b> | <b><i>Attribute confidentiality profiles.....</i></b>  | <b><i>73</i></b>  |
| 7.1      | De-identification .....                                | 73                |
| <b>8</b> | <b><i>Security .....</i></b>                           | <b><i>80</i></b>  |
| 8.1      | Security Profiles .....                                | 80                |
| 8.1.1    | Time Synchronization Profiles .....                    | 80                |
| 8.1.2    | Basic TLS Secure Transport Connection Profile .....    | 80                |
| 8.2      | Association Level Security .....                       | 80                |
| 8.3      | Application Level Security .....                       | 80                |
| <b>9</b> | <b><i>Annexes.....</i></b>                             | <b><i>81</i></b>  |
| 9.1      | IOD Contents .....                                     | 81                |
| 9.1.1    | Enhanced MR Image IOD.....                             | 81                |
| 9.1.2    | Enhanced MR Color Image IOD .....                      | 90                |
| 9.1.3    | MR Spectroscopy IOD .....                              | 90                |
| 9.1.4    | MR Image .....   | 98                |
| 9.1.5    | Raw Data .....   | 101               |
| 9.1.6    | Multi-frame Grayscale Word SC Image.....               | 101               |
| 9.1.7    | Secondary Capture Image.....                           | 104               |
| 9.1.8    | Grayscale Softcopy Presentation State.....             | 105               |
| 9.1.9    | Evidence Documents.....                                | 111               |
| 9.2      | Data Dictionary of Private Attributes .....            | 111               |
| 9.3      | Grayscale Image Consistency .....                      | 112               |
| 9.4      | DICOM Print SCU – detailed status displays .....       | 113               |
| 9.4.1    | Common Status Information .....                        | 113               |
| 9.4.1    | Additional DICOM Execution Status Information.....     | 116               |
|          | <b><i>Annex A: Index of Tables.....</i></b>            | <b><i>117</i></b> |
|          | <b><i>Annex B: Table of Figures .....</i></b>          | <b><i>119</i></b> |

## 3 Introduction

### 3.1 Scope and Field

This DICOM Conformance Statement refers to SIEMENS MR Products using software *syngo®* MR XA30A. The following table relates *syngo®* MR XA30A software versions to SIEMENS *syngo®* MR products.

| Software Name  | SIEMENS MR Product |
|----------------|--------------------|
| syngo MR XA30A | MAGNETOM Aera      |
| syngo MR XA30A | MAGNETOM Skyra     |
| syngo MR XA30A | MAGNETOM Prisma    |

The syngo® MR product is a “syngo®-speaking<sup>a</sup>” Imaging Modality or workstation. The syngo® MR product is designed to be integrated into an environment of medical DICOM-based devices.

### 3.2 Audience

This document is intended for hospital staff, health system integrators, software designers or implementers. It is assumed that the reader has a working understanding of DICOM.

### 3.3 Remarks

The scope of this DICOM Conformance Statement is to facilitate integration between syngo® MR and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [1]. DICOM by itself does not guarantee interoperability.

The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

The comparison of conformance statements is the first step towards assessing interconnectivity and interoperability between syngo® MR and other DICOM conformant equipment.

Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens representative for the most recent product information.

### 3.4 Definitions, Terms and Abbreviations

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM standard.

Additional Abbreviations and terms are as follows:

|       |  |
|-------|--|
| AE    | DICOM Application Entity                           |
| AET   | Application Entity Title                           |
| ASCII | American Standard Code for Information Interchange |
| CSE   | Customer Service Engineer                          |

<sup>a</sup> *syngo* is a registered trademark of Siemens Healthcare GmbH



|       |  |
|-------|--|
| DCS   | DICOM Conformance Statement                    |
| DICOM | Digital Imaging and Communications in Medicine |
| FSC   | File Set Creator                               |
| FSR   | File Set Reader                                |
| FSU   | File Set Updater                               |
| GSDF  | Grayscale Standard Display Function            |
| IOD   | DICOM Information Object Definition            |
| ISO   | International Standard Organization            |
| n. a. | not applicable                                 |
| NEMA  | National Electrical Manufacturers Association  |
| O     | Optional Key Attribute                         |
| PDU   | DICOM Protocol Data Unit                       |
| R     | Required Key Attribute                         |
| SCU   | DICOM Service Class User (DICOM client)        |
| SCP   | DICOM Service Class Provider (DICOM Server)    |
| SOP   | DICOM Service-Object Pair                      |
| SPS   | Scheduled Procedure Step                       |
| SR    | Structured Report                              |
| TFT   | Thin Film Transistor (Display)                 |
| TID   | Template ID                                    |
| U     | Unique Key Attribute                           |
| UID   | Unique Identifier                              |
| UTF-8 | Unicode Transformation Format-8                |
| VR    | Value Representation                           |

## 3.5 References

- [1] Digital Imaging and Communications in Medicine (DICOM PS3.1-3.20 2016a), National Electrical Manufacturers Association (NEMA), <http://medical.nema.org>
- [2] Integrating the Healthcare Enterprise – IHE Radiology Technical Framework – <http://www.ihe.net>

## **4      Networking**

### **4.1      Implementation Model**

syngo® MR supports storing DICOM images to remote nodes like workstations or Archiving Systems. Using the Storage Commitment Service it can request safe keeping of previously stored instances from an Archiving system. Additionally, the syngo® MR can query remote notes, retrieve and store selected instances from that node. Using the Modality Worklist service the syngo® MR can query a HIS/RIS for scheduled procedures. Performed procedure status and other procedure data can be returned to the HIS/RIS using the Modality Performed Procedure Step (MPPS) Service. Furthermore, printing of color and grayscale images is supported.

#### **4.1.1      Application Data Flow**

The following figure provides a functional overview of the syngo® MR Application Entities (AE). Relationships are shown between user-invoked activities (in the circles at the left of the AEs) and the associated real-world activities provided by DICOM service providers (in the circles at the right of the AEs)




Figure 1: syngo® MR DICOM Data Flow diagram – Acquisition Workflow




Figure 2: syngo® MR DICOM Data Flow diagram - Printing

## 4.1.2 Functional Definitions of Application Entities

The SCP components of the Application Entities of the syngo® MR operate as background server processes. They exist as soon as the system is powered up and wait for association requests. Upon accepting an association with a negotiated Presentation Context they start to receive and process the requests described in the following sections.

The SCU components of the Application Entity are invoked upon requests from the user interface or indirect by trigger from internal processes.

### 4.1.2.1 Functional Definition of Verification AE

syngo® MR supports the Verification service as a SCP and SCU. As a SCU, Verification can be activated from the Administrator Portal during system configuration.

As a SCP of the Verification Service the syngo® MR processes and responds to incoming verification requests using the C-ECHO-Response.

### 4.1.2.2 Functional Definition of Storage AE

The syngo® MR Storage SCU is invoked either directly by the user, by an auto-archive trigger or internally by the Query/Retrieve Application Entity that is responsible for processing retrieve requests. The job consists of data describing the composite objects selected for storage and the destination Application Entity Title. An association is negotiated with the destination Application Entity and the image data is transferred using the DIMSE C-STORE -Request. The transfer status is reported to the initiator of the Storage request.

The Storage SCP component of syngo® MR starts to receive the Composite Image Objects and import them into the database after accepting an association with a negotiated Presentation Context. The system can be configured in such a way, that Responses to the Storage Request are sent immediately after reception of the data, after persistent storage on the hard disc or after storage and indexing in the local database.

### 4.1.2.3 Functional Definition of the Storage Commitment AE

If configured, syngo® MR can serve as a SCU for the DICOM Storage Commitment service. Upon successful completion of a storage job, the system uses the DIMSE N-ACTION Request to request storage commitment from a remote DICOM Storage Commitment SCP. This can either be the same as the storage destination or a different system depending on the system configuration. Storage Commitment Requests are sent after a configurable delay after storing the objects. syngo® MR can receive the N-EVENT-REPORT-Request on the same or a different association.

### 4.1.2.4 Functional Definition of Query/Retrieve AE

syngo® MR supports DICOM Query/Retrieve as a SCU: The user can initiate a query to a remote node using the C-Find Request. After matching the specified keys, the remote Query /Retrieve SCP uses the C-Find Response to return the results of its search, which will be displayed to the user. Depending on user action the syngo® MR Query/Retrieve DICOM SCU sends a C-MOVE Request to initiate a C-STORE sub-operation on the SCP to start an image transfer from remote Storage SCU (running on Query/Retrieve SCP) to the system's Storage SCP.

syngo® MR supports the following query models:

- Study Root Query Model.
- Patient Root Query Model
- Patient/Study Only Query Model

Furthermore, the SCU services may issue relational queries, if supported by the remote Query/Retrieve SCP node and required by the querying Application.

The syngo® MR DICOM Query/Retrieve SCP accepts C-Find Request, queries the local database based on the provided matching keys and returns the matches using the C-FIND Response. Depending on further request from the remote Query/Retrieve SCU, syngo® MR responds to C-MOVE Requests by initiating a C-STORE sub-operation to send image objects to the Storage SCP of the querying system.

#### **4.1.2.5 Functional Definition of Modality Worklist AE**

The syngo® MR Modality Worklist SCU issues DICOM Modality Worklist requests using C-FIND Requests. The results in the C-FIND Response are stored in internal database. The provided Patient and Procedure information is used for patient registration prior to starting an exam.

#### **4.1.2.6 Functional Definition of Modality Performed Procedure Step SCU AE**

The syngo® MR MPPS SCU uses the N-CREATE Request to inform an Information System that a procedure step is IN PROGRESS.

syngo® MR MPPS SCU uses the N-SET Request to inform the Information System about the finalization of the Procedure Step, using either a status of COMPLETED or DISCONTINUED.


#### **4.1.2.7 Functional Definition of Print AE**

The Print SCU of syngo® MR is invoked by the user interface to setup film-sheet layout and whenever an image is ready to be printed on film. The Print SCU will hold and maintain all data needed to compile a complete film-sheet from the data (images, layout, configuration) received. Whenever a film-sheet is ready to print the related data is used to supply the Information to the SOP Classes of the Print Management Service Class. A queue is maintained, in order to immediately store several film-sheets in case of resource problems on printer. The SCU will only supply and require the mandatory SOP Classes of the Print Management Service Class.

### **4.1.3 Sequencing of Activities**


This section describes the sequencing of Real-World Activities performed by syngo® MR Entities using a UML sequence diagram. Real-World Activities are depicted as vertical bars and arrows show the events exchanged between them.

#### 4.1.3.1 System Configuration



**Figure 3: Sequence Diagram for Real World Activities - System Configuration**


### 4.1.3.2 Acquisition Workflow



**Figure 4: Sequence Diagram for Real World Activities - Acquisition workflow**



### 4.1.3.3 Printing Workflow



**Figure 5: Sequence Diagram for Real World Activities - Printing**

## 4.2 Application Entity Specification

This section outlines the specifications for each of the Application Entities that are part of syngo® MR.

### 4.2.1 Verification AE Specification

#### 4.2.1.1 SOP Classes

The Verification AE of syngo® MR provides standard conformance to the Verification SOP Class listed in "Table 1: Network Services" section "Verification" in the ["Conformance Statement Overview"](#).

#### 4.2.1.2 Association Policy

The syngo® MR Admin Portal attempts to open an association for verification request whenever the "Echo" function is activated.

**Table 4: Association Policies**

|  |                       |
|--|-----------------------|
| <b>Application Context Name</b>  | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>  | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

##### 4.2.1.2.1 Asynchronous Nature

syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any size is supported.

**Table 5: Asynchronous Nature as an Association Initiator**

|  |    |
|--|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|--|----|

##### 4.2.1.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to "Table 3 - Implementation Identifying Information" in the ["Conformance Statement Overview"](#).

#### 4.2.1.3 Association Initiation Policy

##### 4.2.1.3.1 Activity – "Send Verification" Request

###### 4.2.1.3.1.1 Description and Sequencing of Activity

syngo® MR serves as a SCU of the Verification Service Class. A C-ECHO Request is initiated by the Administrator Portal whenever "Verification" is requested. If an association to a remote Application Entity is successfully established, Verification with the configured AET is requested via the open association. If the C-ECHO Response from the remote Application contains a status other than "Success" this will be indicated to the user and the association is closed.

<sup>1</sup> Default, the value is configurable

#### 4.2.1.3.1.2 Proposed Presentation Contexts

Table 6 - Presentation Context Table "Verification" below lists the supported presentation contexts for verification requests.

**Table 6 - Presentation Context Table "Verification"**

| Presentation Context Table – "Verification" |                   |                           |                   |      |                      |
|---|-------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax                             |                   | Transfer Syntax           |                   | Role | Extended Negotiation |
| Name  | UID               | Name List                 | UID List          |      |                      |
| Verification                                | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |

#### 4.2.1.3.1.3 SOP Specific Conformance – Verification SCU

The application conforms to the definition of the Verification SCU in accordance to the DICOM Standard.

#### 4.2.1.4 Association Acceptance Policy

##### 4.2.1.4.1 Activity – "Receive Verification Request"

##### 4.2.1.4.1.1 Description and Sequencing of Activity

syngo® MR serves as a SCP of the Verification Service Class. If the Verification SCP accepts an association, it will respond to C-ECHOC Requests. If the Called AE Title does not match any pre-configured AE Title shared by SCP, the association will be rejected.

##### 4.2.1.4.1.2 Accepted Presentation Contexts

The syngo® MR DICOM application will accept Presentation Contexts as shown in the following table:

**Table 7 - Presentation Context Table "Verification"**

| Presentation Context Table – "Verification" |                   |                           |                   |      |                      |
|---|-------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax                             |                   | Transfer Syntax           |                   | Role | Extended Negotiation |
| Name  | UID               | Name List                 | UID List          |      |                      |
| Verification                                | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP  | None                 |

##### 4.2.1.4.1.3 SOP Specific Conformance – Verification SCP

The ECHO-SCP provides standard conformance to the Verification Service Class.

## 4.2.2 Storage AE Specification

### 4.2.2.1 SOP Classes

The Storage AE provides Standard Conformance to the the SOP Classes listed in “Table 1: Network Services” section “SOP Classes Created by syngo® MR” and “SOP Classes Managed by syngo® MR” in the [“Conformance Statement Overview”](#).

### 4.2.2.2 Association Policy

**Table 8: Association Policies**

|  |                       |
|--|-----------------------|
| <b>Application Context Name</b>  | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>  | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.2.2.1 Asynchronous Nature

syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 9: Asynchronous Nature as an Association Initiator**

|  |    |
|--|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|--|----|

#### 4.2.2.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

### 4.2.2.3 Association Initiation Policy

#### 4.2.2.3.1 Activity – “Send Storage Request”

##### 4.2.2.3.1.1 Description and Sequencing of Activities

syngo® MR serves as a SCU of the Storage Service Class. The Storage SCU is triggered by the transfer job queue or by an external retrieve request. An association request is sent to the destination AE. Upon successful negotiation of a Presentation Context, the transfer is started. Objects will be transferred sequentially on the same open association

<sup>1</sup> Default, the value is configurable

syngo® MR does not provide any automated retry mechanism.

#### 4.2.2.3.1.2 Proposed Presentation Contexts

For all Image Objects listed in Table 1 in the Conformance Statement Overview the Transfer Syntaxes marked with “yes” in the Image Objects Column of the table below are supported.

For all Non-Image Objects listed in Table 1 in the Conformance Statement Overview the Transfer Syntaxes marked with “yes” in the Non-Image Objects Column of the table below are supported.

For a distinction between Image and Non-Image Objects please refer to the DICOM Standard PS3.3 Section A.1.4 “Overview of the Composite IOD Module Content”.

**Table 10: Proposed Presentation Contexts for Storage**

| UID value              | Transfer Syntax  | Image Objects | Non-Image Objects |
|------------------------|--|---------------|-------------------|
| 1.2.840.10008.1.2      | Implicit Value Representation Little Endian native                                       | yes           | yes               |
| 1.2.840.10008.1.2.1    | Explicit Value Representation Little Endian native                                       | yes           | yes               |
| 1.2.840.10008.1.2.2    | Explicit Value Representation Big Endian   | yes           | yes               |
| 1.2.840.10008.1.2.4.50 | JPEG Baseline (Process 1) lossy compressed   | yes           | no                |
| 1.2.840.10008.1.2.4.51 | JPEG Extended (Process 2 & 4) lossy compressed   | yes           | no                |
| 1.2.840.10008.1.2.4.70 | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) lossless compressed | yes           | no                |
| 1.2.840.10008.1.2.4.90 | JPEG 2000 Image Compression (Lossless Only) compressed                                   | yes           | no                |
| 1.2.840.10008.1.2.4.91 | JPEG 2000 Image Compression lossy compressed   | yes           | no                |
| 1.2.840.10008.1.2.5    | RLE Lossless compressed  | yes           | no                |

Depending on the configuration, the Storage SCU will choose a compressed or uncompressed Transfer Syntax among those accepted by the SCP. The Transfer Syntax chosen is the preferred one among the compressed and uncompressed ones. The preference order is the order of occurrence in the configuration. It is possible to configure for a specific node, which Transfer Syntax shall be used and which one shall be excluded. The configuration can even be extended, based on the combination of SOP Classes and supported Transfer Syntaxes.

An instance will be JPEG lossless (Process 1 and Process 2+4) compressed only if it fulfills the following criteria:

- Is an image and not already compressed
- Photometric Interpretation (0028,0004) is MONOCHROME2 or RGB or YBR\_FULL or YBR\_FULL\_422
- Bits Allocated (0028,0100) equal to 16'D or 8'D
- Bits Stored (0028,0101) is equal to 12'D or 8'D
- High Bit (0028,0102) equal to Bits Stored (0028,0101) - 1
- Pixel Representation (0028,0103) equal to 0'D

An instance will be JPEG lossy compressed only if it fulfills the following criteria:

- Is an image and not already compressed
- Photometric Interpretation (0028,0004) is MONOCHROME2 or RGB
- Bits Allocated (0028,0100) equal to 16'D or 8'D
- Bits Stored (0028,0101) equal to 12'D or 8'D
- High Bit (0028,0102) equal to Bits Stored (0028,0101) - 1
- Pixel Representation (0028,0103) equal to 0'D
- Only lossy transfer syntaxes are supported (Implicit Little Endian is not supported) at the remote side

An instance will be JPEG 2000 lossless compressed only if it fulfills the following criteria:

- Is an image and not already compressed
- Photometric interpretation (0028,0004) neither MONOCHROME, RGB, YBR\_FULL nor YBR\_FULL\_422
- Bits Allocated (0028,0100) neither 16'D nor 8'D

An instance will be RLE lossless compressed only if it fulfills the following criteria:

- Is an image and not already compressed
- Photometric Interpretation (0028,0004) is MONOCHROME, RGB, YBR\_FULL or RGB
- Bits Allocated (0028,0100) neither 16'D nor 8'D

An instance will be JPEG 2000 lossy compressed only if it fulfills the following criteria:

- Is an image and not already compressed
- Photometric interpretation (0028,0004) is MONOCHROME or RGB
- Bits Stored (0028,0101) equal to 12'D or 8'D
- Only lossy transfer syntaxes are supported (Implicit Little Endian is not supported) at the remote side

There is no extended negotiation as an SCU.

#### 4.2.2.3.1.3 SOP specific Conformance for SOP classes

syngo® MR will not add or change private attributes by default, even in case objects are compressed or the image header is updated according to the IHE Patient Information Reconciliation Profile.

The behavior of syngo® MR when encountering status codes in a C-STORE response is summarized in Table 11:

**Table 11: DICOM Command Response Status Handling Behavior**

| Service Status | Further Meaning              | Error Code | Behavior   |
|----------------|------------------------------|------------|--|
| Error          | Any other DIMSE Error Status | 0xXXXX     | Send is continued till the end. Log message is created.                          |
| Success        | Image is successfully stored | 0000       | If configured, Storage Commitment is requested for successfully stored instances |

Table 12 below indicates the behavior if exceptions occur:

**Table 12: DICOM Command Communication Failure Behavior**

| Exception           | Behavior   |
|---------------------|--|
| Timeout             | Log message is created (Timeout configurable; default 30s) |
| Association Aborted | Send is failed. Log message is created.                    |

#### 4.2.2.3.1.4 Correction and Rearrangement

When a Study is moved to a different:

- Procedure received through a DICOM Modality Worklist, the Study Instance UID is overwritten with the Study Instance UID and Accession Number of the Procedure.
- Patient, the system generates a new Study Instance UID.

The system will not update references to the changed Study Instance UIDs, therefore it is possible that there will be broken links between Studies after such move operations.

In case of Patient Merge and Correction no UIDs are changed, therefore it is advised to delete any corrected or rearranged objects from the PACS before attempting to archive them again, to ensure that the PACS system can store them successfully.

When the Patient Position (0018,5100) attribute is corrected, the following attributes are recalculated by the system (UIDs are changed only, if the option "DICOM UID change" is enabled in Administrator Portal):

- 1) Image Position (0020,0032)
- 2) Image Orientation (0020,0037)
- 3) Patient Orientation (0020,0020)

Also the value of the Slice Location (0020,1041) attribute is emptied and a new Frame of Reference UID (0020,0052) is generated for the corrected series.

When the Patient Birth Date or the Study Date is corrected, the system recalculates the Patient Age.

A new item containing attributes that were removed or replaced by other values is added to the Original Attribute Sequence (0040,0561).

## 4.2.2.4 Association Acceptance Policy

### 4.2.2.4.1 Activity – “Receive Storage Request”

#### 4.2.2.4.1.1 Description and Sequencing of Activities

syngo® MR serves as a SCP of the Storage Service Class. The storage SCP accepts incoming C-Store Request from any configured AE Title, receives supported objects transmitted on that association and stores them in the local database.

#### 4.2.2.4.1.2 Accepted Presentation Contexts

For all supported Transfer Objects (see “Table 1: Network Services” section “SOP Classes Created by syngo® MR” and “SOP Classes Managed by syngo® MR” in the [“Conformance Statement Overview”](#).) the appropriate Transfer Syntaxes are supported.

Generally, all Presentation Contexts are accepted as long as they contain at least one suitable Transfer Syntax. All other Presentation Contexts are rejected.

There is no Extended Negotiation as an SCP

#### 4.2.2.4.1.3 SOP-specific Conformance Statement for Storage SOP classes

syngo® MR conforms to the Full Storage Class at Level 2.

In case of a successful C-STORE operation, the image has successfully been written on disk either in Explicit Little-Endian format or in the compression format received.

The Storage AE of syngo® MR returns the status “success” when the data is stored to disk and a minimal image header validation has been performed.

The following header attributes must be available and filled:

- SOP Class UID,
- Study Instance UID,
- Series Instance UID and
- SOP Instance UID.

Table 13 below list the status codes that syngo® MR can return:

**Table 13: Storage C-STORE Response Status**

| Service Status | Further Meaning                   | Error Code | Reason   |
|----------------|-----------------------------------|------------|--|
| Success        | Success                           | 0000       | Image received correctly (success notification is done after receiving, before indexing and storing) |
| Failure        | Out-of-resource                   | A700       | No resource left in the Short Term Storage   |
| Failure        | Unable to Process                 | Cxxx       | Error during instance reception  |
| Failure        | Data set does not match SOP Class | A9xx       | The data set is not conform to the SOP Class contained in the resource.                              |



Restriction: successful operation does not guarantee storage on disk and storage of header data in the database.

#### 4.2.2.4.1.4 Other SOP specific behavior

If an image is received that is already stored in the database - identified by the SOP Instance UID  
- the new image will be ignored. The existing instance is not superseded.

## 4.2.3 Storage Commitment AE Specification

### 4.2.3.1 SOP Classes

The Storage Commitment AE of syngo® MR provides standard conformance to the SOP Class listed in "Table 1: Network Services" section "Storage Commitment" in the ["Conformance Statement Overview"](#).

### 4.2.3.2 Association Policy

**Table 14: Association Policies**

|  |                       |
|--|-----------------------|
| <b>Application Context Name</b>  | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>  | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.3.2.1 Asynchronous Nature

syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 15: Asynchronous Nature as an Association Initiator**

|  |    |
|--|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|--|----|

#### 4.2.3.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to "Table 3 - Implementation Identifying Information" in the ["Conformance Statement Overview"](#).

### 4.2.3.3 Association Initiation Policy

#### 4.2.3.3.1 Activity "Send Initial Storage Commitment"

##### 4.2.3.3.1.1 Description and Sequencing of Activities

syngo® MR serves as a SCU of the Storage Commitment Service Class. After successful transfer of Imaging Objects to a configured Archive, the Storage Commitment SCU initiates an N-Action Request, if Storage Commitment is configured. This request will be sent on a different association than the storage request.

<sup>1</sup> Default, the value is configurable

The Storage Commitment Request will be sent out with a delay, in order to ensure that the remote node properly indexes received instances. The delay time is configurable with a default delay of 10 minutes.

The system may issue one N-ACTION Request for a complete set (bundle) of instances or issue one N-ACTION-Request per instance. This behavior is configurable; the default value is "bundled".

syngo® MR will accept the N-Event-Report-Request on the same association if sent immediately after the N-ACTION-Response. However, it will not wait for it. The association is closed after three seconds.

#### 4.2.3.3.1.2 Proposed Presentation Contexts

Table 16 below lists the supported presentation contexts for Storage Commitment Service Class.

**Table 16: Proposed Presentation Contexts for Storage Commitment**

| Presentation Context Table    |                      |                           |                     |      |           |
|-------------------------------|----------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax               |                      | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                          | UID                  | Name List                 | UID List            |      |           |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                               |                      | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                               |                      | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

#### 4.2.3.3.1.3 SOP specific Conformance for SOP classes

The behavior of syngo® MR when encountering status codes in an N-ACTION response is summarized in Table 17:

**Table 17: DICOM Command Response Status Handling Behavior**

| Service Status | Further Meaning                                | Error Code         | Behavior  |
|----------------|--|--------------------|---|
| Error          | Any failure that occurs                        | Any none null Code | Failure reported to user; corresponding object(s) will be marked as "Archived failed"                             |
| Success        | All Instances are available on the remote node | 0000               | Success reported to user; in case failures exist, the corresponding instances will be marked as "Archived failed" |

Table 18 below indicates the behavior if exceptions occur:

**Table 18: DICOM Command Communication Failure Behavior**

| Exception           | Behavior  |
|---------------------|---|
| Timeout             | Failure reported to user (Timeout configurable; default 30s); the request will be retried |
| Association Aborted | Failure reported to user; the request will be retried                                     |

#### 4.2.3.4 Association Acceptance Policy

##### 4.2.3.4.1 Activity “Receive Reply to Initial Storage Commitment”

###### 4.2.3.4.1.1 Description and Sequencing of Activities

syngo® MR supports the reverse role negotiation of the Storage Commitment Service Class as the SCU. It accepts incoming N-EVENT-REPORT Request, if they do not arrive on the same association as the N-ACTION-Request.

###### 4.2.3.4.1.2 Accepted Presentation Contexts

The syngo® MR DICOM application supports the presentation contexts listed in the following table for the Storage Commitment Service Class.

**Table 19 - Presentation Context Table "Update Flag Information"**

| Presentation Context Table – “Update Flag Information” |                      |                           |                     |      |           |
|--|----------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax  |                      | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name   | UID                  | Name List                 | UID List            |      |           |
| Storage Commitment Push Model                          | 1.2.840.10008.1.20.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCP  | None      |
|  |                      | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |
|  |                      | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |           |

###### 4.2.3.4.1.3 SOP-specific Conformance Statement for Storage SOP classes

The behavior of syngo® MR when encountering status codes in an N-EVENT-REPORT response is summarized in the following table:

**Table 20: DICOM Command Response Status Handling Behavior**

| Service Status | Further Meaning                   | Error Code         | Behavior                             |
|----------------|-----------------------------------|--------------------|--------------------------------------|
| Error          | Storage Commitment Reply ignored. | Any none null Code | Storage Commitment will be repeated. |
| Success        | Storage Commitment Reply noticed. | 0000               | Success reported to user.            |

## 4.2.4 Query/Retrieve AE Specification

### 4.2.4.1 SOP Classes

The Query/Retrieve AE provides Standard Conformance to the the SOP Classes listed in “Table 1: Network Services” section “Query/Retrieve” in the [“Conformance Statement Overview”](#).

### 4.2.4.2 Association Policy

**Table 21: Association Policies**

|  |                       |
|--|-----------------------|
| <b>Application Context Name</b>  | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>  | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.4.2.1 Asynchronous Nature

syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 22: Asynchronous Nature as an Association Initiator**

|  |    |
|--|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|--|----|

#### 4.2.4.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

### 4.2.4.3 Association Initiation Policy

#### 4.2.4.3.1 Activity “Querying a Remote Node” for Instances

##### 4.2.4.3.1.1 Description and Sequencing of Activities

syngo® MR serves as a SCU for the following SOP Classes

<sup>1</sup> Default, the value is configurable

- Patient Root Q/R Information Model - FIND SOP Class
- Study Root Q/R Information Model –FIND SOP Class
- Patient/Study only Q/R Information Model – FIND SOP Class.

Using the attributes specified by the user as Query Keys (in accordance with the query model) the Query SCU initiates a C-FIND Request and displays the responses to the user.

#### 4.2.4.3.1.2 Proposed Presentation Contexts

syngo® MR will propose Presentation Contexts as shown in the following table:

**Table 23: Proposed Presentation Contexts for Query**

| Presentation Context Table                                  |                             |                           |                     |      |           |
|---|-----------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax   |                             | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name  | UID                         | Name List                 | UID List            |      |           |
| Patient Root Query/Retrieve Information Model – FIND        | 1.2.840.10008.5.1.4.1.2.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | Yes       |
|   |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|   |                             | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |
| Study Root Query/ Retrieve Information Model – FIND         | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | Yes       |
|   |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|   |                             | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |
| Patient/Study Only Query/ Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.3.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | No        |
|   |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|   |                             | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

**Table 24: Extended Negotiation as an SCU**

| Name   | UID                         | Extended Negotiation  |
|--|-----------------------------|---|
| Patient Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Relational Query will be negotiated if necessary as defined in DICOM PS3.4. |
| Study Root Query/ Retrieve Information Model – FIND  | 1.2.840.10008.5.1.4.1.2.2.1 | Relational Query will be negotiated if necessary as defined in DICOM PS3.4. |

#### 4.2.4.3.1.3 SOP Specific Conformance Statement to Query SOP classes

syngo® MR checks for the following status codes in the Query SCP's C-FIND-Response:

**Table 25: DICOM Command Response Status Handling Behavior**

| Service Status | Further Meaning  | Error Code         | Behavior                           |
|----------------|--|--------------------|------------------------------------|
| Failure        | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process                  | Any none null Code | Failure reported to user           |
| Pending        | All optional keys are supported the same manner as Required Keys.  | FF00               | Pending state is indicated to user |
|                | Matching Operation continues; some of the optional keys were not supported the same way as the required keys | FF01               | Pending state is indicated to user |
| Success        | Query has been performed successfully.   | 0000               | Success reported to user           |

**Table 26: DICOM Command Communication Failure Behavior**

| exception           | Behavior  |
|---------------------|---|
| Timeout             | Failure reported to user<br>(Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user  |

syngo® MR supports the following query levels:

- Study
- Series

Matching Keys on Instance Level is not supported by the syngo® MR SCU.

The following table lists the various attributes at Study and Series levels, which can be used for hierarchical queries as well as return values for display. The display capabilities are highly configurable and “yes” indicates that it is possible to configure display of the data:

A “yes” in the UI column will indicate that the attribute may be visualized when browsing the Query results with the Browser. The Browser display is additionally influenced by the related Browser configuration

**Table 27: Attributes supported for instance Query - SCU**

| Attribute Name             | Tag         | Type | User input  | UI  |
|----------------------------|-------------|------|-------------|-----|
| <b>Study Level</b>         |             |      |             |     |
| Patient's Name             | (0010,0010) | O    | enter value | yes |
| Patient ID                 | (0010,0020) | O    | enter value | yes |
| Issuer of Patient ID       | (0010,0021) | O    | enter value | yes |
| Patient's Birth Date       | (0010,0030) | O    | enter value | yes |
| Patient's Birth Time       | (0010,0032) | O    | enter value | yes |
| Patient's Sex              | (0010,0040) | O    | enter value | yes |
| Accession Number           | (0008,0050) | O    | enter value | yes |
| Study ID                   | (0020,0010) | O    | enter value | yes |
| Study Instance UID         | (0020,000D) | U    | enter value | yes |
| Study Date                 | (0008,0020) | O    | enter value | yes |
| Study Time                 | (0008,0030) | O    | enter value | yes |
| Referring Physician's Name | (0008,0090) | O    | enter value | yes |

| Attribute Name  | Tag                       | Type | User input  | UI  |
|---|---------------------------|------|-------------|-----|
| Study Description   | (0008,1030)               | O    | enter value | yes |
| Number of Study related Instances                         | (0020,1208)               | O    | -           | yes |
| Modalities in Study                                       | (0008,0061)               | O    | enter value | yes |
| Number of Study Related Series                            | (0020,1206)               | O    | -           | yes |
| <b>Series Level</b>                                       |                           |      |             |     |
| Modality  | (0008,0060)               | O    | enter value | yes |
| Series Date   | (0008,0021)               | O    | enter value | yes |
| Series Time   | (0008,0031)               | O    | enter value | yes |
| Number of Series related Instances                        | (0020,1209)               | O    | -           | yes |
| Series Number   | (0020,0011)               | O    | enter value | yes |
| Series Description  | (0008,103E)               | O    | enter value | yes |
| Request Attributes Sequence \ Requested Procedure ID      | (0040,0275) \ (0040,1001) | O    | enter value | yes |
| Request Attributes Sequence \ Scheduled Procedure Step ID | (0040,0275) \ (0040,0009) | O    | enter value | yes |
| Performed Procedure Step Start Date                       | (0040,0244)               | O    | enter value | yes |
| Performed Procedure Step Start Time                       | (0040,0245)               | O    | enter value | yes |
| Series Instance UID                                       | (0020,000E)               | U    | -           | yes |

#### 4.2.4.3.1 Activity “Retrieve Instances from a remote node”

##### 4.2.4.3.1.1 Description and Sequencing of Activities

syngo® MR serves as a SCU for the following SOP Classes

- Patient Root Q/R Information Model - MOVE SOP Class
- Study Root Q/R Information Model – MOVE SOP Class
- Patient/Study only Q/R Information Model – MOVE SOP Class.

The C-MOVE-Request is used to retrieve the selected imaging objects. The Retrieve AE supports the query model Patient Root, Study Root, Patient/Study Root.

##### 4.2.4.3.1.2 Proposed Presentation Contexts

syngo® MR proposes Presentation Contexts as shown in the following table:

**Table 28: Proposed Presentation Contexts for Retrieve and Activity “MOVE SCU”**

| Presentation Context Table               |                             |                           |                     |      |           |
|--|-----------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax                          |                             | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                                     | UID                         | Name List                 | UID List            |      |           |
| Patient Root Query/Retrieve Model – MOVE | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | No        |
|  |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|  |                             | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |
|  |                             | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | No        |



|   |                                    |                           |                     |     |    |
|---|------------------------------------|---------------------------|---------------------|-----|----|
| Study Root Query/Retrieve Model – MOVE                | 1.2.840.10008.5.1.4.1.2.2.2        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |     |    |
|   |                                    | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |     |    |
| <b>Patient/Study Root Query/Retrieve Model – MOVE</b> | <b>1.2.840.10008.5.1.4.1.2.2.3</b> | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU | No |
|   |                                    | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |     |    |
|   |                                    | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |     |    |

#### 4.2.4.3.1.3 SOP Specific Conformance Statement for Move SCU Classes

The presentation context is negotiated at association establishment time. When the C-MOVE-Request is processed, the Move Destination attribute (receiver of images) is ignored. However, the Move Destination AE must conform to the DICOM conventions (value representation AE).

The behavior of syngo® MR when encountering status codes in a C-MOVE response is summarized in Table 29

**Table 29: DICOM Command Response Status Handling Behavior**

| Service Status | Further Meaning   | Error Code         | Behavior                          |
|----------------|---|--------------------|-----------------------------------|
| Error          | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process; Move destination unknown | Any none null Code | Failure reported to user          |
| Pending        | Move Operation continues  | FF00               | Operation continues in background |
| Success        | Move has been performed successfully.   | 0000               | Success reported to user          |

Table 30 below indicates the behavior if exceptions occur:

**Table 30: DICOM Command Communication Failure Behavior**

| Exception           | Behavior   |
|---------------------|--|
| Timeout             | Failure reported to user (Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user                                     |

#### 4.2.4.4 Association Acceptance Policy

syngo® MR provides SCU and SCP functionality.

## 4.2.5 Modality Worklist AE Specification

### 4.2.5.1 SOP Classes

The Modality Worklist AE provides Standard Conformance to the the SOP Classes listed in “Table 1: Network Services” section “Worklist Management” in the [“Conformance Statement Overview”](#).

### 4.2.5.2 Association Policy

**Table 31: Association Policies**

|  |                       |
|--|-----------------------|
| <b>Application Context Name</b>  | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>  | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.5.2.1 Asynchronous Nature

syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 32: Asynchronous Nature as an Association Initiator**

|  |    |
|--|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|--|----|

#### 4.2.5.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

### 4.2.5.3 Association Initiation Policy

#### 4.2.5.3.1 Activity “Querying a Remote Node” for Modality Worklist

##### 4.2.5.3.1.1 Description and Sequencing of Activities

syngo® MR serves as a SCU of the Modality Worklist service. It performs worklist queries by issuing a C-FIND request at regular intervals. In addition, a worklist request can be triggered manually.

<sup>1</sup> Default, the value is configurable

#### 4.2.5.3.1.2 Proposed Presentation Contexts

syngo® MR will propose Presentation Contexts as shown in the following table:

**Table 33: Proposed Presentation Contexts for Worklist**

| Presentation Context Table |                        |                           |                     |      |           |
|----------------------------|------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax            |                        | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                       | UID                    | Name List                 | UID List            |      |           |
| Modality Worklist-FIND     | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | No        |
|                            |                        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                            |                        | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

#### 4.2.5.3.1.3 SOP Specific Conformance for SOP Classes

##### Search Key Attributes of the Worklist C-FIND

syngo® MR Modality Worklist SCU supports "broad worklist queries" with all required search keys. The following tables describe the "broad query" search keys that the SCU supports. The list is configurable in 'DICOM Modality Worklist Query'.

**Table 34: Broad Query search keys**

| Attribute Name                         | Tag         | Matching Key Type | Query Value                  |
|--|-------------|-------------------|------------------------------|
| Scheduled Procedure Step               |             |                   |                              |
| Scheduled Procedure Step Sequence      | (0040,0100) | R                 |                              |
| >Modality                              | (0008,0060) | R                 | "*" or <configured Modality> |
| >Scheduled Station AE Title            | (0040,0001) | R                 | <own AET> or "*"b            |
| >Scheduled Procedure Step Start Date   | (0040,0002) | R                 | Range from UIc               |
| >Scheduled Procedure Step Description  | (0040,0007) | O                 |                              |
| >Scheduled Station Name                | (0040,0010) | O                 |                              |
| >Scheduled Procedure Step Location     | (0040,0011) | O                 |                              |
| >Scheduled Procedure Step Status       | (0040,0020) | O                 |                              |
| >Scheduled Performing Physician's Name | (0040,0006) | O                 |                              |
| >Scheduled Protocol Code Sequence      | (0040,0008) | O                 |                              |
| >>Code Value                           | (0008,0100) | O                 |                              |
| Requested Procedure Description        | (0032,1060) | O                 |                              |

<sup>b</sup> This depends on user configuration (Administration Portal->Technical Configuration->DICOM Nodes->Local DICOM Node->Worklist) if the "own AET" is provided or not.

<sup>c</sup> A time window can be configured by defining how many days to look into the past and into the future (Administration Portal-> Technical Configuration->DICOM Nodes->Local DICOM Node->Worklist).

| Attribute Name                    | Tag          | Matching Key Type | Query Value |
|-----------------------------------|--------------|-------------------|-------------|
| Requested Procedure Priority      | (0040,1003)  | O                 |             |
| Patient Transport Arrangements    | (0040,1004)  | O                 |             |
| Requested Procedure Comments      | (0040,1400)  | O                 |             |
| Requested Procedure Code Sequence | (0032,1064)  | O                 |             |
| >Code Value                       | (0008,0100)  | O                 |             |
| Requesting Physician              | (0032,1032)  | O                 |             |
| Referring Physicians Name         | (0008,0090)  | O                 |             |
| Current Patient Location          | (0038,0300)  | O                 |             |
| Pregnancy Status                  | (0010, 21C0) | O                 |             |
| Medical Alerts                    | (0010,2000)  | O                 |             |
| Allergies                         | (0010,2110)  | O                 |             |

### Return Key Attributes of the Modality Worklist C-FIND

The syngo® MR Modality Worklist SCU supports worklist queries with return key attributes of all types. The following tables describe the return keys that the SCU supports.

An “x” in the **UI** column indicates that the attribute may be visualized when browsing the Worklist results with the Browser. The Browser display is additionally influenced by the related Browser configuration.

**Table 35: Modality Worklist C-Find Return keys**

| Attribute Name                         | Tag         | Return Key Type | UI | Notes   |
|--|-------------|-----------------|----|---|
| <b>SOP Common</b>                      |             |                 |    |   |
| Specific Character Set                 | (0008,0005) | 1C              | -  |   |
| <b>Scheduled Procedure Step</b>        |             |                 |    |   |
| Scheduled Procedure Step Sequence      | (0040,0100) | 1               |    |   |
| >Modality                              | (0008,0060) | 1               | x  |   |
| >Scheduled Station AE Title            | (0040,0001) | 1               |    | “Scheduled Station AE Title” is taken as default for “Performed Station AE Title”                     |
| >Scheduled Procedure Step Start Date   | (0040,0002) | 1               | -  |   |
| >Scheduled Procedure Step Start Time   | (0040,0003) | 1               | -  |   |
| >Scheduled Procedure Step End Date     | (0040,0004) | 3               | -  |   |
| >Scheduled Procedure Step End Time     | (0040,0005) | 3               | -  |   |
| >Scheduled Performing Physician's Name | (0040,0006) | 1               | x  | “Scheduled Performing Physician's Name” is taken as default for “Performing Physician's Name”         |
| >Scheduled Procedure Step Description  | (0040,0007) | 1C              | x  | “Scheduled Procedure Step Description” is taken as default for “Performed Procedure Step Description” |
| >Scheduled Protocol Code Sequence **   | (0040,0008) | 1C              | -  | Uses universal sequence match   |

| Attribute Name                            | Tag         | Return Key Type | UI | Notes  |
|---|-------------|-----------------|----|--|
|   |             |                 |    | "Scheduled Protocol Code Sequence" is taken as default for "Performed Protocol Code Sequence"                          |
| >>Code Value                              | (0008,0100) | 1C              | -  |  |
| >>Coding Scheme Designator                | (0008,0102) | 1C              | -  |  |
| >>Coding Scheme Version                   | (0008,0103) | 3               | -  |  |
| >>Code Meaning                            | (0008,0104) | 3               | -  |  |
| >Scheduled Procedure Step ID              | (0040,0009) | 1               | x  | "Scheduled Procedure Step ID" is taken as default for "Performed Procedure Step ID"                                    |
| >Scheduled Station Name                   | (0040,0010) | 2               | x  |  |
| >Scheduled Procedure Step Location        | (0040,0011) | 2               | -  | "Scheduled Procedure Step Location" is taken as default for "Performed Location"                                       |
| >Scheduled Procedure Step Status          | (0040,0020) | 3               | -  |  |
| >Comments on the Scheduled Procedure Step | (0040,0400) | 3               | -  |  |
| <b>Requested Procedure</b>                |             |                 |    |  |
| Study Date                                | (0008,0020) | 3               | x  |  |
| Study Time                                | (0008,0030) | 3               | x  |  |
| Referenced Study Sequence **              | (0008,1110) | 2               | -  | Uses universal sequence match  |
| >Referenced SOP Class UID                 | (0008,1150) | 1C              | -  |  |
| >Referenced SOP Instance UID              | (0008,1155) | 1C              | -  |  |
| Study Instance UID                        | (0020,000D) | 1               | -  |  |
| Requested Procedure Description           | (0032,1060) | 1C              | x  |  |
| Requested Procedure Code Sequence **      | (0032,1064) | 1C              | -  | Uses universal sequence match<br>"Requested Procedure Code Sequence" is taken as default for "Procedure Code Sequence" |
| >Code Value                               | (0008,0100) | 1C              | -  |  |
| >Coding Scheme Designator                 | (0008,0102) | 1C              | -  |  |
| >Coding Scheme Version                    | (0008,0103) | 3               | -  |  |
| >Code Meaning                             | (0008,0104) | 3               | -  |  |
| Requested Procedure ID                    | (0040,1001) | 1               | x  | "Requested Procedure ID" is taken as default for "Study ID"  |
| Reason for the Requested Procedure        | (0040,1002) | 3               | -  |  |
| Requested Procedure Priority              | (0040,1003) | 2               | x  |  |
| Patient Transport Arrangements            | (0040,1004) | 2               | -  |  |
| Confidentiality Code                      | (0040,1008) | 3               | -  |  |
| Reporting Priority                        | (0040,1009) | 3               | x  |  |
| Names of intended Recipients of Results   | (0040,1010) | 3               | -  |  |
| Requested Procedure Comments              | (0040,1400) | 3               | -  |  |
| <b>Imaging Service Request</b>            |             |                 |    |  |
| Accession Number                          | (0008,0050) | 2               | x  |  |
| Referring Physician's Name                | (0008,0090) | 2               | x  |  |
| Requesting Physician                      | (0032,1032) | 2               | x  |  |
| Requesting Service                        | (0032,1033) | 3               | -  |  |
| Issuing Date of Imaging Service Request   | (0040,2004) | 3               | -  |  |

| Attribute Name                                  | Tag         | Return Key Type | UI | Notes  |
|---|-------------|-----------------|----|--|
| Issuing Time of Imaging Service Request         | (0040,2005) | 3               | -  |  |
| Placer Order Number / Imaging Service Request * | (0040,2016) | 3               | -  | Old tag (0040,2006) is retired and not used. |
| Filler Order Number / Imaging Service Request * | (0040,2017) | 3               | -  | Old tag (0040,2007) is retired and not used. |
| Order entered by ...                            | (0040,2008) | 3               | -  |  |
| Order Enterer's location                        | (0040,2009) | 3               | -  |  |
| Order Callback Phone Number                     | (0040,2010) | 3               | -  |  |
| Imaging Service Request Comments                | (0040,2400) | 3               | -  |  |
| <b>Visit Identification</b>                     |             |                 |    |  |
| Admission ID                                    | (0038,0010) | 2               | x  |  |
| Issuer of Admission ID                          | (0038,0011) | 3               | -  |  |
| <b>Visit Status</b>                             |             |                 |    |  |
| Current Patient Location                        | (0038,0300) | 2               | x  |  |
| <b>Visit Admission</b>                          |             |                 |    |  |
| Admitting Diagnosis Description                 | (0008,1080) | 3               | x  |  |
| Admitting Date                                  | (0038,0020) | 3               | -  |  |
| <b>Patient Identification</b>                   |             |                 |    |  |
| Patient's Name                                  | (0010,0010) | 1               | x  |  |
| Patient ID                                      | (0010,0020) | 1               | x  |  |
| Issuer of Patient ID                            | (0010,0021) | 3               | x  |  |
| Other Patient IDs                               | (0010,1000) | 3               | x  |  |
| Other Patient Names                             | (0010,1001) | 3               | x  |  |
| Patient's Birth Name                            | (0010,1005) | 3               | -  |  |
| <b>Patient Demographic</b>                      |             |                 |    |  |
| Patient's Birth Date                            | (0010,0030) | 2               | x  |  |
| Patient's Birth Time                            | (0010,0032) | 3               | x  |  |
| Patient's Sex                                   | (0010,0040) | 2               | x  |  |
| Patient's Insurance Plan Code Sequence **       | (0010,0050) | 3               | -  | Uses universal sequence match                |
| >Code Value                                     | (0008,0100) | 1C              | -  |  |
| >Coding Scheme Designator                       | (0008,0102) | 1C              | -  |  |
| >Coding Scheme Version                          | (0008,0103) | 3               | -  |  |
| >Code Meaning                                   | (0008,0104) | 3               | -  |  |
| Patient's Age                                   | (0010,1010) | 3               | -  |  |
| Patient's Size                                  | (0010,1020) | 3               | x  |  |
| Patient's Weight                                | (0010,1030) | 2               | x  |  |
| Patient's Address                               | (0010,1040) | 3               | x  |  |
| Military Rank                                   | (0010,1080) | 3               | x  |  |
| Branch of Service                               | (0010,1081) | 3               | -  |  |
| Ethnic Group                                    | (0010,2160) | 3               | x  |  |
| Patient Comments                                | (0010,4000) | 3               | x  |  |
| <b>Patient Medical</b>                          |             |                 |    |  |
| Medical Alerts                                  | (0010,2000) | 2               | x  |  |
| Allergies                                       | (0010,2110) | 2               | x  |  |
| Pregnancy Status                                | (0010,21C0) | 2               | x  |  |
| Smoking Status                                  | (0010,21A0) | 3               |    |  |
| Last Menstrual Date                             | (0010,21D0) | 3               |    |  |
| Additional Patient History                      | (0010,21B0) | 3               |    |  |
| Special Needs                                   | (0038,0050) | 2               | x  |  |

syngo® MR only supports a one-to-one relationship between Requested Procedure and Scheduled Procedure Steps. If multiple Schedule Procedure Steps are scheduled for a procedure, they will result in one Performed Procedure Step.

The behavior of syngo® MR when encountering status codes in a C-FIND response is summarized in Table 36:

**Table 36: DICOM Command Response Status Handling Behavior**

| Service Status | Further Meaning  | Error Code         | Behavior                           |
|----------------|--|--------------------|------------------------------------|
| Error          | e.g. Out of Resources; Cancellation; Identifier does not match SOP Class; Unable to process                  | Any none null Code | Failure reported to user           |
| Pending        | All optional keys are supported the same manner as Required Keys.  | FF00               | Pending state is indicated to user |
|                | Matching Operation continues; some of the optional keys were not supported the same way as the required keys | FF01               | Pending state is indicated to user |
| Success        | Query has been performed successfully.   | 0000               | Success reported to user           |

Table 37 below indicates the behavior if exceptions occur:

**Table 37: DICOM Command Communication Failure Behavior**

| exception           | Behavior  |
|---------------------|---|
| Timeout             | Failure reported to user<br>(Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user  |

#### 4.2.5.4 Association Acceptance Policy

syngo® MR does not provide the functionality of a SCP of the Modality Worklist – Find SOP Class.

## 4.2.6 Modality Performed Procedure Step AE Specification

### 4.2.6.1 SOP Classes

The Modality Performed Procedure Step AE provides Standard Conformance to the the SOP Classes listed in “Table 1: Network Services” section “Worklist Management” in the [“Conformance Statement Overview”](#).

### 4.2.6.2 Association Policy

**Table 38: Association Policies**

|  |                       |
|--|-----------------------|
| <b>Application Context Name</b>  | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>  | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.6.2.1 Asynchronous Nature

syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 39: Asynchronous Nature as an Association Initiator**

|  |    |
|--|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|--|----|

#### 4.2.6.2.2 Implementation Identifying Information

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

### 4.2.6.3 Association Initiation Policy

#### 4.2.6.3.1 Activity “Create Modality Performed Procedure Step”

##### 4.2.6.3.1.1 Description and Sequencing of Activities

syngo® MR serves as a SCU of the Modality Performed Procedure Step SOP Class. It sends N-CREATE request to inform the Information System that a Procedure Step has been started.

<sup>1</sup> Default, the value is configurable



#### 4.2.6.3.1.2 Accepted Presentation Contexts

syngo® MR proposes Presentation Contexts as shown in the following table:

**Table 40: Acceptable Presentation Contexts Activity “Create MPPS”**

| Presentation Context Table        |                         |                           |                     |      |           |
|-----------------------------------|-------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax                   |                         | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                              | UID                     | Name List                 | UID List            |      |           |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                                   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                                   |                         | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

#### 4.2.6.3.1.3 SOP specific Conformance for MPPS SOP class

The behavior of syngo® MR when encountering status codes in an N-CREATE-RSP response is summarized in Table 41:

**Table 41: MPPS N-CREATE Response Status Handling Behavior**

| Service Status | Further Meaning                               | Error Code         | Behavior             |
|----------------|---|--------------------|----------------------|
| Error          | MPPS creation request could not be processed. | Any none null Code | MPPS is not created. |
| Success        | MPPS creation request processed successfully. | 0000               | MPPS is created.     |

#### 4.2.6.3.2 Activity “Update Modality Performed Procedure Step”

##### 4.2.6.3.2.1 Description and Sequencing of Activities

When the procedure step has been finished, syngo® MR sends N-SET request to inform the Information System about the finalization of the procedure step (COMPLETED or DISCONTINUED).

##### 4.2.6.3.2.2 Proposed Presentation Contexts

syngo® MR proposes Presentation Contexts as shown in the following table:

**Table 42: Acceptable Presentation Contexts Activity “Update MPPS”**

| Presentation Context Table        |                         |                           |                     |      |           |
|-----------------------------------|-------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax                   |                         | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                              | UID                     | Name List                 | UID List            |      |           |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                                   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                                   |                         | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

#### 4.2.6.3.2.3 SOP specific Conformance for MPPS SOP class

The behavior of syngo® MR when encountering status codes in an N-SET-RSP response is summarized in Table 43:

**Table 43: MPPS N-SET Response Status Handling Behavior**

| Service Status | Further Meaning                                   | Error Code         | Behavior             |
|----------------|---|--------------------|----------------------|
| Error          | MPPS update request could not be processed.       | Any none null Code | MPPS is not updated. |
| Success        | MPPS update request could processed successfully. | 0000               | MPPS is updated.     |

#### 4.2.6.4 Association Acceptance Policy

syngo® MR does not provide the functionality of a SCP of the Modality Performed Procedure Step SOP Class.

## 4.2.7 Print AE Specification

### 4.2.7.1 SOP Classes

The Print AE provides Standard Conformance to the the SOP Classes listed in “Table 1: Network Services”section “Print Management” in the [“Conformance Statement Overview”](#).

### 4.2.7.2 Association Policy

**Table 44: Association Policies**

|  |                       |
|--|-----------------------|
| <b>Application Context Name</b>  | 1.2.840.10008.3.1.1.1 |
| <b>PDU size</b>  | 32 kB <sup>1</sup>    |
| <b>Maximum number of simultaneous associations as an association acceptor</b>  | 12 <sup>1</sup>       |
| <b>Maximum number of simultaneous associations as an association initiator</b> | unlimited             |

syngo® MR contains a limitation of 512 kB for the maximum PDU size. By default, the maximum PDU size is set to 32kB.

The maximum number of simultaneous receiving associations (SCP) is configurable at run time, based on the system resources available. By default, the maximum number of associations is set to 12.

There is no inherent limit to the number of outgoing associations (SCU), other than limits imposed by the computer operating system.

#### 4.2.7.2.1 Asynchronous Nature

syngo® MR supports asynchronous communication (multiple outstanding transactions over a single association). On the SCU side the Window size proposed is infinite. On the SCP Side any non-infinite maximum size will be accepted.

**Table 45: Asynchronous Nature as an Association Initiator**

|  |    |
|--|----|
| <b>Maximum number of outstanding asynchronous transactions</b> | 10 |
|--|----|

### 4.2.7.3 Association Initiation Policy

#### 4.2.7.3.1 Activity Print Film

##### 4.2.7.3.1.1 Description and Sequencing of Activities

Whenever a film-sheet is prepared by the user, it is forwarded to the Printer Job queue. As soon as the associated Printer device is available the job is activated, and an association is established.

After the film sheet is internally processed, converted to a Standard/1,1 layout and the page image is sent to the printer, the status is controlled by awaiting any N-EVENT-REPORT message throughout the transfer until the last image or film-sheet is sent.

If the response from the remote application contains a status other than Success or Warning the printing is stopped and the job status is set to Aborted.

<sup>1</sup> Default, the value is configurable

#### 4.2.7.3.1.2 Proposed Presentation Context

syngo® MR proposes Presentation Contexts as shown in the following table:

**Table 46: Presentation Contexts for the Activity “Print Film”**

| Presentation Context Table                      |                         |                           |                   |      |           |
|---|-------------------------|---------------------------|-------------------|------|-----------|
| Abstract Syntax                                 |                         | Transfer Syntax           |                   | Role | Ext. Neg. |
| Name  | UID                     | Name List                 | UID List          |      |           |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9   | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Color Print Management Meta SOP Class     | 1.2.840.10008.5.1.1.18  | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Film Sesssion SOP Class                   | 1.2.840.10008.5.1.1.1   | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Film Box SOP Class                        | 1.2.840.10008.5.1.1.2   | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Grayscale Image Box SOP Class             | 1.2.840.10008.5.1.1.4   | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Basic Color Image Box SOP SOP Class             | 1.2.840.10008.5.1.1.4.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Printer SOP Class                               | 1.2.840.10008.5.1.1.16  | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Print Job SOP Class                             | 1.2.840.10008.5.1.1.14  | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |
| Presentation LUT SOP Class                      | 1.2.840.10008.5.1.1.23  | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |

#### 4.2.7.3.1.3 SOP Specific Conformance

syngo® MR Print SCU conforms to the DICOM Basic Grayscale Print Management Meta SOP Class and the Basic Color Print Management Meta SOP Class.

The application uses a configuration platform to define the properties of the connected DICOM SCP, e.g.:

- supported film sizes of the connected DICOM SCP
- supported film formats of the DICOM SCP

The printing is only suspended in the case of a failure return status of the SCP.

#### Basic Film Session SOP Class

The Basic Film Session information object definition describes all the user-defined parameters, which are common for all the films of a film session. The Basic Film Session refers to one or more Basic Film Boxes that are printed on one hardcopy printer.

syngo® MR Print Management SCU supports the following DIMSE Service elements for the Basic Film Session SOP Class as SCU:

- N-CREATE
- N-DELETE

The Basic Film Session SOP Class N-CREATE-RQ (SCU) uses the attributes listed in the table below:

**Table 47: Attributes for the N-CREATE-RQ of the Basic Film Session**

| Attribute Name   | Tag         | Usage SCU | Supported Values |
|------------------|-------------|-----------|------------------|
| Number of Copies | (2000,0010) | U         | 1                |
| Medium Type      | (2000,0030) | U         | BLUE FILM        |
|                  |             |           | CLEAR FILM       |
|                  |             |           | PAPER            |
|                  |             |           | MAMMO BLUE FILM  |
|                  |             |           | MAMMO CLEAR FILM |

The number of Copies sent to the DICOM Printer is always 1, a number higher than 1 is not supported in this version.

The Affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g. N-DELETE-RQ) on the Basic Film Session.

The Basic Film Session SOP class interprets the status codes (from N-CREATE-RSP messages) listed in the table below:

**Table 48: N-CREATE-RSP Status Handling Behavior for the Basic Film Session**

| Service Status | Further Meaning                   | Error Codes | Behavior                               |
|----------------|-----------------------------------|-------------|--|
| Warning        | Memory Allocation not supported   | B600        | Print job continues, warning is logged |
| Success        | Film session successfully created | 0000        | Print job continues                    |

The N-DELETE-RQ on the Basic Film Session SOP Class is used to remove the complete Basic Film Session SOP Instance hierarchy.

### Basic Film Box SOP Class

The Basic Film Box information object definition describes all user-defined parameters of one film of the film session including presentation parameters, which are common for all images on a given film sheet.

The Basic Film Box refers to one or more Image Boxes.

syngo® MR DICOM Print Management SCU supports the following DIMSE Service elements for the Basic Film Box SOP Class as SCU:

- N-CREATE
- N-ACTION
- N-DELETE

The Basic Film Box SOP Class N-CREATE-RQ message uses the attributes listed below. The actual values for each attribute depend on DICOM printer configuration within the syngo® MR DICOM Print Management SCU:

**Table 49: Attributes for the N-CREATE-RQ of the Basic Film Session**

| Attribute Name                                 | Tag         | Usage SCU | Supported Values   |
|--|-------------|-----------|--|
| Image Display Format                           | (2010,0010) | M         | STANDARD\1,1   |
| Referenced Film Session Sequence               | (2010,0500) | M         |  |
| > Referenced SOP Class UID                     | (0008,1150) | M         | 1.2.840.10008.5.1.1.1  |
| > Referenced SOP Instance UID                  | (0008,1155) | M         |  |
| Film Orientation                               | (2010,0040) | M         | PORTRAIT, LANDSCAPE  |
| Film Size ID                                   | (2010,0050) | M         | 8INX10IN, 10INX12IN, 10INX14IN, 11INX14IN,, 14INX14IN, 14INX17IN, 24CMX24CM, 24CMX30CM |
| Magnification Type                             | (2010,0060) | M         | BILINEAR, CUBIC, NONE, REPLICATE   |
| Border Density                                 | (2010,0100) | U         | BLACK, WHITE   |
| Max Density                                    | (2010,0130) | U         | 0 < Value  |
| Min Density                                    | (2010,0120) | U         | 0 < Value < 50   |
| <b>Required if Presentation LUT is present</b> |             |           |  |
| Reflective Ambient Light                       | (2010,0160) | U         | 0 < Value  |
| Illumination                                   | (2010,015E) | U         | 0 < Value  |
| Referenced Presentation LUT Sequence           | (2050,0500) | U         |  |

For Page Mode printing, the Image Display format used is Standard\1,1.

The N-CREATE-RSP message from the Print SCP includes the Referenced Image Box Sequence with SOP Class/Instance UID pairs which will be kept internally to be further used for the subsequent Basic Image Box SOP Class N-SET-RQ messages.

When all Image Boxes (including parameters) for the film-sheet have been set, the syngo® MR print manager will issue an N-ACTION-RQ message with the SOP Instance UID of the Basic Film Box and the Action Type ID of 1.

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g. N-DELETE-RQ) on the Basic Film Box.

The Basic Film Box SOP class interprets the status codes listed in the tables below:

**Table 50: N-CREATE-RSP Status Handling Behavior for Basic Film Box**

| Service Status | Meaning   | Error Codes | Behavior   |
|----------------|---|-------------|--|
| Failure        | There is an existing Film Box that has not been printed and N-ACTION at the Film Session level is not supported. A new Film Box will not be created when a previous Film Box has not been printed | C616        | Print job is marked as failed and the reason is logged |
| Warning        | Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead   | B605        | Print job continues and warning is logged              |
| Success        | Film Box successfully created   | 0000        | Print job continues                                    |

**Table 51: N-ACTION-RSP Status Handling Behavior for Basic Film Box**

| Service Status | Meaning   | Error Codes | Behavior   |
|----------------|---|-------------|--|
| Failure        | Unable to create print job, print queue is full   | C602        | Print job is marked as failed and the reason is logged |
|                | Image size is larger than images box size   | C603        | Print job is marked as failed and the reason is logged |
|                | Combined Print Image size is larger than the Image Box size   | C613        | Print job is marked as failed and the reason is logged |
| Warning        | Film box does not contain image box (empty page)  | B603        | Print job continues and warning is logged              |
|                | Image size is larger than image box size, the image has been demagnified  | B604        | Print job continues and warning is logged              |
|                | Image size is larger than the Image Box size. The Image has been cropped to fit.  | B609        | Print job continues and warning is logged              |
|                | Image size or Combined Print Image size is larger than the Image Box size. Image or Combined Print Image has been decimated to fit. | B60A        | Print job continues and warning is logged              |
| Success        | Film accepted for printing  | 0000        | Print job continues                                    |

## Basic Grayscale Image Box SOP Class

The Basic Grayscale Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Image Box information describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Grayscale Image Box SOP Class uses only the N-SET-RQ with the attributes listed in the table below:

**Table 52: Attributes for N-SET-RQ of Basic Grayscale Image Box**

| Attribute Name                 | Tag         | Usage SCU | Supported Values |
|--------------------------------|-------------|-----------|------------------|
| Image Position                 | (2020,0010) | M         | 1                |
| Basic Grayscale Image Sequence | (2020,0110) | M         |                  |
| > Samples per Pixel            | (0028,0002) | M         | 1                |

| Attribute Name               | Tag         | Usage SCU | Supported Values |
|------------------------------|-------------|-----------|------------------|
| > Photometric Interpretation | (0028,0004) | M         | MONOCHROME2      |
| > Rows                       | (0028,0010) | M         |                  |
| > Columns                    | (0028,0011) | M         |                  |
| > Pixel Aspect Ratio         | (0028,0034) | M         |                  |
| > Bits Allocated             | (0028,0100) | M         | 8,16             |
| > Bits Stored                | (0028,0101) | M         | 8,12             |
| > High Bit                   | (0028,0102) | M         | 7,11             |
| > Pixel Representation       | (0028,0103) | M         | 0                |
| > Pixel Data                 | (7FE0,0010) | M         |                  |

The Basic Grayscale Image Box SOP class interprets the status codes as listed below:

**Table 53: N-SET-RSP Status Handling Behavior for the Basic Grayscale Image Box SOP Class**

| Service Status | Further Meaning   | Error Codes | Behavior   |
|----------------|---|-------------|--|
| Failure        | Image contains more pixel than printer can print in Image Box   | C603        | Print job is marked as failed and the reason is logged |
|                | Insufficient memory in printer to store the image   | C605        | Print job is marked as failed and the reason is logged |
|                | Combined Print Image size is larger than the Image Box size   | C613        | Print job is marked as failed and the reason is logged |
| Warning        | Image size is larger than image box size, the image has been demagnified.   | B604        | Print job continues and the reason is logged           |
|                | Requested MinDensity or MaxDensity outside of Printer's operating range   | B605        | Print job continues and the reason is logged           |
|                | Image size is larger than the Image Box size. The Image has been cropped to fit.  | B609        | Print job continues and warning is logged              |
|                | Image size or Combined Print Image size is larger than the Image Box size. Image or Combined Print Image has been decimated to fit. | B60A        | Print job continues and warning is logged              |
| Success        | Image successfully stored in Image Box  | 0000        | Print job continues                                    |

## Basic Color Image Box SOP Class

The Basic Color Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Image Box information describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Color Image Box SOP Class uses only the N-SET-RQ with the attributes listed below:

**Table 54: Attributes for N-SET-RQ of Basic Color Image Box**



| Attribute Name               | Tag         | Usage SCU | Supported Values |
|------------------------------|-------------|-----------|------------------|
| Image Position               | (2020,0010) | M         | 1                |
| BASIC Color Image Sequence   | (2020,0111) | M         |                  |
| > Samples per Pixel          | (0028,0002) | M         | 3                |
| > Photometric Interpretation | (0028,0004) | M         | RGB              |
| > Planar Configuration       | (0028,0006) | M         | 0                |
| > Rows                       | (0028,0010) | M         |                  |
| > Columns                    | (0028,0011) | M         |                  |
| > Pixel Aspect Ratio         | (0028,0034) | M         |                  |
| > Bits Allocated             | (0028,0100) | M         | 8                |
| > Bits Stored                | (0028,0101) | M         | 8                |
| > High Bit                   | (0028,0102) | M         | 7                |
| > Pixel Representation       | (0028,0103) | M         | 0                |
| > Pixel Data                 | (7FE0,0010) | M         |                  |

The Color Image Box SOP class interprets the status codes listed below:

**Table 55: N-SET-RSP Status Handling Behavior for the Color Grayscale Image Box**

| Service Status | Meaning   | Error Codes | Behavior   |
|----------------|---|-------------|--|
| Failure        | Image contains more pixel than printer can print in Image Box   | C603        | Print job is marked as failed and the reason is logged |
|                | Insufficient memory in printer to store the image   | C605        | Print job is marked as failed and the reason is logged |
|                | Combined Print Image size is larger than the Image Box size   | C613        | Print job is marked as failed and the reason is logged |
| Warning        | Image size larger than image box size   | B604        | Print job continues and the reason is logged           |
|                | Image size is larger than the Image Box size. The Image has been cropped to fit.  | B609        | Print job continues and warning is logged              |
|                | Image size or Combined Print Image size is larger than the Image Box size. Image or Combined Print Image has been decimated to fit. | B60A        | Print job continues and warning is logged              |
| Success        | Image successfully stored in Image Box  | 0000        | Print job continues                                    |

## Presentation LUT SOP Class

The objective of the Presentation LUT is to realize image hardcopy printing tailored for specific modalities, applications and user preferences.

The output of the Presentation LUT is Presentation Values (P-Values). P-Values are approximately related to human perceptual response. They are intended to facilitate common input for hardcopy. P-Values are intended to be independent of the specific class or characteristics of the hardcopy device.

The syngo® MR DICOM Print Management SCU supports the following DIMSE Service elements for the Basic Film Session SOP Class as SCU:

- N-CREATE
- N-DELETE

The Presentation LUT SOP Class uses only the N-CREATE-RQ with the attributes listed below:

**Table 56: Attributes for N-CREATE-RQ of Presentation LUT SOP Class**

| Attribute Name         | Tag         | Usage SCU | Supported Values |
|------------------------|-------------|-----------|------------------|
| Presentation LUT Shape | (2050,0020) | U         | IDENTITY         |

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and is used for later requests on the Basic Film Box (N-CREATE-RQ) and on the Presentation LUT (N-DELETE-RQ).

The Presentation LUT SOP class interprets the status codes listed below:

**Table 57: N-CREATE-RSP Status Handling Behavior for the Presentation LUT SOP Class**

| Service Status | Further Meaning  | Error Codes | Behavior                                     |
|----------------|--|-------------|--|
| Warning        | Requested MinDensity or MaxDensity outside of HCD's operating range. HCD will use its respective minimum or maximum density value instead. | B605        | Print job continues and the reason is logged |
| Success        | Presentation LUT successfully created  | 0000        | Print job continues                          |

## Printer SOP Class

The Printer SOP Class provides the possibility to monitor the status of the hardcopy printer in a synchronous and in an asynchronous way.

When used synchronously the syngo® MR Print SCU uses the N-GET-RQ to request information about the printer status. It uses the attributes listed in the table below.

**Table 58: Attributes for N-GET-RQ of the Printer SOP Class**

| Attribute Name      | Tag         | Usage SCP | Supported Values                              |
|---------------------|-------------|-----------|---|
| Printer Status      | (2110,0010) | M         | NORMAL, FAILURE, WARNING                      |
| Printer Status Info | (2110,0020) | M         | See table in chapter 9.4 for possible values. |

The command communication failure behavior listed below applies to all SOP classes used for the "Print Film" activity:

**Table 59: DICOM Command Communication Failure Behavior**

| Exception           | Behavior  |
|---------------------|---|
| Timeout             | Failure reported to user<br>(Timeout configurable; default 30s) |
| Association Aborted | Failure reported to user  |

#### 4.2.7.4 Association Acceptance Policy

##### 4.2.7.4.1 Activity Print Film

###### 4.2.7.4.1.1 Description and Sequencing of Activities

syngo® MR supports the reverse role negotiation of the Printer SOP Class. Receiving the N-EVENT-REPORT-RQ from a printer syngo® MR is asynchronously informed about changes of the printer status.

###### 4.2.7.4.1.2 Accepted Presentation Context

syngo® MR accepts Presentation Contexts as shown in the following table:

**Table 60: Presentation Contexts for the Activity “Print Film”**

| Presentation Context Table |                         |                           |                   |      |           |
|----------------------------|-------------------------|---------------------------|-------------------|------|-----------|
| Abstract Syntax            |                         | Transfer Syntax           |                   | Role | Ext. Neg. |
| Name                       | UID                     | Name List                 | UID List          |      |           |
| Printer SOP Class          | 1.2.840.10008.5.1.1.1.6 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |

###### 4.2.7.4.1.3 SOP Specific Conformance

The arguments of the N-EVENT-REPORT-RQ are defined in the table below:

**Table 61: Attributes for the N-EVENT-REPORT-RQ of the Printer SOP Class**

| Event-type Name | Event | Attributes          | Tag         | Usage SCU |
|-----------------|-------|---------------------|-------------|-----------|
| Normal          | 1     | -                   | -           | -         |
| Warning         | 2     | Printer Status Info | (2110,0020) | U         |
| Failure         | 3     | Printer Status Info | (2110,0020) | U         |

##### 4.2.7.4.2 Activity Print Management

###### 4.2.7.4.2.1 Description and Sequencing of Activities

syngo® MR supports the reverse role negotiation of the Print Job SOP Class. Receiving the N-EVENT-REPORT-RQ from a printer syngo® MR is asynchronously informed about the status of a print job for monitoring its progress.

#### 4.2.7.4.2.2 Accepted Presentation Context

syngo® MR accepts Presentation Contexts as shown in the following table:

**Table 62: Presentation Contexts for the Activity “Print Management”**

| Presentation Context Table |                         |                           |                   |      |           |
|----------------------------|-------------------------|---------------------------|-------------------|------|-----------|
| Abstract Syntax            |                         | Transfer Syntax           |                   | Role | Ext. Neg. |
| Name                       | UID                     | Name List                 | UID List          |      |           |
| Print Job SOP Class        | 1.2.840.10008.5.1.1.1.4 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None      |

#### 4.2.7.4.2.3 SOP Specific Conformance

Attributes that can be handled by the Print AE of the product are listed in the table below.

**Table 63: Attributes for the N-EVENT-REPORT-RQ of the Print Job SOP Class**

| Event-type Name | Event | Attributes            | Tag         | Usage SCU  |
|-----------------|-------|-----------------------|-------------|--|
| Pending         | 1     | Execution Status Info | (2100,0030) | U  |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U  |
|                 |       | Printer Name          | (2110,0030) | U  |
| Printing        | 2     | Execution Status Info | (2100,0030) | U  |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U  |
|                 |       | Printer Name          | (2110,0030) | U  |
| Done            | 3     | Execution Status Info | (2100,0030) | U  |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U  |
|                 |       | Printer Name          | (2110,0030) | U  |
| Failure         | 4     | Execution Status Info | (2100,0030) | U  |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U  |
|                 |       | Printer Name          | (2110,0030) | U  |

## **4.2.8 Print SCU Specification**

### **4.2.8.1 SOP Classes**

For SOP Classes supported, please refer to “Table 1: Network Services” section “Print Management” in the [“Conformance Statement Overview”](#).

### **4.2.8.2 Association Policies**

#### **4.2.8.2.1 General**

Whenever a film-sheet is completely set up and printed by command or automated rule, the job is prepared for processing. As soon as the queue is ready to process the job, it is activated and worked according to the processing data. The Print application will initiate an association to the print destination and process the printing.

The default PDU size used will be used.

#### **4.2.8.2.2 Number of Associations**

The syngo® MR DICOM application initiates one association at a time for each different print device configured.

#### **4.2.8.2.3 Asynchronous Nature**

The syngo® MR DICOM print application does not support asynchronous communication (multiple outstanding transactions over a single association).

#### **4.2.8.2.4 Implementation Identifying Information**

For Implementation Identifying Information please refer to “Table 3 - Implementation Identifying Information” in the [“Conformance Statement Overview”](#).

### **4.2.8.3 Association Initiation Policy**

Triggered by the Print job queue the Print Management SCU establishes an association by using the DICOM association services. An N-GET request determines the printer status prior to printing. If the printer status is “normal”, the print job is started.

After the last film is printed from queue, the Print application will leave open the association for another 60 seconds. If a new film job is ready for printing within this time-limit, the job will be immediately processed over the still open association. If there is no new job, the association is closed.

During the “idle-time” (no open association to printer) the Print application will issue a cyclic camera status request (using N-GET of the Printer SOP Class) every 5 minutes.

#### **4.2.8.3.1 Activity - Print Film**

##### **4.2.8.3.1.1 Description and Sequencing of Activity**

The film sheet is internally processed, converted to a Standard/1-1 page and then the page image is sent. Status is controlled by awaiting any N-EVENT message through the transfer until the last image or film-sheet is sent.

If the response from the remote application contains a status other than Success or Warning the association is aborted.

#### 4.2.8.3.1.2 Proposed Presentation Context

The syngo® MR DICOM application will propose Presentation Contexts as shown in the following table:

**Table 64 - Presentation Context Table "Print Film"**

| Presentation Context Table                      |                        |  |  |      |           |
|---|------------------------|--|--|------|-----------|
| Abstract Syntax                                 |                        | Transfer Syntax  |  | Role | Ext. Neg. |
| Name  | UID                    | Name List  | UID List                                 |      |           |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9  | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None      |
| Basic Color Print Management Meta SOP Class     | 1.2.840.10008.5.1.1.18 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None      |
| Print Job SOP Class                             | 1.2.840.10008.5.1.1.14 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None      |
| Presentation LUT SOP Class                      | 1.2.840.10008.5.1.1.23 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None      |

#### 4.2.8.3.1.3 SOP Specific Conformance

The syngo® MR DICOM print management SCU conforms to the DICOM Basic Grayscale Print Management Meta SOP Class and to Basic Color Print Management Meta SOP Class.

The application uses a setting platform to define the properties of the connected DICOM SCP, e.g.:

- supported film sizes of the connected DICOM SCP
- supported film formats of the DICOM SCP
- LUT type to be attached.

The printing is suspended in the case of a failure return status of the SCP or when the user cancels the job.

#### Basic Film Session SOP Class

The Basic Film Session information object definition describes all the user-defined parameters, which are common for all the films of a film session. The Basic Film Session refers to one or more Basic Film Boxes that are printed on one hardcopy printer.

The syngo® MR DICOM print management SCU supports the following DIMSE Service elements for the Basic Film Session SOP Class as SCU:

- N-CREATE
- N-DELETE

The Basic Film Session SOP Class N-CREATE-RQ (SCU) uses the following attributes:

**Table 65 - Basic Film Session N-CREATE attributes**

| Attribute Name   | Tag         | Usage SCU | Supported Values                 |
|------------------|-------------|-----------|----------------------------------|
| Number of Copies | (2000,0010) | U         | "1"                              |
| Medium Type      | (2000,0030) | U         | BLUE FILM<br>CLEAR FILM<br>PAPER |

U = User Option

The number of copies sent to the DICOM Printer is always 1, the job is sent n times for n copies.

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (e.g. N-DELETE-RQ) on the Basic Film Session – see below:

| Attribute Name             | Tag                         | Source of Information   |
|----------------------------|-----------------------------|---|
| Requested SOP Instance UID | (0000,1000)<br>→(0000,1001) | Affected SOP Instance UID of N-CREATE-RSP on Basic Film Session |

The N-DELETE-RQ on the Basic Film Session SOP Class is used to remove the complete Basic Film Session SOP Instance hierarchy.

The Basic Film Session SOP Class interprets the following status codes (from N-CREATE-RSP, N-DELETE-RSP messages):

**Table 66 - Basic Film Session Status Codes**

| Service Status | Meaning  | Error Codes |
|----------------|--|-------------|
| Failure        | Film session SOP instances hierarchy does not contain film box SOP instances | C600        |
|                | Unable to create print job, print queue is full                              | C601        |
|                | Image size is larger than images box size                                    | C603        |
| Warning        | Memory allocation not supported  | B600        |
|                | Film session printing is not supported                                       | B601        |
|                | Film box does not contain image box (empty page)                             | B602        |
| Success        | Film belonging to the film session are accepted for printing                 | 0000        |

## Basic Film Box SOP Class

The Basic Film Box information object definition describes all the user-defined parameter of one film of the film session. The Basic Film Box information description defines the presentation parameters, which are common for all images on a given sheet of film.

The Basic Film Box refers to one or more Image Boxes.

Supported Service Elements as SCU are:

- N-CREATE
- N-ACTION
- N-DELETE

The Basic Film Box SOP Class N-CREATE-RQ message uses the following attributes (the actual values for each attribute depend on DICOM printer configuration within the syngo® MR DICOM print management SCU):

**Table 67 - Basic Film Box N-CREATE attributes**

| Attribute Name                       | Tag         | Usage SCU | Supported Values   |
|--------------------------------------|-------------|-----------|--|
| Image Display Format                 | (2010,0010) | M         | STANDARD\1,1   |
| Referenced Film Session Sequence     | (2010,0500) | M         | n. a.  |
| > Referenced SOP Class UID           | (0008,1150) | M         | 1.2.840.10008.5.1.1.1  |
| > Referenced SOP Instance UID        | (0008,1155) | M         |  |
| Film Orientation                     | (2010,0040) | M         | PORTRAIT   |
| Film Size ID                         | (2010,0050) | M         | 8INX10IN, 10INX12IN, 10INX14IN, 11INX14IN,, 14INX14IN, 14INX17IN, 24CMX24CM, 24CMX30CM |
| Magnification Type                   | (2010,0060) | M         | BILINEAR, CUBIC, NONE, REPLICATE   |
| Border Density                       | (2010,0100) | U         | BLACK, WHITE   |
| Max Density                          | (2010,0130) | U         | > 0  |
| Min Density                          | (2010,0120) | U         | 50 > value > 0   |
| Illumination                         | (2010,015E) | U         | > 0<br>Required if Presentation LUT is present.  |
| Reflective Ambient Light             | (2010,0160) | U         | > 0<br>Required if Presentation LUT is present.  |
| Referenced Presentation LUT Sequence | (2050,0500) | U         |  |

**M** = Mandatory, **U** = User Option

The N-CREATE-RSP message from the Print SCP includes the Referenced Image Box Sequence with SOP Class/Instance UID pairs which will be kept internally and used for the subsequent Basic Image Box SOP Class N-SET-RQ messages.

When all Image Boxes (including parameters) for the film-sheet have been set, the DICOM print manager will issue a N-ACTION-RQ message with the SOP Instance UID of the Basic Film Box and the Action Type ID of 1.

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and used for later requests (N-ACTION-RQ and N-DELETE-RQ) on the Basic Film Box - see below:

| Attribute Name             | Tag                         | Source of Information                                       |
|----------------------------|-----------------------------|---|
| Requested SOP Instance UID | (0000,1000)<br>→(0000,1001) | Affected SOP Instance UID of N-CREATE-RSP on Basic Film Box |

The Basic Film Box SOP Class interprets the following status codes:

**Table 68 - Basic Film Box Status Codes**

| Service Status | Meaning   | Error Codes |
|----------------|---|-------------|
| Failure        | Unable to create print job, print queue is full                         | C601        |
|                | Image size is larger than images box size                               | C603        |
| Warning        | Film box does not contain image box (empty page)                        | B603        |
|                | Requested MinDensity or MaxDensity outside of Printer's operating range | B605        |
| Success        | Film accepted for printing  | 0000        |

## Basic Grayscale Image Box SOP Class

The Basic Grayscale Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Grayscale Image Box information



describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Grayscale Image Box SOP Class uses only the N-SET-RQ with the following attributes:

**Table 69 - Basic Grayscale Image Box N-SET attributes**

| Attribute Name                 | Tag         | Usage SCU | Supported Values      |
|--------------------------------|-------------|-----------|-----------------------|
| Image Position                 | (2020,0010) | M         | 1                     |
| Basic Grayscale Image Sequence | (2020,0110) | M         | n.a.                  |
| > Samples per Pixel            | (0028,0002) | M         | 1                     |
| > Photometric Interpretation   | (0028,0004) | M         | MONOCHROME2           |
| > Rows                         | (0028,0010) | M         | <Printer/Film config> |
| > Columns                      | (0028,0011) | M         | <Printer/Film config> |
| > Pixel Aspect Ratio           | (0028,0034) | M         | (1:1)                 |
| > Bits Allocated               | (0028,0100) | M         | 8, 16                 |
| > Bits Stored                  | (0028,0101) | M         | 8, 12                 |
| > High Bit                     | (0028,0102) | M         | 7, 11                 |
| > Pixel Representation         | (0028,0103) | M         | 0                     |
| > Pixel Data                   | (7FE0,0010) | M         |                       |

**M** = Mandatory

The Grayscale Image Box SOP Class interprets the following status codes:

**Table 70 - Basic Grayscale Image Box Status Codes**

| Service Status | Meaning   | Error Codes |
|----------------|---|-------------|
| Failure        | Image contains more pixel than printer can print in Image Box           | C603        |
|                | Insufficient memory in printer to store the image                       | C605        |
| Warning        | Requested MinDensity or MaxDensity outside of Printer's operating range | B605        |
| Success        |   | 0000        |

## Basic Color Image Box SOP Class

The Basic Color Image Box information object definition is the presentation of an image and image related data in the image area of a film. The Basic Color Image Box information describes the presentation parameters and image pixel data, which apply to a single image of a sheet of film.

The Color Image Box SOP Class uses only the N-SET-RQ with the following attributes:

**Table 71 - Basic Color Image Box N-SET attributes**

| Attribute Name               | Tag         | Usage SCU | Supported Values      |
|------------------------------|-------------|-----------|-----------------------|
| Image Position               | (2020,0010) | M         | 1                     |
| Basic Color Image Sequence   | (2020,0111) | M         | n.a.                  |
| > Samples per Pixel          | (0028,0002) | M         | 3                     |
| > Photometric Interpretation | (0028,0004) | M         | RGB                   |
| > Planar Configuration       | (0028,0006) | M         | 1                     |
| > Rows                       | (0028,0010) | M         | <Printer/Film config> |
| > Columns                    | (0028,0011) | M         | <Printer/Film config> |
| > Pixel Aspect Ratio         | (0028,0034) | M         | (1:1)                 |
| > Bits Allocated             | (0028,0100) | M         | 8                     |
| > Bits Stored                | (0028,0101) | M         | 8,                    |
| > High Bit                   | (0028,0102) | M         | 7                     |
| > Pixel Representation       | (0028,0103) | M         | 0                     |
| > Pixel Data                 | (7FE0,0010) | M         |                       |

**M** = Mandatory

The Color Image Box SOP Class interprets the following status codes:

**Table 72 - Basic Color Image Box Status Codes**

| Service Status | Meaning   | Error Codes |
|----------------|---|-------------|
| Failure        | Image contains more pixel than printer can print in Image Box           | C603        |
|                | Insufficient memory in printer to store the image                       | C605        |
| Warning        | Requested MinDensity or MaxDensity outside of Printer's operating range | B605        |
| Success        |   | 0000        |

## Presentation LUT SOP Class

The Presentation LUT tailors image hardcopy printing for specific modalities, applications and user preferences.

The output of the Presentation LUT is Presentation Values (P-Values). P-Values are approximately related to human perceptual response. They are intended to facilitate common input for hardcopy. P-Values are intended to be independent of the specific class or characteristics of the hardcopy device.

The Presentation LUT SOP Class uses only the N-CREATE-RQ with the following attributes:

**Table 73 - Presentation LUT N-CREATE attribute**

| Attribute Name         | Tag         | Usage SCU | Supported Values |
|------------------------|-------------|-----------|------------------|
| Presentation LUT Shape | (2050,0020) | U         | IDENTITY         |

**U** = User Option

The affected SOP Instance UID received with N-CREATE-RSP message will be kept internally and is used for later requests on the Basic Film Box (N-CREATE-RQ) and on the Presentation LUT (N-DELETE-RQ) - see below:

| Attribute Name             | Tag                         | Source of Information   |
|----------------------------|-----------------------------|---|
| Requested SOP Instance UID | (0000,1000)<br>→(0000,1001) | Affected SOP Instance UID of N-CREATE-RSP on Presentation LUT |

The Presentation LUT SOP Class interprets the following status codes:

**Table 74 - Presentation LUT Status Codes**

| Service Status | Meaning  | Codes |
|----------------|--|-------|
| Success        | Presentation LUT successfully created  | 0000  |
| Warning        | Requested MinDensity or MaxDensity outside of printer's operating range. The printer will use its respective minimum or maximum density value instead. | B605  |

## Printer SOP Class

The Printer SOP Class allows to monitor the status of the hardcopy printer in a synchronous and an asynchronous way.

The SCU uses the mandatory N-EVENT Report DIMSE service to monitor the changes of the printer status in an asynchronous way.

The following returned information is supported:

**Table 75 - Used Printer N-EVENT Report attributes**

| Event-type Name | Event | Attributes          | Tag         | Usage SCU |
|-----------------|-------|---------------------|-------------|-----------|
| Normal          | 1     |                     |             |           |
| Warning         | 2     | Printer Status Info | (2110,0020) | U         |
| Failure         | 3     | Printer Status Info | (2110,0020) | U         |

**U** = User Option

**Table 76 - Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes**

| Attribute Name      | Tag         | Usage SCP | Supported Values                 |
|---------------------|-------------|-----------|----------------------------------|
| Printer Status      | (2110,0010) | M         | NORMAL, FAILURE, WARNING         |
| Printer Status Info | (2110,0020) | M         | See tables in Annex for details. |

**M** = Mandatory

**Note:** For a detailed description on how syngo® MR reacts on different printer status messages, please refer to the appropriate Annex section.

## Printer Job SOP Class

The Print Job SOP Class allows to monitor the execution of the print process.

The syngo® MR DICOM Print Management application supports the optional N-EVENT-REPORT DIMSE Service to receive the changes of the Print Job Status in an asynchronous way.

It can receive Events from the Print SCP asynchronously:

- N-EVENT-REPORT

**Note:** The underlying syngo DICOM Print AE does not support receiving of N-EVENT-REPORT messages from camera during open print sessions. This is typically configurable in the camera setup.

The following information is supported:

**Table 77 - Used Print Job N-EVENT Report attributes**

| Event-type Name | Event | Attributes            | Tag         | Usage SCU  |
|-----------------|-------|-----------------------|-------------|--|
| Normal          | 1     | Execution Status Info | (2100,0030) | U  |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U  |
|                 |       | Printer Name          | (2110,0030) | U  |
| Printing        | 2     | Execution Status Info | (2100,0030) | U  |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U  |
|                 |       | Printer Name          | (2110,0030) | U  |
| Done            | 3     | Execution Status Info | (2100,0030) | U  |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |

| Event-type Name | Event | Attributes            | Tag         | Usage SCU  |
|-----------------|-------|-----------------------|-------------|--|
| Failure         | 4     | Film Session Label    | (2000,0050) | U  |
|                 |       | Printer Name          | (2110,0030) | U  |
|                 |       | Execution Status Info | (2100,0030) | U  |
|                 |       | Print Job ID          | (2100,0010) | --<br>(Print Queue Management SOP Class not supported) |
|                 |       | Film Session Label    | (2000,0050) | U  |
|                 |       | Printer Name          | (2110,0030) | U  |

U = User Option

**Note:** For a detailed description on how syngo® MR reacts on different printer status messages, please refer to the appropriate Annex section".

### 4.2.8.3.2 Activity - Show Device Status

#### 4.2.8.3.2.1 Description and Sequencing of Activity

With no printing activity ongoing ("idle time"), the syngo® MR DICOM Print SCU application will cyclically request the printer status to update the related printer state in the Printing UI.

#### 4.2.8.3.2.2 Proposed Presentation Context

The syngo® MR DICOM application will propose Presentation Contexts as shown in the following table:

**Table 78 - Presentation Context Table "Show Device Status"**

| Presentation Context Table |                        |  |   |      |           |
|----------------------------|------------------------|--|---|------|-----------|
| Abstract Syntax            |                        | Transfer Syntax  |   | Role | Ext. Neg. |
| Name                       | UID                    | Name List  | UID List  |      |           |
| Printer SOP Class          | 1.2.840.10008.5.1.1.16 | Implicit VR Little Endian<br>Explicit VR Little Endian<br>Explicit VR Big Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1<br>1.2.840.10008.1.2.2 | SCU  | None      |

#### 4.2.8.3.2.3 SOP Specific Conformance

The Printer SOP Class allows to monitor the status of the hardcopy printer in a synchronous and an asynchronous way.

The Print SCU AE application will cyclically "ask" the Printer (SCP) for its status synchronously:

- N-GET as SCU

The following information is supported:

**Table 79 - Used Printer N-EVENT Report attributes**

| Event-type Name | Event | Attributes          | Tag         | Usage SCU |
|-----------------|-------|---------------------|-------------|-----------|
| Normal          | 1     |                     |             |           |
| Warning         | 2     | Printer Status Info | (2110,0020) | U         |
| Failure         | 3     | Printer Status Info | (2110,0020) | U         |

**U** = User Option

<modify If detailed status displays of the DICOM Print SCU are described in the Annex, add a reference to the Annex section below

**Table 80 - Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes**

| Attribute Name      | Tag         | Usage SCP | Supported Values                 |
|---------------------|-------------|-----------|----------------------------------|
| Printer Status      | (2110,0010) | M         | NORMAL, FAILURE, WARNING         |
| Printer Status Info | (2110,0020) | M         | See tables in Annex for details. |

**M** = Mandatory

**Note:** For a detailed description on how syngo® MR reacts on different printer status messages, please refer to the Annex section....

>

#### 4.2.8.4 Association Acceptance Policy

The syngo® MR DICOM application does not support Print Management Services as an SCP.

## 4.3 Network Interfaces

### 4.3.1 Physical Network Interface

syngo® MR provides DICOM 3.0 TCP/IP network communication support as defined in Part 8 of the DICOM Standard. The network communication is independent from the physical medium over which TCP/IP executes; it inherits this from the Windows OS system upon which it executes.

### 4.3.2 Additional Protocols

none

### 4.3.3 IPv4 and IPv6 Support

IPv4 and IPv6 are supported. Regarding IPv6 please note, that the complete networking infrastructure in the hospital (firewalls, DNS-Servers ...) must support IPv6 in order to get a functioning communication.

## 4.4 Configuration

### 4.4.1 AE Title/Presentation Address Mapping

AE Titles shall be unique within the hospital. A common way to achieve that is to use the host-name as part of the AE Titles. The string can be up to 16 characters and must not contain any extended characters. Only 7-bit ASCII characters (excluding Control Characters) are allowed according to the DICOM Standard.

#### 4.4.1.1 Local AE Titles

syngo® MR allows configuring AETitles, Ports and Services in any wished way. Default delivery is that all services are using the same AE title and only one port number. In case the connected systems cannot handle this default, the customer service engineer is able to configure for each service its own AE title and Port number.

| Parameter           | Configurable | Default Value  |
|---------------------|--------------|--|
| Default AE title    | Yes          | hostname in uppercase characters; limited to 16 characters |
| Default Port        | Yes          | 104  |
| Default Secure Port | Yes          | 2762   |

#### 4.4.1.2 Remote AE Title/Presentation Address Mapping

##### 4.4.1.2.1 Remote Association Initiators

All relevant remote applications that may setup DICOM associations towards syngo® MR need to be configured in syngo® MR, before the association can be established. This behavior is configurable, but it is recommended, not to change this behavior.

The mapping of external AE Titles to TCP/IP addresses and ports is configurable and initially set at the time of installation by Installation Personnel. Changes can later on also be performed by the local system administrator. The Application Entity Titles and supported transfer syntaxes need to be known for configuration.

To enable a fast and efficient configuration possibility Siemens will deliver templates for known configuration examples, so that the behavior (usage of one AE title, default port numbers, supported services) is determined already through the template.

Remote Application Entities can be configured without restarting the process.

#### 4.4.1.2.2 Remote Association Acceptors

For remote applications that shall be able to accept DICOM associations from syngo® MR, the following information needs to be available:

- Application Entity Title
- Host Name / IP address on which the remote application service runs
- Port number on which the remote application accepts association requests.

The remote system will be indicated in the UI of syngo® MR with a logical name, that is also entered when configuring the node in the administration UI.

To enable a fast and efficient configuration possibility Siemens will deliver templates for known configuration examples, so that the behavior (usage of one AE title, default port numbers, supported services) is determined already through the template.

Remote Application Entities can be configured without restarting the process.

#### 4.4.1.3 Secure DICOM Communication

The system supports configuring the DICOM communication to use secure channel (TLS) between syngo® MR and configured remote nodes. As a security measure the certificate thumbprint or certificate trust chain of the remote nodes shall be added (pinned) to the syngo® MR system to authorize the incoming connection.

Detailed instructions how to set up secure DICOM communication are available in the Administrator Online Help.

Note: The default DICOM port will change to 2762.

If the certificate of remote node contains Enhanced Key Usage (Extended Key Usage) field, then:

- If the remote node acts as DICOM SCP it shall contain Server Authentication (1.3.6.1.5.5.7.3.1)
- If the remote node acts as a DICOM SCU it shall contain Client Authentication (1.3.6.1.5.5.7.3.2)

Otherwise syngo® MR will not accept the certificate.

#### 4.4.2 Parameters

The next table lists configuration parameters, which are true for all Application Entities.

**Table 81: Parameter List**

| Parameter  | Configurable | Default Value |
|--|--------------|---------------|
| max PDU size   | Yes          | 32768 Bytes   |
| time-out for accepting/rejecting an association request      | Yes          | 30 s          |
| time-out for responding to an association open/close request | Yes          | 30 s          |


| Parameter  | Configu-<br>rable | Default Value |
|--|-------------------|---------------|
| time-out for accepting a message over network  | Yes               | 30 s          |
| time-out for waiting for data between TCP/IP-packets   | Yes               | 5 s           |
| time-outs for waiting for a Service Request/Response message from the remote node (Storage SCP/SCU)        | Yes               | 30 s          |
| time-outs for waiting for a Service Request/Response message from the remote node (Query/Retrieve SCP/SCU) | Yes               | 30 s          |
| time-out for waiting for a C-MOVE-RSP  | No                | 1200 s        |
| number of image collection before saving to database   | Yes               | 20            |
| max matches query limit  | Yes               | 100           |
| max number of parallel receiving associations  | Yes               | 12            |



## 5 Media Interchange

### 5.1 Implementation Model

#### 5.1.1 Application Data Flow Diagram



**Figure 6: Media Interchange Application Data Flow Diagram**

syngo® MR provides the functionality to Import or Export DICOM Instances from and to the File System. During export, a DicomDIR may also be generated (user selection). A complete ISO Image ready-to-burn can also be generated. All SOP Classes defined in Table 1 are supported for the Import/Export functionality.


## 5.1.2 Functional definitions of AEs

The syngo® MR application is capable of

- creating a new File-set in the File System (Export to ...)
- importing SOP Instances from the medium onto local storage
- writing the File-sets DICOMDIR information into the file system and joining it to an ISO image.

## 5.1.3 Sequencing of Real-World Activities

Whenever data is written to an external media, syngo® MR creates a DICOMDIR from the selected data and creates an ISO image of the selected data on the local hard disk. Depending on the selected data and options (selected media size, with or without compression) either General Purpose CD profile or DVD-J2K profile are used.



**Figure 7: Sequence diagram – Media creation**

### 5.1.4 File Meta Information for Implementation Class and Version

This section describes the values assigned to the File Meta Information attributes (see part PS 3.10) that pertain to the Implementation Class and Version. The Implementation Class UID and the Implementation Version Name in the File Meta Header are the same as the values specified for networking.

**Table 82: Implementation Class/Version Name - Media Interchange**

|                                      |                   |
|--------------------------------------|-------------------|
| <b>File Meta Information Version</b> | 0001              |
| <b>Implementation Class UID</b>      | 1.3.12.2.1107.5.2 |
| <b>Implementation Version Name</b>   | SYNGO_MR_XA30A    |

## 5.2 AE SPECIFICATIONS

### 5.2.1 Media Storage AE – Specification

syngo® MR provides conformance to the following Application Profiles as an FSC as well as an FSR. The FSU role is only supported on a non-optical storage device (e.g. USB stick).

In addition, augmented conformance is provided to store extra data objects important for the full feature support of the *syngo®*-based products. Details are listed below:

**Table 83: Media - Application Profiles and Real-World Activities**

| Application Profiles Supported | Real-World Activity  | Role             | Service Class Option |
|--------------------------------|--|------------------|----------------------|
| AUG-GEN-CD                     | Browse Directory Information<br>Import into Application<br>Export to local Archive Media | FSR, FSC,<br>FSU | Interchange          |
| AUG-GEN-DVD                    |  |                  |                      |
| AUG- GEN-DVD-J2K               |  |                  |                      |
| AUG- GEN-USB-J2K               |  |                  |                      |
| STD-GEN-CD                     | Browse Directory Information<br>Import into Application<br>Export to local Archive Media | FSR, FSC,<br>FSU | Interchange          |
| STD-GEN-DVD                    |  |                  |                      |
| STD-GEN-DVD-J2K                |  |                  |                      |
| STD-GEN-USB-J2K                |  |                  |                      |

#### 5.2.1.1 Real-World Activities

##### 5.2.1.1.1 Activity “Browse Directory Information”

syngo® MR acts as FSR using the interchange option when requested to read the media directory.

syngo® MR will read the DICOMDIR and insert those directory entries that are valid for the application profiles supported, into a local database. The database then is used for browsing media contents.

**Note:** The "Icon Image Sequence" is also supported in DICOMDIR. But only those Icon Images with "Bits Allocated" (0028,0100) equal to 8 and size of 64x64 or 128x128 pixels are imported into database and are visible in the Browser.

#### 5.2.1.1.2 Real World Activity "Import into Application"

syngo® MR application acts as FSR using the interchange option when requested to read SOP Instances from the medium into the application.

The SOP Instance selected from the media directory will be copied into the running Application. Only SOP Instances, that are valid for the application profile supported and supported by syngo® MR can be retrieved from media.

#### 5.2.1.1.3 Real-World Activity "Export to local Archive Media"

The syngo® MR application acts as FSU (for media with existing DICOM file-set) or FSC (media not initialized) using the interchange option when requested to copy SOP Instances from the local storage to local Archive Medium. The activity as FSU is only possible as long as the local burning SW of syngo® MR has not already processed the generated ISO file.

The syngo® MR application will receive a list of SOP Instances to be copied to the local archive medium. Depending on the profile selected (Standard: uncompressed, with DICOMDIR; Patient: compressed with DICOMDIR) the SOP Instances will be taken, and an ISO file is being generated that includes the DICOMDIR and the corresponding objects.

It is then up to the syngo® MR local configuration (if equipped with a local media burner) to burn the ISO file to the appropriate media.

#### 5.2.1.2 SOP Classes and Transfer Syntaxes

These Application Profiles are based on the Media Storage Service Class with the Interchange Option. syngo® MR provides Standard Conformance to the SOP Classes listed in "Table 1: Network Services" section "SOP Classes Created by syngo® MR" and "SOP Classes Managed by syngo® MR" in the ["Conformance Statement Overview"](#).

Using the Application Profiles supporting compression (AUG- GEN-DVD-J2K, AUG- GEN-USB-J2K, STD-GEN-DVD-J2K, STD-GEN-USB-J2K) the following Transfer Syntaxes are supported:

**Table 84: Transfer Syntaxes for STD-GEN-DVD-J2K and STD-GEN-USB-J2K**

| UID value              | Transfer Syntax  | Image Objects | Non-Image Objects |
|------------------------|--|---------------|-------------------|
| 1.2.840.10008.1.2.1    | Explicit Value Representation Little Endian native                                       | yes           | yes               |
| 1.2.840.10008.1.2.4.50 | JPEG Baseline (Process 1) lossy compressed   | yes           | no                |
| 1.2.840.10008.1.2.4.51 | JPEG Extended (Process 2 & 4) lossy compressed   | yes           | no                |
| 1.2.840.10008.1.2.4.70 | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) lossless compressed | yes           | no                |

| UID value              | Transfer Syntax  | Image Objects | Non-Image Objects |
|------------------------|--|---------------|-------------------|
| 1.2.840.10008.1.2.4.90 | JPEG 2000 Image Compression (Lossless Only) compressed | yes           | no                |
| 1.2.840.10008.1.2.4.91 | JPEG 2000 Image Compression lossy compressed           | yes           | no                |
| 1.2.840.10008.1.2.5    | RLE Lossless compressed                                | yes           | no                |

Using the Application Profiles that do not support compression (AUG- GEN-DVD, AUG- GEN-USB, STD-GEN-DVD, STD-GEN-USB) only Explicit Value Representation Little Endian (1.2.840.10008.1.2.1) is supported.

## 5.3 AUGMENTED AND PRIVATE APPLICATION PROFILES

### 5.3.1 Augmented Application Profiles

The standard application profiles are augmented with private object Siemens CSA Non-Image.

**Table 85: Private SOP Classes and Transfer Syntaxes for Augmented Media Profiles**

| Information Object Definition | SOP Class UID       | Transfer Syntax UID   | FSC | FSR |
|-------------------------------|---------------------|---|-----|-----|
| CSA Non-Image Storage         | 1.3.12.2.1107.5.9.1 | Explicit VR Little Endian Uncompressed<br>1.2.840.10008.1.2.1 | O   | M   |

The Siemens non-image is typically used for raw data and 3D private data.

## 5.4 MEDIA CONFIGURATION

none

## 6 Support of Extended Character Sets

The syngo® MR DICOM application supports the following character sets as defined in the four tables below:

**Table 86: Single-Byte Character Sets without Code Extension**

| Character Set Description | Defined Term | ISO registration number | Character Set        |
|---------------------------|--------------|-------------------------|----------------------|
| Default repertoire        | None         | ISO_IR 6                | ISO 646:             |
| Latin alphabet No. 1      | ISO_IR 100   | ISO_IR 100              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646:             |
| Latin alphabet No. 2      | ISO_IR 101   | ISO_IR 101              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Latin alphabet No. 3      | ISO_IR 109   | ISO_IR 109              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Latin alphabet No. 4      | ISO_IR 110   | ISO_IR 110              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Cyrillic                  | ISO_IR 144   | ISO_IR 144              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Arabic                    | ISO_IR 127   | ISO_IR 127              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Greek                     | ISO_IR 126   | ISO_IR 126              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Hebrew                    | ISO_IR 138   | ISO_IR 138              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Latin alphabet No. 5      | ISO_IR 148   | ISO_IR 148              | Supplementary set    |
|                           |              | ISO_IR 6                | ISO 646              |
| Japanese                  | ISO_IR 13    | ISO_IR 13               | JIS X 0201: Katakana |
|                           |              | ISO_IR 14               | JIS X 0201: Romaji   |
| Thai                      | ISO_IR166    | ISO_IR166               | TIS 620-253 (1990)   |
|                           |              | ISO_IR 6                | ISO 646              |

**Table 87: Single-Byte Characters Sets with Code Extension**

| Character Set Description | Defined Term    | Standard for Code Extension | ESC sequence    | ISO registration number | Character Set     |
|---------------------------|-----------------|-----------------------------|-----------------|-------------------------|-------------------|
| Default repertoire        | ISO 2022 IR 6   | ISO 2022                    | ESC 02/08 04/02 | ISO-IR 6                | ISO 646           |
| Latin alphabet No.1       | ISO 2022 IR 100 | ISO 2022                    | ESC 02/13 04/01 | ISO-IR 100              | Supplementary set |
|                           |                 | ISO 2022                    | ESC 02/08 04/02 | ISO-IR 6                | ISO 646           |
| Latin alphabet No.2       | ISO 2022 IR 101 | ISO 2022                    | ESC 02/13 04/02 | ISO-IR 101              | Supplementary set |
|                           |                 | ISO 2022                    | ESC 02/08 04/02 | ISO-IR 6                | ISO 646           |
| Latin alphabet No.3       | ISO 2022 IR 109 | ISO 2022                    | ESC 02/13 04/03 | ISO-IR 109              | Supplementary set |

|                                     |          |                 |            |                   |
|-------------------------------------|----------|-----------------|------------|-------------------|
|                                     | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6   | ISO 646           |
| Latin alphabet No.4 ISO 2022 IR 110 | ISO 2022 | ESC 02/13 04/04 | ISO-IR 110 | Supplementary set |
|                                     | ISO 2022 | ESC 02/08 04/02 | ISO-IR 6   | ISO 646           |

#### Multi-Byte Character Sets without Code Extension

**Table 88: Multi-Byte Character Sets without Code Extension**

| Character Set Description | Defined Term | ISO registration number | Character Set   |
|---------------------------|--------------|-------------------------|---|
| Unicode                   | ISO_IR 192   | ISO 10646               | Unicode in UTF-8                                      |
| Chinese                   | GB18030      | GB18030                 | GB 18030-2000 (China Association for Standardization) |

**Table 89: Multi-Byte Character Sets with Code Extension**

| Character Set Description | Defined Term    | Standard for Code Extension | ESC sequence          | ISO registration number | Character Set                       |
|---------------------------|-----------------|-----------------------------|-----------------------|-------------------------|-------------------------------------|
| Japanese                  | ISO 2022 IR 159 | ISO 2022                    | ESC 02/04 02/08 04/04 | ISO-IR 159<br>ISO-IR 87 | JIS X 0212: Supplementary Kanji set |
| Korean                    | ISO 2022 IR 149 | ISO 2022                    | ESC 02/04 02/09 04/03 | ISO-IR 149              | KS X 1001: Hangul and Hanja         |

All Special Character Sets (SCS) listed above are supported for incoming Data. When creating new Instances, the system will use the default SCS (or SCS List) configured on the machine.

When there is a mismatch between the given character set in attribute (0008,0005) and the characters in an IOD received by the system, then the following measures are taken to make the characters DICOM conform:

- Convert each illegal character to a '?'.

There are three categories of character sets which have to be differentiated because of their different encoding formats:

- Conventional ISO character sets: ISO\_IR 6, ISO 2022 IR 6, ISO\_IR 100, etc.  
→ encoded in ISO 2022
- ISO\_IR 192 → encoded in UTF-8
- GB18030 → encoded in GB18030

It is not possible to recognize the following mismatches automatically on receiving or importing:

- An attribute value is encoded in ISO\_IR 192 ↔ (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in GB18030 ↔ (0008,0005) contains a conventional ISO character set as primary character set
- An attribute value is encoded in ISO 2022 ↔ (0008,0005) contains ISO\_IR 192
- An attribute value is encoded in ISO 2022 ↔ (0008,0005) contains GB18030

An IOD that contains one of the above mentioned inconsistencies is not DICOM conform. As these kinds of inconsistencies cannot be recognized by the system, the IOD will not be rejected but the character data might be corrupted.

syngo® MR supports Kanji characters in the byte zones after 74 (79, 7A, 7B and 7C).



## 7 Attribute confidentiality profiles

### 7.1 De-identification

The syngo® MR application can de-identify attributes, when exporting to Media. Three different levels of de-identification are supported:

- Full de-identification
- Reduced de-identification
- Service de-identification

The user needs to select the appropriate de-identification level during export.

For full and reduced de-identification private attributes are not included in anonymized studies. For service de-identification all private attributes are included in anonymized studies.

Note: reduced anonymization applies the following options: Retain UIDs , Patient Chars, Long, Full Dates

In the following table for attributes marked with:

- 'Yes' - data are anonymized
- 'No' - data are kept

**Table 90: Application Level Confidentiality Profile Attributes (standard tags)**

| DICOM Tag   | Attribute Name                                | Full | Re-duced | Service |
|-------------|---|------|----------|---------|
| (0000,1000) | Affected SOP Instance UID                     | Yes  | No       | No      |
| (0000,1001) | Requested SOP Instance UID                    | Yes  | No       | No      |
| (0002,0003) | Media Storage SOP Instance UID                | Yes  | No       | No      |
| (0004,1511) | Referenced SOP Instance UID in File           | Yes  | No       | No      |
| (0008,0014) | Instance Creator UID                          | Yes  | No       | No      |
| (0008,0015) | Instance Coercion DateTime                    | Yes  | No       | No      |
| (0008,0018) | SOP Instance UID                              | Yes  | No       | No      |
| (0008,0020) | Study Date                                    | Yes  | No       | No      |
| (0008,0021) | Series Date                                   | Yes  | No       | No      |
| (0008,0022) | Acquisition Date                              | Yes  | No       | No      |
| (0008,0023) | Content Date                                  | Yes  | No       | No      |
| (0008,0024) | Overlay Date                                  | Yes  | No       | No      |
| (0008,0025) | Curve Date                                    | Yes  | No       | No      |
| (0008,002A) | Acquisition DateTime                          | Yes  | No       | No      |
| (0008,0030) | Study Time                                    | Yes  | No       | No      |
| (0008,0031) | Series Time                                   | Yes  | No       | No      |
| (0008,0032) | Acquisition Time                              | Yes  | No       | No      |
| (0008,0033) | Content Time                                  | Yes  | No       | No      |
| (0008,0034) | Overlay Time                                  | Yes  | No       | No      |
| (0008,0035) | Curve Time                                    | Yes  | No       | No      |
| (0008,0050) | Accession Number                              | Yes  | Yes      | No      |
| (0008,0058) | Failed SOP Instance UID List                  | Yes  | No       | No      |
| (0008,0080) | Institution Name                              | Yes  | Yes      | No      |
| (0008,0081) | Institution Address                           | Yes  | Yes      | No      |
| (0008,0082) | Institution Code Sequence                     | Yes  | Yes      | No      |
| (0008,0090) | Referring Physician's Name                    | Yes  | Yes      | Yes     |
| (0008,0092) | Referring Physician's Address                 | Yes  | Yes      | Yes     |
| (0008,0094) | Referring Physician's Telephone Numbers       | Yes  | Yes      | Yes     |
| (0008,0096) | Referring Physician's Identification Sequence | Yes  | Yes      | No      |
| (0008,010D) | Context Group Extension Creator UID           | Yes  | No       | No      |

| DICOM Tag   | Attribute Name                                    | Full | Re-duced | Service |
|-------------|---|------|----------|---------|
| (0008,0201) | Timezone Offset From UTC                          | Yes  | No       | No      |
| (0008,1010) | Station Name                                      | Yes  | Yes      | Yes     |
| (0008,1030) | Study Description                                 | Yes  | Yes      | No      |
| (0008,103E) | Series Description                                | Yes  | Yes      | No      |
| (0008,1040) | Institutional Department Name                     | Yes  | Yes      | No      |
| (0008,1048) | Physician(s) of Record                            | Yes  | Yes      | Yes     |
| (0008,1049) | Physician(s) of Record Identification Sequence    | Yes  | Yes      | No      |
| (0008,1050) | Performing Physicians' Name                       | Yes  | Yes      | Yes     |
| (0008,1052) | Performing Physicians' Identification Sequence    | Yes  | Yes      | No      |
| (0008,1060) | Name of Physician(s) Reading Study                | Yes  | Yes      | Yes     |
| (0008,1062) | Physician Reading Study Identification Sequence   | Yes  | Yes      | No      |
| (0008,1070) | Operators' Name                                   | Yes  | Yes      | Yes     |
| (0008,1072) | Operators' Identification Sequence                | Yes  | Yes      | No      |
| (0008,1080) | Admitting Diagnoses Description                   | Yes  | Yes      | No      |
| (0008,1084) | Admitting Diagnoses Code Sequence                 | Yes  | Yes      | No      |
| (0008,1110) | Referenced Study Sequence                         | Yes  | No       | No      |
| (0008,1111) | Referenced Performed Procedure Step Sequence      | Yes  | No       | No      |
| (0008,1120) | Referenced Patient Sequence                       | Yes  | Yes      | No      |
| (0008,1140) | Referenced Image Sequence                         | Yes  | No       | No      |
| (0008,1155) | Referenced SOP Instance UID                       | Yes  | No       | No      |
| (0008,1195) | Transaction UID                                   | Yes  | No       | No      |
| (0008,2111) | Derivation Description                            | Yes  | No       | No      |
| (0008,2112) | Source Image Sequence                             | Yes  | No       | No      |
| (0008,3010) | Irradiation Event UID                             | Yes  | No       | No      |
| (0008,4000) | Identifying Comments                              | Yes  | Yes      | No      |
| (0008,9123) | Creator Version UID                               | Yes  | No       | No      |
| (0010,0010) | Patient's Name                                    | Yes  | Yes      | Yes     |
| (0010,0020) | Patient ID  | Yes  | Yes      | Yes     |
| (0010,0021) | Issuer of Patient ID                              | Yes  | Yes      | No      |
| (0010,0030) | Patient's Birth Date                              | Yes  | Yes      | Yes     |
| (0010,0032) | Patient's Birth Time                              | Yes  | Yes      | No      |
| (0010,0040) | Patient's Sex                                     | Yes  | No       | No      |
| (0010,0050) | Patient's Insurance Plan Code Sequence            | Yes  | Yes      | Yes     |
| (0010,0101) | Patient's Primary Language Code Sequence          | Yes  | Yes      | Yes     |
| (0010,0102) | Patient's Primary Language Modifier Code Sequence | Yes  | Yes      | Yes     |
| (0010,1000) | Other Patient IDs                                 | Yes  | Yes      | Yes     |
| (0010,1001) | Other Patient Names                               | Yes  | Yes      | Yes     |
| (0010,1002) | Other Patient IDs Sequence                        | Yes  | Yes      | Yes     |
| (0010,1005) | Patient's Birth Name                              | Yes  | Yes      | Yes     |
| (0010,1010) | Patient's Age                                     | Yes  | No       | No      |
| (0010,1020) | Patient's Size                                    | Yes  | No       | No      |
| (0010,1030) | Patient's Weight                                  | Yes  | No       | No      |
| (0010,1040) | Patient Address                                   | Yes  | Yes      | Yes     |
| (0010,1050) | Insurance Plan Identification                     | Yes  | Yes      | No      |
| (0010,1060) | Patient's Mother's Birth Name                     | Yes  | Yes      | Yes     |
| (0010,1080) | Military Rank                                     | Yes  | Yes      | No      |
| (0010,1081) | Branch of Service                                 | Yes  | Yes      | No      |
| (0010,1090) | Medical Record Locator                            | Yes  | Yes      | No      |
| (0010,1100) | Referenced Patient Photo Sequence                 | Yes  | Yes      | No      |
| (0010,2000) | Medical Alerts                                    | Yes  | Yes      | No      |
| (0010,2110) | Allergies   | Yes  | Yes      | No      |
| (0010,2150) | Country of Residence                              | Yes  | Yes      | No      |
| (0010,2152) | Region of Residence                               | Yes  | Yes      | No      |
| (0010,2154) | Patient's Telephone Number                        | Yes  | Yes      | Yes     |
| (0010,2160) | Ethnic Group                                      | Yes  | No       | No      |
| (0010,2180) | Occupation  | Yes  | Yes      | No      |

| DICOM Tag   | Attribute Name                            | Full | Re-duced | Service |
|-------------|---|------|----------|---------|
| (0010,21A0) | Smoking Status                            | Yes  | No       | No      |
| (0010,21B0) | Additional Patient's History              | Yes  | Yes      | Yes     |
| (0010,21C0) | Pregnancy Status                          | Yes  | No       | No      |
| (0010,21D0) | Last Menstrual Date                       | Yes  | No       | No      |
| (0010,21F0) | Patient's Religious Preference            | Yes  | Yes      | No      |
| (0010,2203) | Patient Sex Neutered                      | Yes  | No       | No      |
| (0010,2297) | Responsible Person                        | Yes  | Yes      | No      |
| (0010,2299) | Responsible Organization                  | Yes  | Yes      | No      |
| (0010,4000) | Patient Comments                          | Yes  | Yes      | Yes     |
| (0018,0010) | Contrast Bolus Agent                      | Yes  | Yes      | No      |
| (0018,1000) | Device Serial Number                      | Yes  | Yes      | No      |
| (0018,1002) | Device UID                                | Yes  | No       | No      |
| (0018,1004) | Plate ID                                  | Yes  | Yes      | No      |
| (0018,1005) | Generator ID                              | Yes  | Yes      | No      |
| (0018,1007) | Cassette ID                               | Yes  | Yes      | No      |
| (0018,1008) | Gantry ID                                 | Yes  | Yes      | No      |
| (0018,1030) | Protocol Name                             | Yes  | Yes      | No      |
| (0018,1400) | Acquisition Device Processing Description | Yes  | Yes      | No      |
| (0018,2042) | Target UID                                | Yes  | No       | No      |
| (0018,4000) | Acquisition Comments                      | Yes  | Yes      | No      |
| (0018,700A) | Detector ID                               | Yes  | Yes      | No      |
| (0018,9424) | Acquisition Protocol Description          | Yes  | Yes      | No      |
| (0018,9516) | Start Acquisition DateTime                | Yes  | No       | No      |
| (0018,9517) | End Acquisition DateTime                  | Yes  | No       | No      |
| (0018,A003) | Contribution Description                  | Yes  | Yes      | No      |
| (0020,000D) | Study Instance UID                        | Yes  | No       | No      |
| (0020,000E) | Series Instance UID                       | Yes  | No       | No      |
| (0020,0010) | Study ID                                  | Yes  | Yes      | No      |
| (0020,0052) | Frame of Reference UID                    | Yes  | No       | No      |
| (0020,0200) | Synchronization Frame of Reference UID    | Yes  | No       | No      |
| (0020,3401) | Modifying Device ID                       | Yes  | Yes      | No      |
| (0020,3404) | Modifying Device Manufacturer             | Yes  | Yes      | No      |
| (0020,3406) | Modified Image Description                | Yes  | Yes      | No      |
| (0020,4000) | Image Comments                            | Yes  | Yes      | No      |
| (0020,9158) | Frame Comments                            | Yes  | Yes      | No      |
| (0020,9161) | Concatenation UID                         | Yes  | No       | No      |
| (0020,9164) | Dimension Organization UID                | Yes  | No       | No      |
| (0028,1199) | Palette Color Lookup Table UID            | Yes  | No       | No      |
| (0028,1214) | Large Palette Color Lookup Table UID      | Yes  | No       | No      |
| (0028,4000) | Image Presentation Comments               | Yes  | Yes      | No      |
| (0032,0012) | Study ID Issuer                           | Yes  | Yes      | No      |
| (0032,1020) | Scheduled Study Location                  | Yes  | Yes      | No      |
| (0032,1021) | Scheduled Study Location AE Title         | Yes  | Yes      | No      |
| (0032,1030) | Reason for Study                          | Yes  | Yes      | No      |
| (0032,1032) | Requesting Physician                      | Yes  | Yes      | No      |
| (0032,1033) | Requesting Service                        | Yes  | Yes      | No      |
| (0032,1060) | Requested Procedure Description           | Yes  | Yes      | No      |
| (0032,1070) | Requested Contrast Agent                  | Yes  | Yes      | No      |
| (0032,4000) | Study Comments                            | Yes  | Yes      | No      |
| (0038,0004) | Referenced Patient Alias Sequence         | Yes  | Yes      | No      |
| (0038,0010) | Admission ID                              | Yes  | Yes      | No      |
| (0038,0011) | Issuer of Admission ID                    | Yes  | Yes      | No      |
| (0038,001E) | Scheduled Patient Institution Residence   | Yes  | Yes      | No      |
| (0038,0020) | Admitting Date                            | Yes  | No       | No      |
| (0038,0021) | Admitting Time                            | Yes  | No       | No      |
| (0038,0040) | Discharge Diagnosis Description           | Yes  | Yes      | No      |

| DICOM Tag   | Attribute Name  | Full | Re-duced | Service |
|-------------|---|------|----------|---------|
| (0038,0050) | Special Needs   | Yes  | Yes      | No      |
| (0038,0060) | Service Episode ID  | Yes  | Yes      | No      |
| (0038,0061) | Issuer of Service Episode ID  | Yes  | Yes      | No      |
| (0038,0062) | Service Episode Description   | Yes  | Yes      | No      |
| (0038,0300) | Current Patient Location  | Yes  | Yes      | No      |
| (0038,0400) | Patient's Institution Residence                                     | Yes  | Yes      | No      |
| (0038,0500) | Patient State   | Yes  | Yes      | No      |
| (0038,4000) | Visit Comments  | Yes  | Yes      | No      |
| (0040,0001) | Scheduled Station AE Title  | Yes  | Yes      | No      |
| (0040,0002) | Scheduled Procedure Step Start Date                                 | Yes  | No       | No      |
| (0040,0003) | Scheduled Procedure Step Start Time                                 | Yes  | No       | No      |
| (0040,0004) | Scheduled Procedure Step End Date                                   | Yes  | No       | No      |
| (0040,0005) | Scheduled Procedure Step End Time                                   | Yes  | No       | No      |
| (0040,0006) | Scheduled Performing Physician Name                                 | Yes  | Yes      | No      |
| (0040,0007) | Scheduled Procedure Step Description                                | Yes  | Yes      | No      |
| (0040,000B) | Scheduled Performing Physician Identification Sequence              | Yes  | Yes      | No      |
| (0040,0010) | Scheduled Station Name  | Yes  | Yes      | No      |
| (0040,0011) | Scheduled Procedure Step Location                                   | Yes  | Yes      | No      |
| (0040,0012) | Pre-Medication  | Yes  | Yes      | No      |
| (0040,0241) | Performed Station AE Title  | Yes  | Yes      | No      |
| (0040,0242) | Performed Station Name  | Yes  | Yes      | No      |
| (0040,0243) | Performed Location  | Yes  | Yes      | No      |
| (0040,0244) | Performed Procedure Step Start Date                                 | Yes  | No       | No      |
| (0040,0245) | Performed Procedure Step Start Time                                 | Yes  | No       | No      |
| (0040,0250) | Performed Procedure Step End Date                                   | Yes  | No       | No      |
| (0040,0251) | Performed Procedure Step End Time                                   | Yes  | No       | No      |
| (0040,0253) | Performed Procedure Step ID   | Yes  | Yes      | No      |
| (0040,0254) | Performed Procedure Step Description                                | Yes  | Yes      | No      |
| (0040,0275) | Request Attributes Sequence   | Yes  | Yes      | No      |
| (0040,0280) | Comments on Performed Procedure Step                                | Yes  | Yes      | No      |
| (0040,0555) | Acquisition Context Sequence  | Yes  | Yes      | No      |
| (0040,1001) | Requested Procedure ID  | Yes  | Yes      | No      |
| (0040,1004) | Patient Transport Arrangements                                      | Yes  | Yes      | No      |
| (0040,1005) | Requested Procedure Location  | Yes  | Yes      | No      |
| (0040,1010) | Names of Intended Recipient of Results                              | Yes  | Yes      | No      |
| (0040,1011) | Intended Recipients of Results Identification Sequence              | Yes  | Yes      | No      |
| (0040,1101) | Person Identification Code Sequence                                 | Yes  | Yes      | No      |
| (0040,1102) | Person Address  | Yes  | Yes      | No      |
| (0040,1103) | Person Telephone Numbers  | Yes  | Yes      | No      |
| (0040,1400) | Requested Procedure Comments  | Yes  | Yes      | No      |
| (0040,2001) | Reason for Imaging Service Request                                  | Yes  | Yes      | No      |
| (0040,2008) | Order Entered By  | Yes  | Yes      | No      |
| (0040,2009) | Order Enterer Location  | Yes  | Yes      | No      |
| (0040,2010) | Order Callback Phone Number   | Yes  | Yes      | No      |
| (0040,2016) | Placer Order Number of Imaging Service Request                      | Yes  | Yes      | No      |
| (0040,2017) | Filler Order Number of Imaging Service Request                      | Yes  | Yes      | No      |
| (0040,2400) | Imaging Service Request Comments                                    | Yes  | Yes      | No      |
| (0040,3001) | Confidentiality Constraint on Patient Data Description              | Yes  | Yes      | No      |
| (0040,4005) | Scheduled Procedure Step Start DateTime                             | Yes  | No       | No      |
| (0040,4010) | Scheduled Procedure Step Modification DateTime                      | Yes  | No       | No      |
| (0040,4011) | Expected Completion Date Time                                       | Yes  | No       | No      |
| (0040,4023) | Referenced General Purpose Scheduled Procedure Step Transaction UID | Yes  | No       | No      |
| (0040,4025) | Scheduled Station Name Code Sequence                                | Yes  | Yes      | No      |
| (0040,4027) | Scheduled Station Geographic Location Code Sequence                 | Yes  | Yes      | No      |
| (0040,4028) | Performed Station Name Code Sequence                                | Yes  | Yes      | No      |

| DICOM Tag   | Attribute Name                                      | Full | Re-duced | Service |
|-------------|---|------|----------|---------|
| (0040,4030) | Performed Station Geographic Location Code Sequence | Yes  | Yes      | No      |
| (0040,4034) | Scheduled Human Performers Sequence                 | Yes  | Yes      | No      |
| (0040,4035) | Actual Human Performers Sequence                    | Yes  | Yes      | No      |
| (0040,4036) | Human Performers Organization                       | Yes  | Yes      | No      |
| (0040,4037) | Human Performers Name                               | Yes  | Yes      | No      |
| (0040,4050) | Performed Procedure Step Start DateTime             | Yes  | No       | No      |
| (0040,4051) | Performed Procedure Step End DateTime               | Yes  | No       | No      |
| (0040,4052) | Procedure Step Cancellation DateTime                | Yes  | No       | No      |
| (0040,A027) | Verifying Organization                              | Yes  | Yes      | No      |
| (0040,A073) | Verifying Observer Sequence                         | Yes  | Yes      | No      |
| (0040,A075) | Verifying Observer Name                             | Yes  | Yes      | No      |
| (0040,A078) | Author Observer Sequence                            | Yes  | Yes      | No      |
| (0040,A07A) | Participant Sequence                                | Yes  | Yes      | No      |
| (0040,A07C) | Custodial Organization Sequence                     | Yes  | Yes      | No      |
| (0040,A088) | Verifying Observer Identification Code Sequence     | Yes  | Yes      | No      |
| (0040,A123) | Person Name   | Yes  | Yes      | No      |
| (0040,A124) | UID   | Yes  | Yes      | No      |
| (0040,A171) | Observation UID                                     | Yes  | No       | No      |
| (0040,A172) | Referenced Observation UID (Trial)                  | Yes  | No       | No      |
| (0040,A192) | Observation Date (Trial)                            | Yes  | No       | No      |
| (0040,A193) | Observation Time (Trial)                            | Yes  | No       | No      |
| (0040,A307) | Current Observer (Trial)                            | Yes  | Yes      | No      |
| (0040,A352) | Verbal Source (Trial)                               | Yes  | Yes      | No      |
| (0040,A353) | Address (Trial)                                     | Yes  | Yes      | No      |
| (0040,A354) | Telephone Number (Trial)                            | Yes  | Yes      | Yes     |
| (0040,A358) | Verbal Source Identifier Code Sequence (Trial)      | Yes  | Yes      | No      |
| (0040,A402) | Observation Subject UID (Trial)                     | Yes  | No       | No      |
| (0040,A730) | Content Sequence                                    | Yes  | Yes      | No      |
| (0040,DB0C) | Template Extension Organization UID                 | Yes  | No       | No      |
| (0040,DB0D) | Template Extension Creator UID                      | Yes  | No       | No      |
| (0070,0001) | Graphic Annotation Sequence                         | Yes  | Yes      | No      |
| (0070,0084) | Content Creator's Name                              | Yes  | Yes      | No      |
| (0070,0086) | Content Creator's Identification Code Sequence      | Yes  | Yes      | No      |
| (0070,031A) | Fiducial UID  | Yes  | No       | No      |
| (0088,0140) | Storage Media Fileset UID                           | Yes  | No       | No      |
| (0088,0200) | Icon Image Sequence                                 | Yes  | Yes      | No      |
| (0088,0904) | Topic Title   | Yes  | Yes      | No      |
| (0088,0906) | Topic Subject                                       | Yes  | Yes      | No      |
| (0088,0910) | Topic Author  | Yes  | Yes      | No      |
| (0088,0912) | Topic Keywords                                      | Yes  | Yes      | No      |
| (0400,0100) | Digital Signature UID                               | Yes  | Yes      | No      |
| (0400,0402) | Referenced Digital Signature Sequence               | Yes  | Yes      | No      |
| (0400,0403) | Referenced SOP Instance MAC Sequence                | Yes  | Yes      | No      |
| (0400,0404) | MAC   | Yes  | Yes      | No      |
| (0400,0550) | Modified Attributes Sequence                        | Yes  | Yes      | No      |
| (0400,0561) | Original Attributes Sequence                        | Yes  | Yes      | No      |
| (2030,0020) | Text String   | Yes  | Yes      | No      |
| (3006,0024) | Referenced Frame of Reference UID                   | Yes  | No       | No      |
| (3006,00C2) | Related Frame of Reference UID                      | Yes  | No       | No      |
| (3008,0105) | Source Serial Number                                | Yes  | No       | No      |
| (300A,0013) | Dose Reference UID                                  | Yes  | No       | No      |
| (300E,0008) | Reviewer Name                                       | Yes  | Yes      | No      |
| (4000,0010) | Arbitrary   | Yes  | Yes      | No      |
| (4000,4000) | Text Comments                                       | Yes  | Yes      | No      |
| (4008,0042) | Results ID Issuer                                   | Yes  | Yes      | No      |
| (4008,0102) | Interpretation Recorder                             | Yes  | Yes      | No      |

| DICOM Tag   | Attribute Name                       | Full | Re-duced | Service |
|-------------|--------------------------------------|------|----------|---------|
| (4008,010A) | Interpretation Transcriber           | Yes  | Yes      | No      |
| (4008,010B) | Interpretation Text                  | Yes  | Yes      | No      |
| (4008,010C) | Interpretation Author                | Yes  | Yes      | No      |
| (4008,0111) | Interpretation Approver Sequence     | Yes  | Yes      | No      |
| (4008,0114) | Physician Approving Interpretation   | Yes  | Yes      | No      |
| (4008,0115) | Interpretation Diagnosis Description | Yes  | Yes      | No      |
| (4008,0118) | Results Distribution List Sequence   | Yes  | Yes      | No      |
| (4008,0119) | Distribution Name                    | Yes  | Yes      | No      |
| (4008,011A) | Distribution Address                 | Yes  | Yes      | No      |
| (4008,0202) | Interpretation ID Issuer             | Yes  | Yes      | No      |
| (4008,0300) | Impressions                          | Yes  | Yes      | No      |
| (4008,4000) | Results Comments                     | Yes  | Yes      | No      |
| (50xx,xxxx) | Curve Data                           | Yes  | Yes      | No      |
| (60xx,0100) | Overlay Bits Allocated               | Yes  | Yes      | No      |
| (60xx,0102) | Overlay Bit Position                 | Yes  | Yes      | No      |
| (60xx,3000) | Overlay Data                         | Yes  | Yes      | No      |
| (60xx,4000) | Overlay Comments                     | Yes  | Yes      | No      |
| (FFFA,FFFA) | Digital Signatures Sequence          | Yes  | Yes      | Yes     |
| (FFFC,FFFC) | Data Set Trailing Padding            | Yes  | Yes      | Yes     |

**Table 91: Application Level Confidentiality Profile Attributes (private tags)**

| DICOM Tag                              | Attribute Name                   | Full | Reduced | Service |
|--|----------------------------------|------|---------|---------|
| (0019, SIEMENS CT VA0 COAD, 90)        | Osteo offset                     | Yes  | No      | No      |
| (0019, SIEMENS CT VA0 COAD, 92)        | Osteo Regression Line Slope      | Yes  | No      | No      |
| (0019, SIEMENS CT VA0 COAD, 93)        | Osteo Regression Line Inter-cept | Yes  | No      | No      |
| (0019, SIEMENS CT VA0 COAD, 96)        | Osteo Phantom Number             | Yes  | No      | No      |
| (0043, GEMS_PARM_01, 1E)               | GE Delta Start Time              | Yes  | No      | No      |
| (0029, SIEMENS CT EXAM IMAGE, 49)      | Metal Artifact Reduction Type    | Yes  | No      | No      |
| (0029, SIEMENS CSA ENVELOPE, 10)       | Syngo Report Data                | Yes  | No      | No      |
| (0029, SIEMENS CSA ENVELOPE, 11)       | Syngo Report Presentation        | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 08)         | Modality Image Header Type       | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 09)         | Modality Image Header Version    | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 10)         | Modality Image Header Info       | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 18)         | Modality Series Header Type      | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 19)         | Modality Series Header Version   | Yes  | No      | No      |
| (0029, SIEMENS CSA HEADER, 20)         | Modality Series Header Info      | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 40)      | Application Header Sequence      | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 41)      | Application Header Type          | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 42)      | Application Header ID            | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 43)      | Application Header Version       | Yes  | No      | No      |
| (0029, SIEMENS MEDCOM HEADER, 44)      | Application Header Info          | Yes  | No      | No      |
| (0029, SIEMENS CT APPL DATASET, 00)    | Dual Energy Algorithm Parameters | Yes  | No      | No      |
| (0029, SIEMENS CT APPL ALG PARAMS, 20) | Perfusion Result Set Id          | Yes  | No      | No      |
| (0021, SIEMENS MR SDS 01, 0C)          | Positive PCS Directions          | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDS 01, 5E)          | Field Of View Text               | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDS 01, 5F)          | Relative Table Position Text     | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDS 01, FE)          | Series Data Sequence             | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDI 02, 4F)          | Coil String                      | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDI 02, 56)          | PAT Mode Text                    | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDI 02, 58)          | Acquisition Matrix Text          | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDI 02, 88)          | Slice Position                   | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDI 02, 89)          | Slice Position Text              | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDI 02, FE)          | Image Data Sequence              | Yes  | Yes     | No      |

| DICOM Tag                     | Attribute Name          | Full | Reduced | Service |
|-------------------------------|-------------------------|------|---------|---------|
| (0021, SIEMENS MR SDR 01, 01) | Creator Identifier      | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDR 01, 02) | Application Identifier  | Yes  | Yes     | No      |
| (0021, SIEMENS MR SDR 01, 03) | Cause Identifier        | Yes  | Yes     | No      |
| (0051, SIEMENS MR HEADER, 0A) | Meas Duration           | Yes  | Yes     | No      |
| (0051, SIEMENS MR HEADER, 0C) | Field Of View           | Yes  | Yes     | No      |
| (0051, SIEMENS MR HEADER, 0D) | Slice Position          | Yes  | Yes     | No      |
| (0051, SIEMENS MR HEADER, 12) | Rel Table Position      | Yes  | Yes     | No      |
| (0051, SIEMENS MR HEADER, 13) | Positive PCS Directions | Yes  | Yes     | No      |
| (7FE1, SIEMENS MR IMA, 10)    | Raw Data                | Yes  | Yes     | No      |

## 8 Security

### 8.1 Security Profiles

#### 8.1.1 Time Synchronization Profiles

Time Synchronization Profiles: syngo® MR acts as an NTP Client in the Maintain Time Transaction.

#### 8.1.2 Basic TLS Secure Transport Connection Profile

Basic TLS Secure Transport Connection Profile supports TLS version 1.0, 1.1 and 1.2 protocols with the following features:

| Supported TLS Feature      | Mechanism              |
|----------------------------|------------------------|
| Entity Authentication      | RSA based certificates |
| Exchange of Master Secrets | RSA                    |
| Data Integrity             | SHA                    |
| Privacy                    | Triple DES EDE, CBC    |

The default secure DICOM port is 2762 (can be reconfigured).

### 8.2 Association Level Security

It is possible to configure whether the SCP will only answer to known AETs or to any AET.

### 8.3 Application Level Security

- User must login with own password
- For configuration and Maintenance, Service Technician must login with a separate password.



## 9 Annexes

### 9.1 IOD Contents

#### 9.1.1 Enhanced MR Image IOD

This chapter describes the DICOM attributes of Enhanced MR Image Instances performed by the MR acquisition.

##### 9.1.1.1 Patient Module

| Attribute Name           | Tag         | Supported Values              |
|--------------------------|-------------|-------------------------------|
| Patient's Name           | (0010,0010) | RIS defined or set by creator |
| Patient ID               | (0010,0020) | RIS defined or set by creator |
| Issuer of Patient ID     | (0010,0021) | RIS defined                   |
| Patient's Birth Date     | (0010,0030) | RIS defined or set by creator |
| Patient's Birth Time     | (0010,0032) | RIS defined or set by creator |
| Patient's Sex            | (0010,0040) | RIS defined or set by creator |
| Other Patient IDs        | (0010,1000) | RIS defined or set by creator |
| Other Patient Names      | (0010,1001) | RIS defined or set by creator |
| Ethnic Group             | (0010,2160) | RIS defined or set by creator |
| Patient Comments         | (0010,4000) | RIS defined or set by creator |
| Other Patient Names      | (0010,1001) | RIS defined or set by creator |
| Patient Identity Removed | (0012,0062) | set by creator                |

##### 9.1.1.2 General Study Module

| Attribute Name             | Tag         | Supported Values              |
|----------------------------|-------------|-------------------------------|
| Study Instance UID         | (0020,000D) | RIS defined or set by creator |
| Study Date                 | (0008,0020) | set by creator                |
| Study Time                 | (0008,0030) | set by creator                |
| Referring Physician's Name | (0008,0090) | RIS defined or set by creator |
| Requesting Physician       | (0032,1032) | RIS defined or set by creator |
| Study ID                   | (0020,0010) | set by creator                |
| Accession Number           | (0008,0050) | RIS defined or set by creator |
| Study Description          | (0008,1030) | set by creator                |
| Procedure Code Sequence    | (0008,1032) | RIS defined                   |
| >Code Value                | (0008,0100) | RIS defined                   |
| >Code Scheme Designator    | (0008,0102) | RIS defined                   |
| >Code Meaning              | (0008,0104) | RIS defined                   |

##### 9.1.1.3 Patient Study Module

| Attribute Name   | Tag         | Supported Values              |
|------------------|-------------|-------------------------------|
| Patient's Age    | (0010,1010) | RIS defined or set by creator |
| Patient's Size   | (0010,1020) | RIS defined or set by creator |
| Patient's Weight | (0010,1030) | RIS defined or set by creator |

#### 9.1.1.4 General Series Module

| Attribute Name                               | Tag         | Supported Values              |
|--|-------------|-------------------------------|
| Modality                                     | (0008,0060) | MR                            |
| Series Instance UID                          | (0020,000E) | set by creator                |
| Series Number                                | (0020,0011) | set by creator                |
| Series Date                                  | (0008,0021) | set by creator                |
| Series Time                                  | (0008,0031) | set by creator                |
| Performing Physicians' Name                  | (0008,1050) | RIS defined or set by creator |
| Protocol Name                                | (0018,1030) | set by creator                |
| Series Description                           | (0008,103E) | set by creator                |
| Referenced Performed Procedure Step Sequence | (0008,1111) | set by creator                |
| >Referenced SOP Class UID                    | (0008,1150) | set by creator                |
| >Referenced SOP Instance UID                 | (0008,1155) | set by creator                |
| Body Part Examined                           | (0018,0015) | set by creator                |
| Patient Position                             | (0018,5100) | set by creator                |
| Request Attributes Sequence                  | (0040,0275) | RIS defined                   |
| >Requested Procedure ID                      | (0040,1001) | RIS defined                   |
| >Accession Number                            | (0008,0050) | RIS defined                   |
| >Study Instance UID                          | (0020,000D) | RIS defined                   |
| >Requested Procedure Description             | (0032,1060) | RIS defined                   |
| >Scheduled Procedure Step ID                 | (0040,0009) | RIS defined                   |
| >Scheduled Procedure Step Description        | (0040,0007) | RIS defined                   |
| > Referenced Study Sequence                  | (0008,1110) | RIS defined                   |
| >> Referenced SOP Class UID                  | (0008,1150) | RIS defined                   |
| >> Referenced SOP Instance UID               | (0008,1155) | RIS defined                   |
| > Requested Procedure Code Sequence          | (0032,1064) | RIS defined                   |
| >> Code Value                                | (0008,0100) | RIS defined                   |
| >> Code Scheme Designator                    | (0008,0102) | RIS defined                   |
| >> Code Meaning                              | (0008,0104) | RIS defined                   |
| > Scheduled Protocol Code Sequence           | (0040,0008) | RIS defined                   |
| >> Code Value                                | (0008,0100) | RIS defined                   |
| >> Code Scheme Designator                    | (0008,0102) | RIS defined                   |
| >> Code Meaning                              | (0008,0104) | RIS defined                   |
| Performed Procedure Step ID                  | (0040,0253) | set by creator                |
| Performed Procedure Step Start Date          | (0040,0244) | set by creator                |
| Performed Procedure Step Start Time          | (0040,0245) | set by creator                |
| Performed Procedure Step Description         | (0040,0254) | set by creator                |

#### 9.1.1.5 MR Series Module

| Attribute Name | Tag         | Supported Values |
|----------------|-------------|------------------|
| Modality       | (0008,0060) | MR               |

#### 9.1.1.6 Frame of Reference Module

| Attribute Name               | Tag         | Supported Values |
|------------------------------|-------------|------------------|
| Frame of Reference UID       | (0020,0052) | set by creator   |
| Position Reference Indicator | (0020,1040) | empty            |

### 9.1.1.7 General Equipment Module

| Attribute Name            | Tag         | Supported Values              |
|---------------------------|-------------|-------------------------------|
| Manufacturer              | (0008,0070) | Siemens                       |
| Institution Name          | (0008,0080) | RIS defined or set by creator |
| Station Name              | (0008,1010) | set by creator                |
| Institution Address       | (0008,0081) | RIS defined set by creator    |
| Manufacturer's Model Name | (0008,1090) | set by creator                |
| Device Serial Number      | (0018,1000) | set by creator                |
| Software Versions         | (0018,1020) | set by creator                |

### 9.1.1.8 Enhanced General Equipment Module

| Attribute Name | Tag         |         |
|----------------|-------------|---------|
| Manufacturer   | (0008,0070) | Siemens |

### 9.1.1.9 Image Pixel Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Samples per Pixel          | (0028,0002) | 1                |
| Photometric Interpretation | (0028,0004) | MONOCHROME2      |
| Rows                       | (0028,0010) | set by creator   |
| Columns                    | (0028,0011) | set by creator   |
| Bits Allocated             | (0028,0100) | 16               |
| Bits Stored                | (0028,0101) | 12 or 16         |
| High Bit                   | (0028,0102) | 11 or 15         |
| Pixel Representation       | (0028,0103) | 0                |
| Pixel Data                 | (7FE0,0010) | set by creator   |
| Smallest Image Pixel Value | (0028,0106) | set by creator   |
| Largest Image Pixel Value  | (0028,0107) | set by creator   |

### 9.1.1.10 Enhanced Contrast/Bolus Module

| Attribute Name                                | Tag         | Supported Values |
|---|-------------|------------------|
| Contrast/Bolus Agent Sequence                 | (0018,0012) | set by creator   |
| >Code Value                                   | (0008,0100) | set by creator   |
| >Code Scheme Designator                       | (0008,0102) | set by creator   |
| >Code Meaning                                 | (0008,0104) | set by creator   |
| >Contrast/Bolus Agent Number                  | (0018,9337) | 1                |
| >Contrast/Bolus Administration Route Sequence | (0018,0014) | set by creator   |
| >>Code Value                                  | (0008,0100) | set by creator   |
| >>Code Scheme Designator                      | (0008,0102) | set by creator   |
| >>Code Meaning                                | (0008,0104) | set by creator   |
| >Contrast/Bolus Ingredient Code Sequence      | (0018,9338) | set by creator   |
| >>Code Value                                  | (0008,0100) | set by creator   |
| >>Code Scheme Designator                      | (0008,0102) | set by creator   |
| >>Code Meaning                                | (0008,0104) | set by creator   |
| >Contrast/Bolus Volume                        | (0018,1041) | set by creator   |
| >Contrast/Bolus Ingredient Concentration      | (0018,1049) | set by creator   |

### 9.1.1.11 Multi-frame Functional Groups Module

| Attribute Name   | Tag         | Supported Values |
|------------------|-------------|------------------|
| Instance Number  | (0020,0013) | set by creator   |
| Content Date     | (0008,0023) | set by creator   |
| Content Time     | (0008,0033) | set by creator   |
| Number of Frames | (0028,0008) | set by creator   |

### 9.1.1.12 Multi-frame Dimension Module

| Attribute Name                  | Tag         | Supported Values |
|---------------------------------|-------------|------------------|
| Dimension Organization Sequence | (0020,9221) | set by creator   |
| Dimension Index Sequence        | (0020,9222) | set by creator   |

### 9.1.1.13 Enhanced MR Image Functional Groups

| Functional Group Macro | Attribute Name                       | Tag         | Supported Values |
|------------------------|--------------------------------------|-------------|------------------|
| Pixel Measures         | Pixel Measures Sequence              | (0028,9110) | set by creator   |
|                        | >Pixel Spacing                       | (0018,0030) | set by creator   |
|                        | >Slice Thickness                     | (0028,0050) | set by creator   |
| Frame Content          | Frame Content Sequence               | (0020,9111) | set by creator   |
|                        | >Frame Acquisition Number            | (0020,9156) | set by creator   |
|                        | >Frame Reference Date Time           | (0018,9151) | set by creator   |
|                        | >Frame Acquisition Date Time         | (0018,9074) | set by creator   |
|                        | >Frame Acquisition Duration          | (0018,9220) | set by creator   |
|                        | >Dimension Index Values              | (0020,9157) | set by creator   |
|                        | >Temporal Position Index             | (0020,9128) | set by creator   |
|                        | >Stack ID                            | (0020,9056) | set by creator   |
|                        | >In-Stack Position Number            | (0020,9057) | set by creator   |
|                        | >Frame Comments                      | (0020,9158) | set by creator   |
| Plane Position         | Plane Position Sequence              | (0020,9113) | set by creator   |
|                        | >Image Position (Patient)            | (0020,0032) | set by creator   |
| Plane Orientation      | Plane Orientation Sequence           | (0020,9116) | set by creator   |
|                        | >Image Orientation (Patient)         | (0020,0037) | set by creator   |
| Referenced Image       | Referenced Image Sequence            | (0008,1140) | set by creator   |
|                        | >Referenced SOP Class UID            | (0008,1150) | set by creator   |
|                        | >Referenced SOP Instance UID         | (0008,1155) | set by creator   |
|                        | >Referenced Frame Number             | (0008,1160) | set by creator   |
|                        | >Purpose of Referenced Code Sequence | (0040,A170) | set by creator   |
|                        | >>Code Value                         | (0008,0100) | e.g. 121311      |
|                        | >>Coding Scheme Designator           | (0008,0102) | e.g. DCM         |
|                        | >>Code Meaning                       | (0008,0104) | e.g. Localizer   |
| Derivation Image       | Derivation Image Sequence            | (0008,9124) | set by creator   |
|                        | >Derivation Code Sequence            | (0008,9215) | set by creator   |
|                        | >>Code Value                         | (0008,0100) | set by creator   |
|                        | >>Coding Scheme Designator           | (0008,0102) | set by creator   |
|                        | >>Code Meaning                       | (0008,0104) | set by creator   |
|                        | >Source Image Sequence               | (0008,2112) | set by creator   |
|                        | >>Referenced SOP Class UID           | (0008,1150) | set by creator   |
|                        | >>Referenced SOP Instance UID        | (0008,1155) | set by creator   |

|                                  |  |  |  |
|----------------------------------|--|--|--|
|                                  | >>Referenced Frame Number<br>>>Purpose of Referenced Code Sequence<br>>>>Code Value<br>>>>Coding Scheme Designator<br>>>>Code Meaning  | (0008,1160)<br>(0040,A170)<br>(0008,0100)<br>(0008,0102)<br>(0008,0104)  | set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator   |
| Cardiac Synchronization          | Cardiac Synchronization Sequence<br>>Nominal Percentage of Cardiac Phase<br>>Nominal Cardiac Trigger Delay Time<br>>Actual Cardiac Trigger Delay Time<br>>Nominal Cardiac Trigger Time Prior to R-peak<br>>Actual Cardiac Trigger Time Prior to R-peak<br>>Intervals Acquired<br>>Intervals Rejected<br>>Heart Rate<br>>R-R Interval Time Nominal<br>>Low R-R Value<br>>High R-R Value | (0018,9118)<br>(0020,9241)<br>(0020,9153)<br>(0020,9252)<br>(0020,9154)<br>(0020,9155)<br>(0018,1083)<br>(0018,1084)<br>(0018,1088)<br>(0020,9251)<br>(0018,1081)<br>(0018,1082) | set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator |
| Frame Anatomy                    | Frame Anatomy Sequence<br>>Frame Laterality<br>>Anatomic Region Sequence<br>>>Code Value<br>>>Coding Scheme Designator<br>>>Code Meaning   | (0020,9071)<br>(0020,9072)<br>(0008,2218)<br>(0008,0100)<br>(0008,0102)<br>(0008,0104)   | set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator   |
| Pixel Value Transformation       | Pixel Value Transformation Sequence<br>>Rescale Intercept<br>>Rescale Slope<br>>Rescale Type   | (0028,9145)<br>(0028,1052)<br>(0028,1053)<br>(0028,1054)   | set by creator<br>set by creator<br>set by creator<br>set by creator   |
| Frame VOI LUT                    | Frame VOI LUT Sequence<br>>Window Center<br>>Window Width  | (0028,9132)<br>(0028,1050)<br>(0028,1051)  | set by creator<br>set by creator<br>set by creator   |
| Real World Value Mapping         | Real World Value Mapping Sequence<br>>Real World Value Intercept<br>>Real World Value Slope<br>>Measurement Units Code Sequence<br>>>Code Value<br>>>Coding Scheme Designator<br>>>Code Meaning  | (0040,9096)<br>(0040,9224)<br>(0040,9225)<br>(0040,08EA)<br>(0008,0100)<br>(0008,0102)<br>(0008,0104)  | set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator   |
| Contrast/Bolus Usage             | Contrast/Bolus Usage Sequence<br>>Contrast/Bolus Agent Number<br>>Contrast/Bolus Agent Administered<br>>Contrast/Bolus Agent Detected<br>>Contrast/Bolus Agent Phase   | (0018,9341)<br>(0018,9337)<br>(0018,9342)<br>(0018,9343)<br>(0018,9344)  | set by creator<br>1<br>set by creator<br>set by creator<br>set by creator  |
| MR Image Frame Type              | MR Image Frame Type Sequence<br>>Frame Type<br>>Pixel Presentation<br>>Volumetric Properties<br>>Volume Based Calculation Technique<br>>Complex Image Component<br>>Acquisition Contrast<br>>Functional Settling Phase Frames Present  | (0018,9226)<br>(0008,9007)<br>(0008,9205)<br>(0008,9206)<br>(0008,9207)<br>(0008,9208)<br>(0008,9209)<br>(0018,9622)   | set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator   |
| MR Timing and Related Parameters | MR Timing and Related Parameters Sequence<br>>Repetition Time<br>>Flip Angle   | (0018,9112)<br>(0018,0080)<br>(0018,1314)  | set by creator<br>set by creator<br>set by creator   |

|                     |   |   |  |
|---------------------|---|---|--|
|                     | >Echo Train Length<br>>RF Echo Train Length<br>>Gradient Echo Train Length<br>>Specific Absorption Rate Sequence <sup>d</sup><br>>>Specific Absorption Rate Definition <sup>d</sup><br>>>Specific Absorption Rate Value <sup>d</sup><br>>Gradient Output Type <sup>d</sup><br>>Gradient Output <sup>d</sup><br>>Operation Mode Sequence <sup>d</sup><br>>>Operating Mode Type <sup>d</sup><br>>>Operating Mode <sup>d</sup> | (0018,0091)<br>(0018,9240)<br>(0018,9241)<br>(0018,9239)<br>(0018,9179)<br>(0018,9181)<br>(0018,9180)<br>(0018,9182)<br>(0018,9176)<br>(0018,9177)<br>(0018,9178)   | set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator   |
| MR FOV/Geometry     | MR FOV/Geometry Sequence<br>>In-plane Phase Encoding Direction<br>>MR Acquisition Frequency Encoding Steps<br>>MR Acquisition Phase Encoding Steps in-plane<br>>MR Acquisition Phase Encoding Steps out-of-plane<br>>Percent Sampling<br>>Percent Phase Field of View   | (0018,9125)<br>(0018,1312)<br>(0018,9058)<br>(0018,9231)<br>(0018,9232)<br>(0018,0093)<br>(0018,0094)   | set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator   |
| MR Echo             | MR Echo Sequence<br>>Effective Echo Time  | (0018,9114)<br>(0018,9082)  | set by creator<br>set by creator   |
| MR Modifier         | MR Modifier Sequence<br>>Inversion Recovery<br>>Inversion Times<br>>Flow Compensation<br>>Flow Compensation Direction<br>>Spoiling<br>>T2 Preparation<br>>Spectrally Selected Excitation<br>>Spatial Pre-saturation<br>>Partial Fourier<br>>Partial Fourier Direction<br>>Parallel Acquisition<br>>Parallel Acquisition Technique<br>>Parallel Reduction Factor In-plane<br>>Parallel Reduction Factor out-of-plane         | (0018,9115)<br>(0018,9009)<br>(0018,9079)<br>(0018,9010)<br>(0018,9183)<br>(0018,9016)<br>(0018,9021)<br>(0018,9026)<br>(0018,9027)<br>(0018,9081)<br>(0018,9036)<br>(0018,9077)<br>(0018,9078)<br>(0018,9069)<br>(0018,9155) | set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator |
| MR Imaging Modifier | MR Imaging Modifier Sequence<br>>Magnetization Transfer<br>>Blood Signal Nulling<br>>Tagging<br>>Tag Spacing First Dimension<br>>Tag Spacing Second Dimension<br>>Tag Angle First Axis<br>>Tag Angle Second Axis<br>>Tag Thickness<br>>Tagging Delay<br>>Transmitter Frequency<br>>Pixel Bandwidth  | (0018,9006)<br>(0018,9020)<br>(0018,9022)<br>(0018,9028)<br>(0018,9030)<br>(0018,9218)<br>(0018,9019)<br>(0018,9219)<br>(0018,9035)<br>(0018,9184)<br>(0018,9098)<br>(0018,0095)  | set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator<br>set by creator   |
| MR Receive Coil     | MR Receive Coil Sequence<br>>Receive Coil Name  | (0018,9042)<br>(0018,1250)  | set by creator<br>set by creator   |

<sup>d</sup> Specific Absorption Rate Sequence, Gradient Output Type, Gradient Output, and Operating Mode Sequence are not set for GRASP reconstructed images due to technical reasons. The values can be found in the related preview images of the GRASP acquisition.

|                           |  |             |                |
|---------------------------|--|-------------|----------------|
|                           | >Receive Coil Manufacturer Name        | (0018,9041) | set by creator |
|                           | >Receive Coil Type                     | (0018,9043) | set by creator |
|                           | >Quadrature Receive Coil               | (0018,9044) | set by creator |
|                           | >Multi-Coil Definition Sequence        | (0018,9045) | set by creator |
|                           | >>Multi-Coil Element Name              | (0018,9047) | set by creator |
|                           | >>Multi-Coil Element Used              | (0018,9048) | set by creator |
| MR Transmit Coil          | MR Transmit Coil Sequence              | (0018,9049) | set by creator |
|                           | >Transmit Coil Name                    | (0018,1251) | set by creator |
|                           | >Transmit Coil Manufacturer Name       | (0018,9050) | set by creator |
|                           | >Transmit Coil Type                    | (0018,9051) | set by creator |
| MR Diffusion              | MR Diffusion Sequence                  | (0018,9117) | set by creator |
|                           | >Diffusion b-value                     | (0018,9087) | set by creator |
|                           | >Diffusion Directionality              | (0018,9075) | set by creator |
|                           | >Diffusion Gradient Direction Sequence | (0018,9076) | set by creator |
|                           | >>Diffusion Gradient Orientation       | (0018,9089) | set by creator |
|                           | >Diffusion b-matrix Sequence           | (0018,9601) | set by creator |
|                           | >>Diffusion b-value XX                 | (0018,9602) | set by creator |
|                           | >>Diffusion b-value XY                 | (0018,9603) | set by creator |
|                           | >>Diffusion b-value XZ                 | (0018,9604) | set by creator |
|                           | >>Diffusion b-value YY                 | (0018,9605) | set by creator |
|                           | >>Diffusion b-value YZ                 | (0018,9606) | set by creator |
|                           | >>Diffusion b-value ZZ                 | (0018,9607) | set by creator |
| MR Averages               | MR Averages Sequence                   | (0018,9119) | set by creator |
|                           | >Number of Averages                    | (0018,0083) | set by creator |
| MR Velocity Encoding      | MR Velocity Encoding Sequence          | (0018,9197) | set by creator |
|                           | >Velocity Encoding Direction           | (0018,9090) | set by creator |
|                           | >Velocity Encoding Minimum Value       | (0018,9091) | set by creator |
|                           | >Velocity Encoding Maximum Value       | (0018,9217) | set by creator |
| MR Arterial Spin Labeling | MR Arterial Spin Labeling Sequence     | (0018,9251) | set by creator |
|                           | >ASL Technique Description             | (0018,9252) | set by creator |
|                           | >ASL Context                           | (0018,9257) | set by creator |
|                           | >ASL Slab Sequence                     | (0018,9260) | set by creator |
|                           | >>ASL Slab Number                      | (0018,9253) | set by creator |
|                           | >>ASL Slab Thickness                   | (0018,9254) | set by creator |
|                           | >>ASL Slab Orientation                 | (0018,9255) | set by creator |
|                           | >>ASL Mid Slab Position                | (0018,9256) | set by creator |
|                           | >>ASL Pulse Train Duration             | (0018,9258) | set by creator |
|                           | >ASL Crusher Flag                      | (0018,9259) | set by creator |
|                           | >ASL Crusher Flow Limit                | (0018,925A) | set by creator |
|                           | >ASL Crusher Description               | (0018,925B) | set by creator |
|                           | >ASLBolus Cut-off Flag                 | (0018,925C) | set by creator |
|                           | >ASLBolus Cut-off Timing Sequence      | (0018,925D) | set by creator |
|                           | >>ASL Bolus Cut-off Delay Time         | (0018,925F) | set by creator |
|                           | >>ASL Bolus Cut-off Technique          | (0018,925E) | set by creator |

#### 9.1.1.14 Cardiac Synchronization Module

| Attribute Name                    | Tag         | Supported Values |
|-----------------------------------|-------------|------------------|
| Cardiac Synchronization Technique | (0018,9037) | set by creator   |
| Cardiac Signal Source             | (0018,9085) | set by creator   |
| Cardiac RR Interval Specified     | (0018,9070) | set by creator   |
| Low R-R Value                     | (0018,1081) | set by creator   |

|                    |             |                |
|--------------------|-------------|----------------|
| High R-R Value     | (0018,1082) | set by creator |
| Intervals Acquired | (0018,1083) | set by creator |
| Intervals Rejected | (0018,1084) | set by creator |

#### 9.1.1.15 Respiratory Synchronization Module

| Attribute Name                            | Tag         | Supported Values |
|---|-------------|------------------|
| Respiratory Motion Compensation Technique | (0018,9170) | set by creator   |

#### 9.1.1.16 Bulk Motion Synchronization Module

| Attribute Name                     | Tag         | Supported Values |
|------------------------------------|-------------|------------------|
| Bulk Motion Compensation Technique | (0018,9172) | set by creator   |
| Bulk Motion Signal Source          | (0018,9173) | set by creator   |

#### 9.1.1.17 Supplemental Palette Color Lookup Table Module

| Attribute Name                              | Tag         | Supported Values |
|---|-------------|------------------|
| Red Palette Color Lookup Table Descriptor   | (0028,1101) | set by creator   |
| Green Palette Color Lookup Table Descriptor | (0028,1102) | set by creator   |
| Blue Palette Color Lookup Table Descriptor  | (0028,1103) | set by creator   |
| Red Palette Color Lookup Table Data         | (0028,1201) | set by creator   |
| Green Palette Color Lookup Table Data       | (0028,1202) | set by creator   |
| Blue Palette Color Lookup Table Data        | (0028,1203) | set by creator   |

#### 9.1.1.18 Acquisition Context Module

| Attribute Name               | Tag         | Supported Values |
|------------------------------|-------------|------------------|
| Acquisition Context Sequence | (0040,0555) | TID 99_3100      |

##### 9.1.1.18.1 MR Acquisition Context TID 99\_3100

The MR Acquisition Context TID 99\_3100 is used in the Acquisition Context Sequence (0040,0555) of Enhanced MR Images.

Type: Extensible

|   | NL | VT   | Concept Name                          | VM | Req<br>Type | Condition | Value Set Constraint |
|---|----|------|---------------------------------------|----|-------------|-----------|----------------------|
| 1 |    | CODE | EV (A-52, 99SMS_CTMR, "MR Technique") | 1  | M           |           | DCID(A-200)          |

##### 9.1.1.18.2 MR Technique Context Group A-200

Type: Extensible

| Coding Scheme Designator (0008,0102) | Coding Scheme Version (0008,0103) | Code Value (0008,0100) | Code Meaning (0008,0104) |
|--------------------------------------|-----------------------------------|------------------------|--------------------------|
| 99SMS_CTMR                           | 1.0                               | A-200                  | GRASP                    |



|            |     |       |           |
|------------|-----|-------|-----------|
| 99SMS_CTMR | 1.0 | A-202 | FastView  |
| 99SMS_CTMR | 1.0 | A-204 | Angio     |
| 99SMS_CTMR | 1.0 | A-206 | BOLD      |
| 99SMS_CTMR | 1.0 | A-208 | Diffusion |
| 99SMS_CTMR | 1.0 | A-210 | DTI       |

### 9.1.1.19 Enhanced MR Image Module

| Attribute Name                     | Tag         | Supported Values |
|------------------------------------|-------------|------------------|
| Acquisition Number                 | (0020,0012) | set by creator   |
| Acquisition Date Time              | (0008,002A) | set by creator   |
| Acquisition Duration               | (0018,9073) | set by creator   |
| Content Qualification              | (0018,9004) | set by creator   |
| Resonant Nucleus                   | (0018,9100) | set by creator   |
| k-space Filtering                  | (0018,9064) | set by creator   |
| Magnetic Field Strength            | (0018,0087) | set by creator   |
| Applicable Safety Standard Agency  | (0018,9174) | set by creator   |
| Image Comments                     | (0020,4000) | set by creator   |
| Image Type                         | (0008,0008) | set by creator   |
| Pixel Presentation                 | (0008,9205) | set by creator   |
| Volumetric Properties              | (0008,9206) | set by creator   |
| Volume Based Calculation Technique | (0008,9207) | set by creator   |
| Complex Image Component            | (0008,9208) | set by creator   |
| Acquisition Contrast               | (0008,9209) | set by creator   |
| Samples per Pixel                  | (0028,0002) | 1                |
| Photometric Interpretation         | (0028,0004) | MONOCHROME2      |
| Bits Allocated                     | (0028,0100) | 16               |
| Bits Stored                        | (0028,0101) | 12 or 16         |
| High Bit                           | (0028,0102) | 11 or 15         |
| Pixel Representation               | (0028,0103) | 0                |
| Burned In Annotation               | (0028,0301) | NO               |
| Lossy Image Compression            | (0028,2110) | 00               |
| Presentation LUT Shape             | (2050,0020) | IDENTITY         |

### 9.1.1.20 MR Pulse Sequence Module

| Attribute Name                         | Tag         | Supported Values |
|--|-------------|------------------|
| Pulse Sequence Name                    | (0018,9005) | set by creator   |
| MR Acquisition Type                    | (0018,0023) | set by creator   |
| Echo Pulse Sequence                    | (0018,9008) | set by creator   |
| Multiple Spin Echo                     | (0018,9011) | set by creator   |
| Multi-planar Excitation                | (0018,9012) | set by creator   |
| Phase Contrast                         | (0018,9014) | set by creator   |
| Velocity Encoding Acquisition Sequence | (0018,9092) | set by creator   |
| > Velocity Encoding Direction          | (0018,9090) | set by creator   |
| Time of Flight Contrast                | (0018,9015) | set by creator   |
| Arterial Spin Labeling Contrast        | (0018,9250) | set by creator   |
| Steady State Pulse Sequence            | (0018,9017) | set by creator   |
| Echo Planar Pulse Sequence             | (0018,9018) | set by creator   |
| Saturation Recovery                    | (0018,9024) | set by creator   |

|                                     |             |                |
|-------------------------------------|-------------|----------------|
| Spectral Selected Suppression       | (0018,9025) | set by creator |
| Oversampling Phase                  | (0018,9029) | set by creator |
| Geometry of k-Space Traversal       | (0018,9032) | set by creator |
| Rectilinear Phase Encode Reordering | (0018,9034) | set by creator |
| Segmented k-Space Traversal         | (0018,9033) | set by creator |
| Coverage of k-Space                 | (0018,9094) | set by creator |
| Number of k-Space Trajectories      | (0018,9093) | set by creator |

### 9.1.1.21 SOP Common Module

| Attribute Name         | Tag         | Supported Values                  |
|------------------------|-------------|-----------------------------------|
| SOP Class UID          | (0008,0016) | Enhanced MR Storage SOP Class UID |
| SOP Instance UID       | (0008,0018) | set by creator                    |
| Specific Character Set | (0008,0005) | set by creator                    |
| Instance Creation Date | (0008,0012) | date the SOP instance was created |
| Instance Creation Time | (0008,0013) | time the SOP instance was created |

## 9.1.2 Enhanced MR Color Image IOD

The Enhanced MR Color Image IOD supports the same DICOM attributes as the Enhanced MR Image IOD with the exception of the Photometric Interpretation which is RGB. Also, Supplemental Palette Color Lookup Table module, Pixel Value Transformation Macro, Frame VOI LUT Macro, and Real World Value Mapping Macro are not part of this IOD.

### 9.1.2.1 Image Pixel Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Samples per Pixel          | (0028,0002) | 3                |
| Photometric Interpretation | (0028,0004) | RGB              |
| Rows                       | (0028,0010) | set by creator   |
| Columns                    | (0028,0011) | set by creator   |
| Bits Allocated             | (0028,0100) | 8                |
| Bits Stored                | (0028,0101) | 8                |
| High Bit                   | (0028,0102) | 7                |
| Pixel Representation       | (0028,0103) | 0                |

### 9.1.2.2 Enhanced MR Image Module

| Attribute Name     | Tag         | Supported Values |
|--------------------|-------------|------------------|
| Pixel Presentation | (0008,9205) | TRUE COLOR       |

## 9.1.3 MR Spectroscopy IOD

This chapter describes the DICOM attributes of MR Spectroscopy Instances performed by the MR acquisition.

### 9.1.3.1 Patient Module

| Attribute Name | Tag | Supported Values |
|----------------|-----|------------------|
|----------------|-----|------------------|

|                          |             |                               |
|--------------------------|-------------|-------------------------------|
| Patient's Name           | (0010,0010) | RIS defined or set by creator |
| Patient ID               | (0010,0020) | RIS defined or set by creator |
| Issuer of Patient ID     | (0010,0021) | RIS defined                   |
| Patient's Birth Date     | (0010,0030) | RIS defined or set by creator |
| Patient's Birth Time     | (0010,0032) | RIS defined or set by creator |
| Patient's Sex            | (0010,0040) | RIS defined or set by creator |
| Other Patient IDs        | (0010,1000) | RIS defined or set by creator |
| Other Patient Names      | (0010,1001) | RIS defined or set by creator |
| Ethnic Group             | (0010,2160) | RIS defined or set by creator |
| Patient Comments         | (0010,4000) | RIS defined or set by creator |
| Other Patient Names      | (0010,1001) | RIS defined or set by creator |
| Patient Identity Removed | (0012,0062) | set by creator                |

### 9.1.3.2 General Study Module

| Attribute Name             | Tag         | Supported Values              |
|----------------------------|-------------|-------------------------------|
| Study Instance UID         | (0020,000D) | RIS defined or set by creator |
| Study Date                 | (0008,0020) | set by creator                |
| Study Time                 | (0008,0030) | set by creator                |
| Referring Physician's Name | (0008,0090) | RIS defined or set by creator |
| Study ID                   | (0020,0010) | set by creator                |
| Accession Number           | (0008,0050) | RIS defined or set by creator |
| Study Description          | (0008,1030) | set by creator                |
| Procedure Code Sequence    | (0008,1032) | RIS defined                   |
| >Code Value                | (0008,0100) | RIS defined                   |
| >Code Scheme Designator    | (0008,0102) | RIS defined                   |
| >Code Meaning              | (0008,0104) | RIS defined                   |

### 9.1.3.3 Patient Study Module

| Attribute Name   | Tag         | Supported Values              |
|------------------|-------------|-------------------------------|
| Patient's Age    | (0010,1010) | RIS defined or set by creator |
| Patient's Size   | (0010,1020) | RIS defined or set by creator |
| Patient's Weight | (0010,1030) | RIS defined or set by creator |

### 9.1.3.4 General Series Module

| Attribute Name                               | Tag         | Supported Values              |
|--|-------------|-------------------------------|
| Modality                                     | (0008,0060) | MR                            |
| Series Instance UID                          | (0020,000E) | set by creator                |
| Series Number                                | (0020,0011) | set by creator                |
| Series Date                                  | (0008,0021) | set by creator                |
| Series Time                                  | (0008,0031) | set by creator                |
| Performing Physicians' Name                  | (0008,1050) | RIS defined or set by creator |
| Protocol Name                                | (0018,1030) | set by creator                |
| Series Description                           | (0008,103E) | set by creator                |
| Referenced Performed Procedure Step Sequence | (0008,1111) | set by creator                |
| >Referenced SOP Class UID                    | (0008,1150) | set by creator                |
| >Referenced SOP Instance UID                 | (0008,1155) | set by creator                |
| Body Part Examined                           | (0018,0015) | set by creator                |

|                                       |             |                |
|---------------------------------------|-------------|----------------|
| Patient Position                      | (0018,5100) | set by creator |
| Request Attributes Sequence           | (0040,0275) | RIS defined    |
| >Requested Procedure ID               | (0040,1001) | RIS defined    |
| >Accession Number                     | (0008,0050) | RIS defined    |
| >Study Instance UID                   | (0020,000D) | RIS defined    |
| >Requested Procedure Description      | (0032,1060) | RIS defined    |
| >Scheduled Procedure Step ID          | (0040,0009) | RIS defined    |
| >Scheduled Procedure Step Description | (0040,0007) | RIS defined    |
| > Referenced Study Sequence           | (0008,1110) | RIS defined    |
| >> Referenced SOP Class UID           | (0008,1150) | RIS defined    |
| >> Referenced SOP Instance UID        | (0008,1155) | RIS defined    |
| > Requested Procedure Code Sequence   | (0032,1064) | RIS defined    |
| >> Code Value                         | (0008,0100) | RIS defined    |
| >> Code Scheme Designator             | (0008,0102) | RIS defined    |
| >> Code Meaning                       | (0008,0104) | RIS defined    |
| > Scheduled Protocol Code Sequence    | (0040,0008) | RIS defined    |
| >> Code Value                         | (0008,0100) | RIS defined    |
| >> Code Scheme Designator             | (0008,0102) | RIS defined    |
| >> Code Meaning                       | (0008,0104) | RIS defined    |
| Performed Procedure Step ID           | (0040,0253) | set by creator |
| Performed Procedure Step Start Date   | (0040,0244) | set by creator |
| Performed Procedure Step Start Time   | (0040,0245) | set by creator |
| Performed Procedure Step Description  | (0040,0254) | set by creator |

### 9.1.3.5 MR Series Module

| Attribute Name | Tag         | Supported Values |
|----------------|-------------|------------------|
| Modality       | (0008,0060) | MR               |

### 9.1.3.6 Frame of Reference Module

| Attribute Name               | Tag         | Supported Values |
|------------------------------|-------------|------------------|
| Frame of Reference UID       | (0020,0052) | set by creator   |
| Position Reference Indicator | (0020,1040) | empty            |

### 9.1.3.7 General Equipment Module

| Attribute Name            | Tag         | Supported Values |
|---------------------------|-------------|------------------|
| Manufacturer              | (0008,0070) | Siemens          |
| Institution Name          | (0008,0080) | set by creator   |
| Institution Address       | (0008,0081) | set by creator   |
| Manufacturer's Model Name | (0008,1090) | set by creator   |
| Device Serial Number      | (0018,1000) | set by creator   |
| Software Versions         | (0018,1020) | set by creator   |

### 9.1.3.8 Enhanced General Equipment Module

| Attribute Name | Tag         | Supported Values |
|----------------|-------------|------------------|
| Manufacturer   | (0008,0070) | Siemens          |

### 9.1.3.9 Enhanced Contrast/Bolus Module

| Attribute Name                                | Tag         | Supported Values |
|---|-------------|------------------|
| Contrast/Bolus Agent Sequence                 | (0018,0012) | set by creator   |
| >Code Value                                   | (0008,0100) | set by creator   |
| >Code Scheme Designator                       | (0008,0102) | set by creator   |
| >Code Meaning                                 | (0008,0104) | set by creator   |
| >Contrast/Bolus Agent Number                  | (0018,9337) | 1                |
| >Contrast/Bolus Administration Route Sequence | (0018,0014) | set by creator   |
| >>Code Value                                  | (0008,0100) | set by creator   |
| >>Code Scheme Designator                      | (0008,0102) | set by creator   |
| >>Code Meaning                                | (0008,0104) | set by creator   |
| >Contrast/Bolus Ingredient Code Sequence      | (0018,9338) | set by creator   |
| >>Code Value                                  | (0008,0100) | set by creator   |
| >>Code Scheme Designator                      | (0008,0102) | set by creator   |
| >>Code Meaning                                | (0008,0104) | set by creator   |
| >Contrast/Bolus Volume                        | (0018,1041) | set by creator   |
| >Contrast/Bolus Ingredient Concentration      | (0018,1049) | set by creator   |

### 9.1.3.10 Multi-frame Functional Groups Module

| Attribute Name  | Tag         | Supported Values |
|-----------------|-------------|------------------|
| Instance Number | (0020,0013) | set by creator   |
| Content Date    | (0008,0023) | set by creator   |
| Content Time    | (0008,0033) | set by creator   |

### 9.1.3.11 Multi-frame Dimension Module

| Attribute Name                  | Tag         | Supported Values |
|---------------------------------|-------------|------------------|
| Dimension Organization Sequence | (0020,9221) | set by creator   |
| Dimension Index Sequence        | (0020,9222) | set by creator   |

### 9.1.3.12 MR Spectroscopy Functional Groups

| Functional Group<br>Macro | Attribute Name               | Tag         | Supported Values |
|---------------------------|------------------------------|-------------|------------------|
| Pixel Measures            | Pixel Measures Sequence      | (0028,9110) | set by creator   |
|                           | >Pixel Spacing               | (0018,0030) | set by creator   |
|                           | >Slice Thickness             | (0028,0050) | set by creator   |
| Frame Content             | Frame Content Sequence       | (0020,9111) | set by creator   |
|                           | >Frame Acquisition Number    | (0020,9156) | set by creator   |
|                           | >Frame Reference Date Time   | (0018,9151) | set by creator   |
|                           | >Frame Acquisition Date Time | (0018,9074) | set by creator   |
|                           | >Frame Acquisition Duration  | (0018,9220) | set by creator   |
|                           | >Dimension Index Values      | (0020,9157) | set by creator   |
|                           | >Temporal Position Index     | (0020,9128) | set by creator   |
|                           | >Stack ID                    | (0020,9056) | set by creator   |
|                           | >In-Stack Position Number    | (0020,9057) | set by creator   |
| Plane Position            | >Frame Comments              | (0020,9158) | set by creator   |
|                           | Plane Position Sequence      | (0020,9113) | set by creator   |
|                           | >Image Position (Patient)    | (0020,0032) | set by creator   |

|                                  |   |             |                |
|----------------------------------|---|-------------|----------------|
| Plane Orientation                | Plane Orientation Sequence                | (0020,9116) | set by creator |
|                                  | >Image Orientation (Patient)              | (0020,0037) | set by creator |
| Referenced Image                 | Referenced Image Sequence                 | (0008,1140) | set by creator |
|                                  | >Referenced SOP Class UID                 | (0008,1150) | set by creator |
|                                  | >Referenced SOP Instance UID              | (0008,1155) | set by creator |
|                                  | >Referenced Frame Number                  | (0008,1160) | set by creator |
|                                  | >Purpose of Referenced Code Sequence      | (0040,A170) | set by creator |
|                                  | >>Code Value                              | (0008,0100) | set by creator |
|                                  | >>Coding Scheme Designator                | (0008,0102) | set by creator |
|                                  | >>Code Meaning                            | (0008,0104) | set by creator |
| Derivation Image                 | Derivation Image Sequence                 | (0008,9124) | set by creator |
|                                  | >Derivation Code Sequence                 | (0008,9215) | set by creator |
|                                  | >>Code Value                              | (0008,0100) | set by creator |
|                                  | >>Coding Scheme Designator                | (0008,0102) | set by creator |
|                                  | >>Code Meaning                            | (0008,0104) | set by creator |
|                                  | >Source Image Sequence                    | (0008,2112) | set by creator |
|                                  | >>Referenced SOP Class UID                | (0008,1150) | set by creator |
|                                  | >>Referenced SOP Instance UID             | (0008,1155) | set by creator |
|                                  | >>Referenced Frame Number                 | (0008,1160) | set by creator |
|                                  | >>Purpose of Referenced Code Sequence     | (0040,A170) | set by creator |
|                                  | >>>Code Value                             | (0008,0100) | set by creator |
|                                  | >>>Coding Scheme Designator               | (0008,0102) | set by creator |
|                                  | >>>Code Meaning                           | (0008,0104) | set by creator |
| Frame Anatomy                    | Frame Anatomy Sequence                    | (0020,9071) | set by creator |
|                                  | >Frame Laterality                         | (0020,9072) | set by creator |
|                                  | >Anatomic Region Sequence                 | (0008,2218) | set by creator |
|                                  | >>Code Value                              | (0008,0100) | set by creator |
|                                  | >>Coding Scheme Designator                | (0008,0102) | set by creator |
|                                  | >>Code Meaning                            | (0008,0104) | set by creator |
| Contrast/Bolus Usage             | Contrast/Bolus Usage Sequence             | (0018,9341) | set by creator |
|                                  | >Contrast/Bolus Agent Number              | (0018,9337) | 1              |
|                                  | >Contrast/Bolus Agent Administered        | (0018,9342) | set by creator |
|                                  | >Contrast/Bolus Agent Detected            | (0018,9343) | set by creator |
|                                  | >Contrast/Bolus Agent Phase               | (0018,9344) | set by creator |
| MR Spectroscopy Frame Type       | MR Spectroscopy Frame Type Sequence       | (0018,9227) | set by creator |
|                                  | >Frame Type                               | (0008,9007) | set by creator |
|                                  | >Volumetric Properties                    | (0008,9206) | set by creator |
|                                  | >Volume Based Calculation Technique       | (0008,9207) | set by creator |
|                                  | >Complex Image Component                  | (0008,9208) | set by creator |
|                                  | >Acquisition Contrast                     | (0008,9209) | set by creator |
| MR Timing and Related Parameters | MR Timing and Related Parameters Sequence | (0018,9112) | set by creator |
|                                  | >Repetition Time                          | (0018,0080) | set by creator |
|                                  | >Flip Angle                               | (0018,1314) | set by creator |
|                                  | >Echo Train Length                        | (0018,0091) | set by creator |
|                                  | >RF Echo Train Length                     | (0018,9240) | set by creator |
|                                  | >Gradient Echo Train Length               | (0018,9241) | set by creator |
|                                  | >Specific Absorption Rate Sequence        | (0018,9239) | set by creator |
|                                  | >>Specific Absorption Rate Definition     | (0018,9179) | set by creator |
|                                  | >>Specific Absorption Rate Value          | (0018,9181) | set by creator |
|                                  | >Gradient Output Type                     | (0018,9180) | set by creator |
|                                  | >Gradient Output                          | (0018,9182) | set by creator |
|                                  | >Operation Mode Sequence                  | (0018,9176) | set by creator |
|                                  | >>Operating Mode Type                     | (0018,9177) | set by creator |

|                              |  |             |                |
|------------------------------|--|-------------|----------------|
|                              | >>Operating Mode                                   | (0018,9178) | set by creator |
| MR Spectroscopy FOV/Geometry | MR Spectroscopy FOV/Geometry Sequence              | (0018,9103) | set by creator |
|                              | >Spectroscopy Acquisition Data Columns             | (0018,9127) | set by creator |
|                              | >Spectroscopy Acquisition Phase Rows               | (0018,9095) | set by creator |
|                              | >Spectroscopy Acquisition Phase Columns            | (0018,9234) | set by creator |
|                              | >Spectroscopy Acquisition Out-of-Plane Phase Steps | (0018,9159) | set by creator |
|                              | >Percent Sampling                                  | (0018,0093) | set by creator |
|                              | >Percent Phase Field of View                       | (0018,0094) | set by creator |
| MR Echo                      | MR Echo Sequence                                   | (0018,9114) | set by creator |
|                              | >Effective Echo Time                               | (0018,9082) | set by creator |
| MR Modifier                  | MR Modifier Sequence                               | (0018,9115) | set by creator |
|                              | >Inversion Recovery                                | (0018,9009) | set by creator |
|                              | >Inversion Times                                   | (0018,9079) | set by creator |
|                              | >Flow Compensation                                 | (0018,9010) | set by creator |
|                              | >Flow Compensation Direction                       | (0018,9183) | set by creator |
|                              | >Spoiling  | (0018,9016) | set by creator |
|                              | >T2 Preparation                                    | (0018,9021) | set by creator |
|                              | >Spectrally Selected Excitation                    | (0018,9026) | set by creator |
|                              | >Spatial Pre-saturation                            | (0018,9027) | set by creator |
|                              | >Partial Fourier                                   | (0018,9081) | set by creator |
|                              | >Partial Fourier Direction                         | (0018,9036) | set by creator |
|                              | >Parallel Acquisition                              | (0018,9077) | set by creator |
|                              | >Parallel Acquisition Technique                    | (0018,9078) | set by creator |
|                              | >Parallel Reduction Factor In-plane                | (0018,9069) | set by creator |
|                              | >Parallel Reduction Factor out-of-plane            | (0018,9155) | set by creator |
| MR Receive Coil              | MR Receive Coil Sequence                           | (0018,9042) | set by creator |
|                              | >Receive Coil Name                                 | (0018,1250) | set by creator |
|                              | >Receive Coil Manufacturer Name                    | (0018,9041) | set by creator |
|                              | >Receive Coil Type                                 | (0018,9043) | set by creator |
|                              | >Quadrature Receive Coil                           | (0018,9044) | set by creator |
|                              | >Multi-Coil Definition Sequence                    | (0018,9045) | set by creator |
|                              | >>Multi-Coil Element Name                          | (0018,9047) | set by creator |
|                              | >>Multi-Coil Element Used                          | (0018,9048) | set by creator |
| MR Transmit Coil             | MR Transmit Coil Sequence                          | (0018,9049) | set by creator |
|                              | >Transmit Coil Name                                | (0018,1251) | set by creator |
|                              | >Transmit Coil Manufacturer Name                   | (0018,9050) | set by creator |
|                              | >Transmit Coil Type                                | (0018,9051) | set by creator |
| MR Averages                  | MR Averages Sequence                               | (0018,9119) | set by creator |
|                              | >Number of Averages                                | (0018,0083) | set by creator |

### 9.1.3.13 Cardiac Synchronization Module

| Attribute Name                    | Tag         | Supported Values |
|-----------------------------------|-------------|------------------|
| Cardiac Synchronization Technique | (0018,9037) | set by creator   |
| Cardiac Signal Source             | (0018,9085) | set by creator   |
| Cardiac RR Interval Specified     | (0018,9070) | set by creator   |
| Low R-R Value                     | (0018,1081) | set by creator   |
| High R-R Value                    | (0018,1082) | set by creator   |
| Intervals Acquired                | (0018,1083) | set by creator   |
| Intervals Rejected                | (0018,1084) | set by creator   |

### 9.1.3.14 Respiratory Synchronization Module

| Attribute Name                            | Tag         | Supported Values |
|---|-------------|------------------|
| Respiratory Motion Compensation Technique | (0018,9170) | set by creator   |

### 9.1.3.15 Bulk Motion Synchronization Module

| Attribute Name                     | Tag         | Supported Values |
|------------------------------------|-------------|------------------|
| Bulk Motion Compensation Technique | (0018,9172) | set by creator   |
| Bulk Motion Signal Source          | (0018,9173) | set by creator   |

### 9.1.3.16 Acquisition Context Module

| Attribute Name               | Tag         | Supported Values |
|------------------------------|-------------|------------------|
| Acquisition Context Sequence | (0040,0555) | empty            |

### 9.1.3.17 MR Spectroscopy Module

| Attribute Name                       | Tag         | Supported Values |
|--------------------------------------|-------------|------------------|
| Acquisition Number                   | (0020,0012) | set by creator   |
| Acquisition Date Time                | (0018,002A) | set by creator   |
| Acquisition Duration                 | (0018,9073) | set by creator   |
| Referenced Raw Data Sequence         | (0008,9121) | set by creator   |
| Content Qualification                | (0018,9004) | PRODUCT          |
| Resonant Nucleus                     | (0018,9100) | set by creator   |
| k-space Filtering                    | (0018,9064) | set by creator   |
| Magnetic Field Strength              | (0018,0087) | set by creator   |
| Applicable Safety Standard Agency    | (0018,9174) | set by creator   |
| Image Comments                       | (0020,4000) | set by creator   |
| Image Type                           | (0008,0008) | set by creator   |
| Volumetric Properties                | (0008,9206) | set by creator   |
| Volume Based Calculation Technique   | (0008,9207) | set by creator   |
| Complex Image Component              | (0008,9208) | set by creator   |
| Acquisition Contrast                 | (0008,9209) | set by creator   |
| Transmitter Frequency                | (0018,9098) | set by creator   |
| Spectral Width                       | (0018,9052) | set by creator   |
| Chemical Shift Reference             | (0018,9053) | set by creator   |
| Volume Localisation Technique        | (0018,9054) | set by creator   |
| Volume Localization Sequence         | (0018,9126) | set by creator   |
| >Slab Thickness                      | (0018,9104) | set by creator   |
| >Slab Orientation                    | (0018,9105) | set by creator   |
| >Mid SlabPosition                    | (0018,9106) | set by creator   |
| De-coupling                          | (0018,9059) | set by creator   |
| De-coupling Nucleus                  | (0018,9060) | set by creator   |
| De-coupling Frequency                | (0018,9061) | set by creator   |
| De-coupling Method                   | (0018,9062) | set by creator   |
| De-coupling Chemical Shift Reference | (0018,9063) | set by creator   |
| Time Domain Filtering                | (0018,9065) | set by creator   |
| Number Of Zero Fills                 | (0018,9066) | set by creator   |
| Baseline Correction                  | (0018,9067) | set by creator   |



|                                     |             |                |
|-------------------------------------|-------------|----------------|
| Frequency Correction                | (0018,9101) | set by creator |
| First Order Phase Correction        | (0018,9198) | set by creator |
| Water Referenced Phase Correction   | (0018,9199) | set by creator |
| Water Reference Acquisition         | (0018,9297) | set by creator |
| Referenced Instance Sequence        | (0008,114A) | set by creator |
| >Purpose of Reference Code Sequence | (0040,A170) | set by creator |

### 9.1.3.18 MR Spectroscopy Pulse Sequence Module

| Attribute Name                      | Tag         | Supported Values |
|-------------------------------------|-------------|------------------|
| Pulse Sequence Name                 | (0018,9005) | set by creator   |
| MR Spectroscopy Acquisition Type    | (0018,9200) | set by creator   |
| Echo Pulse Sequence                 | (0018,9008) | set by creator   |
| Multi Spin Echo                     | (0018,9011) | set by creator   |
| Multi-planar Excitation             | (0018,9012) | set by creator   |
| Steady State Pulse Sequence         | (0018,9017) | set by creator   |
| Echo Planar Pulse Sequence          | (0018,9018) | set by creator   |
| Spectrally Selected Suppression     | (0018,9025) | set by creator   |
| Geometry of k-Space Traversal       | (0018,9032) | set by creator   |
| Rectilinear Phase Encode Reordering | (0018,9034) | set by creator   |
| Segmented k-Space Traversal         | (0018,9033) | set by creator   |
| Coverage of k-Space                 | (0018,9094) | set by creator   |
| Number of k-Space Trajectories      | (0018,9093) | set by creator   |
| Echo Peak Position                  | (0018,9298) | set by creator   |

### 9.1.3.19 MR Spectroscopy Data Module

| Attribute Name                     | Tag         | Supported Values |
|------------------------------------|-------------|------------------|
| Rows                               | (0028,0010) | set by creator   |
| Columns                            | (0028,0011) | set by creator   |
| Data Point Rows                    | (0028,9001) | set by creator   |
| Data Point Columns                 | (0028,9002) | set by creator   |
| Data Representation                | (0028,9108) | set by creator   |
| Signal Domain Columns              | (0028,9003) | set by creator   |
| Signal Domains Rows                | (0028,9235) | set by creator   |
| First Order Phase Correction Angle | (5600,0010) | set by creator   |
| Spectroscopy Data                  | (5600,0020) | set by creator   |

### 9.1.3.20 SOP Common Module

| Attribute Name         | Tag         | Supported Values                  |
|------------------------|-------------|-----------------------------------|
| SOP Class UID          | (0008,0016) | MR Spectroscopy SOP Class UID     |
| SOP Instance UID       | (0008,0018) | set by creator                    |
| Specific Character Set | (0008,0005) | set by creator                    |
| Instance Creation Date | (0008,0012) | date the SOP instance was created |
| Instance Creation Time | (0008,0013) | time the SOP instance was created |

## 9.1.4 MR Image

This chapter describes the DICOM attributes of MR Image instances created in interoperability mode. The following table lists the modules created by syngo® MR.

| Module                    | Reference                   |
|---------------------------|-----------------------------|
| Patient                   | see chapter 9.1.1.1         |
| General Study             | see chapter 9.1.1.2         |
| Patient Study             | see chapter 9.1.1.3         |
| General Series            | see chapter 9.1.1.4         |
| Frame of Reference        | see chapter 9.1.1.6         |
| General Equipment         | see chapter 9.1.1.7         |
| General Image             | see chapter 0               |
| Image Plane               | see chapter 9.1.4.2         |
| Image Pixel               | see chapter 9.1.1.9         |
| Contrast/Bolus            | see chapter 9.1.4.3         |
| MR Image                  | see chapter 9.1.4.4         |
| SOP Common                | see chapter 9.1.4.59.1.1.21 |
| Common Instance Reference | see chapter 9.1.4.6         |

### 9.1.4.1 General Image Module

| Attribute Name         | Tag         | Supported Values |
|------------------------|-------------|------------------|
| Instance Number        | (0020,0013) | set by creator   |
| Content Date           | (0008,0023) | set by creator   |
| Content Time           | (0008,0033) | set by creator   |
| Image Type             | (0008,0008) | set by creator   |
| Acquisition Number     | (0020,0012) | set by creator   |
| Acquisition Date       | (0008,0022) | set by creator   |
| Acquisition Time       | (0008,0032) | set by creator   |
| Acquisition DateTime   | (0008,002A) | set by creator   |
| Instance Number        | (0020,0013) | set by creator   |
| Patient Orientation    | (0020,0020) | set by creator   |
| Content Date           | (0008,0023) | set by creator   |
| Content Time           | (0008,0033) | set by creator   |
| Image Type             | (0008,0008) | set by creator   |
| Acquisition Number     | (0020,0012) | set by creator   |
| Acquisition Date       | (0008,0022) | set by creator   |
| Acquisition Time       | (0008,0032) | set by creator   |
| Acquisition DateTime   | (0008,002A) | set by creator   |
| Image Comments         | (0020,4000) | set by creator   |
| Presentation LUT Shape | (2050,0020) | set by creator   |

### 9.1.4.2 Image Plane Module

| Attribute Name              | Tag         | Supported Values |
|-----------------------------|-------------|------------------|
| Pixel Spacing               | (0028,0030) | set by creator   |
| Image Orientation (Patient) | (0020,0037) | set by creator   |
| Image Position (Patient)    | (0020,0032) | set by creator   |
| Slice Thickness             | (0018,0050) | set by creator   |
| Slice Location              | (0020,1041) | set by creator   |

### 9.1.4.3 Contrast/Bolus Module

| Attribute Name                               | Tag         | Supported Values |
|--|-------------|------------------|
| Contrast/Bolus Agent                         | (0018,0010) | set by creator   |
| Contrast/Bolus Agent Sequence                | (0018,0012) | set by creator   |
| >Code Value                                  | (0008,0100) | set by creator   |
| >Coding Scheme Designator                    | (0008,0102) | set by creator   |
| >Coding Scheme Version                       | (0008,0103) | set by creator   |
| >Code Meaning                                | (0008,0104) | set by creator   |
| Contrast/Bolus Route                         | (0018,1040) | set by creator   |
| Contrast/Bolus Administration Route Sequence | (0018,0014) | set by creator   |
| >Code Value                                  | (0008,0100) | set by creator   |
| >Coding Scheme Designator                    | (0008,0102) | set by creator   |
| >Coding Scheme Version                       | (0008,0103) | set by creator   |
| >Code Meaning                                | (0008,0104) | set by creator   |
| Contrast/Bolus Volume                        | (0018,1041) | set by creator   |
| Contrast/Bolus Start Time                    | (0018,1042) | set by creator   |
| Contrast/Bolus Stop Time                     | (0018,1043) | set by creator   |
| Contrast/Bolus Total Dose                    | (0018,1044) | set by creator   |
| Contrast Flow Rate                           | (0018,1046) | set by creator   |
| Contrast Flow Duration                       | (0018,1047) | set by creator   |
| Contrast/Bolus Ingredient                    | (0018,1048) | set by creator   |
| Contrast/Bolus Ingredient Concentration      | (0018,1049) | set by creator   |

### 9.1.4.4 MR Image Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Image Type                 | (0008,0008) | set by creator   |
| Samples per Pixel          | (0028,0002) | set by creator   |
| Photometric Interpretation | (0028,0004) | set by creator   |
| Bits Allocated             | (0028,0100) | set by creator   |
| Scanning Sequence          | (0018,0020) | set by creator   |
| Sequence Variant           | (0018,0021) | set by creator   |
| Scan Options               | (0018,0022) | set by creator   |
| MR Acquisition Type        | (0018,0023) | set by creator   |
| Repetition Time            | (0018,0080) | set by creator   |
| Echo Time                  | (0018,0081) | set by creator   |
| Echo Train Length          | (0018,0091) | set by creator   |

|                                   |             |                |
|-----------------------------------|-------------|----------------|
| Inversion Time                    | (0018,0082) | set by creator |
| Trigger Time                      | (0018,1060) | set by creator |
| Scan Options                      | (0018,0022) | set by creator |
| MR Acquisition Type               | (0018,0023) | set by creator |
| Repetition Time                   | (0018,0080) | set by creator |
| Echo Time                         | (0018,0081) | set by creator |
| Echo Train Length                 | (0018,0091) | set by creator |
| Inversion Time                    | (0018,0082) | set by creator |
| Trigger Time                      | (0018,1060) | set by creator |
| Sequence Name                     | (0018,0024) | set by creator |
| Angio Flag                        | (0018,0025) | set by creator |
| Number of Averages                | (0018,0083) | set by creator |
| Imaging Frequency                 | (0018,0084) | set by creator |
| Imaged Nucleus                    | (0018,0085) | set by creator |
| Echo Number(s)                    | (0018,0086) | set by creator |
| Magnetic Field Strength           | (0018,0087) | set by creator |
| Spacing Between Slices            | (0018,0088) | set by creator |
| Number of Phase Encoding Steps    | (0018,0089) | set by creator |
| Percent Sampling                  | (0018,0093) | set by creator |
| Percent Phase Field of View       | (0018,0094) | set by creator |
| Pixel Bandwidth                   | (0018,0095) | set by creator |
| Nominal Interval                  | (0018,1062) | set by creator |
| Beat Rejection Flag               | (0018,1080) | set by creator |
| Low R-R Value                     | (0018,1081) | set by creator |
| High R-R Value                    | (0018,1082) | set by creator |
| Intervals Acquired                | (0018,1083) | set by creator |
| Intervals Rejected                | (0018,1084) | set by creator |
| PVC Rejection                     | (0018,1085) | set by creator |
| Skip Beats                        | (0018,1086) | set by creator |
| Heart Rate                        | (0018,1088) | set by creator |
| Cardiac Number of Images          | (0018,1090) | set by creator |
| Trigger Window                    | (0018,1094) | set by creator |
| Reconstruction Diameter           | (0018,1100) | set by creator |
| Receive Coil Name                 | (0018,1250) | set by creator |
| Transmit Coil Name                | (0018,1251) | set by creator |
| Acquisition Matrix                | (0018,1310) | set by creator |
| In-plane Phase Encoding Direction | (0018,1312) | set by creator |
| Flip Angle                        | (0018,1314) | set by creator |
| SAR                               | (0018,1316) | set by creator |
| Variable Flip Angle Flag          | (0018,1315) | set by creator |
| dB/dt                             | (0018,1318) | set by creator |
| Temporal Position Identifier      | (0020,0100) | set by creator |
| Number of Temporal Positions      | (0020,0105) | set by creator |
| Temporal Resolution               | (0020,0110) | set by creator |

#### 9.1.4.5 SOP Common Module

| Attribute Name | Tag         | Supported Values       |
|----------------|-------------|------------------------|
| SOP Class UID  | (0008,0016) | MR Image SOP Class UID |

|                        |             |                                   |
|------------------------|-------------|-----------------------------------|
| SOP Instance UID       | (0008,0018) | set by creator                    |
| Specific Character Set | (0008,0005) | set by creator                    |
| Instance Creation Date | (0008,0012) | date the SOP instance was created |
| Instance Creation Time | (0008,0013) | time the SOP instance was created |

#### 9.1.4.6 Common Instance Reference Module

| Attribute Name                 | Tag         | Supported Values |
|--------------------------------|-------------|------------------|
| Referenced Series Sequence     | (0008,1115) | set by creator   |
| >Series Instance UID           | (0020,000E) | set by creator   |
| > Referenced Instance Sequence | (0008,114A) | set by creator   |
| >>Referenced SOP Class UID     | (0008,1150) | set by creator   |
| >> Referenced SOP Instance UID | (0008,1155) | set by creator   |

#### 9.1.5 Raw Data

Neuro Diffusion Tensor Imaging data are encoded in the Raw Data SOP Class.

#### 9.1.6 Multi-frame Grayscale Word SC Image

For some usecases Multi-frame Grayscale Word SC Image instances are created for documentation and reporting purposes. The following table lists the modules created by syngo® MR.

| Module                        | Reference           |
|-------------------------------|---------------------|
| Patient                       | see chapter 9.1.1.1 |
| General Study                 | see chapter 9.1.1.2 |
| Patient Study                 | see chapter 9.1.1.3 |
| General Series                | see chapter 9.1.1.4 |
| General Equipment             | see chapter 9.1.1.7 |
| SC Equipment                  | see chapter 9.1.6.1 |
| Frame of Reference            | see chapter 9.1.1.6 |
| General Image                 | see chapter 0       |
| Image Pixel                   | see chapter 9.1.1.9 |
| Cine                          | see chapter 9.1.6.2 |
| Multi-frame                   | see chapter 9.1.6.3 |
| Frame Pointers                | see chapter 9.1.6.4 |
| Multi-frame Functional Groups | see chapter 9.1.6.5 |
| SC Image                      | see chapter 9.1.6.7 |
| SC Multi-frame Image          | see chapter 9.1.6.8 |
| SC Multi-frame Vector         | see chapter 9.1.6.9 |

|                           |                             |
|---------------------------|-----------------------------|
| VOI LUT                   | see chapter 9.1.6.10        |
| SOP Common                | see chapter 9.1.4.59.1.1.21 |
| Common Instance Reference | see chapter 9.1.4.6         |

### 9.1.6.1 SC Equipment Module

| Attribute Name  | Tag         | Supported Values |
|-----------------|-------------|------------------|
| Conversion Type | (0008,0064) | set by creator   |
| Modality        | (0008,0060) | MR               |

### 9.1.6.2 Cine Module

| Attribute Name                 | Tag         | Supported Values |
|--------------------------------|-------------|------------------|
| Preferred Playback Sequencing  | (0018,1244) | set by creator   |
| Frame Time                     | (0018,1063) | set by creator   |
| Frame Time Vector              | (0018,1065) | set by creator   |
| Start Trim                     | (0008,2142) | set by creator   |
| Stop Trim                      | (0008,2143) | set by creator   |
| Recommended Display Frame Rate | (0008,2144) | set by creator   |
| Cine Rate                      | (0018,0040) | set by creator   |
| Frame Delay                    | (0018,1066) | set by creator   |
| Image Trigger Delay            | (0018,1067) | set by creator   |
| Effective Duration             | (0018,0072) | set by creator   |
| Actual Frame Duration          | (0018,1242) | set by creator   |

### 9.1.6.3 Multi-frame Module

| Attribute Name          | Tag         | Supported Values |
|-------------------------|-------------|------------------|
| Number of Frames        | (0028,0008) | set by creator   |
| Frame Increment Pointer | (0028,0009) | set by creator   |

### 9.1.6.4 Frame Pointers Module

| Attribute Name                  | Tag         | Supported Values |
|---------------------------------|-------------|------------------|
| Representative Frame Number     | (0028,6010) | set by creator   |
| Frame Numbers of Interest (FOI) | (0028,6020) | set by creator   |
| Frame of Interest Description   | (0028,6022) | set by creator   |
| Frame of Interest Type          | (0028,6023) | set by creator   |

### 9.1.6.5 Multi-frame Functional Groups Module

| Attribute Name | Tag | Supported Values |
|----------------|-----|------------------|
|----------------|-----|------------------|

|                                      |             |                |
|--------------------------------------|-------------|----------------|
| Shared Functional Groups Sequence    | (5200,9229) | set by creator |
| Per-frame Functional Groups Sequence | (5200,9230) | set by creator |
| Instance Number                      | (0020,0013) | set by creator |
| Content Date                         | (0008,0023) | set by creator |
| Content Time                         | (0008,0033) | set by creator |
| Number of Fr                         | (0028,0008) | set by creator |
| Representative Frame Number          | (0028,6010) | set by creator |

### 9.1.6.6 Multi-frame Grayscale Word SC Image Functional Group Macros

| Functional Group Macro   | Attribute Name               | Tag         | Supported Values |
|--|------------------------------|-------------|------------------|
| Pixel Measures<br>(if Plane Position (Patient) or Plane Orientation<br>(Patient) Macros Present) | Pixel Measures Sequence      | (0028,9110) | set by creator   |
|  | >Pixel Spacing               | (0018,0030) | set by creator   |
|  | >Slice Thickness             | (0028,0050) | set by creator   |
| Plane Position (Patient)<br>(if Pixel Measures or Plane Orientation (Patient)<br>Macros Present) | Plane Position Sequence      | (0020,9113) | set by creator   |
|  | >Image Position (Patient)    | (0020,0032) | set by creator   |
| Plane Orientation (Patient)<br>(if Pixel Measures or Plane Position (Patient)<br>Macros Present) | Plane Orientation Sequence   | (0020,9116) | set by creator   |
|  | >Image Orientation (Patient) | (0020,0037) | set by creator   |

### 9.1.6.7 SC Image Module

| Attribute Name                        | Tag         | Supported Values |
|---------------------------------------|-------------|------------------|
| Date of Secondary Capture             | (0018,1012) | set by creator   |
| Time of Secondary Capture             | (0018,1014) | set by creator   |
| Nominal Scanned Pixel Spa             | (0018,2010) | set by creator   |
| Pixel Spacing                         | (0028,0030) | set by creator   |
| Pixel Spacing Calibration Type        | (0028,0A02) | set by creator   |
| Pixel Spacing Calibration Description | (0018,1012) | set by creator   |

### 9.1.6.8 SC Multi-frame Image Module

| Attribute Name                | Tag         | Supported Values |
|-------------------------------|-------------|------------------|
| Burned In Annotation          | (0028,0301) | set by creator   |
| Recognizable Visual Features  | (0028,0302) | set by creator   |
| Presentation LUT Shape        | (2050,0020) | set by creator   |
| Illumination                  | (2010,015E) | set by creator   |
| Reflected Ambient Light       | (2010,0160) | set by creator   |
| Rescale Intercept             | (0028,1052) | set by creator   |
| Rescale Slope                 | (0028,1053) | set by creator   |
| Rescale Type                  | (0028,1054) | set by creator   |
| Frame Increment Pointer       | (0028,0009) | set by creator   |
| Nominal Scanned Pixel Spacing | (0018,2010) | set by creator   |
| Pixel Spacing                 | (0028,0030) | set by creator   |

|                                       |             |                |
|---------------------------------------|-------------|----------------|
| Pixel Spacing Calibration Type        | (0028,0A02) | set by creator |
| Pixel Spacing Calibration Description | (0028,0A04) | set by creator |
| Digitizing Device Transport Direction | (0018,2020) | set by creator |
| Rotation of Scanned Film              | (0018,2030) | set by creator |

#### 9.1.6.9 SC Multi-frame Vector Module

| Attribute Name               | Tag         | Supported Values |
|------------------------------|-------------|------------------|
| Frame Time Vector            | (0018,1065) | set by creator   |
| Page Number Vector           | (0018,2001) | set by creator   |
| Frame Label Vector           | (0018,2002) | set by creator   |
| Frame Primary Angle Vector   | (0018,2003) | set by creator   |
| Frame Secondary Angle Vector | (0018,2004) | set by creator   |
| Slice Location Vector        | (0018,2005) | set by creator   |
| Display Window Label Vector  | (0018,2006) | set by creator   |

#### 9.1.6.10 VOI LUT Module

| Attribute Name                      | Tag         | Supported Values |
|-------------------------------------|-------------|------------------|
| VOI LUT Sequence                    | (0028,3010) | set by creator   |
| >LUT Descriptor                     | (0028,3002) | set by creator   |
| >LUT Explanation                    | (0028,3003) | set by creator   |
| >LUT Data                           | (0028,3006) | set by creator   |
| Window Center                       | (0028,1050) | set by creator   |
| Window Width                        | (0028,1051) | set by creator   |
| Window Center and Width Explanation | (0028,1055) | set by creator   |
| VOI LUT Function                    | (0028,1056) | set by creator   |

### 9.1.7 Secondary Capture Image

This chapter describes the DICOM attributes of Secondary Capture Image instances created in interoperability mode. The following table lists the modules created by syngo® MR.

| Module            | Reference           |
|-------------------|---------------------|
| Patient           | see chapter 9.1.1.1 |
| General Study     | see chapter 9.1.1.2 |
| Patient Study     | see chapter 9.1.1.3 |
| General Series    | see chapter 9.1.1.4 |
| General Equipment | see chapter 9.1.1.7 |
| SC Equipment      | see chapter 9.1.6.1 |
| General Image     | see chapter 0       |
| Image Pixel       | see chapter 9.1.1.9 |



|                           |                             |
|---------------------------|-----------------------------|
| SC Image                  | see chapter 9.1.6.7         |
| Overlay Plane             | see chapter 9.1.7.1         |
| Modality LUT              | see chapter 9.1.7.2         |
| VOI LUT                   | see chapter 9.1.6.10        |
| SOP Common                | see chapter 9.1.4.59.1.1.21 |
| Common Instance Reference | see chapter 9.1.4.6         |

### 9.1.7.1 Overlay Plane Module

| Attribute Name         | Tag         | Supported Values |
|------------------------|-------------|------------------|
| Overlay Rows           | (60xx,0010) | set by creator   |
| Overlay Columns        | (60xx,0011) | set by creator   |
| Overlay Type           | (60xx,0040) | set by creator   |
| Overlay Origin         | (60xx,0050) | set by creator   |
| Overlay Bits Allocated | (60xx,0100) | set by creator   |
| Overlay Bit Position   | (60xx,0102) | set by creator   |
| Overlay Data           | (60xx,3000) | set by creator   |
| Overlay Description    | (60xx,0022) | set by creator   |
| Overlay Subtype        | (60xx,0045) | set by creator   |
| Overlay Label          | (60xx,1500) | set by creator   |
| ROI Area               | (60xx,1301) | set by creator   |
| ROI Mean               | (60xx,1302) | set by creator   |
| ROI Standard Deviation | (60xx,1303) | set by creator   |

### 9.1.7.2 Modality LUT Module

| Attribute Name        | Tag         | Supported Values |
|-----------------------|-------------|------------------|
| Modality LUT Sequence | (0028,3000) | set by creator   |
| >LUT Descriptor       | (0028,3002) | set by creator   |
| >LUT Explanation      | (0028,3003) | set by creator   |
| >Modality LUT Type    | (0028,3004) | set by creator   |
| >LUT Data             | (0028,3006) | set by creator   |
| Rescale Intercept     | (0028,1052) | set by creator   |
| Rescale Slope         | (0028,1053) | set by creator   |
| Rescale Type          | (0028,1054) | set by creator   |

## 9.1.8 Grayscale Softcopy Presentation State

For some usecases Grayscale Softcopy Presentation State instances are created referencing Enhanced MR Image or Multi-frame Grayscale Word SC Image instances. The following table lists the modules created by syngo® MR.

| Module | Reference |
|--------|-----------|
|--------|-----------|

|                                   |                             |
|-----------------------------------|-----------------------------|
| Patient                           | see chapter 9.1.1.1         |
| General Study                     | see chapter 9.1.1.2         |
| Patient Study                     | see chapter 9.1.1.3         |
| General Series                    | see chapter 9.1.1.4         |
| Presentation Series               | see chapter 9.1.8.1         |
| General Equipment                 | see chapter 9.1.1.7         |
| Presentation State Identification | see chapter 9.1.8.2         |
| Presentation State Relationship;  | see chapter 9.1.8.3         |
| Displayed Area                    | see chapter 9.1.8.4         |
| Graphic Annotation                | see chapter 9.1.8.5         |
| Spatial Transformation            | see chapter 9.1.8.6         |
| Graphic Layer                     | see chapter 9.1.8.7         |
| Graphic Group                     | see chapter 9.1.8.8         |
| Modality LUT                      | see chapter 9.1.7.2         |
| Softcopy VOI LUT                  | see chapter 9.1.8.9         |
| Softcopy Presentation LUT         | see chapter 9.1.8.10        |
| SOP Common                        | see chapter 9.1.4.59.1.1.21 |

### 9.1.8.1 Presentation Series Module

| Attribute Name | Tag         | Supported Values |
|----------------|-------------|------------------|
| Modality       | (0008,0060) | PR               |

### 9.1.8.2 Presentation State Identification Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Presentation Creation Date | (0070,0082) | set by creator   |
| Presentation Creation Time | (0070,0083) | set by creator   |
| Instance Number            | (0020,0013) | set by creator   |
| Content Label              | (0070,0080) | set by creator   |
| Content Description        | (0070,0081) | set by creator   |
| Content Creator's Name     | (0070,0084) | set by creator   |

### 9.1.8.3 Presentation State Relationship Module

| Attribute Name             | Tag         | Supported Values |
|----------------------------|-------------|------------------|
| Referenced Series Sequence | (0008,1115) | set by creator   |
| >Series Instance UID       | (0020,000E) | set by creator   |

|                               |             |                |
|-------------------------------|-------------|----------------|
| >Referenced Image Sequence    | (0008,1140) | set by creator |
| >>Referenced SOP Class UID    | (0008,1150) | set by creator |
| >>Referenced SOP Instance UID | (0008,1155) | set by creator |
| >>Referenced Frame Number     | (0008,1160) | set by creator |

#### 9.1.8.4 Displayed Area Module

| Attribute Name                           | Tag         | Supported Values |
|--|-------------|------------------|
| Displayed Area Selection Sequence        | (0070,005A) | set by creator   |
| >Referenced Image Sequence               | (0008,1140) | set by creator   |
| >>Referenced SOP Class UID               | (0008,1150) | set by creator   |
| >>Referenced SOP Instance UID            | (0008,1155) | set by creator   |
| >>Referenced Frame Number                | (0008,1160) | set by creator   |
| >Pixel Origin Interpretation             | (0048,0301) | set by creator   |
| >Displayed Area Top Left Hand Corner     | (0070,0052) | set by creator   |
| >Displayed Area Bottom Right Hand Corner | (0070,0053) | set by creator   |
| >Presentation Size Mode                  | (0070,0100) | set by creator   |
| >Presentation Pixel Spacing              | (0070,0101) | set by creator   |
| >Presentation Pixel Aspect Ratio         | (0070,0102) | set by creator   |
| >Presentation Pixel Magnification Ratio  | (0070,0103) | set by creator   |

#### 9.1.8.5 Graphic Annotation Module

| Attribute Name                  | Tag         | Supported Values |
|---------------------------------|-------------|------------------|
| Graphic Annotation Sequence     | (0070,0001) | set by creator   |
| >Referenced Image Sequence      | (0008,1140) | set by creator   |
| >>Referenced SOP Class UID      | (0008,1150) | set by creator   |
| >>Referenced SOP Instance UID   | (0008,1155) | set by creator   |
| >>Referenced Frame Number       | (0008,1160) | set by creator   |
| >Graphic Layer                  | (0070,0002) | set by creator   |
| >Text Object Sequence           | (0070,0008) | set by creator   |
| >>Bounding Box Annotation Units | (0070,0003) | set by creator   |
| >>Anchor Point Annotation Units | (0070,0004) | set by creator   |
| >>Unformatted Text Value        | (0070,0006) | set by creator   |
| >>Text Style Sequence           | (0070,0231) | set by creator   |
| >>>Font Name                    | (0070,0227) | set by creator   |
| >>>Font Name Type               | (0070,0228) | set by creator   |
| >>>CSS Font Name                | (0070,0229) | set by creator   |
| >>>Text Color CIE Lab Value     | (0070,0241) | set by creator   |
| >>>Horizontal Alignment         | (0070,0242) | set by creator   |
| >>>Vertical Alignment           | (0070,0243) | set by creator   |
| >>>Shadow Style                 | (0070,0244) | set by creator   |
| >>>Shadow Offset X              | (0070,0245) | set by creator   |
| >>>Shadow Offset Y              | (0070,0246) | set by creator   |
| >>>Shadow Color CIE Lab Val     | (0070,0247) | set by creator   |
| >>>Shadow Opacity               | (0070,0258) | set by creator   |
| >>>Underlined                   | (0070,0248) | set by creator   |

|  |             |                |
|--|-------------|----------------|
| >>>Bold                                      | (0070,0249) | set by creator |
| >>>Italic                                    | (0070,0250) | set by creator |
| >>Bounding Box Top Left Hand Corner          | (0070,0010) | set by creator |
| >>Bounding Box Bottom Right Hand Corner      | (0070,0011) | set by creator |
| >>Bounding Box Text Horizontal Justification | (0070,0012) | set by creator |
| >>Anchor Point                               | (0070,0014) | set by creator |
| >>Anchor Point Visibility                    | (0070,0015) | set by creator |
| >>Compound Graphic Instance ID               | (0070,0226) | set by creator |
| >>Graphic Group ID                           | (0070,0295) | set by creator |
| >Graphic Object Sequence                     | (0070,0009) | set by creator |
| >>Graphic Annotation Units                   | (0070,0005) | set by creator |
| >>Graphic Dimensions                         | (0070,0020) | set by creator |
| >>Number of Graphic Points                   | (0070,0021) | set by creator |
| >>Graphic Data                               | (0070,0022) | set by creator |
| >>Graphic Type                               | (0070,0023) | set by creator |
| >>Line Style Sequence                        | (0070,0232) | set by creator |
| >>>Pattern On Color CIELab Value;            | (0070,0251) | set by creator |
| >>>Pattern Off Color CIELab Value            | (0070,0252) | set by creator |
| >>>Pattern On Opacity                        | (0070,0284) | set by creator |
| >>>Pattern Off Opacity                       | (0070,0285) | set by creator |
| >>>Line Thickness                            | (0070,0253) | set by creator |
| >>>Line Dashing Style                        | (0070,0254) | set by creator |
| >>>Line Pattern                              | (0070,0255) | set by creator |
| >>>Shadow Style                              | (0070,0244) | set by creator |
| >>>Shadow Offset X                           | (0070,0245) | set by creator |
| >>>Shadow Offset Y                           | (0070,0246) | set by creator |
| >>>Shadow Color CIELab Val                   | (0070,0247) | set by creator |
| >>>Shadow Opacity                            | (0070,0258) | set by creator |
| >>Graphic Filled                             | (0070,0024) | set by creator |
| >>Fill Style Sequence                        | (0070,0233) | set by creator |
| >>>Pattern On Color CIELab Value;            | (0070,0251) | set by creator |
| >>>Pattern Off Color CIELab Value            | (0070,0252) | set by creator |
| >>>Pattern On Opacity                        | (0070,0284) | set by creator |
| >>>Pattern Off Opacity                       | (0070,0285) | set by creator |
| >>>Fill Mode                                 | (0070,0257) | set by creator |
| >>>Fill Pattern                              | (0070,0256) | set by creator |
| >>Compound Graphic Instance ID               | (0070,0226) | set by creator |
| >>Graphic Group ID                           | (0070,0295) | set by creator |
| >Compound Graphic Sequence                   | (0070,0209) | set by creator |
| >>Compound Graphic Instance                  | (0070,0226) | set by creator |
| >>Compound Graphic Units                     | (0070,0282) | set by creator |
| >>Graphic Dimensions                         | (0070,0020) | set by creator |
| >>Number of Graphic Points                   | (0070,0021) | set by creator |
| >>Graphic Data                               | (0070,0022) | set by creator |
| >>Compound Graphic Type                      | (0070,0294) | set by creator |
| >Text Style Sequence                         | (0070,0231) | set by creator |
| >>>Font Name                                 | (0070,0227) | set by creator |
| >>>Font Name Type                            | (0070,0228) | set by creator |
| >>>CSS Font Name                             | (0070,0229) | set by creator |
| >>>Text Color CIELab Value                   | (0070,0241) | set by creator |

|                                   |             |                |
|-----------------------------------|-------------|----------------|
| >>>Horizontal Alignment           | (0070,0242) | set by creator |
| >>>Vertical Alignment             | (0070,0243) | set by creator |
| >>>Shadow Style                   | (0070,0244) | set by creator |
| >>>Shadow Offset X                | (0070,0245) | set by creator |
| >>>Shadow Offset Y                | (0070,0246) | set by creator |
| >>>Shadow Color CIELab Val        | (0070,0247) | set by creator |
| >>>Shadow Opacity                 | (0070,0258) | set by creator |
| >>>Underlined                     | (0070,0248) | set by creator |
| >>>Bold                           | (0070,0249) | set by creator |
| >>>Italic                         | (0070,0250) | set by creator |
| >>Line Style Sequence             | (0070,0232) | set by creator |
| >>>Pattern On Color CIELab Value; | (0070,0251) | set by creator |
| >>>Pattern Off Color CIELab Value | (0070,0252) | set by creator |
| >>>Pattern On Opacity             | (0070,0284) | set by creator |
| >>>Pattern Off Opacity            | (0070,0285) | set by creator |
| >>>Line Thickness                 | (0070,0253) | set by creator |
| >>>Line Dashing Style             | (0070,0254) | set by creator |
| >>>Line Pattern                   | (0070,0255) | set by creator |
| >>>Shadow Style                   | (0070,0244) | set by creator |
| >>>Shadow Offset X                | (0070,0245) | set by creator |
| >>>Shadow Offset Y                | (0070,0246) | set by creator |
| >>>Shadow Color CIELab Val        | (0070,0247) | set by creator |
| >>>Shadow Opacity                 | (0070,0258) | set by creator |
| >>Rotation Angle                  | (0070,0230) | set by creator |
| >>Rotation Point                  | (0070,0273) | set by creator |
| >>Gap Length                      | (0070,0261) | set by creator |
| >>Diameter of Visibility          | (0070,0262) | set by creator |
| >>Major Ticks Sequence            | (0070,0287) | set by creator |
| >>>Tick Position                  | (0070,0288) | set by creator |
| >>>Tick Label                     | (0070,0289) | set by creator |
| >>Tick Alignment                  | (0070,0274) | set by creator |
| >>Tick Label Alignment            | (0070,0279) | set by creator |
| >>Show Tick Label                 | (0070,0278) | set by creator |
| >>Graphic Filled                  | (0070,0024) | set by creator |
| >>Fill Style Sequence             | (0070,0233) | set by creator |
| >>>Pattern On Color CIELab Value; | (0070,0251) | set by creator |
| >>>Pattern Off Color CIELab Value | (0070,0252) | set by creator |
| >>>Pattern On Opacity             | (0070,0284) | set by creator |
| >>>Pattern Off Opacity            | (0070,0285) | set by creator |
| >>>Fill Mode                      | (0070,0257) | set by creator |
| >>>Fill Pattern                   | (0070,0256) | set by creator |
| >>Graphic Group ID                | (0070,0295) | set by creator |

### 9.1.8.6 Spatial Transformation Module

| Attribute Name        | Tag         | Supported Values |
|-----------------------|-------------|------------------|
| Image Rotation        | (0070,0042) | set by creator   |
| Image Horizontal Flip | (0070,0041) | set by creator   |

### 9.1.8.7 Graphic Layer Module

| Attribute Name                                     | Tag         | Supported Values |
|--|-------------|------------------|
| Graphic Layer Sequence                             | (0070,0060) | set by creator   |
| >Graphic Layer                                     | (0070,0002) | set by creator   |
| >Graphic Layer Order                               | (0070,0062) | set by creator   |
| >Graphic Layer Recommended Display Grayscale Value | (0070,0066) | set by creator   |
| >Graphic Layer Recommended Display CIELab Value    | (0070,0401) | set by creator   |
| >Graphic Layer Description                         | (0070,0068) | set by creator   |

### 9.1.8.8 Graphic Group Module

| Attribute Name              | Tag         | Supported Values |
|-----------------------------|-------------|------------------|
| Graphic Group Sequence      | (0070,0234) | set by creator   |
| > Graphic Group ID          | (0070,0295) | set by creator   |
| > Graphic Group Label       | (0070,0207) | set by creator   |
| > Graphic Group Description | (0070,0208) | set by creator   |

### 9.1.8.9 Softcopy VOI LUT Module

| Attribute Name                       | Tag         | Supported Values |
|--------------------------------------|-------------|------------------|
| Softcopy VOI LUT Sequence            | (0028,3110) | set by creator   |
| >Referenced Image Sequence           | (0008,1140) | set by creator   |
| >>Referenced SOP Class UID           | (0008,1150) | set by creator   |
| >>Referenced SOP Instance UID        | (0008,1155) | set by creator   |
| >>Referenced Frame Number            | (0008,1160) | set by creator   |
| >VOI LUT Sequence                    | (0028,3010) | set by creator   |
| >>LUT Descriptor                     | (0028,3002) | set by creator   |
| >>LUT Explanation                    | (0028,3003) | set by creator   |
| >>LUT Data                           | (0028,3006) | set by creator   |
| >Window Center                       | (0028,1050) | set by creator   |
| >Window Width                        | (0028,1051) | set by creator   |
| >Window Center and Width Explanation | (0028,1055) | set by creator   |
| >VOI LUT Function                    | (0028,1056) | set by creator   |

### 9.1.8.10 Softcopy Presentation LUT Module

| Attribute Name            | Tag         | Supported Values |
|---------------------------|-------------|------------------|
| Presentation LUT Sequence | (2050,0010) | set by creator   |
| >LUT Descriptor           | (0028,3002) | set by creator   |
| >LUT Explanation          | (0028,3003) | set by creator   |
| >LUT Data                 | (0028,3006) | set by creator   |
| Presentation LUT Shape    | (2050,0020) | set by creator   |

## 9.1.9 Evidence Documents

Evidence Documents will be created by applications e.g. cardiac analysis to store evaluated results. The Evidence Documents are encoded in the SR Enhanced SOP Class.

### 9.1.9.1 Evidence Document Templates

The finding related results of applications are stored into the Content Sequence of a Structured Evidence Document. The Content Items of a Content Sequence are specified in Structured Reporting Templates.

Examples of Structured Reporting Templates are:

- BI-RADS reporting
- Cardiac reporting
- PI-RADS 2 reporting
- PhoenixZip documentation.

## 9.2 Data Dictionary of Private Attributes

The following table Table 92: Private Data Element Dictionary lists private attributes created by syngo® MR which may be included in the generated instances. These private attributes may be deprecated or replaced with standard attributes in the future.

**Table 92: Private Data Element Dictionary**

| DICOM Tag                                     | Name                            | VR | VM |
|---|---------------------------------|----|----|
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,01) | Business Unit Code              | CS | 1  |
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,02) | Application Type                | LO | 1  |
| (0027,SIEMENS SYNGO ENHANCED IDATASET API,03) | Application Attributes Sequence | SQ | 1  |
| (0029,SIEMENS SYNGO FUNCTION ASSIGNMENT,01)   | Data Reference                  | LO | 1  |
| (0009,SIEMENS SYNGO INDEX SERVICE,20)         | Object Insertion Date           | DA | 1  |
| (0009,SIEMENS SYNGO INDEX SERVICE,A0)         | Sender System Device Name       | LO | 1  |
| (0029,SIEMENS SYNGO VOLUME,12)                | Slices                          | US | 1  |
| (0029,SIEMENS SYNGO VOLUME,14)                | Volume Histogram                | OB | 1  |
| (0029,SIEMENS SYNGO VOLUME,18)                | Volume Level                    | IS | 1  |
| (0029,SIEMENS SYNGO VOLUME,30)                | Voxel Spacing                   | DS | 3  |
| (0029,SIEMENS SYNGO VOLUME,32)                | Volume Position (Patient)       | DS | 3  |
| (0029,SIEMENS SYNGO VOLUME,37)                | Volume Orientation (Patient)    | DS | 9  |
| (0029,SIEMENS SYNGO VOLUME,40)                | Resampling Flag                 | CS | 1  |
| (0029,SIEMENS SYNGO VOLUME,42)                | Normalization Flag              | CS | 1  |

| DICOM Tag                            | Name                           | VR | VM  |
|--------------------------------------|--------------------------------|----|-----|
| (0029,SIEMENS SYNGO VOLUME,44)       | SubVolume Sequence             | SQ | 1-n |
| (0071,SIEMENS SYNGO REGISTRATION,20) | Registered Image Sequence      | SQ | 1   |
| (0071,SIEMENS SYNGO REGISTRATION,21) | Registration Is Validated Flag | CS | 1   |
| (7FDF,SIEMENS SYNGO DATA PADDING,FC) | Pixel Data Leading Padding     | OB | 1   |
| (0021, SIEMENS MR SDS 01, 0C)        | Positive PCS Directions        | SH | 1   |
| (0021, SIEMENS MR SDS 01, 5E)        | Field Of View Text             | LO | 1   |
| (0021, SIEMENS MR SDS 01, 5F)        | Relative Table Position Text   | SH | 1   |
| (0021, SIEMENS MR SDS 01, FE)        | Series Data Sequence           | SQ | 1   |
| (0021, SIEMENS MR SDI 02, 4F)        | Coil String                    | LO | 1   |
| (0021, SIEMENS MR SDI 02, 56)        | PAT Mode Text                  | LO | 1   |
| (0021, SIEMENS MR SDI 02, 58)        | Acquisition Matrix Text        | SH | 1   |
| (0021, SIEMENS MR SDI 02, 88)        | Slice Position                 | DS | 1   |
| (0021, SIEMENS MR SDI 02, 89)        | Slice Position Text            | SH | 1   |
| (0021, SIEMENS MR SDI 02, FE)        | Image Data Sequence            | SQ | 1   |
| (0021, SIEMENS MR SDR 01, 01)        | Creator Identifier             | LO | 1   |
| (0021, SIEMENS MR SDR 01, 02)        | Application Identifier         | LO | 1   |
| (0021, SIEMENS MR SDR 01, 03)        | Cause Identifier               | LO | 1   |
| (0051, SIEMENS MR HEADER, 0A)        | Meas Duration                  | LO | 1   |
| (0051, SIEMENS MR HEADER, 0C)        | Field Of View                  | LO | 1   |
| (0051, SIEMENS MR HEADER, 0D)        | Slice Position                 | SH | 1   |
| (0051, SIEMENS MR HEADER, 12)        | Rel Table Position             | SH | 1   |
| (0051, SIEMENS MR HEADER, 13)        | Positive PCS Directions        | SH | 1   |
| (7FE1, SIEMENS MR IMA, 10)           | Raw Data                       | OB | 1   |

Interpretation of the DICOM Tags from the above table:

(gggg, pp,ee) -> (gggg, ppee)

gggg - odd group number

pp - private creator identification code

ee - private element

## 9.3 Grayscale Image Consistency

The high resolution TFT display monitor option of syngo® MR comes with a DICOM Grayscale Standard Display Function (GSDF) compliant factory pre-setting. A typical working environment setup is assumed for ambient light.



## 9.4 DICOM Print SCU – detailed status displays

The following tables document the behavior of the syngo® MR DICOM Print AE in response to messages received for the printer SOP class and the print job SOP class.

Definitions of camera symbols:

- Idle: Camera is installed and ready; idle icon is displayed.
- Interact: The user has to react in near future, but not immediately.  
Example: A camera was low in 8x10 clear sheets: LOW 8x10 CLR was sent by N-EVENT-REPORT.
- Queue Stopped: The user has to react immediately. Either the camera needs immediate interaction or a job has been aborted.  
Example: A camera is out of 8x10 clear sheets, or camera is down, or a film job is aborted.

Note: different camera symbols are displayed according to the Printer Status Info.

### 9.4.1 Common Status Information

**Table 93: Common Printer Status Information**

| Printer Status Info/<br>Execution Status Info | Description  | Message string visible<br>in 'Status Bar'                 | Other action for UI/<br>'camera symbol'              |
|---|--|---|--|
| NORMAL  | Camera is ready  | Camera is ready   | <None>/idle  |
| BAD RECEIVE MGZ                               | There is a problem with the film receive magazine. Films from the printer cannot be transported into the magazine.                 | Problem with receive magazine.                            | <None>/interact                                      |
| BAD SUPPLY MGZ                                | There is a problem with the film supply magazine. Films from this magazine cannot be transported into the printer.                 | Problem with supply magazine.                             | <None>/interact                                      |
| CALIBRATING                                   | Printer is performing self calibration, it is expected to be available for normal operation shortly.                               | Self calibration. Please wait.                            | <None>/idle  |
| CALIBRATION ERR                               | An error in the printer calibration has been detected, quality of processed films may not be optimal.                              | Problem in calibration. Film quality may not be optimal.  | <None>/interact                                      |
| CHECK CHEMISTRY                               | A problem with the processor chemicals has been detected, quality of processed films may not be optimal.                           | Problem with chemistry. Film quality may not be optimal.  | <None>/interact                                      |
| CHECK SORTER                                  | There is an error in the film sorter   | Error in film sorter.                                     | <None>/interact                                      |
| CHEMICALS EMPTY                               | There are no processing chemicals in the processor, films will not be printed and processed until the processor is back to normal. | Camera chemistry empty. Please check.                     | <None>/interact                                      |
| CHEMICALS LOW                                 | The chemical level in the processor is low, if not corrected, it will probably shut down soon.                                     | Camera chemistry low. Please check.                       | <None>/interact                                      |
| COVER OPEN                                    | One or more printer or processor covers, drawers, doors are open.  | Camera cover, drawer or door open.                        | <None>/interact                                      |
| ELEC CONFIG ERR                               | Printer configured improperly for this job.  | Camera configured improperly for this job. Queue stopped. | Queue for this camera will be STOPPED/ Queue stopped |
| ELEC DOWN                                     | Printer is not operating due to some unspecified electrical hardware problem.  | Camera electrical hardware Problem.                       | <None>/interact                                      |
| ELEC SW ERROR                                 | Printer not operating for some unspecified software error.   | Camera software problem. Queue stopped.                   | Queue for this camera will be STOPPED/ Queue stopped |
| EMPTY 8X10                                    | The 8x10 inch film supply magazine is empty.   | 8x10 film supply empty.                                   | <None>/interact                                      |
| EMPTY 8X10 BLUE                               | The 8x10 inch blue film supply magazine is empty.  | 8x10 blue film supply empty.                              | <None>/interact                                      |
| EMPTY 8X10 CLR                                | The 8x10 inch clear film supply magazine is empty.   | 8x10 clear film supply empty.                             | <None>/interact                                      |
| EMPTY 8X10 PAPER                              | The 8x10 inch paper supply magazine is empty.  | 8x10 paper supply empty.                                  | <None>/interact                                      |

| Printer Status Info/<br>Execution Status Info | Description   | Message string visible<br>in 'Status Bar' | Other action for UI/<br>'camera symbol' |
|---|---|---|---|
| EMPTY 10X12                                   | The 10x12 inch film supply magazine is empty.   | 10x12 film supply empty.                  | <None>/interact                         |
| EMPTY 10X12 BLUE                              | The 10x12 inch blue film supply magazine is empty.                                    | 10x12 blue film supply empty.             | <None>/interact                         |
| EMPTY 10X12 CLR                               | The 10x12 inch clear film supply magazine is empty.                                   | 10x12 clear film supply empty.            | <None>/interact                         |
| EMPTY 10X12 PAPR                              | The 10x12 inch paper supply magazine is empty.  | 10x12 paper supply empty.                 | <None>/interact                         |
| EMPTY 10X14                                   | The 10x14 inch film supply magazine is empty.   | 10x14 film supply empty.                  | <None>/interact                         |
| EMPTY 10X14 BLUE                              | The 10x14 inch blue film supply magazine is empty.                                    | 10x14 blue film supply empty.             | <None>/interact                         |
| EMPTY 10X14 CLR                               | The 10x14 inch clear film supply magazine is empty.                                   | 10x14 clear film supply empty.            | <None>/interact                         |
| EMPTY 10X14 PAPR                              | The 10x14 inch paper supply magazine is empty.  | 10x14 paper supply empty.                 | <None>/interact                         |
| EMPTY 11X14                                   | The 11x14 inch film supply magazine is empty.   | 11x14 film supply empty.                  | <None>/interact                         |
| EMPTY 11X14 BLUE                              | The 11x14 inch blue film supply magazine is empty.                                    | 11x14 blue film supply empty.             | <None>/interact                         |
| EMPTY 11X14 CLR                               | The 11x14 inch clear film supply magazine is empty.                                   | 11x14 clear film supply empty.            | <None>/interact                         |
| EMPTY 11X14 PAPR                              | The 11x14 inch paper supply magazine is empty.  | 11x14 paper supply empty.                 | <None>/interact                         |
| EMPTY 14X14                                   | The 14x14 inch film supply magazine is empty.   | 14x14 film supply empty.                  | <None>/interact                         |
| EMPTY 14X14 BLUE                              | The 14x14 inch blue film supply magazine is empty.                                    | 14x14 blue film supply empty.             | <None>/interact                         |
| EMPTY 14X14 CLR                               | The 14x14 inch clear film supply magazine is empty.                                   | 14x14 clear film supply empty.            | <None>/interact                         |
| EMPTY 14X14 PAPR                              | The 14x14 inch paper supply magazine is empty.  | 14x14 paper supply empty.                 | <None>/interact                         |
| EMPTY 14X17                                   | The 14x17 inch film supply magazine is empty.   | 14x17 film supply empty.                  | <None>/interact                         |
| EMPTY 14X17 BLUE                              | The 14x17 inch blue film supply magazine is empty.                                    | 14x17 blue film supply empty.             | <None>/interact                         |
| EMPTY 14X17 CLR                               | The 14x17 inch clear film supply magazine is empty.                                   | 14x17 clear film supply empty.            | <None>/interact                         |
| EMPTY 14X17 PAPR                              | The 14x17 inch paper supply magazine is empty.  | 14x17 paper supply empty.                 | <None>/interact                         |
| EMPTY 24X24                                   | The 24x24 inch film supply magazine is empty.   | 24x24 film supply empty.                  | <None>/interact                         |
| EMPTY 24X24 BLUE                              | The 24x24 inch blue film supply magazine is empty.                                    | 24x24 blue film supply empty.             | <None>/interact                         |
| EMPTY 24X24 CLR                               | The 24x24 inch clear film supply magazine is empty.                                   | 24x24 clear film supply empty.            | <None>/interact                         |
| EMPTY 24X24 PAPR                              | The 24x24 inch paper supply magazine is empty.  | 24x24 paper supply empty.                 | <None>/interact                         |
| EMPTY 24X30                                   | The 24x30 inch film supply magazine is empty.   | 24x30 film supply empty.                  | <None>/interact                         |
| EMPTY 24X30 BLUE                              | The 24x30 inch blue film supply magazine is empty.                                    | 24x30 blue film supply empty.             | <None>/interact                         |
| EMPTY 24X30 CLR                               | The 24x30 inch clear film supply magazine is empty.                                   | 24x30 clear film supply empty.            | <None>/interact                         |
| EMPTY 24X30 PAPR                              | The 24x30 inch paper supply magazine is empty.  | 24x30 paper supply empty.                 | <None>/interact                         |
| EMPTY A4 PAPR                                 | The A4 paper supply magazine is empty.  | A4 paper supply empty                     | <None>/interact                         |
| EMPTY A4 TRANS                                | The A4 transparency supply magazine is empty.   | A4 transparency supply empty.             | <None>/interact                         |
| EXPOSURE FAILURE                              | The exposure device has failed due to some unspecified reason.                        | Exposure device has failed.               | <None>/interact                         |
| FILM JAM                                      | A film transport error has occurred and a film is jammed in the printer or processor. | Film jam.                                 | <None>/interact                         |
| FILM TRANSP ERR                               | There is a malfunction with the film transport, there may or may not be a film jam.   | Film transport problem.                   | <None>/interact                         |
| FINISHER EMPTY                                | The finisher is empty.  | Finisher is empty.                        | <None>/interact                         |
| FINISHER ERROR                                | The finisher is not operating due to some unspecified reason                          | Finisher problem.                         | <None>/interact                         |

| Printer Status Info/<br>Execution Status Info | Description                                       | Message string visible<br>in 'Status Bar' | Other action for UI/<br>'camera symbol' |
|---|---|---|---|
| FINISHER LOW                                  | The finisher is low on supplies.                  | Finisher low.                             | <None>/interact                         |
| LOW 8X10                                      | The 8x10 inch film supply magazine is low.        | 8x10 film supply low.                     | <None>/interact                         |
| LOW 8X10 BLUE                                 | The 8x10 inch blue film supply magazine is low.   | 8x10 blue film supply low.                | <None>/interact                         |
| LOW 8X10 CLR                                  | The 8x10 inch clear film supply magazine is low.  | 8x10 clear film supply low.               | <None>/interact                         |
| LOW 8X10 PAPR                                 | The 8x10 inch paper supply magazine is low.       | 8x10 paper supply low.                    | <None>/interact                         |
| LOW 10X12                                     | The 10x12 inch film supply magazine is low.       | 10x12 film supply low.                    | <None>/interact                         |
| LOW 10X12 BLUE                                | The 10x12 inch blue film supply magazine is low.  | 10x12 blue film supply low.               | <None>/interact                         |
| LOW 10X12 CLR                                 | The 10x12 inch clear film supply magazine is low. | 10x12 clear film supply low.              | <None>/interact                         |
| LOW 10X12 PAPR                                | The 10x12 inch paper supply magazine is low.      | 10x12 paper supply low.                   | <None>/interact                         |
| LOW 10X14                                     | The 10x14 inch film supply magazine is low.       | 10x14 film supply low.                    | <None>/interact                         |
| LOW 10X14 BLUE                                | The 10x14 inch blue film supply magazine is low.  | 10x14 blue film supply low.               | <None>/interact                         |
| LOW 10X14 CLR                                 | The 10x14 inch clear film supply magazine is low. | 10x14 clear film supply low.              | <None>/interact                         |
| LOW 10X14 PAPR                                | The 10x14 inch paper supply magazine is low.      | 10x14 paper supply low.                   | <None>/interact                         |
| LOW 11X14                                     | The 11x14 inch film supply magazine is low.       | 11x14 film supply low.                    | <None>/interact                         |
| LOW 11X14 BLUE                                | The 11x14 inch blue film supply magazine is low.  | 11x14 blue film supply low.               | <None>/interact                         |
| LOW 11X14 CLR                                 | The 11x14 inch clear film supply magazine is low. | 11x14 clear film supply low.              | <None>/interact                         |
| LOW 11X14 PAPR                                | The 11x14 inch paper supply magazine is low.      | 11x14 paper supply low.                   | <None>/interact                         |
| LOW 14X14                                     | The 14x14 inch film supply magazine is low.       | 14x14 film supply low.                    | <None>/interact                         |
| LOW 14X14 BLUE                                | The 14x14 inch blue film supply magazine is low.  | 14x14 blue film supply low.               | <None>/interact                         |
| LOW 14X14 CLR                                 | The 14x14 inch clear film supply magazine is low. | 14x14 clear film supply low.              | <None>/interact                         |
| LOW 14X14 PAPR                                | The 14x14 inch paper supply magazine is low.      | 14x14 paper supply low.                   | <None>/interact                         |
| LOW 14X17                                     | The 14x17 inch film supply magazine is low.       | 14x17 film supply low.                    | <None>/interact                         |
| LOW 14X17 BLUE                                | The 14x17 inch blue film supply magazine is low.  | 14x17 blue film supply low.               | <None>/interact                         |
| LOW 14X17 CLR                                 | The 14x17 inch clear film supply magazine is low. | 14x17 clear film supply low.              | <None>/interact                         |
| LOW 14X17 PAPR                                | The 14x17 inch paper supply magazine is low.      | 14x17 paper supply low.                   | <None>/interact                         |
| LOW 24X24                                     | The 24x24 inch film supply magazine is low.       | 24x24 film supply low.                    | <None>/interact                         |
| LOW 24X24 BLUE                                | The 24x24 inch blue film supply magazine is low.  | 24x24 blue film supply low.               | <None>/interact                         |
| LOW 24X24 CLR                                 | The 24x24 inch clear film supply magazine is low. | 24x24 clear film supply low.              | <None>/interact                         |
| LOW 24X24 PAPR                                | The 24x24 inch paper supply magazine is low.      | 24x24 paper supply low.                   | <None>/interact                         |
| LOW 24X30                                     | The 24x30 inch film supply magazine is low.       | 24x30 film supply low.                    | <None>/interact                         |
| LOW 24X30 BLUE                                | The 24x30 inch blue film supply magazine is low.  | 24x30 blue film supply low.               | <None>/interact                         |
| LOW 24X30 CLR                                 | The 24x30 inch clear film supply magazine is low. | 24x30 clear film supply low.              | <None>/interact                         |
| LOW 24X30 PAPR                                | The 24x30 inch paper supply magazine is low.      | 24x30 paper supply low.                   | <None>/interact                         |
| LOW A4 PAPR                                   | The A4 paper supply magazine is low.              | A4 paper supply low.                      | <None>/interact                         |
| LOW A4 TRANS                                  | The A4 transparency supply magazine is low.       | A4 transparency supply low.               | <None>/interact                         |
| NO RECEIVE MGZ                                | The film receive magazine is not available.       | Film receiver not available.              | <None>/interact                         |

| Printer Status Info/<br>Execution Status Info | Description   | Message string visible<br>in 'Status Bar'                                | Other action for UI/<br>'camera symbol'              |
|---|---|--|--|
| NO RIBBON                                     | The ribbon cartridge needs to be re-placed.   | Replace ribbon cartridge.  | <None>/interact                                      |
| NO SUPPLY MGZ                                 | The film supply magazine is not available.  | Film supply not available.   | <None>/interact                                      |
| CHECK PRINTER                                 | The printer is not ready at this time, operator intervention is required to make the printer available.   | Check camera.  | <None>/interact                                      |
| CHECK PROC                                    | The processor is not ready at this time, operator intervention is required to make the printer available.   | Check processor.   | <None>/interact                                      |
| PRINTER DOWN                                  | The printer is not operating due to some unspecified reason.  | Camera down.   | <None>/interact                                      |
| PRINTER INIT                                  | The printer is not ready at this time, it is expected to become available without intervention. For example, it may be in a normal warm-up state.   | Camera initializing.   | <None>/Idle  |
| PRINTER OFFLINE                               | The printer has been disabled by an operator or service person.   | Camera off-line.   | <None>/interact                                      |
| PROC DOWN                                     | The processor is not operating due to some unspecified reason.  | Processor down.  | <None>/interact                                      |
| PROC INIT                                     | The processor is not ready at this time, it is expected to become available without intervention. For example, it may be in a normal warm-up state. | Processor initializing.  | <None>/Idle  |
| PROC OVERFLOW FL                              | Processor chemicals are approaching the overflow full mark.   | Processor chemicals near overflow.                                       | <None>/interact                                      |
| PROC OVERFLOW HI                              | Processor chemicals have reached the overflow full mark.  | Processor chemicals overflow.  | <None>/interact                                      |
| QUEUED  | Print job in Queue  | --   | <None>/Idle  |
| RECEIVER FULL                                 | The film receive magazine is full.  | Receiver full.   | <None>/interact                                      |
| REQ MED NOT INST                              | The requested film, paper, or other media supply magazine is installed in the printer, but may be available with operator intervention.             | Install media supply.  | <None>/interact                                      |
| REQ MED NOT AVAI                              | The requested film, paper, or other media requested is not available on this printer.   | Media supply not available on this camera. Queue stopped. Change camera. | Queue for this camera will be STOPPED/ Queue stopped |
| RIBBON ERROR                                  | There is an unspecified problem with the print ribbon.  | Error with print ribbon.   | <None>/interact                                      |
| SUPPLY EMPTY                                  | The printer is out of film.   | Camera out of film.  | <None>/interact                                      |
| SUPPLY LOW                                    | The film supply is low.   | Film supply low.   | <None>/interact                                      |
| UNKNOWN                                       | There is an unspecified problem.  | Unspecified problem with camera.   | <None>/interact                                      |

## 9.4.1 Additional DICOM Execution Status Information

Printer Status Info and Execution Status Info are defined terms and can therefore be extended or reduced by camera manufacturers. Therefore VA30A shall be flexible.

If any other printer status info or execution status info is received (as described in Table 9.4.1, VA30A will react as shown in the following table:

**Table 94: Additional Printer Status Information**

| Printer Status /<br>Execution | Printer / Execution Sta-<br>tus Info | Description               | Message string<br>visible in the Job<br>status bar | Other action for<br>syngo / camera<br>symbol                    |
|-------------------------------|--------------------------------------|---------------------------|--|---|
| WARNING                       | <any other>                          | <not defined status info> | Camera info: <sta-<br>tus info>                    | <None>/Interact   |
| FAILURE                       | <any other>                          | <not defined status info> | Camera info: <sta-<br>tus info><br>Queue stopped.  | Queue for this cam-<br>era will be<br>STOPPED/<br>Queue stopped |

## Annex A: Index of Tables

|   |    |
|---|----|
| Table 1: Network Services.....  | 2  |
| Table 2 - Media Services.....   | 4  |
| Table 3 - Implementation Identifying Information.....   | 5  |
| Table 4: Association Policies.....  | 18 |
| Table 5: Asynchronous Nature as an Association Initiator.....                                 | 18 |
| Table 6 - Presentation Context Table "Verification".....                                      | 19 |
| Table 7 - Presentation Context Table "Verification".....                                      | 19 |
| Table 8: Association Policies.....  | 20 |
| Table 9: Asynchronous Nature as an Association Initiator.....                                 | 20 |
| Table 10: Proposed Presentation Contexts for Storage.....                                     | 21 |
| Table 11: DICOM Command Response Status Handling Behavior.....                                | 23 |
| Table 12: DICOM Command Communication Failure Behavior.....                                   | 23 |
| Table 13: Storage C-STORE Response Status.....  | 24 |
| Table 14: Association Policies.....   | 26 |
| Table 15: Asynchronous Nature as an Association Initiator.....                                | 26 |
| Table 16: Proposed Presentation Contexts for Storage Commitment.....                          | 27 |
| Table 17: DICOM Command Response Status Handling Behavior.....                                | 27 |
| Table 18: DICOM Command Communication Failure Behavior.....                                   | 27 |
| Table 19 - Presentation Context Table "Update Flag Information".....                          | 28 |
| Table 20: DICOM Command Response Status Handling Behavior.....                                | 28 |
| Table 21: Association Policies.....   | 29 |
| Table 22: Asynchronous Nature as an Association Initiator.....                                | 29 |
| Table 23: Proposed Presentation Contexts for Query.....                                       | 30 |
| Table 24: Extended Negotiation as an SCU.....   | 30 |
| Table 25: DICOM Command Response Status Handling Behavior.....                                | 30 |
| Table 26: DICOM Command Communication Failure Behavior.....                                   | 31 |
| Table 27: Attributes supported for instance Query - SCU.....                                  | 31 |
| Table 28: Proposed Presentation Contexts for Retrieve and Activity "MOVE SCU".....            | 32 |
| Table 29: DICOM Command Response Status Handling Behavior.....                                | 33 |
| Table 30: DICOM Command Communication Failure Behavior.....                                   | 33 |
| Table 31: Association Policies.....   | 34 |
| Table 32: Asynchronous Nature as an Association Initiator.....                                | 34 |
| Table 33: Proposed Presentation Contexts for Worklist.....                                    | 35 |
| Table 34: Broad Query search keys.....  | 35 |
| Table 35: Modality Worklist C-Find Return keys.....   | 36 |
| Table 36: DICOM Command Response Status Handling Behavior.....                                | 39 |
| Table 37: DICOM Command Communication Failure Behavior.....                                   | 39 |
| Table 38: Association Policies.....   | 40 |
| Table 39: Asynchronous Nature as an Association Initiator.....                                | 40 |
| Table 40: Acceptable Presentation Contexts Activity "Create MPPS".....                        | 41 |
| Table 41: MPPS N-CREATE Response Status Handling Behavior.....                                | 41 |
| Table 42: Acceptable Presentation Contexts Activity "Update MPPS".....                        | 41 |
| Table 43: MPPS N-SET Response Status Handling Behavior.....                                   | 42 |
| Table 44: Association Policies.....   | 43 |
| Table 45: Asynchronous Nature as an Association Initiator.....                                | 43 |
| Table 46: Presentation Contexts for the Activity "Print Film".....                            | 44 |
| Table 47: Attributes for the N-CREATE-RQ of the Basic Film Session.....                       | 45 |
| Table 48: N-CREATE-RSP Status Handling Behavior for the Basic Film Session.....               | 45 |
| Table 49: Attributes for the N-CREATE-RQ of the Basic Film Session.....                       | 46 |
| Table 50: N-CREATE-RSP Status Handling Behavior for Basic Film Box.....                       | 47 |
| Table 51: N-ACTION-RSP Status Handling Behavior for Basic Film Box.....                       | 47 |
| Table 52: Attributes for N-SET-RQ of Basic Grayscale Image Box.....                           | 47 |
| Table 53: N-SET-RSP Status Handling Behavior for the Basic Grayscale Image Box SOP Class..... | 48 |
| Table 54: Attributes for N-SET-RQ of Basic Color Image Box.....                               | 48 |
| Table 55: N-SET-RSP Status Handling Behavior for the Color Grayscale Image Box.....           | 49 |
| Table 56: Attributes for N-CREATE-RQ of Presentation LUT SOP Class.....                       | 50 |

|  |     |
|--|-----|
| Table 57: N-CREATE-RSP Status Handling Behavior for the Presentation LUT SOP Class ..... | 50  |
| Table 58: Attributes for N-GET-RQ of the Printer SOP Class .....                         | 50  |
| Table 59: DICOM Command Communication Failure Behavior .....                             | 50  |
| Table 60: Presentation Contexts for the Activity "Print Film" .....                      | 51  |
| Table 61: Attributes for the N-EVENT-REPORT-RQ of the Printer SOP Class .....            | 51  |
| Table 62: Presentation Contexts for the Activity "Print Management" .....                | 52  |
| Table 63: Attributes for the N-EVENT-REPORT-RQ of the Print Job SOP Class .....          | 52  |
| Table 64 - Presentation Context Table "Print Film" .....                                 | 54  |
| Table 65 - Basic Film Session N-CREATE attributes .....                                  | 55  |
| Table 66 - Basic Film Session Status Codes .....   | 55  |
| Table 67 - Basic Film Box N-CREATE attributes .....                                      | 55  |
| Table 68 - Basic Film Box Status Codes .....   | 56  |
| Table 69 - Basic Grayscale Image Box N-SET attributes .....                              | 57  |
| Table 70 - Basic Grayscale Image Box Status Codes .....                                  | 57  |
| Table 71 - Basic Color Image Box N-SET attributes .....                                  | 57  |
| Table 72 - Basic Color Image Box Status Codes .....                                      | 58  |
| Table 73 - Presentation LUT N-CREATE attribute .....                                     | 58  |
| Table 74 - Presentation LUT Status Codes .....   | 58  |
| Table 75 - Used Printer N-EVENT Report attributes .....                                  | 59  |
| Table 76 - Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes .....               | 59  |
| Table 77 - Used Print Job N-EVENT Report attributes .....                                | 59  |
| Table 78 - Presentation Context Table "Show Device Status" .....                         | 60  |
| Table 79 - Used Printer N-EVENT Report attributes .....                                  | 60  |
| Table 80 - Mandatory Printer N-GET-RSP, N-EVENT-REPORT-RQ attributes .....               | 61  |
| Table 81: Parameter List .....   | 63  |
| Table 82: Implementation Class/Version Name - Media Interchange .....                    | 67  |
| Table 83: Media - Application Profiles and Real-World Activities .....                   | 67  |
| Table 84: Transfer Syntaxes for STD-GEN-DVD-J2K and STD-GEN-USB-J2K .....                | 68  |
| Table 85: Private SOP Classes and Transfer Syntaxes for Augmented Media Profiles .....   | 69  |
| Table 86: Single-Byte Character Sets without Code Extension .....                        | 70  |
| Table 87: Single-Byte Characters Sets with Code Extension .....                          | 70  |
| Table 88: Multi-Byte Character Sets without Code Extension .....                         | 71  |
| Table 89: Multi-Byte Character Sets with Code Extension .....                            | 71  |
| Table 90: Application Level Confidentiality Profile Attributes (standard tags) .....     | 73  |
| Table 91: Application Level Confidentiality Profile Attributes (private tags) .....      | 78  |
| Table 92: Private Data Element Dictionary .....  | 111 |
| Table 93: Common Printer Status Information .....  | 113 |
| Table 94: Additional Printer Status Information .....                                    | 116 |

## Annex B: Table of Figures

|   |           |
|---|-----------|
| <b>Figure 1: syngo® MR DICOM Data Flow diagram – Acquisition Workflow .....</b>   | <b>11</b> |
| <b>Figure 2: syngo® MR DICOM Data Flow diagram - Printing .....</b>               | <b>12</b> |
| Figure 3: Sequence Diagram for Real World Activities - System Configuration ..... | 15        |
| Figure 4: Sequence Diagram for Real World Activities - Acquisition workflow ..... | 16        |
| <b>Figure 5: Sequence Diagram for Real World Activities - Printing .....</b>      | <b>17</b> |
| Figure 6: Media Interchange Application Data Flow Diagram.....                    | 65        |
| Figure 7: Sequence diagram – Media creation .....                                 | 66        |

**Siemens Healthcare Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone +49 9131 84-0  
[siemens.com/healthcare](https://www.siemens.com/healthcare)

HRT HOOD05162003085912 | © Siemens Healthcare GmbH, 05. 2020

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.