Siemens Healthcare Diagnostics, a global leader in clinical diagnostics, provides healthcare professionals in hospital, reference, and physician office laboratories and point-of-care settings with the vital information required to accurately diagnose, treat, and monitor patients. Our innovative portfolio of performance-driven solutions and personalized customer care combine to streamline workflow, enhance operational efficiency, and support improved patient outcomes.

CA 19-9, CA 125II, and CA 15-3 are trademarks of Fujirebio.

ADVIA Centaur, Dimension Vista, LOCI, V-LYTE, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Order No. A91DX-CAI-121224-GC1-4A00 10-2012 | All rights reserved © 2012 Siemens Healthcare Diagnostics Inc.

Global Siemens Headquarters Siemens AG Wittelsbacherplatz 2 80333 Muenchen Germany Global Siemens Healthcare Headquarters Siemens AG

Siemens AG Healthcare Sector Henkestrasse 127 91052 Erlangen Telephone: +49 9131 84-0

Germany

www.siemens.com/healthcare

Global Division
Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
www.siemens.com/diagnostics

www.siemens.com/diagnostics

Dimension Vista LOCI CA 125II,
CA 15-3, and CA 19-9

Assay Specifications

SIEMENS

www.siemens.com/diagnostics Answers for life.





Dimension Vista LOCI CA 125II, CA 15-3, and CA 19-9 Assay Specifications

The Siemens Healthcare Diagnostics Dimension Vista® LOCI® CA 125II,™ CA 15-3, and CA 19-9 assays are homogeneous, sandwich chemiluminescent immunoassays based on LOCI advanced technology.

Outstanding Assay Performance

- Excellent precision to ensure accurate monitoring
 (LOCI CA 125II: 2.2%-3.8% CV; LOCI CA 15-3: 2.4%-3.5% CV; LOCI CA 19-9: 2.4%-8.9% CV)
- Broad dynamic assay range
- (LOCI CA 125II: 1.5-1,000 U/mL; LOCI CA 15-3: 1.0-300 U/mL; LOCI CA 19-9: 2.0-1,000 U/mL)
- Rapid assay kinetics
- (LOCI CA 125II: 21 minutes; LOCI CA 15-3: 16 minutes; LOCI CA 19-9: 10 minutes)

Clinical Utility

- CA 125II: as an aid in monitoring disease progress or response to therapy or for recurrent or residual disease for patients with epithelial ovarian cancer
- CA 15-3: as an aid in the management of previously treated stage II and III breast cancer patients and for monitoring response to therapy in metastatic breast cancer patients
- CA 19-9: as an aid in managing patients diagnosed with cancers of the exocrine pancreas

Dimension Vista System—Intelligence at Work

- Ultra-integration—four technologies, including best-in-class LOCI advanced chemiluminescence, photometry, nephelometry, and V-LYTE® integrated multisensor technology
- LOCI advanced chemiluminescence—the only homogeneous chemiluminescent technology
- Onboard automation—increased efficiency, simplicity, and convenience for your laboratory

LOCI CA 125II

Changes in CA 125II concentrations and in disease status were analyzed on a per visit basis. Patients were categorized as Active/Progressive, Responding, Stable, or No Evidence of Disease (NED) by the attending physician based on the clinical information (medical imaging, physical examination, and other clinical investigations). All 75 patient sets were analyzed to determine the change in disease status per sequential pair (n = 255). Table 1 shows the distribution of results when compared to the disease status.


Table 1. Disease State Frequency Using the Dimension Vista LOCI CA 125II Assay

| | Change in Disease State | | | | | | |
|-----------------------|-------------------------|------------|------------|------------|--------------|--|--|
| Change in CA 125II | Responding n (%) | , , , | | | | | |
| >70.1% Increase | 3 (1.2%) | 17 (6.8%) | 4 (1.6%) | 42 (16.5%) | 66 (25.9%) | | |
| No Significant Change | 22 (8.6%) | 45 (17.7%) | 54 (21.2%) | 29 (11.4%) | 150 (58.8%) | | |
| >70.1% Decrease | 11 (4.3%) | 13 (5.1%) | 12 (4.7%) | 3 (1.2%) | 39 (15.3%) | | |
| Total | 36 (14.1%) | 75 (29.4%) | 70 (27.5%) | 74 (29.0%) | 255 (100.0%) | | |

Per patient visit clinical performance results for the Dimension Vista LOCI CA 125II assay and predicate devices are given in Table 2. In this evaluation, disease status was classified as "Progression" and "No Progression" with "No Progression" consisting of responding, stable, and no evidence of disease.

Table 2. Dimension Vista LOCI CA 125II Value vs. Disease Progression

| | Progression | No Progression | Total | |
|----------------------|-------------|----------------|----------------------------------|--|
| >70.1% Increase | 42 | 24 | 66 | |
| ≤70.1% Increase | 32 | 157 | 189 | |
| Total | otal 74 | | 255 | |
| | | Estimate | Exact 95% Confidence Interval | |
| Total Concordance | | 78.0% | (72.5%-83.0%) | |
| Positive Cor | ncordance | 56.8% | (44.7%-68.2%) | |
| Negative Concordance | | 86.7% | (80.9%-91.3%) | |



LOCI CA 15-3

Changes in CA 15-3 concentrations and in disease status were analyzed on a per visit basis. Patients were categorized as Active/Progressive, Responding, Stable, or No Evidence of Disease (NED) by the attending physician based on the clinical information (medical imaging, physical examination, and other clinical investigations). All 75 patient sets were analyzed to determine the change in disease status per sequential pair (n = 258). Table 3 shows the distribution of results when compared to the disease status.


Table 3. Disease State Frequency Using the Dimension Vista LOCI CA 15-3 Assay

| | Change in Disease State | | | | |
|-----------------------|-------------------------|-----------------|------------------------------------|----------------------|--------------|
| Change in CA 15-3 | Responding n (%) | Stable n (%) | No Evidence of Disease n (%) | Progression n (%) | Total |
| > 27.9% Increase | 4 (1.5%) | 27 (10.5%) | 3 (1.2%) | 64 (26.8%) | 98 (38%) |
| No Significant Change | 10 (3.9%) | 42 (16.3%) | 37 (14.3%) | 17 (6.6%) | 106 (41%) |
| >27.9% Decrease | 8 (3.1%) | 28 (10.9%) | 3 (1.2%) | 15 (5.8%) | 54 (21%) |
| Total | 22 (8.5%) | 97 (37.7%) | 43 (16.7%) | 96 (37.2%) | 258 (100.0%) |

Per patient visit clinical performance results for the Dimension Vista LOCI CA 15-3 assay and predicate devices are given in Table 4. In this evaluation, disease status was classified as "Progression" and "No Progression" with "No Progression" consisting of responding, stable, and no evidence of disease.

Table 4. Dimension Vista LOCI CA 15-3 Value vs. Disease Progression

| | Progression | No Progression | Total | |
|------------------|---------------------|----------------|----------------------------------|--|
| > 27.9% Increase | > 27.9% Increase 64 | | 98 | |
| ≤ 27.9% Increase | 32 | 128 | 160 | |
| Total 96 | | 162 | 258 | |
| | | Estimate | Exact 95% Confidence Interval | |
| Total Conco | rdance | 74.4% | (68.6%-79.6%) | |
| Positive Cor | ncordance | 66.7% | (56.3%-76.0%) | |
| Negative Co | oncordance | 79.0% | (71.9%-85.0%) | |



 $\mathbf{4}$

LOCI CA 19-9

Changes in CA 19-9 concentrations and in disease status were analyzed on a per visit basis. Patients were categorized as Active/Progressive, Responding, Stable, or No Evidence of Disease (NED) by the attending physician based on the clinical information (medical imaging, physical examination, and other clinical investigations). All 72 patient sets were analyzed to determine the change in disease status per sequential pair (n = 189). Table 5 shows the distribution of results when compared to the disease status.


Table 5. Disease State Frequency Using the Dimension Vista LOCI CA 19-9 Assay

| | Change in Disease State | | | | |
|-----------------------|-------------------------|-----------------|------------------------------------|----------------------|--------------|
| Change in CA 19-9 | Responding n (%) | Stable n (%) | No Evidence of Disease n (%) | Progression n (%) | Total |
| >84.7% Increase | 3 (1.6%) | 10 (5.3%) | 0 (0.0%) | 14 (7.4%) | 27 (14.3%) |
| No Significant Change | 44 (23.3%) | 50 (26.5%) | 8 (4.2%) | 59 (31.2%) | 161 (81.5%) |
| >84.7% Decrease | 0 (0.0%) | 0 (0.0%) | 1 (0.5%) | 0 (0.0%) | 1 (3.2%) |
| Total | 47 (24.9%) | 60 (31.8%) | 9 (4.7%) | 73 (38.6%) | 189 (100.0%) |

Per patient visit clinical performance results for the Dimension Vista LOCI CA 19-9 Assay and predicate devices are given in Table 6. In this evaluation, disease status was classified as "Progression" and "No Progression" with "No Progression" consisting of responding, stable, and no evidence of disease.

Table 6. Dimension Vista LOCI CA 19-9 Value vs. Disease Progression

| | Progression | No Progression | Total | |
|----------------------|-------------|----------------|--------------------------------|--|
| > 84.7% Increase | 14 | 13 | 27 | |
| ≤ 84.7% Increase | 59 | 103 | 162 | |
| Total | 73 | 116 | 189 | |
| | | Estimate | Exact 95% Confidence Limits | |
| Total Concordance | | 61.9% | (54.6%-68.9%) | |
| Positive Cor | ncordance | 19.2% | (10.9%-30.1%) | |
| Negative Concordance | | 88.8% | (81.6%-93.3%) | |



Performance Summary

LOCI CA Assays Performance Summary

| | Sample Type | Sample Volume | Assay Range | Limit of Detection | Cutoff | Calibration Interval | Onboard Stability |
|---------------|----------------|------------------|-------------------|-----------------------|---------|-------------------------|----------------------|
| LOCI CA 125II | Serum/Plasma | 5 μL | 1.5–1,000 U/mL | 1.5 U/mL | 35 U/mL | 21 days | 30 days |
| LOCI CA 15-3 | Serum/Plasma | 1 μL | 1.0-300 U/mL | 1.0 U/mL | 35 U/mL | 30 days | 30 days |
| LOCI CA 19-9 | Serum/Plasma | 4 μL | 2.0-1,000 U/mL | 2.0 U/mL | 37 U/mL | 21 days | 30 days |

| Ordering Information | | | | |
|----------------------|--------------------------------|--------------|--|--|
| Catalog No. | Description | Contents | | |
| K6455 | CA 125 Flex® reagent cartridge | 160 tests | | |
| K6456 | CA 15-3 Flex reagent cartridge | 160 tests | | |
| K6457 | CA 19-9 Flex reagent cartridge | 160 tests | | |
| KC604 | LOCI 6 CAL | 2 x 5 levels | | |
| KC605 | LOCI 7 CAL | 2 x 5 levels | | |

LOCI CA 125II requires LOCI 6 CAL. LOCI CA 15-3 and LOCI CA 19-9 require LOCI 7 CAL.



6