


Characterization of the solitary pulmonary nodule with PET imaging

While conventional imaging establishes the presence of solitary pulmonary nodules (SPNs), invasive procedures that involve risk may be required to characterize the lesions.



PET/CT utilization as a diagnostic tool could reduce the number of unnecessary biopsies or thoracotomies on benign SPNs.¹

“Nodules that were classified as indeterminate on CT were correctly characterized on PET in over 80% of the cases.”

“PET should be obtained in the diagnostic work-up of patients with SPN.” — Fletcher et al.^{2,3}

PET/CT characterizes SPNs with diagnostic accuracy.

Diagnostic accuracy of PET in SPN characterization^{4,5}

	Sensitivity	Specificity	NPV	PPV	Accuracy
PET/CT	97%	85%	92%	93%	92%

PET enables the physician to make a more informed clinical decision.

American College of Chest Physicians recommends: In patients with low-to-moderate pretest probability of malignancy (5 to 65%) and an indeterminate SPN that measures at least 8 mm in diameter, we suggest that functional imaging, preferably with positron emission tomography (PET), should be performed to characterize the nodule.⁶

PET/CT not only has the potential for uncovering malignancy, but also has other roles in patient management, such as staging, evaluation of a response to therapy and detection of recurrence.

Medicare recognizes the utility of PET and PET/CT in solitary pulmonary nodule (SPN).⁷

Initial Treatment Strategy

Medicare covers PET scans for characterization of a single pulmonary nodule under Initial Treatment Strategy of Lung Cancer.

- To determine whether or not the beneficiary is an appropriate candidate for an invasive diagnostic or therapeutic procedure
- To determine the optimal anatomic location for an invasive procedure
- To determine the anatomic extent of tumor when the recommended anti-tumor treatment reasonably depends on the extent of the tumor.


The ordering physician is responsible for documenting the medical necessity of the PET scan and that it meets these criteria.

Society of Nuclear Medicine and Molecular Imaging (SNMMI) Center of Excellence PET PROS Diagnosis of Pulmonary Nodules:

CMS Coverage Policy

PET and PET/CT are approved by the Centers for Medicare and Medicaid Services (CMS) for characterization of solitary pulmonary nodules not exceeding 4 cm to determine the likelihood of malignancy. Claims should include evidence of the initial detection of a primary lung nodule, usually by computed tomography.^{7,8}

PET/CT helps differentiate between benign and malignant solitary pulmonary nodules



Data courtesy of University of Tennessee Medical Center, Knoxville, Tennessee, USA

History

A 58-year-old patient developed bronchitis with onset of shortness of breath. Non-smoker, history of secondhand smoke exposure from mother, who died of small cell lung cancer.

A chest X-ray revealed abnormality in the right apex that warranted a CT scan.

CT demonstrated 2.0 x 1.5 cm RUL nodule.

PET/CT findings and pathology

Further evaluation with PET/CT revealed a nodule as hypermetabolic concerning for malignancy, no evidence of distant disease.

A biopsy was then performed resulting in confirmation of suspicion of malignancy, adenocarcinoma in origin.

Overall evaluation

PET/CT helps in proper clinical management. Results of PET/CT scan revealed patient at high risk of malignancy and therefore required pathological evaluation. Biopsy of nodule proved to be malignant through histopathology and consistent with PET findings. SPN was then surgically resected.

The PET/CT scan also revealed no distant disease eliminating the need for invasive regimen of chemo and/or radiation therapy.

¹ Sim YT, et al. Imaging of solitary pulmonary nodule—a clinical review. *Quant Imaging Med Surg*. 2013 Dec;3(6): 316-326.

² Fletcher JW, et al. A comparison of the diagnostic accuracy of ¹⁸F-FDG PET and CT in the characterization of solitary pulmonary nodules. *J Nucl Med*. 2008 Feb;49(2): 179-85.

³ Fletcher JW, et al. Recommendations on the use of ¹⁸F-FDG PET in oncology. *J Nucl Med*. 2008 Mar;49(3): 480-508.

⁴ Ambrosini V, et al. PET/CT imaging in different types of lung cancer: an overview. *Eur J Radiol*. 2012 May;81(5): 988–1001.

⁵ Kim SK, et al. Accuracy of PET/CT in characterization of solitary pulmonary lesions. *J Nucl Med*. 2007 Feb;48(2): 214-20.

⁶ Gould MK, et al. Evaluation of individuals with pulmonary nodules: when is it lung cancer? Diagnosis and management of lung cancer, 3rd ed: American College of Chest Physicians evidence-based clinical practice guidelines. *Chest*. 2013 May;143(5 Suppl):e93S-120S.

⁷ CMS Publication 100-03, Medicare National Coverage Determinations Manual, Chapter 1, Part 4, Section 220.6). Available at http://www.cms.hhs.gov/manuals/downloads/ncd103c1_part4.pdf

⁸ SNM PET PROS. Available at http://interactive.snm.org/docs/PET_PROS/DiagnosisPulmonaryNodules.pdf

Trademarks and service marks used in this material are property of Siemens Healthcare GmbH. All other company, brand, product, and service names may be trademarks or registered trademarks of their respective holders.

Please contact your local Siemens Healthineers sales representative for the most current information or contact one of the addresses listed below. Note: Original images always lose a certain amount of detail when reproduced.

“Siemens Healthineers” is considered a brand name. Its use is not intended to represent the legal entity to which this product is registered. Please contact your local Siemens organization for further details.

All photographs © 2021 Siemens Healthcare GmbH. All rights reserved.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Published by

PETNET Solutions Inc.
810 Innovation Drive
Knoxville, TN 37932
Phone:+1 877 473 8638
petnetsolutions.com