MAGNETOM Flash

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vol.29

[siemens.com/magnetom-world](https://www.siemens.com/magnetom-world)

Page 3

MAGNETOM Flow を選択

Henrik Michaely, M.D.

Page 8

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| |
|---------------------------------------|
| Tricia Thompson |
| Siemens Healthcare Pty Ltd. Australia |

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| |
|--|
| Xinyu Wang, M.D. ¹ ; Yuehong Liu, M.D. ¹ ; Yichen Tang, M.D. ¹ ; Chen Zhang, Ph.D. ² ; Alto Stemmer ³ ; Qi Yang, M.D., Ph.D. ^{1,4,5} |
| ¹ Department of Radiology, Beijing Chao-yang Hospital, Capital Medical University, Beijing, China |
| ² MR Research Collaboration, Siemens Healthineers, Beijing, China |
| ³ Siemens Healthineers, Erlangen, Germany |
| ⁴ Laboratory for Clinical Medicine, Capital Medical University, Beijing, China |
| ⁵ Key Laboratory of Medical Engineering for Cardiovascular Disease, Ministry of Education, Beijing, China |

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| |
|--|
| Hanns-Christian Breit, M.D. ¹ ; Jan Vosschenrich, M.D. ¹ ; Martin Clauss, M.D. ^{2,3} ; Markus M. Obmann, M.D. ¹ ; Michael Bach, Ph.D. ¹ ; Dorothee Harder, M.D. ¹ ; Ricardo Donners, M.D. ¹ |
| ¹ Department of Radiology, University Hospital Basel, University of Basel, Switzerland |
| ² Center for Musculoskeletal Infections, University Hospital Basel, University of Basel, Switzerland |
| ³ Department for Orthopedics and Trauma Surgery, University Hospital Basel, University of Basel, Switzerland |

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on-time performance と patient experience

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Contact

Associate Professor Susie Y. Huang, M.D., Ph.D.
Radiologist in the Division of Neuroradiology
Department of Radiology
Massachusetts General Hospital
Harvard Medical School
55 Fruit St
Boston MA 02114
USA
susie.huang@mgh.harvard.edu

Alexander Johannes Herold, M.D.
Research Fellow in Radiology
Department of Radiology
Massachusetts General Hospital
Harvard Medical School
55 Fruit St
Boston MA 02114
USA
AHEROLD@mgh.harvard.edu

References

1 Smith-Bindman R, Kwan ML, Marlow EC, Theis MK, Bolch W, Cheng SY, et al. Trends in Use of Medical Imaging in US Health Care Systems and in Ontario, Canada, 2000-2016. JAMA.2019;322(9):843–856.

2 Lin DJ, Doshi AM, Fritz J, Recht MP. Designing Clinical MRI for Enhanced Workflow and Value. J Magn Reson Imaging. 2024;60(1):29–39.

3 Recht MP, Block KT, Chandarana H, Friedland J, Mullholland T, Teahan D, et al. Optimization of MRI Turnaround Times Through the Use of Dockable Tables and Innovative Architectural Design Strategies. AJR Am J Roentgenol. 2019;212(4):855–858.

4 Lang M, Lo WC, Sharp A, Hartmann SP, Pianykh OS, Melski LM, et al. Improving Workflow Efficiency at an Outpatient MRI Imaging Facility: A Case Study. J Am Coll Radiol. 2024:S1546–1440(24)00531–3.

5 Del Grande F, Rashidi A, Luna R, Delcogliano M, Stern SE, Dalili D, et al. Five-Minute Five-Sequence Knee MRI Using Combined Simultaneous Multislice and Parallel Imaging Acceleration: Comparison with 10-Minute Parallel Imaging Knee MRI. Radiology. 2021;299(3):635–646.

6 Fritz J, Fritz B, Zhang J, Thawait GK, Joshi DH, Pan L, et al. Simultaneous Multislice Accelerated Turbo Spin Echo Magnetic Resonance Imaging: Comparison and Combination With In-Plane Parallel Imaging Acceleration for High-Resolution Magnetic Resonance Imaging of the Knee. Invest Radiol. 2017;52(9):529–537.

7 Hammernik K, Klatzer T, Kobler E, Recht MP, Sodickson DK, Pock T, et al. Learning a variational network for reconstruction of accelerated MRI data. Magn Reson Med. 2018;79(6):3055–3071.

8 Lang M, Cartmell S, Tabari A, Briggs D, Pianykh O, Kirsch J, et al. Evaluation of the Aggregated Time Savings in Adopting Fast Brain MRI Techniques for Outpatient Brain MRI. Acad Radiol. 2023;30(2):341–348.

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¹Department of Radiology, University Hospital Basel, University of Basel, Switzerland
²Center for Musculoskeletal Infections, University Hospital Basel, University of Basel, Switzerland
³Department for Orthopedics and Trauma Surgery, University Hospital Basel, University of Basel, Switzerland

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References

- 1 Runge VM, Heverhagen JT. Advocating the Development of Next-Generation, Advanced-Design Low-Field Magnetic Resonance Systems. Invest Radiol. 2020;55(12):747–753.
- 2 Vosschenrich J, Breit HC, Bach M, Merkle EM. Ökonomische Aspekte der Niederfeld-Magnetresonanztomographie : Anschaffung, Installation und Unterhaltskosten von 0,55 T-Geräten [Economic aspects of low-field magnetic resonance imaging: Acquisition, installation, and maintenance costs of 0.55 T systems]. Radiologe. 2022;62(5):400–404. German. Epub Mar 29
- 3 Rusche T, Vosschenrich J, Winkel DJ, Donners R, Segeroth M, Bach M, et al. More Space, Less Noise—New-generation Low-Field Magnetic Resonance Imaging Systems Can Improve Patient Comfort: A Prospective 0.55 T–1.5 T-Scanner Comparison. J Clin Med. 2022;11(22):6705.
- 4 Farahani K, Sinha U, Sinha S, Chiu LC, Lufkin RB. Effect of field strength on susceptibility artifacts in magnetic resonance imaging. Comput Med Imaging Graph. 1990;14(6):409–413.
- 5 Bernstein MA, Huston III J, Ward HA. Imaging artifacts at 3.0T. J Magn Reson Imaging. 2006;24(4):735–746.
- 6 Abdelaal MS, Restrepo C, Sharkey PF. Global Perspectives on Arthroplasty of Hip and Knee Joints. Orthop Clin North Am. 2020;51(2):169–176.
- 7 Jungmann PM, Agten CA, Pfirrmann CW, Sutter R. Advances in MRI around metal. J Magn Reson Imaging. 2017;46(4):972–991.
- 8 Khodarahmi I, Brinkmann IM, Lin DJ, Bruno M, Johnson PM, Knoll F, et al. New-generation low-field magnetic resonance imaging of hip arthroplasty implants using slice encoding for metal artifact correction: first in vitro experience at 0.55 T and comparison with 1.5 T. Invest Radiol. 2022;57(8):517–526.

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References

- Ogawa R, Kido T, Shiraishi Y, Yagi Y, Yoon S, Wetzl J, et al. Neural network-based fully automated cardiac resting phase detection algorithm compared with manual detection in patients. Acta Radiologica Open. 2022;11(10):20584601221137772.
- Wood G, Pedersen A, Kunze K, Neji R, Hajhosseiny R, Wetzl J, et al. Automated detection of cardiac rest period for trigger delay calculation for image-based navigator coronary magnetic resonance angiography. J Cardiovasc Magn Reson. 2023;25(1):52.
- McDermott S, Wetzl J, Schmidt M, Benbow M, Yoon S, Geppert C, et al. AI-based Cardiac Scan Automation: A Prospective Comparison of Highly Automated Scan Workflows in 32 Patients. J Cardiovasc Magn Reson. 2024;26(1):100674.
- Haenel A, Wetzl J, Yoon S, Schmidt M, Weißgerber AB, Barkhausen J, et al. Analysis of an AI-based prototype for user-independent automated cardiac MRI (CMR) scanning. Presented at: European Congress of Radiology; March 1–5, 2023; Vienna, Austria.
- Glessgen C, Crowe L, Wetzl J, Schmidt M, Yoon S, Poletti PA, et al. Fully Automated vs. Manual Cardiac MRI: a Single-Center Prospective Evaluation of Quality and Scan Times. European Radiology [in revision].
- Yoon S, Preuhs E, Schmidt M, Forman C, Chitiboi T, Sharma P, et al. Automated Cardiac Resting Phase Detection Targeted on the Right Coronary Artery. Machine Learning for Biomedical Imaging. arXiv. 2023;2:2109.02342.
- Yoon S, Schmidt M, Rick M, Chitiboi T, Sharma P, Emrich T, et al. Validation of a deep learning based automated myocardial inversion time selection for late gadolinium enhancement imaging in a prospective study. Presentation at International Society for Magnetic Resonance in Medicine (ISMRM) 29th Annual Meeting & Exhibition, Virtual.

From left to right: Jens Wetzl, Manuela Rick, Seung Su Yoon, Michaela Schmidt, Miriam Van de Stadt-Lagemaat, Dominik Rupp

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¹Department of Radiology, Beijing Chao-yang Hospital, Capital Medical University, Beijing, China
²MR Research Collaboration, Siemens Healthineers, Beijing, China
³Siemens Healthineers, Erlangen, Germany
⁴Laboratory for Clinical Medicine, Capital Medical University, Beijing, China
⁵Key Laboratory of Medical Engineering for Cardiovascular Disease, Ministry of Education, Beijing, China

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
| The following scan parameters were used: | |
|--|---|
| Repetition time (TR) | 700 ms |
| Echo time (TE) | 11 ms |
| FOV | 400 × 400 mm ² |
| Matrix | 352 × 387 |
| Slices per slab | 256 |
| Voxel size | 0.6 × 0.6 × 0.7 mm ³ |
| Blood suppression | DANTE |
| Fat suppression | Water excitation |
| Acceleration mode | Generalized autocalibrating partially parallel acquisition (GRAPPA) |
| Acceleration factor | 4 |
| Bandwidth | 1184 Hz/pixel |
| Number of DANTE pulses | 200 |
| DANTE gradient duration | 1 ms |
| DANTE gradient amplitude | 17 mT/m |
| DANTE radiofrequency pulse flip angle | 15° |
| Scan time for each station | 4 minutes |
| Total scan time | 12 minutes |

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
References

1 Mendichovszky IA, Priest AN, Bowden DJ, Hunter S, Joubert I, Hilborne S, et al. Combined MR direct thrombus imaging and non-contrast magnetic resonance venography reveal the evolution of deep vein thrombosis: a feasibility study. Eur Radiol. 2017;27(6):2326–2332.

2 Zhuang G, Tang C, He X, Liang J, He Z, Ye Y, et al. DANTE-SPACE: a new technical tool for DVT on 1.5T MRI. Int J Cardiovasc Imaging. 2019;35(12):2231–2237.

3 Chen H, He X, Xie G, Liang J, Ye Y, Deng W, et al. Cardiovascular magnetic resonance black-blood thrombus imaging for the diagnosis of acute deep vein thrombosis at 1.5 Tesla. J Cardiovasc Magn Reson. 2018;20(1):42.

4 Treitl KM, Treitl M, Kooijman-Kurfuerst H, Kammer NN, Coppenrath E, Suderland E, et al. Three-dimensional black-blood T1-weighted turbo spin-echo techniques for the diagnosis of deep vein thrombosis in comparison with contrast-enhanced magnetic

resonance imaging: a pilot study. Invest Radiol. 2015;50(6):401–8.

5 Li L, Miller KL, Jezzard P. DANTE-Prepared Pulse Trains: A Novel Approach to Motion-Sensitized and Motion-Suppressed Quantitative Magnetic Resonance Imaging. Magn Reson Med. 2012;68(5):1423–38.

6 Jeong MJ, Kwon H, Noh M, Ko GY, Gwon DI, Lee JS, et al. Relationship of Lower-extremity Deep Venous Thrombosis Density at CT Venography to Acute Pulmonary Embolism and the Risk of Postthrombotic Syndrome. Radiology. 2019;293(3):687–694.

7 Xie G, Chen H, He X, Liang J, Deng W, He Z, et al. Black-blood thrombus imaging (BTI): a contrast-free cardiovascular magnetic resonance approach for the diagnosis of non-acute deep vein thrombosis. J Cardiovasc Magn Reson. 2017;19(1):4.

8 Saha P, Andia ME, Modarai B, Blume U, Humphries J, Patel AS, et al. Magnetic resonance T1 relaxation time of venous thrombus is determined by iron processing and predicts susceptibility to lysis. Circulation. 2013;128(7):729–736.

Contact

Professor Qi Yang, M.D., Ph.D.
Chair, Department of Radiology
Beijing Chaoyang Hospital
Capital Medical University
No. 8 Gongti South Road
Chaoyang District
Beijing 100020
China
yangyangqiqi@gmail.com