



The INNOVANCE Anti-Xa Assay*

Increasing Safety and Efficiency in Heparin Therapy Management

siemens-healthineers.com

*Valid also for INNOVANCE Heparin assay as both products are seen to be equivalent in respect of performance of the reagents itself as well as the determination of unfractionated heparin (UFH) and low molecular weight heparin (LMWH).

The INNOVANCE Anti-Xa Assay—Increasing Safety and Efficiency in Heparin Therapy Management

Heparin and its clinical applications

Heparin is a highly sulfated mucopolysaccharide typically prepared from porcine or bovine gut mucosa and used clinically as an anticoagulant for intravenous or subcutaneous application. Heparins are heterogeneous with respect to molecular size, anticoagulant activity, and pharmacokinetic properties.

Heparin is available for clinical use in two common forms: natural, unfractionated heparin (UFH), with a mean molecular size of about 15,000 kDa, and low-molecular weight heparin (LMWH), with a typical mean molecular size of 5000 kDa. LMWH can be derived from UFH by chemical or enzymatic depolymerization.


Both UFH and LMWH require antithrombin as a cofactor for their anticoagulant activity. Both heparin types accelerate FXa (and FIXa, FXIa, and FXIIa) inactivation, while UFH also inactivates thrombin (anti-IIa activity).

The minimum length of the polysaccharide chain required for anti-Xa activity is 5 units, whereas anti-IIa activity requires polysaccharide chains of at least 18 units. Such longer polysaccharide chains are present in significant amounts only in UFH.

Despite being more expensive, LMWH has replaced UFH in many clinical applications because it offers certain advantages:

- Less heparin-induced thrombocytopenia (HIT), a severe complication of heparin therapy with a high mortality rate
- Therapeutic use by subcutaneous rather than intravenous application
- No routine monitoring required
- Less bone destruction with long-term use


Nevertheless, in certain indications, UFH is the preferred heparin, e.g., for hemodialysis or extra-corporal circulation, or patients with severe kidney disease.



Mechanism of action of unfractionated and low-molecular-weight heparin

To inactivate thrombin, UFH forms a ternary complex with antithrombin and thrombin. Because of their lower molecular weight, LMWH species are unable to form the ternary complexes with antithrombin and thrombin. Thus, LMWHs produce their anticoagulant effect mainly by inhibiting factor Xa.

Molecular weight distribution of LMWH and UFH



Characteristics of LMWH and UFH

	LMWH (e.g., enoxaparin, dalteparin, nadroparin)	UFH
Molecular size	Mean: 5000 kDa Range: 2000–9000 kDa	Mean: 15,000 kDa Range: 3000–90,000 kDa
Molecular targets	Reduced anti-IIa activity relative to anti-Xa activity	Anti-IIa and anti-Xa activity
Plasma half-life	About 3–6 hours	Highly variable and dose-dependent 30–150 min
Route of application	Subcutaneous (once or twice daily)	Intravenous (continuous infusion) or subcutaneous (twice daily)
Monitoring	More predictable anticoagulant response; no routine monitoring required Anti-Xa assay only	Requires monitoring in therapeutic concentration APTT or anti-Xa assay
Clearance	Renal clearance	Larger molecules are cleared more rapidly than smaller; only partial renal clearance
Immunogenicity	Reduced complex formation, resulting in less HIT antibodies	Complex formation with PF4 (platelet factor 4), which may induce HIT antibodies
Side effects	Lower risk of osteoporosis	Long-term use associated with bone loss

Heparin therapy modalities

UFH is the preferred drug in patients with severe kidney disease, and for situations where very high anticoagulant levels are required. Neutralization of UFH is possible by complexation with protamine. Applications for UFH include:

- Prophylactic regimen for prevention of venous thromboembolism (VTE), although UFH has been nearly completely replaced by LMWH for this purpose. Typically administered twice daily subcutaneously. No monitoring is required.
- Therapeutic regimen for therapy of VTE or special situations with high coagulation triggering activity, such as hemodialysis or extra-corporal circulation. Administered by continuous intravenous infusion. Routine monitoring by APTT or anti-Xa assay is required.

LMWH is the drug of choice for most heparin-therapy applications. There is a risk of accumulation in kidney-disease patients, and no real antidote is available; protamine is only partially effective. Applications for LMWH include:

- Prophylactic regimen for prevention of VTE. LMWH has become the standard routine prophylaxis for VTE. Typically administered once daily subcutaneously. No monitoring is required.
- Therapeutic dosing for therapy of VTE. Dosing must be adjusted for body weight. Administered once or twice daily subcutaneously. Routine monitoring is not required, but for specific indications, such as obese patients, infants, pregnancy, or patients with renal impairment; monitoring is recommended and can be performed via anti-Xa assay only.

A very small and special heparinoid is fondaparinux. This synthetic pentasaccharide analog is the smallest polysaccharide with anticoagulant activity. Fondaparinux acts only through anti-Xa activity potentiation, has an increased half-life versus LMWH, and is exclusively cleared by the kidneys. Monitoring is not recommended on a routine basis, but, if indicated, a fondaparinux-specific anti-Xa assay is the method of choice. When monitoring is indicated, fondaparinux is required for assay calibration.

Target population for anti-Xa testing

1. All patients under LMWH therapy who should be monitored, such as:
 - Pregnant women
 - Infants and children
 - Patients with renal impairment
 - Patients with high risk for bleeding or recurrences
 - Very obese and very underweight patients circulation, or patients with severe kidney disease.
2. Patients under UFH therapy who show an abnormal baseline APTT, including:
 - Patients with liver cirrhosis
 - Patients in the bridging phase from coumadin therapy
 - Patients with lupus anticoagulants
 - Patients under or having completed thrombolytic therapy
 - Patients with congenital or acquired factor deficiencies
 - Critically ill patients
 - Patients with APTT prolongation due to high C reactive protein (CRP) levels

Recommended therapeutic target ranges¹⁻³

APTT therapeutic ranges must be calibrated versus anti-Xa activity as specified below:

	Target Range for Anti-Xa Activity	Timing for Testing
LMWH	0.6–1.0 IU/mL	At peak level 4 hours after application
UFH	0.3–0.7 IU/mL	4–6 hours after initiation or dose adjustment

For more-specific information, please refer to the prescribing information provided by the heparin manufacturers.

The anti-Xa assay: A better choice for UFH monitoring

The APTT test is affected by many biological factors unrelated to heparin concentration and is not a standardized test. The historical goal of 1.5–2.5 times prolongation of the APTT is considered an antiquated and unsafe practice.

The full effect of the preanalytical, analytical, and biological variables became apparent when it was demonstrated that less than half of the variation in APTT values in patients receiving heparin was explained by differences in heparin concentrations.

Clinical data from the last 10–20 years support a conversion from APTT to anti-Xa monitoring. Advantages of anti-Xa monitoring include:


- A smoother dose-response curve
- More-stable heparin levels during therapy
- Fewer blood samples required
- Fewer dosage adjustments

Increasing data suggest that fewer overall laboratory tests are necessary when the anti-Xa assay is used, due to its smoother dose-response curve. Consequently, although the cost per anti-Xa assay is higher than that of the APTT test, overall institutional costs may be nearly equivalent due to the fewer monitoring tests and dosage adjustments necessary when the anti-Xa assay is used.⁴

Several studies support the use of anti-Xa assays for UFH monitoring:

- Vandiver⁵ reports for the use of an anti-Xa assay-based UFH-monitoring protocol a higher percentage of within-range blood plasma heparin monitoring tests, fewer monitoring tests for the patient to achieve blood plasma monitoring tests within goal range, and fewer dose adjustments compared with a protocol based on blood plasma monitoring using the APTT.
- Liveris⁶ found the anti-Xa assay to correlate better with heparin dosing than ACT or APTT. In pediatric extracorporeal membrane oxygenation, the anti-Xa assay was a more valuable monitor of heparin administration.
- Adatya⁷ observed a high rate (74%) of discordant results between APTT and anti-Xa in patients with left ventricular assist devices. Many of the discordant results were related either to bridging from warfarin therapy or presence of hemolysis, which may falsely elevate APTT.
- Van Roessel⁸ investigated critically ill patients; paired analysis of APTT and anti-Xa was discordant in 63% of the samples. Anti-Xa heparin levels were considered the gold standard; APTT was neither sensitive nor specific to detect under- or overdosing.

Distribution of APTT and anti-Xa levels in critically ill patients⁸



INNOVANCE Anti-Xa assay: Single calibration improves reliability and efficiency

An important feature of the INNOVANCE® Anti-Xa assay from Siemens Healthineers is that only one reference curve and a universal calibrator set are required for measurement of both UFH and LMWH samples. The universal calibrator set is traceable to the international reference preparations for both UFH and LMWH.

Current guidelines and recommendations point to the need for different calibrators for UFH and LMWH. However, with the INNOVANCE Anti-Xa assay, equal units of UFH and LMWH correspond to the same measuring signal, i.e., the same absorbance signal is obtained, irrespective of whether 0.5 IU/mL of UFH or LMWH is present in the sample. This is made possible by the specific, sophisticated assay conditions of the INNOVANCE Anti-Xa assay, allowing the use of a single calibration curve for samples containing either UFH or LMWH.


Note: The guidelines do not include assays that utilize single calibration because no commercial assays with this feature were available at the time the guidelines were established.

Single calibration saves time, eliminates potential errors, and improves patient safety

In many cases, the lab does not receive information specifying which type of heparin a specific patient is receiving. Incorrectly allocating the sample to either a UFH-calibrated or LMWH-calibrated heparin assay most likely results in an erroneous heparin concentration, with all the consequences that may arise from such an invalid result. To prevent this from happening, labs often spend considerable time and effort to get this important information.

Instead, with the INNOVANCE Anti-Xa assay, heparin samples can be tested immediately, without having to verify which type of heparin the patient is receiving. Benefits to the lab include the ability to run fewer heparin assays overall and lower costs for calibration materials, as well as substantial time savings for lab personnel and the ability to deliver heparin test results more quickly. But the most important benefit is the improved patient safety that results from eliminating interpretation errors and getting the correct result regardless of which type of heparin is used.

INNOVANCE Anti-Xa Assay: 2nd International Standard for LMWH compared to the 6th International Standard for UFH⁹





References:

1. Garcia DA, Baglin TP, Weitz JI, Samama MM; American College of Chest Physicians. Current international guideline on parenteral anticoagulant. Parenteral anticoagulants: Antithrombotic Therapy and Prevention of Thrombosis, 9th ed: American College of Chest Physicians Evidence-Based Clinical Practice Guidelines. Chest. 2012;141 (2 Suppl):e24S-43S.
2. Kitchen S, Gray E, Mackie I, Baglin T, Makris M; BCSH committee. British guideline on non-coumadin anticoagulants, including heparins: Measurement of non-coumadin anticoagulants and their effects on tests of Haemostasis: Guidance from the British Committee for Standards in Haematology. Br J Haematol. 2014;166(6): 830-41. Free access: <http://onlinelibrary.wiley.com/doi/10.1111/bjh.12975/pdf>
3. Kitchen S. Problems in laboratory monitoring of heparin dosage. Br J Haematol. 2000;111(2):397-406.
4. Vandiver JW, Vondracek TG. Antifactor Xa levels versus activated partial thromboplastin time for monitoring unfractionated heparin. Pharmacotherapy. 2012;32(6):546-58.
5. Vandiver JW, Vondracek TG. A comparative trial of anti-factor Xa levels versus the activated partial thromboplastin time for heparin monitoring. Hosp Pract. 2013;41(2):16-24.
6. Liveris A, Bello RA, Friedmann P, Duffy MA, Manwani D, Killinger JS, Rodriquez D, Weinstein S. Anti-factor Xa assay is a superior correlate of heparin dose than activated partial thromboplastin time or activated clotting time in pediatric extracorporeal membrane oxygenation. Pediatr Crit Care Med. 2014;15(2):e72-9.
7. Adaya S, Uriel N, Yarmohammadi H, Holley CT, Feng A, Roy SS, Reding MT, John R, Eckman P, Zantek ND. Antifactor Xa and activated partial thromboplastin time measurements for heparin monitoring in mechanical circulatory support. JACC Heart Fail. 2015;3(4):314-22.
8. van Roessel S, Middeldorp S, Cheung YW, Zwinderman AH, de Pont AC. Accuracy of aPTT monitoring in critically ill patients treated with unfractionated heparin. Neth J Med. 2014;72(6):305-10. Free access: <http://www.njmonline.nl/getpdf.php?id=1463>
9. Wilkens M, Becker-Scheidemann B, Borchert A, et al. Preliminary performance data of a new assay for the quantitative determination of the activity of unfractionated heparin (UFH) and low molecular weight heparin (LMWH) using a single calibration curve for both types of heparin. Poster at ISTH 2015, Toronto.

At Siemens Healthineers, we pioneer breakthroughs in healthcare. For everyone. Everywhere. By constantly bringing breakthrough innovations to market, we enable healthcare professionals to deliver high-quality care, leading to the best possible outcome for patients.

Our portfolio, spanning from in-vitro and in-vivo diagnostics to image-guided therapy and innovative cancer care, is crucial for clinical decision-making and treatment pathways. With our strengths in patient twinning, precision therapy, as well as digital, data, and artificial intelligence (AI), we are well positioned to take on the biggest challenges in healthcare. We will continue to build on these strengths to help fight the world's most threatening diseases, improving the quality of outcomes, and enabling access to care.

We are a team of 66,000 highly dedicated employees across more than 70 countries passionately pushing the boundaries of what's possible in healthcare to help improve people's lives around the world.

INNOVANCE and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

This white paper is intended for use outside U.S. only (OUS). Display or promotion to U.S. audience is prohibited.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
Phone: +1 914-631-8000