

Non-Contrast Portal Vein Imaging Using *syngo* NATIVE MR Angiography

Marcelo Fernandes Arêas


Siemens Healthineers, Erlangen, Germany

Evaluation of the portal venous system is an important part of abdominal magnetic resonance imaging (MRI) studies, particularly in patients with chronic liver disease, suspected portal vein thrombosis, or portal hypertension, or before interventional and surgical procedures. Contrast-enhanced MR angiography is commonly used for this purpose; however, in daily practice a non-contrast technique is desirable in many situations. These include impaired renal function, repeated follow-up examinations, limited intravenous access, or even allergies. Ultrasound is often the first-line imaging method, but MRI offers a larger field of view and a better anatomical overview [1, 2].

The aim of our protocol is to provide a reproducible, patient-friendly, and clinically useful examination of the

portal venous system without the need for gadolinium. The protocol can be done using a 1.5T or 3T system. The examples below are from two 1.5T systems (MAGNETOM Sola and MAGNETOM Flow) and one 3T system (MAGNETOM Lumina), respectively. We can use either the Contour or Body Matrix coils together with the Spine coil integrated into the table.


We can benefit from BioMatrix Technology integrated into the Spine coil and use the respiratory signal to gate and trigger the data acquisition. You can take the native_trufi3d_tra_resp_trig protocol directly from the Siemens Healthineers protocol tree (Fig. 1). You still have the possibility to use the Resp. control on the Physio card and PACE (Fig. 2).



1 Where to find the native_trufi3d_tra_resp_trig protocol.

Routine	Contrast	Resolution	Geometry	System	Physio	Inline	Sequence
Signal	Cardiac	PACE					
Resp. Control	Gate Gate & Follow Gate Monitor Only Off						
Scout Mode							
Scout Duration	19 s						
Scout TR	100 ms						
Accept Window $\pm$	4,00 mm						
Accept Position (green)	128,00 mm						
Search Window $\pm$	32,00 mm						
Search Position (red)	128,00 mm						
Tracking Factor	0,60						
Chronological Position	Before Echo						
RF Pulse Type	Crossed Pair						
Resp. Motion Adaptation	☒						
Concatenations	1						


2 Respiratory gating options.




3 How to position the slices and the inversion graphics.

Some sequence adjustments are needed. Change the acquisition orientation to coronal. Increasing the slice thickness to 1.5 mm helps to get more signal. The inversion graphics and the inversion time (TI) values are important. We used two inversion graphics positioned diagonally to the liver and perpendicular to the portal

vein (Fig. 3). At 1.5T, we set the TI value at 1500 ms. Depending on the size of the patient, you can increase the slice oversampling. The sequence can be modified to be bright blood with Trufi 3D (Fig. 4), or dark blood with GRE_3d (Fig. 5). Changing the sequence type is easy to do on the sequence tab (Fig. 6).




4 The sequence can be modified to be bright blood with Trufi 3D. Images acquired on a 1.5T MAGNETOM Sola.





- 5 The sequence can be modified to be dark blood with GRE_3d MR. Images acquired on a 1.5T MAGNETOM Sola (5A) and MAGNETOM Flow (5B).

- 6 Changing the sequence between Trufi and GRE.

Routine	Contrast	Resolution	Geometry	System	Physio	Inline	Sequence
Part 1	Part 2	Assistant					
Sequence Name		tfl		Bandwidth		781 Hz/Px	
Dimension		3D		Echo Spacing		3,48 ms	
Sequence Type		Trufi		Asymmetric Echo		Weak	
Excitation		<Gre>		Optimization		Min. TE	
RF Pulse Type		Fast		Define		Shots	
Gradient Mode		Normal		Shots per Slice		1	
Reordering		Linear		Segments		188	



7 The MAGNETOM Lumina 3T protocol.



Study ID: 2aaca4901

Conclusion

Native non-contrast MRI of the portal vein is a useful addition to abdominal MRI at 1.5T. With careful planning, patient coaching, and source-image review, the technique can provide clinically relevant information about portal venous patency and anatomy.

Contact

Marcelo Fernandes Arêas
Siemens Healthineers
SHS DI MR M&S S CVP
Allee am Röthelheimpark 2
91052 Erlangen
Germany
marcelo.arêas@siemens-healthineers.com



References

- Ohno T, Isoda H, Furuta A, Togashi K. Non-contrast-enhanced MR portography and hepatic venography with time-spatial labeling inversion pulses: comparison at 1.5 Tesla and 3 Tesla. Acta Radiol Open [Internet]. 2015;4(5):2058460115584110. Available from: <http://dx.doi.org/10.1177/2058460115584110>
- Lee CH. Portal vein thrombosis on unenhanced MRI: a case series. BJR Case Rep [Internet]. 2022;9(1):20220059. Available from: <http://dx.doi.org/10.1259/bjrcr.20220059>