MAGNETOM Flash

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vol.30

siemens.com/magnetom-world

Page 4

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¹Medical Radiation Physics, Clinical Sciences Lund, Lund University, Lund, Sweden

²Cardiff University Brain Research Imaging Centre (CUBRIC), Cardiff University, Cardiff, United Kingdom

³Siemens Healthcare Ltd, Camberly, United Kingdom

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Bowen Wu¹, Yitong Lu¹, Yueluan Jiang², Xiaoye Wang³, Liang Zhu¹

¹Department of Radiology, Peking Union Medical College Hospital, Beijing, China

²MR Research Collaboration Team, Siemens Healthineers Ltd., Beijing, China

³MR Clinical Marketing, Siemens Healthineers Ltd., Beijing, China

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¹Department of Diagnostic Imaging and Nuclear Medicine, Graduate School of Medicine, Kyoto University, Kyoto, Japan

²Siemens Healthineers, MR Research and Development, Portland, OR, USA

³Siemens Healthineers, MR Research and Development, San Francisco, CA, USA

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¹Department of Radiology, Zhongshan Hospital, Fudan University, Shanghai Institute of Medical Imaging, Shanghai, China

²MR Application Development, Siemens (Shenzhen) Magnetic Resonance Ltd., Shenzhen, China

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¹Department of Radiology, Balgrist University Hospital, Faculty of Medicine, University of Zurich, Switzerland

²Swiss Innovation Hub, Siemens Healthineers AG, Zurich, Switzerland

³Swiss Center for Musculoskeletal Imaging (SCMI), Balgrist Campus, Zurich, Switzerland

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¹School of Biomedical Engineering and Imaging Sciences, King's College London, United Kingdom

²MR Predevelopment, Siemens Healthineers, Erlangen, Germany

³Department of Diagnostic and Interventional Radiology, University Hospital Frankfurt, Frankfurt am Main, Germany

⁴Centre for Medical Physics and Biomedical Engineering, Medical University of Vienna, Austria

⁵Department of Neuroradiology, Heidelberg University Hospital, Heidelberg, Germany

⁶School of Cancer & Pharmaceutical Sciences, King's College London, United Kingdom

⁷Department of Pathology, Hôpital Beaujon, Clichy, France

⁸Department of Radiology, Hôpital Beaujon, Clichy, France

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¹⁰Laboratory for Vascular Translational Science (LVTS), INSERM U1148, Paris, France

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References

- Dyba T, Randi G, Bray F, Martos C, Giusti F, Nicholson N, et al. The European cancer burden in 2020: Incidence and mortality estimates for 40 countries and 25 major cancers. *Eur J Cancer*. 2021;157:308–347.
- Siegel RL, Giaquinto AN, Jemal A. Cancer statistics, 2024. *CA Cancer J Clin*. 2024 Jan-Feb;74(1):12–49. Erratum in: *CA Cancer J Clin*. 2024 Mar-Apr;74(2):203.
- WHO Classification of Tumours Editorial Board. Urinary and Male Genital Tumours. 5th ed, Vol. 8. International Agency for Research on Cancer; 2022.
- James ND, Tannock I, N'Dow J, Feng F, Gillissen S, Ali SA, et al. The Lancet Commission on prostate cancer: planning for the surge in cases. *Lancet*. 2024;403(10437):1683–1722. Erratum in: *Lancet*. 2024;403(10437):1634.
- Weinreb JC, Barentsz JO, Choyke PL, Cornud F, Haider MA, Macura KJ, et al. PI-RADS Prostate Imaging - Reporting and Data System: 2015, Version 2. *Eur Urol*. 2016;69(1):16–40.
- Le Bihan D, Breton E, Lallemand D, Grenier P, Cabanis E, Laval-Jeantet M. MR imaging of intravoxel incoherent motions: application to diffusion and perfusion in neurologic disorders. *Radiology*. 1986;161(2):401–407.
- Weinreb JC, Barentsz JO, Choyke PL, Cornud F, Haider MA, Macura KJ, et al. PI-RADS Prostate Imaging - Reporting and Data System: 2015, Version 2. *Eur Urol*. 2016;69(1):16–40.
- Ahmed HU, El-Shater Bosaily A, Brown LC, Gabe R, Kaplan R, Parmar MK, et al. Diagnostic accuracy of multi-parametric MRI and TRUS biopsy in prostate cancer (PROMIS): a paired validating confirmatory study. *Lancet*. 2017;389(10071):815–822.
- Thompson JE, van Leeuwen PJ, Moses D, Shnier R, Brenner P, Delprado W, et al. The Diagnostic Performance of Multiparametric Magnetic Resonance Imaging to Detect Significant Prostate Cancer. *J Urol*. 2016;195(5):1428–1435.
- Ahmed HU, Arya M, Freeman A, Emberton M. Do low-grade and low-volume prostate cancers bear the hallmarks of malignancy? *Lancet Oncol*. 2012;13(11):e509–e517.
- Johnson DC, Raman SS, Mirak SA, Kwan L, Bajgiran AM, Hsu W, et al. Detection of Individual Prostate Cancer Foci via Multiparametric Magnetic Resonance Imaging. *Eur Urol*. 2019;75(5):712–720.
- Lee MS, Moon MH, Kim YA, Sung CK, Woo H, Jeong H, et al. Is Prostate Imaging Reporting and Data System Version 2 Sufficiently Discovering Clinically Significant Prostate Cancer? Per-Lesion Radiology-Pathology Correlation Study. *AJR Am J Roentgenol*. 2018;211(1):114–120.
- Borofsky S, George AK, Gaur S, Bernardo M, Greer MD, Mertan FV, et al. What Are We Missing? False-Negative Cancers at Multiparametric MR Imaging of the Prostate. *Radiology*. 2018;286(1):186–195.
- Bratan F, Melodelima C, Souchon R, Hoang Dinh A, Mège-Lechevallier F, Cruzet S, et al. How accurate is multiparametric MR imaging in evaluation of prostate cancer volume? *Radiology*. 2015;275(1):144–154.
- Van Houdt PJ, Ghobadi G, Schoots IG, Heijmink SWTPJ, de Jong J, van der Poel HG, et al. Histopathological Features of MRI-Invisible Regions of Prostate Cancer Lesions. *J Magn Reson Imaging*. 2020;51(4):1235–1246.
- Hasler SW, Kallehaug JF, Hansen RH, Samsøe E, Arp DT, Nissen HD, et al. Geometric distortions in clinical MRI sequences for radiotherapy: insights gained from a multicenter investigation. *Acta Oncol*. 2023;62(11):1551–1560.
- Weygand J, Fuller CD, Ibbott GS, Mohamed AS, Ding Y, Yang J, et al. Spatial Precision in Magnetic Resonance Imaging-Guided Radiation Therapy: The Role of Geometric Distortion. *Int J Radiat Oncol Biol Phys*. 2016;95(4):1304–1316.
- Panagiotaki E, Chan RW, Dikaio N, Ahmed HU, O'Callaghan J, Freeman A, et al. Microstructural characterization of normal and malignant human prostate tissue with vascular, extracellular, and restricted diffusion for cytometry in tumours magnetic resonance imaging. *Invest Radiol*. 2015;50(4):218–227.

- Singh S, Rogers H, Kanber B, Clemente J, Pye H, Johnston EW, et al. Avoiding Unnecessary Biopsy after Multiparametric Prostate MRI with VERDICT Analysis: The INNOVATE Study. *Radiology*. 2022;305(3):623–630.
- Hahn EL. Spin echoes. *Physical review*, 1950;80(4):580.
- Stejskal EO, Tanner JE. Spin diffusion measurements: spin echoes in the presence of a time-dependent field gradient. *J Chem Phys*. 1965;42(1):288–292.
- Stehling MK, Turner R, Mansfield P. Echo-planar imaging: magnetic resonance imaging in a fraction of a second. *Science*. 1991;254(5028):43–50.
- Luna A, Ribes R, Soto JA. Diffusion MRI Outside the Brain: A Case-Based Review and Clinical Applications. Springer Berlin, Heidelberg; 2011.
- Gaddamanugu S, Shafaat O, Sotoudeh H, Sarrami AH, Rezaei A, Saadatpour Z, et al. Clinical applications of diffusion-weighted sequence in brain imaging: beyond stroke. *Neuroradiology*. 2022;64(1):15–30.
- Kimmlingen R. Pioneers of connectome gradients. *MAGNETOM Flash*. 2017;68(2):122–136.
- Setsonpop K, Kimmlingen R, Eberlein E, Witzel T, Cohen-Adad J, McNab JA, et al. Pushing the limits of in vivo diffusion MRI for the Human Connectome Project. *Neuroimage*. 2013;80:220–233.
- Huang SY, Witzel T, Keil B, Scholz A, Davids M, Dietz P, et al. Connectome 2.0: Developing the next-generation ultra-high gradient strength human MRI scanner for bridging studies of the micro-, meso- and macro-connectome. *Neuroimage*. 2021;243:118530.
- Siemens Healthineers. MAGNETOM Cima.X. Available at <https://www.siemens-healthineers.com/magnetic-resonance-imaging/3t-mri-scanner/magnetom-cima-x>. Cited Jan 16, 2025.
- Foo T, Graziani D, Hua Y, Bai Y, Ricci J, Vasil C, et al. MAGNUS: An ultra-high efficiency head-only gradient coil for imaging the brain microstructure. *Proc Intl Soc Mag Reson Med*. 2018;26:0839.
- Foo TKF, Tan ET, Vermilyea ME, Hua Y, Fiveland EW, Piel JE, et al. Highly efficient head-only magnetic field insert gradient coil for achieving simultaneous high gradient amplitude and slew rate at 3.0T (MAGNUS) for brain microstructure imaging. *Magn Reson Med*. 2020;83(6):2356–2369.
- Jones DK, Alexander DC, Bowtell R, Cercignani M, Dell'Acqua F, McHugh DJ, et al. Microstructural imaging of the human brain with a 'super-scanner': 10 key advantages of ultra-strong gradients for diffusion MRI. *Neuroimage*. 2018;182:8–38.
- Szczepankiewicz F, Westin CF, Nilsson M. Gradient waveform design for tensor-valued encoding in diffusion MRI. *J Neurosci Methods*. 2021;348:109007.
- Jelescu IO, Veraart J, Fieremans E, Novikov DS. Degeneracy in model parameter estimation for multi-compartmental diffusion in neuronal tissue. *NMR Biomed*. 2016;29(1):33–47.
- Novikov DS, Fieremans E, Jespersen SN, Kiselev VG. Quantifying brain microstructure with diffusion MRI: Theory and parameter estimation. *NMR Biomed*. 2019;32(4):e3998.
- McNab JA, Edlow BL, Witzel T, Huang SY, Bhat H, Heberlein K, et al. The Human Connectome Project and beyond: initial applications of 300 mT/m gradients. *Neuroimage*. 2013;80:234–45.
- Fan Q, Eichner C, Afzali M, Mueller L, Tax CMW, Davids M, et al. Mapping the human connectome using diffusion MRI at 300 mT/m gradient strength: Methodological advances and scientific impact. *Neuroimage*. 2022;254:118958.
- Genc S, Raven EP, Drakesmith M, Blakemore SJ, Jones DK. Novel insights into axon diameter and myelin content in late childhood and adolescence. *Cereb Cortex*. 2023;33(10):6435–6448.
- Kleban E, Tax CMW, Rudrapatna US, Jones DK, Bowtell R. Strong diffusion gradients allow the separation of intra- and extra-axonal gradient-echo signals in the human brain. *Neuroimage*. 2020;217:116793.
- Tax CMW, Genc S, MacIver CL, Nilsson M, Wardle M, Szczepankiewicz F, et al. Ultra-strong diffusion-weighted MRI reveals cerebellar grey matter abnormalities in movement disorders. *Neuroimage Clin*. 2023;38:103419.

- 40 Huang SY, Tian Q, Fan Q, Witzel T, Wichtmann B, McNab JA, et al. High-gradient diffusion MRI reveals distinct estimates of axon diameter index within different white matter tracts in the in vivo human brain. *Brain Struct Funct*. 2020;225(4):1277–1291.
- 41 Nilsson M, Lasič S, Drobnjak I, Topgaard D, Westin CF. Resolution limit of cylinder diameter estimation by diffusion MRI: The impact of gradient waveform and orientation dispersion. *NMR Biomed*. 2017;30(7):e3711.
- 42 Fan F, Beare R, Genc S, Shapiro JS, Takagi M, Hearps SJC, et al. White matter fiber morphology in persisting postconcussive symptoms and posttraumatic headache after pediatric concussion: a fixel-based analysis. *J Neurosurg Pediatr*. 2024;34(4):373–383.
- 43 Meyer CH, Hu BS, Nishimura DG, Macovski A. Fast spiral coronary artery imaging. *Magn Reson Med*. 1992;28(2):202–213.
- 44 Glover GH, Pauly JM. Projection reconstruction techniques for reduction of motion effects in MRI. *Magn Reson Med*. 1992;28(2):275–89.
- 45 Wilm BJ, Hennel F, Roessler MB, Weiger M, Pruessmann KP. Minimizing the echo time in diffusion imaging using spiral readouts and a head gradient system. *Magn Reson Med*. 2020;84(6):3117–3127.
- 46 Lee Y, Wilm BJ, Brunner DO, Gross S, Schmid T, Nagy Z, et al. On the signal-to-noise ratio benefit of spiral acquisition in diffusion MRI. *Magn Reson Med*. 2021;85(4):1924–1937.
- 47 Feizollah S, Tardif CL. High-resolution diffusion-weighted imaging at 7 Tesla: Single-shot readout trajectories and their impact on signal-to-noise ratio, spatial resolution and accuracy. *Neuroimage*. 2023;274:120159.
- 48 Maki JH, MacFall JR, Johnson GA. The use of gradient flow compensation to separate diffusion and microcirculatory flow in MRI. *Magn Reson Med*. 1991;17(1):95–107.
- 49 Wetscherek A, Stieltjes B, Laun FB. Flow-compensated intravoxel incoherent motion diffusion imaging. *Magn Reson Med*. 2015;74(2):410–419.
- 50 Sjölund J, Szczepankiewicz F, Nilsson M, Topgaard D, Westin CF, Knutsson H. Constrained optimization of gradient waveforms for generalized diffusion encoding. *J Magn Reson*. 2015;261:157–168.
- 51 Westin CF, Knutsson H, Pasternak O, Szczepankiewicz F, Özarslan E, van Westen D, et al. Q-space trajectory imaging for multidimensional diffusion MRI of the human brain. *Neuroimage*. 2016;135:345–62.
- 52 Jones DK, Cercignani M. Twenty-five pitfalls in the analysis of diffusion MRI data. *NMR Biomed*. 2010 Aug;23(7):803–20.
- 53 Boesch C, Gruetter R, Martin E. Temporal and spatial analysis of fields generated by eddy currents in superconducting magnets: optimization of corrections and quantitative characterization of magnet/gradient systems. *Magn Reson Med*. 1991;20(2):268–84.
- 54 Aldefeld B, Börnert P. Effects of gradient anisotropy in MRI. *Magn Reson Med*. 1998;39(4):606–14.
- 55 Vannesjo SJ, Haerberlin M, Kasper L, Pavan M, Wilm BJ, Barmet C, et al. Gradient system characterization by impulse response measurements with a dynamic field camera. *Magn Reson Med*. 2013;69(2):583–93.
- 56 Bernstein MA, Zhou XJ, Polzin JA, King KF, Ganin A, Pelc NJ, et al. Concomitant gradient terms in phase contrast MR: analysis and correction. *Magn Reson Med*. 1998;39(2):300–8.
- 57 Dietrich BE, Brunner DO, Wilm BJ, Barmet C, Gross S, Kasper L, et al. A field camera for MR sequence monitoring and system analysis. *Magn Reson Med*. 2016;75(4):1831–40.
- 58 Wilm BJ, Barmet C, Pavan M, Pruessmann KP. Higher order reconstruction for MRI in the presence of spatiotemporal field perturbations. *Magn Reson Med*. 2011;65(6):1690–701.
- 59 Wilm BJ, Nagy Z, Barmet C, Vannesjo SJ, Kasper L, Haerberlin M, et al. Diffusion MRI with concurrent magnetic field monitoring. *Magn Reson Med*. 2015;74(4):925–33.
- 60 Molendowska M, Palombo M, Foley KG, Narahari K, Fasano F, Jones DK, et al. Diffusion MRI in prostate cancer with ultra-strong whole-body gradients. *NMR Biomed*. 2024;37(12):e5229.
- 61 Zhu A, Tarasek M, Hua Y, Fiveland E, Maier SE, Mazaheri Y, et al. Human prostate MRI at ultrahigh-performance gradient: A feasibility study. *Magn Reson Med*. 2024;91(2):640–648.
- 62 Afzali M, Mueller L, Coveney S, Fasano F, Evans CJ, Engel M, et al. In vivo diffusion MRI of the human heart using a 300 mT/m gradient system. *Magn Reson Med*. 2024;92(3):1022–1034.
- 63 Molendowska M, Mueller L, Fasano F, Jones DK, Tax CMW, Engel M. Giving the prostate the boost it needs: Spiral diffusion MRI using a high-performance whole-body gradient system for high b-values at short echo times. *Magn Reson Med*. 2025;93(3):1256–1272.
- 64 De Santis S, Assaf Y, Jones D. The influence of T2 relaxation in measuring the restricted volume fraction in diffusion MRI. *Proc Intl Soc Mag Reson Med*. 2016;24:1998.
- 65 Hompland T, Fjeldbo CS, Lyng H. Tumor Hypoxia as a Barrier in Cancer Therapy: Why Levels Matter. *Cancers (Basel)*. 2021;13(3):499.
- 66 Hompland T, Hole KH, Ragnum HB, Aarnes EK, Vlatkovic L, Lie AK, et al. Combined MR Imaging of Oxygen Consumption and Supply Reveals Tumor Hypoxia and Aggressiveness in Prostate Cancer Patients. *Cancer Res*. 2018;78(16):4774–4785.
- 67 Le Bihan D. What can we see with IVIM MRI? *Neuroimage*. 2019;187:56–67.
- 68 Iima M. Perfusion-driven Intravoxel Incoherent Motion (IVIM) MRI in Oncology: Applications, Challenges, and Future Trends. *Magn Reson Med Sci*. 2021;20(2):125–138.
- 69 While PT. A comparative simulation study of bayesian fitting approaches to intravoxel incoherent motion modeling in diffusion-weighted MRI. *Magn Reson Med*. 2017;78(6):2373–2387.
- 70 Novikov DS, Kiselev VG, Jespersen SN. On modeling. *Magn Reson Med*. 2018;79(6):3172–3193.
- 71 Ahlgren A, Knutsson L, Wirestam R, Nilsson M, Ståhlberg F, Topgaard D, et al. Quantification of microcirculatory parameters by joint analysis of flow-compensated and non-flow-compensated intravoxel incoherent motion (IVIM) data. *NMR Biomed*. 2016;29(5):640–649.
- 72 Molendowska M, Rashid IA, Foley K, Mueller L, Olsson LE, Brynolfsson P, et al. State-of-the-art perfusion fraction mapping in prostate cancer: intravoxel incoherent motion (IVIM) with velocity-compensated encoding and spiral readout. Poster session presented at the European Society for Magnetic Resonance in Medicine and Biology (ESMRMB); 2024 October 2–5; Barcelona, Spain.
- 73 Molendowska M, et al. Perfusion Fraction Mapping with IVIM in Prostate: The Benefit of Velocity-Compensated Encoding. Paper presented at: ISMRM Workshop on 40 Years of Diffusion: Past, Present & Future Perspectives. 2025 February 16–20. Kyoto, Japan.
- 74 Petersson E, Sohlén M, Ljungberg M, Jalnefjord O. Improving accuracy and precision of intra-voxel incoherent motion parameter estimates with b-value optimization and flow compensation. Poster session presented at the European Society for Magnetic Resonance in Medicine and Biology (ESMRMB); 2024 October 2–5; Barcelona, Spain.
- 75 Molendowska M, et al. A stimulation-adapted gradient design for velocity-compensated IVIM demonstrated in prostate at 200 mT/m. Paper presented at: ISMRM Workshop on 40 Years of Diffusion: Past, Present & Future Perspectives. 2025 February 16–20. Kyoto, Japan.
- 76 Cory DG, Garroway AN, Miller JB. Applications of spin transport as a probe of local geometry. In: Abstracts of Papers of the American Chemical Society. 1990;199:105–POLYAMER CHEMICAL SOC 1155 16TH ST, NW, WASHINGTON, DC 20036.
- 77 Mori S, van Zijl PC. Diffusion weighting by the trace of the diffusion tensor within a single scan. *Magn Reson Med*. 1995;33(1):41–52.
- 78 Heid O, Weber J. Diffusion tensor trace pulse sequence. *Proceedings of 5th ISMRM*. 1997:224.
- 79 Topgaard D. Isotropic diffusion weighting in PGSE NMR: Numerical optimization of the q-MAS PGSE sequence. *Microporous Mesoporous Mater*. 2013;178:60–63.
- 80 Szczepankiewicz F, Lasič S, van Westen D, Sundgren PC, Englund E, Westin CF, et al. Quantification of microscopic diffusion anisotropy disentangles effects of orientation dispersion from microstructure: applications in healthy volunteers and in brain tumors. *Neuroimage*. 2015;104:241–52.
- 81 Lasič S, Szczepankiewicz F, Eriksson S, Nilsson M, Topgaard D. Microanisotropy imaging: quantification of microscopic diffusion anisotropy and orientational order parameter by diffusion MRI with magic-angle spinning of the q-vector. *Front Phys*. 2014;2:11.
- 82 Jespersen SN, Lundell H, Sønderby CK, Dyrby TB. Orientationally invariant metrics of apparent compartment eccentricity from double pulsed field gradient diffusion experiments. *NMR Biomed*. 2013;26(12):1647–62.
- 83 Topgaard D. Multidimensional diffusion MRI. *J Magn Reson*. 2017;275:98–113.
- 84 Szczepankiewicz F, van Westen D, Englund E, Westin CF, Ståhlberg F, Lätt J, et al. The link between diffusion MRI and tumor heterogeneity: Mapping cell eccentricity and density by diffusional variance decomposition (DIVIDE). *Neuroimage*. 2016;142:522–532.
- 85 Langbein BJ, Szczepankiewicz F, Westin CF, Bay C, Maier SE, Kibel AS, et al. A Pilot Study of Multidimensional Diffusion MRI for Assessment of Tissue Heterogeneity in Prostate Cancer. *Invest Radiol*. 2021;56(12):845–853.
- 86 Nilsson M, Eklund G, Szczepankiewicz F, Skorpil M, Bryskhe K, Westin CF, et al. Mapping prostatic microscopic anisotropy using linear and spherical b-tensor encoding: A preliminary study. *Magn Reson Med*. 2021;86(4):2025–2033.
- 87 Molendowska M. Diffusion MRI of the prostate at 300 mT/m [thesis]. Cardiff, Wales. Cardiff University; 2023.
- 88 Gore JC, Xu J, Colvin DC, Yankeelov TE, Parsons EC, Does MD. Characterization of tissue structure at varying length scales using temporal diffusion spectroscopy. *NMR Biomed*. 2010;23(7):745–56.
- 89 Novikov DS, Jensen JH, Helpert JA, Fieremans E. Revealing mesoscopic structural universality with diffusion. *Proc Natl Acad Sci U S A*. 2014;111(14):5088–5093.
- 90 Xu J, Does MD, Gore JC. Quantitative characterization of tissue microstructure with temporal diffusion spectroscopy. *J Magn Reson*. 2009;200(2):189–197.
- 91 Jiang X, Li H, Xie J, McKinley ET, Zhao P, Gore JC, et al. In vivo imaging of cancer cell size and cellularity using temporal diffusion spectroscopy. *Magn Reson Med*. 2017;78(1):156–164.
- 92 Wu D, Jiang K, Li H, Zhang Z, Ba R, Zhang Y, et al. Time-Dependent Diffusion MRI for Quantitative Microstructural Mapping of Prostate Cancer. *Radiology*. 2022;303(3):578–587.
- 93 Molendowska M, et al. A novel framework for restriction-weighted q-space trajectory imaging (resQ) demonstrated in prostate cancer. Paper presented at: ISMRM Workshop on 40 Years of Diffusion: Past, Present & Future Perspectives. 2025 February 16–20. Kyoto, Japan.
- 94 Szczepankiewicz F, Safi M, Ceberg C, Gottschalk M, Sereti E, Bjartell A, et al. Multidimensional diffusion MRI for monitoring radiotherapy response in human prostate cancer xenografts in mice: A longitudinal pilot study. *Proc Intl Soc Magn Reson Med*. 2023;Abstract #4157.
- 95 Jespersen SN, Olesen JL, İlanuş A, Shemesh N. Effects of non-gaussian diffusion on "isotropic diffusion" measurements: An ex-vivo microimaging and simulation study. *J Magn Reson*. 2019;300:84–94.
- 96 de Swiet TM, Mitra PP. Possible Systematic Errors in Single-Shot Measurements of the Trace of the Diffusion Tensor. *J Magn Reson B*. 1996;111(1):15–22.
- 97 Beaulieu C, Fenrich FR, Allen PS. Multicomponent water proton transverse relaxation and T2-discriminated water diffusion in myelinated and nonmyelinated nerve. *Magn Reson Imaging*. 1998;16(10):1201–1210.
- 98 Does MD, Gore JC. Compartmental study of diffusion and relaxation measured in vivo in normal and ischemic rat brain and trigeminal nerve. *Magn Reson Med*. 2000;43(6):837–844.
- 99 de Almeida Martins JP, Topgaard D. Multidimensional correlation of nuclear relaxation rates and diffusion tensors for model-free investigations of heterogeneous anisotropic porous materials. *Sci Rep*. 2018;8(1):2488.
- 100 Ning L, Gagoski B, Szczepankiewicz F, Westin CF, Rathi Y. Joint RELaxation-Diffusion Imaging Moments to Probe Neurite Microstructure. *IEEE Trans Med Imaging*. 2020;39(3):668–677.
- 101 Benjamini D, Iacono D, Komlos ME, Perl DP, Brody DL, Basser PJ. Diffuse axonal injury has a characteristic multidimensional MRI signature in the human brain. *Brain*. 2021;144(3):800–816.
- 102 Benjamini D, Priemer DS, Perl DP, Brody DL, Basser PJ. Mapping astrogliosis in the individual human brain using multidimensional MRI. *Brain*. 2023;146(3):1212–1226.
- 103 Chatterjee A, Bourne RM, Wang S, Devaraj A, Gallan AJ, Antic T, et al. Diagnosis of Prostate Cancer with Noninvasive Estimation of Prostate Tissue Composition by Using Hybrid Multidimensional MR Imaging: A Feasibility Study. *Radiology*. 2018;287(3):864–873.
- 104 Palombo M, Valindria V, Singh S, Chiou E, Giganti F, Pye H, et al. Joint estimation of relaxation and diffusion tissue parameters for prostate cancer with relaxation-VERDICT MRI. *Sci Rep*. 2023;13(1):2999.
- 105 Chatterjee A, Antic T, Gallan AJ, Paner GP, Lin LI, Karczmar GS, et al. Histological validation of prostate tissue composition measurement using hybrid multi-dimensional MRI: agreement with pathologists' measures. *Abdom Radiol (NY)*. 2022;47(2):801–813.
- 106 Lemberskiy G, Fieremans E, Veraart J, Deng FM, Rosenkrantz AB, Novikov DS. Characterization of prostate microstructure using water diffusion and NMR relaxation. *Front Phys*. 2018 Sep;6:91.

Contact

Malwina Molendowska, Ph.D.
Department of Medical Radiation Physics
Lund University
Barrgatan 4, Skånes Universitetssjukhus
22185, Lund
Sweden
malwina.molendowska@med.lu.se

Fabrizio Fasano
Senior Collaboration Scientist
Siemens Healthineers
Camberley, Surrey, GU15 3YL
United Kingdom
fabrizio.fasano@siemens-healthineers.com

Derek K Jones, Professor
Cardiff University Brain Research Imaging
Centre (CUBRIC), Cardiff University
Maindy Rd, CF24 4HQ Cardiff, Wales
United Kingdom
jonesd27@cardiff.ac.uk

Filip Szczepankiewicz
Associate Professor
Department of Medical Radiation Physics
Lund University
Barrgatan 4, Skånes Universitetssjukhus
22185, Lund
Sweden
filip.szczepankiewicz@med.lu.se

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

Professor Xianlin Han
Department of General Surgery, Peking Union Medical College Hospital, Beijing, China

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
Con = conventional; HR = high resolution; DL = deep learning

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Bowen Wu, Radiographer Yitong Lu, Radiographer Professor Liang Zhu, M.D.

Contact

Liang Zhu, M.D.
Associate Professor, Department of Radiology
Peking Union Medical College Hospital
Shuaifuyuan No. 1, Dongcheng District
Beijing, 100730
China
zhuliang_pumc@163.com

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
	Con T2WI cor	HR T2WI cor	Con T2WI ax	HR T2WI ax	Conv DWI	HR DWI
Field of view (mm ³)	360 × 360	400 × 400	360 × 360	380 × 309	380 × 309	300 × 208
Matrix	320 × 320	448 × 358	320 × 320	384 × 269	128 × 128	150 × 150
Spatial resolution (mm ³)	1.1 × 1.1 × 4.0	0.45 × 0.45 × 4.0	1.1 × 1.1 × 6.0	0.5 × 0.5 × 3.0	1.5 × 1.5 × 4.0	1.0 × 1.0 × 3.0
Slices	32	32	25	48	38	48
Echo time (ms)	92	87	85	89	56	50
Repetition time (ms)	2000	567	2390	666	2000	2400
Averages	1	1	1	1	1, 4	1, 8
Diffusion gradients	–	–	–	–	3	3
B values (s/mm ²)	–	–	–	–	50, 800	50, 800
Bandwidth (Hz/pixel)	710	657	710	420	2298	1960
Acquisition time (min:s)	1:30	0:18	3:00	0:42	2:04	3:20
Fat suppression	–	–	SPAIR	SPAIR	SPAIR	Fat saturation
DL reconstruction	–	On	–	On	–	On
Sequence	HASTE	HASTE	BLADE	HASTE	SE-EPI	ZOOMit SE-EPI
Scanner	MAGNETOM Skyra	MAGNETOM Cima.X	MAGNETOM Skyra	MAGNETOM Cima.X	MAGNETOM Skyra	MAGNETOM Cima.X
Interpolation	Off	On	Off	On	On	On
Resp. control	Trigger	Breath-hold	Trigger	Breath-hold	Trigger	Trigger
Acceleration factor	3	4	2	4	2	2
Acceleration mode	GRAPPA	GRAPPA	GRAPPA	GRAPPA	SMS, GRAPPA	GRAPPA

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
- Mandell DM, Mossa-Basha M, Qiao Y, Hess CP, Hui F, Matouk C, et al. Intracranial Vessel Wall MRI: Principles and Expert Consensus Recommendations of the American Society of Neuroradiology. *AJNR Am J Neuroradiol.* 2017;38(2):218–229.
- Li L, Miller KL, Jezzard P. DANTE-prepared pulse trains: a novel approach to motion-sensitized and motion-suppressed quantitative magnetic resonance imaging. *Magn Reson Med.* 2012;68(5):1423–38.
- Li L, Chai JT, Biasioli L, Robson MD, Choudhury RP, Handa AI, et al. Black-blood multicontrast imaging of carotid arteries with DANTE-prepared 2D and 3D MR imaging. *Radiology.* 2014;273(2):560–9.
- Fushimi Y, Yoshida K, Okawa M, Maki T, Nakajima S, Sakata A, et al. Vessel wall MR imaging in neuroradiology. *Radiol Med.* 2022;127(9):1032–1045.
- Okuchi S, Fushimi Y, Okada T, Yamamoto A, Okada T, Kikuchi T, et al. Visualization of carotid vessel wall and atherosclerotic plaque: T1-SPACE vs. compressed sensing T1-SPACE. *Eur Radiol.* 2019;29(8):4114–4122.
- Takayama N, Maki T, Fushimi Y, Okawa M, Mineharu Y, Yamada K, et al. Association between ischemic stroke and hyperintense plaques detected by high-resolution vessel wall MRI in Japanese patients with intracranial atherosclerotic disease. *Front Stroke.* 2025; 4: 2025.
- Oichi Y, Mineharu Y, Agawa Y, Morimoto T, Funaki T, Fushimi Y, et al. Plasticity of the bony carotid canal and its clinical use for assessing negative remodeling of the internal carotid artery. *PLoS One.* 2021;16(12):e0261235.
- Matsukawa S, Ishii A, Fushimi Y, Grinstead J, Ahn S, Kikuchi T, et al. Efficacy of high-resolution vessel wall MRI in the postoperative assessment of intracranial aneurysms following flow diversion treatment. *J Neurosurg.* 2024;142(1):88–97.
- Matsukawa S, Ishii A, Fushimi Y, Terada Y, Natsuhara H, Kikuchi T, et al. Ruptured thrombosed vertebral artery dissecting aneurysm treated with staged flow diverter after prediction of the rupture point by vessel wall MRI. *Neuroradiol J.* 2024;37(6):772–775.
- Tokunaga S, Yamao Y, Maki T, Ishii A, Miyake T, Yasuda K, et al. Vessel wall imaging and carotid artery stenting for recurrent cervical internal carotid artery vasospasm syndrome: illustrative case. *J Neurosurg Case Lessons.* 2023;5(23):CASE2341.
- Tagawa H, Fushimi Y, Funaki T, Nakajima S, Sakata A, Okuchi S, et al. Vessel wall MRI in moyamoya disease: arterial wall enhancement varies depending on age, arteries, and disease progression. *Eur Radiol.* 2024;34(4):2183–2194.
- Komatsu K, Takagi Y, Ishii A, Kikuchi T, Yamao Y, Fushimi Y, et al. Ruptured intranidal aneurysm of an arteriovenous malformation diagnosed by delay alternating with nutation for tailored excitation (DANTE)-prepared contrast-enhanced magnetic resonance imaging. *Acta Neurochir (Wien).* 2018;160(12):2435–2438.
- Shimada A, Suda K, Nakano E, Tagawa M, Miyata M, Kashii S, et al. Accuracy of Diagnosing Optic Neuritis Using DANTE T1-SPACE Imaging. *Eye Brain.* 2024;16:65–73.
- Tsubouchi R, Ohira J, Sawamura M, Fushimi Y, Grinstead J, Ahn S, et al. Multiple cranial neuritis depicted with DANTE-prepared contrast-enhanced MRI. *Neurology and Clinical Neuroscience.* 2020;8:220–1.
- Oshima S, Fushimi Y, Okada T, Nakajima S, Yokota Y, Shima A, et al. Neuromelanin-Sensitive Magnetic Resonance Imaging Using DANTE Pulse. *Mov Disord.* 2021;36(4):874–882.
- Park S, Sung YH, Kim WR, Noh Y, Kim EY. Correlation Between Neuromelanin-Sensitive MRI and ¹⁸F-FP-CIT PET in Early-Stage Parkinson's Disease: Utility of a Voxel-Wise Analysis by Using High-Spatial-Resolution MRI. *J Clin Neurol.* 2023;19(2):156–164.
- Oshima S, Fushimi Y, Miyake KK, Nakajima S, Sakata A, Okuchi S, et al. Denoising approach with deep learning-based reconstruction for neuromelanin-sensitive MRI: image quality and diagnostic performance. *Jpn J Radiol.* 2023;41(11):1216–1225.

Contact

Associate Professor Yasutaka Fushimi, M.D., Ph.D.
Department of Diagnostic Imaging
and Nuclear Medicine
Graduate School of Medicine, Kyoto University
54 Shōgoin Kawaharachō, Sakyo-ku
Kyoto 606-8507
Japan
yfushimi@kuhp.kyoto-u.ac.jp

From left to right:
Dr. Sakata, Dr. Fushimi, Dr. Okuchi

Front row left to right: Dr. Fushimi and Dr. Okuchi
Back row left to right: Dr. Yamamoto, Dr. Nakajima,
and colleagues

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Contact

Professor Mengsu Zeng, M.D.
Director, Radiology Department
Zhongshan Hospital
Fudan University
No. 180, Fenglin Road, Xuhui
Shanghai
China
zeng.mengsu@zs-hospital.sh.cn

Professor Shengxiang Rao, M.D.
Vice Director, Radiology Department
Zhongshan Hospital
Fudan University
No. 180, Fenglin Road, Xuhui
Shanghai
China
rao.shengxiang@zs-hospital.sh.cn

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References

- Chandra SS, Bran Lorenzana M, Liu X, Liu S, Bollmann S, Crozier S. Deep learning in magnetic resonance image reconstruction. J Med Imaging Radiat Oncol. 2021 Aug;65(5):564–577.
- Herrmann J, Lingg A, Kündel M, Gassenmaier S, Afat S, Al-Mansour H, et al. Clinical Implementation of Deep Learning-Accelerated HASTE and TSE. MAGNETOM Flash. 2021;79(2):23–27.
- Herrmann J, Nickel D, Mugler JP 3rd, Arberet S, Gassenmaier S, Afat S, et al. Development and Evaluation of Deep Learning-Accelerated Single-Breath-Hold Abdominal HASTE at 3 T Using Variable Refocusing Flip Angles. Invest Radiol. 2021;56(10):645–652.
- Herrmann J, Gassenmaier S, Nickel D, Arberet S, Afat S, Lingg A, et al. Diagnostic Confidence and Feasibility of a Deep Learning Accelerated HASTE Sequence of the Abdomen in a Single Breath-Hold. Invest Radiol. 2021;56(5):313–319.
- Liu K, Li Q, Wang X, Fu C, Sun H, Chen C, et al. Feasibility of deep learning-reconstructed thin-slice single-breath-hold HASTE for detecting pancreatic lesions: A comparison with two conventional T2-weighted imaging sequences. Res Diagn Interv Imaging. 2024;9:100038.
- Park MJ, Kim YK, Lee MW, Lee WJ, Kim YS, Kim SH, et al. Small hepatocellular carcinomas: improved sensitivity by combining gadoxetic acid-enhanced and diffusion-weighted MR imaging patterns. Radiology. 2012;264(3):761–70.
- Chung WS, Kim MJ, Chung YE, Kim YE, Park MS, Choi JY, et al. Comparison of gadoxetic acid-enhanced dynamic imaging and diffusion-weighted imaging for the preoperative evaluation of colorectal liver metastases. J Magn Reson Imaging. 2011;34(2):345–53.

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Sequence	CE T1w FS	FLAIR FS
Echo time (ms)	10	91
Repetition time (ms)	755	6070
Slice thickness (mm)	2.5	2.5
Acquisition matrix	448 × 314	204 × 274
Imaged field of view (mm)	149 × 149	180 × 180
Readout bandwidth (Hz/pixel)	180	289
Acquisition time (min:sec)	4:50	1:38
iPAT factor	2	4
Inversion time (ms)	N/A	1855

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(Image quality and diagnostic confidence)

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



The study is available as an open access article here:

<https://doi.org/10.1097/rli.0000000000001065>


References

- Filippucci E, Di Geso L, Grassi W. Progress in imaging in rheumatology. Nat Rev Rheumatol. 2014;10(10):628-34.
- Roemer FW, Kassim Javid M, Guermazi A, Thomas M, Kiran A, Keen R, et al. Anatomical distribution of synovitis in knee osteoarthritis and its association with joint effusion assessed on non-enhanced and contrast-enhanced MRI. Osteoarthritis Cartilage. 2010;18(10):1269-74.
- Ramalho J, Semelka RC, Ramalho M, Nunes RH, AIObaidy M, Castillo M. Gadolinium-Based Contrast Agent Accumulation and Toxicity: An Update. AJNR Am J Neuroradiol. 2016;37(7):1192-8. Orthop J Sports Med. 2020;8(3):2325967120909090.
- McDonald RJ, McDonald JS, Kallmes DF, Jentoft ME, Murray DL, Thielen KR, et al. Intracranial Gadolinium Deposition after Contrast-enhanced MR Imaging. Radiology. 2015;275(3):772-82.
- Yoo HJ, Hong SH, Oh HY, Choi JY, Chae HD, Ahn JM, et al. Diagnostic Accuracy of a Fluid-attenuated Inversion-Recovery Sequence with Fat Suppression for Assessment of Peripatellar Synovitis: Preliminary Results and Comparison with Contrast-enhanced MR Imaging. Radiology. 2017;283(3):769-778.
- Johnson PM, Lin DJ, Zbontar J, Zitnick CL, Sriram A, Muckley M, et al. Deep Learning Reconstruction Enables Prospectively Accelerated Clinical Knee MRI. Radiology. 2023;307(2):e220425.
- Feuerriegel GC, Goller SS, von Deuster C, Sutter R. Inflammatory Knee Synovitis: Evaluation of an Accelerated FLAIR Sequence Compared With Standard Contrast-Enhanced Imaging. Invest Radiol. 2024;59(8):599-604.
- Hammernik K, Klatzer T, Kobler E, Recht MP, Sodickson DK, Pock T, et al. Learning a variational network for reconstruction of accelerated MRI data. Magn Reson Med. 2018;79(6):3055-3071.

Contact

Georg C. Feuerriegel, M.D.
Fellow
Department of Radiology
Balgrist University Hospital
Forchstrasse 340
8008 Zurich
Switzerland
georg.feuerriegel@balgrist.ch

Reto Sutter, M.D.
Professor of Radiology
Head of Department
Department of Radiology
Balgrist University Hospital
Forchstrasse 340
8008 Zurich
Switzerland
reto.sutter@balgrist.ch

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References


- Ophir J, Céspedes I, Ponnekanti H, Yazdi Y, Li X. Elastography: a quantitative method for imaging the elasticity of biological tissues. *Ultrason Imaging*. 1991;13(2):111–134.
- Lewa GJ. Elastic properties imaging by periodical displacement NMR measurements (EPMRI). *Proceedings of the Ultrasonics Symposium IEEE*, 1994;2:691–694.
- Muthupillai R, Lomas DJ, Rossman PJ, Greenleaf JF, Manduca A, Ehman RL. Magnetic resonance elastography by direct visualization of propagating acoustic strain waves. *Science*. 1995;269(5232):1854–1857.
- Sinkus R, Lorenzen J, Schrader D, Lorenzen M, Dargatz M, Holz D. High-resolution tensor MR elastography for breast tumour detection. *Phys Med Biol*. 2000;45(6):1649–1664.
- Allen AM, Shah VH, Therneau TM, Venkatesh SK, Mounajjed T, Larson JJ, et al. The Role of Three-Dimensional Magnetic Resonance Elastography in the Diagnosis of Nonalcoholic Steatohepatitis in Obese Patients Undergoing Bariatric Surgery. *Hepatology*. 2020;71(2):510–521.
- Sinkus R, Lambert S, Abd-Elmoniem KZ, Morse C, Heller T, Guenther C, et al. Rheological determinants for simultaneous staging of hepatic fibrosis and inflammation in patients with chronic liver disease. *NMR Biomed*. 2018;31(10):e3956.
- McKnight AL, Kugel JL, Rossman PJ, Manduca A, Hartmann LC, Ehman RL. MR elastography of breast cancer: preliminary results. *AJR Am J Roentgenol*. 2002;178(6):1411–1417.
- Bunevicius A, Schregel K, Sinkus R, Golby A, Patz S. REVIEW: MR elastography of brain tumors. *Neuroimage Clin*. 2020;25:102109.
- Flogstad Svensson S, Fuster-Garcia E, Latysheva A, Fraser-Green J, Nordhøy W, Isam Darwish O, et al. Decreased tissue stiffness in glioblastoma by MR elastography is associated with increased cerebral blood flow. *Eur J Radiol*. 2022;147:110136.
- Glaser KJ, Manduca A, Ehman RL. Review of MR elastography applications and recent developments. *J Magn Reson Imaging*. 2012;36(4):757–774.
- Rump J, Klatt D, Braun J, Warmuth C, Sack I. Fractional encoding of harmonic motions in MR elastography. *Magn Reson Med*. 2007 Feb;57(2):388–395.
- Guenther C, Sethi S, Troelstra M, Dokumaci AS, Sinkus R, Kozerke S. Ristretto MRE: A generalized multi-shot GRE-MRE sequence. *NMR Biomed*. 2019 May;32(5):e4049.
- Muthupillai R, Rossman PJ, Lomas DJ, Greenleaf JF, Riederer SJ, Ehman RL. Magnetic resonance imaging of transverse acoustic strain waves. *Magn Reson Med*. 1996 Aug;36(2):266–274.
- Bercoff J, Tanter M, Fink M. Supersonic shear imaging: a new technique for soft tissue elasticity mapping. *IEEE Trans Ultrason Ferroelectr Freq Control*. 2004;51(4):396–409.
- Kruse SA, Smith JA, Lawrence AJ, Dresner MA, Manduca A, Greenleaf JF, et al. Tissue characterization using magnetic resonance elastography: preliminary results. *Phys Med Biol*. 2000;45(6):1579–1590.
- Van Houten EE, Paulsen KD, Miga MI, Kennedy FE, Weaver JB. An overlapping subzone technique for MR-based elastic property reconstruction. *Magn Reson Med*. 1999 Oct;42(4):779–786.
- Sinkus R, Siegmann K, Xydeas T, Tanter M, Claussen C, Fink M. MR elastography of breast lesions: understanding the solid/liquid duality can improve the specificity of contrast-enhanced MR mammography. *Magn Reson Med*. 2007;58(6):1135–1144.
- Manduca A, Lake DS, Kruse SA, Ehman RL. Spatio-temporal directional filtering for improved inversion of MR elastography images. *Med Image Anal*. 2003;7(4):465–473.
- Uffmann K, Ladd M. Actuation systems for MR elastography: design and applications. *IEEE Eng Med Biol Mag*. 2008;27(3):28–34.
- Van Houten EE, Doyley MM, Kennedy FE, Weaver JB, Paulsen KD. Initial in vivo experience with steady-state subzone-based MR elastography of the human breast. *J Magn Reson Imaging*. 2003;17(1):72–85.
- Yin M, Talwalkar JA, Glaser KJ, Manduca A, Grimm RC, Rossman PJ, et al. Assessment of hepatic fibrosis with magnetic resonance elastography. *Clin Gastroenterol Hepatol*. 2007;5(10):1207–1213.e2.
- Dittmann F, Tzschätzsch H, Hirsch S, Barnhill E, Braun J, Sack I, et al. Tomoelastography of the abdomen: Tissue mechanical properties of the liver, spleen, kidney, and pancreas from single MR elastography scans at different hydration states. *Magn Reson Med*. 2017;78(3):976–983.
- Runge JH, Hoelzl SH, Sudakova J, Dokumaci AS, Nelissen JL, Guenther C, et al. A novel magnetic resonance elastography transducer concept based on a rotational eccentric mass: preliminary experiences with the gravitational transducer. *Phys Med Biol*. 2019;64(4):045007.
- Fovargue D, Kozerke S, Sinkus R, Nordsletten D. Robust MR elastography stiffness quantification using a localized divergence free finite element reconstruction. *Med Image Anal*. 2018;44:126–142.
- Wochner P, Schneider T, Stockmann J, Lee J, Sinkus R. Diffusion phase-imaging in anisotropic media using non-linear gradients for diffusion encoding. *PLoS One*. 2023;18(3):e0281332.
- Venkatesh SK, Yin M, Ehman RL. Magnetic resonance elastography of liver: technique, analysis, and clinical applications.

- J Magn Reson Imaging. 2013;37(3):544–555.
- 27 Guenthner C, Runge JH, Sinkus R, Kozerke S. Hadamard encoding for magnetic resonance elastography. In: Proceedings of the 25th Annual Meeting ISMRM. 2017. Honolulu.
- 28 Garteiser P, Sahebjavaher RS, Ter Beek LC, Salcudean S, Vilgrain V, Van Beers BE, et al. Rapid acquisition of multifrequency, multislice and multidirectional MR elastography data with a fractionally encoded gradient echo sequence. NMR Biomed. 2013;26(10):1326–1335.
- 29 Barth M, Breuer F, Koopmans PJ, Norris DG, Poser BA. Simultaneous multislice (SMS) imaging techniques. Magn Reson Med. 2016;75(1):63–81.
- 30 Darwish OI, Gharib AM, Jeljeli S, Metwalli NS, Feeley J, Rotman Y, et al. Single Breath-Hold 3-Dimensional Magnetic Resonance Elastography Depicts Liver Fibrosis and Inflammation in Obese Patients. Invest Radiol. 2023;58(6):413–419.
- 31 Johnson CL, McGarry MD, Van Houten EE, Weaver JB, Paulsen KD, Sutton BP, et al. Magnetic resonance elastography of the brain using multishot spiral readouts with self-navigated motion correction. Magn Reson Med. 2013;70(2):404–412.
- 32 Einstein A, Podolsky B, Rosen N. Can Quantum-Mechanical Description of Physical Reality Be Considered Complete? Phys Rev. 1935;47(10):777–780.
- 33 Landau LD, Lifshitz EM. Course of Theoretical Physics: Theory of Elasticity. 2nd English ed. Oxford: Pergamon Press; 1970. p. 165.
- 34 Bhatia H, Norgard G, Pascucci V, Bremer P T. The Helmholtz-Hodge Decomposition – A Survey. IEEE Transactions on Visualization and Computer Graphics. 2013;19(8):1386–1404.
- 35 Van Houten EE, Miga MI, Weaver JB, Kennedy FE, Paulsen KD. Three-dimensional subzone-based reconstruction algorithm for MR elastography. Magn Reson Med. 2001;45(5):827–837.
- 36 Verdier C. Rheological Properties of Living Materials. From Cells to Tissues. Journal of Theoretical Medicine. 2003;5(2):67–91.
- 37 Schiessel H, Blumen A. Mesoscopic Pictures of the Sol-Gel Transition: Ladder Models and Fractal Networks. Macromolecules. 1995;28(11):4013–4019.
- 38 Holm S. Waves with Power-Law Attenuation. New York: Springer, 2019.
- 39 Lambert SA, Näsholm SP, Nordsletten D, Michler C, Juge L, Serfaty JM, et al. Bridging Three Orders of Magnitude: Multiple Scattered Waves Sense Fractal Microscopic Structures via Dispersion. Phys Rev Lett. 2015;115(9):094301.
- 40 Gidener T, Ahmed OT, Larson JJ, Mara KC, Therneau TM, Venkatesh SK, et al. Liver Stiffness by Magnetic Resonance Elastography Predicts Future Cirrhosis, Decompensation, and Death in NAFLD. Clin Gastroenterol Hepatol. 2021;19(9):1915–1924.e6.
- 41 Sheka AC, Adeyi O, Thompson J, Hameed B, Crawford PA, Ikramuddin S. Nonalcoholic Steatohepatitis: A Review. JAMA. 2020;323(12):1175–1183.
- 42 Rada P, González-Rodríguez Á, García-Monzón C, Valverde ÁM. Understanding lipotoxicity in NAFLD pathogenesis: is CD36 a key driver? Cell Death Dis. 2020;11(9):802.
- 43 Loomba R, Adams LA. The 20% Rule of NASH Progression: The Natural History of Advanced Fibrosis and Cirrhosis Caused by NASH. Hepatology. 2019;70(6):1885–1888.
- 44 Pepin KM, McGee KP, Arani A, Lake DS, Glaser KJ, Manduca A, et al. MR Elastography Analysis of Glioma Stiffness and IDH1-Mutation Status. AJNR Am J Neuroradiol. 2018;39(1):31–36.
- 45 Svensson SF, Halldórsson S, Latysheva A, Fuster-Garcia E, Hjørnevik T, Fraser-Green J, et al. MR elastography identifies regions of extracellular matrix reorganization associated with shorter survival in glioblastoma patients. Neurooncol Adv. 2023;5(1):vdad021.
- 46 Hiscox LV, Schwarb H, McGarry MDJ, Johnson CL. Aging brain mechanics: Progress and promise of magnetic resonance elastography. Neuroimage. 2021;232:117889.
- 47 Wuerfel J, Paul F, Beierbach B, Hamhaber U, Klatt D, Papazoglou S, et al. MR-elastography reveals degradation of tissue integrity in multiple sclerosis. Neuroimage. 2010;49(3):2520–2525.
- 48 Bohte AE, Nelissen JL, Runge JH, Holub O, Lambert SA, de Graaf L, et al. Breast magnetic resonance elastography: a review of clinical work and future perspectives. NMR Biomed. 2018;31(10):e3932.
- 49 Jurrius P, Darwish O, Shifa B, Bell J, Spence J, Fovargue D, et al. Measuring the tumour IFP using MRE: a mechanical biomarker for the prediction of metastatic potential in women with invasive breast cancer. Abstract #1436 in: Proceedings of the 30th Annual Meeting ISMRM. 2021.
- 50 Jurrius P, Darwish O, Fovargue D, Shifa B, Welsh K, Bell J, et al. Magnetic Resonance Elastography for Early Assessment of Response to Neoadjuvant Chemotherapy in Women with Breast Cancer. Abstract #1591 in: Proceedings of the 31st Annual Meeting ISMRM. 2022.

Contact

Ralph Sinkus, Ph.D.
Professor of Biomedical Engineering
at King's College London, UK,
and Research Director CNRS
at INSERM, Paris, France
King's College London
School of Biomedical Engineering &
Imaging Sciences
Rayne Institute, 4th Floor, Lambeth Wing
St. Thomas' Hospital
Westminster Bridge Road
London SE1 7EH, United Kingdom
ralph.sinkus@kcl.ac.uk

Omar Darwish
Siemens Healthineers
SHS DI MR RCT CLS BODY
Allee am Roethelheimpark 2
91052 Erlangen
Germany
omar.darwish@siemens-healthineers.com