


Clinical Value of Cios Spin and Intuitive Ion

Interpretation of outcomes using robotic-assisted bronchoscopy with integrated cone beam CT

[siemens-healthineers.com](https://www.siemens-healthineers.com)



Results from a high volume center

A recent prospective cohort study from the University Hospital Zürich (Universitätsspital Zürich) evaluated the impact of robotic-assisted bronchoscopy (RAB) using the Ion endoluminal system combined with integrated Cios Spin cone beam CT (CBCT) on the diagnostic yield and workflow in patients with peripheral pulmonary nodules (PPNs).

Key clinical outcomes¹

89% diagnostic yield


(ATS 2024 strict criteria):
median lesion size 10 mm,
15% with positive bronchus sign

1.25 spins/lesion

CBCT-integration into RAB
allows for a low number
of spins

37 min. procedure time²

CBCT+RAB procedure times are
similar to conventional bronchoscopy
with 2D fluoroscopy (34 Min)



The team's learning curve extended over a period of up to six months, during which a procedure time reduction (scope-in to scope-out) of 40–60 % can be expected. Regression analysis demonstrated significant improvements over this period, including a 26 % reduction in radiation dose (Gy), a 28 % decrease in fluoroscopy time, and 12 % fewer cone-beam CT spins (all $p < 0.01$), each correlating with the procedural learning curve. In contrast, key clinical parameters such as lesion size, complication rates, and diagnostic yield remained stable throughout the learning phase.

¹ First 6 months of program with 141 patients and 213 PPLs

² 44% reduction in procedure time, from initial 67 min to 37 min, defined by anesthesia (scope in to scope out).

The information shown herein refers to a product of 3rd party manufacturer Intuitive and thus is in their regulatory responsibility. Please contact the 3rd party manufacturer for further information and availability.



“With the integration of cone-beam CT into robotic-assisted bronchoscopy, our team adapted quickly – procedure times now average 37 minutes, only slightly longer than the 34 minutes it takes us to perform conventional bronchoscopy under 2D fluoroscopy. Also, radiation exposure was lowered while increasing patient throughput, improving safety and access.”

Thomas Gaisl, MD, PhD

Managing consultant, Department of Pulmonology University Hospital Zurich

The integration of robotic-assisted bronchoscopy (RAB) with cone beam CT (Cios Spin) from Siemens Healthineers at University Hospital Zürich transformed lung cancer detection.



Dr. Gaisl’s research shows a remarkable surge in early-stage (1A) detections. This significant stage shift enables earlier intervention when treatments can be more cost-effective.³

³ Earlier detection also has substantial economic implications. According to current data, monthly excess expenditures for lung cancer treatment increase significantly with disease progression, ranging from approximately €3,115 for Stage IA cancer to €8,156 for Stage IVB cancer [de Nijs K, de Koning HJ, van der Aalst C, Ten Haaf K. Medical costs of lung cancer by stage, histology and first-line treatment modality in the Netherlands (2012-2021). *Eur J Cancer*. 2024;208:114231.].

Publication: *American Journal of Respiratory and Critical Care Medicine*, Vol 211, Abstract A5512, DOI: 10.1164/ajrccm.2025.211.Abstracts.A5512 Author: Dr. Thomas Gaisl, MD, MPH, PhD (thomas.gaisl@usz.ch) Disclaimer: This flyer summarizes data from one institution’s clinical experience and interpretation of outcomes. Results may vary. Device performance and workflow outcomes depend on multiple factors, including operator experience and institutional protocols. This communication is intended for informational use only and does not constitute a claim of clinical performance. Refer to official labeling for intended use, contraindications, and safety information.

The outcomes and statements provided by customers of Siemens Healthineers are unique to each customer’s setting. Since there is no “typical” hospital and many variables exist (e.g., hospital size, case mix, and level of service/technology adoption), there can be no guarantee that others will achieve the same results.

Dr. Gaisl is employed by an institution that receives financial support from Siemens Healthineers for collaborations.

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
siemens-healthineers.com