Vitamin D Total Assay

The only CDC certified method with defined pediatric expected values











Ranges established for pediatrics

- Easily adopt Atellica® IM and ADVIA Centaur® Vitamin D Total assay pediatric reference ranges*
- Optimize clinical operations by reducing lab costs and resources associated with establishing ranges.*
- Reduce the potential for over-recovery and risk of missed deficiency, with only 1% cross-reactivity to 3-epi-25(OH)vitamin D3
- Prevent underestimation of results with equimolar detection of vitamin D2 (104.5%) and D3 (100.7%)

Based on the 95% confidence interval, the following pediatric values were obtained:^{5,6}

	Specimen type	Observed values, pediatric (12 months up to 21 years)
Median 25(OH) vitamin D	Serum	23.8 ng/mL (59.5 nmol/L)
Observed range 2.5th to 97.5th percentile	Serum	11.4–45.8 ng/mL (28.5–114.5 nmol/L)

^{*}Pediatric reference ranges are available for Atellica IM and ADVIA Centaur. As with all in vitro diagnostic assays, each laboratory should determine its own reference interval for the diagnostic evaluation of patient results. Consider these values as guidance only. Refer to CLSI EP28-A3c for recommendations on establishing and verifying ranges.

Atellica IM and ADVIA Centaur Vitamin D Total reference values were established in accordance with CLSI guideline EP28-A3c using adult and pediatric serum samples.

The pediatric serum samples used to establish the pediatric reference values included:5,6

- 237 apparently healthy male and female subjects of light and dark skin types:
 - 32 between the ages of 1-3 years
 - 114 from 3-12 years
 - 91 from 13-21 years
- Samples from different seasons
- Samples from different geographical regions of the United States
- Samples with normal values for PTH and TSH

Expected pediatric values based on literature reviewed by Siemens Healthineers:

Vitamin D Status	Range, Pediatric
Deficiency	< 15 ng/mL (37.5 nmol/L)
Insufficiency	15- < 20 ng/mL (37.5- < 50 nmol/L)
Sufficiency	20-100 ng/mL (50-250 nmol/L)


siemens-healthineers.us/Vitamin-D



CDC Certification (VDSCP*)1,2,3

- Ensures accuracy and reliability of vitamin D testing using well established procedures
- Certifies those methods that have a certain accuracy and precision:
 - Maximum allowable bias of ≤5%
 - Maximum allowable imprecision of <10%
- Program participation requires quarterly performance assessments for re-certification

According to the CDC, there is an improvement in *Mean Bias among VDSCP participants:1



^{*}The difference between the average of measurements of vitamin D made by VDSCP participants and the CDC-assigned reference value.

The CDC has certified the following assays from Siemens Healthineers:3,4

- The ADVIA Centaur® Vitamin D Total assay
- The Atellica® IM Vitamin D Total assay
- The Dimension® EXL™ LOCI Vitamin D Total assay

A complete list of certified methods/vendors are posted and continually updated on the CDC website.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH Henkestr. 127 91052 Erlangen, Germany Phone: +49 9131 84-0 siemens-healthineers.com

USA

Siemens Healthcare Diagnostics Inc. Laboratory Diagnostics 511 Benedict Avenue Tarrytown, NY 10591-5005, USA Phone: +1-914-631-8000 siemens-healthineers.us On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. For the most current information, please contact your local sales representative from Siemens Healthineers.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

^{*}Vitamin D Standardization-Certification Program

¹https://www.cdc.gov/labstandards/vdscp.html

²https://www.cdc.gov/labstandards/vdscp_procedures.html

³https://www.cdc.gov/labstandards/vdscp_participants.html

⁴https://www.cdc.gov/labstandards/pdf/hs/CDC_Certified_Vitamin_D_Procedures.pdf

⁵Atellica IM Analyzer Vitamin D Total (VitD)

⁶ADVIA Centaur, ADVIA Centaur XP, ADVIA XPT Immunoassay Systems Vitamin D Total (VitD)