



syngo[®] Dynamics

Online Help - HL7 Interface Specifications Manual
(Standalone)

VA41

This page is left blank intentionally.

1	Intended Use/Intended Purpose	5
	Intended Use/Indications for Use Statement	5
	Intended Purpose Statement	5
	Indications	5
	Contraindications	6
2	Responsibility Statement	7
3	Introduction	9
4	Legal Notes	11
	User Profile	11
	Training	13
5	Documentation Overview	15
	Conventions Used in this Documentation	15
	Labeling of Safety Advises	16
	Safety Advisory	16
	About this Manual	16
6	HL7 Overview	17
	Purpose	17
	Scope	22
	Abbreviations	22
7	HL7 Study Reconciliation of PHI Data	27
8	Pre-procedural Documentation	29
9	Report Upload to OPENLink Option	31
	HL7 Message Details Message Specifications	31
	syngo Dynamics Report Format	31
	Default Segment Listing	33
	Report Upload Message Examples	38
10	IHE DRPT Report Creator	41
	IHE DRPT Message Examples	42
11	Detailed Financial Transaction (DFT) Messages	43
	Default Segment Listing	46
	MSH Segment	46
	EVN Segment	47
	PID Segment	47
	PV1 Segment (Optional)*	48
	OBR Segment (Optional)	50
	FT1 Segment	50

	PR1 Segment (Optional)*	52
12	syngo Dynamics Reporting	55
	Allergy Information	55
	Read-only PHI Fields	55
	Order Messages (ORM) Mapping	57
13	Message Segments.	61
	Message Segment: MSH	61
	Message Segment: EVN	62
	Message Segment: PID	62
	Message Segment: PV1	65
	Message Segment: ORC	66
	Message Segment: MRG	69
	Message Segment: DG1	69
	Message Segment: AL1	69
	Message Segment: NK1	70
	Message Segment: ZDS (Customized Segment)	71
	Message Segment: ZSU	71
14	Appendix A: Sample Messages/Generic HL7 Information	73
	ORM-O01 Message	73
	ADT-A01 Message	73
	ADT-A40 Message (Merge Identifier)	74
	Message with Z Segment	74
	ZSU Segment in ORM Message	74
	ACK Messages	75
	ACK Scenarios Currently Handled	75
15	Appendix B: Configuration Files for ADT/ORM Option	79

1 Intended Use/Intended Purpose

1.1 Intended Use/Indications for Use Statement

syngo Dynamics is a multimodality, vendor agnostic Cardiology image and information system intended for medical image management and processing that provides capabilities relating to the review and digital processing of medical images.

syngo Dynamics supports clinicians by providing image post-processing functions for image manipulation, and/or quantification that are intended for use in the interpretation and analysis of medical images for disease detection, diagnosis, and/or patient management within the healthcare institution's network.

syngo Dynamics is not intended to be used for display or diagnosis of digital mammography images in the U.S.

1.2 Intended Purpose Statement

Intended Purpose Statement under the European Medical Device Regulation 2017/745:

Software solutions intended to process, communicate, display, read, and archive medical data for informing and driving clinical management.

1.3 Indications

syngo Dynamics supports clinicians by providing image post-processing functions for image manipulation, and/or quantification that are intended for use in the interpretation and analysis of medical images for disease detection, diagnosis, and/or patient management within the healthcare institution's network.

syngo Dynamics is not indicated for any specific disease, condition, patient population.

1.4 Contraindications

syngo Dynamics is not intended to be used for display or diagnosis of digital mammography images in the U.S.

There are no known clinical contraindications for *syngo* Dynamics.

2 Responsibility Statement

The evidence-based reporting in *syngo* Dynamics is intended to be a supporting tool for the diagnosing physician. Physicians should look at all relevant images and measurements prior to concluding on a diagnosis. While Siemens Healthineers may provide sample templates for user customization, the accuracy and appropriateness of clinical phrases, clinical observations, and user-defined measurements and calculations in *syngo* Dynamics worksheets and reports are the responsibility of the diagnosing physician and/or facility.

All quantitative data ranges are derived from the clinical experience of echocardiology laboratories and are included in observation libraries for *syngo* Dynamics users. Siemens Healthineers strongly recommends that clinicians review these ranges with their individual diagnostic needs in mind prior to using *syngo* Dynamics for clinical reporting.



Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

This page is left blank intentionally.

3 Introduction

syngo Dynamics is aimed at providing physicians and staff with fast, easy access to important patient data. *syngo* Dynamics has made the all-digital cardiology department a practical reality, incorporating images and information from echocardiography, cardiac catheterization, nuclear medicine, OB, and ultrasound departments. Ease of use, customizable reporting templates, and clinical tools all combine to provide our customers with improved workflow efficiencies and help provide optimal patient care. *syngo* Dynamics supports a streamlined, rapid and efficient decision-making process based on better and more accessible information.

For better readability, we refer to the user in the masculine form.

All parameters and images shown in this document are examples. Only the parameters displayed in your system are definite. All names and data of patients and institutions that are used in this document are entirely fictional. Any resemblance to names of existing people or organizations past or present is entirely coincidental.

Configuration-dependent designations are names of drives, network nodes, and databases. The names that are used in this document are usually not the same as the designations to be found on a particular installation of the system in a clinical environment. For confidentiality reasons, some of the data may have been masked in the images.

This page is left blank intentionally.

4 Legal Notes

This document is confidential, proprietary to Siemens Healthineers, protected by copyright laws in the U.S. and abroad, and licensed for use by customers only in strict accordance with the license agreement governing its use.

Any reports or other figures that appear in this document are merely illustrative and do not contain names or data of real people.

Any similarity in names of people, living or dead, or in data is strictly coincidental and is expressly disclaimed.

Siemens Healthineers does not warrant that the material contained in its documentation is error-free.

Documentation supplied to Siemens Healthineers by third parties and included with this documentation is not warranted for accuracy or completeness. The information contained in this document is subject to change. Revisions and updates will be issued from time to time to document changes and/or additions.

4.1 User Profile

Please note that the following profiles may vary in practice depending on (hospital) organization, qualification, and personal responsibilities and can only be considered as a general guide.

The following profiles may vary depending on the site, qualifications, and personal responsibilities. They should be considered general guides.

Technologist

Cath Lab Technologist: The cardiac technician, RCVT or RN assists in performing the cath lab exam. The tech is usually responsible for documenting the procedures performed in the cath lab, such as techniques and equipment used. He enters patient information associated with the report, records measurements, documents hemodynamics, and produces the procedural report.

Cardiac Sonographer: The sonographer performing the echo exam (US workflow model). The sonographer enters patient information associated with the report, performs measurements on images, and produces a preliminary report.

Vascular Technologist: The sonographer performing the exam. The sonographer enters patient information associated with the report, performs measurements on images, and produces a preliminary report.

Radiologic Technologist: The radiology technician captures the images for the exam and satisfies the procedure step for the patient.

OB/GYN and MFM Sonographer: The OB/GYN and MFM (Maternal Fetal Medicine) sonographer performs the exam and, in some cases, assists the physician with creating the reports.

- Physician**
- Cath Lab Cardiologist:** The cardiologist performing the cath lab exam and producing reports of the procedure for their records, for the patient, and for the referring physician.
 - Cardiac Surgeon:** The cardiac surgeon who reviews images and/or reports before surgery.
 - Echo Cardiologist:** The cardiologist performing the echo exam and producing reports of the procedure for their records, for the patient, and for the referring physician.
 - Vascular Lab Physician:** The physician performing the exam and producing reports of the procedure for their records, for the patient, and for the referring physician.
 - Radiologist:** The radiologist reads the images captured by the technician and produces reports for the referring physician and for their own records.
 - OB/GYN and MFM Physician:** The OB/GYN and MFM physician reads the exam and creates the reports.

- Clinical Administrator**
- The clinical administrator configures application settings, DICOM nodes, archiving, and manages short term storage and licenses, performs patient data administration, and is the first support contact for reading physicians.
- The clinical administrator is an application specialist who has work experience, product experience, and knowledge of the clinical workflow and is English speaking.

- IT Administrator**
- The IT administrator is responsible for data security, data protection, user management, configuration of DICOM nodes, as well as backup and archiving. The IT administrator manages the IT infrastructure of the clinical network and is responsible for the administration and configuration of IT components.
- The IT administrator must have expert knowledge of networks, operating systems, user administration, and knowledge of the workflow and is English speaking.

- Siemens Healthineers Service**
- The Siemens Healthineers Service user role is responsible to restore the normal system operation as fast as possible while minimizing adverse impact on the customer's business operation.

They are responsible for trouble shooting hardware and software issues identified by end customers.

They must have expert knowledge of networks, operating systems, user administration, and knowledge of the workflow.

They must have Product Experience, Professional Education or work experience and Product training certificate or comparable.

4.2 Training

syngo Dynamics can be used intuitively. Nevertheless, a training course tailored to your specific user situation is mandatory.

To enable users to routinely operate the *syngo* Dynamics system and give them the necessary theoretical and practical experience for the efficient usage, different kinds of training are available:

- User documentation packaged as part of the *syngo* Dynamics software package
- Integrated Help, which is part of the respective application
- Classroom trainings which take place at the Siemens Healthineers Training Centers
- On-site trainings which provide the most intensive training for you. The trainings are conducted by trainers from the local Siemens Healthineers organization. On-site trainings can be optionally extended depending on your specific needs.
- Healthcare professional will get access to PEPconnect – the Siemens Healthineers learning platform – so that they can expand and improve their knowledge at their own pace, anytime, anywhere. See pep.siemens-info.com, and search for *syngo* Dynamics to see relevant e-learning contents.
- See www.siemens-healthineers.com. Select your country and then select the category **Education & Training**.

If you require an additional training, please contact your local Siemens Healthineers Sales Representative.

This page is left blank intentionally.


5 Documentation Overview

The *syngo* Dynamics software offers several levels of user assistance, beginning with tooltips on the screen, and accompanied by the basic Online Help and Operator Manual.

The Online Help documents are available only in English and applicable for the tools and applications designed to be used with/for the *syngo* Dynamics software. The Operator Manuals, if available, may be localized and contain safety advisories and information on frequently-used user interaction.

5.1 Conventions Used in this Documentation

To help you find and understand the appropriate information, the following visual orientation aids are used.

- ◆ One-step operation
- 1 Operating sequence
- List item
- Guidelines
- Bold letters** Any user interface item, such as Window titles, menu items, function names, field names, and buttons are written in bold letters
-  Note for optimal handling
- Screenshots** Screenshots that are intended to illustrate general principles are shown in a neutral form for greater simplicity. Screenshots that illustrate instructions contain text to show the context within the application.

5.2 Labeling of Safety Advices

Cautions are specially marked in the documentation.

The content of a caution is structured in three different sections:

- Source of danger
- Consequence
- Countermeasure

Caution CAUTION indicates potential risk that may result in minor physical injury or material damage.



CAUTION

First, the source of danger is addressed!

Then, possible consequences are described.

◆ Finally, measures are given to prevent a dangerous situation.

5.3 Safety Advisory



All safety advisories that users need to consider for the use of *syngo* Dynamics are noted in the *syngo* Dynamics Safety Advisory manual. Refer to this manual specific to the software version installed at your hospital.

Additionally, all safety advisories applicable in an application workflow are indicated in the respective manual where the workflow is described.

5.4 About this Manual

This manual contains information on HL7 specifications required for interfacing with *syngo* Dynamics.

6 HL7 Overview

6.1 Purpose

This document provides the requirements to interoperate with *syngo* Dynamics using the Health Level Seven (HL7) standard. As per STIG V-76761, it is recommended that you disable TLS 1.1 and lower versions. It is also recommended that you enable TLS 1.2. All inbound HL7 communication is supported on TLS 1.2.

The objective of this document is to provide a standard set of Hospital Information System (HIS) requirements to support the following features:

- Report Upload

The *syngo* Dynamics system components communicate observation reporting information to the Cerner OPENLink Engine via a proprietary communications protocol. The Cerner OPENLink Engine constructs the HL7 message from information supplied by the *syngo* Dynamics workstation.

- ADT/Order Messages

This specification provides the *syngo* Dynamics baseline requirements for a Hospital Information System (HIS) to effectively send ADT/Order messages.

- HL7 Study Reconciliation of PHI (Protected Health Information) Data

When *syngo* Dynamics is integrated with Electronic Health Record System (EHR) and the necessary configurations exist in the **Server Configuration Tool (HL7 interface and Enable EHR are selected)**, Protected Health Information (PHI) data pertaining to studies are synchronized in both systems whenever there are updates to the study data.

On arrival of a HL7 message to *syngo* Dynamics Server from the Electronic Health Record System (EHR), the HL7 message is cached and queued. The queue is continuously processed for reconciliation.

ORM reconciliation

For an ORM message, the images are sent to *syngo* Dynamics from EHR PACS server. The patient information from the images is associated and reconciled to a matching order (ORM) based on the Study Instance UID or ORM matching rules configuration in the Site IDs tab in SysAdmin. For order association, a study must undergo successful demographic reconciliation. Note that the order association is based on a study date range which is configured on the server for all departments. Contact Siemens Healthineers Customer Care for study date range configuration. When the study type is modified on the worksheet and the study has had demographic reconciliation, the system updates the order association based on the updated study type. If configured, the user is prevented from verifying a report and/or marking the study as read when an order cancellation has been received for the order associated with the study.

Procedure type to department mapping

Mapping procedure type to department is a prerequisite for the studies (pre-procedure or new study from EHR) to be associated to a *syngo* Dynamics department. The ORM message contains the procedure type information. The study is assigned to the configured department based on the procedure type to department mapping in the **Server Configuration Tool**. If no mapping exists, the study is assigned to Default department. Contact your system administrator if this is not done. Your system administrator must contact Siemens Healthineers Customer Care to ensure all "default" studies are assigned to a specific department.

Study correction

If the study is associated with an incorrect order that needs to be corrected, *syngo* Dynamics supports study UID correction via HL7 ORM 001 messages sent by EHR. The ZSU segment contains details of the old and new study UID. The study UID correction is done on the study with the old study UID. For studies that went through study UID correction, the studies are associated with a new order ORM message. The match is based on the new study UID or Accession number present in the ORM message sent for the study correction. If there are multiple matching ORM messages, then the latest ORM is associated with the study. The studies are then reconciled with the new order ORM message, which in turn might lead to patient-level demographics update.

If the new study UID from the study UID ORM message exists in *syngo* Dynamics, the old study is merged to the new study. On successful study merge, all information (image data, image series data, and report) pertaining to the source study is moved to the target study. The merged studies are updated with the target study details.

Audit messages are logged when study UID correction is done and when a study is merged.

In the integrated environment, the following messages are displayed if a study is open with read-write access.

Study cannot be edited as study is merged to Patient ID <number> - Indicates that merge is in progress. The open source study changes to read-only mode.

Study merge has occurred, please reload the study - Indicates that the target study is open with read-write access, and merge has occurred.

In DICOM Image Review, you shall receive the same message for similar situations, but with the **OK** button.



If the study UID correction at your site is via HL7 ORM message, ensure that EHRHL7BasedStudyUpdate is set to **TRUE** in HL7CacheConfiguration.xml. This .xml file is present in

<Server>\F:\Acuson\Configuration\HL7Cache\HL7CacheConfiguration.xml.

By default, the EHRHL7BasedStudyUpdate is set to FALSE. This can be retained to be FALSE if the study UID correction is via DICOM.

Patient merge

In situations where the synchronization warrants patient data to be merged with existing patient, the ADT HL7 message when sent with merge segment is processed only if a matching study is found on the *syngo* Dynamics server.

The matching study is identified based on the matching rules specified in the HL7PatientIdentification.xml on the server in the location <Server>\\F:\Acuson\Configuration. By default, **Patient ID** and **Site ID** are used for Patient Identification. The System Administrator can modify the HL7PatientIdentification.xml to modify the default values to contain only **Patient ID**, as applicable for the hospital site. When a match is found, all information (image data, image series data, and report) pertaining to the source patient is moved to the target patient. The merged patient demographics are updated with the target patient demographics. Audit messages are logged when patient information is merged. When studies being merged are already open in *syngo* Dynamics, the following message is displayed indicating about the merge – **Patient demographics updated and able to continue with the reading of the study**. Patient level demographics are updated for both current and prior studies. Audit messages are logged when patients are merged.

Errors during merge

If the merge fails (report merge or image merge), the source study is retained. If report merge succeeds and image merge fails, the target study contains the merged reports along with the images (partial set). Proper logs are maintained for such merge failures.

Patient Demographics Update

PHI data of the studies in *syngo* Dynamics is automatically reconciled with the corresponding procedure of the EHR through a HL7 ADT message from EHR. Automatic reconciliation is done only for the EHR linked studies.

The HL7 message is processed only if a matching study is found on the *syngo* Dynamics server. The matching study is identified based on the patient identification fields configured in the HL7PatientIdentification.xml on the server in the location <Server>\F:\Acuson\Configuration. Matching criteria can be **Patient ID** and **Site ID/Patient ID** only. The **Site ID** should be configured for a department in SysAdmin. By default, **Patient ID** and **Site ID** is used for patient identification. On successful processing of the HL7 message, PHI data of the study automatically reconciles with the updates from the HL7 message.

During reconciliation, patient level demographics are updated for both current and prior studies. When reports are opened in Reporting component, demographic updates can be seen on reports if the reports are in unlocked state. For a locked report, users must unlock the report to see the demographics updates.

The following fields are auto reconciled as part of HL7 ADT update: **Last Name, First Name, Middle Name, Prefix, Suffix, Gender, Date of Birth.**

When a study being updated with the PHI is already opened, Reporting is automatically refreshed with the updated PHI data and the user is notified with the message - **Patient Demographics Updated.**



The demographic data is not updated for the following conditions.

- Study is marked as Read with lock enabled.
- Report is locked and Verified.

Using the HL7 Reconciliation Tool

Use the HL7 Reconciliation Tool to reconcile ORM and ADT messages for non-EHR studies. You may also use the HL7 Reconciliation Tool to manually reconcile ORM/ADT message updates for EHR studies for which reconciliation has failed.

6.2 Scope

This document specifically refers to *syngo* Dynamics VA41 server and workstations and the supported Cerner OPENLink component version.

6.3 Abbreviations

(DR)	Short for Demographic Reconciliation in this document.
(OA)	Short for Order Association in this document.
ACK	Acknowledgment
ADT	Patient administration messages within the HL7 standard. Typically initiated by the HIS to inform ancillary systems that a patient has been admitted, discharged, transferred, merged, that other demographic data about the patient has changed (name, insurance, next of kin, etc.) or that some visit information has changed (patient location, attending physician, etc.).
ADT^A01	An HL7 Admit/Discharge/Transfer message for Patient Admit. This occurs when the patient is admitted to the hospital.
ADT^A04	An HL7 Admit/Discharge/Transfer message for Patient Registered. This occurs when the patient is registered at the hospital.
ADT^A08, ADT^A31	An HL7 Admit/Discharge/Transfer message for Patient Update. This occurs when the patient's demographic information has been updated from what was previously sent when the patient was admitted. An A31 event can be used to update person information on an MPI. It is similar to an A08 (update patient information) event, but an A08 (update patient information) event should be used to update patient information for a current episode. A31 can also be used for back-loading MPI information for the patient or for back-loading patient and historical information.
ADT^A18, A34, A40	An HL7 Admit/Discharge/Transfer message for Patient Merge. These are used when two different identifications for the same patient are found and need to be merged into one identification.
AL1	Allergy Information
All servers	For the purposes of this document, the term "all servers" refers to all <i>syngo</i> Dynamics servers that are included in the multi-server login configuration of the current server.
CA	The order control code for the order that is canceled. It is used in ORC segment of HL7 ORM Message
DDI	Data Dictionary Item

DG1	Diagnosis
DOB	Patient date of birth
Enterprise	<p>For the purposes of this document, the term “enterprise” refers to the collection of all departments on all servers that are included in the multi-server login configuration of the current server.</p> <p>The scope of “enterprise” is not defined as all departments of all <i>syngo</i> Dynamics servers that are physically located in an institution.</p> <p>Each <i>syngo</i> Dynamics server within an institution server pool has its own enterprise scope, as the scope is defined by the multi-server login configuration of that individual server.</p>
FON	Short for Filler Order Number in this document - The order number assigned to the Service Request by the party filling the order.
HIS	Health Information System
HIS Verification	A method used by the <i>syngo</i> Dynamics Server’s Vault services to reconcile patient/study information by using a DICOM Modality Worklist Server to confirm scheduled study data.
HL7	Specific message standards created by the Health Level Seven organization for the exchange, management and integration of electronic healthcare information.
IE	Interface Engine
IP	Internet Protocol
M	Mandatory
MDM	<p>Medical Document Management</p> <p>Helps managing medical records by transmitting new or updated documents, or by transmitting important status information and/or updates for the record.</p>
Master Patient Index	<p>Master Patient Index (MPI) is an index which references all existing patients registered at a health care organization in such manner that each patient is logically represented only once and can be uniquely identified. Use of MPI allows clinical and demographic information to be cross-referenced between different facilities within the enterprise of a health care organization.</p> <p>MPI can be a specific SID/PID pair among all SID/PID pairs that are assigned to a patient. For example, for a health care organization that defines “EPI” as the MPI SID, the value “123456789^^^EPI” indicates it is a MPI and “123456789” is the MPI PID.</p>

MPI	<p>Master Patient Index (MPI) is an index which references all existing patients registered at a health care organization in such manner that each patient is logically represented only once and can be uniquely identified. Use of MPI allows clinical and demographic information to be cross-referenced between different facilities within the enterprise of a health care organization.</p> <p>MPI can be a specific SID/PID pair among all SID/PID pairs that are assigned to a patient. For example, for a health care organization that defines "EPI" as the MPI SID, the value "123456789^^^EPI" indicates it is a MPI and "123456789" is the MPI PID.</p>
MRG	Merge
MRN	Medical Record Number in HL7 (Same as Patient ID in DICOM)
MSH	Message Header
NACK	Negative Acknowledgment
NK1	Next of Kin
NW	The order control code for new order that is used in ORC segment of HL7 ORM Message
O	Optional
OBR	The Observation Request segment that HL7 ORM message uses to transmit information specific to an order for a diagnostic study or observation, physical exam, or assessment.
ORC	The Common Order segment that HL7 ORM message uses to transmit fields that are common to all orders
ORM^O01	An HL7 General Order Message. Also referred to as an "Order" or "HL7 Order" in this document.
ORU	Observation Result that provides clinical observations in response to an order. ORU messages provide structured patient-oriented clinical data between systems.
PID	Patient ID or refers to the patient identification segment in a HL7 message.
PON	Short for Placer Order Number in this document - the order number assigned to the Service Request by the party placing the order.
PV1	Patient visit segment in a HL7 message which includes the patient class.
RP-RO	The order control code used for the order that is replacing the previous one. It is used in ORC segment of HL7 ORM Message. For each replacement order, ORC segment use ORC-1-order control value of RP, followed by another ORC segment with ORC-1-order control value of RO.
SID	Site ID. Another name for an Assigning Authority (e.g. an Information System) from the HL7 Standard.
SID/PID	A Site ID and Patient ID combination.

SIU	Scheduling Information Unsolicited
UI	User Interface

This page is left blank intentionally.

7 HL7 Study Reconciliation of PHI Data

Access the **HL7 Reconciliation Tool** to associate orders with studies and reconcile studies. For complete information on how to perform the HL7 study reconciliation of PHI (Protected Health Information) data, refer to the *syngo* Dynamics HL7 Reconciliation Tool Manual specific to the software version installed at your hospital.

These order associations and demographic reconciliations can also be broken.

Break order association

In *syngo* Dynamics, an associated order with a study is broken in the following situations.

- An ORM “cancel” order is received and the **BreakOrderReconOnCancelOrder** variable is set to **True** in `sdConfig.xml` file.
- Study type is modified.
However, if a new ORM is found that matches the study type, the study is re-associated automatically with the new order.
- Break the order association manually from Study List. For more information, refer to the *syngo* Dynamics Workplace Manual specific to the software version installed at your hospital.

Break reconciliation

In *syngo* Dynamics, a demographic reconciliation that was done in error can be broken manually from Study List. When a reconciliation is broken, all order associations with the study will be broken too. For more information, refer to the *syngo* Dynamics Workplace Manual specific to the software version installed at your hospital.



CAUTION

Demographic data which was stored in a *syngo* Dynamics study during reconciliation is replaced with the demographic data from Sensis when a *syngo* Dynamics study that is linked with Sensis is opened at the Workplace. Incorrect prior studies may be identified when the study is displayed in *syngo* Dynamics if Sensis and *syngo* Dynamics do not receive HL7 messages from the EHR at the same time.

Misdiagnosis: Incorrect prior studies are displayed for a sensis study which is opened in sD when the Sensis changes the demographic data in the study.

- ◆ To prevent the display of incorrect priors and potential misdiagnosis, mark the study as read in *syngo* Dynamics to persist the Sensis data in the *syngo* Dynamics database and then use the **Break Reconciliation** feature and the HL7 Reconciliation Tool to correct the patient demographic data.

8 Pre-procedural Documentation

Pre-procedural documentation help you document patient-related data ahead of the planned procedure. Pre-procedural documentation cannot be done for existing studies.

When *syngo* Dynamics is integrated with an Electronic Health Record System (EHR) and the necessary configurations exist in the **Server Configuration Tool**, you can create pre-procedural documentation for the matching inbound ORM messages. Information from the ORM messages are mapped to the *syngo* Dynamics database and the corresponding values updated in the *syngo* Dynamics worksheet and report. The matching is based on Study Instance UID or Accession Number. This helps in pre-procedural reporting.

Additionally, in an EHR-integrated environment, if UTC (Coordinated Universal Time) is included in the Observation Date/Time field in the HL7 message, the time is converted per the Server time and associated with the pre-procedural documentation.

Refer to the following:

- For further information on order mapping, refer to “Order Messages (ORM) Mapping” on page 57.
- For more information on pre-procedural documentation, refer to *syngo* Dynamics System Administration Manual specific to the software version installed at your hospital.
- For information on entering pre-procedural documentation for unscheduled procedures, refer to *syngo* Dynamics Workplace Manual specific to the software version installed at your hospital.

This page is left blank intentionally.

9 Report Upload to OPENLink Option

9.1 HL7 Message Details Message Specifications

The HL7 outgoing observation reporting interface on the Cerner OPENLink Engine allows the *syngo* Dynamics system to transmit structured patient-oriented clinical data to external systems via unsolicited messages. The standard used by the Cerner OPENLink Engine is a subset of the HL7 Standard, version 2.3.1.

The Observation Reporting transaction set is used to transmit results of diagnostic studies from the *syngo* Dynamics system to interested information systems. Transactions are initiated by the filler application using the ORU message to provide the recorded observation. A single trigger event is used by the Observation Reporting transaction set.

R01: Unsolicited Transmission of Observation

This message is generated by the Cerner OPENLink Engine to provide reports from the *syngo* Dynamics system to the Information System. The *syngo* Dynamics workstation provides the report information to the Cerner OPENLink Engine using a proprietary communication interface.

The following table lists the ORU Message Type Events.

Event Code	Event Description	Message
R01	Unsolicited Transmission of Observation Report (Preliminary, Final, or Corrected)	MSH PID OBR {OBX}

9.2 syngo Dynamics Report Format

The *syngo* Dynamics Report is generated in an XML format prior to being sent to the Cerner OPENLink Engine. Tabs are converted to spaces (0x20) in the XML output. The linefeed character (0x0A) is replaced with a value of \.br\ by the OPENLink Engine to indicate new lines within the formatted report for the HL7 output.

The syngo Dynamics Report may contain HL7 message delimiters in the report text. To ensure an HL7 compliant character set is sent by the Cerner OPENLink Engine, the delimiter characters will be replaced using the table below. The default values are shown. Site-defined substituted character strings are supported.

The following table contains information on the HL7 Delimiter Substitution.

Message Delimiter	ASCII Character	Substituted Character String
Field Separator		\F\
Component Separator	^	\S\
Subcomponent Separator	&	\T\
Repetition Separator	~	\R\
Escape Character	\	\E\

The entire report text is inserted into a single OBX segment. This must be supported by the receiving system.

9.3 Default Segment Listing

The following tables list the HL7 segments and their corresponding data elements that are used by this Report Upload Interface.

Message Header Segment (MSH):

Seq	Type	Element Name	<i>syngo</i> Dynamics Element	Notes
1	ST	Field Separator	NA	Default Value:
2	ST	Encoding Characters	NA	Default Value: ^~\&
3	HD	Sending Application	NA	"KinetDx"
4	HD	Sending Facility	KinetDxHospitalName	Configured <i>syngo</i> Dynamics Hospital Name*
5	HD	Receiving Application	NA	Optional
6	HD	Receiving Facility	NA	Optional
7	TS	Date/Time of Message	NA	NA
9	CM	Message Type	NA	ORU^R01
10	ST	Message Control ID	NA	Unique number
11	PT	Processing ID	NA	P sent as default
12	ID	Version ID	NA	2.4 - for ORU 2.5 - for MDM
18	ID	Character Set	NA	Default Value: ISO_IR 100

*This field is currently populated with the study institution name.

Patient Identification Segment (PID):

Seq	Type	Element Name	<i>syngo</i> Dynamics Element	Notes
1	SI	Set ID	NA	1
3	CX	Patient ID (external)	Patientid	DICOM Patient ID
5	XP	Patient Name	PatientName	DICOM Patient Name*
7	TS	Date/time of birth	PatientBirthDate	
8	IS	Sex	PatientSex	M, F, or O

*This field is currently sent with the DICOM PN encoding rules.

Observation Request Segment (OBR):

Seq	Type	Element Name	<i>syngo</i> Dynamics Element	Notes
1	SI	Set ID	NA	1
3	EI	Filler Order Number	AccessionNumber	DICOM Accession Number
4	CE	Universal Service ID	'^' + StudyDescription	DICOM Scheduled Procedure Step Description. OPENLink to concatenate and add SubComponent Separator

Seq	Type	Element Name	<i>syngo</i> Dynamics Element	Notes
7	TS	Observation Date/Time	StudyDate + StudyTime	DICOM Study Date + Study Time (Concatenate Values)
22	TS	Results rpt/status change date/time	ReportDate + ReportTime	Date/Time Report uploaded from Workstation (Concatenate Values)
25	ID	Results Status	ReportStatus	F – If <i>syngo</i> Dynamics status is Read P – If <i>syngo</i> Dynamics status is Completed C – If <i>syngo</i> Dynamics report has been amended or corrected
32	CM	Principal result interpreter	PhysicianReadingStudy	<i>syngo</i> Dynamics Diagnosing Physician Name*

*This field is currently sent with the DICOM PN encoding rules.

OBX – Observation/Result Segments:**OBX1 – Impressions**

Seq	Type	Element Name	<i>syngo</i> Dynamics Element	Notes
1	SI	Set ID	NA	1
2	ID	Value Type	CE	Default value depending on the report section CE = Impression
3	CE	Observation Identifier	AccessionNumber&IMP	OPENLink to concatenate and add SubComponent Separator
5	*	Observation Value	Text	Contains <i>syngo</i> Dynamics Impressions Text. OPENLink to concatenate and add SubComponent Separator. Text value should be "HL7 Safe"
11	ID	Observation Result Status	ReportStatus	Same as OBR-25 (F, P, C)

OBX2 – Report Text

Seq	Type	Element Name	<i>syngo</i> Dynamics Element	Notes
1	SI	Set ID	NA	2
2	ID	Value Type	FT	Default value depending on the report section FT = Report or Addendum

Seq	Type	Element Name	<i>syngo</i> Dynamics Element	Notes
3	CE	Observation Identifier	AccessionNumber&GDT	OPENLink to concatenate and add SubComponent Separator
5	*	Observation Value	Text	Contains <i>syngo</i> Dynamics Report Text. Text value should be "HL7 Safe"
11	ID	Observation Result Status	ReportStatus	Same as OBR-25 (F, P, C)

OBX3 – Addendum

Seq	Type	Element Name	Notes
2	ID	Value Type	FT
3	CE	Observation Identifier	DICOM Accession Number&ADT
5	*	Observation Value	Contains <i>syngo</i> Dynamics Addendum Text
11	ID	Observation Result Status	Same as OBR-25 (F, P, C)

9.4 Report Upload Message Examples

The following examples contain information on syngo Dynamics Report with standard linefeed substitution.

Final Report MSH|^~|&|KinetDx|Ann Arbor Lab|Broker|Connectivity
Lab|20030305105716||ORU^R01|00002455|P|2.2|||||ISO_IR 100
PID||M555||SUMNER^ALEX^WILLIAM||19680212|M|
OBR|1||275|^Adult
Heart|||20030305120744|||||||||20030305162300|||F|||||
|Welby^Marcus^T
OBX|1|CE|275&IMP||^Acute MI unspecified site, initial episode of care-
410.91|||||F
OBX|2|FT|275&GDT||Ann Arbor Medical Center|.br\Echocardiography
Laboratory|.br\Telephone: (734)-555-1234 Fax: (734)-555-
1235|.br|.br\ECHO REPORT|.br\Patient Name: ALEX WILLIAM SUMNER Patient
ID: M555|.br\Exam Date: March 5, 2003 Patient DOB: February 12,
1968|.br\Date of Report: March 5, 2003 Height: \.br\Referring Phys. YEVES
ANDREW. TIESSEN Weight: \.br\Referring Diag. Acute MI unspecified site, initial
episode of care-410.91 \.br\Resting BP: / \.br\Indication: Kinet, Acute MI
anterolateral, subsequent episode of care-410.02|.br|.br\IM-mode
Data|.br\LVDD: AO: \.br\LVSD: EPSS: \.br\Fract. Shortn: PW: \.br\RVDD: EDV
Teichholz: \.br\IVSTD: ESV Teichholz: \.br\IVSTS: SV: \.br\Fract. Change: EF:
\.br\LAD: \.br|.br\2D Study|.br\There is mild left atrial enlargement. \.br\There
is apical akinesis with wall thinning compatible with scar. The|.br\distribution
suggests a distal left anterior descending coronary artery|.br\distribution.
\.br|.br|.br|.br\Doppler|.br\There is trivial aortic regurgitation present.
\.br|.br|.br|.br|.br|_____|.br\Marcus T
Welby|.br|.br|F

Corrected Report MSH|^~|&|KinetDx|Ann Arbor Medical Center|Broker|Connectivity
Lab|20030305105716||ORU^R01|00002455|P|2.2|||||ISo_IR 100
PID||M555||SUMNER^ALEX^WILLIAM||19680212|M|
OBR|1||275|^Adult
Heart|||20030305120744|||||||||20030305162300|||C|||||
|Welby^Marcus^T
OBX|1|CE|275&IMP||^Acute MI unspecified site, initial episode of care-
410.91|||||C

OBX|2|FT|275&GDT||Ann Arbor Medical Center.l.br|Echocardiography
Laboratory.l.br|Telephone: (734)-555-1234 Fax: (734)-555-
1235.l.br|.br|ECHO REPORT.l.br|Patient Name: ALEX WILLIAM SUMNER Patient
ID: M555.l.br|Exam Date: March 5, 2003 Patient DOB: February 12,
1968.l.br|Date of Report: March 5, 2003 Height: .l.br|Referring Phys. YEVES
ANDREW. TIESSEN Weight: .l.br|Referring Diag. Acute MI unspecified site, initial
episode of care-410.91 .l.br|Resting BP: / .l.br|Indication: Kinet, Acute MI
anterolateral, subsequent episode of care-410.02.l.br|.br|M-mode
Data.l.br|LVDD: AO: .l.br|LVSD: EPSS: .l.br|Fract. Shortn: PW: .l.br|RVDD: EDV
Teichholz: .l.br|IVSTD: ESV Teichholz: .l.br|IVSTS: SV: .l.br|Fract. Change: EF:
.l.br|LAD: .l.br|.br|2D Study.l.br|There is mild left atrial enlargement. .l.br|There
is apical akinesis with wall thinning compatible with scar. The.l.br|distribution
suggests a distal left anterior descending coronary artery.l.br|distribution.
.l.br|.br|.br|.br|Doppler.l.br|There is trivial aortic regurgitation present.
.l.br|.br|.br|.br|.br|.br|_____l.br|Marcus T
Welby.l.br|.br|C

OBX|3|FT|275&ADT||April 23, 2003 TIESEN, YEVE\$|.br|Further
examination of this study indicates the possibility for some|.br|additional
aortic regurgitation|.br|C


This page is left blank intentionally.

10 IHE DRPT Report Creator

The Displayable Reports Profile specifies transactions supporting the creation, revision, intra/inter-department transmission, and reading of display-ready clinical reports. In the imaging procedure context, it provides linkage between the report, the imaging and other evidence of the procedure. The report is provided to users outside the department for broad distribution (e.g., using other profiles).

The DRPT Profile requires the use of the Portable Document Format (PDF), which has emerged as a means of encoding documents ready for presentation, including graphical content. For reporting on imaging procedures, especially on cardiac procedures, PDF is able to present the full range of documentation generated by a wide variety of reporting packages. Furthermore, PDF allows the reporting physician to control the “look” of the report, which is important for both clinical and business reasons.

syngo Dynamics is a Report Creator – A system that generates and transmits clinical reports.



Complete details of the IHE DRPT implementation requirements can be found at: http://www.ihe.net/Technical_Framework/index.cfm#cardiology.

10.1 IHE DRPT Message Examples

Below is a sequence sample of a test study from preliminary to final to amended in *syngo* Dynamics and the corresponding MDM T2 and MDM T10 messages.

```
MSH|^~\&|syngo Dynamics|Hospital Name|Receiving Application|Receiving
Facility|201011151311||MDM^T02|6|P|2.3.1|||||8859/1
EVN||201011151311 PID|||drpt2||sample^drpt^^^mr OBR|1|^syngo
Dynamics||^Adult Cardiac|||20101115131051|||||||||||||R
TXA|1001|DI|AP|||20101115131124.431||||1.3.12.2.1107.5.8.11.105.1777
6024256.374677893148344||||PA
OBX|1|HD|^113014^DICOMStudy^DCM||1.3.12.2.1107.5.8.11.105.177760
24256.374677874100265||||O
OBX|2|ED|^Physician||^Application^PDF^Base64^JVBeriOxLjENCg0KMyAwIG
9iaG0CAglCAglCAglCA0YNCg==||||R
```

```
MSH|^~\&|syngo Dynamics|Hospital Name|Receiving Application|Receiving
Facility|201011151311||MDM^T10|9|P|2.3.1|||||8859/1
EVN||201011151311 PID|||drpt2||sample^drpt^^^mr OBR|1|^syngo
Dynamics||^Adult Cardiac|||20101115131051|||||||||||||F
TXA|1002|DI|AP|||20101115131147.228||||^dehar^todd|1.3.12.2.1107.5.8.
11.105.17776024256.374677909041117|1.3.12.2.1107.5.8.11.105.17776
024256.374677893148344||||AU||||^dehar^todd^^^^^^^^^^^^^20101115
131147.228
OBX|1|HD|^113014^DICOMStudy^DCM||1.3.12.2.1107.5.8.11.105.177760
24256.374677874100265||||O
OBX|2|ED|^Physician||^Application^PDF^Base64^JVBeriOxLjENCg0KMyAwIG
9iaCg==||||F
```

```
MSH|^~\&|syngo Dynamics|Hospital Name|Receiving Application|Receiving
Facility|201011151312||MDM^T10|12|P|2.3.1|||||8859/1
EVN||201011151312 PID|||drpt2||sample^drpt^^^mr OBR|1|^syngo
Dynamics||^Adult Cardiac|||20101115131051|||||||||||||C
TXA|1003|DI|AP|||20101115131147.228||||^dehar^todd|1.3.12.2.1107.5.8.
11.105.17776024256.374677944757801|1.3.12.2.1107.5.8.11.105.17776
024256.374677909041117||||AU||||^dehar^todd^^^^^^^^^^^^^20101115
131147.228
OBX|1|HD|^113014^DICOMStudy^DCM||1.3.12.2.1107.5.8.11.105.177760
24256.374677874100265||||O
OBX|2|ED|^Physician||^Application^PDF^Base64^JVBeriOxLjENCg0KMTEgM
CBvYmoQo=||||C
```

11 Detailed Financial Transaction (DFT) Messages

As a patient is treated, the clinical data is transferred to the appropriate departments throughout the hospital, including billing. The main HL7 messages transferred to billing are as follows:

- ADT - Message containing patient demographics
- DFT - Message containing detailed financial transactions
- ORU - Message containing unsolicited observation result

Every healthcare environment bills for services provided. The goal of billing is to increase reimbursement in a timely fashion. Automating clinical data delivery through HL7 interfaces provides accurate information quickly to the billing department. Delivery of billing information needs to happen periodically – most often on a daily basis – and can be accomplished through an off-hours batch method. This is done using the DFT message.

The DFT message describes a financial transaction that is sent to a billing system and is used for patient accounting purposes. This message is used to send charges from *syngo* Dynamics to a Hospital Information System (HIS) or a Radiology Information System (RIS). This message may include ancillary charges or patient deposits, and is sent between the DSS/Order Filler and the Charge Processor. The DSS/Order Filler then verifies whether the procedure is complete. Credit, void, or other adjustment transactions are not supported.

The segments and groups of segments in a DFT-P03 message include the following:

Segment/ Group	Name	Optional/Repeatable	Remarks
MSH	Message header	Required	NA
EVN	Event type	Required	NA
PID	Patient identification	Required	NA
PID1	Patient demographics	Optional	NA
PV1	Patient visit	Optional	NA
PV2	Patient visit – additional information	Optional	NA

Segment/ Group	Name	Optional/Repeatable	Remarks
DB1	Disability segment	Optional, Repeatable	NA
OBX	Observation	Optional, Repeatable	NA
FT1	Financial Transaction	Required	Transaction group
PR1	Procedures	Required	Procedure group, part of Transaction group
ROL	Role	Optional, Repeatable	Procedure group, part of Transaction group
IN1	Insurance	Repeatable	Insurance group
IN2	Insurance additional info	Optional	Insurance group
IN3	Insurance additional info – certification	Optional	Insurance group
DG1	Diagnosis	Optional, Repeatable	Additional segment, not part of a group
DRG	Diagnosis Related Group	Optional	Additional segment, not part of a group
GT1	Guarantor	Optional, Repeatable	Additional segment, not part of a group
ACC	Accident	Optional	Additional segment, not part of a group

The key segments containing data are as follows:

- PID: Patient Information
- PV1: Visit Information (DOS)
- FT1: Financial Transaction
- IN1: Insurance Information
- IN2: Additional Insurance Information
- IN3: Additional Insurance Info, Certification
- GT1: Guarantor
- AUT: Authorization Information

Message Processing

syngo Dynamics generates a charge file for every verified report in each department that has the billing option enabled. The files are created once a day, typically at day-end, and are sent to the billing upload folder for retrieval by the OPENLink DFT interface instance on *syngo* Dynamics or Magic Link A. The OPENLink interface then parses the charge files and converts them into HL7 DFT transactions. The DFT transactions are transmitted outbound by OPENLink via a TCP/IP socket connection using HL7 Minimal Lower Level Protocol (MLLP) transaction framing. Once the DFT transactions are submitted to the destination system, the source charge files are deleted.

Example: DFT Message

The following message shows the addition of a test, Lipid Panel, to a patient's bill.

```
MSH|^~\&|FrapLab|StJohn|HIS|StJohn|20071217094822||DFT^P03|MSGIDv|P
|2.3
PID||3|82828||Simpson^Margaret^^^Mrs.||19650525|F|||12 Maple
St.^Springfield^OH^21003^USA
PV1|4|^22^1|||2360^England^Mikey|||IP|||||4|||||||20071217
094755|20071217094813
FT1|1|6|4|20071217094821||Credit|303756^Lipid Panel^L|||2|115
```

The following message shows a single FT1 segment.

```
MSH|^~\&|SYNGODYNAMICS|Hospital
Name|||201506172300||DFT^P03|15289507743831|P|2.4
EVN|P03|20150617230008
PID|1||124789||SAMPLE^TESTING||19541215|F|||||||78451239
OBR|1||8574963|||20150617162026
FT1|1|8574963||20150617||CG|456|||1|||||||123^EXPORT^DOC||||99307^
Echocardiography, transthoracic
```

The following message shows multiple FT1 segments.

```
MSH|^~\&|SYNGODYNAMICS|Hospital
Name|||201506182300||DFT^P03|15289519973811|P|2.4
EVN|P03|20150618230047

PID|1||2424||TEST3^TESTING||19750909|M|||||||654654
OBR|1||5467|||||20150618104603
FT1|1|5467||20150618||CG|456|||1|||||||123^EXPORT^DOC|44^REQUESTIN
G^DOC||||93306^Echocardiography, transthoracic real time image
FT1|2|5467||20150618||CG|456|||1|||||||123^EXPORT^DOC|44^REQUESTIN
G^DOC||||93307^Echocardiography, transthoracic
FT1|3|5467||20150618||CG|456|||1|||||||123^EXPORT^DOC|44^REQUESTIN
G^DOC||||92982^Percutaneous transluminal coronary angioplasty (PTCA),
single vessel
```

11.1 Default Segment Listing

The tables listed below describe the required and optional HL7 data elements, the corresponding *syngo* Dynamics database fields, and the default values used by this interface. Elements that have an “N” in the “OPT HL7/*syngo* Dynamics” column are not supported in this interface by *syngo* Dynamics.

11.1.1 MSH Segment

Seq	Len	DT	OPT HL7/ <i>syngo</i> Dynamics	Element Name	DB Field in <i>syngo</i> Dynamics/ Comments
1	1	ST	R	Field Separator	
2	4	ST	R	Encoding Characters	^~\&
3	180	HD	O	Sending Application	Default: SYNGODYNAMICS
4	180	HD	O	Sending Facility	[Hospital Name]
5	180	HD	O	Receiving Application	User defined
6	180	HD	O	Receiving Facility	User defined
7	26	TS	R	Date/Time of Message	Date/Time of Message
8	40	ST	O	Security	NA
9	13	CM	R	Message Type	DFT^P03

Seq	Len	DT	OPT HL7/ <i>syngo</i> Dynamics	Element Name	DB Field in <i>syngo</i> Dynamics/ Comments
10	20	ST	R	Message Control ID	Unique ID
11	3	PT	R	Processing ID	P (configurable)
12	60	VID	R	Version ID	HL7 Version (2.4)

11.1.2 EVN Segment

Seq	Len	DT	OPT HL7/ <i>syngo</i> Dynamics	Element Name	DB Field in <i>syngo</i> Dynamics/ Comments
1	3	ID	B	Event Type Code	P03
2	26	TS	R	Recorded Date/Time	Date/Time of Message

11.1.3 PID Segment

Seq	Len	DT	OPT HL7/ <i>syngo</i> Dynamics	Element Name	DB Field in <i>syngo</i> Dynamics/ Comments
1	4	SI	O	Set ID	1
2	20	CX	B	Patient ID	NA
3	250	CX	R	Patient ID List	Patient ID only
4	20	CX	B	Alternate Patient ID	NA
5	250	XP	R	Patient Name	Patient Name
6	250	XP	O/N	Mother's Maiden Name	NA
7	26	TS	O	Date/Time of Birth	Date of Birth
8	1	IS	O	Administrative Sex	M, F, or O
9	250	XP	O/N	Patient Alias	NA
10	250	CE	O/N	Race	NA
11	250	XAD	O/N	Patient Address	NA

Seq	Len	DT	OPT HL7/ syngo Dynamics	Element Name	DB Field in syngo Dynamics/ Comments
12	4	IS	O/N	Country Code	NA
13	250	XTN	O/N	Phone – Home	NA
14	250	XTN	O/N	Phone – Business	NA
15	250	CE	O/N	Primary Language	NA
16	250	CE	O/N	Marital Status	NA
17	250	CE	O/N	Religion	NA
18	250	CX	O/R	Patient Account Number	Admission ID (Can also be sent in PV1-19)

11.1.4 PV1 Segment (Optional)*

Seq	Len	DT	OPT HL7/ syngo Dynamics	Element Name	DB Field in syngo Dynamics/ Comments
1	4	SI	O	Set ID	1
2	1	IS	R/O	Patient Class	Default “U”
3	80	PL	O	Patient Location	NA
4	2	IS	O/N	Admission Type	NA
5	250	CX	O/N	Preadmit Number	NA
6	80	PL	O/N	Prior Patient Location	NA
7	250	XCN	O/N	Attending Doctor	NA
8	250	XCN	O	Referring Doctor	Referring Physician
9	250	XCN	O/N	Consulting Doctor	NA
10	3	IS	O/N	Hospital Service	NA
11	80	PL	O/N	Temporary Location	NA
12	2	IS	O/N	Preadmit Test Indicator	NA
13	2	IS	O/N	Re-admission Indicator	NA
14	6	IS	O/N	Admit Source	NA

Seq	Len	DT	OPT HL7/ <i>syngo</i> Dynamics	Element Name	DB Field in <i>syngo</i> Dynamics/ Comments
15	2	IS	O/N	Ambulatory Status	NA
16	2	IS	O/N	VIP Indicator	NA
17	250	XCN	O/N	Admitting Doctor	NA
18	2	IS	O/N	Patient Type	NA
19	250	CX	O	Visit Number	Admission ID (Can also be sent in PID-18)

* This segment is not enabled by default.

11.1.5 OBR Segment (Optional)

Seq	Len	DT	OPT HL7/ <i>syngo</i> Dynamics	Element Name	DB Field in <i>syngo</i> Dynamics/ Comments
1	4	SI	O	Set ID	1
2	22	EI	C/O	Placer Order Number	NA
3	22	EI	C/O	Filler Order Number	NA
4	250	CE	R/O	Universal Service ID	External charge code linked to <i>syngo</i> Dynamics study type
5	2	ID	B/N	Priority	NA
6	26	TS	B/N	Requested Date/Time	NA
7	16	TS	C/O	Observation Date/Time	Study Date/Time
8	26	TS	O/N	Observation End Date/Time	NA
9	20	CQ	O/N	Collection Volume	NA
10	250	XCN	O/N	Collector Identifier	NA
11	1	ID	O/N	Specimen Action Code	NA
12	250	CE	O/N	Danger Code	NA
13	300	ST	O/N	Relevant Clinical Info	NA
14	26	TS	C/N	Specimen Received Date/Time	NA
15	300	CM	O/N	Specimen Source	NA
16	250	XCN	O	Ordering Provider	Requesting physician

11.1.6 FT1 Segment

Seq	Len	DT	OPT HL7/ <i>syngo</i> Dynamics	Element Name	DB Field in <i>syngo</i> Dynamics/ Comments
1	4	SI	O	Set ID	FT1 segment sequence number
2	12	ST	O	Transaction ID	Transaction ID

Seq	Len	DT	OPT HL7/ <i>syngo</i> Dynamics	Element Name	DB Field in <i>syngo</i> Dynamics/ Comments
3	10	ST	O/N	Transaction Batch ID	NA
4	26	TS	R	Transaction Date	NA
5	26	TS	O	Transaction Posting Date	NA
6	8	IS	R	Transaction Type	CG (Charge); credits and other adjustments not supported
7	250	CE	R	Transaction Code	External charge code linked to <i>syngo</i> Dynamics Study Type (e.g. Universal Service ID)
8	40	ST	B	Transaction Description	NA
9	40	ST	B/N	Transactions Description – Alternate	NA
10	6	NM	O	Transaction Quantity	Default “1”
11	12	CP	O/N	Transaction Amount – Extended	NA
12	12	CP	O/N	Transaction Amount – Unit	NA
13	250	CE	O/N	Department Code	NA
14	250	CE	O/N	Insurance Plan ID	NA
15	12	CP	O/N	Insurance Amount	NA
16	80	PL	O/N	Assigned Patient Location	NA
17	1	IS	O/N	Fee Schedule	NA
18	2	IS	O/N	Patient Type	NA
19	250	CE	O	Diagnosis Code	From <i>syngo</i> Dynamics Diagnosis library
20	250	XCN	O	Performed By Code	Billing ID & name of physician reading study
21	250	XCN	O	Ordered By Code	Requesting physician
22	12	CP	O/N	Unit Cost	NA
23	22	EI	O	Filler Order Number	NA
24	250	XCN	O	Entered By Code	NA

Seq	Len	DT	OPT HL7/ syngo Dynamics	Element Name	DB Field in syngo Dynamics/ Comments
25	250	CE	O	Procedure Code	From syngo Dynamics Procedure library
26	250	CE	O	Procedure Code Modifier	NA

11.1.7 PR1 Segment (Optional)*

Seq	Len	DT	OPT HL7/ syngo Dynamics	Element Name	DB Field in syngo Dynamics/ Comments
1	4	SI	R	Set ID	PR1 segment sequence number
2	3	IS	B	Procedure Coding Method	NA
3	250	CE	R	Procedure Code	From syngo Dynamics Procedure library
4	40	ST	B	Procedure Description	NA
5	26	TS	R	Procedure Date/Time	Report verification date/time
6	2	IS	O/N	Procedure Functional Type	NA
7	4	NM	O/N	Procedure Minutes	NA
8	250	XCN	B/N	Anesthesiologist	NA
9	2	IS	O/N	Anesthesia Code	NA
10	4	NM	O/N	Anesthesia Minutes	NA
11	250	XCN	B/N	Surgeon	NA
12	250	XCN	B	Procedure Practitioner	NA
13	250	CE	O/N	Consent Code	NA
14	2	ID	O/N	Procedure Priority	NA
15	250	CE	O	Associated Diagnosis Code	NA
16	250	CE	O	Procedure Code Modifier	NA
17	20	IS	O/N	Procedure DRG Type	NA
18	250	CE	O/N	Tissue Type Code	NA
19	427	EI	C	Procedure Identifier	NA

Seq	Len	DT	OPT HL7/ <i>syngo</i> Dynamics	Element Name	DB Field in <i>syngo</i> Dynamics/ Comments
20	1	ID	C	Procedure Action Code	NA

* This segment is not enabled by default.

This page is left blank intentionally.

12 syngo Dynamics Reporting

12.1 Allergy Information

Allergy information is communicated in the optional Allergy segment, AL1.

The AL1 segment contains information on a single allergy and is repeated for each allergy that is recorded for the patient. Within the AL1 segment, the allergy type and allergy description may be repeated, as required. While the AL1 segment can appear in ADT^A01 (Admit), ADT^A31 (Update Person), RDE (Pharmacy Encoded Order - Perfected), and ORM (Pharmacy Prescription Order - non-perfected), allergy information is imported only from ORM messages.

AL1 segment contains the allergy type and allergy description. Refer to "Message Segment: AL1" on page 69 for detailed information.

Allergy information from HL7 ADT and ORM messages, imported into study data, are displayed on the syngo Dynamics worksheet and report. Allergy information is imported into a study when the order is associated with the study. The allergy information is purged when the order association is broken.

The variables that contain the allergy information are as follows:

- Sensis: ALLER2 Group (ALLER2)
The allergy information is assigned to Sensis from the Study context and not from data contexts.
- Philips: PATN ALLERGIES (in 3rd Party Hemo) – C-priority, not for Phase Zero

12.2 Read-only PHI Fields

In the syngo Dynamics worksheet, Protected Health Information (PHI) data fields updated through HL7 orders in Electronic Health Record System (EHR) appear as read-only. This means that users cannot modify data for certain PHI fields that are present in syngo Dynamics worksheet and marked as EHR studies in the syngo Dynamics database.

Refer to the following PHI fields that are marked as read-only in the syngo Dynamics worksheet. If you do not want a certain field to appear read-only in the worksheet, remove the field from the EHRReadOnlyFields.txt file available on the syngo Dynamics server in <Server>\IF:\Acuson\Configuration.

- PatientID
- PatientsName

- PatientsBirthDate
- OtherPatientIDs
- AccessionNumber
- PatientSex
- PatientsSize
- PatientsWeight
- ReferringPhysiciansName
- MasterPatientID
- PatientsAddress1
- PatientsAddress2
- PatientsCity
- PatientsCountry
- PatientsOfficePhoneNumber
- PatientsEmailAddress
- PatientsHomePhoneNumber
- PatientsAddressState
- PatientsZipCode
- PatientsMobilePhoneNumber
- NameIDConflictExists - (The **Resolve Patient/ID Conflict** button will be in disabled state.)
- CurrentPatientLocation
- MultiMRNSiteName

12.3 Order Messages (ORM) Mapping

When *syngo* Dynamics is integrated with an Electronic Health Record System (EHR), the following information is extracted from the inbound HL7 ORM messages: Patient demographics, Visit details, Allergy details. The details are mapped in the *syngo* Dynamics database and the values pertaining to these are updated in the *syngo* Dynamics worksheet and report.

This information is required to complete pre-procedural documentation and reporting. For more information on pre-procedural information, refer to "Pre-procedural Documentation" on page 29.

Prerequisites to ensure that the mapping works are given as follows:

- HL7 interface (ADT/ORM) is enabled in the **Server Configuration Tool**.
- ORM matching rule (Filler Order Number vs Study Accession Number) is configured in SysAdmin.
- Site ID & Patient ID is configured in the ADT matching rules in SysAdmin.
- Fields mentioned in the table below are configured in the report and worksheet template.

The following table contains details about the HL7 field information and the corresponding *syngo* Dynamics database field name and the *syngo* Dynamics table containing the field. These HL7 fields are predefined in **HL7SegmentMapping.xml** present on the server at **F:\Acuson\Configuration\HL7Cache**.

HL7 Field	<i>syngo</i> Dynamics Database Field (Reporting Model)	Context Name	Report/Worksheet
PID 3.1 ID	PATIENT_ID	STUDY	Pt. ID
PID -3.4 Assigning Authority	SiteID/MultiMRNSiteID	STUDY	Site ID
PID -5 Patient Name	PATIENT_NAME	STUDY	Pt. Name
PID -7 Date of Birth	PATIENT_BIRTH_DATE	STUDY	Pt. Birth Date
PID-8 Sex	PATIENT_SEX	STUDY	Pt. Sex
PID 18.1 Patient Account Number	Department_Visit_Number	Study	AdmissionId
ORC-1 Order Control	Used For Processing*	NA	NA
ORC-2 Placer Order Number	PlacerOrderNumberImagin gServiceRequestLO	STUDY	Placer Order Number

HL7 Field	syngo Dynamics Database Field (Reporting Model)	Context Name	Report/Worksheet
ORC-3 Filler Order Number	ACCESSION_NUMBER	STUDY	Accession Number
ORC-5 Order Status	Used For Processing*	NA	NA
ORC-12 Ordering Provider	OrderingProviderName	STUDY	Ordering Provider Name
OBR-2 Placer Order Number	PlacerOrderNumberImagingServiceRequestLO	STUDY	Placer Order Number
OBR-3 Filler Order Number	ACCESSION_NUMBER	STUDY	Accession Number
OBR-4 Universal Service Identifier	UniversalServiceId	STUDY	Universal Service Identifier
OBR-5 Priority	Used For Processing*	NA	NA
OBR-7 Observation Date/Time	Study Date/Time(for pre-procedure)	STUDY	Study Date / Study Time
OBR-16 Ordering Provider	OrderingProviderName	STUDY	Ordering Provider Name
OBR-20 Filler Field 1	ACCESSION_NUMBER, Used for Processing*(HNA_ACCN)	STUDY, STUDY	Accession Number
OBR-31.2 Reason For Study	Indications	STUDY	Indications
PV1-1.2 Patient Class	PatientClass	STUDY	Patient Class
PV1-1.8 Referring Doctor	REFERRING_PHYSICIAN	STUDY	Referring Physician
MRG-1.1 ID	Used For Processing*	NA	NA
MRG-1.4 Assigning Authority	Used For Processing*	NA	NA
AL1-2 Allergy Type	Allergy_Type_obs	Allergies_ctx	Allergy Type
AL1-3 Allergy Code/Mnemonic/Description Can be customized to AL1-3.2	Allergy_Description_txt	Allergies_ctx	Allergy Description
AL1-4 Allergy Severity	Allergy_Severity_obs	Allergies_ctx	Allergy Severity
AL1-5 Allergy Reaction	Allergy_Reaction_txt	Allergies_ctx	Allergy Reaction

HL7 Field	syngo Dynamics Database Field (Reporting Model)	Context Name	Report/Worksheet
AL1-6 Identification Date	Allergy_Identification_Date_dt	Allergies_ctx	Allergy Identification Date

* Values are not stored in the database. The values are used for processing the inbound HL7 message.

System administrators can also configure additional HL7 fields required for processing HL7 ORM messages. For more information on how to configure additional HL7 fields, refer to the *syngo Dynamics System Administration Manual*.

This page is left blank intentionally.

13 Message Segments

The list specifies the information at each field level. The Fields are categorized into Field, Field Name, Type and Notes. In the Field Type “M” means Mandatory & “O” means optional. All of the “Not Supported Fields” and “save and forward option available” in the below tables are intended to be supported in future versions and currently they are supported as save and forward option only. The save and forward option in the ORM will store all mandatory fields for MWL providers. Currently, the EVN Segment is optional.



All of the “Non Supported Fields” can be configured to be supported by configuring them as additional HL7 fields. For more information on configuring additional HL7 fields, refer to syngo Dynamics System Administration Manual specific to the application version installed at your hospital.

13.1 Message Segment: MSH

Field	Field Name	Field Type	Notes
MSH 1	Separator	NA	
MSH 2	Encoding Characters	NA	^~\& Supported Currently
MSH 3	Sending Application	O	NA
MSH 4	Sending Facility	M	The MSH4 is Sending Facility (Configuration) becomes MSH6 Receiving Facility as per Inbound. For outbound ORM it will be vice versa.
MSH 5	Receiving Application	O	Configurable
MSH 6	Receiving Facility	O	NA
MSH 7	Date and Time of Message	O	NA
MSH 8	Security	O	NA
MSH 9.1	Message Type	M	E.g., ADT A01
MSH 9.2	Trigger Event		

Field	Field Name	Field Type	Notes
MSH 10	Message Control ID	O	Order ID is used as Message control ID in Outbound ORMs
MSH 11	Processing ID	O	No Validation
MSH 12	Version ID	M	2.3.1 version only supported. Any other formats will have to convert to the required format via standard interface engines.

13.2 Message Segment: EVN

Field	Field Name	Field Type	Notes
EVN 1	Event type Code	O	NA
EVN 2	Recorded Date and Time	O	Date and Time of the Event

13.3 Message Segment: PID

Field	Field Name	Field Type	Notes
PID 2	Patient ID	O	Save and Forward Only.
PID 3.1	Patient Identifier	M	Patient ID –Unique ID assigned by the HIS. If multiple PIDS (MRNs) are supported, they are separated by commas. Values Currently Accepted are MRN and MPI.(EMPI or EE SS Needs to mapped to either MR or MPI Only)
PID 3.4	Assigning Authority	M	PID assigner if Type MR. For MPI it can be optional or blank. Note: - If this value is blank for MR TYPE then NACKS is sent to IE.
PID 3.5	Identifier Type	M	MR or MPI .The type needs to be configured under Tab Patient ID mapping in data porting tool Siemens.Cvis.Utilities.DataPorting which will be located in the server path :- C:\Program Files\Siemens\Cvis\bin
PID 4.1	Alternate Patient ID	O	Save and Forward Option available

Field	Field Name	Field Type	Notes
PID 5.1	Family +Last Name	O	Anyone Name (PID 5.1 or 5.2) Should be available else the message is Error Logged. For PV1-2(Patient Class:E) The first name or last name is not mandatory
PID 5.2	Given Name	O	Anyone Name(PID 5.1 or 5.2) Should be available else the message is error logged
PID 5.3	Middle Name	O	Save and Forward Option available
PID 5.4	Suffix	O	Save and Forward Option available
PID 5.5	Prefix	O	Save and Forward Option available
PID 6	Mother's Maiden Name	O	Save and Forward Option available
PID 7.1	Date of Birth	O	Not mandatory field. If PID 7.1 is null, default value 1900/01/01) will be the DOB of the patient (YYMMDD)
PID 8	Gender	O	If no value is sent, the gender is updated as unknown. Supported are M,F,U and B .B is not recognized as standard HL7 defined Gender and it will updated as "Unknown "or "ambiguous" in UI
PID 9.1	Alias Last Name	O	UI will display the entire alias name in a single field separated by Comma.(ln,mn,fn)
PID 9.2	Alias Given Name	O	UI will display the entire alias name in a single field separated by Comma.
PID 9.3	Alias Middle Name	O	UI will display the entire alias name in a single field separated by Comma.
PID 10.1	Race	O	All Supported DDI Should be Part of the Configuration .Users is Advised to create the Same Abbreviation as used in the IE. Also there is an option to edit the Factory Default Abbreviation
PID 11.1 to 6	Patient Address	O	Street Address, Other Designation, City, State, ZIP and Country are Supported Each Value needs to separated by ^ and "&" needs to be replace by IT\Escape)
PID 13.1	Home Tel No	O	Save and Forward Option Available
PID 14.1	Business Tel no	O	Save and Forward Option Available

Field	Field Name	Field Type	Notes
PID 13.5,6,7	Telephone No -HOME	O	Supports up to 20 Char. Checks first this field and then 13.1
PID 14.5,6,7	Telephone No -Office	O	Supports up to 20 Char Checks first this field and then 14.1
PID 15.1	Language	O	All Supported DDI Should be Part of the Configuration .Users are advised to create the Same Abbreviation as used in the IE. Additionally users are allowed to edit DDI abbreviation.
PID 16.1	Marital status	O	All Supported DDI Should be Part of the Configuration .Users are Advised to create the Same Abbreviation as used in the IE Also there is an option to edit the Factory Default Abbreviation
PID 17.1	Religion	O	All Supported DDI Should be Part of the Configuration .Users are Advised to create the Same Abbreviation as used in the IE Also there is an option to edit the Factory Default Abbreviation
PID 18.1	Account Number	O	Unique Account Number
PID 19	SSN Number	O	Alternate ID Supported
PID 20	Driving License	O	Alternate ID Supported
PID 21	Mothers Identifier	O	Supported as Save and Forward Option.
PID 22	Ethnic Group	O	All Supported DDI Should be Part of the Configuration .Users are Advised to create the Same Abbreviation as used in the IE. Also there is an option to edit the Factory Default Abbreviation
PID 23	Birth Place	O	Save and Forward Option available
PID 24	Multiple Birth Indicator	O	Save and Forward Option available
PID 25	Birth Order	O	Save and Forward Option available
PID 26	Citizenship	O	Save and Forward Option available
PID 27	Veterans Military Status	O	Save and Forward Option available

Field	Field Name	Field Type	Notes
PID 28	Nationality	O	Save and Forward Option available
PID 29	Patient Death Date and Time	O	Save and Forward Option available
PID 30	Patient Death Indicator	O	Save and Forward Option available

13.4 Message Segment: PV1

Field	Field Name	Field Type	Notes
PV1-2	Patient Class	M	Only I, O, P and E Classes are Supported. Other Classes are also honored and treated as a separate class "Other". All outbound ORMs will be sent with the same class as received.
PV1-3.1,3.2,3.3	Assigned Patient Location	O	POC, Room and Bed Format Supported 3.1 and 3.2 are both mandatory fields required to process anyvalue in PV1-3 segment
PV1-4	Admission Type	O	Expects the Patient Class only and not the type .Save and Forward Option Available
PV1 -5	Pre-admit No	O	Save and Forward Option Available
PV1-6.1,6.2,6.3	Prior Patient Location	O	POC, Room and Bed Supported –Will be replaced by PV1-3.1 for A02 message
PV1-7	Attending Doctor	O	Save and forward option available
PV1-8	Referring Doctor	O	Tool available to import physicians. Performing and admitting doctors should be part of referring doctors import list.
PV1 -9	Consulting Doctor	O	Save and Forward Option Available
PV1-10	Hospital Service	O	Departments Like Cardiology, EP, Cath.
PV1-14	Admit Source	O	Save and Forward Option Available
PV1-15.1	Ambulatory Status	O	All Supported DDI Should be Part of the Configuration .Users are Advised to create the Same Abbreviation as used in the IE
PV1-16	VIP Indicator	O	Y (1),N (0) Supported

Field	Field Name	Field Type	Notes
PV1-17.1,17.2	Admitting Doctor	O	Supported
PV1-18	Patient Type	O	Save and Forward Option Available
PV1-19	Visit Number	M	<p>Unique Visit ID as received by the Placer Application (HIS).</p> <p>19.1 Visit Number retained and the same is sent as outbound.</p> <p>For 19.2 –Default value is 1</p> <p>19.3 –AA Should be same as PID AA or Blank if different the procedure cannot scheduled in WB.</p>
PV1-19.4	Visit Assigning Authority	M	Value is expected to be similar to PID AA and if different, needs to be negotiated. Save and Forward Option Available.
PV1-20.1	Financial Class	O	Save and Forward Option Available
PV1 -44	Admit Date	M	Date of Admission of the Patient Supported. If the Value is absent then Error is Logged
PV1-45	Discharge Date	O	Mandatory only for A03 (Discharge) Message ,else will log error and send positive Ack.

13.5 Message Segment: ORC

Field	Field Name	Field Type	Notes
ORC 1	Order Control	M	<p>Only RO, NW, CA, RP and DC, SC Supported. Any other Order control needs to converted or mapped to the Above types by the interface engine.</p> <p>RO modifies the order. CA, RP, and DC remove the order and are treated like a "cancel" in syngo Dynamics.</p>

Field	Field Name	Field Type	Notes
ORC 2	Placer Order Number (PON)	M	The Placer Order Number is referenced by a unique Filling Order Number (FON) Which is directly mapped to the Study UID generated by syngo Dynamics. The Same Filler ID is used as PON for Outbound ORMS. Mandatory only for Inbound.
ORC 3	Filler Order ID	M	For Outbound Only. Retain the FON if send by the placer and send as it is in outbound ORM for IMC and Deprecated for GA. If no FON is sent, new FON number is updated.
ORC 5	Order Status	O	Save and Forward Option Available
ORC 7	Quantity and Timing		Day only is supported in ORC 7.4 and in Outbound ORM the Exact day and time is sent. ORC 7.6 priority Supported as save and forward option available
ORC 12	Order Provider Name	O	Check OBR 16 first and if no value then the Value in ORC 12.
ORC 15	Order Effective Date/Time	O	Contains the date/time that the changes to the request took effect or are supposed to take effect.
OBR -1	Observation Request ID	O	Save and Forward Option Available
OBR-2.1	Placer Order Number	M	Unique Order ID generated by HIS (Placer Application). OBR2 or ORC2 either one is mandatory, not both fields.OBR2 is checked first and then ORC 2.If the Value is absent then error will be logged and send Pos ACK.
OBR -3.1	Filler Order Number	O	Retain id (accession number) sent by HIS/ RIS/inbound. Same value will be sent back IN ORM Outbound. (16 Character Restriction)

Field	Field Name	Field Type	Notes
OBR -4 4.1 – Identifier 4.2 –Text 4.3 – Name of the Coding System	Universal Service ID	M (Only 4.1)	Procedure types are sent here. Only the procedure types which are mapped to factory defaults via procedure mapping tool will be processed.b If not Pos Ack is sent and error is logged. Any new Procedure created needs to be uploaded via procedure mapping too 4.2 = HIS Procedure Code
OBR 5	Priority	O	If no priority is sent then the value is set to elective. No NACKS will be sent
OBR 6	Requested Date and time	O	If no date and time is sent then the value is set to current day and time No NACKS will be sent. Request date should be within next 12 months.
OBR 7	Observation Date and Time	M	If requested time is null, the observation time is taken. No NACKS will be sent
OBR 16	Order Provider Name	O	Check OBR 16 first and if no value then the Value in ORC 12.
OBR 18	Placer Field 1	O	Custom Field-Not Supported
OBR 19	Placer field 2	O	Custom Field-Not Supported
OBR 27.4	Start Date and Time	O	Not Supported. Available only for outbound ORM
OBR 27.6	Priority	O	Save and Forward Option Available
OBR 44	Procedure Code	O	Save and Forward Option Available
OBR 45	Procedure Code Modifier	O	Save and Forward Option Available

13.6 Message Segment: MRG

Field	Field Name	Field Type	Notes
MRG 1.1	Source Patient ID	M	Destination Patient Info is same as PID segment
MRG 1.4	Assigning Authority	M	PID assigning Authority. Mandatory only if the MRG1.5 value is 'MR'.
MRG 1.5	Identifier Code	M	MR or MPI
MRG 1.6	Assigning Facility	O	NA
MRG 3.1	Prior Patient Account Number	M	Only FOR A41
MRG 5.1	Prior Patient Visit Number	M	Only for A42

13.7 Message Segment: DG1

Field	Field Name	Field Type	Notes
DG 1.3.2	Diagnosis Code Text	O	Only the Description or Text is supported and other attributes of diagnosis are ignored.
DG 1.4	Diagnosis Description	O	Supported

13.8 Message Segment: AL1

Field	Field Name	Field Type	Notes
AL1 3	Allergy Code ,Mnemonic or Description	O	Currently only Text is supported and other attributes are ignored

One AL1 segment is sent for each separate patient allergy. Therefore a series of (none, 1 or more) AL1 segment(s) may be included in ADT messages, or in pharmacy order (RDE) messages. Refer to the table below for detailed information on the contents of the AL1 segment.

Sequence	Length	Format	Field Type	Element Name
0	3		R	Segment ID = "AL1"
1	4	SI	O	Set ID - Internal
2	2	IS	O	Allergy type The values and description are as follows: DA - Drug Allergy FA - Food Allergy MA - Miscellaneous Allergy
3	60	CE	R	Allergy description This field consists of several components as follows: <Allergy identifier>^<Text>^<Coding system>
4	2	IS	O	Allergy severity
5	15	ST	O	Allergy reaction
6	8	DT	O	Identification date

13.9 Message Segment: NK1

Field	Field Name	Field Type	Notes
NK 1.1	ID	O	Supported
NK 1.2	Name	O	Supported
NK 1.3	Relationship	O	Supported
NK 1.4	Address	O	Supported
NK 1.5	Phone No	O	Supported
NK 1.6	Business Phone No	O	Supported

13.10 Message Segment: ZDS (Customized Segment)

Field	Field Name	Field Type	Notes
ZDS 1	Study ID	O	Only for ORM outbound if it is Scheduled

If the Mandatory Fields are absent in the message then the system sends NACKS or Positive ACKs based on the sanity check. It is advised to configure the interface engine appropriately to receive or to ignore NACKS.

13.11 Message Segment: ZSU

Field	Field Name	Field Type	Notes
ZSU 1	Study instance ID (old)	R	ORM message containing the old study UID for the study that warrants a study UID correction
ZSU 1	Study instance ID (new)	R	ORM message containing the new study UID for the study that warrants a study UID correction

This page is left blank intentionally.

14 Appendix A: Sample Messages/Generic HL7 Information

14.1 ORM-O01 Message

```
MSH|^~\&|7EDIT|HL7EDIT|SDIS|SIEMENS|20110826160426||ORM^O01|MSGI
D200205122028|P|2.3.1
PID|||PS4234241^^^CVIS^MR||Donald^Duck^^^Mr.||20110826160411|M|||
123 West St.^MICHIGAN^CO^80020^USA||2498652892|2304972374|||
PV1||O|OP^20A^|||2342^Jones^Bob|||CAR|||||||239668^^^CVIS|||||||
|||||20110826160426|
ORC|NW|2010825122028
OBR|1|20110825160426||003038^AdultEcho^L|Elective|20110826160358|
20110826160404|||||||
```

14.2 ADT-A01 Message

```
MSH|^~\&|7edit|1|||20110826110607||ADT^A01|599102|P|2.3.1|||ER|AL|20
110825095518||||
EVN|A01|200708181123|
PID|1||10016604^^^CVIS^MR^1||Dido^Bear^D||19391010|M||1|111 Helm
ST^^NoMansland^CA^54321^^P|1|88853451212|8885551212|1|2||40477
16^^^7edit^VN^1|123561234|||||||N
PV1|1||PREOP^101^1^1^^^S|3|||37^Mickey Mouse
(AT)^^^^^^7edit^^^^CI|33337^ Mickey Mouse (R)||CAR|||1|||37^ Micky
Mouse
(ADM))^^^^^^7edit^^^^CI|2|41047334^^^CVIS^VN|4|||||||1||G||20
110826115152|20110826115156||||
DG1|1||9|71596^CMP ^I9|CMP ||A|
```

14.3 ADT–A40 Message (Merge Identifier)

```
MSH|^~\&||Test lab^00D0000000^GUID|
ELR|DOH|20061101102700||ADT^A40|2006280000002|T|2.3.1|
EVN|A40|20061101102700|
PID|||PID0700^5^M11^ADT1^MR^GOOD HEALTH
HOSPITAL|Road^^City^State^D|^M||^123^5551212|^123^5551
213||||345678901|||U|PV1||N||U||20090223150551|5101^KELL^DON
^P^DR|A|PUL||A|R|4|B4~A8|A|A|V105-
1^^^ADT1|A|A|A|A|A|A|A|A|A|A|04|A|A|A|H|A||200008161300||||TAX|V|
A|
MRG|PID0701^5^M11^ADT1^MR^GOOD HEALTH HOSPITAL
```

14.4 Message with Z Segment

The ZDS segment is used to send the study UID in all of the outbound ORMs. The study UID is required for the Modality Work List.

Inbound Z segments are not supported currently. If any other miscellaneous info needs to be supported, then a mutual understanding is needed.

```
MSH|^~\&|sDIS|sDIS||20110906143321||ORM^O01^ORM_O01|ORD000000
0003691|P|2.3.1
PID||PAT0000000003379^^^MR||Donald^Duck||190101010012|U
PV1||E||||||||||||VIS0000000003390
ORC|RP||sDIS_00000003720^^sDIS_00000003720||SC||^201109062130
^^Emergency||||^^administra||201109062130
OBR||sDIS_00000003720^^sDIS_00000003720|Biopsy US^Biopsy
US|Emergency||201109062130|||||||||US|||||^201109062130^^Emerg
ency
ZDS|1.3.12.2.1107.5.8.11.104110162960070568.634509161808692986
```

14.5 ZSU Segment in ORM Message

When a study UID correction is warranted, the ZSU segment contains the old and the new study UID.

```
MSH|^~\&|^ASYNQ_QUEUEING_SERVICE||||20170728133912.132-
0500||ORM^O01^ZSU_ORM_O01|1|P|S|T|2.6
```

```

PID||347823^^^&10.30.20.50.45.23|||APACHETEST^CRAWLER^MIDDLE||198
20508000000-0500|F
ZSU|2.16.840.1.113669.632.2.2626886709563223956.756970231717774
8753|1.2.840.113970.1.2.840.113970.2280106.20170727.1124721ORC|O
P|||||||0170728125959.19-0400
OBR|||00000CV170000169|ECHO TEE 3RD PARTY

```

14.6 ACK Messages

ACK –ACK messages are sent when the message is successfully received. An ACK does not guarantee that a message will be fully processed (e.g. a duplicate ADT^A01 may return an ACK but may not be fully processed).

```

MSH|^~\&|syngo Dynamics server|||||ACK|T|2.3.1|||||||
MSA|AA|20180322132200|||||

```

NACK message is also displayed in ACK message format.

14.7 ACK Scenarios Currently Handled

NACK message is also displayed in ACK message format.

Message	Segment and Field	Message Issue/ Scenario	Expected Behavior	Reason for NACK/ Additional Information
General	NA	HL7 Version Number Mismatch	NACK	Sanity Check
General	NA	HL7 Format is not proper	NACK	Sanity Check
All	PID 3.5	PID 3 segment does not contain MR or MPI number	NACK	Cannot identify/create the patient
All	PID 3	Patient Identifier(PID 3) is missing	NACK	Cannot identify/create the patient
ADT	PV1.19	Visit Number is missing	NACK	Cannot identify/create the Visit uniquely

Message	Segment and Field	Message Issue/ Scenario	Expected Behavior	Reason for NACK/ Additional Information
ADT	PV1.44	Date of Admission is missing	NACK	Not able to create the Visit without the Admission Date
ADT	PV1.2	Visit type is not mentioned or not handled	NACK	I, E, O, P is handled
ADT	PID 5	Both Patient Last name and First name is empty	NACK (Only for Non-Emergency Patient)	For Non emergency Patient, at least one of the names need to be present there. For emergency patient both names are optional
ADT	PID 7	Date of Birth Cannot be in future	NACK	NA
ORM (All)	OBR2	Order Placer Number is missing	NACK	Cannot create/identify order without the ID
ORM (All)	OBR 4.2	Requested Procedure is not supported	ACK	Unknown procedure s cannot be scheduled in WB. If the Procedure is not available in mapping table then positive ACK is sent and the error is logged in the system.
ORM (All)	ORC 1	Requested Order Control is Not supported	NACK	Order Control is not RP, NW, cancel, DL. No business logic for handling other than RP, NW, cancel, DL
ORM (All)	OBR 4.2	Procedure Type is Missing	NACK	Cannot create Procedure without procedure Type
Merge - All	PID 3	Patient ID type is neither MPI nor the MRN	NACK	Not supporting other ID types
Merge - A40	MRG1	Patient ID type is neither MPI nor the MRN	NACK	Not supporting other ID types

Message	Segment and Field	Message Issue/ Scenario	Expected Behavior	Reason for NACK/ Additional Information
Merge - A41	PID 18	Target patient Account Number is missing	NACK	Not able to identify the account to do the account merge
Merge - A41	MRG 3	Source Patient Account Number is missing	NACK	Not able to identify the account to do the account merge
Merge - A42	PV1 19	Target patient Visit Number is missing	NACK	Not able to identify the visit to do the visit merge
Merge - A42	MRG 5	Source Patient Visit Number is missing	NACK	Not able to identify the visit to do the visit merge
ORM	Enable /Disable	Option to Configure Outbound ORM	Y/N	
ORM (RP, Cancel)	OBR 2	Order with the given Placer Number is not found in the system (RP)	ACK, Do not Process; log error	Send ACK & log. Do not process
ORM (All)	OBR 6	Orders can be placed only for 12 months in future	ACK	Send ACK; Do not process; Log error if Order date is beyond 12 months.
ADT	PV1.19.4	Assigning Authority is not provided for the Visit	ACK	AA is expected when type is MRN. SD integration requires AA for creating Study (Under discussion)
ADT	PID 7.1	DOB missing	ACK	Accept & Process even if the DOB is missing
A03	PV1 19	Visit ID missing	ACK	Accept and Process and log error

Message	Segment and Field	Message Issue/ Scenario	Expected Behavior	Reason for NACK/ Additional Information
ORM –O01	NA	Repeat ORMs	ACK	Process and Log Error

15 Appendix B: Configuration Files for ADT/ORM Option

HL7CacheConfiguration.xml and HL7SegmentMapping.xml files are only used for the ADT or ORM reconciliation options. ADT 31 is supported now. These two xml configuration files can be located on the syngo Dynamics server at F:\Acuson\Configuration\HL7Cache

The first file HL7CacheConfiguration.xml contains the list of ADT messages we will listen for. To prevent syngo Dynamics from receiving one of these message types you must remove its entry from this file and restart the ADT Cache Manager service.

```
<HL7CacheConfig>
<Version>4.0</Version>
<Enable>True</Enable>
<EnterpriseAssigningAuthority>MPI</EnterpriseAssigningAuthority>
<PIDTypeBlackList>SS</PIDTypeBlackList>
<ADTCacheConfig>
<Enable>True</Enable>
<TCPPort>8000</TCPPort>
<PurgeTimerInterval>86400000</PurgeTimerInterval>
<CacheDays>3</CacheDays>
<SupportedADTMessages>
<ADTMessage>ADT^A01</ADTMessage>
<ADTMessage>ADT^A02</ADTMessage>
<ADTMessage>ADT^A04</ADTMessage>
<ADTMessage>ADT^A06</ADTMessage>
<ADTMessage>ADT^A07</ADTMessage>
<ADTMessage>ADT^A08</ADTMessage>
<ADTMessage>ADT^A18</ADTMessage>
<ADTMessage>ADT^A34</ADTMessage>
<ADTMessage>ADT^A36</ADTMessage>
<ADTMessage>ADT^A40</ADTMessage>
```

```

<ADTMessage>ADT^A47</ADTMessage>
</SupportedADTMessages>
<ADTReconciliationRules>
  <ADTReconciliationRule Matches="Site ID,Patient ID,First Name,Last
Name,Date of Birth" Enabled="true" />
  <ADTReconciliationRule Matches="Site ID,Patient ID,First Initial,Last Name"
Enabled="true" />
  <ADTReconciliationRule Matches="Site ID,Patient ID,Date of Birth"
Enabled="true" />
  <ADTReconciliationRule Matches="Site ID,Patient ID" Enabled="true" />
  <ADTReconciliationRule Matches="Patient ID,First Name,Last Name,Date of
Birth" Enabled="true" />
  <ADTReconciliationRule Matches="First Initial,Last Name,Date of Birth,"
Enabled="true" />
  <ADTReconciliationRule Matches="First Initial,Last Name" Enabled="true" />
  <ADTReconciliationRule Matches="Master Patient Index must exist in
successful match" Enabled="true" />
</ADTReconciliationRules>
</ADTCacheConfig>
<ORMCacheConfig>
</HL7CacheConfig>

```

The second file is the HL7SegmentMapping.xml file which contains the mapping info for each item we accept. This file is not intended to be modified. This information is presented for your reference only.

```

- <HL7MessageMapping xmlns:xs="http://www.w3.org/2001/XMLSchema"
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance" Version="1.0"
xsi:schemaLocation="syngoDynamics
Z:\Server\Implementation\HL7MessageParser\Schemas\MessageDefinition.xs
d" xmlns="syngoDynamics">
  <Platform />
  <PlatformVersions />
  <ObservationIDFieldPath>3</ObservationIDFieldPath>
- <SegmentDefinitions>
- <NonOBXSegment>
  <SegmentName>PID</SegmentName>

```



```

- <MappedFields>
  <ValueFieldPath>3</ValueFieldPath>
  <sDxName>PatientIDList</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>MRN_PatientIDList</Description>
</MappedFields>
- <MappedFields>
  <ValueFieldPath>3.1</ValueFieldPath>
  <sDxName>PatientID</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>MRN_PatientID</Description>
</MappedFields>
- <MappedFields>
  <ValueFieldPath>3.4</ValueFieldPath>
  <sDxName>SiteID</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>MRN_SiteID</Description>
</MappedFields>
- <MappedFields>
  <ValueFieldPath>3.5</ValueFieldPath>
  <sDxName>PIDType</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>MRN_Type</Description>
</MappedFields>

```

```

- <MappedFields>
  <ValueFieldPath>5</ValueFieldPath>
  <sDxName>PatientsName</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>PersonName</DataType>
  <Description>Patient Name</Description>
</MappedFields>
- <MappedFields>
  <ValueFieldPath>7</ValueFieldPath>
  <sDxName>PatientsBirthDate</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Date</DataType>
  <Description>Patient Birth Date</Description>
</MappedFields>
- <MappedFields>
  <ValueFieldPath>8</ValueFieldPath>
  <sDxName>PatientsGender</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>Patient Gender</Description>
</MappedFields>
</NonOBXSegment>
- <NonOBXSegment>
  <SegmentName>PV1</SegmentName>
- <MappedFields>
  <ValueFieldPath>2</ValueFieldPath>
  <sDxName>PatientClass</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>

```

```

<DataType>Text</DataType>
<Description>Patient Class</Description>
</MappedFields>
</NonOBXSegment>
- <NonOBXSegment>
  <SegmentName>MRG</SegmentName>
- <MappedFields>
  <ValueFieldPath>1</ValueFieldPath>
  <sDxName>PatientIDList</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>MRN_PatientIDList</Description>
  </MappedFields>
- <MappedFields>
  <ValueFieldPath>1.1</ValueFieldPath>
  <sDxName>PatientID</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>MRN_PatientID</Description>
  </MappedFields>
- <MappedFields>
  <ValueFieldPath>1.4</ValueFieldPath>
  <sDxName>SiteID</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>MRN_SiteID</Description>
  </MappedFields>
- <MappedFields>
  <ValueFieldPath>1.5</ValueFieldPath>

```

```

<sDxName>PIDType</sDxName>
<IsTableValue>>false</IsTableValue>
<DefaultUnits>unitless</DefaultUnits>
<DataType>Text</DataType>
<Description>MRN_Type</Description>
</MappedFields>
</NonOBXSegment>
- <NonOBXSegment>
  <SegmentName>ORC</SegmentName>
- <MappedFields>
  <ValueFieldPath>1</ValueFieldPath>
  <sDxName>OrderControl</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>OrderControl</Description>
  </MappedFields>
- <MappedFields>
  <ValueFieldPath>2</ValueFieldPath>
  <sDxName>PlacerOrderNumber</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>PlacerOrderNumber</Description>
  </MappedFields>
- <MappedFields>
  <ValueFieldPath>3</ValueFieldPath>
  <sDxName>FillerOrderNumber</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>FillerOrderNumber</Description>

```

```

    </MappedFields>
- <MappedFields>
    <ValueFieldPath>5</ValueFieldPath>
    <sDxName>OrderStatus</sDxName>
    <IsTableValue>>false</IsTableValue>
    <DefaultUnits>unitless</DefaultUnits>
    <DataType>Text</DataType>
    <Description>OrderStatus</Description>
    </MappedFields>
- <MappedFields>
    <ValueFieldPath>8</ValueFieldPath>
    <sDxName>Parent</sDxName>
    <IsTableValue>>false</IsTableValue>
    <DefaultUnits>unitless</DefaultUnits>
    <DataType>Text</DataType>
    <Description>Parent</Description>
    </MappedFields>
- <MappedFields>
    <ValueFieldPath>12</ValueFieldPath>
    <sDxName>OrderingProvider</sDxName>
    <IsTableValue>>false</IsTableValue>
    <DefaultUnits>unitless</DefaultUnits>
    <DataType>Text</DataType>
    <Description>OrderingProvider</Description>
    </MappedFields>
    </NonOBXSegment>
- <NonOBXSegment>
    <SegmentName>OBR</SegmentName>
- <MappedFields>
    <ValueFieldPath>2</ValueFieldPath>
    <sDxName>PlacerOrderNumber</sDxName>
    <IsTableValue>>false</IsTableValue>

```


```

<DefaultUnits>unitless</DefaultUnits>
<DataType>Text</DataType>
<Description>PlacerOrderNumber</Description>
</MappedFields>
- <MappedFields>
  <ValueFieldPath>3</ValueFieldPath>
  <sDxName>FillerOrderNumber</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>FillerOrderNumber</Description>
  </MappedFields>
- <MappedFields>
  <ValueFieldPath>4</ValueFieldPath>
  <sDxName>UniversalServiceId</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Text</DataType>
  <Description>UniversalServiceId</Description>
  </MappedFields>
- <MappedFields>
  <ValueFieldPath>7</ValueFieldPath>
  <sDxName>ObservationDate</sDxName>
  <IsTableValue>>false</IsTableValue>
  <DefaultUnits>unitless</DefaultUnits>
  <DataType>Date</DataType>
  <Description>ObservationDate</Description>
  </MappedFields>
- <MappedFields>
  <ValueFieldPath>16</ValueFieldPath>
  <sDxName>OrderingProvider</sDxName>
  <IsTableValue>>false</IsTableValue>

```

```
<DefaultUnits>unitless</DefaultUnits>  
<DataType>Text</DataType>  
<Description>OrderingProvider</Description>  
</MappedFields>  
</NonOBXSegment>  
</SegmentDefinitions>  
</HL7MessageMapping>
```

This page is left blank intentionally.



This page is left blank intentionally.



Manufacturer's note:

All product designations and company names are trademarks or registered trademarks of the corresponding companies.

Siemens Healthineers reserves the right to modify the design and specifications contained herein without prior notice.

Some of the specifications described herein may not be currently available in all countries. Please contact your local Siemens Healthineers Sales representative for the most current information.

Caution: US federal law restricts the herein described devices to sale by or on the order of a physician.

The original language of this document is English.

Origin US

Local Contact Information in the USA

Siemens Medical Solutions USA, Inc.
40 Liberty Boulevard
Malvern, PA, 19355-9998
U.S.A
Telephone:
+1-610-448-4500
1-888-826-9702

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Manufacturer

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany