



Improving the efficacy of cerebrospinal fluid testing for the diagnosis of multiple sclerosis

by a hyperbolic reference range according to Reiber (Reibergram) using the N Latex FLC kappa Assay

Dr. Andreas Rechner, Siemens Healthcare Diagnostics Products GmbH
Prof. Dr. Thomas Skripuletz, Department of Neurology, Hannover Medical School

[siemens-healthineers.com](https://www.siemens-healthineers.com)



Multiple sclerosis (MS) is an autoimmune inflammatory neurological disease with an incidence of approximately 6560 new cases per year and a prevalence of 143,000 to 199,500 patients in Germany alone.¹ In 2013, 2.3 million people worldwide were estimated to have MS, with a much higher prevalence in North America and most European countries.² Early diagnosis and treatment can significantly improve management and slow the progression of this not-yet-curable disease.³ In the latest revision of the McDonald criteria for the diagnosis of multiple sclerosis,⁴ the qualitative detection of oligoclonal bands (OCBs) in cerebrospinal fluid (CSF) allows the diagnosis of MS in selected patients as a substitute for fulfilling the requirement of dissemination in time (DIT). The presence of OCBs in CSF is an indicator of an intrathecal humoral immune response. The consensus recommendation in the McDonald criteria defines the presence of CSF OCBs as evidence of DIT, which is characterized by the development of an inflammatory activity in the central nervous system at different times.

According to latest findings, the qualitative detection of OCBs in CSF has a diagnostic sensitivity of up to 98% in MS patients.^{5,6} However, it is a laborious method in the routine laboratory that requires skilled and experienced staff. The fully automated, nephelometric determination of kappa free light chain (κ FLC) levels in CSF has been shown to be a suitable screening test to detect an intrathecal humoral immune response. Siemens Healthineers was first to introduce a CE-marked assay application for FLC testing in CSF samples.

Several studies used the Siemens Healthineers N Latex FLC kappa Assay to assess the presence of MS, applying different approaches and algorithms.⁷⁻⁹ For example, Valencia-Vera et al. concluded: “ κ FLC determination is rapid and automatized, but it has no higher sensitivity and specificity than OCB in MS diagnosis. Nevertheless, when used in screening, it could reduce the number of manual OCB tests.”

All approaches aim to compensate for the large inter-individual variation in blood-brain barrier function as well as the serum levels of κ FLC. Reiber et al.¹⁰ now propose the use of a hyperbolic reference range in CSF/serum quotient diagrams to further improve and standardize result interpretation of κ FLC levels in CSF. For decades, similar diagrams have been used for the assessment of inflammatory processes in the central nervous system (CNS) and the detection of intrathecal synthesis of immunoglobulins. In this study, the N Latex FLC kappa Assay (site 1) and FREELITE assay (site 2) were used to derive a hyperbolic reference range from a total of 433 defined controls by fitting a hyperbolic function to the CSF/serum quotients of albumin plotted against the respective CSF/serum quotients of κ FLC (Figure 1, Reibergram). A method-independent hyperbolic reference range was then established by adding ± 3 times the mean coefficient of variation (CV).




Figure 1. κ FLC control data in a double logarithmic quotient diagram (Reibergram).

The hyperbolic functions are shown up to $QAlb = 150 \times 10^{-3}$. Q_{mean} (dashed line) is shown to characterize the optical asymmetry of Q_{mean} in the double log diagram. The intrathecal fraction, IF, is represented either with reference to Q_{lim} for diagnostic purposes as a relative fraction (KIF in %) or with reference to Q_{mean} for statistical purposes, KIF(mean), eventually calculated as quantitative value $Kloc(mean)$ in mg/L.

Applying the Reibergram for the diagnosis of MS appeared to have a comparable sensitivity to the reference method of OCBs in CSF. The evaluation of 95 MS patients (45 MS patients and 50 clinically isolated syndrome patients, later confirmed as MS patients) showed a sensitivity of 93% when using the suggested κ FLC analytic approach (Figure 2).

“Reiber’s κ FLC diagram shows a great diagnostic performance to detect an intrathecal κ FLC production in patients with MS.”




Figure 2. MS and CV patient data in the κ FLC Reibergram.

Filled circles are multiple sclerosis patients ($n = 45$), and open circles represent clinically isolated syndrome patients, later found as MS ($n = 50$). FLC-K analysis detects 93% of a total of 95 MS patients.


Another study used the N Latex FLC kappa Assay for the quantification of κ FLC in CSF and serum to investigate this new algorithm by Reiber, which was able to identify 98 of 100 confirmed MS patients, just one patient less (99/100) than with OCB analysis.¹¹ In comparison, Reiber's approach appeared to be superior to other approaches for assessing MS using the κ FLC results, as summarized in Table 1.

MS (McDonald 2017)			
Characteristics	MS, n = 100	CIS to MS, n = 24	Stable CIS, n = 44
Oligoclonal bands, n (%)	99/100 (99%)	21/24 (88%)	11/44 (25%)
Reiber's FLC diagram, n (%)	98/100 (98%)	21/24 (88%)	9/44 (20%)
Presslauer's FLC curve, n (%)	96/100 (96%)	19/24 (79%)	9/44 (20%)
Senel's FLC curve, n (%)	96/100 (96%)	19/24 (79%)	8/44 (18%)
FLC index >5.9, n (%)	96/100 (96%)	19/24 (79%)	7/44 (16%)
Intrathecal IgG synthesis, n (%)	59/100 (59%)	18/24 (75%)	2/44 (5%)
Intrathecal IgM synthesis, n (%)	33/100 (33%)	4/24 (17%)	2/44 (5%)
Intrathecal IgA synthesis, n (%)	11/100 (11%)	2/24 (8%)	0/44 (0%)
Age, median (range)	32 (16–73)	30.5 (15–73)	36 (16–53)
Females, n (%)	73/100 (73%)	15/24 (63%)	28/44 (64%)

Table 1. Demographic and laboratory characteristics of patients diagnosed with MS according to the McDonald criteria of 2017, patients clinically isolated syndrome (CIS) converted to MS during follow-up, and patients with stable CIS. Laboratory characteristics include the determination of oligoclonal bands, the kappa free light chains index, and the proportion of an intrathecal synthesis of FLC using the methods of Reiber,¹⁰ Presslauer,¹² and Senel.⁷

The authors of this direct-comparison study consequently concluded that "Reiber's κ FLC diagram shows a great diagnostic performance to detect an intrathecal κ FLC production in patients with MS."

Although additional studies need to be conducted to confirm these results, the fully automated and quantitative determination of κ FLC in CSF and serum could significantly reduce the number of time-consuming OCB analyses in CSF by providing evidence of dissemination in time for MS diagnostics according to the McDonald criteria of 2017. One possibility might be the following two-stage approach intended to improve the workflow in the laboratory:



*Perform OCB testing if patients are profoundly suspicious for MS.

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all made possible by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 120 years of experience and 18,000 patents globally. Through the dedication of more than 50,000 colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

This white paper is intended for use outside U.S. only (OUS). Display or promotion to U.S. audience is prohibited.

References:

1. Kip M, Zimmerman A, Bless HH. Epidemiologie der Multiplen Sklerose. In: Kip M, et al., Hrsg. Weißbuch Multiple Sklerose. 2016.
2. Browne P, Chandraratna D, Angood C, Tremlett H, Baker C, Taylor BV, Thompson AJ. Atlas of multiple sclerosis 2013: a growing global problem with widespread inequity. *Neurol.* 2014;83(11):1022-4.
3. Cerqueira JJ, Compston DAS, Ghera R, Rosa MM, Schmierer K, Thompson A, Tinelli M, Palace J. Time matters in multiple sclerosis: can early treatment and long-term follow-up ensure everyone benefits from the latest advances in multiple sclerosis? *J Neurol Neurosurg Psychiatry.* 2018;89:844-50.
4. Thompson AJ, Banwell BL, Barkhof F, Carroll WM, Coetzee T, Comi G, Correale J, Fazekas F, Filippi M, Freedman MS, Fujihara K, Galetta SL, Hartung HP, Kappos L, Lublin FD, Marrie RA, Miller AE, Miller DH, Montalban X, Mowry EM, Sorensen PS, Tintoré M, Traboulsee AL, Trojano M, Uitdehaag BMJ, Vukusic S, Waubant E, Weinshenker BG, Reingold SC, Cohen JA. Diagnosis of multiple sclerosis: 2017 revisions of the McDonald criteria. *Lancet Neurol.* 2018;17(2):162-73.
5. Reiber H, Teut M, Pohl D, Rostasy KM, Hanefeld F. Paediatric and adult multiple sclerosis: age related differences and time course of the neuroimmunological response in cerebrospinal fluid. *Mult Scler.* 2009;15:1466-80.
6. Schwenkenbecher P, Wurster U, Konen FF, Ginge S, Sühs KW, Wattjes MP, Stangel M, Skripuletz T. Impact of the McDonald criteria 2017 on early diagnosis of relapsing-remitting multiple sclerosis. *Front Neurol.* 2019;10:188.
7. Senel M, Mojib-Yezdani F, Braisch U, Bachhuber F, Lewerenz J, Ludolph AC, Otto M, Tumani H. CSF free light chains as a marker of intrathecal immunoglobulin synthesis in multiple sclerosis: a blood-CSF barrier related evaluation in a large cohort. *Front Immunol.* 2019;10:641.
8. Valencia-Vera E, Martinez-Escribano Garcia-Ripoll A, Enguix A, Abalos-Garcia C, Segovia-Cuevas MJ. Application of κ free light chains in cerebrospinal fluid as a biomarker in multiple sclerosis diagnosis: development of a diagnosis algorithm. *Clin Chem Lab Med.* 2018;56(4):609-13.
9. Vasilj M, Kes VB, Vrkic N, Vukasovic I. Relevance of κ FLC quantification to differentiate clinically isolated syndrome from multiple sclerosis at clinical onset. *Clin Neurol Neurosurg.* 2018;174:220-9.
10. Reiber H, Zeman D, Kušnierová P, Mundwiler E, Bernasconi L. Diagnostic relevance of free light chains in cerebrospinal fluid - the hyperbolic reference range for reliable data interpretation in quotient diagrams. *Clin Chim Acta.* 2019;497:153-62.
11. Schwenkenbecher P, Konen FF, Wurster U, Witte T, Ginge S, Sühs KW, Stangel M, Skripuletz T. Reiber's diagram for kappa free light chains: the new standard for assessing intrathecal synthesis? *Diagnostics.* 2019;9:194.
12. Presslauer S, Milosavljevic D, Huebl W, Parigger S, Schneider-Koch G, Bruecke T. Kappa free light chains: diagnostic and prognostic relevance in MS and CIS. *PLoS ONE.* 2014;9:e89945.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
Laboratory Diagnostics
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
Phone: +1 914-631-8000