# MAGNETOM Flash

The Magazine of MR 翻訳版

SIEMENS

## Research

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## References

- 1 Jackson A, Buckley DL, Parker GJ. Dynamic Contrast-Enhanced Magnetic Resonance Imaging in Oncology. Springer, 2004.
- 2 Tofts PS. Quantitative MRI of the brain: measuring changes caused by disease. John Wiley, 2003.
- 3 Tofts PS, Brix G, Buckley DL, Evelhoch JL, Henderson E, Knopp MV, Larsson HB, Lee TY, Mayr NA, Parker GJ, Port RE, Taylor J, Weisskoff RM. Estimating kinetic parameters from dynamic contrast-enhanced T(1)-weighted MRI of a diffusable tracer: standardized quantities and symbols. J Magn Reson Imaging 1999;10:223-232.
- 4 Leach MO, Brindle KM, Evelhoch JL, Griffiths JR, Horsman MR, Jackson A, Jayson GC, Judson IR, Knopp MV, Maxwell RJ, McIntyre D, Padhani AR, Price P, Rathbone R, Rustin GJ, Tofts PS, Tozer GM, Vennart W, Waterton JC, Williams SR, Workman P. The assessment of antiangiogenic and antivascular therapies in early-stage clinical trials using magnetic resonance imaging: issues and recommendations. Br J Cancer 2005;92:1599-1610.
- 5 Tofts PS, Kermode AG. Measurement of the blood-brain barrier permeability and leakage space using dynamic MR imaging. 1. Fundamental concepts. Magn Reson Med 1991;17:357-367.

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Quantity	symbol	units	type
flip angle <sup>a</sup>	FA	degrees	fixed
haematocrit	Hct	%	fixed (42%)
onset time	t <sub>onset</sub>	s	free
rate constant <sup>b</sup>	k <sub>ep</sub>	min <sup>-1</sup>	free
transfer constant	K <sup>trans</sup>	min <sup>-1</sup>	free
T <sub>1</sub> relaxivity	r <sub>1</sub>	s <sup>-1</sup> mM <sup>-1</sup>	fixed (4.5 s <sup>-1</sup> mM <sup>-1</sup> )
T <sub>1</sub> of blood	T <sub>10</sub> <sup>blood</sup>	s	fixed (1.4 s)
T <sub>1</sub> of tissue	T <sub>10</sub>	s	fixed
TR	TR	s	fixed
fractional volume of EES <sup>c</sup>	v <sub>e</sub>	0<v <sub>e</sub> <100%	free
fractional volume of blood plasma in tissue	v <sub>p</sub>	0<v <sub>p</sub> <100%	free

<sup>a</sup>θ is the flip angle in radians, <sup>b</sup>k<sub>ep</sub> = K<sup>trans</sup>/v<sub>e</sub>, <sup>c</sup>Extravascular Extracellular Space

- 6 Kety SS. The theory and applications of the exchange of inert gas at the lungs and tissues. *Pharmacol Rev* 1951;3:1-41.
- 7 Tofts PS. Modeling tracer kinetics in dynamic Gd-DTPA MR imaging. *J Magn Reson Imaging* 1997;7:91-101.
- 8 St Lawrence KS, Lee TY. An adiabatic approximation to the tissue homogeneity model for water exchange in the brain: I. Theoretical derivation. *J Cereb Blood Flow Metab* 1998;18:1365-1377.
- 9 Donaldson SB, West CM, Davidson SE, Carrington BM, Hutchison G, Jones AP, Sourbron SP, Buckley DL. A comparison of tracer kinetic models for T<sub>1</sub>-weighted dynamic contrast-enhanced MRI: application in carcinoma of the cervix. *Magn Reson Med* 2010;63:691-700.
- 10 Tofts PS, Cutajar M, Mendichovszky IA, Gordon I. Accurate and precise measurement of renal filtration and vascular parameters using DCE-MRI and a 3-compartment model. *Proc Intl Soc Mag Reson Med*, 18th annual meeting, Stockholm.2010; 326.
- 11 Parker GJ, Roberts C, Macdonald A, Buonaccorsi GA, Cheung S, Buckley DL, Jackson A, Watson Y, Davies K, Jayson GC. Experimentally-derived functional form for a population-averaged hightemporal-resolution arterial input function for dynamic contrast-enhanced MRI. *Magn Reson Med* 2006;56:993-1000.
- 12 Dowell NG, Tofts PS. Fast, accurate, and precise mapping of the RF field in vivo using the 180 degrees signal null. *Magn Reson Med* 2007;58:622-630.
- 13 Tofts PS, Berkowitz B, Schnall MD. Quantitative analysis of dynamic Gd-DTPA enhancement in breast tumors using a permeability model. *Magn Reson Med* 1995;33:564-568.
- 14 Fritz-Hansen T, Rostrup E, Larsson HB, Sondergaard L, Ring P, Henriksen O. Measurement of the arterial concentration of Gd-DTPA using MRI: a step toward quantitative perfusion imaging. *Magn Reson Med* 1996;36:225-231.
- 15 Tofts PS, Davies GR, Dehmshki J. Histograms: measuring subtle diffuse disease (chapter 18). In: Paul Tofts, editor. *Quantitative MRI of the brain: measuring changes caused by disease*. Chichester: John Wiley, 2003: 581-610.
- 16 Tofts PS, Steens SC, Cercignani M, Admiraal-Behloul F, Hofman PA, van Osch MJ, Teeuwisse WM, Tozer DJ, van Waesberghe JH, Yeung R, Barker GJ, van Buchem MA. Sources of variation in multi-centre brain MTR histogram studies: body-coil transmission eliminates inter-centre differences. *Magn Reson Mater Phy* 2006;19:209-222.
- 17 Tofts PS, Benton CE, Weil RS, Tozer DJ, Altmann DR, Jager HR, Waldman AD, Rees JH. Quantitative analysis of whole-tumor Gd enhancement histograms predicts malignant transformation in low-grade gliomas. *J Magn Reson Imaging* 2007;25:208-214.
- 18 Donaldson SB, Buckley DL, O'Connor JP, Davidson SE, Carrington BM, Jones AP, West CM. Enhancing fraction measured using dynamic contrast-enhanced MRI predicts disease-free survival in patients with carcinoma of the cervix. *Br J Cancer* 2010;102:23-26.
- 19 Dehmshki J, Ruto AC, Arridge S, Silver NC, Miller DH, Tofts PS. Analysis of MTR histograms in multiple sclerosis using principal components and multiple discriminant analysis. *Magn Reson Med* 2001;46:600-609.
- 20 Spees WM, Yablonskiy DA, Oswood MC, Ackerman JJ. Water proton MR properties of human blood at 1.5 Tesla: magnetic susceptibility, T(1), T(2), T\*(2), and non-Lorentzian signal behavior. *Magn Reson Med* 2001;45:533-542.

#### Contact

Professor Paul Tofts  
 Clinical Imaging Sciences Centre  
 Brighton and Sussex Medical School  
 University of Sussex  
 Brighton  
 BN1 9PX  
 United Kingdom  
 Website: [www.paul-tofts-phd.org.uk](http://www.paul-tofts-phd.org.uk)  
[dce@paul-tofts.org.uk](mailto:dce@paul-tofts.org.uk)

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#### References

- Ross R, Harisinghani M. Prostate cancer imaging – what the urologic oncologist needs to know. *Radiol Clin North Am* 2006; 44:711-722, viii.
- Somford DM, Futterer JJ, Hambrock T, Barentsz JO. Diffusion and perfusion MR imaging of the prostate. *Magn Reson Imaging Clin N Am* 2008; 16:685-695, ix.
- Jung JA, Coakley FV, Vigneron DB, et al. Prostate depiction at endorectal MR spectroscopic imaging: investigation of a standardized evaluation system. *Radiology* 2004; 233:701-708.
- Futterer JJ, Heijmink SW, Scheenen TW, et al. Prostate cancer localization with dynamic contrast-enhanced MR imaging and proton MR spectroscopic imaging. *Radiology* 2006; 241:449-458.
- Bloch BN, Furman-Haran E, Helbich TH, et al. Prostate cancer: accurate determination of extracapsular extension with high-spatial-resolution dynamic contrast-enhanced and T2-weighted MR imaging – initial results. *Radiology* 2007; 245:176-185.
- Weinreb JC, Blume JD, Coakley FV, et al. Prostate cancer: sextant localization at MR imaging and MR spectroscopic imaging before prostatectomy – results of ACRIN prospective multi-institutional clinicopathologic study. *Radiology* 2009; 251:122-133.
- Langer DL, van der Kwast TH, Evans AJ, Plotkin A, Trachtenberg J, Wilson BC, Haider MA. Prostate tissue composition and MR measurements: investigating the relationships between ADC, T2,  $K^{trans}$ ,  $v(e)$ , and corresponding histologic features. *Radiology*. 2010 May;255(2):485-94.
- Turkbey B, Pinto PA, Mani H, Bernardo M, Pang Y, McKinney YL, Khurana K, Ravizzini GC, Albert PS, Merino MJ, Choyke PL. Prostate cancer: value of multiparametric MR imaging at 3 T for detection – histopathologic correlation. *Radiology*. 2010 Apr;255(1):89-99.
- Groenendaal G, Moman MR, Korporaal JG, van Diest PJ, van Vulpen M, Philippens ME, van der Heide UA. Validation of functional imaging with pathology for tumor delineation in the prostate. *Radiother Oncol*. 2010 Feb;94(2):145-50.
- Langer DL, van der Kwast TH, Evans AJ, Trachtenberg J, Wilson BC, Haider MA. Prostate cancer detection with multi-parametric MRI: logistic regression analysis of quantitative T2, diffusion-weighted imaging, and dynamic contrast-enhanced MRI. *J Magn Reson Imaging*. 2009 Aug;30(2):327-34.
- Franiel T, Lüdemann L, Rudolph B, Rehbein H, Stephan C, Taupitz M, Beyersdorff D. Prostate MR imaging: tissue characterization with pharmacokinetic volume and blood flow parameters and correlation with histologic parameters. *Radiology*. 2009 Jul;252(1):101-8.
- Ren J, Huan Y, Wang H, Ge Y, Chang Y, Yin H, Sun L. Seminal vesicle invasion in prostate cancer: prediction with combined T2-weighted and diffusion-weighted MR imaging. *Eur Radiol*. 2009 Oct;19(10):2481-6.
- Futterer JJ, Barentsz JO, Heijmink SW. Value of 3-T magnetic resonance imaging in local staging of prostate cancer. *Top Magn Reson Imaging*. 2008 Dec;19(6):285-9.
- McMahon CJ, Bloch BN, Lenkinski RE, Rofsky NM. Dynamic contrast-enhanced MR imaging in the evaluation of patients with prostate cancer. *Magn Reson Imaging Clin N Am*. 2009 May;17(2):363-83.

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## References

- Houssami, N., et al., Accuracy and surgical impact of magnetic resonance imaging in breast cancer staging: systematic review and meta-analysis in detection of multifocal and multicentric cancer. *J Clin Oncol*, 2008. 26(19): p. 3248-58.
- Peters, N.H., et al., Meta-analysis of MR imaging in the diagnosis of breast lesions. *Radiology*, 2008. 246(1): p. 116-24.
- Warner, E., et al., Surveillance of BRCA1 and BRCA2 mutation carriers with magnetic resonance imaging, ultrasound, mammography, and clinical breast examination. *Jama*, 2004. 292(11): p. 1317-25.
- Kuhl, C.K., et al., MRI for diagnosis of pure ductal carcinoma in situ: a prospective observational study. *Lancet*, 2007. 370(9586): p. 485-92.
- Vag, T., et al., Diagnosis of ductal carcinoma in situ using contrast-enhanced magnetic resonance mammography compared with conventional mammography. *Clin Imaging*, 2008. 32(6): p. 438-42.
- Ikeda, D.M., et al., MRI Breast Imaging Reporting And Data System Atlas. 2003: 1 ed. Reston: American College of Radiology.
- Kaiser, W., Signs in MR-Mammography. Vol. Auflage: 1. 2007, Berlin: Springer.
- Malich, A., et al., Potential MRI interpretation model: differentiation of benign from malignant breast masses. *AJR Am J Roentgenol*, 2005. 185(4): p. 964-70.
- Schnall, M.D., et al., Diagnostic architectural and dynamic features at breast MR imaging: multicenter study. *Radiology*, 2006. 238(1): p. 42-53.

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