

# Pediatric radius head fracture

Multitom Rax Real3D<sup>1</sup> clinical case  
Children's and Youth Hospital "Auf der Bult" Hannover



<sup>1</sup>Option

Results from case studies are not predictive of results in other cases. Results in other cases may vary.



Study ID 5aad369

# Clinical background and indication for Multitom Rax Real3D<sup>1</sup> examination

## Patient

Male | Age range 10 - 15 years

## Anamnesis


Patient jumped on the trampoline and fell poorly. Significant pain in the proximal forearm. Unclear position of the fracture on conventional X-ray examination.

## Indication for Real3D<sup>1</sup> examination

Better assessment of the fragment position for classification needed.



*AP and lateral view  
Conventional X-ray examination*



*Position the patient with head first for left elbow.*

*Fixate the arm of the patient.*

*The elbow must be covered by the light field.*

*Adapt lateral collimation and adjust table height if necessary.*

The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

<sup>1</sup> Option

# Multitom Rax Real3D<sup>1</sup>

## Settings



Study ID 5aad369

### Settings for tableside scan with Standard Protocol

<i>Tube voltage</i>	81 kV
<i>Current time product</i>	157 mAs
<i>Dose area product</i>	113.7 $\mu\text{Gy}\cdot\text{m}^2$
<i>Calculated value for <math>\text{CTDI}_{\text{vol},32}</math></i>	1.6 mGy
<i>Scan time</i>	12 sec
<i>Number of projections</i>	314

### Reconstruction settings for sectional views

<i>Pixel size</i>	0.4 mm
<i>Reconstruction kernel</i>	sharp (equivalent to Br69)
<i>Slice thickness</i>	2 mm

The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

<sup>1</sup> Option

# Multitom Rax Real3D<sup>1</sup>

## Diagnostic findings

Assessment and position of radius head fracture on Real3D views possible.

Further treatment and therapy for patient with physical rest and immobilization of the affected extremity and an X-ray control of the position in the course.



*Real3D sagittal view on  
2 mm slice thickness*



*Real3D sagittal view on  
0.5 mm slice thickness*

Study ID 5aad369

The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

<sup>1</sup> Option



*“Image contrast and spatial resolution in normal scan mode of Real3D are comparable to conventional multidetector CT examinations. Bony structures and fragments are detectable and assessable.”<sup>1</sup>*

**Jürgen Weidemann, MD**

Chief Physician Pediatric Radiology, Children's and Youth Hospital "Auf der Bult" Hannover

<sup>1</sup>The statements by Siemens Healthineers customers described herein are based on results that were achieved in the customer's unique setting. Since there is no “typical” hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.





The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

Results from case studies are not predictive of results in other cases. Results in other cases may vary.

Dr. Jürgen Weidemann is employed by an institution that receives financial support from Siemens Healthineers for collaborations.