

DICOM Conformance Statement

MULTIX Impact

VA10, VA11, VA20 and higher

[siemens-healthineers.com](https://www.siemens-healthineers.com)



Table of Contents

| | |
|---|----|
| Table of Contents | 2 |
| 1 Conformance Statement Overview | 4 |
| 2 Introductions | 6 |
| 2.1 Intended Audience..... | 6 |
| 2.2 SCOPE AND FIELD OF APPLICATION | 6 |
| 2.3 Remarks..... | 6 |
| 2.4 Terms and Definitions | 7 |
| 2.5 Abbreviations | 10 |
| 2.6 Reference | 10 |
| 3 Networking..... | 11 |
| 3.1 Implementation Model | 11 |
| 3.1.1 Application Data Flow..... | 11 |
| 3.1.2 Functional Definition of AEs..... | 12 |
| 3.1.3 Sequencing of Real-World Activities | 13 |
| 3.2 AE Specifications | 14 |
| 3.2.1 Modality Worklist Application Entity Specification | 14 |
| 3.2.2 Modality Performed Procedure Step Application Entity Specifications | 19 |
| 3.2.3 Storage Application Entity Specification..... | 25 |
| 3.2.4 Storage Commitment Application Entity Specification..... | 30 |
| 3.2.5 Query Application Entity Specification | 33 |
| 3.2.6 Print Application Entity Specification..... | 36 |
| 3.2.7 Verification Application Entity Specification | 43 |
| 3.3 Network Interfaces | 44 |
| 3.3.1 Physical Network Interface..... | 44 |
| 3.3.2 Additional Protocols | 45 |
| 3.3.3 IPv4 and IPv6 Support..... | 46 |
| 3.4 Configuration | 46 |
| 3.4.1 AE Title/Presentation Address Mapping | 46 |
| 3.4.2 Parameters | 47 |
| 4 Media Interchange | 49 |
| 4.1 Implementation Model | 49 |
| 4.1.1 Application Data Flow..... | 49 |
| 4.1.2 Functional Definition of AEs..... | 49 |
| 4.1.3 Sequencing of Real-World Activities | 49 |
| 4.1.4 File Meta Information Options | 49 |
| 4.2 AE Specifications | 50 |

| | |
|---|-----------|
| 4.2.1 Offline-Media Application Entity Specification | 50 |
| 4.3 Augmented and Private Application Profiles | 51 |
| 4.4 Media Configuration | 51 |
| 5 Support of Character Sets | 52 |
| 6 Security | 53 |
| 7 Annexes | 54 |
| 7.1 IOD Contents | 54 |
| 7.1.1 Created SOP Instances | 54 |
| 7.1.2 Used Fields in Received IOD by Application | 67 |
| 7.1.3 Attribute Mapping | 67 |
| 7.1.4 Coerced/Modified Fields | 68 |
| 7.2 Coded Terminology and Templates | 68 |
| 7.3 Grayscale Image Consistency | 68 |
| 7.4 Standard Extended / Specialized / Private SOP Classes | 68 |
| 7.4.1 X-Ray Digital Image Storage SOP Class | 68 |
| 7.5 Private Transfer Syntaxes | 68 |
| 8 List of Figures | 69 |
| 9 List of Tables | 70 |

1 Conformance Statement Overview

This document is the DICOM Conformance Statement for MULTIX Impact Image System (IS). Note that the format of this document strictly follows the format defined in DICOM Standard PS 3.2 (Conformance). Please refer to that part of the standard while reading this document.

Image System is an integrated digital radiography operating console to support the necessary DICOM Services to allow a smooth integration into the clinical network.

With DICOM services to export/import images to/from remote workstations:

The DICOM service "Storage" indicates established image storage in the DICOM image archive; and the DICOM service "Storage Commitment (Push Model)" can request storage commitment from the remote node.

With "DICOM Basic Worklist Management Service Class" the worklist containing patient data can be retrieved from the remote hospital or radiology department information system and afterwards examination data can be sent back to the information system by using "Modality Performed Procedure Step Service Class".

Finally, IS supports the "Print Management Service Class" and the "Media Storage Service Class".

Table 1-1 provides an overview of the network services supported by IS.

Table 1-1 Network Services

| SOP Classes | User of Service (SCU) | Provider of Service (SCP) |
|--|-----------------------|---------------------------|
| Transfer | | |
| Computed Radiography Image Storage | Yes | Yes |
| Digital X-Ray Image Storage - For Presentation | Yes | Yes |
| Digital X-Ray Image Storage - For Processing | Yes | Yes |
| Secondary Capture Image Storage | Yes | Yes |
| X-Ray Radiation Dose SR | Yes | No |
| Query | | |
| Study Root Query Information Model - FIND | Yes | No |
| Study Root Query Information Model - MOVE | Yes | No |
| Workflow Management | | |
| Modality Worklist Information Model - FIND | Yes | No |
| Storage Commitment Push Model | Yes | No |
| Modality Performed Procedure Step | Yes | No |
| Print Management | | |
| Basic Grayscale Print Management Meta | Yes | No |
| Basic Film Session | Yes | No |
| Basic Film Box | Yes | No |
| Basic Grayscale Image Box | Yes | No |
| Printer | Yes | No |

Table 1-2 provides an overview of the Media Storage Application Profiles supported by IS.

Table 1-2 Media Services

| Media Storage Application Profile | Write Files (FSC or FSU) | Read Files (FSR) |
|-----------------------------------|--------------------------|------------------|
| General Purpose CD-R | Yes | Yes |

| Media Storage Application Profile | Write Files (FSC or FSU) | Read Files (FSR) |
|---|--------------------------|------------------|
| General Purpose DVD-RAM | Yes | Yes |
| General Purpose USB Media Interchange with JPEG | Yes | Yes |

2 Introductions

2.1 Intended Audience

The reader of this document is concerned with software design and/or system integration issues. It is assumed that the reader of this document is familiar with the DICOM standard and with the terminology and concepts which are used in those Standards.

2.2 SCOPE AND FIELD OF APPLICATION

This document intends to provide an unambiguous specification for implementations. This specification, called a Conformance Statement, includes a DICOM standard Conformance Statement and is necessary to ensure proper processing and interpretation of medical data exchanged using DICOM standard. The Conformance Statements are available to the public.

The reader of this DICOM Conformance Statement should be aware that different devices can use different Information Object Definitions. For example, a IS may send images using the DX Information Object, CR Information Object, etc.

Included in this DICOM Conformance Statement are the Module Definitions which define all data elements used by this implementation. If the user encounters unspecified private data elements while parsing Data Set, the user is well advised to ignore those data elements (per the DICOM standard). Unspecified private data element information is subject to change without notice. If, however, the device is acting as a "full fidelity storage device", it should retain and re-transmit all the private data elements which are sent by devices.

2.3 Remarks

The use of these DICOM Conformance Statements, in conjunction with the DICOM standard, is intended to facilitate communication with imaging equipment. However, by itself, it is not sufficient to ensure that inter-operation will be successful. The user (or user's agent) needs to proceed with caution and address at least four issues:

- Integration - The integration of any device into an overall system of interconnected devices goes beyond the scope of DICOM standards, and of this introduction and associated DICOM Conformance Statements when interoperability with non-equipment is desired. The responsibility to analyse the applications requirements and to design a solution that integrates imaging equipment with non-manufacture systems is the user's responsibility and should not be underestimated. The user is strongly advised to ensure that such an integration analysis is correctly performed.
- Validation - Testing the complete range of possible interactions between any device and non-manufacture devices, before the connection is declared operational, should not be overlooked. Therefore, the user should ensure that any other provider accepts full responsibility for all validation required for their connection with our devices. This includes the accuracy of the image data once it has crossed the interface between our imaging equipment and the other device and the stability of the image data for the intended applications. Such a validation is required before any clinical use (diagnosis and/or treatment) is performed. It applies when images acquired on our imaging equipment are processed/displayed on somebody else's device, as well as when images acquired on somebody else's equipment is processed/displayed on our console or workstation.
- Future Evolution - The DICOM Standard will evolve to meet the user's growing requirements. DICOM standard will incorporate new features and technologies and we may follow the evolution of the Standard. The protocol is based on DICOM standard as specified in each DICOM Conformance Statement. Evolution of the Standard may require changes to devices which have implemented DICOM standard. In addition, we reserve the right to discontinue or make changes to the support of communications features (on its products) reflected on by these DICOM Conformance Statements. The user should ensure that any other provider, which connects with our devices, also plans evolution of the DICOM Standard. Failures to do so will likely result in the loss of function and/or connectivity as the DICOM Standard changes and

Products are enhanced to support these changes.

- **Interaction** - It is the sole responsibility of the other provider to ensure that communication with the interfaced equipment does not cause degradation of our imaging equipment performance and/or function.

2.4 Terms and Definitions

A

Abstract Syntax: A DICOM term which is identical to a DICOM SOP Class; it identifies a set of SOPs which, when taken together, represent a logical grouping. An Abstract Syntax identifies one SOP Class or Meta SOP Class.

ACR: American College of Radiology.

Annotation Box: A DICOM name for annotation text printed on the film or other media.

ANSI: American National Standards Institute.

Application Entity (AE): A DICOM term for defining a particular user at an IP address.

Association: A DICOM term for a communication context which is used by two Application Entities that communicate to one another.

Association Negotiation: The software handshaking that occurs between two DICOM Application Entities to set up an Association.

Attribute: Each DICOM information object has its own set of characteristics or attributes. Each attribute has a name and may have a value (see IOD), depending on its category.

B

Big Endian: A term for encoding data where the most-significant byte appears first and remaining bytes follow in descending order of significance; sometimes known as "Motorola" format (see Little Endian). (The term is used because of an analogy with the story Gulliver's Travels, in which Jonathan Swift imagined a never-ending fight between the kingdoms of the Big-Endian and the Little-Endian, whose only difference is in where they crack open a hard-boiled egg.)

C

Calling (Requesting) AE Title: The name used by the receiver in a DICOM Association to indicate which Application Entity it received the data from. It is the AE Title of the AE that is initiating the transfer.

Called (Receiving) AE Title: The name used by the sender in a DICOM Association to indicate which Application Entity it wants to transmit its data to. It is the AE Title of the AE that is receiving the transfer.

Command Element: An encoding of a parameter of a command which conveys this parameter's value.

Command Stream: The result of encoding a set of DICOM Command Elements using the DICOM encoding scheme.

Composite Information Object: A DICOM information object (see IOD) whose attributes contain multiple real world objects.

Conformance: Conformance in the DICOM sense means to be in compliance with the parts of the DICOM Standard.

Conformance Statement: A document whose organization and content are mandated by the DICOM Standard, which allows users to communicate how they have chosen to comply with the Standard in their implementations (see Section 8).

Combined Print Image: a pixel matrix created by superimposing an image and an overlay, the size of which is defined by the smallest rectangle enclosing the superimposed image and overlay.

D

Data Dictionary: A registry of DICOM Data Elements which assigns a unique tag, a name, value characteristics, and semantics to each Data Element (see the DICOM Data Element Dictionary in DICOM PS 3.6-2004).

Data Element: A unit of information as defined by a single entry in the data dictionary. An encoded Information Object Definition (IOD) Attribute that is composed of, at a minimum, three fields: a Data Element Tag, a Value Length, and a Value Field. For some specific Transfer Syntaxes, a Data Element also contains a VR Field where the Value Representation of that Data Element is specified explicitly.

Data Set: Exchanged information consisting of a structured set of Attribute values directly or indirectly related to Information Objects. The value of each Attribute in a Data Set is expressed as a Data Element.

Data Stream: The result of encoding a Data Set using the DICOM encoding scheme (Data Element Numbers and representations as specified by the Data Dictionary).

DICOM: Digital Imaging and Communications in Medicine.

DICOM File: A DICOM File is a file with a content formatted according to the requirements of DICOM.

DICOM File Format: The DICOM File Format provides a means to encapsulate in a File the Data Set representing a SOP Instance related to a DICOM Information Object.

DIMSE: DICOM Message Service Element. This represents an abstraction of a common set of things that a user would do to a data element, would likely use over and over, and would appear in various different contexts.

DIMSE-C: DICOM Message Service Element—Composite.

DIMSE-C services: A subset of the DIMSE services which supports operations on Composite SOP Instances related to composite Information Object Definitions with peer DIMSE-service-users.

DIMSE-N: DICOM Message Service Element—Normalized.

DIMSE-N services: A subset of the DIMSE services which supports operations and notifications on Normalized SOP Instances related to Normalized Information Object Definitions with peer DIMSE service-users.

E, F

Film Box: A Normalized Information Object which is the DICOM name for the equivalent of a sheet of physical film.

Film Session: A Normalized Information Object which is the DICOM name for the equivalent of a typical “study” or “series”.

G, H, I

HIS: Hospital Information System.

IE: Information Entity.

Image Box: A Normalized Information Object which is the DICOM name for the equivalent of a typical “frame” or “image”.

Information Object Class or

Information Object [Definition] (IOD): A software representation of a real object (e.g., CT Image, Study, etc.). An Information Object is generally a list of characteristics (Attributes) which completely describe the object as far as the software is concerned. The formal description of an Information Object generally includes a description of its purpose and the Attributes it possesses.

Information Object Instance or

Instance (of an IOD): A software representation of a specific occurrence of a real object or entity, including values for the Attributes of the Information Object Class to which the entity belongs..

J, K, L

Little Endian: A term for encoding data where the least-significant byte appears first and remaining bytes follow in ascending order of significance; sometimes known as “Intel” format (see Big Endian).

LUT: Lookup Table.

M

Message: A data unit of the Message Exchange Protocol exchanged between two cooperating DICOM Application Entities. A Message is composed of a Command Stream followed by an optional Data Stream.

Meta SOP Class: A collection or group of related SOP Classes identified by a single Abstract Syntax UID, which, when taken together, represent a logical grouping and which are used together to provide a high-level functionality, e.g., for the purpose of negotiating the use of the set with a single item.

Module: A logical group of the valid attributes of DICOM information objects.

N

NEMA: National Electrical Manufacturers Association.

Normalized Information Object: A DICOM Information Object (see IOD) whose attributes contain a single real

Copyright © SSME 2021. All rights reserved. Restricted.

world object. Note: the differentiation of normalized versus composite information object definitions is not strongly enforced in DICOM 3.0.

O, P

Presentation Context: A Presentation Context consists of an Abstract Syntax plus a list of acceptable Transfer Syntaxes. The Presentation Context defines both what data will be sent (Abstract Syntax) and how the data are encoded to be sent (Transfer Syntax).

Print Job SOP Class: A DICOM representation of a Print Job which consists of a set of IODs which describe a Print Job and a set of services which can be performed on those IODs.

Print Management Service Class or Print Service Class (PSC): A DICOM term for a logical grouping of Service Classes which all involve printing, also referred to as Print Management Service Class (an example of a Meta SOP Class).

Printer SOP Class: A DICOM representation of a Printer which consists of a set of IODs which describe a Printer and a set of services which can be performed on those IODs.

Protocol Data Unit (PDU): A data object which is exchanged by software protocol devices (entities, machines) within a given layer of the protocol stack.

Q, R

RIS: Radiology Information System.

Real-World Activity: Something which exists in the real world and which pertains to specific area of information processing within the area of interest of the DICOM Standard. A Real-World Activity may be represented by one or more SOP Classes.

Real-World Object: Something which exists in the real world and upon which operations may be performed which are within the area of interest of the DICOM Standard. A Real-World Object may be represented through a SOP Instance.

S

Service Class: A group of operations that a user might want to perform on particular Information Objects. Formally, a structured description of a service which is supported by cooperating DICOM Application Entities using specific DICOM Commands acting on a specific class of Information Object.

Service Class Provider (SCP, Provider, Server): A device which provides the services of a DICOM Service Class or Classes which are utilized by another device (SCU) and which performs operations and invokes notifications on a specific Association.

Service Class User (SCU, User, Client): A device which utilizes the DICOM Service Class or Classes which are provided by another device (SCP) and which invokes operations and performs notifications on a specific Association.

Service-Object Pair (SOP): The combination of a DICOM Information Object and the Service Class which operates upon that object.

SOP Class: A DICOM term which is identical to an Abstract Syntax; it identifies a set of SOPs which, when taken together, represent a logical grouping (see Meta SOP Class).

Storage Service Class (SSC): A DICOM term for a logical grouping of Service Classes which all involve storage of images.

T

Tag: A unique identifier for an element of information composed of an ordered pair of numbers (a Group Number followed by an Element Number), which is used to identify Attributes and corresponding Data Elements.

TCP/IP: Transmission Control Protocol / Internet Protocol.

Transfer Syntax: A part of the DICOM Presentation Context which specifies a set of encoding rules that allow Application Entities to unambiguously negotiate the encoding techniques (e.g., Data Element structure, byte ordering, compression) they are able to support, thereby allowing these Application Entities to communicate.

U

Unique Identifier (UID): A globally unique identifier (based on the structure defined by ISO 8824 for OSI Object

Identifiers) which is assigned to every DICOM information object as specified by the DICOM Standard (see Section 2.1.1.4) and which guarantees global unique identification for objects across multiple countries, sites, vendors and equipment.

V

Value Representation (VR): A VR is the defined format of a particular data element.

W, X, Y, Z

2.5 Abbreviations

| | |
|-----------|--|
| ACC | American College of Cardiology |
| ACR | American College of Radiology |
| ASCII | American Standard Code for Information Interchange |
| AE | Application Entity |
| ANSI | American National Standards Institute |
| CEN TC251 | Comite Europeen de Normalisation – Technical Committee 251 – Medical Informatics |
| DICOM | Digital Imaging and Communications in Medicine |
| DIMSE | DICOM Message Service Element |
| DIMSE-C | DICOM Message Service Element - Composite |
| DIMSE-N | DICOM Message Service Element – Normalized |
| HIS | Hospital Information System |
| HL7 | Health Level 7 |
| IE | Information Entity |
| IOD | Information Object Definition |
| ISO | International Standard Organization |
| NEMA | National Electrical Manufacturers Association |
| OSI | Open Systems Interconnection |
| PDU | Protocol Data Unit |
| RIS | Radiology Information System |
| SCP | Service Class Provider |
| SCU | Service Class User |
| SOP | Service-Object Pair |
| TCP/IP | Transmission Control Protocol/Internet Protocol |
| UID | Unique Identifier |

2.6 Reference

- [1] [DICOM] Digital Imaging and Communications in Medicine (DICOM), Part 1 – 16
- [2] NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <http://medical.nema.org/>

3 Networking

3.1 Implementation Model

3.1.1 Application Data Flow




Figure 1 Application Data Flow Diagram

- The Modality Worklist Application Entity receives Worklist information from a remote AE. It is associated with the local real-world activities "Worklist Update" and "Acquire Images". When the "Worklist Update" local real-world activity is performed the Modality Worklist Application Entity queries a remote AE for worklist items and provides the set of worklist items matching the query request. "Worklist Update" is performed as a result of an operator request or can be performed automatically at specific time intervals.
- MPPS Application Entity creates and updates Modality Performed Procedure Step instances (sends MPPS information to a remote AE).
- The Storage Application Entity sends images to a remote AE. It is associated with the local real-world activity "Send Images". "Send Images" is performed upon user request for each study completed or for specific images selected. When activated by user's settings (auto-send), each marked set of images and RDSR can be immediately stored to a preferred destination whenever a Patient/Study is closed by the user.
- If the remote AE is configured as an archive device the Storage AE will request Storage Commitment and if a commitment is successfully obtained will record this information in the local database.
- The Print Application Entity prints images on a remote AE (Printer). It is associated with the local real-world activity "Print Images". "Print Images" creates a print-job within the print queue containing one or more virtual film sheets composed from images selected by the user.

3.1.2 Functional Definition of AEs

3.1.2.1 Functional Definition of Modality Worklist Application Entity

Worklist Update attempts to download a Worklist from a remote node. If the Modality Worklist AE establishes an Association to a remote AE, it will transfer all worklist items via the open Association. During receiving the worklist response items are counted and the results will be displayed in a separate list, which will be cleared with the next Worklist Update.

3.1.2.2 Functional Definition of MPPS Application Entity

The MPPS AE performs the creation of a MPPS Instance automatically. Further updates on the MPPS data can be performed interactively from the related MPPS user interface. The MPPS "Complete" or "Discontinued" states can only be set from the user interface.

3.1.2.3 Functional Definition of Storage Application Entity

The existence of a send-job queue entry with associated network destination will activate the Storage AE. An association request is sent to the destination AE and upon successful negotiation of a Presentation Context the image transfer is started. If the association cannot be opened, the related send-job is set to an error state and can be restarted by the user via job control interface. By default, the Storage AE will not try to initiate another association for this send-job automatically. However, an automatic retry (retry-timer, retry-count) can be configured.

3.1.2.4 Functional Definition of Storage Commitment Application Entity

The Storage Commitment AE accepts an association for Storage Commitment notification (N-EVENT-REPORT) only as a SCU. The Storage Commitment Provider initiating the association must use the role selection negotiation.

3.1.2.5 Functional Definition of Query Application Entity

IS DICOM query SCU requests the remote query SCP to perform a search (C-FIND) and match to the keys specified in the request in order to display the results in the user interface. Depending on user action (query data are input) the IS DICOM query SCU sends a C-MOVE DIMSE service to initiate a C-STORE sub-operation on the SCP to start an image transfer from remote Storage SCU (running on Query SCP) to local Storage SCP. The user has the possibility to cancel a running import by an "Import Cancel" function. IS DICOM application will perform only one running import operation at a time.

3.1.2.6 Functional Definition of Print Application Entity

The existence of a print-job in the print queue will activate the Print AE. An association is established with the printer and the printer's status determined. If the printer is operating normally, the film sheets described within

the print-job will be printed. Changes in printer status will be detected and reported to the user. If the printer is not operating normally, the print-job will set to an error state and can be restarted by the user via the job control interface.

3.1.2.7 Functional Definition of Verifications Application Entity

IS DICOM application opens an association when a "C-ECHO" of a remote application is requested from the Service User Interface. This can be done to verify the correct setup of a remote destinations configuration data.

3.1.3 Sequencing of Real-World Activities




Figure 2 Sequencing Constraints

Under normal scheduled Modality Worklist conditions the sequencing constraints illustrated in Figure 2 apply:

1. Query Worklist
2. Receive Worklist
3. Select Work item from Worklist
4. Start acquisition and create MPPS
5. Acquire Images
6. Print acquired images (optional step)
7. Complete acquisition and finalize MPPS
8. Store acquired images
9. If the Storage Provider is configured as an archive device the Storage AE will request Storage Commitment for the images

Other workflow situations (e.g., unscheduled procedure steps) will have other sequencing constraints. Printing could equally take place after the acquired images have been stored. Printing could be omitted completely if no printer is connected, or hard copies are not required.

3.2 AE Specifications

3.2.1 Modality Worklist Application Entity Specification

3.2.1.1 SOP Classes

IS provides Standard Conformance to the following SOP Classes as SCU:

Table 3-1 SOP Classes for AE Modality Worklist

| SOP Class Name | SOP Class UID |
|--|------------------------|
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 |

3.2.1.2 Association Policies

3.2.1.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 3-2 DICOM Application Context

| | |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

3.2.1.2.2 Number of Associations

IS initiates one Association at a time for a Worklist request.

3.2.1.2.3 Asynchronous Nature

IS does not support asynchronous communication (multiple outstanding transactions over a single Association).

3.2.1.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 3-3 DICOM Implementation Identifying Information

| | |
|-----------------------------|----------------------|
| Implementation Class UID | 1.2.840.69677977.1.0 |
| Implementation Version Name | EPACS20001.0 |

3.2.1.3 Association Initiation Policy

3.2.1.3.1 Activity - Worklist Update

3.2.1.3.1.1 Description and Sequencing of Activities

The request for a Worklist Update is initiated by user interaction, pressing the buttons "Query" on the user interface. The Modality Worklist AE will then initiate an association with the remote AE in order to query for the worklist. A user can configure a number of parameters which directly control the worklist query request. The user can request worklist items that are intended for the system the user is working at, all items that apply to the modality of the system the user is working at or all worklist items available. These selections and their effects on worklist query parameters are given below:

The interactive Patient Worklist Query will display a dialog for entering data as search criteria. When the Query is started on user request, only the data from the dialog will be inserted as matching keys into the query.

With automated worklist queries (including "Worklist Update") IS always requests all items for a Scheduled Procedure Step Start Date (actual date), Modality (DX/CR) and Scheduled Station AE Title. Query for the Scheduled Station AE Title is configurable by a Service Engineer.

Upon initiation of the request, IS will build an Identifier for the C-FIND request, will initiate an Association to send the request and will wait for Worklist responses. After retrieval of all responses, IS will access the local database to add or update patient demographic data.

IS will initiate an Association in order to issue a C-FIND request according to the Modality Worklist Information Model.




Figure 3 Sequencing of Activity - Worklist Update

The sequence of interactions between the Modality Worklist AE and a MWL/MPPS Provider (e.g., a device such as a RIS or HIS that supports the Modality Worklist SOP Class as an SCP) is illustrated in the Figure above:

1. The Modality Worklist AE opens an association with the MWL/MPPS Provider
2. The Modality Worklist AE sends a C-FIND request to the MWL/MPPS Provider containing the Worklist Query attributes.
3. The MWL/MPPS Provider returns a C-FIND response containing the requested attributes of the first matching Worklist Item.
4. The MWL/MPPS Provider returns another C-FIND response containing the requested attributes of the second matching Worklist Item.
5. The MWL/MPPS Provider returns another C-FIND response with status Success indicating that no further matching Worklist Items exist. This example assumes that only 2 Worklist items match the Worklist Query.
6. The Modality Worklist AE closes the association with the MWL/MPPS Provider.
7. The user selects a Worklist Item from the Worklist and prepares to acquire new images.

3.2.1.3.1.2 Proposed Presentation Contexts

IS will propose Presentation Contexts as shown in the following table:

Table 3-4 Proposed Presentation Contexts for Activity Worklist Update

| Presentation Context Table | | | | | |
|---------------------------------------|------------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Modality Worklist Information Model - | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

| Presentation Context Table | | | | | |
|----------------------------|-----|-----------------|----------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| FIND | | | | | |

3.2.1.3.1.3 SOP Specific Conformance for Modality Worklist

The behavior of IS when encountering status codes in a Modality Worklist C-FIND response is summarized in the table below. If any other SCP response status than "Success" or "Pending" is received by IS, a message "query failed" will appear on the user interface.

Table 3-5 Modality Worklist C-FIND Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|--|------------------------|--|
| Success | Matching is complete | 0000 | The SCP has completed the matches. Worklist items are available for display or further processing. |
| Refused | Out of Resources | A700 | The Association is aborted and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |
| Failed | Identifier does not match SOP Class | A900 | The Association is aborted and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |
| Failed | Unable to Process | C000 - CFFF | The Association is aborted and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |
| Cancel | Matching terminated due to Cancel request | FE00 | If the query was cancelled due to too many worklist items then the SCP has completed the matches. Worklist items are available for display or further processing. Otherwise, The Association is aborted and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. |
| Pending | Matches are continuing | FF00 | The worklist item contained in the Identifier is collected for later display or further processing. |
| Pending | Matches are continuing - Warning that one or more Optional Keys were not supported | FF01 | The worklist item contained in the Identifier is collected for later display or further processing. The status meaning is logged only once for each C-FIND operation. |
| * | * | Any other status code. | The Association is aborted and the worklist is marked as failed. The status |

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|-----------------|------------|--|
| | | | meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |

The behavior of IS during communication failure is summarized in the table below.

Table 3-6 Modality Worklist Communication Failure Behavior

| Exception | Behavior |
|--|--|
| Association aborted by the SCP or network layers | The worklist query is marked as failed. The reason is logged and reported to the user if an interactive query. |

Acquired images will always use the Study Instance UID specified for the Scheduled Procedure Step (if available). If an acquisition is unscheduled, a Study Instance UID will be generated locally.

The table below provides a description of the IS Worklist Request Identifier and specifies the attributes that are copied into the images. Unexpected attributes returned in a C-FIND response are ignored.

Requested return attributes not supported by the SCP are set to have no value. Non-matching responses returned by the SCP due to unsupported optional matching keys are ignored. No attempt is made to filter out possible duplicate entries.

Table 3-7 Worklist Request Identifier

| Module Name Attribute Name | Tag | VR | M | R | Q | D | IOD |
|--|-------------|----|---|---|---|---|-----|
| SOP Common | | | | | | | |
| Specific Character Set | (0008,0005) | CS | | x | | | |
| Scheduled Procedure Step | | | | | | | |
| Scheduled Procedure Step Sequence | (0040,0100) | SQ | | x | | | |
| >Scheduled Station AE Title | (0040,0001) | AE | S | x | | x | |
| >Scheduled Procedure Step Start Date | (0040,0002) | DA | S | x | | x | |
| >Scheduled Procedure Step Start Time | (0040,0003) | TM | S | x | | x | |
| >Pre-Medication | (0040,0012) | LO | | x | | x | |
| >Modality | (0008,0060) | CS | S | x | | | |
| >Scheduled Performing Physician's Name | (0040,0006) | PN | | x | | | x |
| >Scheduled Protocol Code Sequence | (0040,0008) | SQ | | x | | | x |
| >> Code Value | (0008,0100) | SH | | x | | | x |
| >> Coding Scheme Designator | (0008,0102) | SH | | x | | | x |
| >> Coding Scheme Version | (0008,0103) | SH | | x | | | x |
| >> Code Meaning | (0008,0104) | LO | | x | | | x |
| >Scheduled Procedure Step ID | (0040,0009) | SH | | x | | x | x |
| >Scheduled Procedure Step Description | (0040,0007) | LO | | x | | x | x |
| Requested Procedure | | | | | | | |
| Requested Procedure Comments | (0040,1400) | LT | | | | | |
| Requested Procedure ID | (0040,1001) | SH | | x | | x | x |
| Requested Procedure Description | (0032,1060) | LO | | x | | x | x |
| Requested Procedure Code Sequence | (0032,1064) | SQ | | x | | | x |

| Module Name Attribute Name | Tag | VR | M | R | Q | D | IOD |
|--|-------------|----|---|---|---|---|-----|
| > Code Value | (0008,0100) | SH | | x | | | x |
| > Coding Scheme Designator | (0008,0102) | SH | | x | | | x |
| > Coding Scheme Version | (0008,0103) | SH | | x | | | x |
| > Code Meaning | (0008,0104) | LO | | x | | | x |
| Names of Intended Recipients of Results | (0040,1010) | PN | | | | | |
| Study Instance UID | (0020,000D) | UI | | x | | | x |
| Requested Procedure Priority | (0040,1003) | SH | | x | | | x |
| Referenced Study Sequence | (0008,1110) | SQ | | x | | | x |
| > Referenced SOP Class UID | (0008,1150) | UI | | x | | | x |
| > Referenced SOP Instance UID | (0008,1155) | UI | | x | | | x |
| Imaging Service Request | | | | | | | |
| Accession Number | (0008,0050) | SH | | x | x | x | x |
| Requesting Physician | (0032,1032) | PN | | x | | | x |
| Referring Physician's Name | (0008,0090) | PN | | x | | | x |
| Requesting Service | (0032,1033) | LO | | x | | | |
| Imaging Service Request Comments | (0040,2400) | LT | | x | | | |
| Visit Identification | | | | | | | |
| Admission ID | (0038,0010) | LO | | x | | | |
| Institution Name | (0008,0080) | LO | | x | | x | x |
| Visit Status | | | | | | | |
| Current Patient Location | (0038,0300) | LO | | x | | | |
| Visit Relationship | | | | | | | |
| Referenced Patient Sequence | (0008,1120) | SQ | | x | | x | x |
| > Referenced SOP Class UID | (0008,1150) | UI | | x | | x | x |
| > Referenced SOP Instance UID | (0008,1155) | UI | | x | | x | x |
| Visit Admission | | | | | | | |
| Admitting Diagnosis Description | (0008,1080) | LO | | x | | x | |
| Patient Identification | | | | | | | |
| Patient Name | (0010,0010) | PN | | x | x | x | x |
| Patient ID | (0010,0020) | LO | | x | x | x | x |
| Other Patient Ids | (0010,1000) | LO | | x | | | x |
| Patient Demographic | | | | | | | |
| Patient's Birth Date | (0010,0030) | DA | | x | x | x | x |
| Patient's Sex | (0010,0040) | CS | | x | x | x | x |
| Confidentiality constraint on patient data | (0040,3001) | LO | | x | | | |
| Ethnic Group | (0010,2160) | SH | | x | | | x |
| Patient Comments | (0010,4000) | LT | | x | | | x |

| Module Name Attribute Name | Tag | VR | M | R | Q | D | IOD |
|-------------------------------|-------------|----|---|---|---|---|-----|
| Patient Address | (0010,1040) | LO | | x | | | |
| Patient Medical | | | | | | | |
| Patient State | (0038,0500) | LO | | x | | | |
| Pregnancy Status | (0010,21C0) | US | | x | | | |
| Medical Alerts | (0010,2000) | LO | | x | | | |
| Allergies | (0010,2110) | LO | | x | | | |
| Special Needs | (0038,0050) | LO | | x | | | |
| Additional Patient History | (0010,21B0) | LT | | x | | | |

The above table should be read as follows:

Module Name The name of the associated module for supported worklist attributes.

Attribute Name Attributes supported to build an EXAMPLEINTEGRATED-MODALITY Worklist Request Identifier.

Tag DICOM tag for this attribute.

VR DICOM VR for this attribute.

M Matching keys for (automatic) Worklist Update. A "S" will indicate that IS will supply an attribute value for Single Value Matching, a "R" will indicate Range Matching and a "*" will denote wild card matching. It can be configured if "Scheduled Station AE Title" is additionally supplied "(S)" and if Modality is set to DX or CR.

R Return keys. An "x" will indicate that IS will supply this attribute as Return Key with zero length for Universal Matching. IS will support retired date format (yyyy.mm.dd) for "Patient's Birth Date" and "Scheduled Procedure Step Start Date" in the response identifiers. For "Scheduled Procedure Step Start Time" also retired time format as well as unspecified time components are supported.

Q Interactive Query Key. An "x" will indicate that IS will supply this attribute as matching key, if entered in the Query Patient Worklist dialog. For example, the Patient Name can be entered thereby restricting Worklist responses to Procedure Steps scheduled for the patient.

D Displayed keys. An "x" indicates that this worklist attribute is displayed to the user during a patient registration dialog. For example, Patient Name will be displayed when registering the patient prior to an examination.

IOD An "x" indicates that this Worklist attribute is included into all Object Instances created during performance of the related Procedure Step.

The default Query Configuration is set to "Modality" and "Date". Optionally, matching for the own "Modality", "AE Title" and "Date" is configurable.

For the Modality Type attribute (0008,0060) of the DICOM Modality Worklist query any combination of the following values may be configured:

- CR
- DX

For "Date" one of the following settings could be configured:

- "1 Day"
- "2 Days"
- "1 Week"
- "All The Time"

3.2.1.4 Association Acceptance Policy

The Modality Worklist Application Entity does not accept Associations.

3.2.2 Modality Performed Procedure Step Application Entity Specifications

3.2.2.1 SOP Classes

IS provides Standard Conformance to the following SOP Classes as SCU:

Table 3-8 SOP Classes for AE MPPS

| | |
|-----------------------------------|-------------------------|
| SOP Class Name | SOP Class UID |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 |

3.2.2.2 Association Policies

3.2.2.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 3-9 DICOM Application Context

| | |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

3.2.2.2.2 Number of Associations

IS initiates one Association at a time for a communication of MPPS information.

3.2.2.2.3 Asynchronous Nature

IS does not support asynchronous communication (multiple outstanding transactions over a single Association).

3.2.2.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 3-10 DICOM Implementation Class and Version

| | |
|-----------------------------|----------------------|
| Implementation Class UID | 1.2.840.69677977.1.0 |
| Implementation Version Name | EPACS20001.0 |

3.2.2.3 Association Initiation Policy

3.2.2.3.1 Activity - Acquire Images

3.2.2.3.1.1 Description and Sequencing of Activities

After Patient registration, IS is awaiting the 1st application of X-Ray Dose to the patient. The trigger to create a MPPS SOP Instance is derived from this event. An Association to the configured MPPS SCP system is established immediately and the related MPPS SOP Instance will be created.

A manual update can be performed with the MPPS user interface where it is possible to set the final state of the MPPS to "COMPLETED" or "DISCONTINUED". In the "Discontinued" case the user can also select the discontinuation reason from a list corresponding to CID 9300 "Procedure Discontinuation Reasons". A MPPS Instance that has been sent with a state of "COMPLETED" or "DISCONTINUED" can no longer be updated.

IS will support creation of "unscheduled cases" by allowing MPPS Instances to be communicated for locally registered Patients.

IS will initiate an Association to issue an:

- N-CREATE request according to the CREATE Modality Performed Procedure Step SOP Instance operation or a
- N-SET request to update the contents and state of the MPPS according to the SET Modality Performed Procedure Step Information operation.




Figure 4 Sequencing of Activity - Acquire Images

A possible sequence of interactions between the MPPS AE and a MWL/MPPS Provider (e.g., a device such as a RIS or HIS that supports the MPPS SOP Class as an SCP) is illustrated in Figure 4:

1. The MPPS AE opens an association with the MWL/MPPS Provider.
2. The MPPS AE sends an N-CREATE request to the MWL/MPPS Provider to create an MPPS instance with status of "IN PROGRESS" and create all necessary attributes. The MWL/MPPS Provider acknowledges the MPPS creation with an N-CREATE response (status success).
3. The MPPS AE closes the association with the MWL/MPPS Provider.
4. All images are acquired and stored in the local database.
5. The MPPS AE opens an association with the MWL/MPPS Provider.
6. The MPPS AE sends an N-SET request to the MWL/MPPS Provider to update the MPPS instance with status of "COMPLETED" and set all necessary attributes. The MWL/MPPS Provider acknowledges the MPPS update with an N-SET response (status success).
7. The MPPS AE closes the association with the MWL/MPPS Provider.

3.2.2.3.1.2 Proposed Presentation Contexts

IS will propose Presentation Contexts as shown in the following table:

Table 3-11 Proposed Presentation Contexts for Real-World Activity Acquire Images

| Presentation Context Table | | | | | |
|-----------------------------------|-------------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

3.2.2.3.1.3 SOP Specific Conformance for MPPS

The behavior of IS when encountering status codes in an MPPS N-CREATE or N-SET response is summarized in Table 4.2-26. If any other SCP response status than "Success" or "Warning" is received by EXAMPLE-INTEGRATED-MODALITY, a message "MPPS update failed" will appear on the user interface.

Table 3-12 MPPS N-CREATE / N-SET Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|---|------------|--|
| Success | Success | 0000 | The SCP has completed the operation successfully. |
| Failure | Processing Failure - Performed Procedure Step Object may no longer be updated | 0110 | The Association is aborted and the MPPS is marked as failed. The status meaning is logged and reported to the user. Additional information in the Response will be logged (i.e., Error Comment and Error ID). |
| Warning | Attribute Value Out of Range | 0116H | The MPPS operation is considered successful but the status meaning is logged. Additional information in the Response identifying the attributes out of range will be logged (i.e., Elements in the Modification List/Attribute List) |

The behavior of IS during communication failure is summarized in the table below:

Table 3-13 MPPS Communication Failure Behavior

| Exception | Behavior |
|--|--|
| Timeout | The Association is aborted and MPPS marked as failed. The reason is logged and reported to the user. |
| Association aborted by the SCP or network layers | The MPPS is marked as failed. The reason is logged and reported to the user. |

Table 3-14 provides a description of the MPPS N-CREATE and N-SET request identifiers sent by IS. Empty cells in the N-CREATE and N-SET columns indicate that the attribute is not sent. An "x" indicates that an appropriate value will be sent. A "Zero length" attribute will be sent with zero length.

Table 3-14 MPPS N-CREATE / N-SET Request Identifier

| Attribute Name | Tag | VR | N-CREATE | N-SET |
|--|-------------|----|--|-------|
| Specific Character Set | (0008,0005) | CS | X | X |
| Performed Procedure Step Relationship | | | | |
| Scheduled Step Attributes Sequence | (0040,0270) | SQ | X | |
| > Study Instance UID | (0020,000D) | UI | X | |
| > Referenced Study Sequence | (0008,1110) | SQ | Zero length if not received from RIS | |
| >> Referenced SOP Class UID | (0008,1150) | UI | From RIS | |
| >> Referenced SOP Instance UID | (0008,1155) | UI | X | |
| > Accession Number | (0008,0050) | SH | From Modality Worklist or user input. The user can modify values provided via Modality Worklist. | |
| > Requested Procedure ID | (0040,1001) | SH | Zero length if not available | |
| > Requested Procedure Description | (0032,1060) | LO | Zero length if not received from RIS | |
| > Scheduled Procedure Step ID | (0040,0009) | SH | Zero length if not received | |

Copyright © SSME 2021. All rights reserved. Restricted.

| Attribute Name | Tag | VR | N-CREATE | N-SET |
|--|-------------|----|--|--|
| | | | from RIS | |
| > Scheduled Procedure Step Description | (0040,0007) | LO | Zero length if not received from RIS | |
| > Scheduled Protocol Code Sequence | (0040,0008) | SQ | Zero length if not received from RIS | |
| > Code Value | (0008,0100) | | X | |
| > Coding Scheme Designator | (0008,0102) | | X | |
| > Coding Scheme Version | (0008,0103) | | X | |
| > Code Meaning | (0008,0104) | | X | |
| Patient's Name | (0010,0010) | PN | Zero length if not available | |
| Patient ID | (0010,0020) | LO | Zero length if not available | |
| Patient's Birth Date | (0010,0030) | DA | Zero length if not available | |
| Patient's Sex | (0010,0040) | CS | Zero length if not available | |
| Referenced Patient Sequence | (0008,1120) | SQ | Zero length if not received from RIS | |
| >> Referenced SOP Class UID | (0008,1150) | UI | X | |
| >> Referenced SOP Instance UID | (0008,1155) | UI | X | |
| Performed Procedure Step Information | | | | |
| Performed Procedure Step ID | (0040,0253) | SH | copied from (0040,0009) if received from RIS internal created value otherwise | |
| Performed Station AE Title | (0040,0241) | AE | MPPS AE Title | |
| Performed Station Name | (0040,0242) | SH | same as (0008,1010) in storage objects | |
| Performed Location | (0040,0243) | SH | Zero length | Zero Length |
| Performed Procedure Step Start Date | (0040,0244) | DA | = (0008,0020) (IHE-8) | |
| Performed Procedure Step Start Time | (0040,0245) | TM | = (0008,0030) (IHE-8) | |
| Performed Procedure Step End Date | (0040,0250) | DA | Zero length | only in final N-SET |
| Performed Procedure Step End Time | (0040,0251) | TM | Zero length | only in final N-SET |
| Performed Procedure Step Status | (0040,0252) | CS | IN PROGRESS | X |
| Performed Procedure Step Description | (0040,0254) | LO | = (0008,1030) (IHE-8) | = (0008,1030) (IHE-8) or value from UI |
| Performed Procedure Type Description | (0040,0255) | LO | Zero length | |
| Performed Procedure Step | (0040,0281) | SQ | Zero length | |

| Attribute Name | Tag | VR | N-CREATE | N-SET |
|--|-------------|----|---|---|
| Discontinuation Reason Code Sequence | | | | |
| Procedure Code Sequence | (0008,1032) | | Zero length if not received from RIS Requested procedure code otherwise | Zero length if not received from RIS Requested procedure code otherwise |
| Image Acquisition Results | | | | |
| Modality | (0008,0060) | CS | dependent on 1st image of procedure and configuration: "CR" or "DX" | |
| Study ID | (0020,0010) | SH | X | |
| Performed Protocol Code Sequence | (0040,0260) | SQ | Zero length | shall only be present if the performed Exams match to a saved Protocol Codes |
| > Code Value | (0008,0100) | | | X |
| > Coding Scheme Designator | (0008,0102) | | | X |
| > Coding Scheme Version | (0008,0103) | | | X |
| > Code Meaning | (0008,0104) | | | X |
| Performed Series Sequence | (0040,0340) | SQ | Zero length | One item for all performed series, one item for the according DICOM Dose Report |
| > Performing Physician's Name | (0008,1050) | PN | | x |
| > Protocol Name | (0018,1030) | LO | | x |
| > Operator's Name | (0008,1070) | PN | | x |
| > Series Instance UID | (0020,000E) | UI | | x |
| > Series Description | (0008,103E) | LO | | x |
| > Retrieve AE Title | (0008,0054) | AE | | x |
| > Referenced Image Sequence | (0008,1140) | SQ | | One or more items |
| >> Referenced SOP Class UID | (0008,1150) | | | |
| >> Referenced SOP Instance UID | (0008,1155) | | | |
| > Referenced Non-Image Composite SOP Instance Sequence | (0040,0220) | | | Only for dose report |
| >> Referenced SOP Class UID | (0008,1150) | | | X-Ray Radiation Dose SR |
| >> Referenced SOP Instance UID | (0008,1155) | | | UID OF DICOM DOSE SR |
| > Referenced Standalone SOP Instance Seq. | (0040,0220) | SQ | Zero length (SOP classes not supported) | Zero length (SOP classes not supported) |
| Radiation Dose Module | | | | |
| Total Time of Fluoroscopy | (0040,0300) | US | Zero length | Total time |
| Total Number of Exposures | (0040,0301) | US | Zero length | Number of exposures, includes cassette images for series |

| Attribute Name | Tag | VR | N-CREATE | N-SET |
|-----------------------------------|-------------|----|-------------|--|
| | | | | acquisition each frame is counted |
| Distance Source to Detector (SID) | (0018,1110) | DS | Zero length | minimum value |
| Distance Source to entrance | (0040,0306) | | | minimum value |
| Image Area Dose Product | (0018,115E) | DS | Zero length | includes fluoro and cassette images |
| Exposure Dose Sequence | (0040,030E) | | Zero length | |
| Comments on Radiation Dose | (0040,0310) | | Zero length | |
| Entrance Dose in mGy | (0040,8302) | | | -sum of irradiation event entrance dose |
| Billing and Material Code Module | | | | |
| Billing Procedure Step Sequence | (0040,0320) | | Zero length | Zero length |
| Film Consumption Sequence | (0040,0321) | SQ | Zero length | one item for each film size. Shall only be present if Number of Films for this film size > 0 |
| > Film Size ID | (2010,0050) | CS | | X |
| > Number of Films | (2100,0170) | IS | | Expected film number |

3.2.2.4 Association Acceptance Policy

The MMPS Application Entity does not accept Associations.

3.2.3 Storage Application Entity Specification

3.2.3.1 SOP Classes

IS provides Standard Conformance to the following SOP Classes as SCU:

Table 3-15 SOP Classes for AE Storage

| SOP Class Name | SOP Class UID |
|--|-------------------------------|
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 |
| Digital X-Ray Image Storage - For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 |
| Second Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 |
| X-Ray Radiation Dose SR | 1.2.840.10008.5.1.4.1.1.88.67 |

IS provides Standard Conformance to the following SOP Classes as SCP:

Table 3-16 SOP Classes for AE Storage

| SOP Class Name | SOP Class UID |
|--|-----------------------------|
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 |
| Second Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 |

3.2.3.2 Association Policies

3.2.3.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 3-17 DICOM Application Context

| | |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

3.2.3.2.2 Number of Associations

IS accepts one Association at a time for the Storage SCP AE.

3.2.3.2.3 Asynchronous Nature

IS does not support asynchronous communication (multiple outstanding transactions over a single Association).

3.2.3.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 3-18 DICOM Implementation Class and Version

| | |
|-----------------------------|----------------------|
| Implementation Class UID | 1.2.840.69677977.1.0 |
| Implementation Version Name | EPACS20001.0 |

3.2.3.3 Association Initiation Policy

3.2.3.3.1 Activity - Send Images

3.2.3.3.1.1 Description and Sequencing of Activities




Figure 5 Sequencing of Activity - Send Images

Copyright © SSME 2021. All rights reserved. Restricted.

A possible sequence of interactions between the Storage AE and a Storage Provider is illustrated in Figure 5:

1. The Storage AE opens an association with the Storage Provider
2. An acquired image is transmitted to the I Storage Provider using a C-STORE request and the Storage Provider replies with a C-STORE response (status success).
3. Another acquired DX image is transmitted to the Storage Provider using a C-STORE request and the Storage Provider replies with a C-STORE response (status success).
4. An N-ACTION request is transmitted to the Storage Provider to obtain storage commitment of previously transmitted DX images and RDSR. The Storage Provider replies with a N-ACTION response indicating the request has been received and is being processed.
5. The Storage Provider immediately transmits an N-EVENT-REPORT request notifying the Storage AE of the status of the Storage Commitment Request (sent in step 6 using the N-ACTION message). The Storage AE replies with a N-EVENT-REPORT response confirming receipt. The Storage Provider could send this message at any time or omit it entirely in favour of transmitting the N-EVENT-REPORT over a separate dedicated association (see note).
6. The Storage AE closes the association with the Storage Provider.

3.2.3.3.1.2 Proposed Presentation Contexts

IS is capable of proposing the Presentation Contexts shown in the following table:

Table 3-19 Proposed Presentation Contexts for Activity Send Images

| Presentation Context Table | | | | | |
|--|-------------------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Digital X-Ray Image Storage - For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Second Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| X-Ray Radiation Dose SR | 1.2.840.10008.5.1.4.1.1.88.67 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

3.2.3.3.1.3 SOP Specific Conformance Image & RDSR Storage SOP Classes

All Image & RDSR Storage SOP Classes supported by the Storage AE exhibit the same behavior, except where stated, and are described together in this section.

If Digital X-Ray Image Storage SOP Instances are included in the Send Job and a corresponding Presentation Context is not accepted, then the Association is aborted using AP-ABORT and the send job is marked as failed. The job failure is logged and reported to the user via the job control application.

Table 3-20 Storage C-STORE Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|------------------|------------|--|
| Success | Success | 0000 | The SCP has successfully stored the SOP Instance. If all SOP Instances in a send job have status success then the job is marked as complete. |
| Refused | Out of Resources | A700-A7FF | The Association is aborted and the send job |

Copyright © SSME 2021. All rights reserved. Restricted.

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|-----------------------------------|------------|--|
| | | | is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application. This is a transient failure. |
| Error | Data Set does not match SOP Class | A900-A9FF | The Association is aborted and the send job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application. |
| Error | Cannot Understand | C000-CFFF | The Association is aborted and the send job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application. |

The behavior of Storage AE during communication failure is summarized in the table below:

Table 3-21 Storage Communication Failure Behavior

| Exception | Behavior |
|--|--|
| Association aborted by the SCP or network layers | The send job is marked as failed. The reason is logged, and the job failure is reported to the user via the job control application. |

A failed send job can be restarted by user interaction. The system can be configured to automatically resend failed jobs if a transient status code is received. The delay between resending failed jobs and the number of retries is also configurable.

3.2.3.4 Association Acceptance Policy

IS DICOM application attempts to accept a new association for

- DIMSE C-ECHO
- DIMSE C-STORE

service operations. Any Information Objects transmitted on that association will be checked on conformance and stored in database if check was successful.

3.2.3.4.1.1 Proposed Presentation Contexts

IS DICOM application will accept Presentation Contexts as shown in the following table:

Table 3-22 Proposed Presentation Contexts for Receive Images

| Presentation Context Table | | | | | |
|--|-----------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

Copyright © SSME 2021. All rights reserved. Restricted.

| Presentation Context Table | | | | | |
|--|-----------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Digital X-Ray Image Storage - For Processing | 1.2.840.10008.5.1.4.1.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| Second Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |
| | | | | | |
| | | | | | |
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

3.2.3.4.1.2 SOP Specific Conformance Image Storage SOP Classes

Table 3-23 SOP Specific Conformance Image Storage SOP Classes

| Attribute Name | Tag | Accepted Values |
|----------------------------|-------------|------------------------------|
| Image pixel | | |
| Samples per Pixel | (0028,0002) | "1" |
| Photometric Interpretation | (0028,0004) | "MONOCHROME1", "MONOCHROME2" |
| Rows | (0028,0010) | up to 5000 |
| Columns | (0028,0011) | up to 5000 |
| Bits Allocated | (0028,0100) | 8,16 |
| Bits Stored | (0028,0101) | 8,10,12,16 |
| High Bit | (0028,0102) | Bits Stored - 1 |
| Pixel Representation | (0028,0103) | 0 or 1 |
| Overlay Plan | | |
| Overlay Rows | (6000,0010) | equal to 0028,0010 |
| Overlay Columns | (6000,0011) | equal to 0028,0011 |
| Overlay Type | (6000,0040) | "G" |
| Overlay Origin | (6000,0050) | "1\1" |

| Attribute Name | Tag | Accepted Values |
|------------------------|-------------|-------------------------|
| Overlay Bits Allocated | (6000,0100) | "1" or "16" |
| Overlay Bit Position | (6000,0102) | "0" or "12" |
| Overlay Data | (6000,3000) | Only if (6000,0100) = 1 |

If images contain different values, they shall not be imported. An error text shall be written into the network job list and logged too.

3.2.4 Storage Commitment Application Entity Specification

3.2.4.1 SOP Classes

IS provides Standard Conformance to the following SOP Classes as a SCU:

Table 3-24 SOP Classes for AE Storage Commitment

| SOP Class Name | SOP Class UID |
|-------------------------------|----------------------|
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 |

3.2.4.2 Association Policies

3.2.4.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 3-25 DICOM Application Context

| | |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

3.2.4.2.2 Number of Associations

The Storage Commitment AE will initiate only one DICOM association at a time to perform a DICOM storage commitment operation as a SCU to a Remote Host AE.

3.2.4.2.3 Asynchronous Nature

Asynchronous mode is not supported. All operations will be performed synchronously.

3.2.4.2.4 Implementation Identifying Information

Table 3-26 DICOM Implementation Class and Version

| | |
|-----------------------------|----------------------|
| Implementation Class UID | 1.2.840.69677977.1.0 |
| Implementation Version Name | EPACS20001.0 |

3.2.4.3 Association Initiation Policy

3.2.4.3.1 Activity: Send images successfully to a remote host declared as a Storage Commitment Provider

3.2.4.3.1.1 Description and Sequencing of Activities

The Storage AE will request storage commitment for instances of the Digital X-Ray Image Storage SOP Class and X-Ray Radiation Dose Report Storage SOP Class if the Remote AE is configured as an archive device and a presentation context for the Storage Commitment Push Model has been accepted.

When an association is established, an N-ACTION service is sent to request storage commitment for stored images and RDSR. The Storage Commitment AE releases the association immediately after receiving the N-ACTION

response from the remote storage provider.

3.2.4.3.1.2 Proposed Presentation Contexts

Table 3-27 Proposed Presentation Context Table

| Presentation Context Table - Proposed | | | | | |
|---------------------------------------|----------------------|---------------------------|-------------------|------|----------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended |
| Name | UID | Name List | UID List | | |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

3.2.4.3.1.3 SOP Specific DICOM Conformance Statement for the Storage Commitment Push Model SOP Class

The behavior of Storage Commitment AE when encountering status codes in a N-ACTION response is summarized in the table below:

Table 3-28 Storage Commitment N-ACTION Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|-----------------|------------|---|
| Success | Success | 0000 | The request for storage comment is considered successfully sent. A timer is started that will expire if no N-EVENT-REPORT for the Transaction UID is received within a configurable timeout period. |

The behavior of Storage AE during communication failure is summarized in the table below:

Table 3-29 Storage Commitment Communication Failure Behavior

| Exception | Behavior |
|--|---|
| Association aborted by the SCP or network layers | The send job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application. |

3.2.4.4 Association Acceptance Policy

3.2.4.4.1 Activity - Receive Storage Commitment Response

3.2.4.4.1.1 Description and Sequencing of Activities

The Storage commitment AE will accept associations in order to receive responses to a Storage Commitment Request. A possible sequence of interactions between the Storage Commitment AE and a Storage Provider (e.g., a storage or archive device supporting Storage Commitment SOP Classes as an SCP) is:

1. The Storage Provider opens a new association with the Storage AE.
2. The Storage Provider sends an N-EVENT-REPORT request notifying the Storage AE of the status of a previous Storage Commitment Request. The Storage AE replies with a N-EVENT-REPORT response confirming receipt.
3. The Storage Provider closes the association with the Storage AE.

3.2.4.4.1.2 Accepted Presentation Contexts

The Storage Commitment AE will accept Presentation Contexts as shown in the table below.

Table 3-30 Acceptable Presentation Contexts for Activity Receive Storage Commitment Response

| Presentation Context Table | | | |
|----------------------------|--|-----------------|---------------|
| Abstract Syntax | | Transfer Syntax | Role Extended |

| Name | UID | Name List | UID List | | Negotiation |
|-------------------------------|----------------------|---------------------------|---------------------|-----|-------------|
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCP | None |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | |
| | | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | | |

3.2.4.4.1.3 SOP Specific Conformance for Storage Commitment SOP Class

Upon receipt of a N-EVENT-REPORT the timer associated with the Transaction UID will be cancelled.

The behavior of Storage Commitment AE when receiving Event Types within the N-EVENT-REPORT is summarized in Table 3-31.

Table 3-31 Storage Commitment N-EVENT-REPORT Behavior

| Event Type Name | Event Type ID | Behaviour |
|--|---------------|---|
| Storage Commitment Request Successful | 1 | The Referenced SOP Instances under Referenced SOP Sequence (0008,1199) are marked within the database as "Stored & Committed (SC) " to the value of Retrieve AE Title (0008,0054). Successfully committed SOP Instances are candidates for automatic deletion from the local database if local resources become scarce. The conditions under which automatic deletion is initiated and the amount of space freed are site configurable. SOP Instances will not be deleted if they are marked with a lock flag. The least recently accessed SOP Instances are deleted first. |
| Storage Commitment Request Complete - Failures Exist | 2 | The Referenced SOP Instances under Referenced SOP Sequence (0008,1199) are treated in the same way as in the success case (Event Type 1). The Referenced SOP Instances under Failed SOP Sequence (0008,1198) are marked within the database as "Store & Commit Failed (Sf) ". The Failure Reasons are logged and the job failure is reported to the user via the job control application. A send job that failed storage commitment will not be automatically restarted but can be restarted by user interaction. |

The reasons for returning specific status codes in a N-EVENT-REPORT response are summarized in Table 3-32.

Table 3-32 Storage Commitment N-EVENT-REPORT Response Status Reasons

| Service Status | Further Meaning | Error Code | Reasons |
|----------------|------------------------|------------|--|
| Success | Success | 0000 | The storage commitment result has been successfully received. |
| Failure | Unrecognized Operation | 0211H | The Transaction UID in the N-EVENT-REPORT request is not recognized (was never issued within an N-ACTION request). |
| Failure | Resource Limitation | 0213H | The Transaction UID in the N-EVENT-REPORT request has expired (no N-EVENT-REPORT was received within a configurable time limit). |

| Service Status | Further Meaning | Error Code | Reasons |
|----------------|------------------------|------------|---|
| Failure | No Such Event Type | 0113H | An invalid Event Type ID was supplied in the N-EVENT-REPORT request. |
| Failure | Processing Failure | 0110H | An internal error occurred during processing of the N-EVENT-REPORT. A short description of the error will be returned in Error Comment (0000,0902). |
| Failure | Invalid Argument Value | 0115H | One or more SOP Instance UIDs with the Referenced SOP Sequence (0008,1199) or Failed SOP Sequence (0008,1198) was not included in the Storage Commitment Request associated with this Transaction UID. The unrecognized SOP Instance UIDs will be returned within the Event Information of the N-EVENT-REPORT response. |

3.2.4.4.1.4 SOP Specific Conformance for Verification SOP Class

The Storage AE provides standard conformance to the Verification SOP Class as an SCP. If the C-ECHO request was successfully received, a 0000 (Success) status code will be returned in the C-ECHO response. Otherwise, a C000 (Error - Cannot Understand) status code will be returned in the C-ECHO response.

3.2.5 Query Application Entity Specification

3.2.5.1 SOP Classes

IS provides Standard Conformance to the following SOP Classes as a SCU:

Table 3-33 SOP Classes for AE Query/Retrieve

| SOP Class Name | SOP Class UID |
|---|-----------------------------|
| Study Root Query Information Model-Find | 1.2.840.10008.5.1.4.1.2.2.1 |

3.2.5.2 Association Policies

3.2.5.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 3-34 DICOM Application Context

| | |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

3.2.5.2.2 Number of Associations

IS initiates one Association at a time for each query request being processed to a remote node.

3.2.5.2.3 Asynchronous Nature

Asynchronous mode is not supported. All operations will be performed synchronously.

3.2.5.2.4 Implementation Identifying Information

Table 3-35 DICOM Implementation Identifying Information

| | |
|-----------------------------|----------------------|
| Implementation Class UID | 1.2.840.69677977.1.0 |
| Implementation Version Name | EPACS20001.0 |

3.2.5.3 Association Initiation Policy

3.2.5.3.1 Real-World Activity: Query Request (FIND SCU)

3.2.5.3.1.1 Description and Sequencing of Activities

The associated Real-World activity is to fill out a query form with search data and pass it as a query request to the network application that issues one or more C-FIND requests over a previously built association. The remote SCP will respond with related data-entries that will be passed to a browser application. When data transfer is finished, the association is closed.

In case of error, the list is built up to the point where the error occurred. The user is informed about incomplete processing of the query results.

3.2.5.3.1.2 Proposed Presentation Contexts

Table 3-36 Proposed Presentation Contexts for Activity Worklist Update

| Presentation Context Table | | | | | |
|---|-----------------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Study Root Query Information Model-Find | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | | |

3.2.5.3.1.3 SOP Specific Conformance for Query

IS Query SCU supports following DICOM tags are used as return and matching keys:

Table 3-37 return and matching keys:

| Attribute Name | Tag | Matching | UI | Return Key |
|----------------------------------|-------------|----------|-----|------------|
| Patient's Name | (0010,0010) | | Yes | |
| Patient ID | (0010,0020) | | Yes | |
| Patient's Birth Date | (0010,0030) | | | |
| Patient's Sex | (0010,0040) | | | |
| Patient Comments | (0010,4000) | | | |
| Study Instance UID | (0020,000D) | | | Yes |
| Study ID | (0020,0010) | | | |
| Study Date | (0008,0020) | | Yes | Yes |
| Study Time | (0008,0030) | | | Yes |
| Accession Number | (0008,0050) | Wildcard | Yes | |
| Referring Physician's Name | (0008,0090) | | | |
| Study Description | (0008,1030) | | | Yes |
| Number of Study Related Series | (0028,1206) | | | Yes |
| Number of Study Related Instance | (0028,1208) | | | Yes |
| Series Instance UID | (0020,000E) | | | Yes |

| Attribute Name | Tag | Matching | UI | Return Key |
|------------------------------------|-------------|----------|----|------------|
| Series Number | (0020,0011) | | | Yes |
| Modality | (0008,0060) | | | Yes |
| Number of Series Related Instances | (0020,1209) | | | Yes |
| Series Description | (0008,103E) | | | Yes |
| Body Part Examined | (0018,0015) | | | Yes |

Table 3-38 Modality Worklist C-FIND Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behaviour |
|----------------|--|------------------------|--|
| Success | Matching is complete | 0000 | The SCP has completed the matches. Worklist items are available for display or further processing. |
| Refused | Out of Resources | A700 | The Association is aborted and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |
| Failed | Identifier does not match SOP Class | A900 | The Association is aborted and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |
| Failed | Unable to Process | C000 - CFFF | The Association is aborted and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any additional error information in the Response will be logged. |
| Cancel | Matching terminated due to Cancel request | FE00 | If the query was cancelled due to too many worklist items then the SCP has completed the matches. Worklist items are available for display or further processing. Otherwise, The Association is aborted and the worklist query is marked as failed. The status meaning is logged and reported to the user if an interactive query. |
| Pending | Matches are continuing | FF00 | The worklist item contained in the Identifier is collected for later display or further processing. |
| Pending | Matches are continuing - Warning that one or more Optional Keys were not supported | FF01 | The worklist item contained in the Identifier is collected for later display or further processing. The status meaning is logged only once for each C-FIND operation. |
| * | * | Any other status code. | The Association is aborted and the worklist is marked as failed. The status meaning is logged and reported to the user if an interactive query. Any |

Copyright © SSME 2021. All rights reserved. Restricted.

| Service Status | Further Meaning | Error Code | Behaviour |
|----------------|-----------------|------------|--|
| | | | additional error information in the Response will be logged. |

3.2.5.4 Association Acceptance Policy

The Workflow Application Entity does not accept Associations.

3.2.6 Print Application Entity Specification

3.2.6.1 SOP Classes

IS provides Standard Conformance to the following SOP Classes as SCU:

Table 3-39 SOP Classes for AE Print

| SOP Class Name | SOP Class UID |
|---------------------------------------|-----------------------|
| Basic Grayscale Print Management Meta | 1.2.840.10008.5.1.1.9 |

3.2.6.2 Association Policies

3.2.6.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 3-40 DICOM Application Context for AE Print

| | |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

3.2.6.2.2 Number of Associations

IS initiates one Association at a time for each configured Printer. Multiple Printers can be configured.

Table 3-41 Number of Associations Initiated for AE Print

| | |
|---|---------------------------------|
| Maximum number of simultaneous Associations | (number of configured Printers) |
|---|---------------------------------|

3.2.6.2.3 Asynchronous Nature

IS does not support asynchronous communication (multiple outstanding transactions over a single Association). Asynchronous receiving of Printer Status (N_EVENT_REPORT) is supported

3.2.6.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 3-42 DICOM Implementation Class and Version for AE Print

| | |
|-----------------------------|----------------------|
| Implementation Class UID | 1.2.840.69677977.1.0 |
| Implementation Version Name | EPACS20001.0 |

3.2.6.3 Association Initiation Policy

3.2.6.3.1 Activity - Print Images

3.2.6.3.1.1 Description and Sequencing of Activities

A user composes images onto film sheets and requests them to be sent to a specific Printer. The user can select the desired film format and number of copies. Each print-job is forwarded to the job queue and processed individually. The Printer AE is invoked by the job control interface that is responsible for processing network tasks. The job consists of data describing the images and graphics to be printed as well as the requested layout and other parameters. The film sheet is internally processed, converted to a STANDARD/1,1 page and then the page image is sent. If no association to the printer can be established, the print-job is switched to a failed state and the user informed.




Figure 6 Sequencing of Activity - Print Images

A typical sequence of DIMSE messages sent over an association between Printer AE and a Printer is illustrated in Figure 6:

1. Printer AE opens an association with the Printer.
2. N-CREATE on the Film Session SOP Class creates a Film Session.
3. N-CREATE on the Film Box SOP Class creates a Film Box linked to the Film Session. A single Image Box will be created as the result of this operation. (Printer AE only uses the format STANDARD\1,1)
4. N-SET on the Image Box SOP Class transfers the contents of the film sheet to the printer.
5. N-ACTION on the Film Box SOP Class instructs the printer to print the Film Box.
6. The printer prints the requested number of film sheets.
7. The Printer asynchronously reports its status via N-EVENT-REPORT notification (Printer SOP Class). The printer can send this message at any time. Printer AE does not require the N-EVENT-REPORT to be sent. Printer AE is capable of receiving an N-EVENT-REPORT notification at any time during an association. If the Printer reports a status of FAILURE, the print-job is switched to a failed state and the user informed.
8. N-DELETE on the Film Session SOP Class deletes the complete Film Session SOP Instance hierarchy.
9. Printer AE closes the association with the Printer.

Copyright © SSME 2021. All rights reserved. Restricted.

Status of the print-job is reported through the job control interface. Only one job will be active at a time for each separate Printer. If any Response from the remote Application contains a status other than Success or Warning, the Association is aborted and the related Job is switched to a failed state. It can be restarted any time by user interaction or, if configured, by automated retry.

3.2.6.3.1.2 Proposed Presentation Contexts

IS is capable of proposing the Presentation Contexts shown in the table below:

Table 3-43 Proposed Presentation Contexts for Activity Film Images

| Presentation Context Table | | | | | |
|---------------------------------------|-----------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Basic Grayscale Print Management Meta | 1.2.840.10008.5.1.1.9 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

3.2.6.3.1.3 Common SOP Specific Conformance for All Print SOP Classes

The general behavior of Print AE during communication failure is summarized in the table below. This behavior is common for all SOP Classes supported by Print AE.

Table 3-44 Print Communication Failure Behavior

| Exception | Behavior |
|--|--|
| Association aborted by the SCP or network layers | The print-job is marked as failed. The reason is logged and the job failure is reported to the user via the job control application. |

3.2.6.3.1.4 SOP Specific Conformance for the Film Session SOP Class

IS Print AE supports the following DIMSE operations for the Film Session SOP Class:

- N-CREATE
- N-DELETE

Details of the supported attributes and status handling behavior are described in the following subsections.

The attributes supplied in an N-CREATE Request are listed in the table below:

Table 3-45 Film Session SOP Class N-CREATE Request Attributes

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|------------------|-------------|----|--------------------------------|-------------------|--------|
| Number of Copies | (2000,0010) | IS | 1 .. 9 | ALWAYS | User |
| Print Priority | (2000,0020) | | MED | ALWAYS | User |
| Medium Type | (2000,0030) | CS | BLUE FILM, CLEAR FILM or PAPER | ALWAYS | User |
| Film Destination | (2000,0040) | CS | MAGAZINE or PROCESSOR | ALWAYS | User |

The behavior of Print AE when encountering status codes in a N-CREATE response is summarized in the table below:

Table 3-46 Film Session SOP Class N-CREATE Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|------------------------------|------------|--|
| Success | Success | 0000 | The SCP has completed the operation successfully. |
| Warning | Attribute Value Out of Range | 0116H | The N-CREATE operation is considered successful but the status meaning is logged. Additional information in the Response identifying the attributes out of range will be logged (i.e., Elements in the Modification List/Attribute List) |

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|----------------------|------------------------|--|
| Warning | Attribute List Error | 0107H | The N-CREATE operation is considered successful but the status meaning is logged. Additional information in the Response identifying the attributes will be logged (i.e., Elements in the Attribute Identifier List) |
| * | * | Any other status code. | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |

The behavior of Printer AE when encountering status codes in a N-DELETE response is summarized in the table below:

Table 3-47 Printer SOP Class N-DELETE Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|-----------------|------------------------|--|
| Success | Success | 0000 | The SCP has completed the operation successfully. |
| * | * | Any other status code. | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |

3.2.6.3.1.5 SOP Specific Conformance for the Film Box SOP Class

Print AE supports the following DIMSE operations for the Film Box SOP Class:

- N-CREATE
- N-ACTION

Details of the supported attributes and status handling behavior are described in the following subsections.

The attributes supplied in an N-CREATE Request are listed in the table below:

Table 3-48 Film Box SOP Class N-CREATE Request Attributes

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|--------------------------------------|-------------|----|---|-------------------|--------|
| Image Display Format | (2010,0010) | CS | STANDARD\1,1 | ALWAYS | Auto |
| Referenced Film Session Sequence | (2010,0500) | SQ | | ALWAYS | Auto |
| >Referenced SOP Class UID | (0008,1150) | UI | 1.2.840.10008.5.1.1.1 | ALWAYS | Auto |
| >Referenced SOP Instance UID | (0008,1155) | UI | From created Film Session SOP Instance | ALWAYS | Auto |
| Film Orientation | (2010,0040) | CS | PORTRAIT or LANDSCAPE | ALWAYS | User |
| Film Size ID | (2010,0050) | CS | 14INX17IN, 14INX14IN, 11INX14IN, 11INX11IN, 85INX11IN, 8INX10IN | ALWAYS | User |
| Referenced Presentation LUT Sequence | (2050,0500) | SQ | Only sent if Presentation LUT SOP Class has been negotiated. | ANAP | Auto |
| >Referenced SOP Class UID | (0008,1150) | UI | 1.2.840.10008.5.1.1.23 | ALWAYS | Auto |
| >Referenced SOP Instance UID | (0008,1155) | UI | From created Presentation LUT SOP Instance | ALWAYS | Auto |

The behavior of Printer AE when encountering status codes in a N-CREATE response is summarized in the table below:

Table 3-49 Film Box SOP Class N-CREATE Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|---|------------------------|---|
| Success | Success | 0000 | The SCP has completed the operation successfully. |
| Warning | Requested Min Density or Max Density outside of printer's operating range | B605H | The N-CREATE operation is considered successful but the status meaning is logged. |
| * | * | Any other status code. | The Association is aborted, and the print-job is marked as failed. The status meaning is logged and reported to the user. |

An N-ACTION Request is issued to instruct the Print SCP to print the contents of the Film Box. The Action Reply argument in an N-ACTION response is not evaluated.

The behavior of Printer AE when encountering status codes in a N-ACTION response is summarized in the table below:

Table 3-50 Film Box SOP Class N-ACTION Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|---|------------------------|--|
| Success | Success | 0000 | The SCP has completed the operation successfully. The film has been accepted for printing. |
| Warning | Film Box SOP Instance hierarchy does not contain Image Box SOP Instances (empty page) | B603H | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| Warning | Image size is larger than Image Box size. The image has been demagnified. | B604H | The N-ACTION operation is considered successful but the status meaning is logged. |
| Warning | Image size is larger than Image Box size. The image has been cropped to fit. | B609H | The N-ACTION operation is considered successful but the status meaning is logged. |
| Warning | Image size or Combined Print Image Size is larger than Image Box size. The image or combined Print Image has been decimated to fit. | B60AH | The N-ACTION operation is considered successful but the status meaning is logged. |
| Failure | Unable to create Print Job SOP Instance; print queue is full. | C602 | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| Failure | Image size is larger than Image Box size. | C603 | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| Failure | Combined Print Image Size is larger than Image Box size. | C613 | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| * | * | Any other status code. | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |

3.2.6.3.1.6 SOP Specific Conformance for the Image Box SOP Class

Copyright © SSME 2021. All rights reserved. Restricted.

Print AE supports the following DIMSE operations for the Image Box SOP Class:

- N-SET

Details of the supported attributes and status handling behavior are described in the following subsections.

The attributes supplied in an N-SET Request are listed in the table below:

Table 3-51 Image Box SOP Class N-SET Request Attributes

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|--------------------------------|-------------|----|-------------------------------|-------------------|--------|
| Image Position | (2020,0010) | US | 1 | ALWAYS | Auto |
| Basic Grayscale Image Sequence | (2020,0110) | SQ | | ALWAYS | Auto |
| >Samples Per Pixel | (0028,0002) | US | 1 | ALWAYS | Auto |
| >Photometric Interpretation | (0028,0004) | CS | MONOCHROME2 | ALWAYS | Auto |
| >Rows | (0028,0010) | US | Depends on film size | ALWAYS | Auto |
| >Columns | (0028,0011) | US | Depends on film size | ALWAYS | Auto |
| >Pixel Aspect Ratio | (0028,0034) | IS | 1\1 | ALWAYS | Auto |
| >Bits Allocated | (0028,0100) | US | 8 | ALWAYS | Auto |
| >Bits Stored | (0028,0101) | US | 8 | ALWAYS | Auto |
| >High Bit | (0028,0102) | US | 7 | ALWAYS | Auto |
| >Pixel Representation | (0028,0103) | US | 0 | ALWAYS | Auto |
| >Pixel Data | (7FE0,0010) | OB | Pixels of rendered film sheet | ALWAYS | Auto |
| Request Image Size | (2020,0030) | | | | Auto |

The behavior of Printer AE when encountering status codes in a N-SET response is summarized in the table below:

Table 3-52 Image Box SOP Class N-SET Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|---|------------|--|
| Success | Success | 0000 | The SCP has completed the operation successfully. Image successfully stored in Image Box. |
| Warning | Image size is larger than Image Box size. The image has been demagnified. | B604H | The N-SET operation is considered successful but the status meaning is logged. |
| Warning | Requested Min Density or Max Density outside of printer's operating range. | B605H | The N-SET operation is considered successful but the status meaning is logged. |
| Warning | Image size is larger than Image Box size. The image has been cropped to fit. | B609H | The N-SET operation is considered successful but the status meaning is logged. |
| Warning | Image size or Combined Print Image Size is larger than Image Box size. The image or combined Print Image has been decimated to fit. | B60AH | The N-SET operation is considered successful but the status meaning is logged. |
| Failure | Image size is larger than Image Box size. | C603 | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|--|------------------------|--|
| Failure | Insufficient memory in printer to store the image. | C605 | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| Failure | Combined Print Image Size is larger than Image Box size. | C613 | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |
| * | * | Any other status code. | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |

3.2.6.3.1.7 SOP Specific Conformance for the Printer SOP Class

Printer AE supports the following DIMSE operations and notifications for the Printer SOP Class:

- N-GET
- N-EVENT-REPORT

Details of the supported attributes and status handling behavior are described in the following subsections.

Printer AE uses the Printer SOP Class N-GET operation to obtain information about the current printer status. The attributes obtained via N-GET are listed in the table below:

Table 3-53 Printer SOP Class N-GET Request Attributes

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|--------------------------|-------------|----|---------------------|-------------------|---------|
| Printer Status | (2110,0010) | CS | Provided by Printer | ALWAYS | Printer |
| Printer Status Info | (2110,0020) | CS | Provided by Printer | ALWAYS | Printer |
| Printer Name | (2110,0030) | | | | |
| Manufacturer | (0008,0070) | | | | |
| Manufacturer Model Name | (0008,1090) | | | | |
| Device Serial Number | (0018,1000) | | | | |
| Software Versions | (0018,1020) | | | | |
| Date Of Last Calibration | (0018,1200) | | | | |
| Time Of Last Calibration | (0018,1201) | | | | |

The Printer Status information is evaluated as follows:

1. If Printer status (2110,0010) is NORMAL, the print-job continues to be printed.
2. If Printer status (2110,0010) is FAILURE, the print-job is marked as failed. The contents of Printer Status Info (2110,0020) is logged and reported to the user via the job →control application.
3. If Printer status (2110,0010) is WARNING, the print-job continues to be printed. The contents of Printer Status Info (2110,0020) is logged and reported to the user via the job → control application.

The behavior of Printer AE when encountering status codes in a N-GET response is summarized in the table below:

Table 3-54 Printer SOP Class N-GET Response Status Handling Behavior

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|-----------------|------------|--|
| Success | Success | 0000 | The request to get printer status information was success. |

| Service Status | Further Meaning | Error Code | Behavior |
|----------------|-----------------|------------------------|--|
| * | * | Any other status code. | The Association is aborted and the print-job is marked as failed. The status meaning is logged and reported to the user. |

Printer AE is capable of receiving an N-EVENT-REPORT request at any time during an association.

The behavior of Printer AE when receiving Event Types within the N-EVENT-REPORT is summarized in the table below:

Table 3-55 Printer SOP Class N-EVENT-REPORT Behavior

| Event Type Name | Event Type ID | Behavior |
|-----------------|---------------|--|
| Normal | 1 | The print-job continues to be printed. |
| Warning | 2 | The print-job continues to be printed. The contents of Printer Status Info (2110,0020) is logged and reported to the user via the job-control application. |
| Failure | 3 | The print-job is marked as failed. The contents of Printer Status Info (2110,0020) is logged and reported to the user via the job-control application. |
| * | * | An invalid Event Type ID will cause a status code of 0113H to be returned in a N-EVENT-REPORT response. |

The reasons for returning specific status codes in a N-EVENT-REPORT response are summarized in the table below:

Table 3-56 Printer SOP Class N-EVENT-REPORT Response Status Reasons

| Service Status | Further Meaning | Error Code | Reasons |
|----------------|--------------------|------------|---|
| Success | Success | 0000 | The notification event has been successfully received. |
| Failure | No Such Event Type | 0113H | An invalid Event Type ID was supplied in the N-EVENT-REPORT request. |
| Failure | Processing Failure | 0110H | An internal error occurred during processing of the N-EVENT-REPORT. A short description of the error will be returned in Error Comment (0000,0902). |

3.2.6.4 Association Acceptance Policy

The Print Application Entity does not accept Associations.

3.2.7 Verification Application Entity Specification

3.2.7.1 SOP Classes

IS provides Standard Conformance to the following SOP Classes as a SCU:

Table 3-57 SOP Classes for AE Verification

| SOP Class Name | SOP Class UID |
|----------------|-------------------|
| Verification | 1.2.840.10008.1.1 |

IS provides Standard Conformance to the following SOP Classes as a SCP:

Table 3-58 SOP Classes for AE Modality Verification

| SOP Class Name | SOP Class UID |
|----------------|-------------------|
| Verification | 1.2.840.10008.1.1 |

3.2.7.2 Association Establishment Policies

3.2.7.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 3-59 DICOM Application Context

| | |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

3.2.7.2.2 Number of Associations

IS initiates one Association at a time for each configured Printer. Multiple Printers can be configured.

3.2.7.2.3 Asynchronous Nature

IS does not support asynchronous communication (multiple outstanding transactions over a single Association).

3.2.7.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 3-60 DICOM Implementation Class and Version

| | |
|-----------------------------|----------------------|
| Implementation Class UID | 1.2.840.69677977.1.0 |
| Implementation Version Name | EPACS20001.0 |

3.2.7.3 Association Initiation Policy

3.2.7.3.1 Activity – Verification SCU

3.2.7.3.1.1 Description and Sequencing of Activities

The associated Real-World activity is a C-ECHO request initiated by Service User Interface environment whenever a verification of a destination configuration is requested. If an association to a remote Application Entity is successfully established, verification with the configured AET is requested via the open association. If the C-ECHO Response from the remote Application is received, a message will be displayed indicating that the C-Echo was successful, and the association is closed.

3.2.7.3.1.2 Proposed Presentation Contexts

IS is capable of proposing the Presentation Contexts shown in the following table:

Table 3-61 Proposed Presentation Contexts for Activity Verification

| Presentation Context Table | | | | | |
|----------------------------|-------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
| Name | UID | Name List | UID List | | |
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

3.3 Network Interfaces

3.3.1 Physical Network Interface

IS supports a single network interface. One of the following physical network interfaces will be available

Copyright © SSME 2021. All rights reserved. Restricted.

depending on installed hardware options:

Table 3-62 Supported Physical Network Interfaces

| |
|-------------------------------------|
| Ethernet 1000baseT/100baseT/10baseT |
|-------------------------------------|

3.3.2 Additional Protocols

EXAMPLE-INTEGRATED-MODLALITY conforms to the System Management Profiles listed in the table below. All requested transactions for the listed profiles and actors are supported. Support for optional transactions is listed in the table below:

Table 3-63 Supported System Management Profiles

| Profile Name | Actor | Protocols Used | Optional Transactions | Security Support |
|--|-------------|----------------|---------------------------|------------------|
| Network Address Management | DHCP Client | DHCP | N/A | |
| | DNS Client | DNS | N/A | |
| Time Synchronization | NTP Client | NTP | Find NTP Server | |
| | DHCP Client | DHCP | N/A | |
| DICOM Application Configuration Management | LDAP Client | LDAP | Client Update LDAP Server | See Section B.7 |

3.3.2.1 DHCP

DHCP can be used to obtain TCP/IP network configuration information. The network parameters obtainable via DHCP are shown in the table below. The Default Value column of the table shows the default used if the DHCP server does not provide a value. Values for network parameters set in the Service/Installation tool take precedence over values obtained from the DHCP server. Support for DHCP can be configured via the Service/Installation Tool. The Service/Installation tool can be used to configure the machine name. If DHCP is not in use, TCP/IP network configuration information can be manually configured via the Service/Installation Tool.

Table 3-64 Supported DHCP Parameters

| DHCP Parameter | Default Value |
|---------------------|--|
| IP Address | None |
| Hostname | Requested machine name |
| List of NTP servers | Empty list |
| List of DNS servers | Empty list |
| Routers | Empty list |
| Static routes | None |
| Domain name | None |
| Subnet mask | Derived from IP Address (see service manual) |
| Broadcast address | Derived from IP Address (see service manual) |
| Default router | None |
| Time offset | Site configurable (from Time zone) |
| MTU | Network Hardware Dependent |
| Auto-IP permission | No permission |

If the DHCP server refuses to renew a lease on the assigned IP address all active DICOM Associations will be aborted.

3.3.2.2 DNS

DNS can be used for address resolution. If DHCP is not in use or the DHCP server does not return any DNS server addresses, the identity of a DNS server can be configured via the Service/Installation Tool. If a DNS server is not in use, local mapping between hostname and IP address can be manually configured via the Service/Installation Tool.

3.3.2.3 NTP

The NTP client implements the optional Find NTP Server Transaction. The NTP client will issue an NTP broadcast to identify any local NTP servers. If no local servers can be found via NTP broadcast, the NTP Servers identified by DHCP will be used as time references. Additionally, one or more NTP Servers can be configured via the Service/Installation Tool. If no NTP Servers are identified then the local clock will be used as a time reference and a warning written to the system log files.

3.3.3 IPv4 and IPv6 Support

This product only supports IPv4 connections.

3.4 Configuration

3.4.1 AE Title/Presentation Address Mapping

3.4.1.1 Local AE Titles

All local applications use the AE Titles and TCP/IP Ports configured via the Service/Installation Tool. The Field Service Engineer can configure the TCP Port via the Service/Installation Tool. The AE Titles must be configured during installation. The local AE Title used by each individual application can be configured independently of the AE Title used by other local applications. If so configured, all local AEs are capable of using the same AE Title.

Table 3-65 AE Title Configuration Table

| Application Entity | Default AE Title | Default TCP/IP Port |
|------------------------|------------------|---------------------|
| Modality Worklist SCU | No Default | Not Applicable |
| MPPS SCU | No Default | Not Applicable |
| Storage SCU | No Default | Not Applicable |
| Storage SCP | No Default | 104 |
| Storage Commitment SCU | No Default | Not Applicable |
| Query SCU | No Default | Not Applicable |
| Print SCU | No Default | Not Applicable |
| Verification SCU | No Default | Not Applicable |
| Verification SCP | No Default | Not Applicable |

3.4.1.2 Remote AE Title/Presentation Address Mapping

The AE Title, host names and port numbers of remote applications are configured using the IS Service/Installation Tool.

3.4.1.2.1 Modality Worklist

The IS Service/Installation Tool must be used to set the AE Title, port-number, host-name and capabilities of the remote Modality Worklist SCP. Only a single remote Modality Worklist SCP can be defined.

3.4.1.2.2 Modality Performed Procedure Step

The IS Service/Installation Tool must be used to set the AE Title, port-number, host-name and capabilities of the remote MPPS SCP. Only a single remote MPPS SCP can be defined.

3.4.1.2.3 Storage

The IS Service/Installation Tool must be used to set the AE Titles, port-numbers, host-names and capabilities for the remote Storage SCPs. Associations will only be accepted from known AE Titles and associations from unknown AE Titles will be rejected (an AE Title is known if it can be selected within the Service/Installation Tool). Multiple remote Storage SCPs can be defined. Any Storage SCP can be configured to be an "Archive" device causing storage commitment to be requested for images or presentation states transmitted to the device.

3.4.1.2.4 Storage Commitment

If a remote AE is attached to a device containing a DICOM DeviceType attribute with value "ARCHIVE" it will be automatically configured as an "Archive" device provided the AE also supports Storage Commitment as an SCP.

3.4.1.2.5 Query

3.4.1.2.6 Print

The IS Service/Installation Tool must be used to set the AE Titles, port-numbers, host-names, IP→ Addresses and capabilities for the remote Print SCPs.

Multiple remote Print SCPs can be defined.

3.4.1.2.7 Verification

3.4.2 Parameters

A large number of parameters related to acquisition and general operation can be configured using the Service/Installation Tool. The table below only shows those configuration parameters relevant to DICOM communication. See the IS Service Manual for details on general configuration capabilities.

Table 3-66 Configuration Parameters Table

| Parameter | Configurable (Yes/No) | Default Value |
|--|-----------------------|--------------------------------------|
| General Parameters | | |
| Max PDU Receive Size | No | 65536 Bytes(64 kB) |
| Max PDU Send Size(larger PDUs will never be sent, even if the receiver supports a larger Max PDU Receive Size. If the receiver supports a smaller Max PDU Receive Size then the Max PDU Send Size will be reduced accordingly for the duration of the Association. Max PDU Receive Size information is exchanged during DICOM Association Negotiation in the Maximum Length Sub-Item of the A-ASSOCIATION-RQ and A-ASSOCIATE-AC) | Yes | 65536 Bytes(64 kB) |
| Storage Parameters | | |
| Number of times a failed send job may be retried | Yes | 0 (Failed send jobs are not retried) |
| Delay between retrying failed send jobs | Yes | 60 s |
| Maximum number of simultaneously initiated Associations by the Storage AE | Yes | 1 |
| Supported Transfer Syntaxes (separately configurable for each remote AE) | No | Implicit VR Little Endian |

| Parameter | Configurable (Yes/No) | Default Value |
|--|--------------------------|--------------------------------------|
| Storage Commitment Parameters | | |
| Timeout waiting for a Storage Commitment Notification (maximum duration of applicability for a Storage Commitment Transaction UID). | Yes | 24 hours |
| Maximum number of simultaneously accepted Associations by the Storage AE | Yes | 5 |
| Delay association release after sending a Storage Commitment Request (wait for a Storage Commitment Notification over the same association). | Yes | 120 s |
| Modality Worklist Parameters | | |
| Modality Worklist SCU time-out waiting for the final response to a C-FIND-RQ | Yes | 600 s |
| Maximum number of Worklist Items | Yes | 100 |
| Supported Transfer Syntaxes for Modality Worklist | No | Implicit VR Little Endian |
| Delay between automatic Worklist Updates | Yes | 10 mins |
| Query Worklist for specific Scheduled Station AE Title | Yes | EXINTMOD_WFL |
| Query Worklist for specific Modality Value | Yes | DX |
| MPPS Parameters | | |
| Supported Transfer Syntaxes for MPPS | No | Implicit VR Little Endian |
| Print Parameters | | |
| Supported Transfer Syntaxes (separately configurable for each remote printer) | No | Implicit VR Little Endian |
| Number of times a failed print-job may be retried | Yes | 0 (Failed send jobs are not retried) |
| Delay between retrying failed print-jobs | Yes | 60 s |
| Printer correction LUT (separately configurable for each remote printer) | Yes | Identity LUT |

4 Media Interchange

4.1 Implementation Model

4.1.1 Application Data Flow




Figure 7 Application Data Flow Diagram for Media Storage

The Offline-Media Application Entity exports images and Presentation States to DICOM Storage medium. It is associated with the local real-world activity "Export to ..." and "Import to Local Storage".

"Export to ..." is performed upon user request for selected patients, studies, series or instances (images or presentation states).

Import is supported from local media.

4.1.2 Functional Definition of AEs

Functional Definition of Offline-Media Application Entity

IS DICOM offline media storage application entity is capable of:

1. Creating a new File-set onto an unwritten medium (Export to...)
2. Updating an existing File-set by writing new SOP Instances onto the medium (Export to...)
3. Importing SOP Instances from the medium onto local storage

4.1.3 Sequencing of Real-World Activities

- At least one image or presentation state must exist and be selected before the Offline-Media Application Entity can be invoked.
- The operator can insert a new CD-R/DVD or plug in a USB media at any time before or after invocation of the Offline-Media Application Entity.
- The Offline-Media Application Entity will wait indefinitely for a media to be inserted before starting to write to the CD-R/DVD/USB device.
- If no media is available the export job can be cancelled from the job queue.

4.1.4 File Meta Information Options

The implementation information written to the File Meta Header in each file is:

Table 4-1 DICOM Implementation Class and Version for Media Storage

| | |
|-----------------------------|----------------------|
| Implementation Class UID | 1.2.840.69677977.1.0 |
| Implementation Version Name | EPACS20001.0 |

4.2 AE Specifications

4.2.1 Offline-Media Application Entity Specification

The Offline-Media Application Entity provides standard conformance to the Media Storage Service Class. The Application Profiles and roles are listed below:

Table 4-2 Application Profiles, Activities and Roles for Offline-Media

| Application Profiles Supported | Real World Activity | Role |
|--------------------------------|-------------------------------|---------|
| STD-GEN-CD | Export to local Archive Media | FSC/FCU |
| | Import into local Storage | FSR |
| STD-GEN-DVD-JPEG | Export to local Archive Media | FSC/FCU |
| | Import into local Storage | FSR |
| STD-GEN-USB-JPEG | Export to local Archive Media | FSC/FCU |
| | Import into local Storage | FSR |

4.2.1.1 File Meta Information for the Application Entity

The Source Application Entity Title included in the File Meta Header is configurable (see Section 4.4).

4.2.1.2 Real-World Activities

4.2.1.2.1 Activity - Export to ...

The Offline-Media Application Entity acts as an FSC when requested to export SOP Instances from the local database to a local medium.

A dialogue will be presented allowing the user to modify the suggested media label and provides control over the available media capacity. If the contents of the current selection do not fit on a single media an automatic separation into multiple export jobs will be suggested that can be adapted by the user.

The user will be prompted to prepare the local media for each export job. The contents of the export job will be written together with a corresponding DCOMDIR to a single-session CD-R. Writing in multi-session mode is not supported. The user can cancel an export job in the job queue.

The Offline-Media Application Entity supports the SOP Classes and Transfer Syntaxes listed in the table below:

Table 4-3 IODs, SOP Classes and Transfer Syntaxes for Offline Media

| Information Object Definition | SOP Class UID | Transfer Syntax | Transfer Syntax UID |
|--|-------------------------------|---------------------------|---------------------|
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Digital X-Ray Image Storage - For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Second Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |

| Information Object Definition | SOP Class UID | Transfer Syntax | Transfer Syntax UID |
|--|-------------------------------|---------------------------|---------------------|
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Digital X-Ray Image Storage - For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Digital X-Ray Image Storage - For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Second Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| X-Ray Radiation Dose SR | 1.2.840.10008.5.1.4.1.1.88.67 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |

4.3 Augmented and Private Application Profiles

IS does not support any augmented for private application profiles.

4.4 Media Configuration

All local applications use the AE Titles configured via the Service/Installation Tool. The Application Entity Titles configurable for Media Services are listed in the table below:

Table 4-4 AE Title Configuration Table

| Application Entity | Default AE Title |
|--------------------|------------------|
| Offline-Media | NA |

5 Support of Character Sets

IS DICOM applications support the Character Set (0008, 0005) value:

ISO_IR 100: Latin-1 Latin alphabet

ISO_IR 101: Latin-2 Eastern European

ISO_IR 109: Latin alphabet #3

ISO_IR 110: Latin alphabet #4

ISO_IR 144: Cyrillic

ISO_IR 127: Arabic

ISO_IR 126: Greek

ISO_IR 138: Hebrew

ISO_IR 148: Latin alphabet #5

ISO_IR 166: Thai

ISO 2022 IR 13\ISO 2022 IR 87: Japanese

ISO 2022 IR 6\ISO 2022 IR 149: Korean

ISO 2022 IR 6\ISO 2022 IR 58: Chinese

GB18030: Chinese

GBK: Chinese

ISO_IR 192: Unicode

If IS creates Module, the Character Set (0008,0005) value ISO_IR 100, ISO_IR 192 will be supported.

6 Security

IS does not support any specific security measures.

It is assumed that IS is used within a secured environment. It is assumed that a secured environment includes at a minimum:

- a. Firewall or router protections to ensure that only approved external hosts have network access to IS.
- b. Firewall or router protections to ensure that IS only has network access to approved external hosts and services.
- c. Any communication with external hosts and services outside the locally secured environment use appropriate secure network channels (e.g., such as a Virtual Private Network (VPN)).

Other network security procedures such as automated intrusion detection may be appropriate in some environments. Additional security features may be established by the local security policy and are beyond the scope of this conformance statement.

7 Annexes

7.1 IOD Contents

7.1.1 Created SOP Instances

The following tables use a number of abbreviations. The abbreviations used in the "Presence of ..." column are:

VNAP Value Not Always Present (attribute sent zero length if no value is present)

ANAP Attribute Not Always Present

ALWAYS Always Present

EMPTY Attribute is sent without a value

The abbreviations used in the "Source" column:

MWL the attribute value source Modality Worklist

USER the attribute value source is from User input

AUTO the attribute value is generated automatically

MPPS the attribute value is the same as that use for Modality Performed Procedure Step

CONFIG the attribute value source is a configurable parameter

Note:

All dates and times are encoded in the local configured calendar and time. Date, Time and Time zone are configured using the Service/Installation Tool.

7.1.1.1 X-Ray Digital Image IOD

Table 7-1 specifies the attributes of an X-Ray Digital Image transmitted by the IS storage application.

Table 7-2 specifies the attributes of a Grayscale Softcopy Presentation State transmitted by the IS storage application.

Table 7-1 IOD of Created DX SOP Instance

| IE | Module | Reference | Presence of Module |
|-----------|-------------------|-------------------|--------------------|
| Patient | Patient | Table 7-2 | ALWAYS |
| Study | General Study | Table 7-3 | ALWAYS |
| | Patient Study | Table 7-4 | ALWAYS |
| Series | General Series | Table 7-5 | ALWAYS |
| | DX Series | Table 7-16 | ALWAYS |
| Equipment | General Equipment | Table 7-6 | ALWAYS |
| Image | General Image | Table 7-7 | ALWAYS |
| | Image Pixel | Table 7-8 | ALWAYS |
| | X-Ray Acquisition | Table 7-11 | ALWAYS |
| | X-Ray Detector | Table 7-12 | ALWAYS |
| | DX Image | Table 7-13 | ALWAYS |
| | DX Detector | Table 7-14 | ALWAYS |
| | DX Position | Table 7-15 | ALWAYS |
| | DX Anatomy Imaged | Table 7-17 | ALWAYS |
| | X-Ray Dose | Table 7-18 | ALWAYS |

Copyright © SSME 2021. All rights reserved. Restricted.

| IE | Module | Reference | Presence of Module |
|----|------------------|-------------------|---|
| | X-Ray Filtration | Table 7-19 | ALWAYS |
| | X-Ray Grid | Table 7-20 | ALWAYS |
| | Modality LUT | Table 7-21 | Only if Pixel Intensity Relationship (0028,1040) is LOG |
| | SOP Common | Table 7-9 | ALWAYS |

7.1.1.2 Common Modules

Table 7-2 Patient Module of Created SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|-----------------------------|-------------|----|--|-------------------|----------|
| Patient's Name | (0010,0010) | PN | From Modality Worklist or user input. Values supplied via Modality Worklist will be entered as received. Values supplied via user input will contain all 5 components (some possibly empty).. Maximum 64 characters. | ALWAYS | MWL/USER |
| Patient ID | (0010,0020) | LO | From Modality Worklist or user input. Maximum 64 characters. | ALWAYS | MWL/USER |
| Patient's Birth Date | (0010,0030) | DA | From Modality Worklist or user input | ALWAYS | MWL/USER |
| Patient's Sex | (0010,0040) | CS | From Modality Worklist or user input | ALWAYS | MWL/USER |
| Referenced Patient Sequence | (0008,1120) | | From Modality Worklist | VNAP | MWL |
| Referenced SOP Class UID | (0008,1150) | | From Modality Worklist | VNAP | MWL |
| Referenced SOP Instance UID | (0008,1155) | | From Modality Worklist | VNAP | MWL |
| Other Patient IDs | (0010,1000) | | From Modality Worklist | EMPTY | |
| Other Patient Names | (0010,1001) | | From Modality Worklist | VNAP | MWL |
| Ethnic Group | (0010,2160) | | From Modality Worklist or user input. | VNAP | MWL/USER |
| Patient Comments | (0010,4000) | LT | From User Input. Maximum 1024 characters. | VNAP | MWL/USER |

Table 7-3 General Study Module of Created SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|----------------------------|-------------|----|---|-------------------|----------|
| Study Instance UID | (0020,000D) | UI | From Modality Worklist or generated by device | ALWAYS | MWL/AUTO |
| Study Date | (0008,0020) | DA | <yyyymmdd> | ALWAYS | MWL/AUTO |
| Study Time | (0008,0030) | TM | <hhmmss> | ALWAYS | MWL/AUTO |
| Referring Physician's Name | (0008,0090) | PN | From Modality Worklist or user input | VNAP | MWL/USER |
| Study ID | (0020,0010) | SH | Requested Procedure ID from | ALWAYS | MWL/USER |

Copyright © SSME 2021. All rights reserved. Restricted.

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|------------------------------|-------------|----|--|-------------------|----------|
| | | | Worklist or User Input | | |
| Accession Number | (0008,0050) | SH | From Modality Worklist or user input | ALWAYS | MWL/USER |
| Study Description | (0008,1030) | LO | Comment text box in study list. Maximum 1024 characters. | ALWAYS | MWL/AUTO |
| Referenced Study Sequence | (0008,1110) | SQ | From Modality Worklist | VNAP | MWL |
| >Referenced SOP Class UID | (0008,1150) | UI | From Modality Worklist | VNAP | MWL |
| >Referenced SOP Instance UID | (0008,1155) | UI | From Modality Worklist | VNAP | MWL |
| Procedure Code Sequence | (0008,1032) | | From Modality Worklist | VNAP | MWL |
| Code Value | (0008,0100) | | From Modality Worklist | VNAP | MWL |
| Coding Scheme Designator | (0008,0102) | | From Modality Worklist | VNAP | MWL |
| Coding Scheme Version | (0008,0103) | | From Modality Worklist | VNAP | MWL |
| Code Meaning | (0008,0104) | | From Modality Worklist | VNAP | MWL |

Table 7-4 Patient Study Module of Created SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|---------------------------------|-------------|----|--|-------------------|---------------|
| Admitting Diagnosis Description | (0008,1080) | LO | From Modality Worklist | VNAP | MWL |
| Patient's Age | (0010,1010) | AS | Calculated from DoB input on base of actual Date | ALWAYS | MWL/AUTO/USER |
| Patient's Size | (0010,1020) | | From Modality Worklist or user input | VNAP | MWL/USER |
| Patient's Weight | (0010,1030) | DS | From Modality Worklist or user input | VNAP | MWL/USER |
| Additional Patient History | (0010,21B0) | 3 | From Modality Worklist | VNAP | MWL |

Table 7-5 General Series Module of Created SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|-----------------------------|-------------|----|---|-------------------|----------|
| Modality | (0008,0060) | CS | DX | ALWAYS | CONFIG |
| Series Instance UID | (0020,000E) | UI | Generated by device | ALWAYS | AUTO |
| Series Number | (0020,0011) | IS | Generated by device | ALWAYS | AUTO |
| Series Date | (0008,0021) | DA | <yyyymmdd> | ALWAYS | AUTO |
| Series Time | (0008,0031) | TM | <hhmmss> | ALWAYS | AUTO |
| Performing Physician's Name | (0008,1050) | PN | Physician field in Study list. Maximum 64 characters. | VNAP | USER |
| Protocol Name | (0018,1030) | LO | Organ program | ALWAYS | MWL/AUTO |
| Series Description | (0008,103E) | LO | Organ from Study list. Maximum 512 characters. | VNAP | MWL/AUTO |
| Operator's Name | (0008,1070) | PN | Operator field in Study list. Maximum 64 characters. | VNAP | USER |

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|--|-------------|----|--|-------------------|-----------|
| Referenced Performed Procedure Step Sequence | (0008,1111) | SQ | Identifies the MPPS SOP Instance to which this image is related | VNAP | MPPS |
| >Referenced SOP Class UID | (0008,1150) | UI | MPPS SOP Class UID | ALWAYS | MPPS |
| >Referenced SOP Instance UID | (0008,1155) | UI | MPPS SOP Instance UID | ALWAYS | MPPS |
| Body Part Examined | (0018,0015) | CS | Value from Configuration (OGP EDIT) -zero length if not available | ALWAYS | AUTO |
| Request Attributes Sequence | (0040,0275) | SQ | Zero or 1 item will be present | VNAP | AUTO |
| >Requested Procedure ID | (0040,1001) | SH | From Modality Worklist | VNAP | MWL |
| >Scheduled Procedure Step ID | (0040,0009) | SH | From Modality Worklist | VNAP | MWL |
| >Scheduled Procedure Step Description | (0040,0007) | LO | From Modality Worklist | VNAP | MWL |
| >Scheduled Protocol Code Sequence | (0040,0008) | SQ | From Modality Worklist | VNAP | MWL |
| Code Value | (0008,0100) | SH | From Modality Worklist | VNAP | MWL |
| Coding Scheme Designator | (0008,0102) | SH | From Modality Worklist | VNAP | MWL |
| Coding Scheme Version | (0008,0103) | SH | From Modality Worklist | VNAP | MWL |
| Code Meaning | (0008,0104) | LO | From Modality Worklist | VNAP | MWL |
| Performed Procedure Step ID | (0040,0253) | SH | Same as MPPS. | ALWAYS | MPPS/AUTO |
| Performed Procedure Step Start Date | (0040,0244) | DA | Same as MPPS | ALWAYS | MPPS/AUTO |
| Performed Procedure Step Start Time | (0040,0245) | TM | Same as MPPS | ALWAYS | MPPS/AUTO |
| Performed Procedure Step Description | (0040,0254) | LO | Same as MPPS. Maximum 64 characters. | VNAP | MPPS/AUTO |
| Performed Protocol Code Sequence | (0040,0260) | SQ | Same as MPPS | ALWAYS | MPPS/AUTO |

Table 7-6 General Equipment Module of Created SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|---------------------------|-------------|----|------------------------------|-------------------|--------|
| Manufacturer | (0008,0070) | LO | From Configuration | ALWAYS | CONFIG |
| Institution Name | (0008,0080) | LO | From Configuration | ALWAYS | CONFIG |
| Station Name | (0008,1010) | SH | From Configuration/Host name | ALWAYS | CONFIG |
| Manufacturer's Model Name | (0008,1090) | LO | MULTIX Impact | ALWAYS | CONFIG |
| Device Serial Number | (0018,1000) | LO | From Configuration | ALWAYS | CONFIG |
| Software Version | (0018,1020) | LO | From Configuration | ALWAYS | CONFIG |
| Date of Last Calibration | (0018,1200) | DA | Generated by device | VNAP | AUTO |

Table 7-7 General Image Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|----------------|-----|----|-------|-------------------|--------|
|----------------|-----|----|-------|-------------------|--------|

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|---------------------------------|--|----|---|-------------------|----------|
| Instance Number | (0020,0013) | IS | Generated by software | ALWAYS | AUTO |
| Patient Orientation | (0020,0020) | CS | Zero length | VNAP | AUTO |
| Content Date | (0008,0023) | DA | <yyyymmdd> | ALWAYS | AUTO |
| Content Time | (0008,0033) | TM | <hhmmss> | ALWAYS | AUTO |
| Image Type | (0008,0008) | CS | Pixel Data Characteristics and Patient Examination Characteristics | ALWAYS | AUTO |
| Acquisition Number | (0020,0012) | IS | Generated by device | ALWAYS | AUTO |
| Acquisition Date | (0008,0022) | DA | Generated by device | ALWAYS | AUTO |
| Acquisition Time | (0008,0032) | TM | Generated by device | ALWAYS | AUTO |
| Image Comments | (0020,4000) | LT | From user input. Maximum 1024 characters. | VNAP | USER |
| Burned In Annotation | (0028,0301) | CS | NO | ALWAYS | AUTO |
| Lossy Image Compression | (0028,2110) | CS | 00 | ALWAYS | AUTO |
| Presentation LUT Shape | (2050,0020) | CS | IDENTIY | ALWAYS | AUTO |
| Irradiation Event UID | (0008,3010) | UI | UID of the irradiation event associated with the acquisition of the image | ALWAYS | AUTO |
| Rescale Intercept | (0028,1052) | DS | "0" | ALWAYS | AUTO |
| Rescale Slope | (0028,1053) | DS | "1" | ALWAYS | AUTO |
| Rescale Type | (0028,1054) | LO | "US" | ALWAYS | AUTO |
| Anatomic Region Sequence | (0008,2218) | SQ | Zero Length | VNAP | MWL/AUTO |
| > Include 'Code Sequence Macro' | Baseline Context ID is 4009 (see also Section B.8.6) | | | | |

Table 7-8 Image Pixel Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|----------------------------|-------------|----|---------------------------------|-------------------|--------|
| Samples per Pixel | (0028,0002) | US | "1" | ALWAYS | AUTO |
| Photometric Interpretation | (0028,0004) | CS | "MONOCHROME1", "MONOCHROME2" | ALWAYS | AUTO |
| Rows | (0028,0010) | US | Generated by device | ALWAYS | AUTO |
| Columns | (0028,0011) | US | Generated by device | ALWAYS | AUTO |
| Bits Allocated | (0028,0100) | US | Generated by device | ALWAYS | AUTO |
| Bits Stored | (0028,0101) | US | Generated by device | ALWAYS | AUTO |
| High Bit | (0028,0102) | US | Bits Stored-1 | ALWAYS | AUTO |
| Pixel Representation | (0028,0103) | US | 0 or 1 | ALWAYS | AUTO |

Table 7-9 SOP Common Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|------------------------|-------------|----|------------------------------|-------------------|--------|
| Specific Character Set | (0008,0005) | CS | "ISO_IR 100" or "ISO_IR 192" | ALWAYS | CONFIG |
| SOP Class UID | (0008,0016) | UI | 1.2.840.10008.5.1.4.1.1.12.2 | ALWAYS | AUTO |
| SOP Instance UID | (0008,0018) | UI | 1.3.12.2.1107.5.3..... | ALWAYS | AUTO |

Table 7-10 Softcopy Presentation LUT Module of SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|------------------------|-------------|----|----------|-------------------|--------|
| Presentation LUT Shape | (2050,0020) | CS | IDENTITY | ALWAYS | AUTO |

7.1.1.3 X-Ray Digital Image Modules

Table 7-11 X-Ray Acquisition Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|---|-------------|----|-----------------------------|-------------------|--------|
| KVP | (0018,0060) | DS | From Acquisition parameters | ALWAYS | AUTO |
| X-Ray Tube Current | (0018,1151) | IS | From Acquisition parameters | ALWAYS | AUTO |
| X-Ray Tube Current in μ A | (0018,8151) | DS | From Acquisition parameters | ALWAYS | AUTO |
| Exposure Time | (0018,1150) | IS | From Acquisition parameters | ALWAYS | AUTO |
| Exposure | (0018,1152) | IS | From Acquisition parameters | ALWAYS | AUTO |
| Exposure in μ As | (0018,1153) | IS | From Acquisition parameters | ALWAYS | AUTO |
| Image and Fluoroscopy Area Dose Product | (0018,115E) | DS | Generated from DAP | VNAP | AUTO |
| Relative X-Ray Exposure | (0018,1405) | IS | Generated from Calculation | ALWAYS | AUTO |
| Exposure Index | (0018,1411) | DS | Generated from Calculation | ALWAYS | AUTO |
| Target Exposure Index | (0018,1412) | DS | From Configuration | ALWAYS | CONFIG |
| Deviation Index | (0018,1413) | DS | Generated from Calculation | ALWAYS | AUTO |

Table 7-12 X-Ray Detector Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|-----------------------------------|-------------|----|----------------------------|-------------------|--------|
| Detector Type | (0018,7004) | CS | Generated by configuration | ALWAYS | CONFIG |
| Detector Description | (0018,7006) | LT | Generated by configuration | VNAP | CONFIG |
| Detector ID | (0018,700A) | SH | Generated by configuration | VNAP | CONFIG |
| Date of Last Detector Calibration | (0018,700C) | DA | Generated by software | VNAP | AUTO |
| Time of Last Detector Calibration | (0018,700E) | TM | Generated by software | VNAP | AUTO |
| Detector Conditions Nominal Flag | (0018,7000) | CS | Generated by software | VNAP | AUTO |
| Detector Temperature | (0018,7001) | DS | Generated by device | ALWAYS | AUTO |

Table 7-13 DX Image Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|----------------|-----|----|-------|-------------------|--------|
|----------------|-----|----|-------|-------------------|--------|

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|---|-------------|----|---------------------|-------------------|--------|
| Pixel Intensity Relationship | (0028,1040) | CS | "LIN" or "LOG" | ALWAYS | AUTO |
| Acquisition Device Processing Description | (0018,1400) | LO | Generated by device | VNAP | AUTO |
| Acquisition Device Processing Code | (0018,1401) | LO | Generated by device | VNAP | AUTO |

Table 7-14 DX Detector Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|-----------------------------------|-------------|----|---------------------------------|-------------------|--------|
| Detector Type | (0018,7004) | CS | Generated by configuration | VNAP | CONFIG |
| Detector Description | (0018,7006) | LT | Generated by configuration | VNAP | CONFIG |
| Detector ID | (0018,700A) | SH | Generated by configuration | VNAP | CONFIG |
| Date of Last Detector Calibration | (0018,700C) | DA | Generated by software | VNAP | AUTO |
| Time of Last Detector Calibration | (0018,700E) | TM | Generated by software | VNAP | AUTO |
| Detector Conditions Nominal Flag | (0018,7000) | CS | Generated by software | VNAP | AUTO |
| Detector Temperature | (0018,7001) | DS | Generated by device | ALWAYS | AUTO |
| Sensitivity | (0018,6000) | DS | Generated by Configuration | ALWAYS | CONFIG |
| Field of View Shape | (0018,1147) | CS | Generated by device | VNAP | AUTO |
| Field of View Dimension(s) | (0018,1149) | IS | Generated by device or software | ALWAYS | AUTO |
| Field of View Origin | (0018,7030) | DS | Generated by device or software | ALWAYS | AUTO |
| Field of View Rotation | (0018,7032) | DS | Generated by device or software | ALWAYS | AUTO |
| Field of View Horizontal Flip | (0018,7034) | CS | Generated by software | ALWAYS | AUTO |
| Imager Pixel Spacing | (0018,1164) | DS | Generated by device | ALWAYS | AUTO |

Table 7-15 DX Position Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|---|-------------|----|-----------------------|-------------------|--------|
| View Position | (0018,5101) | CS | Generated by software | ALWAYS | AUTO |
| Distance Source to Patient | (0018,1111) | DS | Generated by device | ALWAYS | AUTO |
| Distance Source to Detector | (0018,1110) | DS | Generated by device | ALWAYS | AUTO |
| Estimated Radiographic Magnification Factor | (0018,1114) | DS | Generated by software | VNAP | AUTO |
| Positioner Type | (0018,1508) | CS | Generated by device | "COLUMN" | AUTO |
| Column Angulation | (0018,1450) | DS | Generated by device | VNAP | AUTO |
| Table Type | (0018,113A) | CS | Generated by device | VNAP | AUTO |
| Table Angle | (0018,1138) | DS | Generated by device | VNAP | AUTO |

Table 7-16 DX Series Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|--------------------------|-------------|----|--------------------|-------------------|--------|
| Presentation Intent Type | (0008,0068) | CS | "FOR PRESENTATION" | ALWAYS | AUTO |

Table 7-17 DX Anatomy Imaged Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|--------------------------|-------------|----|----------------------------|-------------------|--------|
| Image Laterality | (0020,0062) | CS | Generated by configuration | VNAP | CONFIG |
| Anatomic Region Sequence | (0008,2218) | SQ | Zero Length | VNAP | AUTO |

Table 7-18 X-Ray Dose Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|---|-------------|----|-----------------------------|-------------------|--------|
| KVP | (0018,0060) | DS | From Acquisition parameters | ALWAYS | AUTO |
| X-Ray Tube Current | (0018,1151) | IS | From Acquisition parameters | ALWAYS | AUTO |
| X-Ray Tube Current in μ A | (0018,8151) | DS | From Acquisition parameters | ALWAYS | AUTO |
| Exposure Time | (0018,1150) | IS | From Acquisition parameters | ALWAYS | AUTO |
| Exposure Time in μ S | (0018,8150) | DS | From Acquisition parameters | ALWAYS | AUTO |
| Exposure | (0018,1152) | IS | Generated from Calculation | ALWAYS | AUTO |
| Exposure in μ As | (0018,1153) | IS | From Acquisition parameters | ALWAYS | AUTO |
| Image and Fluoroscopy Area Dose Product | (0018,115E) | DS | Generated from DAP | ALWAYS | AUTO |
| Relative X-Ray Exposure | (0018,1405) | IS | Generated from Calculation | ALWAYS | AUTO |

Table 7-19 X-Ray Filtration Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|----------------|-------------|----|---------------------|-------------------|--------|
| Filter Type | (0018,1160) | SH | Generated by device | VNAP | AUTO |

Table 7-20 X-Ray Grid Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|---------------------|-------------|----|---------------------|-------------------|--------|
| Grid | (0018,1166) | CS | Generated by device | VNAP | AUTO |
| Grid Focal Distance | (0018,704C) | DS | Generated by device | VNAP | AUTO |

Table 7-21 Modality LUT Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|-----------------------|-------------|----|------------------------------|-------------------|--------|
| Modality LUT Sequence | (0028,3000) | SQ | present if (0028,1040) = LOG | ANAP | AUTO |
| > LUT Descriptor | (0028,3002) | US | <1024,0,16> | ANAP | AUTO |
| > Modality LUT Type | (0028,3004) | LO | US | ANAP | AUTO |
| > LUT Data | (0028,3006) | US | LUT | ANAP | AUTO |

7.1.1.4 Computed Radiography Image IOD

Table 7-22 Computed Radiography Image IOD

| IE | Module | Reference | Presence of Module |
|---------|---------|------------------|--------------------|
| Patient | Patient | Table 7-2 | ALWAYS |

Copyright © SSME 2021. All rights reserved. Restricted.

| IE | Module | Reference | Presence of Module |
|-----------|-------------------|-------------------|--------------------|
| Study | General Study | Table 7-3 | ALWAYS |
| | Patient Study | Table 7-4 | ALWAYS |
| Series | General Series | Table 7-5 | ALWAYS |
| | CR Series | Table 7-23 | ALWAYS |
| Equipment | General Equipment | Table 7-6 | ALWAYS |
| Image | General Image | Table 7-7 | ALWAYS |
| | Image Pixel | Table 7-8 | ALWAYS |
| | CR Image | Table 7-24 | ALWAYS |
| | SOP Common | Table 7-9 | ALWAYS |

7.1.1.5 Computed Radiography Image Modules

Table 7-23 CR Series Module of Created SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|--------------------|-------------|----|--------------------|-------------------|--------|
| Body Part Examined | (0018,0015) | CS | From configuration | ALWAYS | CONFIG |
| View Position | (0018,5101) | CS | From configuration | ALWAYS | CONFIG |

Table 7-24 CR Image Module of Created SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|----------------------------|-------------|----|-----------------------|-------------------|--------|
| Photometric Interpretation | (0028,0004) | CS | Generated by software | ALWAYS | AUTO |

7.1.1.6 X-Ray Radiation Dose SR IOD

Table 7-25 X-Ray Radiation Dose SR IOD

| IE | Module | Reference | Presence of Module |
|-----------|--------------------|-------------------|--------------------|
| Study | General Study | Table 7-3 | ALWAYS |
| | Patient Study | Table 7-4 | ALWAYS |
| Series | General Series | Table 7-5 | ALWAYS |
| | SR Document Series | Table 7-26 | ALWAYS |
| Equipment | General Equipment | Table 7-6 | ALWAYS |
| Document | SOP Common | Table 7-9 | ALWAYS |

7.1.1.7 X-Ray Radiation Dose SR Modules

Table 7-26 SR Document Series Module of Created DX SOP Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|---------------------|-------------|----|---------------------|-------------------|--------|
| Modality | (0008,0060) | CS | SR | ALWAYS | CONFIG |
| Series Instance UID | (0020,000E) | UI | Generated by device | ALWAYS | AUTO |

| Attribute Name | Tag | VR | Value | Presence of Value | Source |
|--------------------|-------------|----|---|-------------------|----------|
| Series Number | (0020,0011) | IS | Generated by device | ALWAYS | AUTO |
| Series Date | (0008,0021) | DA | <yyyymmdd> | ALWAYS | AUTO |
| Series Time | (0008,0031) | TM | <hhmmss> | ALWAYS | AUTO |
| Series Description | (0008,103E) | LO | Organ from Study list. Maximum 512 characters. | VNAP | MWL/AUTO |

Table 7-27 Projection X-Ray Radiation Dose (TID 10001)

| NL | Rel with Parent | VT | Concept Name | Value |
|----|-----------------|-----------|---|--|
| | | CONTAINER | EV (113701, DCM, "X-Ray Radiation Dose Report") | |
| > | HAS CONCEPT MOD | CODE | EV (121058, DCM, "Procedure reported") | (113704, DCM, "Projection X-Ray") |
| >> | HAS CONCEPT MOD | CODE | EV (G-C0E8, SRT, "Has Intent") | (R-408C3, SRT, "Diagnostic Intent") |
| > | | INCLUDE | DTID 1002 "Observer Context" | See Table 7-28 |
| > | HAS OBS CONTEXT | CODE | EV (113705, DCM, "Scope of Accumulation") | (113014, DCM, "Study") |
| >> | HAS PROPERTIES | UIDREF | DCID 10001 "UID Types" | |
| > | CONTAINS | CODE | EV (113945, DCM, "X-Ray Detector Data Available") | (R-0038D, SRT, "Yes") |
| > | CONTAINS | CODE | EV (113943, DCM, "X-Ray Source Data Available") | (R-0038D, SRT, "Yes") |
| > | CONTAINS | CODE | EV (113944, DCM, "X-Ray Mechanical Data Available") | (R-0038D, SRT, "Yes") |
| > | CONTAINS | INCLUDE | DTID 10002 "Accumulated X-Ray Dose" | See Table 7-30 |
| > | CONTAINS | INCLUDE | DTID 10003 "Irradiation Event X-Ray Data" | See Table 7-33 |
| > | CONTAINS | TEXT | EV (121106, DCM, "Comment") | |
| > | CONTAINS | CODE | EV (113854, DCM, "Source of Dose Information") | DCID 10020 "Source of Projection X-Ray Dose Information" (113858, DCM, "MPPS Content") |

Table 7-28 Observer Context (TID 1002)

| NL | Rel with Parent | VT | Concept Name | Value |
|----|-----------------|---------|--|------------------------|
| | HAS OBS CONTEXT | CODE | EV (121005, DCM, "Observer Type") | (121007 DCM, "Device") |
| | HAS OBS CONTEXT | INCLUDE | DTID 1004 "Device Observer Identifying Attributes" | See Table 7-29 |

Table 7-29 Device Observer Identifying Attributes (TID 1004)

| NL | Rel with Parent | VT | Concept Name | Value |
|----|-----------------|--------|---|-----------------------------------|
| | | UIDREF | EV (121012, DCM, "Device Observer UID") | Instance Creator UID value |
| | | TEXT | EV (121013, DCM, "Device Observer | Defaults to value of Station Name |

| NL | Rel with Parent | VT | Concept Name | Value |
|----|-----------------|------|--|--|
| | | | <u>Name</u>) | (0008,1010) in General Equipment Module |
| | | TEXT | EV (<u>121014, DCM, "Device Observer Manufacturer"</u>) | Defaults to value of Manufacturer (0008,0070) in General Equipment Module |
| | | TEXT | EV (<u>121015, DCM, "Device Observer Model Name"</u>) | Defaults to value of Manufacturer's Model Name (0008,1090) in General Equipment Module |
| | | TEXT | EV (<u>121016, DCM, "Device Observer Serial Number"</u>) | Defaults to value of Device Serial Number (0018,1000) in General Equipment Module |

Table 7-30 Accumulated X-Ray Dose (TID 10002)

| NL | Rel with Parent | VT | Concept Name | Value Set Constraint |
|----|-----------------|-----------|---|---------------------------------------|
| | | CONTAINER | EV (113702, DCM, "Accumulated X-Ray Dose Data") | |
| > | HAS CONCEPT MOD | CODE | EV (113764, DCM, "Acquisition Plane") | \$Plane |
| > | CONTAINS | CONTAINER | EV (122505, DCM, "Calibration") | |
| >> | HAS CONCEPT MOD | CODE | EV (113794, DCM, "Dose Measurement Device") | DCID 10010 "Dose Measurement Devices" |
| >> | CONTAINS | DATETIME | EV (113723, DCM, "Calibration Date") | |
| >> | CONTAINS | NUM | EV (122322, DCM, "Calibration Factor") | UNITS = EV (1, UCUM, "no units") |
| >> | CONTAINS | NUM | EV (113763, DCM, "Calibration Uncertainty") | UNITS = EV (% , UCUM, "Percent") |
| >> | CONTAINS | TEXT | EV (113724, DCM, "Calibration Responsible Party") | |
| >> | CONTAINS | TEXT | EV (113720, DCM, "Calibration Protocol") | |

Table 7-31 Accumulated Total Projection Radiography Dose (TID 10007)

| NL | Rel with Parent | VT | Concept Name | Value Set Constraint |
|----|-----------------|------|--|--|
| | | NUM | EV (<u>113722, DCM, "Dose Area Product Total"</u>) | UNITS = EV (Gy.m2, UCUM, "Gy.m2") |
| | | NUM | EV (<u>113725, DCM, "Dose (RP) Total"</u>) | UNITS = EV (Gy, UCUM, "Gy") |
| | | NUM | EV (<u>113737, DCM, "Distance Source to Reference Point"</u>) | UNITS = EV (mm, UCUM, "mm") |
| | | NUM | EV (<u>113731, DCM, "Total Number of Radiographic Frames"</u>) | UNITS = EV (1, UCUM, "no units") |
| | | CODE | EV (<u>113780, DCM, "Reference Point Definition"</u>) | DCID 10025 "Radiation Dose Reference Points" |

Table 7-32 Accumulated Cassette-Based Projection Radiography Dose (TID 10006)

| NL | Rel with Parent | VT | Concept Name | Value Set Constraint |
|----|-----------------|----|--------------|----------------------|
|----|-----------------|----|--------------|----------------------|

| NL | Rel with Parent | VT | Concept Name | Value Set Constraint |
|----|-----------------|-----|---|----------------------------------|
| | | NUM | EV (113731, DCM, "Total Number of Radiographic Frames") | UNITS = EV (1, UCUM, "no units") |

Table 7-33 Irradiation Event X-Ray Data (TID 10003)

| NL | Rel with Parent | VT | Concept Name | Value Set Constraint |
|----|-----------------|-----------|---|---|
| | | CONTAINER | EV (113706, DCM, "Irradiation Event X-Ray Data") | |
| > | HAS CONCEPT MOD | CODE | EV (113764, DCM, "Acquisition Plane") | DCID 10003 "Equipment Plane Identification" |
| > | CONTAINS | UIDREF | EV (113769, DCM, "Irradiation Event UID") | |
| > | CONTAINS | DATETIME | DT (111526, DCM, "DateTime Started") | |
| > | CONTAINS | CODE | EV (113721, DCM, "Irradiation Event Type") | DCID 10002 "Irradiation Event Types" |
| > | CONTAINS | TEXT | EV (125203, DCM, "Acquisition Protocol") | |
| > | CONTAINS | CODE | EV (113745, DCM, "Patient Table Relationship") | DCID 21 "Patient Equipment Relationship" |
| > | CONTAINS | CODE | EV (113743, DCM, "Patient Orientation") | (F-10440, SRT, "erect") |
| >> | HAS CONCEPT MOD | CODE | EV (113744, DCM, "Patient Orientation Modifier") | (F-10320, SRT, "standing") |
| > | CONTAINS | CODE | EV (123014, DCM, "Target Region") | DCID 4031 "Common Anatomic Regions" |
| > | CONTAINS | NUM | EV (122130, DCM, "Dose Area Product") | UNITS = EV (Gy.m2, UCUM, "Gy.m2") |
| > | CONTAINS | NUM | EV (111638, DCM, "Patient Equivalent Thickness") | UNITS = EV (mm, UCUM, "mm") |
| > | CONTAINS | TEXT | EV (113780, DCM, "Reference Point Definition") | |
| > | CONTAINS | TEXT | EV (121106, DCM, "Comment") | |
| > | CONTAINS | INCLUDE | DTID 10003A "Irradiation Event X-Ray Detector Data" | |
| > | CONTAINS | INCLUDE | DTID 10003B "Irradiation Event X-Ray Source Data" | |
| > | CONTAINS | INCLUDE | DTID 10003C "Irradiation Event X-Ray Mechanical Data" | |

Table 7-34 Irradiation Event X-Ray Detector Data (TID 10003A)

| NL | Rel with Parent | VT | Concept Name | Value Set Constraint |
|----|-----------------|-----|---|----------------------------------|
| | | NUM | EV (113845, DCM, "Exposure Index") | UNITS = EV (1, UCUM, "no units") |
| | | NUM | EV (113846, DCM, "Target Exposure Index") | UNITS = EV (1, UCUM, "no units") |
| | | NUM | EV (113847, DCM, "Deviation Index") | UNITS = EV (1, UCUM, "no units") |

| NL | Rel with Parent | VT | Concept Name | Value Set Constraint |
|----|-----------------|---------|------------------------------------|--|
| | | INCLUDE | DTID 1021 "Device Participant" | \$DeviceProcedureRole = EV (113942, DCM, "X-Ray Reading Device") |
| | | IMAGE | EV (113795, DCM, "Acquired Image") | |

Table 7-35 Irradiation Event X-Ray Source Data (TID 10003B)

| NL | Rel with Parent | VT | Concept Name | Value Set Constraint |
|----|-----------------|-----------|--|--|
| | | NUM | EV (113738, DCM, "Dose (RP) ") | UNITS = EV (Gy, UCUM, "Gy") |
| | | TEXT | EV (113780, DCM, "Reference Point Definition") | |
| | | CODE | EV (113780, DCM, "Reference Point Definition") | DCID 10025 "Radiation Dose Reference Points" |
| | | NUM | EV (113742, DCM, "Irradiation Duration") | UNITS = EV (s, UCUM, "s") |
| | | NUM | EV (113733, DCM, "KVP") | UNITS = EV (kV, UCUM, "kV") |
| | | NUM | EV (113734, DCM, "X-Ray Tube Current") | UNITS = EV (mA, UCUM, "mA") |
| | | NUM | EV (113824, DCM, "Exposure Time") | UNITS = EV (ms, UCUM, "ms") |
| | | NUM | EV (113736, DCM, "Exposure") | UNITS = EV (uA.s, UCUM, "uA.s") |
| | | NUM | EV (113766, DCM, "Focal Spot Size") | UNITS = EV (mm, UCUM, "mm") |
| | | CONTAINER | EV (113771, DCM, "X-Ray Filters") | |
| > | CONTAINS | CODE | EV (113772, DCM, "X-Ray Filter Type") | DCID 10007 "X-Ray Filter Types" |
| | | NUM | EV (113790, DCM, "Collimated Field Area") | UNITS = EV (m2, UCUM, "m2") |
| | | NUM | EV (113788, DCM, "Collimated Field Height") | UNITS = EV (mm, UCUM, "mm") |
| | | NUM | EV (113789, DCM, "Collimated Field Width") | UNITS = EV (mm, UCUM, "mm") |
| | | INCLUDE | DTID 1021 "Device Participant" | \$DeviceProcedureRole = EV (113859, DCM, "Irradiating Device") |

Table 7-36 Irradiation Event X-Ray Mechanical Data (TID 10003C)

| NL | Rel with Parent | VT | Concept Name | Value Set Constraint |
|----|-----------------|-----|---------------------------------------|-------------------------------|
| | | NUM | EV (113770, DCM, "Column Angulation") | UNITS = EV (deg, UCUM, "deg") |

Table 7-37 Device Participant (TID 1021)

| NL | Rel with Parent | VT | Concept Name | Value |
|----|-----------------|------|--|--------------------------------------|
| | | CODE | EV (113876, DCM, "Device Role in Procedure") | (113859, DCM, "Irradiating Device ") |
| > | HAS PROPERTIES | TEXT | EV (113877, DCM, "Device Name") | |
| > | HAS PROPERTIES | TEXT | EV (113878, DCM, "Device Manufacturer") | |
| > | HAS PROPERTIES | TEXT | EV (113879, DCM, "Device Model Name") | |
| > | HAS PROPERTIES | TEXT | EV (113880, DCM, "Device Serial Number") | |

| NL | Rel with Parent | VT | Concept Name | Value |
|----|-----------------|--------|---|----------------------------|
| > | HAS PROPERTIES | UIDREF | EV (121012, DCM, "Device Observer UID") | Instance Creator UID Value |

7.1.2 Used Fields in Received IOD by Application

The IS storage application does not receive SOP Instances. The usage of attributes received via Modality Worklist is described in Section 3.2.1 .

7.1.3 Attribute Mapping

The relationships between attributes received via Modality Worklist, stored in acquired images and communicated via MPPS are summarized in Table 7-38. The format and conventions used in DICOM Table 7-38 are the same as the corresponding table in Section J.6 in PS3.17.

Table 7-38 Attribute Mapping Between Modality Worklist, Image and MPPS

| Modality Worklist | Image IOD | MPPS IOD |
|---------------------------------------|--|--|
| Patient Name | Patient Name | Patient Name |
| Patient ID | Patient ID | Patient ID |
| Patient's Birth Date | Patient's Birth Date | Patient's Birth Date |
| Patient's Sex | Patient's Sex | Patient's Sex |
| Patient's Weight | Patient's Weight | |
| Referring Physician's Name | Referring Physician's Name | |
| ---- | ---- | Scheduled Step Attributes Sequence |
| Study Instance UID | Study Instance UID | >Study Instance UID |
| Referenced Study Sequence | Referenced Study Sequence | >Referenced Study Sequence |
| Accession Number | Accession Number | >Accession Number |
| ---- | Request Attributes Sequence | ---- |
| Requested Procedure ID | >Requested Procedure ID | >Requested Procedure ID |
| Requested Procedure Description | | >Requested Procedure Description |
| Scheduled Procedure Step ID | >Scheduled Procedure Step ID | >Scheduled Procedure Step ID |
| Scheduled Procedure Step Description | >Scheduled Procedure Step Description | >Scheduled Procedure Step Description |
| Scheduled Protocol Code Sequence | >Scheduled Protocol Code Sequence | ---- |
| ---- | Performed Protocol Code Sequence | Performed Protocol Code Sequence |
| ---- | Study ID | Study ID |
| ---- | Performed Procedure Step ID | Performed Procedure Step ID |
| ---- | Performed Procedure Step Start Date | Performed Procedure Step Start Date |
| ---- | Performed Procedure Step Start Time | Performed Procedure Step Start Time |
| ---- | Performed Procedure Step Description | Performed Procedure Step Description |
| ---- | Comments on the Performed Procedure Step | Comments on the Performed Procedure Step |
| ---- | ---- | Performed Series Sequence |
| Scheduled Performing Physician's Name | Performing Physician's Name | >Performing Physician's Name |
| Requested Procedure Code | ---- | Procedure Code Sequence |

Copyright © SSME 2021. All rights reserved. Restricted.

| Modality Worklist | Image IOD | MPPS IOD |
|-------------------|-------------------------------------|------------------|
| Sequence | | |
| ---- | Referenced Study Component Sequence | ---- |
| ---- | >Referenced SOP Class UID | SOP Class UID |
| ---- | >Referenced SOP Instance UID | SOP Instance UID |
| ---- | Protocol Name | Protocol Name |

7.1.4 Coerced/Modified Fields

The Modality Worklist AE will truncate attribute values received in the response to a Modality Worklist Query if the value length is longer than the maximum length permitted by the attribute's VR.

7.2 Coded Terminology and Templates

The Workflow AE is capable of supporting arbitrary coding schemes for Procedure and Protocol Codes. The contents of Requested Procedure Code Sequence (0032,1064) and Scheduled Protocol Code Sequence (0040,0008) supplied in Worklist Items will be mapped to Image IOD and MPPS attributes as described in Table 7-38. During installation, a service technician will establish a mapping between the site-specific codes and the Protocol Names used internally to identify acquisition protocols.

The contents of Anatomic Region Sequence (0008,2218) in generated images will be filled with an anatomic code selected by the user from a catalog. The default catalog of anatomic codes corresponds to CID 4009 "DX Anatomy Imaged" but can be extended using the Service/Installation Tool.

The contents of Performed Procedure Step Discontinuation Reason Code Sequence (0040,0281) for a discontinued MPPS will be filled with a code selected by the user from a fixed list corresponding to CID 9300 "Procedure Discontinuation Reasons".

7.3 Grayscale Image Consistency

The high-resolution display monitor attached to the product can be calibrated according to the Grayscale Standard Display Function (GSDf). The Service/Installation Tool is used together with a luminance meter to measure the Characteristic Curve of the display system and the current ambient light. See the IS Service Manual for details on the calibration procedure and supported calibration hardware. The result of the calibration procedure is a Monitor Correction LUT that will be active within the display subsystem after a system reboot.

7.4 Standard Extended / Specialized / Private SOP Classes

No Specialized or Private SOP Classes are supported.

7.4.1 X-Ray Digital Image Storage SOP Class

The X-Ray Digital Image Storage SOP Class is extended to create a Standard Extended SOP Class by addition of standard and private attributes to the created SOP Instances as documented in Section 7.1 .

7.5 Private Transfer Syntaxes

No Private Transfer Syntaxes are supported.

8 List of Figures

| | |
|--|----|
| Figure 1 Application Data Flow Diagram | 11 |
| Figure 2 Sequencing Constraints | 13 |
| Figure 3 Sequencing of Activity - Worklist Update | 15 |
| Figure 4 Sequencing of Activity - Acquire Images | 21 |
| Figure 5 Sequencing of Activity - Send Images | 26 |
| Figure 7 Sequencing of Activity - Print Images..... | 37 |
| Figure 8 Application Data Flow Diagram for Media Storage | 49 |

9 List of Tables

| | |
|--|----|
| Table 1-1 Network Services | 4 |
| Table 1-2 Media Services | 4 |
| Table 3-1 SOP Classes for AE Modality Worklist | 14 |
| Table 3-2 DICOM Application Context | 14 |
| Table 3-3 DICOM Implementation Identifying Information | 14 |
| Table 3-4 Proposed Presentation Contexts for Activity Worklist Update | 15 |
| Table 3-5 Modality Worklist C-FIND Response Status Handling Behavior | 16 |
| Table 3-6 Modality Worklist Communication Failure Behavior | 17 |
| Table 3-7 Worklist Request Identifier | 17 |
| Table 3-8 SOP Classes for AE MPPS | 20 |
| Table 3-9 DICOM Application Context | 20 |
| Table 3-10 DICOM Implementation Class and Version | 20 |
| Table 3-11 Proposed Presentation Contexts for Real-World Activity Acquire Images | 21 |
| Table 3-12 MPPS N-CREATE / N-SET Response Status Handling Behavior | 22 |
| Table 3-13 MPPS Communication Failure Behavior | 22 |
| Table 3-14 MPPS N-CREATE / N-SET Request Identifier | 22 |
| Table 3-15 SOP Classes for AE Storage | 25 |
| Table 3-16 SOP Classes for AE Storage | 25 |
| Table 3-17 DICOM Application Context | 26 |
| Table 3-18 DICOM Implementation Class and Version | 26 |
| Table 3-19 Proposed Presentation Contexts for Activity Send Images | 27 |
| Table 3-20 Storage C-STORE Response Status Handling Behavior | 27 |
| Table 3-21 Storage Communication Failure Behavior | 28 |
| Table 3-22 Proposed Presentation Contexts for Receive Images | 28 |
| Table 3-23 SOP Specific Conformance Image Storage SOP Classes | 29 |
| Table 3-24 SOP Classes for AE Storage Commitment | 30 |
| Table 3-25 DICOM Application Context | 30 |
| Table 3-26 DICOM Implementation Class and Version | 30 |
| Table 3-27 Proposed Presentation Context Table | 31 |
| Table 3-28 Storage Commitment N-ACTION Response Status Handling Behavior | 31 |
| Table 3-29 Storage Commitment Communication Failure Behavior | 31 |
| Table 3-30 Acceptable Presentation Contexts for Activity Receive Storage Commitment Response | 31 |
| Table 3-31 Storage Commitment N-EVENT-REPORT Behavior | 32 |
| Table 3-32 Storage Commitment N-EVENT-REPORT Response Status Reasons | 32 |
| Table 3-33 SOP Classes for AE Query/Retrieve | 33 |
| Table 3-34 DICOM Application Context | 33 |
| Table 3-35 DICOM Implementation Identifying Information | 33 |
| Table 3-36 Proposed Presentation Contexts for Activity Worklist Update | 34 |
| Table 3-37 return and matching keys: | 34 |
| Table 3-38 Modality Worklist C-FIND Response Status Handling Behavior | 35 |

Copyright © SSME 2021. All rights reserved. Restricted.

| | |
|--|----|
| Table 3-41 SOP Classes for AE Print | 36 |
| Table 3-42 DICOM Application Context for AE Print | 36 |
| Table 3-43 Number of Associations Initiated for AE Print | 36 |
| Table 3-44 DICOM Implementation Class and Version for AE Print | 36 |
| Table 3-45 Proposed Presentation Contexts for Activity Film Images..... | 38 |
| Table 3-46 Print Communication Failure Behavior | 38 |
| Table 3-47 Film Session SOP Class N-CREATE Request Attributes..... | 38 |
| Table 3-48 Film Session SOP Class N-CREATE Response Status Handling Behavior | 38 |
| Table 3-49 Printer SOP Class N-DELETE Response Status Handling Behavior | 39 |
| Table 3-50 Film Box SOP Class N-CREATE Request Attributes | 39 |
| Table 3-51 Film Box SOP Class N-CREATE Response Status Handling Behavior | 39 |
| Table 3-52 Film Box SOP Class N-ACTION Response Status Handling Behavior..... | 40 |
| Table 3-53 Image Box SOP Class N-SET Request Attributes..... | 41 |
| Table 3-54 Image Box SOP Class N-SET Response Status Handling Behavior | 41 |
| Table 3-55 Printer SOP Class N-GET Request Attributes..... | 42 |
| Table 3-56 Printer SOP Class N-GET Response Status Handling Behavior | 42 |
| Table 3-57 Printer SOP Class N-EVENT-REPORT Behavior | 43 |
| Table 3-58 Printer SOP Class N-EVENT-REPORT Response Status Reasons | 43 |
| Table 3-59 SOP Classes for AE Verification | 43 |
| Table 3-60 SOP Classes for AE Modality Verification | 43 |
| Table 3-61 DICOM Application Context | 44 |
| Table 3-62 DICOM Implementation Class and Version | 44 |
| Table 3-63 Proposed Presentation Contexts for Activity Verification | 44 |
| Table 3-64 Supported Physical Network Interfaces | 45 |
| Table 3-65 Supported System Management Profiles..... | 45 |
| Table 3-66 Supported DHCP Parameters..... | 45 |
| Table 3-67 AE Title Configuration Table | 46 |
| Table 3-68 Configuration Parameters Table..... | 47 |
| Table 4-1 DICOM Implementation Class and Version for Media Storage | 50 |
| Table 4-2 Application Profiles, Activities and Roles for Offline-Media | 50 |
| Table 4-3 IODs, SOP Classes and Transfer Syntaxes for Offline Media..... | 50 |
| Table 4-4 AE Title Configuration Table | 51 |
| Table 7-1 IOD of Created DX SOP Instance | 54 |
| Table 7-2 Patient Module of Created SOP Instances..... | 55 |
| Table 7-3 General Study Module of Created SOP Instances | 55 |
| Table 7-4 Patient Study Module of Created SOP Instances | 56 |
| Table 7-5 General Series Module of Created SOP Instances..... | 56 |
| Table 7-6 General Equipment Module of Created SOP Instances..... | 57 |
| Table 7-7 General Image Module of Created DX SOP Instances..... | 57 |
| Table 7-8 Image Pixel Module of Created DX SOP Instances..... | 58 |
| Table 7-9 SOP Common Module of Created DX SOP Instances | 59 |
| Table 7-10 Softcopy Presentation LUT Module of SOP Instances | 59 |
| Table 7-11 X-Ray Acquisition Module of Created DX SOP Instances..... | 59 |

Copyright © SSME 2021. All rights reserved. Restricted.

| | |
|---|----|
| Table 7-12 X-Ray Detector Module of Created DX SOP Instances..... | 59 |
| Table 7-13 DX Image Module of Created DX SOP Instances | 59 |
| Table 7-14 DX Detector Module of Created DX SOP Instances..... | 60 |
| Table 7-15 DX Position Module of Created DX SOP Instances..... | 60 |
| Table 7-16 DX Series Module of Created DX SOP Instances..... | 60 |
| Table 7-17 DX Anatomy Imaged Module of Created DX SOP Instances..... | 61 |
| Table 7-18 X-Ray Dose Module of Created DX SOP Instances | 61 |
| Table 7-19 X-Ray Filtration Module of Created DX SOP Instances | 61 |
| Table 7-20 X-Ray Grid Module of Created DX SOP Instances | 61 |
| Table 7-21 Modality LUT Module of Created DX SOP Instances | 61 |
| Table 7-22 Computed Radiography Image IOD..... | 61 |
| Table 7-23 CR Series Module of Created SOP Instances..... | 62 |
| Table 7-24 CR Image Module of Created SOP Instances | 62 |
| Table 7-25 X-Ray Radiation Dose SR IOD..... | 62 |
| Table 7-26 SR Document Series Module of Created DX SOP Instances..... | 62 |
| Table 7-27 Projection X-Ray Radiation Dose (TID 10001) | 63 |
| Table 7-28 Observer Context (TID 1002) | 63 |
| Table 7-29 Device Observer Identifying Attributes (TID 1004) | 63 |
| Table 7-30 Accumulated X-Ray Dose (TID 10002)..... | 64 |
| Table 7-31 Accumulated Total Projection Radiography Dose (TID 10007) | 64 |
| Table 7-32 Accumulated Cassette-Based Projection Radiography Dose (TID 10006) | 64 |
| Table 7-33 Irradiation Event X-Ray Data (TID 10003) | 65 |
| Table 7-34 Irradiation Event X-Ray Detector Data (TID 10003A)..... | 65 |
| Table 7-35 Irradiation Event X-Ray Source Data (TID 10003B)..... | 66 |
| Table 7-36 Irradiation Event X-Ray Mechanical Data (TID 10003C)..... | 66 |
| Table 7-37 Device Participant (TID 1021) | 66 |
| Table 7-38 Attribute Mapping Between Modality Worklist, Image and MPPS..... | 67 |

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone +49 9131 84-0
siemens-healthineers.com

Legal Manufacturer

Siemens Shanghai Medical Equipment
Ltd.
278 Zhou Zhu Road
201318 Shanghai
P.R.China