

Thumb base fracture

Multitom Rax Real3D¹ clinical case
University Hospital Wuerzburg, Germany



¹ Option



Study ID 5aab560

Image reprocessed on *syngo.via* with cinematic VRT.
Cinematic VRT is recommended for communication, education,
and publication purposes and not intended for diagnostic reading.

Clinical background and indication for Multitom Rax Real3D¹ examination

Patient

Female | *1994 | BMI 18.8 kg/m²

Anamnesis

While walking backwards, a horse collided with the patient's right thumb. Immediate pain and swelling were reported by the patient.

Indication for Real3D¹ examination

Radiography revealed a dislocated fracture of the first metacarpal (arrow). Articular involvement was deemed questionable.



Lateral



AP

Conventional X-ray examination

Study ID 5aab938

The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

¹ Option

Multitom Rax Real3D¹

Settings



Study ID 5aab560

Settings for tableside scan with High Quality Protocol

<i>Tube voltage</i>	80.7 kV
<i>Current time product</i>	625 mAs
<i>Dose area product</i>	1440 $\mu\text{Gy}\cdot\text{m}^2$
<i>Calculated value for $\text{CTDI}_{\text{vol},32}$</i>	10.6 mGy
<i>Scan time</i>	12 sec
<i>Number of projections</i>	314

Reconstruction settings for sectional views

<i>Pixel size</i>	0.4 mm
<i>Reconstruction kernel</i>	sharp (equivalent to Br69)
<i>Slice thickness</i>	2 mm


The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

¹ Option

Multitom Rax Real3D¹

Diagnostic findings

The fracture pattern was easily assessable in Real 3D¹ images, allowing for detailed evaluation of fragment dislocation. Coronal images display ulnar-sided displacement of the distal metacarpal I fragment. Image quality is not impaired after cast immobilization. Sagittal reformatting demonstrated articular affliction of the metacarpotrapezial joint (arrows). The fracture pattern could further be visualized in VRT for surgical planning.



Sagittal view



VRT view



Coronal view

Study ID 5aab560

The products/features (mentioned herein) are not commercially available in all countries.

Their future availability cannot be guaranteed.

Cinematic VRT is recommended for communication, education, and publication purposes and not intended for diagnostic reading.

¹ Option



“Spatial resolution and assessment of bone microarchitecture in Real3D images is superior to multidetector CT examinations. Good soft tissue contrast allows for diagnosis of fracture-related findings.”¹

Jan-Peter Grunz, MD
University Hospital Wuerzburg, Germany

¹The statements by Siemens Healthineers customers described herein are based on results that were achieved in the customer’s unique setting. Since there is no “typical” hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.



The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

Results from case studies are not predictive of results in other cases. Results in other cases may vary.

Dr. Jan-Peter Grunz is employed by an institution that receives financial support from Siemens Healthineers for collaborations.