

# Study Protocol Book

for ARTIS systems in  
Interventional Neuroradiology

[siemens-healthineers.com/education/  
clinical-specialty-educational-resources](https://siemens-healthineers.com/education/clinical-specialty-educational-resources)





## Dear Reader,

In neurointerventional radiology, sharing clinical expertise is key to advancing patient care, procedural excellence, and innovation. When experienced specialists exchange insights, practical approaches, and decision-making strategies in complex cerebrovascular cases, the entire neurointerventional community benefits.

Our Study Protocols series has been developed to provide precisely this platform for peer-to-peer knowledge sharing — by experts, for experts. These protocols offer you practical guidance, workflow recommendations, and relevant real-world experience across a range of neurointerventional procedures and imaging strategies.

Our aim is to support your entire team by enhancing daily clinical workflows, strengthening procedural confidence, and developing neurointerventional practice. The protocols also serve as a valuable reference for routine and complex clinical scenarios alike.

We sincerely thank all the contributors who generously shared their protocols, expertise, and clinical experience.

If you'd like to contribute your own protocols or clinical insights, please contact us.

Best wishes,

A handwritten signature in black ink, reading "Angelika Hench". The script is fluid and cursive, with the first name and last name clearly distinguishable.

**Angelika Hench,**

Customer Excellence  
Advanced Therapies

# Thank you

We would like to express our sincere gratitude to the following physicians and hospitals, whose contributions made this Application Protocol Book possible:

**René Chapot, MD,**  
**Hannes Nordmeyer, MD, Ekin Celik, MD,**  
Interventional Neuroradiology,  
Alfried Krupp Hospital, Essen, Germany

**Wan-You Guo, MD,**  
Department of Radiology,  
Taipei Veterans General Hospital, Taiwan

**Kiffon Keigher, APN,**  
**Demetrius K. Lopes, MD,**  
Rush University Medical Center,  
Chicago, USA

**Riitta Rautio, MD,**  
Department of Interventional Radiology,  
TYSK, Turku University Hospital,  
Turku, Finland

**Yuichi Murayama, MD,**  
**Toshihiro Ishibashi, MD, Ichiro Yuki, MD,**  
Department of Neurosurgery,  
The Jikei University School of Medicine,  
Tokyo, Japan

**Andrey Sergeev, MD,**  
Department of Neurosurgery,  
The Federal State Budgetary Institution  
“Almazov National Medical Research Centre”  
of the Ministry of Health, Russia

**PD Marios N. Psychogios, MD,**  
**Ioannis Tsogkas, MD,**  
Department of Diagnostic and  
Interventional Neuroradiology,  
University Medicine Goettingen, Germany

**Naoki Kato, MD,**  
**Yuichi Murayama, MD,**  
Department of Neurosurgery,  
The Jikei University School of Medicine,  
Tokyo, Japan

**Rahul Kumar, MD**  
GS Neuroscience Clinic and  
Research Centre Pvt.Ltd., Patna, Bihar, India

**Tobias Struffert, MD**  
Department of Neurosurgery,  
Universitätsklinikum Giessen, Germany

**Prof. Jan Gralla, MD, Pascal Mordasini, MD,**  
Institute of Diagnostic and  
Interventional Neuroradiology,  
University Hospital -Inselspital- Bern,  
Switzerland

**Martin Skalej, MD,**  
Department of Neuroradiology,  
University Hospital Magdeburg, Germany

**David Niemann, MD,**  
Department of Neurological Surgery,  
University of Wisconsin Hospitals and Clinics,  
USA

# Content

## Stroke treatment

|                                                                                                      |    |
|------------------------------------------------------------------------------------------------------|----|
| Visualization of collateral vessels for stroke thrombectomy                                          | 8  |
| Thrombectomy in stroke using <i>syngo</i> DynaCT Sine Spin                                           | 12 |
| One-stop-shop management for stroke                                                                  | 14 |
| Large vessel occlusion of the left MCA                                                               | 18 |
| Mechanical thrombectomy of left MCA occlusion                                                        | 22 |
| Visualization of pontine hemorrhage in <i>syngo</i> DynaCT Micro                                     | 26 |
| One stop management for aneurysmal subarachnoid hemorrhage                                           | 30 |
| Flow diverter placement in a subarachnoid hemorrhage case supported by 'SuperDyna' imaging technique | 32 |

## Cerebral aneurysms

|                                                                     |    |
|---------------------------------------------------------------------|----|
| Pipeline flow diverting stent follow-up                             | 36 |
| Embolization of an C-ophthalmic artery aneurysm with coil and stent | 38 |
| Flow Diverter Stenting of Large Intracranial Aneurysm               | 42 |

## Cerebral AVM

|                                                                                         |    |
|-----------------------------------------------------------------------------------------|----|
| Large frontal AVM embolization using <i>syngo</i> 3D/3D Fusion technique                | 46 |
| High-flow dural fistula treatment decision supported by <i>syngo</i> iFlow              | 50 |
| Embolization of cerebral AVM                                                            | 56 |
| AVM treatment planning with <i>syngo</i> Dyna4D                                         | 58 |
| Understanding complex angio architectures of Dural AVF supported by <i>syngo</i> Dyna4D | 60 |

## Intravenous injected procedures

|                                                                            |    |
|----------------------------------------------------------------------------|----|
| Follow-up after bilateral carotid stent placement                          | 62 |
| Follow-up after stent-assisted aneurysm coiling                            | 64 |
| Follow-up of aneurysm clipping using <i>syngo</i> DynaCT with IV injection | 66 |
| Visualization of clipped aneurysm                                          | 68 |
| Follow-up of aneurysm treated with a flow diverter                         | 70 |

## Other neurovascular procedures

|                                                                                       |    |
|---------------------------------------------------------------------------------------|----|
| Transradial embolization of the middle meningeal artery for chronic subdural hematoma | 72 |
| <i>syngo</i> DynaCT of petrous bones to visualize cholesteatoma                       | 76 |
| Total ossicular replacement prosthesis                                                | 80 |

Scan it to find more protocols!



---

Clinical  
Specialty  
Educational  
Resources

---



A decorative graphic on the left side of the page, consisting of a vertical line of black circles connected by a grey, wavy line, resembling a spiral binding.

# **Neuro interventions**

# Visualization of collateral vessels for stroke thrombectomy

**Courtesy of** J. Gralla, MD, P. Mordasini, MD, Institute of Diagnostic and Interventional Neuroradiology  
University Hospital – Inselspital – Bern, Switzerland

**Supported by** • syngo DynaCT Multiphase

**System & Software** ARTIS icono VE20, syngo Application Software VE2

## Patient history

56-year-old patient with a wake-up stroke presenting with hemiparesis on the right side and dysarthria, NIHSS 11.

## Diagnosis

M1 occlusion on the left side on the basis of an underlying intracranial atherosclerotic stenosis.

## Procedure description


Mechanical thrombectomy and stent implantation.

## Tips and tricks

syngo DynaCT Multiphase consists of 10 consecutive rotations and results in 10 volumes. Start injection with start of acquisition run. First two volumes required as mask run for potential reconstruction of perfusion maps.\* The other 8 volumes are fill runs.


\* Third party software required.

Graphical representation shows workflow including two mask runs




| Acquisition protocol              | 60sDCT Multiphase                                    |
|-----------------------------------|------------------------------------------------------|
| <b>Injection protocol</b>         |                                                      |
| Catheter position                 | antecubital 18G venous line                          |
| Contrast medium                   | 300 mg/mL iodine                                     |
| Dilution                          | None                                                 |
| Injection volume                  | 50 mL contrast media followed by 40 mL saline chaser |
| Injection rate                    | 5 mL/s CM and saline                                 |
| Duration of injection             | 18 s                                                 |
| X-ray delay                       | 0 s                                                  |
| Power injector used               | Yes                                                  |
| <b>Reconstructions</b>            |                                                      |
| <b>For both mask and fill run</b> |                                                      |
| Name                              | Neuro Perfusion Full HU Auto                         |
| VOI size                          | Full                                                 |
| Slice matrix                      | 512 × 512                                            |
| Kernel type                       | HU                                                   |
| Image characteristics             | Auto                                                 |
| Reconstruction mode               | Nat Fill                                             |
| Viewing preset                    | DynaCT Head                                          |


Initial *syngo* DynaCT Multiphase run before thrombectomy



Fill run 1–8, 10 mm MIP slices




Visualization of M1 occlusion. Fill run 4, MIP 10 mm




Visualization of potential collaterals in later phase. Fill run 6, MIP 10 mm


### Post-thrombectomy and stent placement




MIP 10 mm transversal



MIP 25 mm sagittal



MIP 25 mm coronal



MIP 25 mm transversal

# Thrombectomy in stroke using *syngo* DynaCT Sine Spin

**Courtesy of** René Chapot, MD; Ekin Celik, MD; Interventional Neuroradiology, Alfried Krupp Hospital, Essen, Germany

**Supported by** *syngo* DynaCT Sine Spin

**System & Software** ARTIS icono VE2 with *syngo*, Application Software VE2

## Patient history

83-year-old male patient

## Diagnosis

CTA in another hospital showed an acute occlusion of the left MCA (middle cerebral artery) in the M1 segment.

## Treatment

Intracranial thrombectomy with stent retriever under general anesthesia. Thrombectomy in several maneuvers with Solitaire X stent retriever and Tigertriever to reopen the vessel.

## General comments




Upon arrival of the patient in our department from the other hospital,

the initial CT was 3h old. With a native *syngo* DynaCT Sine Spin, we could check for potential bleeding and demarcation of the infarct before starting thrombectomy. The final *syngo* DynaCT Sine Spin showed a blood-brain barrier disorder, which influences the further treatment of the patient. Severe bleeding would have also been visible.

## Tips and tricks

The correct location of the intubation tube is very important for good image quality. Position the ventilation hose from the mouth towards the foot end.


## Pre-interventional *syngo* DynaCT Sine Spin




MPR 0.5 mm

|                       |                                                        |
|-----------------------|--------------------------------------------------------|
| Acquisition protocol  | 7sDCT Sine Spin - for pre- and post-interventional run |
| Injection protocol    |                                                        |
| Contrast medium (CM)  | N/A                                                    |
| Reconstructions       |                                                        |
| Name                  | DCT Head Clear Nat Fill                                |
| VOI size              | Full                                                   |
| Slice matrix          | 512 × 512                                              |
| Kernel type           | HU                                                     |
| Image characteristics | Smooth                                                 |
| Reconstruction mode   | Nat Fill                                               |
| Viewing preset        | DynaCT Head                                            |


DSA imaging during intervention






MCA occlusion



After first pass



Post-interventional syngo DynaCT Sine Spin



MPR 0.5 mm contrast medium pooling

# One-stop-shop management for stroke

---

|                              |                                                                                                                                                         |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Courtesy of</b>           | Marios N. Psychogios, MD, Ioannis Tsogkas, MD,<br>Department of Diagnostic and Interventional<br>Neuroradiology University Medicine Goettingen, Germany |
| <b>Supported by</b>          | <i>syngo</i> DynaCT                                                                                                                                     |
| <b>System &amp; Software</b> | Artis Q biplane VD11, <i>syngo</i> X Workplace VD10                                                                                                     |

---

## Patient history

56-year-old male. Right hemiparesis and aphasia with 75-minute onset. No intubation or sedation of the patient.

## Diagnosis

Total occlusion of arteria cerebri media.

## Treatment

Thrombectomy with SAVE technique.

## General comments

Direct admission of the patient to the angio suite for diagnostics and therapy (one-stop-shop management). Check for potential bleeding with the help of *syngo* DynaCT. Carotid bifurcation visible. A native FDCT scan was performed within 15 minutes of the patient arriving at the hospital. After IV injection of contrast media, a biphasic FDCTA scan was performed.





We timed the first phase of the FDCTA scan after a bolus-tracked digital subtraction angiography (the contrast agent must be visible in the intracranial internal carotid artery), while the second phase was acquired automatically with a delay of 5 s. The first phase shows the occlusion of the arteria cerebri media. The second phase of the FDCTA scan was used to additionally evaluate late collaterals. An occlusion or relevant stenosis of the proximal internal carotid artery could be excluded with the same acquisition.

## Tips and tricks

Very good stabilization of patient's head possible when a head holder is used, even for non-anesthetized patients.

Protocol – Part 1


|                       |                         |
|-----------------------|-------------------------|
| Acquisition protocol  | 20sDCT Head Clear 109kV |
| Injection protocol    |                         |
| Contrast medium (CM)  | N/A                     |
| Reconstructions       |                         |
| Name                  | DCT Head Clear          |
| VOI size              | Full                    |
| Slice matrix          | 512 × 512               |
| Kernel type           | HU                      |
| Image characteristics | Normal                  |
| Reconstruction mode   | Nat Fill                |
| Viewing preset        | DCT Head                |




DCT native, check for severe bleeding. 1–4 Thin MPRs from skull base to skull cap

## Protocol – Part 2

| Acquisition protocol      | 10sDCT DSA Head 70kV                                                                                                                              |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Injection protocol</b> |                                                                                                                                                   |
| Catheter position         | Intravenous injection                                                                                                                             |
| Contrast medium (CM)      | 400 mg/mL iodine                                                                                                                                  |
| Dilution (CM/Saline)      | No                                                                                                                                                |
| Injection volume          | 60 mL contrast media followed by 60 mL saline chaser                                                                                              |
| Injection rate            | 5 mL/s                                                                                                                                            |
| Duration of injection     | 12 s for CM + 12 s for saline                                                                                                                     |
| X-ray delay               | Bolus watching - start first run when contrast agent is seen in the intracranial internal carotid artery. Second run with automatic delay of 5 s. |
| Power injector used       | Yes                                                                                                                                               |
| <b>Reconstructions</b>    |                                                                                                                                                   |
| <b>Primary</b>            |                                                                                                                                                   |
| Name                      | DCT Head Clear for 1st and 2nd rotation                                                                                                           |
| VOI size                  | Full                                                                                                                                              |
| Slice matrix              | 512 × 512                                                                                                                                         |
| Kernel type               | HU                                                                                                                                                |
| Image characteristics     | Normal                                                                                                                                            |
| Reconstruction mode       | Nat Fill                                                                                                                                          |
| Viewing preset            | DCT Head                                                                                                                                          |




DCT arterial phase. MCA left occlusion.



DCT delayed phase. Collaterals visible.

## Biphasic Flat Detector CT Angiography Protocol



# Large vessel occlusion of the left MCA

---

|                              |                                                                         |
|------------------------------|-------------------------------------------------------------------------|
| <b>Courtesy of</b>           | Rahul Kumar, MD<br>GS Neuroscience Clinic, Patna, Bihar, India          |
| <b>Supported by</b>          | <i>syngo</i> Dyna3D, <i>syngo</i> DualVolume                            |
| <b>System &amp; Software</b> | ARTIS icono biplane VE20 with<br><i>syngo</i> Application Software VE20 |

---

## Patient history

A 48-year-old male presented with transient weakness (complete hemiparesis). Stroke onset 13 hours ago.

## Diagnosis

CT images showed a large vessel occlusion in the left middle cerebral artery (MCA).

## Treatment

The patient was immediately transferred to Angio for a thrombectomy procedure. The occlusion of the left MCA was confirmed in DSA imaging. After acquiring *syngo* Dyna3D runs from both sides, it was found that the MCA occlusion was resolved by collateral filling. The patient's stroke

symptoms improved immediately. It was decided not to proceed with the thrombectomy procedure.

## General comments



*syngo* Dyna3D/3D Fusion can achieve a good overview of the 3D volumes injected from either side.

## Tips and tricks

It is essential to use an X-ray delay to ensure constant filling of the complete vessel structure. Ensure that the patient's head is properly fixed. If the patient is not sedated, inform them about the procedure steps to ensure maximum cooperation and avoid movement during the 3D runs.

| Acquisition protocol  | 4s DSA Dyna3D Head for both sides |                       |
|-----------------------|-----------------------------------|-----------------------|
| Injection protocol    |                                   |                       |
| Catheter position     | Common carotid artery             |                       |
| Contrast medium (CM)  | 300 mg/mL                         |                       |
| Dilution              | 100%                              |                       |
| Injection volume      | 18 mL                             |                       |
| Injection rate        | 3 mL/s                            |                       |
| Duration of injection | 6 s                               |                       |
| X-ray delay           | 2 s                               |                       |
| Power injector used   | Yes                               |                       |
|                       |                                   |                       |
| Reconstructions       | Primary                           | Secondary             |
| Name                  | 3D Head Sub Full                  | 3D Head Nat Fill Full |
| VOI size              | Full                              | Full                  |
| Slice matrix          | 512x512                           | 512x512               |
| Kernel type           | EE                                | HU                    |
| Image characteristics | Auto                              | Auto                  |
| Reconstruction mode   | Sub                               | Nat                   |
| Viewing preset        | 3D Angio                          | 3D Angio              |

DSA




Completely occluded left MCA.


MIP




Injection in right common carotid artery.



Injection in left common carotid artery.



## MIP



Demonstration of the leptomeningeal collaterals.

## VRT



*syngo* DualVolume visualization after  
*syngo* 3D/3D Fusion of both 3D runs.

# Mechanical thrombectomy of left MCA occlusion

---

**Courtesy of** Prof. Wan-You Guo, MD, Department of Radiology, Taipei Veterans General Hospital, Taiwan

**Supported by**

- *syngo* DynaCT
- *syngo* Neuro PBV IR
- *syngo* iFlow

**System & Software** Artis zee biplane VC14, *syngo* X Workplace VB15

---

## Patient history

A 67-year-old man presented with symptoms of acute right hemiplegia of less than 6 hours duration since onset. The patient was taken to CT immediately to exclude hemorrhage with native CT.

## Procedure description

Left middle cerebral artery (MCA) occlusion. Noncontrast CT demonstrated a hyperdense MCA sign and excluded intracranial hemorrhage. Flat-detector CT (FDCT) angiography demonstrated the total occlusion of left MCA. A parenchymal cerebral blood volume (*syngo* Neuro PBV IR) map depicted a large area of hypoperfusion in the corresponding left MCA territory, which was similar to the results of multidetector CT perfusion imaging (not shown). After intraarterial mechanical thrombectomy, recanalization of the left MCA was

demonstrated by FDCT angiography. A *syngo* Neuro PBV IR map depicted the recovery of CBV values in part of the hypoperfused parenchyma after revascularization.

Revascularization with intra-arterial approach. MCA occlusion could be successfully reopened.

## General comments

The individual X-ray delay for *syngo* Neuro PBV IR can be analyzed by *syngo* iFlow.

By combining C-arm rotational acquisition with intraarterial contrast medium injection from the aortic arch, we can obtain

- (1) non-contrast *syngo* DynaCT,
- (2) a three-dimensional volume of intracranial vasculature (*syngo* DynaCT reconstruction of fill run), and
- (3) a *syngo* Neuro PBV IR map.


## Related publication

BioMed Research International Volume 2013, Article ID 873614, 7 pages  
<http://dx.doi.org/10.1155/2013/873614>


| Acquisition protocol  | 8s Neuro PBV IR                                                                                         |
|-----------------------|---------------------------------------------------------------------------------------------------------|
| Injection protocol    |                                                                                                         |
| Catheter position     | Aortic arch 4 F pigtail catheter                                                                        |
| Contrast medium (CM)  | 340 mg iodine/mL                                                                                        |
| Dilution              | 50%                                                                                                     |
| Injection volume      | 85 mL                                                                                                   |
| Injection rate        | 5 mL/s                                                                                                  |
| Duration of injection | 17 s                                                                                                    |
| X-ray delay           | Individually analyzed with syngo iFlow Tmax = 9 seconds for this patient in the superior sagittal sinus |
| Power injector used   | Yes                                                                                                     |

| Reconstructions       | Primary             | Secondary                            |
|-----------------------|---------------------|--------------------------------------|
| Name                  | Neuro PBV HU Smooth | DynaCT Head Nat Fill HU <sup>1</sup> |
| VOI size              | Full                |                                      |
| Slice matrix          | 512 × 512           |                                      |
| Kernel type           | HU                  | HU                                   |
| Image characteristics | Smooth              | Normal                               |
| Reconstruction mode   | Neuro PBV           | Nat Fill                             |
| Viewing preset        | Neuro PBV           | DynaCT Head                          |


<sup>1</sup> Reconstruction of the fill run of Neuro PBV 3D acquisition




Noncontrast CT demonstrated a hyperdense MCA sign (arrow) and excluded intracranial hemorrhage




Thick MIP – FDCT angiography demonstrated the total occlusion of left MCA (arrow)



Thick MIP – After intra-arterial mechanical thrombectomy, recanalization of the left MCA was demonstrated by FDCT



A parenchymal cerebral blood volume (syngo Neuro PBV IR) map depicted a large area of hypoperfusion in the corresponding left MCA territory



A syngo Neuro PBV IR map depicted the recovery of CBV values (circle) in part of the hypoperfused parenchyma after revascularization



# Visualization of pontine hemorrhage in *syngo* DynaCT Micro

---

|                              |                                                                                                |
|------------------------------|------------------------------------------------------------------------------------------------|
| <b>Courtesy of</b>           | Rahul Kumar, MD<br>GS Neuroscience Clinic and Research Centre Pvt. Ltd.<br>Patna, Bihar, India |
| <b>Supported by</b>          | <i>syngo</i> DynaCT Micro                                                                      |
| <b>System &amp; Software</b> | ARTIS icono biplane VE20 with<br><i>syngo</i> Application Software VE20                        |

---

## Patient history

A 32-year-old male patient presented with early acute stroke symptoms of vertigo, dizziness, and severe imbalance. He was to be diagnosed with a central pontine hematoma.

## Diagnosis

The native computed tomography (CT) revealed bleeding in the brain stem. It was decided to perform a DynaCT scan with contrast media injection to better understand the involved vascular structures. The cone beam CT angiography shows the hematoma cavity, which is seen as an avascular zone in the pontine tectum. More interesting is the beading seen on the para-median perforators that are the homologous vessels to the sulco commissurals in the spinal cord.

## Treatment

The patient was found to have elevated homocysteine levels. We decided for conservative treatment and he was medicated appropriately. The patient has responded very well to the medication. A follow-up CT scan showed complete resorption of the hemorrhage.

## General comments

"The *syngo* DynaCT Micro image showing the beading of the vessels is probably the singular most impressive image yet."

## Tips and tricks

Adequate and continuous contrast filling must be maintained throughout the scan. The use of diluted contrast medium is recommended, as high contrast concentrations can generate beam-hardening artifacts, potentially obscuring small vascular structures.


|                             |                           |
|-----------------------------|---------------------------|
| <b>Acquisition protocol</b> | <b>14s DCT Head Micro</b> |
| <b>Injection protocol</b>   |                           |
| Dilution                    | 33–35%                    |
| Injection volume            | 63–84 mL                  |
| Power injector used         | Yes                       |
| Injection rate              | 3–4 mL/s                  |
| Duration of injection       | 21 s                      |
| X-ray delay                 | 7 s                       |
| Catheter position           | Vertebral artery          |
| <b>Reconstructions</b>      |                           |
| <b>Primary</b>              |                           |
| Name                        | DCT Head Clear            |
| VOI size                    | Full                      |
| Slice matrix                | 512 x 512                 |
| Kernel type                 | HU                        |
| Image characteristics       | Normal                    |
| Reconstruction mode         | Fill                      |
| Viewing preset              | DCT Head                  |

CT Head






Hemorrhage is visible.

CT Head follow-up



## DCT Micro



The hematoma cavity is seen as an avascular zone in the pontine tectum. Zoomed image illustrates beading of vessels.

# One stop management for aneurysmal subarachnoid hemorrhage

**Courtesy of** Marios N. Psychogios, MD, Ioannis Tsogkas, MD,  
Department of Diagnostic and Interventional  
Neuroradiology University Medicine Goettingen, Germany

**Supported by** *syngo* DynaCT

**System & Software** Artis Q biplane VD11, *syngo* X Workplace VD10

## Patient history

62-year-old female patient transferred from a peripheral hospital with an aneurysmal subarachnoid hemorrhage. Direct admission to angio suite and acquisition of a native *syngo* DynaCT, followed by an IV injected DynaCT run to assess the progression and the cause of the bleeding.

## Diagnosis

Ruptured aneurysm of the bifurcation of the right middle cerebral artery.

## Treatment

Surgical clipping.


## General comments

The immediate diagnosis, imaging, and treatment of a ruptured aneurysm is very important and prevents rebleeding. Intrahospital time delays can be avoided if a patient from a peripheral hospital with an aneurysmal hemorrhage depicted on NCCT is transferred directly to the angio suite. The ruptured aneurysm can be reliably depicted on an FDCT angiogram and

the DSA can be started without further transfer of the patient. IV *syngo* DynaCT offers an excellent assessment of the aneurysm and adjusted arteries.

## Tips and tricks

Use a headholder to avoid movement when the patient is not sedated.




MPR 3 mm axial slice shows subarachnoid bleeding (minor artifacts due to patient motion of not sedated patient)

| Acquisition protocol      | 10sDCT Head                                          |
|---------------------------|------------------------------------------------------|
| <b>Injection protocol</b> |                                                      |
| Catheter position         | IV injection                                         |
| Contrast medium (CM)      | 400 mg iodine/mL                                     |
| Dilution                  | None                                                 |
| Injection volume          | 60 mL contrast media followed by 60 mL saline chaser |
| Injection rate            | 5.0 mL/s                                             |
| Duration of injection     | 12 s – 12 s                                          |
| X-ray delay               | bolus tracking with digital subtraction angiography  |

| Reconstructions       | Primary        | Secondary      |
|-----------------------|----------------|----------------|
| Name                  | DCT Head Clear | DCT Head Clear |
| VOI size              | Full           | Small          |
| Slice matrix          | 512 × 512      | 512 × 512      |
| Kernel type           | HU             | HU             |
| Image characteristics | Normal         | Normal         |
| Reconstruction mode   | Nat Fill       | Nat Fill       |
| Viewing preset        | DCT Head       | DCT Head       |



MIP 10 mm frontal slice shows aneurysm and surrounding vessel structure (metal artifacts due to dental implants)



MIP 0.3 mm axial slice shows aneurysm in secondary, small VOI reconstruction

# Flow diverter placement in a subarachnoid hemorrhage case supported by 'SuperDyna' imaging technique

---

**Courtesy of** St Vincent's Hospital Radiology Department,  
Melbourne, Australia

**Supported by** *syngo* Dyna3D, *syngo* 3D/3D Fusion, *syngo* DynaCT

**System & Software** ARTIS icono biplane VE40A  
with *syngo* Application Software VE40C

---

## Patient history

60-year-old male patient with a known subarachnoid hemorrhage (grade 2) of unclear origin. Sudden drop in Glasgow Coma Score – query aneurysm/malformation

## Diagnosis

Right internal carotid artery blister aneurysm

## Treatment

A diagnostic cerebral angiogram was performed. A blister aneurysm evident on the right internal carotid artery (ICA) was confirmed by performing a 4 s Dyna3D DSA spin with the injection of contrast media using a catheter in the right ICA. Two flow diverter stents were then inserted to cover the parent vessel of the blister aneurysm (5.0 x 21/14 mm and 5.0 x 15/9 mm). Typically, two or more flow diverters are needed to treat blister aneurysms specifically.


After the deployment of both stents, a selective diagnostic angiogram was performed using a catheter in the right ICA. Then a native 8 s DynaCT run was performed. In a secondary reconstruction this volume was refined with a small VOI to increase spatial resolution. Persistent filling of the right ICA blister aneurysm without significant change in size or appearance could be seen. By fusing the subtracted reconstruction volume of the initial 4 s Dyna3D run, and the refined volume of the 8 s DynaCT run the patency of the two stents was visualized in significant detail, providing a crucial diagnosis.

## Tips and tricks

Proper preparation of the patient is crucial – a head holder and head strap/sponges are essential to prevent patient movement during the procedure. This must be considered prior to creating the sterile procedure field.


|                       |                     |
|-----------------------|---------------------|
| Acquisition protocol  | 4s DSA Dyna3D Head  |
| Injection protocol    |                     |
| Catheter position     | Right ICA           |
| Contrast medium       | 240 mg iodine/mL    |
| Dilution              | 100%                |
| Injection volume      | 16 mL               |
| Injection rate        | 4.0 mL/s            |
| Duration of injection | 4 s                 |
| X-ray delay           | 1 s                 |
| Power injector used   | Yes                 |
| Acquisition protocol  | 8 s DCT Head 109 kV |
| Injection protocol    | NA                  |

| Reconstructions       | Primary           | Secondary         |
|-----------------------|-------------------|-------------------|
| Name                  | DynaCT Head Clear | DynaCT Head Clear |
| VOI size              | Full              | Small             |
| Voxel size            | 0.5 mm            | 0.11 mm           |
| Slice matrix          | 512 x 512         | 512 x 512         |
| Kernel type           | HU                | EE                |
| Image characteristics | Normal            | Sharp             |
| Reconstruction mode   | Nat Fill          | Nat Fill          |
| Viewing preset        | DynaCT Head       | DynaCT Head       |




Subtracted reconstruction of 4 s Dyna3D run visualized in VRT mode, depicting presence of right ICA blister aneurysm (shown with 3D marker)


Flow diverter placement in subarachnoid hemorrhage case supported by 'SuperDyna' imaging technique




MIP thin in alignment with stents



MIP thin in alignment with stents




MIP thin in alignment with stents



Dual volume visualization in VRT mode



Dual volume visualization in VRT mode  
Dual volume visualization with translucent  
color preset



Dual volume visualization in VRT mode  
Dual volume visualization with translucent  
color preset

# Pipeline flow diverting stent follow-up

---

|                              |                                                                                                         |
|------------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Courtesy of</b>           | David Niemann, MD, University of Wisconsin Hospitals and Clinics, USA                                   |
| <b>Supported by</b>          | <ul style="list-style-type: none"><li>• <i>syngo</i> Dyna3D</li><li>• <i>syngo</i> DualVolume</li></ul> |
| <b>System &amp; Software</b> | ARTIS icono biplane VE20 with <i>syngo</i> Applications Software VE2                                    |

---

## Patient history

67-year-old female patient.

## Diagnosis

Ophthalmic aneurysms on right and left side treated with pipeline embolization devices (PEDs)

## Treatment


Longitudinal follow-up of implanted pipeline embolization devices (PEDs) using 3D-DSA acquisition and *syngo* DualVolume visualization on ARTIS icono biplane system. Diagnostic angiogram to follow up right ophthalmic aneurysm treated with three PEDs and left ophthalmic aneurysm treated with a single PED. Volumetric imaging acquired on right ophthalmic aneurysm using ARTIS icono. The

study shows residual contrast inflow into the right ophthalmic aneurysm and opacification within the extruding vessel, as well as good vessel wall apposition of the three PEDs. Downstream vasculature is well opacified and no in-stent stenosis is observed.


## General comments

For the volumetric imaging, we used a 4 s *syngo* Dyna3D run in Twin Spin mode. For the first time, ARTIS icono gives us the ability to perform the 3D run without the need to park the lateral plane.

| Acquisition protocol      | 4sDSA Head              |
|---------------------------|-------------------------|
| <b>Injection protocol</b> |                         |
| Catheter position         | Internal carotid artery |
| Contrast medium (CM)      | 300 mg/mL               |
| Dilution                  | None                    |
| Injection volume          | 15 mL                   |
| Injection rate            | 3 mL/s                  |
| Duration of injection     | 5 s                     |
| X-ray delay               | 2 s                     |
| Power injector used       | Yes                     |
| <b>Reconstructions</b>    |                         |
| <b>Primary</b>            |                         |
| Name                      | Dyna3D DSA Dual         |
| VOI size                  | Medium                  |
| Slice matrix              | 512 × 512               |
| Kernel type               | EE                      |
| Image characteristics     | Normal                  |
| Reconstruction mode       | Dual                    |
| Viewing preset            | 3D Dual Volume          |



Residual aneurysm well visible in VRT, as well as good understanding of stent wall apposition.



In-stent stenosis can be excluded by the MPR image.

# Embolization of an C-ophthalmic artery aneurysm with coil and stent

---

**Courtesy of** Naoki Kato, MD, Yuichi Murayama, MD  
The Jikei University School of Medicine, Tokyo, Japan

**Supported by**

- *syngo* Dyna3D
- *syngo* DynaCT Micro
- *syngo* Dyna4D

**System & Software** ARTIS pheno VE1 with *syngo* Application Software VD2

---

## Patient history

A 54-year-old woman underwent MRI for a headache in 2006. An ophthalmic artery aneurysm was discovered. The size had not changed on followup, but the patient was worried and requested treatment.

## Diagnosis

The maximum diameter of the aneurysm was 5 mm. Based on the location of the aneurysm endovascular treatment was recommended. Diagnostic DSA was performed in 2017 and right ICAG (internal carotid artery angiography) revealed that the aneurysm was located at the origin of the ophthalmic artery at the bifurcation of the C2 segment of the ICA. ECAG (external carotid artery angiography) combined with a balloon occlusion test (BOT) did not show collateral flow (i.e., no choroidal blush). To avoid the complication of loss of vision, it was therefore of utmost importance to preserve the ophthalmic artery.

## Procedure description

Access was established through the right femoral artery and a catheter was advanced until it reached the left ICA. A 4D DSA scan was then taken to allow the relative location of the catheter to be determined before the stent and coil deployment. To perform the embolization, a stent was deployed to cover the neck of the aneurysm and 10 coils were used to fill the aneurysm. A high filling rate of 35.6% was achieved without sacrificing the ophthalmic artery.



## General comments

*syngo* DynaCT provided excellent spatial resolution and thereby helped to determine the border between the aneurysm and ophthalmic artery. This was important for preserving the patient's visual function. The robotic C-arm can be positioned at any angle to optimize visualization during the procedure. In 3D imaging, the robotic C-arm quickly moves to align with any

image position, providing the endovascular surgeons with views from several angles. This is important for the success of the procedure.

**Tips and tricks**  
3D imaging helps to determine the optimal angle of the C-arm for intraprocedural fluoroscopy, ensuring safe deployment of stent and coil placement.


| Acquisition protocol  | 6s Dyna4D       |                                    |
|-----------------------|-----------------|------------------------------------|
| Injection protocol    |                 |                                    |
| Catheter position     | Left ICA        |                                    |
| Contrast medium (CM)  | 320 mg/mL       |                                    |
| Dilution (CM/Saline)  | None            |                                    |
| Injection volume      | 18 mL           |                                    |
| Injection rate        | 3.0 mL/s        |                                    |
| Duration of injection | 6.0 s           |                                    |
| X-ray delay           | 0 s             |                                    |
| Power injector used   | Yes             |                                    |
|                       |                 |                                    |
| Reconstructions       | Primary         | Secondary                          |
| Name                  | Dyna4D arterial |                                    |
| VOI size              | Full            | Small for more detailed resolution |
| Slice matrix          | 512 × 512       |                                    |
| Kernel type           | HU              |                                    |
| Image characteristics | Auto            |                                    |
| Reconstruction mode   | Sub             |                                    |
| Viewing preset        | Dyna4D          |                                    |




Not only does syngo Dyna4D show the morphology from different angles




..but also makes it possible to see it at different timepoints




VRT visualization before intervention



Dual volume visualization after intervention



VRT sub-visualization after intervention



*syngo* DynaCT Micro run, reconstructed with *syngo* DynaCT SMART shows the stent in relation to the coil pack

# Flow Diverter Stenting of Large Intracranial Aneurysm

---

**Courtesy of** Yuichi Murayama, MD, Toshihiro Ishibashi, MD, Ichiro Yuki, MD, Jikeii University School of Medicine, Tokyo, Japan

**Supported by**

- *syngo* Dyna4D
- *syngo* DynaCT Micro

**System & Software** Artis Q biplane VD11, *syngo* X Workplace VD10

---

## Patient history

60-year-old male patient with large intracranial aneurysm.

## Diagnosis

Large intracranial aneurysm located in the cavernous segment of the internal carotid artery. The diameter of the aneurysm at its largest point was 17 mm.

## Procedure description

A flow diverter (Covidien Pipeline®) was inserted under general anesthesia according to the general guidelines for flow diverter (FD) treatment in Japan. After initial deployment of the flow diverter, the following 3D acquisitions were performed: A native 20 s DynaCT Micro run in zoom 3 (22 cm) to visualize the metal implant (flow diverter) followed by a 6 s Dyna4D run to obtain information about the vessels' architecture and the flow dynamics. A secondary reconstruction of the DynaCT Micro run was performed, then this volume was fused with the Dyna4D volume.

Based on these images, the physicians concluded that the apposition of the FD to the vessel wall was sub-optimal (see red arrows in the following image) and could lead to thrombo-embolic events in the future. Therefore, they returned and performed a postdilatation of the flow diverter, followed by a second DynaCT Micro run.

The post-dilatation *syngo* DynaCT image was then fused with the available Dyna4D volume (since the patient was under general anesthesia, it was assumed, observed, and confirmed that there was no head motion between the scans).

By displaying the images side by side, the physicians were able to clearly see the achieved improvement in stent apposition: The *syngo* Dyna4D clearly showed how slow the inflow of blood into the aneurysmal structure had become, confirming the expected effect.

General comments

The physicians were impressed by the image quality, the level of detail achieved, and the accuracy of the fusion functionality. The case was very successful, and the physicians unanimously emphasized that the imaging capabilities of the Artis Q system with the PURE® platform enabled them to better visualize the vessels and the device during the


procedure and to significantly improve their treatment decisions and outcome.

Tips and tricks


Perform a secondary reconstruction of the DynaCT Micro run with “Sharp” image characteristics to improve visibility of stent struts.

|                       |                    |                   |
|-----------------------|--------------------|-------------------|
| Acquisition protocol  | 20s DCT Head Micro |                   |
| Injection protocol    |                    |                   |
| Contrast medium (CM)  | N/A                |                   |
|                       |                    |                   |
| Reconstructions       | Primary            | Secondary         |
| Name                  | DynaCT Head Micro  | DynaCT Head Micro |
| VOI size              | Full               | Full              |
| Slice matrix          | 512 × 512          | 512 × 512         |
| Kernel type           | HU                 | HU                |
| Image characteristics | Normal             | Sharp             |
| Reconstruction mode   | Nat Fill           | Nat Fill          |
| Viewing preset        | DynaCT Head        | DynaCT Head       |


First and second DynaCT Micro run



Fused visualization of secondary reconstruction of first DynaCT Micro run with Dyna4D run



Sub-optimal apposition of FD to the vessel wall (red arrows)



After post-dilatation



Before post-dilatation in orange, after in white, with measurements showing the shortening and the differences in wall location

| Acquisition protocol      | 6s DSA Dyna4D Head      |
|---------------------------|-------------------------|
| <b>Injection protocol</b> |                         |
| Catheter position         | Internal carotid artery |
| Contrast medium (CM)      | 270 mg iodine/mL        |
| Dilution                  | None                    |
| Injection volume          | 18 mL                   |
| Injection rate            | 3 mL/s                  |
| Duration of injection     | 6 s                     |
| X-ray delay               | 0 s                     |
| Power injector used       | Yes                     |

| Reconstructions       | Primary      |
|-----------------------|--------------|
| Name                  | Dyna4D Sub4D |
| VOI size              | Full         |
| Slice matrix          | 512 × 512    |
| Kernel type           | HU           |
| Image characteristics | Auto         |
| Reconstruction mode   | Sub          |
| Viewing preset        | Dyna4D       |



VRT time ranges. Dyna4D volume fused with DynaCT Micro volume, before post-dilatation

# Large frontal AVM embolization using *syngo* 3D/3D Fusion technique

---

**Courtesy of** René Chapot, MD, Interventional Neuroradiology, Alfried Krupp Hospital, Essen, Germany

**Supported by** *syngo* Dyna3D, *syngo* 3D/3D Fusion

**System & Software** ARTIS icono biplane VE30 with *syngo* Application Software VE30

---

## Diagnosis

Large frontal arteriovenous malformation (AVM), fed by the right internal carotid artery (ICA).

## Treatment

Liquid embolization

## General comments

The complicated vessel structure in an AVM is not easy to understand. MPR images make it easier to differentiate arteries from veins and the nidus. We use a technique, which we call “6D”, it consists of two fused color-coded 3D DSA runs. The understanding of the vessels comes from looking at the fused thin-slice MPRs in dedicated color presets, one blue and one orange. The benefit of this fusion technique is the ability to look inside the AVM and see which part comes from which artery (with one 3D volume in blue and one in orange) and where the veins are. Thanks to the color coding and the mix of orange and blue, the veins appear whitish.

## Procedure workflow

In addition to regular DSA scenes, a first 3D DSA run is achieved using an injection from the right ICA. This shows the territory of the middle cerebral artery (MCA) and the anterior cerebral artery (ACA).

Venous filling must be achieved from the beginning of the rotation. This requires a certain X-ray delay.

The second 3D DSA run is done from the left carotid artery and depicts the ACA.

## *syngo* 3D/3D Fusion workflow

- Load the first 3D DSA reconstructed dataset, use the blue color preset.
- Load the second 3D DSA dataset, use the orange color preset (dedicated color presets are required).
- Confirm “Load fused” in the pop-up window.
- Scroll through the MPRs to precisely see which vessel comes from the MCA and which from the ACA.

- Switch from coronal to sagittal view to better understand the main AVM vein. The course can be followed by scrolling through until the sagittal sinus is visible.

### Tips & Tricks

To achieve optimal image quality in the 3D volumes, it is necessary to fix the head (a dedicated head-holder with fixation straps is recommended) and prevent any patient motion.


### Statement

“6D” is a major step forward in understanding the AVM anatomy and defining the treatment approach – whether to start the embolization from the arteries or the veins.

| Acquisition protocol      | 4 s DSA Head Micro (for both runs)                                                   |
|---------------------------|--------------------------------------------------------------------------------------|
| <b>Injection protocol</b> |                                                                                      |
| Catheter position         | 1) right ICA; 2) left ICA                                                            |
| Contrast medium           | 300 mg/mL                                                                            |
| Dilution                  | 70%                                                                                  |
| Injection volume          | 40 mL                                                                                |
| Injection rate            | ~5 mL/s                                                                              |
| Duration of injection     | ~8 s                                                                                 |
| X-ray delay               | ~3 s                                                                                 |
| Power injector used       | Hand injection; start injection during back rotation, when C-arm passes AP position. |


| Reconstructions       | Primary                 | Secondary               |
|-----------------------|-------------------------|-------------------------|
| Name                  | 3D Head Sub Full Normal | 3D Head Sub Full Normal |
| VOI size              | Full                    | Full                    |
| Slice matrix          | 512x512                 | 512x512                 |
| Kernel type           | EE                      | EE                      |
| Image characteristics | Normal                  | Normal                  |
| Reconstruction mode   | Sub                     | Sub                     |
| Viewing preset        | 6D blue                 | 6D orange               |

DSA




DSA image of AVM in lateral view.

VRT




VRT 1: Injected in right ICA




VRT 2: Injected in left ICA

VRT visualization of each 3D run.

## MPR





Arteries coming from the ACA are displayed in orange; those from the MCA are in blue.



As the vein is fed by two compartments (with the vessels from the ACA in orange, and the MCA in blue), it appears whitish.

Coronal MPR slice showing AVM in syngo 3D/3D Fusion.



In the sagittal MPRs, you can easily follow the course of the common vein.

# High-flow dural fistula treatment decision supported by *syngo* iFlow

---

|                              |                                                                                                |
|------------------------------|------------------------------------------------------------------------------------------------|
| <b>Courtesy of</b>           | Rahul Kumar, MD<br>GS Neuroscience Clinic and Research Centre Pvt. Ltd.<br>Patna, Bihar, India |
| <b>Supported by</b>          | <i>syngo</i> Dyna3D, <i>syngo</i> iFlow, <i>syngo</i> DualVolume,<br><i>syngo</i> DynaCT Micro |
| <b>System &amp; Software</b> | ARTIS icono biplane VE21 with<br><i>syngo</i> Application Software VE21                        |

---

## Patient history

A 45-year-old male patient presented with progressively worsening sensorium and obtundation over the preceding 18 months. He was being treated as a case of cortical venous sinus thrombosis. Clinically, he had no lateralizing deficits, but had mild neck stiffness and bilateral papilledema. A lumbar puncture showed increased ICP: 310 mm of water.

## Diagnosis

An MRI was performed and showed features suggestive of intracranial venous hypertension, and multiple flow voids in the region of the left lateral sinus. A conventional catheter angiogram was performed.

DSA revealed a high-flow dural fistula in the region of the left transverse and left sigmoid sinuses, supplied by the left occipital artery, the left posterior cerebral artery, and the left meningo-hypophyseal trunk.

## Treatment

Endovascular embolization was performed, and complete exclusion of the fistula was seen in the immediate postprocedural runs. However, the venous drainage, though better angiographically, was still showing slowing.

The dilemma was whether to continue embolization, do a thrombectomy of the chronically occluded sinuses, or just terminate the procedure and do further medical management.

This is where the *syngo* iFlow application came to our rescue. Based on this analysis, the procedure was completed and no further intervention was done.

Clinically, the patient's condition metamorphosed completely within 72 hours of the procedure.

To check the patient's progress, another angiogram was done two weeks after the procedure.

The *syngo* iFlow application showed further shortening of the transit time across the brain, paralleling the clinical improvement. The crux of this case was the slow flow of the venous blood through the medullary venous system of the brain because of the relatively higher pressure in the superior sagittal sinus, as a result of the dural fistula. Once the fistula is eliminated,

theoretically there should be an increase in the venous flow velocity across the brain parenchyma, thus shortening the transit time of blood through the brain.


**General comments**

The treatment decision and progress was supported strongly by the information obtained from *syngo* iFlow.

**Protocol – Part 1**


| Acquisition protocol      | DSA VFR                                   |
|---------------------------|-------------------------------------------|
| Framerate                 | 2 F/s for 20 frames – 1 F/s for 29 frames |
| Slice matrix              | 2k live imaging                           |
| Zoom                      | 32 cm                                     |
| <b>Injection protocol</b> |                                           |
| Contrast medium (CM)      | 300 mg iodine/mL                          |
| Dilution                  | 100%                                      |
| Injection volume          | n.a.                                      |
| Power injector used       | Hand injection                            |
| Injection rate            | 3 mL/s                                    |
| Duration of injection     | n.a.                                      |
| X-ray delay               | n.a.                                      |
| Catheter position         | Common carotid artery                     |

**syngo iFlow before treatment**




The blood took 19.68 seconds to reach the internal cerebral vein, and 23.94 seconds to reach the superficial middle cerebral vein before the fistula was closed.

**syngo iFlow immediately after treatment**



Post-embolization, the blood took 14.36 seconds to reach the internal cerebral vein, and 18.62 seconds to reach the superficial middle cerebral vein immediately after the fistula was closed.


**syngo iFlow two weeks after treatment**

Post-embolization at the two-week interval, the blood took 12.23 seconds to reach the internal cerebral vein, and 14.36 seconds to reach the superficial middle cerebral vein.


## Protocol – Part 2

| Acquisition protocol      | 14 s DCT Head Micro |
|---------------------------|---------------------|
| <b>Injection protocol</b> |                     |
| Contrast medium (CM)      | 300 mg iodine/mL    |
| Dilution                  | 33–35%              |
| Injection volume          | 63–84 mL            |
| Power injector used       | Yes                 |
| Injection rate            | 3–4 mL/s            |
| Duration of injection     | 21 s                |
| X-ray delay               | 7 s                 |
| <b>Reconstructions</b>    |                     |
| Name                      | DCT Head            |
| VOI size                  | Full                |
| Slice matrix              | 512 x 512           |
| Kernel type               | HU                  |
| Image characteristics     | Sharp               |
| Reconstruction mode       | Nat                 |
| Viewing preset            | DynaCT Head         |

### syngo DynaCT Micro



Frontal view 12 mm MIP



Sagittal view 12 mm MIP


Showing the congested, worm-like structures of the deep medullary veins.

Protocol – Part 3


| Acquisition protocol  | 4 s DSA Head     |
|-----------------------|------------------|
| Injection protocol    |                  |
| Contrast medium (CM)  | 300 mg iodine/mL |
| Dilution              | 100%             |
| Injection volume      | 18 mL            |
| Power injector used   | Yes              |
| Injection rate        | 3 mL/s           |
| Duration of injection | 6 s              |
| X-ray delay           | 2 s              |

| Reconstructions       | Primary          | Secondary         |
|-----------------------|------------------|-------------------|
| Name                  | 3D Head Sub Full | 3D Head Mask Full |
| VOI size              | Full             | Full              |
| Slice matrix          | 512 x 512        | 512 x 512         |
| Kernel type           | EE               | HU                |
| Image characteristics | Normal           | Normal            |
| Reconstruction mode   | Sub              | Nat               |
| Viewing preset        | 3D Head Dual     | 3D Head Dual      |

syngo Dyna3D DualVolume



AP view



Lateral view

Visualization of vessel structure in relation to embolic material in dual volume mode.

# Embolization of cerebral AVM

**Courtesy of** Renè Chapot, MD, Hannes Nordmeyer, MD, Interventional Neuroradiology, Alfried Krupp Hospital, Essen, Germany

**Supported by** *syngo* Dyna4D

**System & Software** Artis zee biplane with PURE VD11, *syngo* X Workplace VD10

## Patient history

47-year-old male with aphasia and motoric deficiency of right hand due to edema.

## Diagnosis


Symptomatic left hemispheric precentral AVM with giant venous aneurysm and brain edema.

## Treatment

Multiple endovascular treatments with embolization (Onyx and PHIL).


## General comments

Normally we make a couple of DSA runs with different angulations of the C-arm to see the feeding vessels in detail and also observe the drainage of the AVM. Thanks to *syngo* Dyna4D this information was readily available after one injection, we could also see the filling and drainage of the AVM in 3D which helped our planning for treatment.



Flow range – *syngo* Dyna4D shows flow patterns in 3D

| Acquisition protocol      | 6s Dyna4D              |
|---------------------------|------------------------|
| <b>Injection protocol</b> |                        |
| Catheter position         | ICA extracranial       |
| Contrast medium (CM)      | 300 mg iodine/mL       |
| Dilution                  | None                   |
| Injection volume          | 20 mL                  |
| Injection rate            | 3.3 mL/s               |
| Duration of injection     | 6 s                    |
| X-ray delay               | 0 s                    |
| Power injector used       | None                   |
| <b>Reconstructions</b>    |                        |
| <b>Primary</b>            |                        |
| Name                      | Dyna4D arterial Sub 4D |
| VOI size                  | Full                   |
| Slice matrix              | 512 × 512              |
| Kernel type               | EE                     |
| Image characteristics     | Auto                   |
| Reconstruction mode       | Sub                    |
| Viewing preset            | Dyna4D                 |



Radial range – syngo Dyna4D shows flow patterns in any angulation

# AVM treatment planning with *syngo* Dyna4D

---

|                              |                                                                                                                                           |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Courtesy of</b>           | Peter Mitchell, MD, Director Neurointervention Service, Department of Radiology, Royal Melbourne Hospital, Parkville, Victoria, Australia |
| <b>Supported by</b>          | <i>syngo</i> Dyna4D                                                                                                                       |
| <b>System &amp; Software</b> | Artis Q biplane with PURE® VD11, <i>syngo</i> X Workplace VD10                                                                            |

---

## Patient history

The patient was referred for a cerebral angiogram to assess an arteriovenous malformation that was discovered as an incidental finding.

## Procedure description

A 2 cm × 2.3 cm arteriovenous malformation demonstrated in the medial aspect of the posterior temporal lobe on the right side. This is supplied by branches arising from the P2 segment of the right PCS, and drains to the right transverse/s-shaped sinus via veins traversing the tentorium cerebelli. There is a 5 mm × 5 mm aneurysm positioned at the posterior aspect of the nidus (likely venous), best appreciated on the *syngo* Dyna4D volume.

The application of *syngo* Dyna4D provided the following additional information: *syngo* Dyna4D clearly demonstrates that one of the two main arteries (both from P2 segment of the rPCA) are not filling the nidus at all, which is not demonstrated on the DSA images.

Out of the one main feeding artery, it is possible to perfectly delineate which of the two “sub-arteries” are supplying the intranidal aneurysm. This allows the interventional team to clearly identify which artery to embolize first.

*syngo* Dyna4D was essential in identifying the treatment plan for this patient. Embolization of the AVM is due to take place later this year.

## General comments


Manipulation of the *syngo* Dyna4D volume allows the clinicians to visualize filling of the AVM in any projection of the 3D model. This results in the most effective planning of the embolization treatment.

| Acquisition protocol      | 12s Dyna4D       |
|---------------------------|------------------|
| <b>Injection protocol</b> |                  |
| Catheter position         | Vertebral artery |
| Contrast medium (CM)      | 240 mg iodine/mL |
| Dilution                  | None             |
| Injection volume          | 20 mL            |
| Injection rate            | 2.5 mL/s         |
| Duration of injection     | 8 s              |
| X-ray delay               | None             |
| Power injector used       | Yes              |

| Reconstructions       | Primary                 |
|-----------------------|-------------------------|
| Name                  | DynaCT Head Nat Fill HU |
| VOI size              | Full                    |
| Slice matrix          | 512 × 512               |
| Kernel type           | HU                      |
| Image characteristics | Normal                  |
| Reconstruction mode   | Nat Fill                |
| Viewing preset        | DynaCT Head             |



syngo Dyna4D is used to “untangle” the AVM to perfectly delineate which of the two subarteries are supplying the intranidal aneurysm



DSA image which demonstrates AVM arising from the right vertebral artery

# Understanding complex angio architectures of Dural AVF supported by *syngo* Dyna4D

---

**Courtesy of** Ichiro Yuki, MD, Yuichi Murayama, MD, Department of Neurosurgery, Jikei University Hospital, Tokyo, Japan

**Supported by** *syngo* Dyna4D

**System & Software** Artis Q biplane VD11, *syngo* X Workplace VD10

---

## Patient history

A 67-year-old male presented with diagnosed progressive visual acuity loss. An MR angiogram (MRA) indicated potential dural arteriovenous fistula (dAVF).

## Procedure description

The patient was referred to the endovascular neurosurgery department and underwent a cerebral angiogram. The angiogram revealed a left transverse sigmoid sinus dAVF with significant venous reflux into the superior sagittal sinus. Normal draining pattern of left transverse sigmoid sinus was impaired, and it was replaced by the high flow shunting between the multiple arterial branches from the left external carotid arteries and the impaired sinus. After completing the conventional 2D DSA, it was not yet clear whether or not the lesion had “cortical venous reflux”, which is decisive in determining the prognosis for the patient. A *syngo* Dyna4D run showed a detailed flow pattern of arteriovenous shunting as

well as a correlation between the complex angio architecture and the cranium. The image enables the 3D structure of cortical veins to be distinguished from overlapped vascular structures. A cortical venous reflux was confirmed. The retrograde filling of the contrast in the cortical vein was clearly visible due to the temporal information provided by *syngo* Dyna4D.


Based on the aforementioned image findings, this dAVF was classified as Cognard IIa + B and endovascular treatment was recommended. Transarterial embolization of the feeding artery followed by a trans-venous approach and occlusion of the left T-S sinus was performed.

## General comments


*syngo* Dyna4D breaks down the complex angio architecture of the shunting disease into different phases of contrast fillings, and allows to see the direction, speed, and amount of flow in one acquisition. This information is crucial to the treatment.

| Acquisition protocol  | 12s Dyna4D          |
|-----------------------|---------------------|
| Injection protocol    |                     |
| Catheter position     | Left carotid artery |
| Contrast medium (CM)  | 270 mg iodine/mL    |
| Dilution              | None                |
| Injection volume      | 21 mL               |
| Injection rate        | 3 mL/s              |
| Duration of injection | 7 s                 |
| X-ray delay           | None                |
| Power injector used   | Yes                 |


| Reconstructions                   | Primary               |
|-----------------------------------|-----------------------|
| Name                              | Sub4D*Full HU<br>Auto |
| VOI size                          | Full                  |
| Slice matrix                      | 512 × 512             |
| Kernel type                       | EE                    |
| Image characteristics             | Auto                  |
| Reconstruction mode               | Sub                   |
| Viewing preset                    | Dyna4D                |
| Sub 4D fused with Native Mask Run |                       |



Yellow arrow: superior temporal artery  
Red arrow: posterior auricular artery



Left external carotid artery angiogram (ECAG) shows the contrast dye welling up in the left sigmoid sinus, then flowing back to the transverse sinus, and finally to the superior sagittal sinus (red arrow)



Cortical venous reflux was seen at the left occipital lobe near the transverse sigmoid junction (yellow arrow)

# Follow-up after bilateral carotid stent placement

---

**Courtesy of** Kiffon Keigher, APN, Demetrius K. Lopes, MD,  
Rush University Medical Center, Chicago, USA

**Supported by** *syngo* DynaCT

**System & Software** Artis zee biplane VC21, *syngo* X Workplace VB21

---

## Patient history

76-year-old female who developed severe bilateral carotid artery stenosis thought to be primarily related to her history of neck radiation for thyroid disease. Re-stenosis of both carotids seen in CTA. Decision for *syngo* DynaCT with IV injection of contrast medium for improved stenosis evaluation and stent visualization.

## Procedure description

A 20sDR DynaCT of the neck region was performed, with a contrast injection through an 18 G IV access in the right antecubital vein, using an X-ray delay of 14 sec. Images were automatically reconstructed on the *syngo* X Workplace and displayed in *syngo* InSpace 3D.

Patient has no new symptoms and no further intervention was recommended with exception of continued dual anticoagulation therapy and stroke risk factor management.

## General comments

At our institution we are using *syngo* DynaCT with IV injection routinely for patients with the following concerns: stent patency, in-stent stenosis, residual filling of aneurysms s/p clipping and/or coiling, and vasospasm.




Coronal MIP showing bilateral stent placement


|                       |                                      |
|-----------------------|--------------------------------------|
| Acquisition protocol  | 20sDCT Head 109kV                    |
| Injection protocol    |                                      |
| Catheter position     | IV injection; right antecubital vein |
| Contrast medium (CM)  | 370 mg iodine/mL                     |
| Dilution              | None                                 |
| Injection volume      | 80 mL                                |
| Injection rate        | 4 mL/s                               |
| Duration of injection | 20 s                                 |
| X-ray delay           | 14 s                                 |
| Power injector used   | Yes                                  |

|                       |                         |                         |
|-----------------------|-------------------------|-------------------------|
| Reconstructions       | Primary                 | Secondary <sup>1</sup>  |
| Name                  | DynaCT Head Nat Fill HU | DynaCT Head Nat Fill HU |
| VOI size              | Full                    | Small                   |
| Slice matrix          | 512 × 512               | 512 × 512               |
| Kernel type           | HU                      | HU                      |
| Image characteristics | Normal                  | Normal                  |
| Reconstruction mode   | Nat Fill                | Nat Fill                |
| Viewing preset        | DynaCT Head             | DynaCT Head             |

<sup>1</sup> In order to improve the visualization of the stent, a 2<sup>nd</sup> reconstruction with a small size VOI was performed, using the same parameter as for the initial reconstruction.



Sagittal MPR demonstrating partial in-stent stenosis (with secondary reconstruction)



Curved MPR demonstrating partial in-stent stenosis (with secondary reconstruction)

# Follow-up after stent-assisted aneurysm coiling

**Courtesy of** Kiffon Keigher, APN, Demetrius K. Lopes, MD,  
Rush University Medical Center, Chicago, USA

**Supported by** *syngo* DynaCT

**System & Software** Artis zee biplane VC21, *syngo* X Workplace VB21

## Patient history

73-year-old female patient. Left middle cerebral artery (MCA) aneurysm that was treated with stent-assisted coil embolization in 2005. Because of her history of contrast allergies, we decided to complete a *syngo* DynaCT run with IV injection of contrast medium in the angiography suite for the purpose of improved visualization and monitoring the patient post IV injection.


## Procedure description

In 3D rotational angiography we identified that the left MCA aneurysm was completely occluded, status post stenting and coiling. We were able to obtain accurate visualization of the parent vessel and verify occlusion of this aneurysm and patency of the stent. There was no evidence of in-stent stenosis.

The quality of the study was very good, with no evidence of any other abnormalities in the circle of Willis. She tolerated this imaging study without incident and was discharged home. Our plan will be to follow up with this patient in 12 months with repeated imaging study.

## General comments

At our institution we are using *syngo* DynaCT runs with IV injection routinely for patients with the following concerns: stent patency, in-stent stenosis, residual filling of aneurysms s/p clipping and/or coiling, and vasospasm.




VRT (Volume Rendering Technique) of *syngo* DynaCT run with IV injection, visualizing the X stent


| Acquisition protocol      | 20sDCT Head 109kV                                            |
|---------------------------|--------------------------------------------------------------|
| <b>Injection protocol</b> |                                                              |
| Catheter position         | IV injection antecubital vein (18 G) allergy prep medication |
| Contrast medium (CM)      | 370 mg iodine/mL                                             |
| Dilution                  | None                                                         |
| Injection volume          | 80 mL                                                        |
| Injection rate            | 4 mL/s                                                       |
| Duration of injection     | 20 s                                                         |
| X-ray delay               | 14 s                                                         |
| Power injector used       | Yes                                                          |

| Reconstructions       | Primary                 | Secondary <sup>1</sup>  |
|-----------------------|-------------------------|-------------------------|
| Name                  | DynaCT Head Nat Fill HU | DynaCT Head Nat Fill HU |
| VOI size              | Full                    | Small                   |
| Slice matrix          | 512 × 512               | 512 × 512               |
| Kernel type           | HU                      | HU                      |
| Image characteristics | Normal                  | Normal                  |
| Reconstruction mode   | Nat Fill                | Nat Fill                |
| Viewing preset        | DynaCT Head             | DynaCT Head             |

<sup>1</sup> In order to improve the visualization of the stent, a 2<sup>nd</sup> reconstruction with a small size VOI was performed, using the same parameter as for the initial reconstruction.



Axial MPR reconstruction outlining full patency of the X stent in the left middle cerebral artery (with secondary reconstruction)



Axial MPR reconstruction showing X stent in the left middle cerebral artery

# Follow-up of aneurysm clipping using *syngo* DynaCT with IV injection

---

**Courtesy of** Peter Mitchell, MD, Department of Radiology,  
The Royal Melbourne Hospital, Melbourne, Australia

**Supported by**

- *syngo* DynaCT DSA
- *syngo* DualVolume

**System & Software** Artis Q biplane VD11, *syngo* X Workplace with  
*syngo* Application Software VD11

---

## Patient history

Past history of anterior communicating artery clipping (13 years ago), with subsequent elective middle cerebral artery embolization. Small bleb-like aneurysm adjacent to tip of aneurysm clips on DSA.

## Diagnosis

Intracranial aneurysms at both the anterior communicating and middle cerebral arteries. Small bleb-like aneurysm for assessment.

## Procedure description

Clip occlusion > 10 years ago, routine follow up showed small new bleb-like aneurysm on DSA. Difficult to identify on MRA. *syngo* DynaCT with intravenous contrast injection in DualVolume mode offered less invasive imaging, able to exclude growth into a true saccular aneurysm, and allowed ongoing surveillance without repeat DSA.

## General comments


When alternative non-invasive arterial imaging is limited, due to dense coil and stent packing, or certain types and numbers of clips, *syngo* DynaCT with intravenous contrast injection in DualVolume mode offers a technique that is less impacted by these artifacts, particularly when the clinical area of suspicion is at the base of the previously treated aneurysm.

## Tips and tricks


- Once the injector is armed, set the 3D protocol to manual mode (injection starts automatically after returning from mask spin)
- Start acquiring the fill run when the carotid artery is displayed on the bolus watch (approx. 8–15 sec)
- Head positioning: symmetrical patient head positioning with strap fixing to reduce patient movement during both spins

| Acquisition protocol      | 10s DSA DCT Head                           |
|---------------------------|--------------------------------------------|
| <b>Injection protocol</b> |                                            |
| Catheter position         | IV injection antecubital vein 18 g cannula |
| Contrast medium (CM)      | 350 mg iodine/mL                           |
| Dilution                  | None                                       |
| Injection volume          | 96 mL                                      |
| Injection rate            | 8 mL/s                                     |
| Duration of injection     | 12 s                                       |
| X-ray delay               | DSA bolus tracking                         |
| Power injector used       | Yes                                        |

| Reconstructions       | Primary                                                                    | Secondary                         |
|-----------------------|----------------------------------------------------------------------------|-----------------------------------|
| Name                  | DynaCT DSA Dual Head<br>(reconstructs Nat Fill +<br>Nat Mask + Sub volume) | DynaCT DSA Dual Head              |
| VOI size              | Full (Nat Fill volume)                                                     | Medium (Nat Mask + Sub<br>volume) |
| Slice matrix          | 512 × 512                                                                  | 512 × 512                         |
| Kernel type           | HU (Nat Fill + Nat Mask<br>volume)                                         | EE (Sub volume)                   |
| Image characteristics | Auto                                                                       | Auto                              |
| Reconstruction mode   | Vasc head                                                                  | Vasc head                         |
| Viewing preset        | Dual Volume 2                                                              | Dual Volume 2                     |



syngo DualVolume visualization of vessels and clips



Transversal MIP 5 mm

# Visualization of clipped aneurysm

**Courtesy of** Andrey Sergeev, MD, Department of Neurosurgery, The Federal State Budgetary Institution “Almazov National Medical Research Centre” of the Ministry of Health, Russia

**Supported by**

- syngo DualVolume
- syngo DynaCT DSA

**System & Software** Artis zeego, VC21, syngo X-Workplace, VB21

## Patient history

Patient hemorrhaged from a giant ACom aneurysm. Urgent admission and surgery was performed the next day.

## Diagnosis

Giant aneurysm of the ACom artery. Acute subarachnoid hemorrhage from the aneurysm, Hunt-Hess II.

## Procedure description

Acute craniotomy with intracranial anastomosis A2-A2 and aneurysm clipping.

Intraoperative imaging of cerebral vessels allows neurosurgeons to visualize the effectiveness of clipped aneurysms and patency of bypasses immediately after the surgery.

The injection of the iodine contrast takes place directly after the mask run. We observe the contrast bolus using DSA. As soon as the carotid artery is filled with contrast, we start the robotic C-arm rotation syngo DynaCT acquired in Zoom 1 (42 cm) delivers better resolution due to the smaller voxel size compared to syngo DynaCT

in Zoom 0 (48 cm), and results in better and more detailed visualization of the clips and vessels.

## General comments


DynaCT with intravenous contrast injection is a fast and easy procedure for noninvasive arterial imaging.



Clipped aneurysm. Intra-intracranial anastomosis A2-A2 with retrograde filling of right A2 vessel (secondary reconstruction)

|                       |                              |
|-----------------------|------------------------------|
| Acquisition protocol  | 10s DSA DCT Head, zoom 42 cm |
| Injection protocol    |                              |
| Catheter position     | Cubital vein 18 g            |
| Contrast medium (CM)  | 370 mg iodine/mL             |
| Dilution              | None                         |
| Injection volume      | 100 mL                       |
| Injection rate        | 5 mL/s                       |
| Duration of injection | 20 s                         |
| X-ray delay           | DSA bolus tracking           |
| Power injector used   | Yes                          |

| Reconstructions       | Primary                    | Secondary                                      |
|-----------------------|----------------------------|------------------------------------------------|
| Name                  | DynaCT Head<br>Nat Fill HU | DynaCT Head<br>Sub MoCo HU                     |
| VOI size              | Full                       | Full                                           |
| Slice matrix          | 512 × 512                  | 512 × 512                                      |
| Kernel type           | HU                         | HU                                             |
| Image characteristics | Auto                       | Very smooth                                    |
| Reconstruction mode   | Nat Fill                   | Sub with MoCo (motion<br>correction algorithm) |
| Viewing preset        | Dual Volume 2              | Dual Volume 2                                  |



3D visualization of the clipped aneurysm  
(secondary reconstruction)



Transversal MIP with vessels and clips (primary  
reconstruction)

# Follow-up of aneurysm treated with a flow diverter

---

**Courtesy of** Martin Skalej, MD, Department of Neuroradiology  
University Hospital Magdeburg, Germany

**Supported by** syngo DynaCT Micro

**System & Software** Artis Q biplane VD11, syngo X Workplace VD10

---

## Patient history

64-year-old male with two innocent paraophthalmic wide neck aneurysms in the left internal carotid artery.

## Treatment


Implantation of a flow diverter (Acandis Derivo 4,5 mm × 20 mm) 6 month ago. syngo DynaCT Micro acquisition for follow-up of flow diverter treated aneurysm to exclude intimal hyperplasia or in-stent stenosis.

## General comments


Cone beam CT with IV injection is a very easy workflow in our daily routine.

The follow up of flow diverter placement is easy in our daily routine. There is no need for catheter placement, only IV contrast injection is required. syngo DynaCT Micro, with its higher spatial resolution, improves the visualization of the stent.

| Acquisition protocol      | 20s DCT Head Micro           |
|---------------------------|------------------------------|
| <b>Injection protocol</b> |                              |
| Catheter position         | IV injection                 |
| Contrast medium (CM)      | 300 mg iodine/mL             |
| Dilution                  | 80%                          |
| Injection volume          | 100 mL                       |
| Injection rate            | 5 mL/s                       |
| Duration of injection     | 20 s                         |
| X-ray delay               | 20 s                         |
| Power injector used       | Yes                          |
| <b>Reconstructions</b>    |                              |
| <b>Primary</b>            |                              |
| Name                      | DynaCT Micro                 |
| VOI size                  | Manual 78 mm × 78 mm × 75 mm |
| Slice matrix              | 512 × 512                    |
| Kernel type               | HU                           |
| Image characteristics     | Normal                       |
| Reconstruction mode       | Nat Fill                     |
| Viewing preset            | DynaCT Head                  |



Implanted flow diverter in the ophthalmic segment of the left ICA lateral: anterior clinoid process, medial: sphenoid sinus, no intimal hyperplasia or in-stent-stenosis



No intimal hyperplasia or in-stent-stenosis

# Transradial embolization of the middle meningeal artery for chronic subdural hematoma

**Courtesy of** Riitta Rautio, MD, Head of Interventional Radiology, TYSK, Turku University Hospital, Turku, Finland

**Supported by** *syngo* DynaCT

**System & Software** ARTIS icono biplane VE30 with *syngo* Application Software VE30

## Patient history

A 79-year-old female patient fell down two weeks ago and thereafter suffered from gradual loss of balance and mild hemiparesis. She currently smokes and has a coronary bypass in her history. Ticagrelor 90 mg x 2 and ASS 100 mg x 1 after recent coronary intervention.

## Diagnosis

Chronic subdural hematoma discovered on the right side in native CT imaging.

## Treatment

Subdural hematoma evacuation was performed with burr hole technique by the Neurosurgeons in the OR theater due to the large mass effect and midline shift. Continuation of antithrombotic medication was required due to the patient's previous coronary intervention. Middle meningeal artery embolization with Onyx was performed two weeks after evacuation.

## General comments


A native flat-detector CT scan is performed before the intervention

to check the size of the existing subdural hematoma (images 1–3).

Aortic arch tortuosity increases with age and can make transfemoral procedures longer and complicated. A transradial procedure allows early ambulation and discharge for suitable patients.

The procedure starts by positioning the patient's arm in the radial support on the angiography table, with the

Figure 1




flat detectors ready for radial approach (image 4). The change to biplane neuro imaging uses a preset system position with a very smooth transition (images below).

**Tips and tricks**

A head holder allows good stabilization of the patient’s head if the patient is awake during the procedure. When using liquid embolics that


contain tantalum for opacity, the empty tantalum roadmap reliably shows the forward penetration of the embolic agent. We are using roadmap presets, which use an X-ray spectrum optimized for the imaging characteristics of the material, in our case tantalum for Onyx, to achieve the best image contrast with the lowest possible dose (Structure Scout).

|                       |                            |
|-----------------------|----------------------------|
| Acquisition protocol  | 7sDCT Sine Spine           |
| Injection protocol    |                            |
| Catheter position     | n.a.                       |
| Reconstructions       |                            |
| Primary               |                            |
| Name                  | DynaCT Head Clear Nat Fill |
| VOI size              | Full                       |
| Slice matrix          | 512x512                    |
| Kernel type           | HU                         |
| Image characteristics | Smooth                     |
| Reconstruction mode   | Nat Fill                   |
| Viewing preset        | DynaCT Head                |




Native syngo DynaCT Sine Spin to visualize the size of the hemorrhage

**Figure 2**




DSA radial approach



Neuro DSA imaging using a variable frame rate protocol (4/2/1)

**Figure 3**




Empty roadmap image (subtracted roadmap with tantalum settings) to watch the inflow of the liquid embolic material



Native single shot after embolization



Radial access setting



With the help of stored system positions, changing from radial access position to biplane imaging position is fast and easy



System position for biplane imaging

# **syngo DynaCT of petrous bones to visualize cholesteatoma**

---

|                              |                                                                |
|------------------------------|----------------------------------------------------------------|
| <b>Courtesy of</b>           | Tobias Struffert, MD,<br>Universitätsklinikum Giessen, Germany |
| <b>Supported by</b>          | syngo DynaCT Head                                              |
| <b>System &amp; Software</b> | ARTIS icono biplane VE2 with<br>syngo Application Software VE2 |

---

## **Patient history**

Male patient, who presented with hearing impairment.

## **Diagnosis**

Cholesteatoma of the left ear. A cholesteatoma is a destructive and expanding growth consisting of keratinizing squamous epithelium in the middle ear and/or mastoid process. The reconstructed *syngo* DynaCT volume shows eroded ossicles, liquid-filled mastoid cells, and swollen intima in the meatus externa.

## **Treatment**

Surgical intervention to remove the cholesteatoma.

## **General comments**

The 8s DCT Head acquisition protocol has the advantage of acquiring both sides at once. Collimate top-down for the desired ROI. A secondary reconstruction of left and right petrous bone separately provides a very high spatial resolution. With this approach, the expansion of the cholesteatoma and related erosions were nicely depicted.


The approach is requested by our neurosurgeons and enhances their confidence during surgery.

## **Tips & Tricks:**


Explain the short procedure to the patient and the need to lie still during the rotation of the system. Ensure proper positioning of the patient and fixation of the head to achieve optimal image quality.

|                       |                                                                  |
|-----------------------|------------------------------------------------------------------|
| Acquisition protocol  | 8sDCT Head<br>Collimated (top-down) for the desired ROI          |
| Injection protocol    |                                                                  |
| Contrast medium (CM)  | no contrast                                                      |
| Test Bolus            | n.a.                                                             |
| Reconstructions       |                                                                  |
| Name                  | DynaCT Head                                                      |
| VOI size              | Manual; adjust VOI size to reach slice thickness of 0.15–0.17 mm |
| Slice matrix          | 512 x 512                                                        |
| Kernel type           | HU                                                               |
| Image characteristics | Sharp                                                            |
| Reconstruction mode   | Nat                                                              |
| Viewing preset        | DynaCT Head                                                      |


Coronal MPR




Petrol arrow shows attic sporn on healthy right petrous bone.



Orange arrow shows eroded attic sporn with thickened intima in meatus externa.




Petrol arrow shows normal anatomy of malleus.




Orange arrow shows destroyed malleus.


**Axial MPR**



Right incus in healthy structure.



Left incus shows erosion due to the cholesteatoma.



Petrol arrow shows normal pneumatized mastoid cells.



Orange arrow shows liquid-filled mastoid cells. Yellow arrow shows thickened intima of the meatus externa.



# Total ossicular replacement prosthesis

The high spatial resolution of *syngo* DynaCT Micro supports diagnosis of complications after TORP implantations

---

**Courtesy of** Martin Skalej Department of Neuroradiology, MD, University Hospital Magdeburg, Germany

**Supported by** *syngo* DynaCT Micro

**System & Software** Artis Q biplane VD11, *syngo* X Workplace VD10

---

## Patient history


35-year-old male with known conductive hearing loss; initially benefits after implantation of a total ossicular replacement prosthesis (TORP). After 3 months, he suddenly suffers from hearing loss again.

## Diagnosis

TORP dislocation


## General comments

*syngo* DynaCT Micro provides higher spatial resolution compared with standard CT images. For the complex anatomy of the inner ear and especially with tiny prosthesis and implants, it helps us to make the right diagnosis.




After tympanoplasty type III with titanium total ossicular replacement prosthesis, the tympanum is filled with soft tissue isodense material; the TORP is dislocated and bent

|                       |                              |
|-----------------------|------------------------------|
| Acquisition protocol  | 20s DCT Head Micro           |
| Injection protocol    |                              |
| Contrast medium (CM)  | N/A                          |
| Reconstructions       |                              |
| Name                  | DynaCT Micro                 |
| VOI size              | Manual 47 mm × 47 mm × 35 mm |
| Slice matrix          | 512 × 512                    |
| Kernel type           | HU                           |
| Image characteristics | Normal                       |
| Reconstruction mode   | Nat Fill                     |
| Viewing preset        | DynaCT Head                  |



Dislocated TORP without contact with the oval window (indicated by arrow)



Sample image of a TORP correctly positioned in the oval window (indicated by arrow)

# Overview of Advanced Applications

## 3D/Advanced 3D

|                                |                                                                                                                                                                                                         |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>syngo</i> Dyna3D            | performs a rotation around the patient, acquiring from different angles.                                                                                                                                |
| <i>syngo</i> Dual Volume       | enables simultaneous visualization of dual-volume imaging (e.g., stent and contrast-enhanced vessels).                                                                                                  |
| <i>syngo</i> DynaCT            | offers dedicated, enhanced cone-beam reconstruction algorithms that reconstruct a low-contrast, three-dimensional image of soft tissue. The C-arm rotates around the area of interest at the isocenter. |
| <i>syngo</i> Dyna4D            | is a rotational acquisition that delivers a virtually unlimited number of DSA runs with no additional radiation dose or contrast medium. Allows you to see flow patterns in 3D.                         |
| <i>syngo</i> DynaCT HighSpeed  | Provides a fast high contrast 3D protocol of less than 3 s, represented with reduced motion artifacts and less contrast media.                                                                          |
| <i>syngo</i> DynaCT Micro      | provides high spatial resolution 3D imaging by using each detector pixel for 40 percent higher resolution compared with a standard <i>syngo</i> DynaCT.                                                 |
| <i>syngo</i> DynaPBV Neuro     | enables 3D functional imaging acquisition protocols that deliver physiological information about the patient intraoperatively.                                                                          |
| <i>syngo</i> DynaCT Sine Spin  | helps neuroradiologists to detect bleeds in the interventional suite, especially when imaging the basal part of the brain and areas close to the skull.                                                 |
| <i>syngo</i> DynaCT Multiphase | enables visualization of collateral vessels with time-resolved DynaCT, depicting up to 10 different timepoints within a period of 60 seconds.                                                           |

## Advanced tools and workflows

|                                         |                                                                                                                                                                                                 |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>syngo</i> 3D Roadmap                 | superimposes 3D objects and/or the 3D color-coded image onto the live X-ray image for real-time guidance.                                                                                       |
| <i>syngo</i> Biplane 3D Roadmap         | superimposes 3D objects and/or the 3D color-coded image onto the live X-ray image in both planes for real-time guidance.                                                                        |
| <i>syngo</i> Fusion package             | allows preprocedural images from CT, MR, and PET-CT to be imported and fused to live fluoroscopy.                                                                                               |
| <i>syngo</i> 2D/3D Fusion               | requires just two fluoro projections to fuse 3D volumes from other imaging modalities for live image guidance.                                                                                  |
| <i>syngo</i> 3D/3D Fusion               | aligns two 3D volumes from the same or a different modality. Any <i>syngo</i> DynaCT or <i>syngo</i> Dyna3D image can be fused with datasets from CT, MR, or PET, for example.                  |
| <i>syngo</i> Toolbox                    | uses points and lines to interactively mark anatomical structures of interest in a 3D volume, e.g., a <i>syngo</i> DynaCT image. It is a generic application.                                   |
| <i>syngo</i> DynaCT SMART               | reduces metal artifacts through secondary reconstruction (streak metal artifacts reduction technology).                                                                                         |
| myNeedle Companion                      | allows the operator to plan a needle-based procedure in a 3D volume by specifying a target and multiple trajectories for needles or other rigid instruments.                                    |
| myEmbolization Guide                    | supports advanced embolization planning and guidance.                                                                                                                                           |
| <i>syngo</i> iFlow                      | allows the visualization of a DSA run in one single image through color-coding of time and intensity of maximum opacification.                                                                  |
| <i>syngo</i> 3D Stenosis Measurement    | allows analysis of a 3D segment within a 3D volume and calculates all data relevant to a stenosis (minimum diameter, maximum diameter, and cross-section of stenosis).                          |
| <i>syngo</i> 3D Aneurysm Guidance Neuro | enables segmentation of a cerebral aneurysm and its parent vessels with three simple mouse clicks, providing all data relevant to an aneurysm (parent vessel dimensions, volume, ostium, etc.). |
| <i>syngo</i> Neuro Virtual              | is a planning tool that provides a 3D virtual intracranial stent.                                                                                                                               |



For more information about ARTIS and Workplace applications in Interventional Neuroradiology visit us online at:

[siemens-healthineers.com/angio/interventional-radiology/clinical-specialities/neuroradiology](https://www.siemens-healthineers.com/angio/interventional-radiology/clinical-specialities/neuroradiology)

## Notes

[illegible]

[illegible]

## Notes

[illegible]

[illegible]

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States or other countries.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications and options described herein without prior notice. Please contact your local Siemens Healthineers sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory new components, we guarantee the quality of these recycled components.

#### **Note**

Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

#### **Caution**

Federal law restricts this device to sale by or on the order of a physician.

The statements by Siemens Healthineers' customers described herein are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.

---

### **Siemens Healthineers Headquarters**

Siemens Healthineers AG  
Siemensstr. 3  
91301 Forchheim, Germany  
Phone: +49 9191 180  
siemens-healthineers.com