

# Secondary Scaphoid Nonunion

Multitom Rax Real3D<sup>1</sup> Hi-Res clinical case  
University Hospital Wuerzburg, Germany



<sup>1</sup>Option



# Clinical background and indication for Multitom Rax Real3D<sup>1</sup> Hi-Res examination

## Patient

Male | \*1972 | BMI 29 kg/m<sup>2</sup>

## Anamnesis

Four months before the present examination, the patient fell on his right wrist suffering an isolated scapholunate ligament (SLL) tear. Due to carpal instability and initial osteoarthritis in the stylo scaphoid joint, SLL reconstruction (using the ECRB tendon) and transfixation were performed.

The patient now reports to the emergency department with increasing pain over the radial side of the proximal carpal row (especially in motion).

## Indication for Real3D<sup>1</sup> Hi-Res examination

Radiography displays dorsal dislocation of the corkscrew anchor in the lateral view (right side; arrow). While transfixation of the scapholunate gap appears satisfactory in AP view (left side), lateral view depicts dorsal extension of the lunate bone suggesting insufficient stability of the scapholunate compartment.



AP

Conventional X-ray examination

Lateral

Study ID 5aab940

The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

<sup>1</sup> Option

# Multitom Rax Real3D<sup>1</sup> Hi-Res Settings



## Settings for tableside scan with dedicated Metal Protocol

<i>Tube voltage</i>	116 kV
<i>Current time product</i>	223 mAs
<i>Dose area product</i>	315 $\mu\text{Gy}\cdot\text{m}^2$
<i>Calculated value for <math>\text{CTDI}_{\text{vol},32}</math></i>	4.7 mGy
<i>Scan time</i>	14 sec
<i>Number of projections</i>	318

## Reconstruction settings for sectional views

<i>Pixel size</i>	0.2 mm
<i>Reconstruction kernel</i>	very sharp (equivalent to Ur77)
<i>Slice thickness</i>	2 mm

The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

<sup>1</sup> Option

# Multitom Rax Real3D<sup>1</sup> Hi-Res

## Diagnostic findings

Rax Real3D images with metal artefact reduction also display widening of the scapholunate gap and dorsal intercalated segment instability (DISI) position of the central carpal column, confirming insufficient transfixation. Coronal image shows additional bone resorption in the scaphoid (arrow).



*Coronal view*



*Sagittal view*

Study ID 5aab562

The products/features (mentioned herein) are not commercially available in all countries.  
Their future availability cannot be guaranteed.


<sup>1</sup> Option

# Multitom Rax Real3D<sup>1</sup> Hi-Res


## Diagnostic findings

Sagittal and axial images display extensive resorption in the scaphoid with displacement of the k-wire (arrows). Image findings suggest a bony avulsion of the corkscrew anchor in the dorsal portion of the scaphoid after transfixation with subsequent nonunion.

Immediate surgical removal of the k-wire was performed with consecutive cast immobilization. With respect to the highly unstable carpus, different salvage procedures were discussed with the patient including proximal row carpectomy and four corner fusion.



*Axial view*



*Sagittal view*

Study ID 5aab562





*“Real3D images possess excellent spatial resolution with easily discernable bone microarchitecture. The metal implants do not cause considerable artefacts. Soft tissue contrast is high.”<sup>1</sup>*

**Jan-Peter Grunz, MD**  
University Hospital Wuerzburg, Germany

<sup>1</sup>The statements by Siemens Healthineers customers described herein are based on results that were achieved in the customer’s unique setting. Since there is no “typical” hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption) there can be no guarantee that other customers will achieve the same results.



The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

Results from case studies are not predictive of results in other cases. Results in other cases may vary.

Dr. Jan-Peter Grunz is employed by an institution that receives financial support from Siemens Healthineers for collaborations.