

### **Accompanying notes for whole-body MRI protocols on Magnetom World.**

- 6-station and 7-station options are included for Aera and Sola protocols.
- 7-station protocol is suitable for coverage from skull vertex to knees in most patients; additional stations may be added for taller patients or if coverage to feet is required.
- 5 stations are usually suitable for coverage from skull vertex to thighs; a sixth station has been included but can be deleted if not required.
- Sites may choose to omit stacked axial T2w Haste sequences in some scenarios (optional in MY-RADS and MET-RADS-P guidelines).
- Sites may choose to omit intermediate b-value ( $b=600 \text{ s/mm}^2$ ) in some scenarios (optional in MY-RADS and MET-RADS-P guidelines).
- Sites may choose to add small field-of-view imaging if required, as per local protocols.