

Table of contents
-------------------

\\USER_MR3
------------

<table border="1"> <tr> <td>MAMMAE</td> </tr> </table>	MAMMAE
MAMMAE	

<table border="1"> <tr> <td> <table border="1"> <tr> <td>MAMMAE</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>MAMMAE</td> </tr> </table>	MAMMAE
<table border="1"> <tr> <td>MAMMAE</td> </tr> </table>	MAMMAE	
MAMMAE		

<table border="1"> <tr> <td> <table border="1"> <tr> <td>KLINIEK/SCREENING</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>KLINIEK/SCREENING</td> </tr> </table>	KLINIEK/SCREENING
<table border="1"> <tr> <td>KLINIEK/SCREENING</td> </tr> </table>	KLINIEK/SCREENING	
KLINIEK/SCREENING		

<table border="1"> <tr> <td> <table border="1"> <tr> <td> <table border="1"> <tr> <td>localizer</td> <td>*</td> </tr> <tr> <td>tirm_tra_prothese_overview</td> <td>*</td> </tr> <tr> <td>tirm_ws_tra_silicone_prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>tirm_tra_NaCl-prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>Cal_t1_fl3d_tra_Dixon</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_pre</td> <td>*</td> </tr> <tr> <td>Perfusie_t1_twist</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t2_tse_dixon_fast_p3</td> <td>*</td> </tr> </table> </td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td> <table border="1"> <tr> <td>localizer</td> <td>*</td> </tr> <tr> <td>tirm_tra_prothese_overview</td> <td>*</td> </tr> <tr> <td>tirm_ws_tra_silicone_prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>tirm_tra_NaCl-prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>Cal_t1_fl3d_tra_Dixon</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_pre</td> <td>*</td> </tr> <tr> <td>Perfusie_t1_twist</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t2_tse_dixon_fast_p3</td> <td>*</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>localizer</td> <td>*</td> </tr> <tr> <td>tirm_tra_prothese_overview</td> <td>*</td> </tr> <tr> <td>tirm_ws_tra_silicone_prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>tirm_tra_NaCl-prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>Cal_t1_fl3d_tra_Dixon</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_pre</td> <td>*</td> </tr> <tr> <td>Perfusie_t1_twist</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t2_tse_dixon_fast_p3</td> <td>*</td> </tr> </table>	localizer	*	tirm_tra_prothese_overview	*	tirm_ws_tra_silicone_prothese	*	resolve_diff_tra_spair_p2	*	tirm_tra_NaCl-prothese	*	resolve_diff_tra_spair_p2	*	resolve_diff_tra_spair_p2	*	Cal_t1_fl3d_tra_Dixon	*	t1_fl3d_tra_Dixon_pre	*	Perfusie_t1_twist	*	t1_fl3d_tra_Dixon_post	*	t1_fl3d_tra_Dixon_post	*	t1_fl3d_tra_Dixon_post	*	t1_fl3d_tra_Dixon_post	*	t2_tse_dixon_fast_p3	*
<table border="1"> <tr> <td> <table border="1"> <tr> <td>localizer</td> <td>*</td> </tr> <tr> <td>tirm_tra_prothese_overview</td> <td>*</td> </tr> <tr> <td>tirm_ws_tra_silicone_prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>tirm_tra_NaCl-prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>Cal_t1_fl3d_tra_Dixon</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_pre</td> <td>*</td> </tr> <tr> <td>Perfusie_t1_twist</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t2_tse_dixon_fast_p3</td> <td>*</td> </tr> </table> </td> </tr> </table>	<table border="1"> <tr> <td>localizer</td> <td>*</td> </tr> <tr> <td>tirm_tra_prothese_overview</td> <td>*</td> </tr> <tr> <td>tirm_ws_tra_silicone_prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>tirm_tra_NaCl-prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>Cal_t1_fl3d_tra_Dixon</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_pre</td> <td>*</td> </tr> <tr> <td>Perfusie_t1_twist</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t2_tse_dixon_fast_p3</td> <td>*</td> </tr> </table>	localizer	*	tirm_tra_prothese_overview	*	tirm_ws_tra_silicone_prothese	*	resolve_diff_tra_spair_p2	*	tirm_tra_NaCl-prothese	*	resolve_diff_tra_spair_p2	*	resolve_diff_tra_spair_p2	*	Cal_t1_fl3d_tra_Dixon	*	t1_fl3d_tra_Dixon_pre	*	Perfusie_t1_twist	*	t1_fl3d_tra_Dixon_post	*	t1_fl3d_tra_Dixon_post	*	t1_fl3d_tra_Dixon_post	*	t1_fl3d_tra_Dixon_post	*	t2_tse_dixon_fast_p3	*	
<table border="1"> <tr> <td>localizer</td> <td>*</td> </tr> <tr> <td>tirm_tra_prothese_overview</td> <td>*</td> </tr> <tr> <td>tirm_ws_tra_silicone_prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>tirm_tra_NaCl-prothese</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>resolve_diff_tra_spair_p2</td> <td>*</td> </tr> <tr> <td>Cal_t1_fl3d_tra_Dixon</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_pre</td> <td>*</td> </tr> <tr> <td>Perfusie_t1_twist</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t1_fl3d_tra_Dixon_post</td> <td>*</td> </tr> <tr> <td>t2_tse_dixon_fast_p3</td> <td>*</td> </tr> </table>	localizer	*	tirm_tra_prothese_overview	*	tirm_ws_tra_silicone_prothese	*	resolve_diff_tra_spair_p2	*	tirm_tra_NaCl-prothese	*	resolve_diff_tra_spair_p2	*	resolve_diff_tra_spair_p2	*	Cal_t1_fl3d_tra_Dixon	*	t1_fl3d_tra_Dixon_pre	*	Perfusie_t1_twist	*	t1_fl3d_tra_Dixon_post	*	t1_fl3d_tra_Dixon_post	*	t1_fl3d_tra_Dixon_post	*	t1_fl3d_tra_Dixon_post	*	t2_tse_dixon_fast_p3	*		
localizer	*																															
tirm_tra_prothese_overview	*																															
tirm_ws_tra_silicone_prothese	*																															
resolve_diff_tra_spair_p2	*																															
tirm_tra_NaCl-prothese	*																															
resolve_diff_tra_spair_p2	*																															
resolve_diff_tra_spair_p2	*																															
Cal_t1_fl3d_tra_Dixon	*																															
t1_fl3d_tra_Dixon_pre	*																															
Perfusie_t1_twist	*																															
t1_fl3d_tra_Dixon_post	*																															
t1_fl3d_tra_Dixon_post	*																															
t1_fl3d_tra_Dixon_post	*																															
t1_fl3d_tra_Dixon_post	*																															
t2_tse_dixon_fast_p3	*																															

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\localizer \*

TA: 8.6 s PM: REF Voxel size: 1.3×1.3×7.0 mmPAT: Off Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	L0.0 A3.9 H12.5 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Slice group	2
Slices	1
Dist. factor	20 %
Position	R90.0 P36.1 H12.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	3
Slices	1
Dist. factor	20 %
Position	L90.0 P36.1 H12.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	4
Slices	1
Dist. factor	20 %
Position	L0.0 P36.1 H12.5 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	400 mm
FoV phase	100,0 %
Slice thickness	7,0 mm
TR	7,4 ms
TE	3,01 ms
Averages	1
Concatenations	4
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Contrast - Common**

TR	7,4 ms
TE	3,01 ms
TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Common**

SWI	Off
-----	-----

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	400 mm
FoV phase	100,0 %
Slice thickness	7,0 mm
Base resolution	320
Phase resolution	75 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	L0.0 A3.9 H12.5 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Slice group	2
Slices	1
Dist. factor	20 %
Position	R90.0 P36.1 H12.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	3
Slices	1
Dist. factor	20 %
Position	L90.0 P36.1 H12.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	4
Slices	1
Dist. factor	20 %
Position	L0.0 P36.1 H12.5 mm
Orientation	Transversal
Phase enc. dir.	A >> P

**Geometry - Common**

FoV read	400 mm
FoV phase	100,0 %
Slice thickness	7,0 mm
TR	7,4 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	4

**Geometry - AutoAlign**

Slice group	1
Position	L0.0 A3.9 H12.5 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Slice group	2
Position	R90.0 P36.1 H12.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	3
Position	L90.0 P36.1 H12.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	4
Position	L0.0 P36.1 H12.5 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	0,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off

**System - Adjustments**

Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0,00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	7,4 ms
Concatenations	4
Segments	1

**Physio - Cardiac**

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	400 mm
FoV phase	100,0 %
Phase resolution	75 %

**Physio - PACE**

Resp. control	Off
Concatenations	4

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off

**Inline - Soft Tissue**

MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	20 deg
Measurements	1
Contrasts	1
TR	7,4 ms
TE	3,01 ms

**Sequence - Part 1**

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Sequential
Bandwidth	320 Hz/Px

**Sequence - Part 2**

Segments	1
Acoustic noise reduction	None
RF pulse type	Fast
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\tirm\_tra\_prothese\_overview \*

TA: 0:30 PM: REF Voxel size: 1.1×1.1×6.0 mmPAT: 2 Rel. SNR: 1.00 : tir

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	5
Dist. factor	100 %
Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	19 %
FoV read	340 mm
FoV phase	100,0 %
Slice thickness	6,0 mm
TR	1800,0 ms
TE	53,0 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Contrast - Common**

TR	1800,0 ms
TE	53,0 ms
MTC	Off
Magn. preparation	Slice-sel. IR
T1	230 ms
Flip angle	120 deg
Fat suppr.	None
Water suppr.	Water sat.
Restore magn.	Off
Freeze suppressed tissue	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	340 mm
FoV phase	100,0 %
Slice thickness	6,0 mm
Base resolution	320
Phase resolution	75 %
Phase partial Fourier	Off
Trajectory	Cartesian

**Resolution - Common**

Interpolation	Off
---------------	-----

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	43
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slice group	1
Slices	5
Dist. factor	100 %
Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
FoV read	340 mm
FoV phase	100,0 %
Slice thickness	6,0 mm
TR	1800,0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	0,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
Water suppr.	Water sat.
Restore magn.	Off
Special sat.	None

**Geometry - Navigator**

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	On
Only after freq. change	Off
Assume Dominant Fat	Off
Assume Silicone	On
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Rotation	90,00 deg
R >> L	340 mm
A >> P	340 mm
F >> H	54 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1800,0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Slice-sel. IR
TI	230 ms
Fat suppr.	None
Dark blood	Off
FoV read	340 mm
FoV phase	100,0 %
Phase resolution	75 %
Trajectory	Cartesian

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Reduce Motion Sens.	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	8,86 ms
Bandwidth	300 Hz/Px

**Sequence - Part 2**

Define	Turbo factor
Echo trains per slice	15
Phase correction	Automatic
Acoustic noise reduction	None
RF pulse type	Low SAR
Gradient mode	Normal
Hyperecho	Off
WARP	Off
Red. EC sensitivity	Off
Turbo factor	11

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\tirm\_ws\_tra\_silicone\_prothese \*

TA: 2:53 PM: FIX Voxel size: 0.9×0.9×4.0 mmPAT: 2 Rel. SNR: 1.00 : tir

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	32
Dist. factor	20 %
Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	29 %
FoV read	340 mm
FoV phase	100,0 %
Slice thickness	4,0 mm
TR	3900,0 ms
TE	64,0 ms
Averages	1
Concatenations	2
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Contrast - Common**

TR	3900,0 ms
TE	64,0 ms
TD	0,0 ms
MTC	Off
Magn. preparation	Slice-sel. IR
T1	230 ms
Flip angle	80 deg
Fat suppr.	None
Water suppr.	Water sat.
Restore magn.	Off
Freeze suppressed tissue	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	340 mm
FoV phase	100,0 %
Slice thickness	4,0 mm
Base resolution	384
Phase resolution	85 %
Phase partial Fourier	Off

**Resolution - Common**

Trajectory	Cartesian
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	39
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slice group	1
Slices	32
Dist. factor	20 %
Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
FoV read	340 mm
FoV phase	100,0 %
Slice thickness	4,0 mm
TR	3900,0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	2

**Geometry - AutoAlign**

Slice group	1
Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	0,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
Water suppr.	Water sat.
Restore magn.	Off
Special sat.	None

**Geometry - Navigator**

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	On
Only after freq. change	Off
Assume Dominant Fat	Off
Assume Silicone	On
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Rotation	90,00 deg
R >> L	340 mm
A >> P	340 mm
F >> H	153 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3900,0 ms
Concatenations	2

**Physio - Cardiac**

Magn. preparation	Slice-sel. IR
TI	230 ms
Fat suppr.	None
Dark blood	Off
FoV read	340 mm
FoV phase	100,0 %
Phase resolution	85 %
Trajectory	Cartesian

**Physio - PACE**

Resp. control	Off
Concatenations	2

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Reduce Motion Sens.	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	9,08 ms
Bandwidth	303 Hz/Px

**Sequence - Part 2**

Define	Turbo factor
Echo trains per slice	21
Phase correction	Automatic
Acoustic noise reduction	None
RF pulse type	Low SAR
Gradient mode	Normal
Hyperecho	On
WARP	Off
Red. EC sensitivity	Off
Turbo factor	11

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s



\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\resolve\_diff\_tra\_spair\_p2 \*

TA: 2:53 PM: FIX Voxel size: 1.5×1.5×4.0 mmPAT: 2 Rel. SNR: 1.00 : resolve

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	47
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off

**Resolution - Filter Rawdata**

Raw filter	Off
------------	-----

**Routine**

Slice group	1
Slices	28
Dist. factor	20 %
Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	340 mm
FoV phase	50,0 %
Slice thickness	4,0 mm
TR	5190 ms
TE 1	46 ms
TE 2	74 ms
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	7BL;7BR;AXL;AXR

**Geometry - Common**

Slice group	1
Slices	28
Dist. factor	20 %
Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	340 mm
FoV phase	50,0 %
Slice thickness	4,0 mm
TR	5190 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R0.8 A2.2 H0.0
R	0,8 mm
A	2,2 mm
H	0,0 mm
Initial Rotation	-0,80 deg
Initial Orientation	Transversal

**Contrast - Common**

TR	5190 ms
TE 1	46 ms
TE 2	74 ms
MTC	Off
Magn. preparation	None
Flip angle	180 deg
Fat suppr.	SPAIR
Fat sat. mode	Strong

**Contrast - Dynamic**

Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Resolution - Common**

FoV read	340 mm
FoV phase	50,0 %
Slice thickness	4,0 mm
Base resolution	232
Phase resolution	80 %
Phase partial Fourier	Off
Readout partial Fourier	5/8
Readout segments	5
Interpolation	Off

**Geometry - Saturation**

Sat. region	1
Thickness	110 mm
Position	R2.5 P121.1 H0.0 mm
Orientation	C > S1.2
Sat. region	2
Thickness	110 mm
Position	L3.2 P180.6 H0.0 mm
Orientation	C > S-1.0
Fat sat. mode	Strong
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	On
Only after freq. change	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Rotation	-0,80 deg
A >> P	170 mm
R >> L	340 mm
F >> H	134 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	8,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	5190 ms
Concatenations	1

**Diff - Neuro**

Diffusion mode	3-Scan Trace
Diff. directions	3
Diffusion Scheme	Monopolar
Diff. weightings	2
b-value 1	0 s/mm <sup>2</sup>
b-value 2	800 s/mm <sup>2</sup>
b-value 1	1
b-value 2	3
Diff. weighted images	Off
Trace weighted images	On
ADC maps	On
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	5

**Diff - Body**

Diffusion mode	3-Scan Trace
Diff. directions	3
Diffusion Scheme	Monopolar
Diff. weightings	2
b-value 1	0 s/mm <sup>2</sup>
b-value 2	800 s/mm <sup>2</sup>
b-value 1	1
b-value 2	3
Diff. weighted images	Off
Trace weighted images	On
ADC maps	On
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	On
Calculated bValue	1400 s/mm <sup>2</sup>
b-Value >=	0 s/mm <sup>2</sup>
Noise level	5

**Diff - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Contrasts	2
Optimization	Min. TE
Multi-slice mode	Interleaved
Echo spacing	0,38 ms
Bandwidth	829 Hz/Px

**Sequence - Part 2**

EPI factor	47
RF pulse type	Low SAR
Gradient mode	Fast
Reacquisition mode	Off

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\tirm\_tra\_NaCl-prothese \*

TA: 2:53 PM: FIX Voxel size: 0.9×0.9×4.0 mmPAT: 2 Rel. SNR: 1.00 : tir

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	32
Dist. factor	20 %
Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	29 %
FoV read	340 mm
FoV phase	100,0 %
Slice thickness	4,0 mm
TR	3900,0 ms
TE	64,0 ms
Averages	1
Concatenations	2
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Contrast - Common**

TR	3900,0 ms
TE	64,0 ms
TD	0,0 ms
MTC	Off
Magn. preparation	Slice-sel. IR
TI	230 ms
Flip angle	80 deg
Fat suppr.	None
Water suppr.	None
Restore magn.	Off
Freeze suppressed tissue	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	340 mm
FoV phase	100,0 %
Slice thickness	4,0 mm
Base resolution	384
Phase resolution	85 %
Phase partial Fourier	Off

**Resolution - Common**

Trajectory	Cartesian
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	39
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slice group	1
Slices	32
Dist. factor	20 %
Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
FoV read	340 mm
FoV phase	100,0 %
Slice thickness	4,0 mm
TR	3900,0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	2

**Geometry - AutoAlign**

Slice group	1
Position	L0.0 P50.0 H0.0 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	0,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
Water suppr.	None
Restore magn.	Off
Special sat.	None

**Geometry - Navigator**

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0,00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3900,0 ms
Concatenations	2

**Physio - Cardiac**

Magn. preparation	Slice-sel. IR
T1	230 ms
Fat suppr.	None
Dark blood	Off
FoV read	340 mm
FoV phase	100,0 %
Phase resolution	85 %
Trajectory	Cartesian

**Physio - PACE**

Resp. control	Off
Concatenations	2

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Reduce Motion Sens.	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	9,08 ms
Bandwidth	303 Hz/Px

**Sequence - Part 2**

Define	Turbo factor
Echo trains per slice	21
Phase correction	Automatic
Acoustic noise reduction	None
RF pulse type	Low SAR
Gradient mode	Normal
Hyperecho	On
WARP	Off
Red. EC sensitivity	Off
Turbo factor	11

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\resolve\_diff\_tra\_spair\_p2 \*

TA: 2:53 PM: FIX Voxel size: 1.5×1.5×4.0 mmPAT: 2 Rel. SNR: 1.00 : resolve

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	47
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off

**Resolution - Filter Rawdata**

Raw filter	Off
------------	-----

**Routine**

Slice group	1
Slices	28
Dist. factor	20 %
Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	340 mm
FoV phase	50,0 %
Slice thickness	4,0 mm
TR	5190 ms
TE 1	46 ms
TE 2	74 ms
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	7BL;7BR;AXL;AXR

**Geometry - Common**

Slice group	1
Slices	28
Dist. factor	20 %
Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	340 mm
FoV phase	50,0 %
Slice thickness	4,0 mm
TR	5190 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R0.8 A2.2 H0.0
R	0,8 mm
A	2,2 mm
H	0,0 mm
Initial Rotation	-0,80 deg
Initial Orientation	Transversal

**Contrast - Common**

TR	5190 ms
TE 1	46 ms
TE 2	74 ms
MTC	Off
Magn. preparation	None
Flip angle	180 deg
Fat suppr.	SPAIR
Fat sat. mode	Strong

**Contrast - Dynamic**

Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Resolution - Common**

FoV read	340 mm
FoV phase	50,0 %
Slice thickness	4,0 mm
Base resolution	232
Phase resolution	80 %
Phase partial Fourier	Off
Readout partial Fourier	5/8
Readout segments	5
Interpolation	Off

**Geometry - Saturation**

Sat. region	1
Thickness	110 mm
Position	R2.5 P121.1 H0.0 mm
Orientation	C > S1.2
Sat. region	2
Thickness	110 mm
Position	L3.2 P180.6 H0.0 mm
Orientation	C > S-1.0
Fat sat. mode	Strong
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	On
Only after freq. change	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Rotation	-0,80 deg
A >> P	170 mm
R >> L	340 mm
F >> H	134 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	8,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	5190 ms
Concatenations	1

**Diff - Neuro**

Diffusion mode	3-Scan Trace
Diff. directions	3
Diffusion Scheme	Monopolar
Diff. weightings	2
b-value 1	0 s/mm <sup>2</sup>
b-value 2	800 s/mm <sup>2</sup>
b-value 1	1
b-value 2	3
Diff. weighted images	Off
Trace weighted images	On
ADC maps	On
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	5

**Diff - Body**

Diffusion mode	3-Scan Trace
Diff. directions	3
Diffusion Scheme	Monopolar
Diff. weightings	2
b-value 1	0 s/mm <sup>2</sup>
b-value 2	800 s/mm <sup>2</sup>
b-value 1	1
b-value 2	3
Diff. weighted images	Off
Trace weighted images	On
ADC maps	On
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	On
Calculated bValue	1400 s/mm <sup>2</sup>
b-Value >=	0 s/mm <sup>2</sup>
Noise level	5

**Diff - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Contrasts	2
Optimization	Min. TE
Multi-slice mode	Interleaved
Echo spacing	0,38 ms
Bandwidth	829 Hz/Px

**Sequence - Part 2**

EPI factor	47
RF pulse type	Low SAR
Gradient mode	Fast
Reacquisition mode	Off

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\resolve\_diff\_tra\_spair\_p2 \*

TA: 4:07 PM: FIX Voxel size: 1.5×1.5×4.0 mmPAT: 2 Rel. SNR: 1.00 : resolve

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	49
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off

**Resolution - Filter Rawdata**

Raw filter	Off
------------	-----

**Routine**

Slice group	1
Slices	28
Dist. factor	20 %
Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	360 mm
FoV phase	40,5 %
Slice thickness	4,0 mm
TR	5110 ms
TE 1	47 ms
TE 2	75 ms
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	7BL;7BR;AXL;AXR

**Geometry - Common**

Slice group	1
Slices	28
Dist. factor	20 %
Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	360 mm
FoV phase	40,5 %
Slice thickness	4,0 mm
TR	5110 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R0.8 A2.2 H0.0
R	0,8 mm
A	2,2 mm
H	0,0 mm
Initial Rotation	-0,80 deg
Initial Orientation	Transversal

**Contrast - Common**

TR	5110 ms
TE 1	47 ms
TE 2	75 ms
MTC	Off
Magn. preparation	None
Flip angle	180 deg
Fat suppr.	SPAIR
Fat sat. mode	Strong

**Contrast - Dynamic**

Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Resolution - Common**

FoV read	360 mm
FoV phase	40,5 %
Slice thickness	4,0 mm
Base resolution	242
Phase resolution	100 %
Phase partial Fourier	Off
Readout partial Fourier	5/8
Readout segments	5
Interpolation	Off

**Geometry - Saturation**

Sat. region	1
Thickness	110 mm
Position	R2.5 P121.1 H0.0 mm
Orientation	C > S1.2
Sat. region	2
Thickness	110 mm
Position	L3.2 P180.6 H0.0 mm
Orientation	C > S-1.0
Fat sat. mode	Strong
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - All

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	On
Only after freq. change	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R0.8 A2.2 H0.0 mm
Orientation	Transversal
Rotation	-0,80 deg
A >> P	146 mm
R >> L	360 mm
F >> H	134 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	8,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	5110 ms
Concatenations	1

**Diff - Neuro**

Diffusion mode	3-Scan Trace
Diff. directions	3
Diffusion Scheme	Monopolar
Diff. weightings	2
b-value 1	50 s/mm <sup>2</sup>
b-value 2	800 s/mm <sup>2</sup>
b-value 1	1
b-value 2	4
Diff. weighted images	Off
Trace weighted images	On
ADC maps	On
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	5

**Diff - Body**

Diffusion mode	3-Scan Trace
Diff. directions	3
Diffusion Scheme	Monopolar
Diff. weightings	2
b-value 1	50 s/mm <sup>2</sup>
b-value 2	800 s/mm <sup>2</sup>
b-value 1	1
b-value 2	4
Diff. weighted images	Off
Trace weighted images	On
ADC maps	On
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	On
Calculated bValue	1400 s/mm <sup>2</sup>
b-Value >=	0 s/mm <sup>2</sup>
Noise level	5

**Diff - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Contrasts	2
Optimization	Min. TE
Multi-slice mode	Interleaved
Echo spacing	0,38 ms
Bandwidth	861 Hz/Px

**Sequence - Part 2**

EPI factor	49
RF pulse type	Low SAR
Gradient mode	Fast
Reacquisition mode	Off

**Sequence - Assistant**

Mode	Off
------	-----



\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\Cal\_t1\_f13d\_tra\_Dixon \*

TA: 1:03 PM: FIX Voxel size: 0.9×0.9×1.0 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Resolution - Common**

Phase partial Fourier	7/8
Slice partial Fourier	6/8
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	20 %
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On
POCS	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Contrast - Common**

TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Flip angle	5,0 deg
Fat suppr.	None
Water suppr.	None
Dixon	On
Dixon evaluation	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
Base resolution	416
Phase resolution	86 %
Slice resolution	55 %

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	90,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
------------	------

**Geometry - Saturation**

Water suppr.	None
Dixon	On
Dixon evaluation	Off
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	90,00 deg
R >> L	360 mm
A >> P	360 mm
F >> H	144 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

View sharing	Off
Flip angle	5,0 deg
Measurements	1

**Inline - Common**

Burn time-to-center	Off
Temporal interpolation	1
3D centric reordering	Off
	null

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	5,0 deg
Measurements	1
Contrasts	2
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reduce Motion Sens.	On
Asymmetric echo	Weak
Contrasts	2
Readout mode	Bipolar
Optimization	In phase
Multi-slice mode	Sequential
Bandwidth 1	860 Hz/Px
Bandwidth 2	860 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast*
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	60 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\t1\_fl3d\_tra\_Dixon\_pre \*

TA: 1:03 PM: FIX Voxel size: 0.9×0.9×1.0 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Resolution - Common**

Phase partial Fourier	7/8
Slice partial Fourier	6/8
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	20 %
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On
POCS	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Contrast - Common**

TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Flip angle	15,0 deg
Fat suppr.	None
Water suppr.	None
Dixon	On
Dixon evaluation	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
Base resolution	416
Phase resolution	86 %
Slice resolution	55 %

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	90,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
------------	------

**Geometry - Saturation**

Water suppr.	None
Dixon	On
Dixon evaluation	Off
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	90,00 deg
R >> L	360 mm
A >> P	360 mm
F >> H	144 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

View sharing	Off
Flip angle	15,0 deg
Measurements	1

**Inline - Common**

Burn time-to-center	Off
Temporal interpolation	1
3D centric reordering	Off
	null

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	15,0 deg
Measurements	1
Contrasts	2
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reduce Motion Sens.	On
Asymmetric echo	Weak
Contrasts	2
Readout mode	Bipolar
Optimization	In phase
Multi-slice mode	Sequential
Bandwidth 1	860 Hz/Px
Bandwidth 2	860 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast*
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	60 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\Perfusie\_t1\_twist \*

TA: 2:03 PM: FIX Voxel size: 0.9×0.9×2.5 mmPAT: 3 Rel. SNR: 1.00 : fldyn

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	On
Auto close inline display	On
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	20 %
Slice oversampling	0,0 %
Slices per slab	60
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	2,50 mm
TR	4,09 ms
TE	2,06 ms
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Contrast - Common**

TR	4,09 ms
TE	2,06 ms
Flip angle	20 deg
Fat suppr.	None

**Contrast - Dynamic**

Reconstruction	Magnitude
Measurements	20
Pause after meas. 1	0,0 s
Pause after meas. 2	0,0 s
Pause after meas. 3	0,0 s
Pause after meas. 4	0,0 s
Pause after meas. 5	0,0 s
Pause after meas. 6	0,0 s
Pause after meas. 7	0,0 s
Pause after meas. 8	0,0 s
Pause after meas. 9	0,0 s
Pause after meas. 10	0,0 s
Pause after meas. 11	0,0 s
Pause after meas. 12	0,0 s
Pause after meas. 13	0,0 s
Pause after meas. 14	0,0 s
Pause after meas. 15	0,0 s
Pause after meas. 16	0,0 s
Pause after meas. 17	0,0 s
Pause after meas. 18	0,0 s
Pause after meas. 19	0,0 s

**Contrast - Dynamic**

Multiple series	Each measurement
-----------------	------------------

**Resolution - Common**

FoV read	360 mm
FoV phase	100,0 %
Slice thickness	2,50 mm
Base resolution	384
Phase resolution	90 %
Slice resolution	77 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
View sharing	TWIST
Central region A	15 %
Sampling density B	10 %
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Slice oversampling	0,0 %
Slices per slab	60
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	2,50 mm
TR	4,09 ms

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm

**Geometry - AutoAlign**

Initial Rotation	90,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Maximum

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0,00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Angio - Common**

View sharing	TWIST
Central region A	15 %
Sampling density B	10 %
Dynamic recon mode	Symmetric Share
Flip angle	20 deg

**Angio - Common**

Temporal resolution	4,57 s
Measurements	20
Pause after meas. 1	0,0 s
Pause after meas. 2	0,0 s
Pause after meas. 3	0,0 s
Pause after meas. 4	0,0 s
Pause after meas. 5	0,0 s
Pause after meas. 6	0,0 s
Pause after meas. 7	0,0 s
Pause after meas. 8	0,0 s
Pause after meas. 9	0,0 s
Pause after meas. 10	0,0 s
Pause after meas. 11	0,0 s
Pause after meas. 12	0,0 s
Pause after meas. 13	0,0 s
Pause after meas. 14	0,0 s
Pause after meas. 15	0,0 s
Pause after meas. 16	0,0 s
Pause after meas. 17	0,0 s
Pause after meas. 18	0,0 s
Pause after meas. 19	0,0 s
Burn time-to-center	On
Temporal interpolation	1
Time to center	18,1 s

**Angio - Inline**

Subtract	On
Save images	Off
Autoscaling	Off
Scaling factor	1
Offset	50
Subtrahend	1
Subtraction indices	
Subtraction group	1
Measurements	20
StdDev	Off
Save original images	On

**Angio - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	On
MIP-Time	On
Save original images	On

**Angio - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Off
Optimization	Min. TE TR
Bandwidth	540 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Phase Enc. Rewinder	On



**Sequence - Assistant**

Mode	Off
Allowed delay	0 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\t1\_fl3d\_tra\_Dixon\_post \*

TA: 1:03 PM: FIX Voxel size: 0.9×0.9×1.0 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Resolution - Common**

Phase partial Fourier	7/8
Slice partial Fourier	6/8
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	20 %
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On
POCS	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Contrast - Common**

TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Flip angle	15,0 deg
Fat suppr.	None
Water suppr.	None
Dixon	On
Dixon evaluation	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
Base resolution	416
Phase resolution	86 %
Slice resolution	55 %

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	90,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
------------	------

**Geometry - Saturation**

Water suppr.	None
Dixon	On
Dixon evaluation	Off
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	90,00 deg
R >> L	360 mm
A >> P	360 mm
F >> H	144 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

View sharing	Off
Flip angle	15,0 deg
Measurements	1

**Inline - Common**

Burn time-to-center	Off
Temporal interpolation	1
3D centric reordering	Off
	null

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	15,0 deg
Measurements	1
Contrasts	2
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reduce Motion Sens.	On
Asymmetric echo	Weak
Contrasts	2
Readout mode	Bipolar
Optimization	In phase
Multi-slice mode	Sequential
Bandwidth 1	860 Hz/Px
Bandwidth 2	860 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast*
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	60 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\t1\_fl3d\_tra\_Dixon\_post \*

TA: 1:03 PM: FIX Voxel size: 0.9×0.9×1.0 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	20 %
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Contrast - Common**

TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Flip angle	15,0 deg
Fat suppr.	None
Water suppr.	None
Dixon	On
Dixon evaluation	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
Base resolution	416
Phase resolution	86 %
Slice resolution	55 %

**Resolution - Common**

Phase partial Fourier	7/8
Slice partial Fourier	6/8
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On
POCS	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	90,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
------------	------

**Geometry - Saturation**

Water suppr.	None
Dixon	On
Dixon evaluation	Off
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	90,00 deg
R >> L	360 mm
A >> P	360 mm
F >> H	144 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

View sharing	Off
Flip angle	15,0 deg
Measurements	1

**Inline - Common**

Burn time-to-center	Off
Temporal interpolation	1
3D centric reordering	Off
	null

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	15,0 deg
Measurements	1
Contrasts	2
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reduce Motion Sens.	On
Asymmetric echo	Weak
Contrasts	2
Readout mode	Bipolar
Optimization	In phase
Multi-slice mode	Sequential
Bandwidth 1	860 Hz/Px
Bandwidth 2	860 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast*
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	60 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\t1\_fl3d\_tra\_Dixon\_post \*

TA: 1:03 PM: FIX Voxel size: 0.9×0.9×1.0 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	20 %
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Contrast - Common**

TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Flip angle	15,0 deg
Fat suppr.	None
Water suppr.	None
Dixon	On
Dixon evaluation	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
Base resolution	416
Phase resolution	86 %
Slice resolution	55 %

**Resolution - Common**

Phase partial Fourier	7/8
Slice partial Fourier	6/8
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On
POCS	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	90,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
------------	------



**Geometry - Saturation**

Water suppr.	None
Dixon	On
Dixon evaluation	Off
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	90,00 deg
R >> L	360 mm
A >> P	360 mm
F >> H	144 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

View sharing	Off
Flip angle	15,0 deg
Measurements	1

**Inline - Common**

Burn time-to-center	Off
Temporal interpolation	1
3D centric reordering	Off
	null

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	15,0 deg
Measurements	1
Contrasts	2
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reduce Motion Sens.	On
Asymmetric echo	Weak
Contrasts	2
Readout mode	Bipolar
Optimization	In phase
Multi-slice mode	Sequential
Bandwidth 1	860 Hz/Px
Bandwidth 2	860 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast*
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	60 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\t1\_fl3d\_tra\_Dixon\_post \*

TA: 1:03 PM: FIX Voxel size: 0.9×0.9×1.0 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	20 %
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	7BL;7BR;AXL;AXR

**Contrast - Common**

TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms
Flip angle	15,0 deg
Fat suppr.	None
Water suppr.	None
Dixon	On
Dixon evaluation	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
Base resolution	416
Phase resolution	86 %
Slice resolution	55 %

**Resolution - Common**

Phase partial Fourier	7/8
Slice partial Fourier	6/8
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On
POCS	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Slice oversampling	22,2 %
Slices per slab	144
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	1,0 mm
TR	5,68 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	90,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
------------	------

**Geometry - Saturation**

Water suppr.	None
Dixon	On
Dixon evaluation	Off
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	90,00 deg
R >> L	360 mm
A >> P	360 mm
F >> H	144 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

View sharing	Off
Flip angle	15,0 deg
Measurements	1

**Inline - Common**

Burn time-to-center	Off
Temporal interpolation	1
3D centric reordering	Off
	null

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	15,0 deg
Measurements	1
Contrasts	2
TR	5,68 ms
TE 1	2,46 ms
TE 2	3,69 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reduce Motion Sens.	On
Asymmetric echo	Weak
Contrasts	2
Readout mode	Bipolar
Optimization	In phase
Multi-slice mode	Sequential
Bandwidth 1	860 Hz/Px
Bandwidth 2	860 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast*
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	60 s

\\USER\_MR3\MAMMAE\MAMMAE\KLINIEK\SCREENING\t2\_tse\_dixon\_fast\_p3 \*

TA: 2:23 PM: FIX Voxel size: 0.9×0.9×2.5 mmPAT: 3 Rel. SNR: 1.00 : tse

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	70 %
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	2,5 mm
TR	5110,0 ms
TE	81 ms
Averages	1
Concatenations	2
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter, Image Filter
Coil elements	7BL;7BR;AXL;AXR

**Contrast - Common**

TR	5110,0 ms
TE	81 ms
TD	0,0 ms
MTC	Off
Magn. preparation	None
Flip angle	120 deg
Fat sat. mode	Strong
Water suppr.	None
Dixon	On
Fast Dixon	On
Restore magn.	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	360 mm
FoV phase	100,0 %
Slice thickness	2,5 mm
Base resolution	384
Phase resolution	80 %

**Resolution - Common**

Phase partial Fourier	Allowed
Trajectory	Cartesian
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	29
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	On
Intensity	Medium
Edge Enhancement	3
Smoothing	3
Unfiltered images	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slice group	1
Slices	60
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
FoV read	360 mm
FoV phase	100,0 %
Slice thickness	2,5 mm
TR	5110,0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	2

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0,0 mm
P	0,0 mm
H	0,0 mm
Initial Rotation	90,00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat sat. mode	Strong
Water suppr.	None
Dixon	On

**Geometry - Saturation**

Fast Dixon	On
Restore magn.	Off
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	On
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	90,00 deg
R >> L	360 mm
A >> P	360 mm
F >> H	150 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123,235025 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1,000
Reset	Off
? Ref. amplitude 1H	0,000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	5110,0 ms
Concatenations	2

**Physio - Cardiac**

Magn. preparation	None
Dark blood	Off

**Physio - Cardiac**

FoV read	360 mm
FoV phase	100,0 %
Phase resolution	80 %
Trajectory	Cartesian

**Physio - PACE**

Resp. control	Off
Concatenations	2

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	Off
Dimension	2D
Reduce Motion Sens.	On
Asymmetric echo	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Echo spacing	11,6 ms
Bandwidth	930 Hz/Px

**Sequence - Part 2**

Define	Turbo factor
Echo trains per slice	13
Phase correction	Automatic
RF pulse type	Normal
Gradient mode	Fast
Hyperecho	On
Red. EC sensitivity	Off
Turbo factor	15

**Sequence - Assistant**

Mode	Off
Allowed delay	30 s