



MOBILETT Elara Max

Stands out from the crowd in mobile X-ray imaging

Member of the MAX family

siemens-healthineers.us/mobilett-elara-max



Using the advantages of mobile digital radiography

The healthcare market is moving away from analog X-ray and computed radiography, toward fully digital X-ray systems (DR). This trend is not just about the many advantages that DR offers in terms of enhanced efficiency, handling, and image quality—it ultimately leads to new requirements and expectations like declining reimbursements for film-based X-rays and computed radiography, as seen, for example, in the U.S.


Embracing digital radiography...

Since 2017, Medicare reimbursements of the technical component have been reduced by 20% to providers submitting claims for film-based X-rays.¹



...and mobile solutions

Medical practitioners are also rapidly shifting their preferences to mobile X-ray equipment because of the advantages it offers over stationary systems.²



Global X-ray equipment market share, by product type (USD Mn)²

¹Massat, MB. *Changes in reimbursement to propel digital radiography*. Available from: <https://www.appliedradiology.com/articles/changes-in-reimbursement-to-propel-digital-radiography> [Accessed August 2018].

²Credence Research. *X-Ray Equipment Market By Product [...], Technology [...], Application [...]* - Growth, Share, Opportunities & Competitive Analysis, 2016–2022. Available from: <http://www.credenceresearch.com/report/x-ray-equipment-market> [Accessed August 2018].



Experience the MAX effect:
Multiply your gains in standardization,
savings, and satisfaction with
MOBILETT Elara Max.



Transform care delivery with MOBILETT Elara Max

Faced with the ongoing trend toward increasing efficiency and growing cost pressure, you as a healthcare provider might be looking for ways to optimize operations and improve productivity—without compromising on quality of care.

MOBILETT Elara Max supports you in this challenge: The mobile X-ray system delivers exceptional results

and offers unique features for enhancing efficiency and quality—and combines this with our extraordinary MAX advantages. Its easy maneuverability and flexible positioning improve productivity and enable consistently high image quality at the patient's bedside for excellent diagnostic certainty. Seamless and secure integration gives you access to relevant data

where and when you need it. All this makes MOBILETT Elara Max a smart, compact, and powerful X-ray system that's always ready to go.

MOBILETT Elara Max helps you to transform care delivery, and makes a real and lasting difference to your entire department.

Discover MOBILETT Elara Max

At a glance	4
Drive efficiency and image quality in your department	6
Accelerate your mobile workflow through seamless and secure integration	8
The MAX effect	10
Clinical images	12
Additional products and services	13
Technical specifications	14

MOBILETT Elara Max

At a glance

Smart design, big effect

Easy maneuverability and flexible positioning allow you to confidently produce high-quality images with ease.

Sleek ergonomic design for easy navigation

The ergonomic design has an antimicrobial coating and fully integrated cables for an unobstructed front view that allows for smooth, safe navigation.

Advanced workflows, greater security

Get access to relevant data, when and where you need it by integrating seamlessly and securely into your IT infrastructure.





The MAX effect

With two or more MAX systems, gain enhanced standardization, savings, and satisfaction with a common user interface and interchangeable detectors.

Keep your IT secure

Prevent cyberattacks and data theft with technology that safeguards your IT environment and protects sensitive patient data.

Personalize it when it matters

Low-dose organ programs allow you to image your youngest patients quickly and easily.

Drive efficiency and image quality in your department

Busy hallways and small patient rooms can be a challenge for mobile X-ray systems. With its exceptional maneuverability and flexible MAXreach arm, MOBILETT Elara Max helps your staff work confidently and with ease so they can produce consistently high image quality—even when all batteries are empty.



Maneuver effortlessly and ergonomically

MOBILETT Elara Max supports your workflow from the very start—on your motor-assisted walk to each patient. The compact, lightweight, and ergonomic design makes moving easy, and the unobstructed view to the front allows you to avoid obstacles on your route. Intelligent solutions such as on-the-go detector charging and long-life batteries help to transform care delivery by minimizing waiting times.

Handle easily at the patient's bedside

With its remarkably compact dimensions, MOBILETT Elara Max fits into very small patient rooms. When you are interacting with patients, the flexible MAXreach arm helps you quickly find the right position for the optimal image in every case: An impressive 180 degrees of lateral arm rotation maximizes your positioning options, while the Micro-Inch Mover buttons make fine adjustments easy. The system can also switch between batteries and main power so that it keeps going as long as you need it to.





The child-friendly giraffe design¹ creates a relaxing environment.

Tailor your imaging to very young patients

MOBILETT Elara Max allows you to image premature babies quickly and easily while limiting the stress placed on these sensitive patients. Special pediatric organ programs take account of higher tissue radiosensitivity and variations in size and weight. Copper filters reduce dose, while DiamondView MAX optimizes image quality to aid in diagnosis. A maximum focus height of 213.5 cm helps you to comply with the minimum required distance between X-ray tube and incubator, and the MAX mini detector fits neatly into incubator trays.

Achieve consistently high image quality

With its powerful 35 kW (450 mA) generator and wireless MAX detectors, MOBILETT Elara Max delivers excellent images fast. Gridless imaging in free examinations using DiamondView MAX postprocessing gives you the clear, detailed, high-contrast images you need for confident diagnostic decisions—at very low dose. At the same time, it reduces error rates and enhances your workflow efficiency.



DiamondView MAX enables high image quality thanks to optimized noise suppression and enhanced contrast.



MAX wi-D detector: 35 cm x 43 cm, 3.3 kg

MAX mini detector: 24 cm x 30 cm, 1.6 kg

¹Option

Accelerate your mobile workflow through seamless and secure integration

If your mobile X-ray system doesn't fit perfectly into your hospital landscape, you risk losing valuable time and facing critical situations. MOBILETT Elara Max integrates seamlessly and securely into your existing IT infrastructure and gives you access to relevant data wherever you need it.




Integrate your digital workflow seamlessly

The common user interface syngo® FLC enables a fully digital, one-stop imaging workflow that covers everything from patient registration to documentation. Wireless connectivity supports convenient and fast transfer of images to the hospital network. MOBILETT Elara Max also features integrated quality management: It provides a comprehensive overview of dose and image quality that allows you to easily monitor your system's performance and determine your staff's training needs.

High efficiency gains in accessing relevant data without having to leave the system¹





**Virtual
workstation
for advanced
use cases**

Access relevant data where and when you need it

Thanks to the virtual workstation¹, you can easily take hospital applications to the patient's bedside. This virtualized Windows 10® environment allows you to install institution-specific third-party software (e.g., an RIS client, a PACS viewer, an EMR client, or the hospital intranet) to efficiently manage operations directly on MOBILETT Elara Max. Since you don't have to leave the system in order to move to a separate workstation, you can spend more time with your patient.

Protect your mobile system and your patient data

MOBILETT Elara Max has high security standards and protective features that give you peace of mind. If you have to leave the system unattended at any time, you can lock it to prevent unauthorized access. In order to safeguard your system and patient health information against cyberattacks, MOBILETT Elara Max offers a comprehensive cybersecurity package including whitelisting, 90-days hotfixes, encrypted data storage as well as an advanced security package¹. This helps you comply with regulations for keeping patient data safe, confidential, and intact.



¹Option

The MAX effect— where gains multiply

Whether healthcare providers have to cope with tighter budgets or want to expand their business: they are forced to do more with less. What if there's a smart way to run a radiology department more efficiently?

Take a look at our high-performance MAX systems. They all share the same usability concept and imaging technology over the broadest range of clinical applications on the market. And if you combine two or more MAX systems, you'll experience the MAX effect: unique standardization across your fleet, so you can save both time and costs and ultimately boost your level of satisfaction.

Multiply the full potential of your department—with a winning combination of MAX systems.

Gain standardization

Achieve the same high-quality results throughout sites, rooms, and personnel—always. One MAX platform allows you to standardize workflows and image impressions the way you want. That means faster working and fewer mistakes. Resulting in better outcomes at lower costs.



Gain savings

Discover the smart way to save with the MAX effect: Save investment costs thanks to detector swapping, reduce the number of rooms required with More-in-1 MAX systems, and keep staff training to a minimum thanks to one usability.

Gain satisfaction

Increase satisfaction of all stakeholders by offering more patient comfort and lower dose, less stress and more confidence for staff, fast results, and high diagnostic accuracy at lower costs. For consistent high-quality care and enhanced reputation.



Multitom Rax
Twin Robotic X-ray

More-in-1
system



MOBILETT Elara Max
High-end digital mobile X-ray imaging



Multix Fusion Max
High-efficiency radiography



Ysio Max
High-performance radiography



Luminos dRF Max
2-in-1 remote fluoroscopy

More-in-1
system



Luminos Agile Max
2-in-1 patient-side fluoroscopy

More-in-1
system

See anatomical details with high image quality

To make fast and confident diagnostic decisions, you need high and reproducible image quality. MOBILETT Elara Max offers high-end imaging, automated workflows, and advanced postprocessing technologies.



Hand, d.v. view¹: Bone structures sharp and clearly visible



Pediatric chest, a.p. view²: Small structures outlined at low dose and reduced noise




Postoperative hip joint with implant, a.p. view¹: Diagnostically valuable results for dense organs

Enhanced contrast and a high level of detail visibility with DiamondView MAX



Without DiamondView MAX With DiamondView MAX



Without DiamondView MAX With DiamondView MAX

¹Courtesy of Centre Hospitalier de Mouscron, Belgium

²The clinical image was acquired with Mobilett Mira Max. Similarities between relevant components and parameters mean that the image impressions here are comparable to those possible with MOBILETT Elara Max.

Service and exchange

Increasing value by partnering throughout the entire equipment life cycle

Equipment Maintenance & Monitoring

Reliably servicing your MOBILETT Elara Max allows you to identify deviations from current norms to maximize equipment availability.

Education Management

Personalized education and training improves your staff's expertise as well as your equipment efficiency.

Fleet Management

A transparent overview of your fleet allows you to manage the performance and maintenance of your Siemens Healthineers equipment, 24/7.

Accessory Solutions

Products from our partner companies complement your use of our equipment in your daily workflow.

Performance Management

An intelligible overview of your radiography performance data helps you make prompt and well-informed decisions.

Asset Management & Planning

Access to innovative medical technology and equipment throughout the lifetime of the entire contract allows you to maximize focus on patient care.

Business Modeling & Financing

Customized business and financial models address your budgetary and enterprise needs, enabling you to remain more competitive.

Departmental Layout Optimization

3D visualization and digital twin analysis create more efficient workflows and a more enjoyable working environment.

Continuously adding value and caring for your equipment, your staff, your fleet, your workflows, your department, and your entire institution.

Technical specifications

Discover how you can stand out from the crowd in mobile X-ray imaging with MOBILETT Elara Max





MAX mini
is ideal for incubators
and extremities

MAX wi-D
is light, thin, fast

MAXswap
for advanced
detector sharing



DiamondView MAX
for enhanced contrast
and detail

syngo FLC
for an easy-to-learn,
common user interface



**Virtual workstation and
wireless connectivity**
for unrestricted data access

127.8 cm (50.3 in)

MOBILETT Elara Max is not commercially available in all countries. Due to regulatory reasons its future availability cannot be guaranteed. Please contact your local Siemens Healthineers organization for further details.

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all enabled by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 170 years of experience and 18,000 patents globally. With more than 48,000 dedicated colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens Healthineers sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/All of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens Healthineers reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. For the most current information, please contact your local sales representative from Siemens Healthineers.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we may recycle certain components where legally permissible. For recycled components we use the same extensive quality assurance measures as for factory-new components.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen, Germany
Phone: +49 9131 84-0
siemens-healthineers.com

USA

Siemens Medical Solutions USA, Inc.
Healthcare
40 Liberty Boulevard
Malvern, PA 19355-9998, USA
Phone: +1-888-826-9702
siemens-healthineers.us