


Study protocol

# Fixation of distal radius fracture

Ortho/trauma surgery



## Case description

### Patient

Male, aged between 46 and 60 years  
BMI range: 40+

### Diagnosis

Fracture of distal radius following  
bicycle accident  
AO classification: 2R3C3  
(Figs. 1–5)

### Surgical procedure

Fixation of distal radius fracture:  
temporary stabilization was  
achieved using K-wires followed by  
definitive fixation. A plate was then  
secured with screws to provide  
stable fixation.  
(Figs. 6–13)

### Benefits of CIARTIC Move

No C-arm technologist was needed,  
the surgeon operated the system  
from within the sterile field using  
**Smart Control**.

With **Position Assist** it is possible to  
store up to 12 procedure-specific  
positions, making surgical workflows  
more efficient.

The dose area product was  
63.02  $\mu\text{Gy}\cdot\text{m}^2$  and the radiation  
time was 141.4 seconds.

**3D imaging** played a vital part  
in confirming the success of the  
surgery.

---

### Courtesy of

BG Trauma Center Ludwigshafen,  
Germany


### System and software

CIARTIC Move, VB10A

---

## Clinical images and illustrations

Fig. 1




Fracture assessment and decision-making

Fig. 2



**Preoperative X-ray imaging:** anterior-posterior view of right distal radius showing fracture with intra-articular components

Fig. 3



**Preoperative X-ray imaging:** lateral view of right distal radius showing ventrally tilted distal radius fracture (flexion fracture)

## Clinical images and illustrations

Fig. 4




**Preoperative CT imaging:** right distal radius showing multifragmentary intra-articular fracture

Fig. 5



**Preoperative CT imaging:** right distal radius showing intra-articular fracture with interposed fragments in joint surface. Sagittal section

Fig. 6



Lateral transfixation with wire and reduction forceps. This was not performed in our case.

Fig. 7




Plate is completely fixed in position with angle-stable locked screws distally for fixation of multifragmentary fracture.

## Clinical images

Fig. 8




**Intraoperative image, anterior-posterior:** for correct positioning, the plate was placed on the radius and initially fixed in the elongated hole with a cortical screw. This image shows the correct position of the plate.

Fig. 9



**Intraoperative image, lateral:** for correct positioning, the plate was placed on the radius and initially fixed in the elongated hole with a cortical screw. This image shows the correct position of the plate and cortical screw.

Fig. 10



**Intraoperative image, anterior-posterior:** temporary fixation of the intermediate fragment using two K-wires. The plate is also already equipped with three additional screws in the shaft area and four screws distally.

Fig. 11




**Intraoperative final image, anterior-posterior:** screws and plate are correctly aligned. There is no intra-articular screw placement. The middle distal screw appears intra-articular due to its length but the lateral view shows that it is entirely extra-articular.

Fig. 12



**Intraoperative final image, lateral:** plate and screws are correctly positioned. Reconstruction of the joint surface is shown.

Fig. 13



**Intraoperative 3D scan:** implant is correctly positioned with no intra-articular screw placement. Joint surface is fully reconstructed.

The information presented in the study protocol is for illustration only and is not intended to be relied upon by the reader for instruction as to the practice of medicine. Any healthcare practitioner reading this information is reminded that they must use their own learning, training, and expertise in dealing with their individual patients. This material does not substitute for that duty and is not intended by Siemens Healthineers to be used for any purpose in that regard.

The statements by customers of Siemens Healthineers presented here are based on results that were achieved in the customer's unique setting. Since there is no "typical" hospital and many variables exist (e.g., hospital size, case mix, level of IT adoption), there can be no guarantee that other customers will achieve the same results.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this case are available throughout the Siemens Healthineers sales organization worldwide.

The product names and/or brands referred to are the property of their respective trademark holders.

All rights reserved.

---

**Siemens Healthineers Headquarters**

Siemens Healthineers AG  
Siemensstr. 3  
91301 Forchheim, Germany  
Phone: +49 9191 18-0  
[siemens-healthineers.com](https://www.siemens-healthineers.com)