

Atellica DCA Analyzer Featuring Auto-Checks Technology

siemens-healthineers.com/poc



Sample Volume and Integrity Check

Mitigates risk of sample collection errors



Optical Clarity Check

Verifies working order of components



Measurement Precision Check

Conducts multiple readings of each measurement to calculate and verify %CV is in range



Thermal Verification

Temperature of entire testing process monitored every second



Cartridge Temperature Check

Verifies cartridge is within appropriate operating range



Environmental Calibration

Adjusts for room temperature and altitude



Air Flow Check

Verifies proper air circulation and filter operation



Internal Bar-code Reader

Automatically reads cartridge lot information, expiration date, and calibration values

Offering lab-quality* results in only 4 minutes, the Atellica® DCA Analyzer is designed to make diabetes testing fast and efficient for patients, providers, and organizations. Employing our trusted, proprietary Auto-Checks™ technology, the Atellica DCA Analyzer runs internal checks and system calibrations before and during every single test. This is diabetes testing, evolved.

*Results comparable or equivalent to those produced by core laboratory equipment.

Published by Siemens Healthcare Diagnostics Inc. · Order No. 40-21-DX-886-76-1 · 06-2021 · ©Siemens Healthcare Diagnostics Inc., 2021

Atellica and Auto-Checks are trademarks of Siemens Healthcare Diagnostics Inc.

Not available for sale in the U.S.A. The products/features (mentioned herein) are not commercially available in all countries. Their future availability cannot be guaranteed.

SIEMENS
Healthineers