


Atellica Connect POC Insights

Device Interface List

siemens-healthineers.com/atellicaconnect


Total number of interfaces: 254

Last updated: June 30, 2024

Manufacturer	Device	Application	Connection	Download operator	Positive patient ID	Configuration
4S Information Systems	DAWN	Supplier software	→			
Abaxis	Piccolo Xpress	Clinical chemistry	→			
Abbott Diabetes Care	FreeStyle Precision Pro (FSPP)	Glucose, Ketone	↔	•	•	•
	FreeStyle Precision Pro (FSPP) Wireless mode	Glucose, Ketone	↔	•	•	•
	PCx Precision	Glucose	→			
	Precision Xceed Pro (PXP)	Glucose, Ketone	↔	•	•	•
Abbott Point of Care	Afinion AS100 UPDATE	HbA1c	↔	•		•
	Afinion 2 UPDATE	HbA1c	↔	•		•
	Alere i (ASTM)	Multitest	→			
	Alere i (POCT1-A)	Multitest	↔	•	•	•
	DIGIVAL	Multitest	→			
	CELL-DYN Emerald	Hematology	→			
	Cholestech LDX	Multitest	→			
	ID NOW	Multitest	↔	•	•	•
	Pima	HIV	→			
	i-STAT 1 (via i-STAT-DE-DMI software) UPDATE	Multitest	↔	•	•	
	i-STAT Alinity	Multitest	↔	•	•	
Acon Laboratories	Mission U120	Urinalysis	→			
Aidian (former Orion Diagnostica)	QuikRead 101	CRP	→			
	QuikRead go (ASTM)	CRP	→			
	QuikRead go (POCT1-A2)	CRP	↔	•	•	•
A. MENARINI Diagnostics	Zenit RA	Immunoassay	→			
Analyticon	CombiScan 100	Urinalysis	→			
	Urilyzer 100 UPDATE	Urinalysis	→			
	Urilyzer 100 Pro UPDATE	Urinalysis	↔		•	
	Urilyzer 100 Pro V2 / DocUReader 2 Pro	Urinalysis	↔	•	•	•
ApexBio	KETOSURE POC	Ketone	↔	•	•	•
Arkray	Adams HA-8160	HbA1c	→			

→ Unidirectional connectivity

↔ Bidirectional connectivity

Manufacturer	Device	Application	Connection	Download operator	Positive patient ID	Configuration
Beckman Coulter	Coulter Ac-T diff	Hematology	→			
	Coulter Ac-T diff 2	Hematology	→			
	Coulter Ac-T 8	Hematology	→			
	Coulter Ac-T 10	Hematology	→			
	DxH 500	Hematology	→			
	DxH 520	Hematology	→			
Becton Dickinson	BD Veritor Plus (via BD Cloud Connectivity Solution)	Multitest	→			
bioMérieux	Biofire Spotfire	Multitest	↔	•		•
Bio-Rad	Bio-Rad Unity Real Time (URT)	Supplier software	→			
Biosynex	BSX Reader	Multitest	→			
biosurfit	Spinit	Multitest	→			
Boditech	AFIAS-1 UPDATE	Immunoassay	→			
	AFIAS-3 UPDATE	Immunoassay	→			
	AFIAS-6 UPDATE	Immunoassay	→			
	AFIAS-10 UPDATE	Immunoassay	→			
Bosch Healthcare	Vivalytic one	Multitest	→			
Boule Medical	Medonic M-series M16/M20	Hematology	→			
	Medonic M-series M16C/M20C	Hematology	→			
	Medonic M-series M16M/M20M	Hematology	→			
	Medonic M-series M16S/M20S	Hematology	→			
	Medonic M-series M32B	Hematology	→			
	Medonic M-series M32C	Hematology	→			
	Medonic M-series M32M	Hematology	→			
	Medonic M-series M32S	Hematology	→			
	Swelab Alfa Basic / Plus Basic	Hematology	→			
	Swelab Alfa Standard / Plus Standard	Hematology	→			
	Swelab Alfa Cap / Plus Cap	Hematology	→			
	Swelab Alfa Sampler / Plus Sampler	Hematology	→			
Caretaker Medical	Caretaker4	Vital Signs Monitor	→			
Cepheid	GeneXpert PC	Multitest	→			
	GeneXpert with Touchscreen	Multitest	↔	•		•
	GeneXpert Infinity	Multitest	→			
	GeneXpert Omni	Multitest	→			
	GeneXpert Xpress	Multitest	↔	•		•
Chirus	DXpress Reader	Multitest	→			
CoaguSense	Coag-Sense PT2	Coagulation	↔	•	•	
concile	Omega 100 (Q100)	Immunoassay	→			
Diagnostics For The Real World (DRW)	Samba II	Multitest	↔		•	
Diatron	Diatron Aquila	Hematology	→			
DnaNudge	NudgeBox	Immunoassay	↔	•		
Dräger	Jaundice Meter JM-105	Bilirubin	→			
Drucker Diagnostics	QBC Star	Hematology	→			

Manufacturer	Device	Application	Connection	Download operator	Positive patient ID	Configuration
EKF Diagnostics	Quo-Lab (via CIB)	HbA1c	↔	•	•	•
	Quo-Test (via CIB)	HbA1c	↔	•	•	•
Enigma Diagnostics	Enigma MiniLab	Multitest	→			
ET Healthcare	Pylon 3D	Multitest	→			
	Pylon IRIS	Multitest	→			
Eurolyser Diagnostica	Eurolyser CUBE	Multitest	→			
	Eurolyser CUBE-S	Multitest	→			
	Eurolyser smart	Multitest	→			
	Eurolyser solo	Multitest	→			
Fuji	Fuji Dri-Chem 4000	Clinical chemistry	→			
	Fuji Dri-Chem NX10N	Clinical chemistry	→			
	Fuji Dri-Chem NX500	Clinical chemistry	→			
	Fuji Dri-Chem NX600	Clinical chemistry	→			
	Fuji Dri-Chem NX700	Clinical chemistry	→			
Genedrive	Genedrive System	Multitest	→			
GenePOC Diagnostics	GenePOC Revogene	Multitest	↔	•		•
Haemonetics	ClotPro	Hemostasis	↔		•	
	TEG 6s (via TEGManager HaemoCommunicator)	Hemostasis	↔		•	
Helena Laboratories	Cascade Abrazo	Coagulation	↔			•
	Cascade POC	Coagulation	↔	•		
HemoCue	Glucose 201 DM UPDATE	Glucose	↔	•	•	•
	Glucose 201 DM RT UPDATE	Glucose	↔	•	•	•
	Hb 201 DM UPDATE	Hemoglobin	↔	•	•	•
	HbA1c 501	Hemoglobin	→			
	WBC Diff UPDATE	Hematology	→			
HemoSonics	Quantra	Hemostasis	↔	•	•	
HLS Therapeutics	CSAN Pronto	Hematology	→			
Hologic	PeriLynx	Immunoassay	→			
	Rapid fFN	Immunoassay	→			
HORIBA	ABX Micros 60	Hematology	→			
	ABX Micros ES 60 UPDATE	Hematology	→			
	ABX Microsemi CRP LC-667G	Hematology	→			
	ABX Microsemi CRP LC-767G	Hematology	→			
	Yumizen H500	Hematology	→			
	Yumizen H550	Hematology	→			
iLine Microsystems	micro INR Expert UPDATE	Coagulation	↔	•	•	
I-SENS	CareSens NEW	Glucose	↔	•	•	•
Lite-On Technology Corp.	skyla HB1 (POCT1-A2)	Clinical chemistry	↔	•		
	skyla HB1 (ASTM)	Clinical chemistry	→			
	skyla Hi	HbA1c	↔	•		
LumiraDx	LumiraDx (Instrument and Host)	Multitest	↔	•	•	•
Macherey-Nagel	Uryxxon Relax (Connect)	Urinalysis	↔	•	•	
Magnasense Technologies	Magnia Reader	Multitest	→			

Manufacturer	Device	Application	Connection	Download operator	Positive patient ID	Configuration
Medcaptain	Immu F6	Multitest	→			
Medtronic	ACT Plus	Coagulation	→			
	HMS Plus	Coagulation	→			
MeMed	MeMed Key	Multitest	→			
Micropoint Biotechnologies	qLabs Q-3 Plus NEW	Coagulation	→			
	qLabs Q-3 Pro NEW	Coagulation	→			
Mitsubishi	PATHFATH	Immunoassay	→			
nal von minden (BioSciTec and Mediwiss)	Rapid Slide Scanner (RSS) III	Drug-Screen	→			
Nihon Kohden	Celltac Alpha (MEK 6500/6510)	Hematology	→			
	Celltac Chemi (CHM-4100)	Hematology	→			
Norma Diagnostika	Norma Icon 3	Hematology	→			
	Norma Icon 5	Hematology	→			
Nova Biomedical	NovaNet*	Supplier software	↔	•		
	Stat Profile pHox Ultra (via NovaNet)	Blood gas	→			
	Stat Profile Prime (via NovaNet)	Blood gas	→			
	Stat Profile Prime Plus (via NovaNet)	Blood gas	↔	•		
	StatSensor Creatinine (via NovaNet)	Creatinine	↔	•		
	StatStrip Glucose (via NovaNet)	Glucose	↔	•		
	StatStrip Glucose/Ketone (via NovaNet)	Glucose, Ketone	↔	•		
	StatStrip Hb/Hct (via NovaNet)	Hemoglobin	↔	•		
	StatStrip HHL (via NovaNet)	Hemoglobin/ Hematocrit/Lactate	↔	•		
	StatStrip Lactate (via NovaNet)	Lactate	↔	•		
	Stat Profile CCX	Blood gas	↔		•	
	Stat Profile pHox / pHox plus	Blood gas	↔		•	
	Stat Profile pHox Ultra	Blood gas	→			
	Stat Profile Prime	Blood gas	→			
	StatSensor Creatinine	Creatinine	↔	•	•	•
	StatStrip Glucose	Glucose	↔	•	•	•
	StatStrip Glucose/Ketone	Glucose, Ketone	↔	•	•	•
	StatStrip Hb/Hct	Hemoglobin	↔	•	•	•
	StatStrip Lactate	Lactate	↔	•	•	•
OPTI Medical	OPTI CCA-TS2	Blood gas	↔	•	•	
Philips	BiliChek	Bilirubin	→			
PixCell Medical	HemoScreen	Hematology	↔	•	•	•
QuidelOrtho	Savanna UPDATE	Multitest	↔	•	•	
	Sofia UPDATE	Immunoassay	→			
	Sofia 2 UPDATE	Immunoassay	→			
	Triage Meter Plus	Cardiac markers	↔	•		•
	Triage Meter Pro	Cardiac markers	↔	•		•

*NovaNet is needed for the StatStrip-Family connections with a firmware version greater than 4.17.8.47.
Other devices can be connected directly or via NovaNet

Manufacturer	Device	Application	Connection	Download operator	Positive patient ID	Configuration
Radiometer	ABL7xx series (ASTM 2.1)	Blood gas	↔		•	
	ABL8xx series / ABL8xx FLEX (ASTM 2.1)	Blood gas	↔		•	
	ABL80 FLEX (ASTM 2.1)	Blood gas	↔		•	
	ABL835 (POCT1-A)	Blood gas	→			
	ABL90 / ABL90 FLEX (ASTM 2.1) UPDATE	Blood gas	↔		•	
	ABL90 / ABL90 FLEX (POCT1-A)	Blood gas	→			
	AQT90 FLEX (ASTM 2.1)	Multitest	↔		•	
	AQURE (ASTM 2.1)	Supplier software	↔	•	•	
	NPT7 (ASTM 2.1)	Blood gas	→			
	RADIANCE (ASTM 2.1)	Supplier software	↔	•	•	
ReLIA Biological Engineering	FY-100	Multitest	→			
	TZ-206	Multitest	→			
	TZ-21	Multitest	→			
	TZ-212	Multitest	→			
	TZ-3000/ TZ-3000Plus	Multitest	→			
	TZ-310	Multitest	→			
	TZ-320	Multitest	→			
	TZ-33	Multitest	→			
Response Biomedical	RAMP 200	Multitest	↔	•		•
Roche Diagnostics	ACCU-CHEK Inform II UPDATE	Glucose	↔	•	•	•
	CoaguChek Pro II UPDATE	Coagulation	↔	•	•	•
	CoaguChek XS Plus UPDATE	Coagulation	↔	•	•	•
	CoaguChek XS Pro UPDATE	Coagulation	↔	•	•	•
	cobas b 101	Hemoglobin	↔	•	•	•
	cobas b 123	Blood gas	↔	•	•	
	cobas c 111	HbA1c	→			
	cobas e 411	Immunoassay	→			
	cobas h 232 UPDATE	Cardiac markers	↔	•	•	•
	cobas 6000 series	Clinical chemistry	→			
	cobas IT 1000 / cobas infinity	Supplier software	↔	•		
	cobas Liat	Multitest	↔	•	•	•
	cobas pulse	Glucose	↔	•	•	•
	cobas u 411	Urinalysis	→			
	Urisys 1100 UPDATE	Urinalysis	↔	•		

Manufacturer	Device	Application	Connection	Download operator	Positive patient ID	Configuration
Samsung	LABGEO HC10	Hematology	↔	•		
	LABGEO IB10	Immunoassay	↔	•		
	LABGEO PT10	Clinical chemistry	↔	•		
SD Biosensor	Standard F200	Multitest	→			
	Standard F2400	Multitest	→			
	Standard M10	Multitest	→			
Siemens Healthcare Diagnostics	ADVIA 360 Hematology System	Hematology	→			
	ADVIA Centaur CP Immunoassay System	Immunoassay	→			
	ADVIA Centaur XP Immunoassay System	Immunoassay	→			
	Atellica DCA Analyzer (Engine and IDM) UPDATE	Hemoglobin	↔	•	•	•
	Atellica 1500 Automated Urinalysis System	Urinalysis	↔			
	Atellica UAS 800 Automated Urinalysis System	Urinalysis	↔			
	Atellica VTLi Immunoassay System UPDATE	Immunoassay	↔	•		
	CLINITEK Advantus Urine Chemistry Analyzer	Urinalysis	→			
	CLINITEK Status Analyzer	Urinalysis	→			
	CLINITEK Status+ Urine Chemistry Analyzer UPDATE	Urinalysis	↔	•		•
	DCA Vantage Analyzer	HbA1c	↔	•		
	epoc Blood Analysis System (Reader and Host) UPDATE	Multitest	↔	•	•	•
	epoc NXS Blood Analysis System (Reader and Host) UPDATE	Multitest	↔	•	•	•
	RAPIDComm Data Management System UPDATE	Supplier software	↔	•	•	
	RAPIDLab 1200 series Blood Gas Systems	Blood gas	↔	•	•	
	RAPIDLab 348EX Blood Gas System	Blood gas	→			
	RAPIDPoint 500 Blood Gas System UPDATE	Blood gas	↔	•	•	
	RAPIDPoint 500e Blood Gas System UPDATE	Blood gas	↔	•	•	•
	Stratus CS Acute Care Diagnostic System	Cardiac markers	→			
	Stratus CS 200 Acute Care Diagnostic System	Cardiac markers	↔	•		

Manufacturer	Device	Application	Connection	Download operator	Positive patient ID	Configuration
Sight Diagnostics	OLO	Hematology	→			
	OLO (POCT1-A)	Hematology	↔	•		•
Sphere Medical	Proxima	Blood gas	→			
Spindiag	Rhonda Player	Multitest	↔	•		
Stat-Diagnostica	DiagCORE	Multitest	→			
Sysmex	Alifax Roller	Hematology	→			
	CA-560	Coagulation	→			
	CA-600	Coagulation	→			
	pocH-100i	Hematology	→			
	UF-500i (via SUIT Manager)	Urinalysis	→			
	UF-1000i (via SUIT Manager)	Urinalysis	→			
	UX-2000 (via SUIT Manager)	Urinalysis	→			
	XE Series (via SUIT Manager)	Hematology	→			
	XN Series (via SUIT Manager)	Hematology	→			
	XN-L Series (via SUIT Manager)	Hematology	→			
	XS Series (via SUIT Manager)	Hematology	→			
	XT Series (via SUIT Manager)	Hematology	→			
	XP-300	Hematology	→			
	XQ-320	Hematology	→			
	XQ-520	Hematology	→			
TestCard	ClearScreen	Multitest	→			
Thermo Fisher Scientific	B-R-A-H-M-S Direct Reader	Immunoassay	→			
Trinity Biotech	Meritas POC Analyzer	Cardiac markers	→			
Universal Biosensors	Xprecia Stride	Coagulation	↔	•		•
Veda.Lab	Easy Reader+	Multitest	↔	•		
Voden	VChemistry S	Clinical chemistry	→			
Welch Allyn	Connex Spot Monitor	Vital Signs Monitor	→			
Werfen	Avoximeter 1000E	Blood gas	→			
	Avoximeter 4000	Blood gas	→			
	BIO-Flash	Immunoassay	→			
	GEM Hemochron 100	Coagulation	↔	•	•	•
	GEM Premier 3000	Blood gas	↔		•	
	GEM Premier 3500	Blood gas	↔		•	
	GEM Premier 4000 UPDATE	Blood gas	↔	•	•	
	GEM Premier 5000 UPDATE	Blood gas	↔	•	•	
	GEM Premier 7000 NEW	Blood gas	↔	•	•	
	GEMweb Plus	Supplier software	↔		•	
	HEMOCHRON Response	Coagulation	→			
	HEMOCHRON Signature Elite UPDATE	Coagulation	↔	•		•
	ROTEM delta (direct or via ROTEM Live or GEMweb Live)	Hemostasis	↔		•	
	ROTEM gamma (direct or via ROTEM Live or GEMweb Live)	Hemostasis	↔		•	
	ROTEM sigma (direct or via ROTEM Live or GEMweb Live)	Hemostasis	↔		•	

At Siemens Healthineers, we pioneer breakthroughs in healthcare. For everyone. Everywhere. Sustainably. As a leader in medical technology, we want to advance a world in which breakthroughs in healthcare create new possibilities with a minimal impact on our planet. By consistently bringing innovations to the market, we enable healthcare professionals to innovate personalized care, achieve operational excellence, and transform the system of care.

Our portfolio, spanning in vitro and in vivo diagnostics to image-guided therapy and cancer care, is crucial for clinical decision-making and treatment pathways. With the unique combination of our strengths in patient twinning,* precision therapy, as well as digital, data, and artificial intelligence (AI), we are well positioned to take on the greatest challenges in healthcare. We will continue to build on these strengths to help overcome the world's most threatening diseases, enable efficient operations, and expand access to care.

We are a team of more than 72,000 Healthineers in over 70 countries passionately pushing the boundaries of what is possible in healthcare to help improve the lives of people around the world.

**Personalization of diagnosis, therapy selection and monitoring, aftercare, and managing health.*

All trademarks are the property of their respective owners. Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim, Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Published by

Siemens Healthcare Diagnostics Inc.
Point of Care Diagnostics
333 Coney Street
East Walpole, MA 02032-1516
USA
Phone: +1 508-668-5000