

# IHE Integration Statement

## Biograph Trinion VK10A

Molecular Imaging

IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens Healthcare GmbH	Biograph Trinion 18cm CT64 Biograph Trinion 24cm CT64 Biograph Trinion 18cm CT128 Biograph Trinion 24cm CT128	VK10A	01/2024
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Acquisition Modality	Patient Based Worklist Query
		Broad Worklist Query
		MPPS - Modality group case → Siemens solution: Multi Study
	Image Display	No options defined
	Evidence Creator	None
Radiation Exposure Monitoring (REM)	Acquisition Modality	None
Patient Information Reconciliation (PIR)	Acquisition Modality	No options defined
Access to Radiology Information (ARI)	Image Display	None
Consistent Presentation of Images (CPI)	Print Composer	No options defined
Consistent Time (CT)	Time Client	No options defined
Portable Data for Imaging (PDI)	Portable Media Creator	No options defined
	Image Display	No options defined
	Print Composer	No options defined
NM Image (NMI)	Acquisition Modality	No options defined
	Image Display	Review Option
Internet address for vendor's IHE information: <a href="https://siemens-healthineers.com/ihe">siemens-healthineers.com/ihe</a>		
Links to Standards Conformance Statements for the Implementation		
HL7	<a href="https://siemens-healthineers.com/hl7">siemens-healthineers.com/hl7</a>	
DICOM	<a href="https://siemens-healthineers.com/dicom">siemens-healthineers.com/dicom</a>	
Links to general information on IHE worldwide: <a href="http://ihe.net/IHE_Worldwide/">http://ihe.net/IHE_Worldwide/</a>		

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries. The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases. Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Please contact your local Siemens sales representative for the most current information. In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components. Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced. Caution: Federal law restricts this device to sale by or on the order of a physician.

---

**Siemens Healthineers Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
[siemens-healthineers.com](http://siemens-healthineers.com)

**Legal Manufacturer**

Siemens Medical Solutions USA, Inc.  
Molecular Imaging  
2501 N. Barrington Rd.  
Hoffman Estates, IL 60192  
USA  
Telephone: +1-888-826-9702  
[www.usa.siemens.com/healthcare](http://www.usa.siemens.com/healthcare)

---

**EU**

EU Authorized Representative is  
Siemens Healthcare GmbH  
Henkestrasse 127  
91052 Erlangen  
Germany