

Siemens AG Healthcare			Date	
	SOMATOM Sensation 16 / Cardiac			
	SOMATOM Sensation 16 / Cardiac		VB42 12/2013	
	SOMATOM Sensation 40 SOMATOM Sensation 64 / Cardia	VB42		
пеаннсате				
	SOMATOM Sensation Open			
This product implements gration Profiles, Actors a	all transactions required in the IHE Techr nd Options listed below:	nical Framework to support t	he IHE Inte-	
Integration Profiles Implemented	Actors Implemented	Options Ir	Options Implemented	
Scheduled Workflow	Acquisition Modality	Acquisition Modality Patient based work		
		Broad worklis	Broad worklist query	
		Modality grou Siemens soluti	ıp case → on: Study split	
	Image Display No options defined in		efined in IHE	
	Evidence Creator none			
Patient Information Reconciliation	Acquisition Modality	No options do	No options defined in IHE	
Access to Radiology nformation	Image Display	none		
Consistent Presentation of Images	Print Composer	No options de	No options defined in IHE	
Radiation Exposure Measurement	Acquisition Modality	none	none	
	Dose Info Consumer (available as a standalone tool by CARE Analytics)		Except transaction "Store dose information"	
nternet address for ve	ndor's IHE information: http://www.sie	emens.com/IHE		
_inks to Standards C	onformance Statements for the Im	plementation		
HL7	n.a.			
DICOM	http://www.siemens.com/dicom			
_inks to general info	_⊥ mation on IHE			

http://www.ihe-europe.net

http://www.ihe-j.org

http://www.ihe.net

## **Global Siemens Headquarters**

Siemens AG Wittelsbacherplatz 2 80333 München Germany

## Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestraße 127
91052 Erlangen
Germany
Phone +49 9131 84-0
www.siemens.com/healthcare

## Legal Manufacturer

Siemens AG Wittelsbacherplatz 2 DE-80333 München Germanv

## **Global Business Unit**

Siemens AG
Medical Solutions
Computed Tomography & Radiation Oncology
Siemensstraße 1
DE-91301 Forchheim
Germany
Phone +49 9191 18-0
www.siemens.com/healthcare

Order No. CT-IHE.640.25.01.02 © 2013, Siemens AG

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries. The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases. Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Please contact your local Siemens sales representative for the most current information. In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components. Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced. Caution: Federal law restricts this device to sale by or on the order of a physician.