

IHE Integration Statement

Product Name: Biograph mCT

Version: VG30A

Date: September 2011

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	Biograph mCT	VG30A	09/13/2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Patient based worklist query	
		Broad worklist query	
	Image Display	No options defined in IHE	
	Evidence Creator	none	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net		In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org