

IHE Integration Statement

Product Name:

syngo Workflow SLR

Version: VA31A

Date: June 2011

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	syngo Workflow SLR	VA31A	7-JUN-2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Information Reconciliation	Order Filler/Department System Scheduler	None	
	Performed Procedure Step Manager	None	
Scheduled Workflow	Order Filler/Department System Scheduler	None	
	Performed Procedure Step Manager	None	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In United States: http://www.iheusa.org	In Japan: http://www.jira-net.or.jp/ihe-j	
In Australia: http://www.ihe.net.au	In Canada: http://www.ihe-canada.com	In Europe: http://www.ihe-europe.net	