

| IHE Integration Statement | | | |
|--|--|--|---|
| Vendor | Product Name | Version | Date |
| Siemens Medical Solutions | syngo® Imaging | V30B | Aug. 2007 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented | |
| Scheduled Workflow (SWF) | Image Archive Image Manager PPS Manager | none | none |
| Patient Information Reconciliation (PIR) | Image Archive Image Manager PPS Manager | none | none |
| Reporting Workflow (RWF) | Image Archive Image Manager PPS Manager | none | none |
| Consistent Presentation of Images (CPI) | Image Archive Image Manager | none | none |
| Key Image Note (KIN) | Image Archive Image Manager | none | none |
| Evidence Documents (ED) | Image Archive Image Manager | none | none |
| Simple Image and Numeric Report (SINR) | Report Repository | none | |
| Access to Radiology Information (ARI) | Image Archive Image Manager Report Repository | none | Multiple Sources Multiple Sources |
| Basic Security ^{*1)} (SEC) | Audit Record Repository Time Server | none | none |
| Nuclear Medicine Imaging (NMI) | Image Archive Image Manager | none | none |
| Portable Data for Imaging (PDI) | Portable Media Creator ^{*2)} Image Display ^{*2)} Print Composer ^{*2)} | none | no options defined no options defined |
| Consistent Time (CT) | Time Client Time Server | none | none |
| Internet address for vendor's IHE information: http://www.siemens.com/ihe | | | |
| Links to Standards Conformance Statements for the Implementation | | | |
| HL7 DICOM | n.a. http://www.siemens.com/dicom | | |
| Links to general information on IHE | | | |
| In North America: http://www.rsna.org/IHE | | In Europe: http://www.ihe-europe.org | In Japan: http://www.jira-net.or.jp/ihe-j |

1) Actors only. No full support of SEC profile.

2) Supported with a syngo® Studio Advanced Workplace