

IHE Integration Statement

Product Name:

syngo Dynamics Workplace

Version: 8.0

Date: December 2010

IHE Integration Statement		Date	01 Dec 2010
Vendor	Product Name	Version	
Siemens Medical Solutions	syngo® Dynamics Workplace	8.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Displayable Reports (DRPT – CARD)	Report Creator	■	
Evidence Documents (ED - CARD)	Evidence Creator	■ Echo Evidence - o Vascular (5100) - o Echocardiography (5200)	
Retrieve ECG for Display (ECG - CARD)	Display	■ Retrieve ECG List (CARD-5) ■ Retrieve ECG Document for Display (CARD-6) o Scalable Vector Graphics	
Retrieve Information for Display (RID - ITI)	Display	■	
Internet address for vendor's IHE information: http://www.siemens.com/IHE			
Links to Standards Conformance Statements for the Implementation			
HL7	N/A		
DICOM	http://www.siemens.com/DICOM		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	