

IHE Integration Statement

Product Name: Luminos Agile
Version: VC10
Date: July 2011

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens AG Agile	Luminos Agile	VC10	July 2011

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
RadioLOGY Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query
		Broad Worklist Query
		PPS Exception Management
		Billing and Material Management
		Assisted Acquisition Protocol Setting
Patient Information Reconciliation	Acquisition Modality	No options defined
Consistent Presentation of Images	Print Composer	No options defined

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org
---	--	---