Clarify vascular enhancement technology:

ACUSON S2000™ ultrasound system: There are no pre-requisites for this software on the ACUSON S2000 system.

ACUSON Antares™ ultrasound system: Systems that shipped before August 1, 2008 must be at revision 3.5 or greater software level.

SONOLINE Antares™ ultrasound system: Revision 2.0 or greater is required for SONOLINE Antares systems.

ACUSON X300™ ultrasound system and ACUSON X300™ ultrasound system, premium edition (PE): Requires 2.0 or greater software level.

ACUSON Sequoia™ ultrasound systems: Requires 8.5s software release or greater.

Please Note: Additional technical pre-requisites may apply. Upon receiving your request, your local Siemens representative will clarify whether your system meets the requirements.