Advanced SieClear spatial compounding:

ACUSON S2000™ ultrasound systems: There are no pre-requisites for Advanced SieClear™ spatial compounding

ACUSON Antares™ ultrasound systems: Systems that shipped before August 1, 2008 must be at revision 5.0 with Flat Panel Monitor.

SONOLINE Antares™ ultrasound systems: Revision 5.0 is required for SONOLINE Antares systems.

ACUSON Sequoia™ ultrasound systems: Requires 12.x software, Progressive Media and Display and HSIX.

Please Note: Additional technical pre-requisites may apply. Upon receiving your request, your local Siemens representative will clarify whether your system meets the requirements.