

IHE Integration Statement

Product Name: Uroskop Omnia
Version: VC10/VD10
Date: August 2012

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens AG Healthcare	Uroskop Omnia	VC10/VD10	August 2012
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
		PPS Exception Management	
		Billing and Material Management	
		Assisted Acquisition Protocol Setting	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Print Composer	No options defined	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	