

IHE Integration Statement

Product Name: *syngo* Workflow SLR

Version: VB10A

Date: 23 Jul 2013

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	<i>syngo</i> Workflow SLR	VB10A	23-JUL-2013

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Patient Information Reconciliation	Order Filler/Department System Scheduler	None
	Performed Procedure Step Manager	None
Scheduled Workflow	Order Filler/Department System Scheduler	None
	Performed Procedure Step Manager	None
Cross-enterprise Document Sharing for Imaging	Document Source	None
	Document Consumer	None
Cross-enterprise Document Sharing	Document Source	None
	Document Consumer	None

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.
DICOM	http://www.siemens.com/dicom

Links to general information on IHE

In North America: http://www.ihe.net	In United States: http://www.iheusa.org	In Japan: http://www.jira-net.or.jp/ihe-j
In Australia: http://www.ihe.net.au	In Canada: http://www.ihe-canada.com	In Europe: http://www.ihe-europe.net