

IHE Integration Statement

syngo.plaza
VB10A
May 2013

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens AG Healthcare Sector	syngo.plaza	VB10A	May 2013

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Image Archive	none
	Image Manager	none
	Image Display	no options defined
	PPS Manager	no options defined
	Evidence Creator	none
Patient Information Reconciliation (PIR)	Image Archive	no options defined
	Image Manager	no options defined
	PPS Manager	no options defined
Portable Data for Imaging (PDI)	Portable Media Creator (Only with separately installed Media writing software)	none
	Portable Media Importer	no options defined
	Image Display	no options defined
	Print Composer	no options defined
Access to Radiology Information (ARI)	Image Archive	no options defined
	Image Manager	no options defined
Consistent Presentation of Images (CPI)	Image Archive	no options defined
	Image Manager	no options defined
	Evidence Creator	no options defined
	Image Display	no options defined
	Print Composer	none
Simple Image and Numeric Report (SINR)	Report Creator	none
Mammography Image Integration Profile	Image Archive	no options defined
	Image Manager	no options defined
	Image Display	none
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer	none
Consistent Time (CT)	Time Client	none
Nuclear Medicine Image (NMI)	Image Manager	none
	Image Archive	none

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	n.a.		
DICOM	http://www.siemens.com/dicom		

Links to general information on IHE

In North America: http://www.rsna.org/IHE	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j
--	--	---