

**SIEMENS**



IHE Integration Statement

# SOMATOM Emotion

Computed Tomography

IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens AG Healthcare	SOMATOM Emotion 6 SOMATOM Emotion 16	VC20	09/2013
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Patient based worklist query	
		Broad worklist query	
		Modality group case → Siemens solution: Study split	
	Image Display	No options defined in IHE	
	Evidence Creator	none	
Patient Information Reconciliation	Acquisition Modality	No options defined in IHE	
Access to Radiology Information	Image Display	none	
Consistent Presentation of Images	Print Composer	No options defined in IHE	
Radiation Exposure Measurement	Acquisition Modality	none	
	Dose Info Consumer (available as a standalone tool by CARE Analytics)	Except transaction “Store dose information”	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.siemens.com/IHE">http://www.siemens.com/IHE</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
HL7	n.a.		
DICOM	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	

**Global Siemens Headquarters**

Siemens AG  
Wittelsbacherplatz 2  
80333 München  
Germany

**Global Siemens Healthcare Headquarters**

Siemens AG  
Healthcare Sector  
Henkestraße 127  
91052 Erlangen  
Germany  
Phone +49 9131 84-0  
[www.siemens.com/healthcare](http://www.siemens.com/healthcare)

**Legal Manufacturer**

Siemens Shanghai Medical Equipment Ltd.  
278 Zhou Zhu Road  
201318 Shanghai  
P.R. China

**EU Authorized Representative**

Siemens AG  
Medical Solutions  
Henkestraße 127  
91052 Erlangen  
Germany

[www.siemens.com/healthcare](http://www.siemens.com/healthcare)

Order No. CT-IHE-C.640.22.01.02

© 2013, Siemens Shanghai Medical Equipment Ltd.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries. The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases. Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Please contact your local Siemens sales representative for the most current information. In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components. Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced. Caution: Federal law restricts this device to sale by or on the order of a physician.

**[www.siemens.com/healthcare](http://www.siemens.com/healthcare)**