

IHE Integration Statement

Product Name: *syngo.via*

Version: VA30A

Date: 2013-12-06

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	syngo.via	VA30A	2013-12-06

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Evidence Creator	None
	Image Display	None
Access to Radiology Information (ARI)	Image Display	Multiple Sources
Portable Data for Imaging (PDI)	Portable Media Creator	None
	Portable Media Importer	No options defined
	Display	No options defined
Evidence Documents(ED)	Evidence Creator	No options defined
Consistent Time (CT)	Time Client	None

Internet address for vendor's IHE information: <http://www.siemens.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	http://www.siemens.com/hl7
DICOM	http://www.siemens.com/dicom

Links to general information on IHE

In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org	In Japan: http://www.ihe-j.org
--	---	--