

IHE Integration Statement

Product Name: ACUSON S Family
Ultrasound System

Version:

ACUSON S1000™ – VC30 and VC31

ACUSON S2000™ – VC30 and VC31

ACUSON S3000™ – VC30 and VC31

Date: September 2014

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	ACUSON S Family Ultrasound System		Sept 2014
	ACUSON S1000	VC30 and VC31	
	ACUSON S2000	VC30 and VC31	
	ACUSON S3000	VC30 and VC31	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology			
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		Patient Based Worklist Query	
Patient Information Reconciliation	Acquisition Modality	None	
Evidence Documents	Acquisition Modality	None	
	Evidence Creator (Template: TID 5000 OB-GYN Ultrasound)	None	
Cardiology			
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient Based Worklist Query	
Evidence Documents	Acquisition Modality	Echo Evidence	
	Evidence Creator (Templates: TID 5100 Vascular Ultrasound TID 5200 Echocardiography)	Echo Evidence	
Internet address for vendor’s IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org	In Japan: http://www.ihe-j.org	