

IHE Integration Statement

Product Name: ACUSON X Family
Ultrasound System

Version:

ACUSON X300™ – 2.5

ACUSON X300PE™ – 7.5

ACUSON X600™ – VA10

ACUSON X700™ – 1.0/1.1

ACUSON X700™ – VB20

Date: September 2014

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	ACUSON X Family Ultrasound System		Sep 2014
	ACUSON X300	2.5	
	ACUSON X300PE	7.5	
	ACUSON X600	VA10	
	ACUSON X700	1.0/1.1	
	ACUSON X700	VB20	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Information Reconciliation	Acquisition Modality	None	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Evidence Documents	Acquisition Modality	OB-GYN Ultrasound Procedure Report (TID 5000)	
		Vascular Ultrasound Procedure Report (TID 5100)	
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Cardiology Evidence Documents	Acquisition Modality	Echo Evidence	
		Adult Echocardiography Procedure Report (TID 5200)	
Internet address for vendor’s IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org	In Japan: http://www.ihe-j.org	