

IHE Integration Statement

Product Name: *syngo* Workflow
SLR

Version: VB20A

Date: 17 Feb 2015

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Medical Solutions	syngo Workflow SLR	VB20A	17-Feb-2015
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Information Reconciliation	Order Filler/Department System Scheduler	None	
	Performed Procedure Step Manager	None	
Scheduled Workflow	Order Filler/Department System Scheduler	None	
	Performed Procedure Step Manager	None	
Cross-enterprise Document Sharing for Imaging	Document Source	None	
	Document Consumer	None	
Cross-enterprise Document Sharing	Document Source	None	
	Document Consumer	None	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	There is no HL7 conformance stantement. The HI7 requirements are described in syngo Workflow SLR Generic_HIS_Interface_Mapping Guide		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In United States: http://www.iheusa.org	In Japan: http://www.jira-net.or.jp/ihe-j	
In Australia: http://www.ihe.net.au	In Canada: http://www.ihe-canada.com	In Europe: http://www.ihe-europe.net	