

# IHE Integration Statement

Product Name: ACUSON OXANA  
Family Ultrasound System

Version:

ACUSON OXANA 1 <sup>TM</sup> – VD10

ACUSON OXANA 2 <sup>TM</sup> – VD10

ACUSON OXANA 3 <sup>TM</sup> – VD10

Date: May 2015

# IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	ACUSON OXANA Family Ultrasound System		May 2015
	ACUSON OXANA 1	VD10	
	ACUSON OXANA 2	VD10	
	ACUSON OXANA 3	VD10	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology			
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		Patient Based Worklist Query	
Patient Information Reconciliation	Acquisition Modality	None	
Evidence Documents	Acquisition Modality	None	
	Evidence Creator (Template: TID 5000 OB-GYN Ultrasound)	None	
Cardiology			
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient Based Worklist Query	
Evidence Documents	Acquisition Modality	Echo Evidence	
	Evidence Creator (Templates: TID 5100 Vascular Ultrasound TID 5200 Echocardiography)	Echo Evidence	
Internet address for vendor’s IHE information: <a href="http://www.siemens.com/ihe">http://www.siemens.com/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	<a href="http://www.siemens.com/dicom">http://www.siemens.com/dicom</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>	