

IHE Integration Statement

Product Name: ACUSON P Family
Ultrasound System

Version:
ACUSON P500™ – VA10

Date: October 2015

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	ACUSON P Family Ultrasound System		Oct 2015
	ACUSON P500	VA10	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Information Reconciliation	Acquisition Modality	None	
Evidence Documents	Acquisition Modality	OB-GYN Ultrasound Procedure Report (TID 5000) Vascular Ultrasound Procedure Report (TID 5100)	
Internet address for vendor’s IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org	In Japan: http://www.ihe-j.org	