

IHE Integration Statement

Product Name: ACUSON NX3
Ultrasound System

Version:

ACUSON NX3 VA10

ACUSON NX3 Elite VA10

Date: August 2015

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare	ACUSON NX3	VA10	Aug 2015
	ACUSON NX3 Elite	VA10	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Patient Information Reconciliation	Acquisition Modality	None	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Evidence Documents	Acquisition Modality	OB-GYN Ultrasound Procedure Report (TID 5000)	
		Vascular Ultrasound Procedure Report (TID 5100)	
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query	
		Patient-based Worklist Query	
Cardiology Evidence Documents	Acquisition Modality	Echo Evidence	
		Adult Echocardiography Procedure Report (TID 5200)	
Internet address for vendor’s IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org	In Japan: http://www.ihe-j.org	