

SIEMENS



IHE Integration Statement

SOMATOM Definition, AS, Edge, Flash

Computed Tomography

IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens AG Healthcare	SOMATOM Definition SOMATOM Definition AS SOMATOM Definition Edge SOMATOM Definition Flash	VA44	07/2012
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Patient based worklist query	
		Broad worklist query	
		Modality group case → Siemens solution: Study split resp. Multi Study	
	Image Display	No options defined in IHE	
	Evidence Creator	none	
Patient Information Reconciliation	Acquisition Modality	No options defined in IHE	
Access to Radiology Information	Image Display	none	
Consistent Presentation of Images	Print Composer	No options defined in IHE	
Radiation Exposure Monitoring	Acquisition Modality	none	
	Dose Info Consumer (available as a standalone tool by CARE Analytics)	Except transaction “Store dose information”	
Internet address for vendor's IHE information: http://www.siemens.com/IHE			
Links to Standards Conformance Statements for the Implementation			
HL7	www.siemens.com/hl7		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	

Global Siemens Headquarters
Siemens AG
Wittelsbacherplatz 2
80333 München
Germany

Global Siemens Healthcare Headquarters
Siemens AG
Healthcare Sector
Henkestraße 127
91052 Erlangen
Germany
Phone +49 9131 84-0
www.siemens.com/healthcare

Legal Manufacturer
Siemens AG
Wittelsbacherplatz 2
DE-80333 München
Germany

Global Business Unit
Siemens AG
Medical Solutions
Computed Tomography & Radiation Oncology
Siemensstraße 1
DE-91301 Forchheim
Germany
Phone +49 9191 18-0
www.siemens.com/computedtomography

Order No. CT-IHE.640.20.02.02
© 2014-2015, Siemens AG

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

Some/All of the features and products described herein may not be available in the United States or other countries. The information in this document contains general technical descriptions of specifications and options as well as standard and optional features that do not always have to be present in individual cases. Siemens reserves the right to modify the design, packaging, specifications and options described herein without prior notice.

Please contact your local Siemens sales representative for the most current information. In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we recycle certain components. Using the same extensive quality assurance measures as for factory-new components, we guarantee the quality of these recycled components. Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced. Caution: Federal law restricts this device to sale by or on the order of a physician.