

A Siemens Multitom Rax robotic X-ray system is shown in a clinical setting. A patient is lying on a table, and the robotic arm is positioned over them. The background features a large, colorful, abstract graphic of a spiral. The Siemens logo is visible in the top left corner.

SIEMENS

IHE Integration Statement

Multitom Rax

Join the new movement with robotic X-ray

[siemens.com/ihe](https://www.siemens.com/ihe)

IHE Integration Statement

Vendor	Product Name	Version	Date
Siemens Healthcare GmbH	Multitom Rax	VE20	Sept. 2015
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Radiology Scheduled Workflow	Acquisition Modality	Patient Based Worklist Query	
		Broad Worklist Query	
		PPS Exception Management	
		Billing and Material Management	
		Assisted Acquisition Protocol Setting	
Patient Information Reconciliation	Acquisition Modality	No options defined	
Consistent Presentation of Images	Print Composer	No options defined	
Internet address for vendor's IHE information: http://www.siemens.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.siemens.com/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org	



Multitom Rax is not commercially available in all countries. Due to regulatory reasons, its future availability cannot be guaranteed. Please contact your local Siemens organization for further information. On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products/services/features included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice. The information in this document contains general descriptions of the technical options available and may not always apply in individual cases. The product names and/or brands referred to are the property of their respective trademark holders.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens sales representative for the most current information. In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we may recycle certain components where legally permissible. For recycled components we use the same extensive quality assurance measures as for factory-new components. Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Siemens Healthcare Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone +49 9131 84-0
siemens.com/healthcare