

IHE Integration Statement

Product Name:
syngo MultiModality
Workplace

Version: VE40F
Date: July 2016

IHE Integration Statement			
Vendor	Product Name	Version	Date
Siemens Healthcare	syngo MultiModality Workplace	VE40F	July 2016
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Display	No options defined	
	Evidence Creator	None	
Consistent Presentation of Images	Print Composer	No options defined	
Simple Image and Numeric Report	Report Creator	None	
	Report Reader	None	
Access to Radiology Information	Image Display	None	
Portable Data for Imaging	Portable Media Creator	None	
	Image Display	No options defined	
	Print Composer	No options defined	
	Report Reader	No options defined	
RT Objects	Contourer	No options defined	
	Geometric Planner	No options defined	
Consistent Time	Time Client	None	
Internet address for vendor's IHE information: http://www.healthcare.siemens.com/services/it-standards/ihe-integrating-the-healthcare-enterprise			
Links to Standards Conformance Statements for the Implementation			
HL7	n.a.		
DICOM	http://www.healthcare.siemens.com/services/it-standards/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net		In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org