

## Ultrasound – Designed for You.

ACUSON S Family™ Ultrasound Systems, HELX™ Evolution with Touch Control



## Workflow

12.1"

**touch display** with Instant Response Technology

- > Easy-to-navigate tab menus
- > Smart-adaptive annotations
- > Dynamic protocol options

33%

fewer clicks coupled with dynamic backlighting for intuitive guidance

44%

fewer tactile keys in home base simplifying


the user interface

fewer software selections through heads-up display for no-look operation





## Imaging Performance



new large-format screen with 60% bigger image area: 21.5-inch, LED-backlit, 1080p LCD monitor more pixel density for greater image resolution provided by SieStream™ HD architecture



2 x

contrast resolution improvements for 6C1 HD & 18L6 HD transducers for improved penetration and skin-line detection

80%

new hardware now supports firmware updates and diagnostics

750/ accuracy with patter

accuracy for eSieFusion™ imaging with pattern recognition software merging CT to ultrasound in real-time



Remote service, auto-calibration and streamlined hardware provide multiple improvements:



Fewer service visits



All systems have Ultrasound System Security for virus protection



> 70% first-time fix rate



> 50% remote software updates



Windows® 7 operating system increasing network security protection



## Sustainability



upgradeability and forward compatibility between systems



recyclable, meeting RoHS standards



integrated **gel** warmer with customizable positions



supports high-speed wireless connectivity for on-the-go DICOM reporting