

SOMATOM

Perspective

Most Economical CT

Healthcare service providers are increasingly facing rising costs and shrinking budgets. At the same time, demand for flawless healthcare has intensified. These two factors make the efficient use of medical devices crucial to success in clinical practice. Siemens Healthineers has the ideal answer to these challenges: SOMATOM® Perspective, a scalable CT scanner applicable for all clinical fields that delivers economic benefits without compromising on patient care.

Improve image quality

SOMATOM Perspective is equipped with Interleaved Volume Reconstruction (IVR), a technology that enhances spatial resolution along the patient's longitudinal axis by increasing the number of data points. For this, IVR generates up to 128 slices during the image reconstruction. This makes it possible to visualize very small structures, such as fine fractures.

Expediting the daily routine

All SOMATOM Perspective configurations feature Fully Assisting Scanner Technologies (FAST) to automate time-consuming processes. FAST Planning enables one-click positioning of the scan and recon range, for a faster, more reproducible and standardized workflow. This reduces the risk of rescanning due to incorrect positioning.

Efficient operation with eCockpit

eCockpit is a bundle of features available for the SOMATOM Perspective that support highly economic scanner usage and reduce overhead costs. eStart gradually warms up the tube after extended periods of non-use (such as overnight), reducing wear and tear associated with cold starts of the tube. eMode automatically helps to enable the optimum balance between dose, efficiency, and image quality—and further extends the scanner's lifespan. eSleep minimizes electricity consumption while the scanner is not in use, for example overnight or during breaks, avoiding increases in operating costs.



Limited-
Time Offer

Save \$20,000
plus your Trade-in value on a
SOMATOM Perspective 16 and 32

OR

Save \$30,000
plus your Trade-in value on a
SOMATOM Perspective
64 and 128

For more information:
Visit usa.healthcare.siemens.com/computed-tomography or call 1-888-826-9702 and ask to speak to your local sales representative.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and is subject to change without prior notice. Some/ALL of the features and products described herein may not be available in the United States.

The information in this document contains general technical descriptions of specifications and options as well as standard and optional features, which do not always have to be present in individual cases.

Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice. For the most current information, please contact your local sales representative from Siemens.

Note: Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthineers

Local Contact Information

Siemens Medical Solutions USA, Inc.
40 Liberty Boulevard
Malvern, PA 19355-9998
USA
Phone: +1-888-826-9702
usa.siemens.com/healthineers